# ClaroCit


#50300113

22.10.2015

Instructions

Anleitung

Mode d'emploi

 **Struers**

# English

## Description

ClaroCit is a clear two component acrylic mounting system, consisting of one powder and one liquid. ClaroCit is a universal mounting material and ideal for filling vias  $\geq 200 \mu\text{m}$ . The system cures at room temperature; to obtain the clearest mount the use of a pressure chamber is recommended.

## Material

### Safety Data

*ClaroCit Powder*

Please read the following carefully:

Contains: Dibenzoyl peroxide.  
May produce an allergic reaction.

*ClaroCit Liquid*

Danger

Contains: Methyl methacrylate.  
H225 Highly flammable liquid and vapour.  
H315 Causes skin irritation.  
H317 May cause an allergic skin reaction.  
H335 May cause respiratory irritation.  
P210 Keep away from heat/sparks/open flames/hot surfaces. - No smoking.  
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.  
P280h Wear protective clothing and gloves.  
P302/352 IF ON SKIN: Wash with plenty of soap and water.  
P333/313 If skin irritation or rash occurs: Get medical advice/attention.

## Mounting of Specimens

### *Measuring of Components*

The ratio for ClaroCit Powder:ClaroCit Liquid is **20:12** by weight and **5:2** by volume. Weighing is recommended to obtain the best mounts and for reproducibility.

#### ***Important!***

Do NOT pour excess liquid back into the bottle as this will significantly reduce the shelf-life.

### *Mixing of Components*

- Measure out the amount of ClaroCit Powder into a mixing cup (max. 60 g / 5 x 20 ml).
- Measure out the amount of ClaroCit Liquid into a mixing cup (max. 36 g / 4 x 10 ml).
- Pour the powder into the liquid while stirring.
- Stir **thoroughly for 1½ min.**, without introducing air bubbles into the mixture.  
Immediately after the components are mixed polymerisation of the resin starts and the viscosity will increase. To get the best filling use the resin as soon as possible after mixing.

### *Selecting the Mounting Cup*

All Struers mounting cups can be used.

### *Pre-treatment of specimens*

Ensure that the specimens are thoroughly cleaned and dry prior to mounting. Each surface must be free from grease, moisture and other contaminants to promote the best possible adhesion of ClaroCit to the specimen.

### *Mounting*

- Place the cleaned specimen in the centre of the mounting cup, mix the necessary amount of resin and carefully pour over the specimen.
- Ensure that no air remains trapped within cracks or cavities in the specimen by moving the specimen or by tapping the mounting cup on the table.

**Curing**

ClaroCit is a room temperature curing system. For a 30 mm dia. mount, the curing time is about 20 minutes. It is recommended to use a pressure chamber (approx. 2 bar) in order to obtain a clear mount with no air trapped within the resin. When using a pressure chamber the curing time is unchanged but the temperature in the chamber will normally rise during curing, especially if many mounts are cured at the same time. We recommend keeping the mount in the mounting cup until the temperature is below 40°C/104°F.

*After Curing*

After removal from the mounting cups, the mounts are immediately ready for preparation.

# Technical Data

	Components	Ratio by weight	Ratio by volume
ClaroCit	ClaroCit Powder	20	5
	ClaroCit Liquid	12	2

**NB** For the best mount use weighing method.

## Density

ClaroCit Powder: 0.66 g/ml

ClaroCit Liquid: 0.95 g/ml

## Curing

Can be prepared (ground and polished) after:	20 min.
Peak temperature:	90°C
Pot life:	1½ min.

30 mm dia. mounting cup without specimen at 21°C / 70°F.

## Properties after Curing

Hardness: 85 Shore D

Density: 1.2 g/ml

## Tips

- ClaroCit can be coloured with EpoDye a fluorescent dye for easily identifying of pores and cracks and AcryDye for colour.
- When mixing very small amounts (< 20 g / 30 ml powder and 12 g / 12 ml liquid) it may be necessary to add a few drops more liquid to ensure a miscible mixture.

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