# **Diamond and CBN Cut-Off Wheels**

Instructions for Use

# **Diamant-und CBN Trennscheiben**

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# Introduction

Struers Diamond and CBN Cut-off Wheels are long-term consumable wheels used for cutting very hard materials. The abrasive grains in the outer rim of the wheels are diamonds or CBN (cubic boron nitride) held in a binder of resin or metal.

# Instructions

## Safety Precautions

Do not exceed the max. rotational speed marked on the cut-off wheel. The workpiece must be securely fixed when cutting.

## WARNING

**Diamond cut-off wheels should not be used on ferrous metals or materials containing Chromium or Vanadium. The diamonds will react chemically with these materials under the high temperature conditions existing in the cutting zone.**

## Mounting of Cut-off Wheels

Make sure that the sides of the cut-off wheel and the flanges are clean. Mount the cut-off wheel between the flanges of the cut-off machine. Do not use cardboard washers. The rotational direction of the cut-off wheel itself is indifferent unless otherwise stated on the cut-off wheel.

## Dressing

New cut-off wheels should always be dressed before use. When cutting, the cutting ability is gradually reduced and the cut-off wheel rim will appear shiny through sample material deposits on the wheel surface and wear on the abrasive. Cutting ability can be restored by dressing. Use the aluminium oxide dressing stick supplied with the cut-off wheel. The dressing stick is also available as a consumable (Cat. No. 40800044).

- Mount the dressing stick like a workpiece. Do not dress manually.
- Cut through the dressing stick using a moderate cutting pressure and plenty of cooling fluid.
- Repeat the treatment if the cut-off wheel does not cut satisfactorily.

**NB** Do not perform more dressing than necessary as this will cause needless wear on the wheel.

## WARNING

**Never use a cut-off wheel with inadequate dressing. Blunt abrasive grains will cause friction and generate heat, which can destroy the cut-off wheel.**

# Cutting

## Cooling

Cutting should always be carried out as wet cutting. Struers Additive for Cooling Fluid, Corrozip, should be added according to the instructions.

**NB** A pH value > 10 might damage the cut-off wheels.

## Oscillating Cutting

Do not use oscillating cutting (e.g. ExciCut). This might damage the cut-off wheel. On precision cut-off machines oscillation or rotation of the sample can be used.

## Feed Speed

Make sure the force between the cut-off wheel and workpiece is limited. Use the Force Limitation feature or Feed Reduction feature if it is available on the cut-off machine. Begin with a low feed speed.

## Spindle Speed

### *Diamond Cut-off Wheels*

Optimal wheel speed is 18 - 24 m/sec.

Wheel speed 18 m/s ~ spindle speed 350,000 / D rpm

Wheel speed 24 m/s ~ spindle speed 450,000 / D rpm

D = the diameter of the cut-off wheel in mm.

### *CBN Cut-off Wheels*

Optimal wheel speed is 25 - 35 m/sec.

Wheel speed 25 m/s ~ spindle speed 480,000 / D rpm

Wheel speed 35 m/s ~ spindle speed 660,000 / D rpm

D = the diameter of the cut-off wheel in mm.

## Optimising the Cutting Results

### *Better surface quality*

Use a low feed speed and no specimen holder rotation (if available).

### *Flatter samples*

Use primarily a low feed speed, maximum wheel speed and the largest flanges. The initial cut is critical. Use a low initial feed speed.

### *Better parallelism*

Use a low feed speed.

### *Faster cutting*

Position the sample so that the wheel will cut the smallest possible cross-section and then increase the feed speed.

### *Composites with metal content*

Composites with metal content should be cut with abrasive cut-off wheels. However, if this causes too much damage to the sample, then a diamond cut-off wheel with resin bonding can be used.

**NB** This could result in excessive wheel wear and damage on the cut-off wheel.

**Return Position**

If a sample is clamped on both sides of the cut-off wheel, the cut-off wheel should not be returned automatically when the cut is completed. This could damage the rim of the cut-off wheel, should the cut be slightly askew. If possible the cut-off wheel should remain in the stop position. Remove the workpiece before returning the cut-off wheel.

**Storage**

Clean and dry the cut-off wheel before storing in its box. If necessary, use ordinary detergents for the cleaning.

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