

# Abrapol-2

## Instruction Manual



Manual No.: 14227001

Date of Release: 18.02.1999



*Abrapol-2*  
*Instruction Manual*

<b>Table of Contents</b>	<b>Page</b>
<b>User's Guide .....</b>	<b>1</b>
<b>Reference Guide .....</b>	<b>15</b>

---

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the *Serial No.* and *Voltage* on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

**Instruction Manuals:** Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

**Service Manuals:** Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2002.

**Struers A/S**  
Valhøjs Allé 176  
DK-2610 Rødovre/Copenhagen  
Denmark  
Telephone +45 36 70 35 00  
Fax +45 38 27 27 01

---



## **Abrapol-2**

### **Safety Precaution Sheet**

#### **To be read carefully before use**

1. The operator should be fully aware of the use of the machine according to the instruction manual.
2. The machine must be placed in an adequate working position.
3. Be sure that the actual voltage corresponds to the voltage stated on the side of the machine. The machine must be earthed.
4. Be sure that the water connections are without leaks. The main water supply should be turned off if you leave the machine unattended.
5. Make sure that the specimens in the specimen holder are securely fixed.
6. Keep clear of the preparation disc when operating.
7. If you observe malfunctions or hear unusual noises - stop the machine and call technical service.

---

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

---

# User's Guide

Table of Contents	Page
<b>1. Installation</b>	
Checking the Contents of Packing.....	3
Abrapol-2 .....	3
Checking the Contents of Packing of Accessories .....	3
Recirculation Cooling Unit.....	3
Dosing Unit (DOTWO).....	3
Solenoid Valve .....	3
Drip Lubricator .....	3
Microstop .....	4
Placing Abrapol-2 .....	4
Getting Acquainted with Abrapol-2.....	5
Front of Abrapol-2.....	5
Supplying Power .....	6
Direction of Rotation.....	6
Supplying Compressed Air .....	6
Mounting a Recirculation Cooling Unit (Accessory) .....	6
Mounting a Solenoid Valve (Accessory).....	7
Mounting a Dosing Unit (Accessory) .....	7
Mounting a Drip Lubricator (Accessory) .....	7
Mounting Microstop (Accessory).....	7
<b>2. Operation</b>	
Using the Controls .....	8
Buttons/Adjusting Knobs.....	8
Front Panel Keys.....	9
Program Cards .....	11
Programming a Preparation Step.....	11
When Grinding with Water .....	11
When Grinding/polishing with Lubricant .....	11
Recommended Dosing Levels.....	12
Editing a Preparation Step .....	12
Transferring a Preparation Method .....	12
Editing a Preparation Method.....	12
Mounting a Preparation Disc.....	12
Inserting the Specimen Holder.....	13
Adjusting the Specimen Holder Position.....	13
Adjusting the Force.....	13
Starting the Preparation Process .....	13
When grinding with water .....	13
When Grinding/polishing with Lubricant .....	13

*Abrapol-2*  
*Instruction Manual*

Stopping the Preparation Process .....	13
Setting Microstop .....	14
Manual Functions .....	14
When Grinding with Water .....	14
When Grinding(polishing with Lubricant .....	14

## **1. Installation**

**Checking the Contents of Packing**  
*Abrapol-2*

In the packing box you should find the following parts:

- 1 Abrapol-2 complete with inlet hose and cable
- 1 Splash ring
- 1 Lid
- 1 Non-return valve for recirculating Cooling Unit
- 1 Outlet hose
- 1 Fitting for outlet hose
- 2 Hose clamps for outlet hose
- 1 Hose for compressed air
- 1 Hose connection for compressed air
- 2 Hose clamps for compressed air
- 1 Hexagon key, 2.5 mm
- 20 Program Cards
- 1 Instruction Manual

**Checking the Contents of Packing of Accessories**  
*Recirculation Cooling Unit*

- 1 Recirculation Cooling Unit
- 1 Drain angle
- 1 Funnel
- 20 Plastic bags

*Dosing Unit (DOTWO)*

- 1 Bottle holder plate
- 1 Nylon Strap
- 4 Screws M5x10
- 1 Screw M4x8
- 4 Washers M5
- 1 Washer M4
- 2 Lubricants
- 3 Suspensions

*Solenoid Valve*

- 1 Solenoid valve with mounting plate
- 1 Clamp for inlet hose
- 2 Washers
- 2 Screws (M5x10)

*Drip Lubricator*

- 1 Lubricant bottle
- 1 Mounting plate with valves and bottle holder
- 2 Screws (M5x10)
- 2 Washers (M5)

*Abrapol-2*  
*Instruction Manual*

*Microstop*

- 1 Microswitch with mounting box
- 2 Screw (M5x10)
- 3 Screws (M5x16)
- 3 Washers (M5)
- 1 Cable clamp
- 1 Washer (M4)
- 1 Screw (M4x8)
- 1 Hexagon key (4 mm)
- 1 Micrometer screw
- 1 Instruction Manual

**Placing Abrapol-2**


Abrapol-2 should be placed on a plane and horizontal floor. The machine must be placed close to the power supply, water mains and water outlet facilities.

*Abrapol-2*  
*Instruction Manual*

**Getting Acquainted with  
 Abrapol-2**

*Front of Abrapol-2*

Take a moment to familiarise yourself with the location and names of the Abrapol-2 components.



- |   |   |
|---|---|
| ① Base unit                                     | ⑨ Front panel                             |
| ② Specimen holder                               | ⑪ Main Switch (left-hand side of cabinet) |
| Ⓜ Specimen holder motor                         | ⑫ Standby                                 |
| Ⓟ Quick-clamping device                         | ⑬ Double start buttons                    |
| ⊗ Quick coupling flange                         | ⑭ Emergency stop                          |
| ⊕ Handle for adjusting specimen holder position | ⑮ Force selector                          |
| ∅ Nozzle block                                  | ⑯ Display                                 |
| ⑧ Cooling water                                 |   |



## Supplying Power

**IMPORTANT**

Check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.

- Mount a plug on the electric cable and connect it as follows:

Yellow/green: earth

Black/brown: phase

## *Direction of Rotation*

Check that the preparation disc rotates counter-clockwise. If not change the two phases.

## Supplying Compressed Air

- Connect the compressed air supply with the inlet on the rear side of the machine by means of the air hose and the hose connection delivered with the machine.
- Fasten the air hose with a hose clamp.

The pressure supply should be 6 bar and should be supplied either from a central compressor, portable compressor with compressed air reservoir or compressed-air bottle. A capacity of 20 l/min at atmospheric pressure is sufficient.

## Mounting a Recirculation Cooling Unit (Accessory)

- Place the recirculation cooling unit on either side of the machine where you find it the most convenient.
- Mount the non-return valve on the pump of the recirculation cooling unit.
- Mount the inlet hose from the machine on the non-return valve.
- Connect the electric cable from the pump to the machine.
- Mount the outlet hose on the underside of the tank.
- Mount the drain angle on the other end of the outlet hose by means of the hose fitting and lead it down into the hole of the recirculation cooling unit.
- Check that there is a steady fall on the whole course of the outlet hose.
- Place a plastic bag in the tank and fill with water and additive.

**Mounting a Solenoid Valve  
(Accessory)**

If no recirculation cooling unit is mounted on Abrapol-2 a solenoid valve must be mounted instead:

- Fasten the solenoid valve on the rear side of the machine by means of two screws with washers.
- Connect the electric cable to the machine.
- Cut the inlet hose at a distance of 100 mm from the machine and connect the part still fixed on the machine to the inlet of the valve and connect one end of the cut off part to the outlet of the valve and the other end to the water tap.

**Mounting a Dosing Unit  
(Accessory)**

- Mount the dosing box on the left side of Abrapol-2 by means of screws.
- Loosen the nozzle plate and screw the nozzle holder firmly on the front plate of the machine. Fasten the nozzle plate again.
- Connect the compressed air hose on the rear side of the machine.
- Connect the plug on the rear side of the machine and fasten the cable with the nylon strap, screw and washer.
- Place the bottle support under the bottle box.
- Place the bottles, fill them, and connect them with the tubing.

Lubricant 1: blue lubricant

Lubricant 2: red lubricant

Suspension 1: DP-Suspension, highest grain size

Suspension 2: DP-Suspension, medium grain size

Suspension 3: DP-Suspension, lowest grain size

**IMPORTANT**

Only use Struers suspension in the dosing unit.

**Mounting a Drip Lubricator  
(Accessory)**

- Fasten the mounting plate on the arm for the specimen holder motor by means of the two screws with washers.
- Connect the plug to the machine.
- Fill the bottle with lubricant and place it in the lubricator.

**Mounting Microstop  
(Accessory)**

- Fasten the microswitch to the left side of the machine by means of screws.
- Connect the plug to the rear side of the machine and fasten the cable with the cable clamp, the screw and washer.
- Clamp the micrometer screw in the holder and mount the holder to the arm for the specimen holder motor by means of screws and washers.

## 2. Operation

### Using the Controls

#### *Buttons/Adjusting Knobs*


See section Getting acquainted with Abrapol-2 for location of buttons.

---



















Name	Key	Function	Name	Key	Function
<b>Buttons/Adjusting Knobs (See Drawing of Abrapol-2)</b>					
⑪		Turns on/off the main power of the machine. The main switch is located on the left-hand side of the machine.	⑭		Interrupts all functions. The specimen holder remains in lower position until the emergency button is released by pulling the button.
MAIN SWITCH			EMERGENCY STOP		
⑫		On/Off switch for daily use.	⑮		For setting the required force. The pre-set force can be read on the manometer.
STAND-BY			FORCE SELECTOR		
⑬		Starts the preparation process. Press the double start button simultaneously until the specimen holder starts rotating.	⑯		For reading the pre-set time of grinding/polishing.
DOUBLE START			DISPLAY		

*Abrapol-2  
Instruction Manual*

*Front Panel Keys*



*Abrapol-2*  
*Instruction Manual*



Name	Key	Function	Name	Key	Function
<b>Data</b> 					
⑦ SPEED		Sets the speed of disc. Choose between 150 rpm and 300 rpm.	⑩ LUBRI-CANT		Activates the dosing of lubricant. Choose between bottle 1-2.
⑧ TIME UP/ DOWN		Increases (⌚▲) or decreases (⌚▼) the preparation time.	⑪ SUSPEN-SION		Activates the dosing of suspension. Choose between bottle 1-3.
⑨ WATER		Activates the water flow. Remember to adjust the flow on the water tap.	⑫ MICRO-STOP		Activates Microstop.
<b>CHOICE</b> 					
⑬ STEP KEYS		8 keys for programming the individual preparation steps. With these keys a complete preparation method can be transferred to the machine memory.	⑭ CLEAR		The present preparation sequence is deleted when AUTO is activated.
<b>MODE</b>					
⑮ CREATE		To program a preparation step CREATE mode should be activated.	⑯ AUTO		When transferring a preparation method or if you want to clear a sequence AUTO mode should be activated.
<b>MANUAL FUNCTIONS</b> 					
⑰ DISC		The disc starts rotating independently of the specimen holder.	⑳ LUBRI-CANT		Activates the dosing of the selected lubricant.
⑱ STOP		The machine is stopped at any time during the preparation process. By restart the machine will start from the point where the preparation was interrupted.	21 SUSPEN-SION		Activates the dosing of the selected suspension.
⑲ WATER		Activates the water supply. Remember to adjust the flow on the water tap.			

## **Program Cards**


The program cards supplied with the machine are used for filling in data about the preparation method of a given material. Data about the single preparation step should be stated on the back of the card and data about the material to be prepared as well as information about the preparation disc and force should be stated on the front of the card. The card should be placed on the front panel so that it covers the Choice functions.

## **Programming a Preparation Step**






8 different preparation steps can be programmed individually in regard of time, speed of disc, water, type and dosing of lubricant and suspension. These 8 steps can be combined in a number of different ways so that different materials can be prepared without having to repeat the programming.

- Set the machine to CREATE mode.
- Press the PREPARATION STEP KEY, 1 to 8.
- Set the time by pressing the key TIME UP  or TIME DOWN . The pre-set time can be read on the display.
- Set the speed of disc to either 150 rpm or 300 rpm by pressing the key SPEED respectively on the left side of the key to select 150 rpm and on the right side to select 300 rpm.  
*NB!*The speed of disc will automatically be pre-set to 150 rpm when an unprogrammed preparation step key is pressed.

### *When Grinding with Water*

- Select water supply by pressing the key WATER .

### *When Grinding/Polishing with Lubricant*

- Select the type of lubricant needed for the step by pressing one of the two LUBRICANT keys Lub.1  or Lub.2 .
- When one of the lubricant keys has been activated the display (left) shows level 0. Set the dosing level (0-10) by pressing the selected lubricant key again. Each time you press, the dosing level increases (see table below for Recommended Dosing Levels).
- Select the suspension needed for the step by pressing one of the three SUSPENSION keys Susp.1 , Susp.2 , or Susp.3 .
- When one of the suspension keys has been activated the display (right) shows level 0. Set the dosing level (0-10) by pressing the selected suspension key. Each time you press, the dosing level increases (see table below for Recommended Dosing Levels).

## Recommended Dosing Levels

The table gives examples of recommended dosing levels for 7 different polishing cloths with different resilience. Remember that these examples are Struers recommendations and therefore adjustments may be needed according to the requirements of the specimens to be prepared.

Name of Cloth	Recommended Dosing Levels	
	Lubricant	Suspensions
DP-Plan	3	5
DP-Pan	3	5
DP-Dur	6	5
DP-Dac	6	5
DP-Plus	8	5
DP-Mol	8	5
DP-Nap	8	5

## Editing a Preparation Step

- Set the machine to CREATE mode.
- Press the key for the preparation step which should be edited.
- Enter the required changes by means of the Data Keys.

## Transferring a Preparation Method

- Place the program card for the material to be prepared on the two pins on the front plate so that it covers the Choice Program. The keys can be activated through the card and a light diode shows when a key is selected.
- Press AUTO mode.
- Press CLEAR to cancel existing step sequence.
- Select a preparation method by pressing the relevant PREPARATION STEP KEYS in the correct sequence. The light diode for the first selected step turns green and the following steps turn red indicating that the step to be activated when the machine is started is the step with the green light. The programmed data for this step can now be read on the data keys.

## Editing a Preparation Method

- Set the machine to AUTO mode.
- Press CLEAR to delete any existing sequence.
- Select a new sequence with the PREPARATION STEP KEYS.

## Mounting a Preparation Disc

- Place the disc on the turntable and move it until the three pins engage with the holes in the turntable.

- Inserting the Specimen Holder**
- Place the specimen holder under the quick coupling.
  - Press and hold the black flange of the column down with your thumbs while guiding the pressure tap of the specimen holder into the coupling.
  - Release the black flange.
  - Turn the specimen holder until the three pins engage with the corresponding holes.
- Adjusting the Specimen Holder Position**
- The position of the specimen holder has to be adjusted correctly in relation to the preparation disc to get the best possible preparation results.
- Loosen the two handles placed on each side of the specimen holder motor.
  - The specimen holder motor can now be moved manually to the sides. The eccentricity can be read on the scale placed on the specimen holder motor.
  - When the correct position is found fasten the handles again.
- Adjusting the Force**
- Turn the adjusting knob (F) to set the required force.
- Starting the Preparation Process**
- Start the preparation step by pressing the two start buttons simultaneously until specimen holder starts rotating.
- When grinding with Water*
- Regulate the water flow on the tap.
- When Grinding/Polishing with Lubricant*
- Adjust the dosing valve of the drip lubricator.
  - The remaining time of preparation can be read on the display.
- Stopping the Preparation Process**
- When the time has elapsed the preparation disc will automatically stop and the specimen holder will return to its initial position. The light diode for the next preparation step in the sequence turns green. The machine is now ready for the next step.
- A preparation step can be stopped at any time during the process by pressing the key STOP (V) in Manual Functions.



## Setting Microstop

- Press the key MICROSTOP  $\frac{\pm}{\mp}$  in the Data field.
- Press both start buttons. The specimen holder is activated independently of the disc and moves downwards without being rotated.
- Turn the micrometer screw clockwise until STOP is read on the display.
- Press the key STOP  $\nabla$ . The specimen holder moves upwards to its initial position.
- Turn the micrometer screw counter-clockwise corresponding to the required removal.

### **IMPORTANT**

When Microstop is activated the machine does not automatically proceed to the next step. Press the key STOP  $\nabla$  and the next step will switch from red to green.

## Manual Functions

Manual operations are possible as the disc can run independently of the specimen holder. The manual functions are found on the lower right corner of the front panel indicated with  $\curvearrowright$ .

- Start the disc separately by pressing the key DISC  $\circ$ . The disc will operate according to the pre-set speed 150 rpm or 300 rpm as indicated with a light diode.

### *When Grinding with Water*

- Start the flow of water independently of a preparation step by pressing the key WATER  $\text{H}_2\text{O}$  in Manual Functions. If lubricant and suspension are selected in the Data Program water cannot be selected.

### *When Grinding/Polishing with Lubricant*

- To be able to give manual dosing of lubricant and suspension you should first choose the type of lubricant and suspension you want to use in the data field and also set the dosing level.
- You can now press the key LUBRICANT  $\text{LUB}$  and the key SUSPENSION  $\text{SUS}$  in Manual Functions. The dosing will be at maximum for as long as the key is pressed.
- Stop the disc and water by pressing the key STOP  $\nabla$  or by pressing the key DISC  $\circ$  and the key WATER  $\text{H}_2\text{O}$ .

# Reference Guide

Table of Contents	Page
<b>1. Accessories and Consumables</b> .....	16
<b>2. Struers Metalog Guide™</b> .....	17
<b>3. Maintenance</b>	
Daily .....	18
Suspension nozzles.....	18
Cleaning of nozzles alone .....	18
Weekly .....	18
Coolant.....	18
Extended Cleaning of Nozzles and Tubes .....	18
Yearly .....	19
V-belts.....	19
General Cleaning .....	19
Cooling Tank .....	19
Painted Surfaces .....	19
Air Filter.....	19
Lubricant Nozzles .....	19
Manual Lubrication .....	19
Lubrication Chart.....	20
<b>4. Machine Adjustments</b>	
Adjustment of counter-pressure .....	21
Adjustment of down-speed.....	21
Adjustment of brake action in bottom position of vertical movement .....	21
Adjustment of quick-coupling .....	21
<b>5. Trouble shooting</b> .....	22
<b>6. Technical Data</b> .....	24

## 1. Accessories and Consumables

<b>Specification</b>	<b>Code</b>
<i>Cooling unit</i> Recirculating cooling unit, complete, 3-phased	TRECI
<i>Additive</i> Bottle with 1 litre additive for cooling fluid	ADDUN
<i>Dosing unit</i> Dosing unit for 2 lubricants and 3 suspensions, complete	DOTWO
<i>Water valve</i> Magnet valve for control of tap water supply, complete	ABVEN
<i>Connection Kit for Multidoser</i> Serial no. > 260 1383	MULAT
<i>Drip lubricator</i> Holder for one lubricant bottle, with valve for ON/OFF and mechanical regulation of quantity	ABPYR
<i>Microstop</i> For stopping the process after removal of the required material	ABKRO
<i>Grinding/polishing disc</i> Wet grinding disc, aluminium, for ø305 mm SiC-paper	ROTAL
<i>DP-polishing disc, aluminium</i> ø230 mm (9") ø250 mm (10") ø300 mm (12")	DETOT DEDIF DEDAL
<i>MD-polishing disc, aluminium</i> ø200 mm (8") ø250 mm (10") ø300 mm (12")	DEMLA DEMIF DEMAL
<i>Petrodisc-M</i> ø290 mm (8") for 200 mm specimen holder ø230 mm (9") for 160 mm specimen holder	PETIL PETHA
<i>Specimen holders</i> For many sizes and shapes of specimens. Ask for separate list.	
<i>Table unit for Abraplan or Abrapol-2</i> With cabinet and trolley for recirculating cooling unit and 1 set of tubes for cooling unit	LABAB
<i>Storage unit</i> For accessories and consumables. The unit consists of 2 chapters, one with 8 drawers, including a drawer with inserts for specimen holders, and another one with 8 shelves.	LABCO
<i>Insert for specimen holders</i> Fits into the specimen holder drawer of the LABCO storage unit	LABOA

***Remember...***

Struers offers a comprehensive range of consumables for grinding and polishing. Please ask for separate leaflets.

## 2. Struers Metalog Guide™



### **Metalog Guide**

In Struers Metalog Guide™ you will find a detailed description of grinding/polishing methods for automated mechanical specimen preparation.

Struers Metalog Guide™ offers preparation methods for the most common materials, based on a simple analysis of two key properties: hardness and ductility. Finding the right method is easy, including choice of consumables. Always consult Struers Metalog Guide™ for the correct preparation method for the actual specimens.

Struers Metalog Guide™ contains 6 useful chapters:

Metalogram: a quick and safe guide to the right preparation method.

Metalog Methods: a complete catalogue of preparation methods, based on Struers' vast store of materialographic experience, and employing Struers' range of consumables.

Preparation Philosophy: the basics of modern specimen preparation, seen from a professional point of view.

Metalog Process: the materialographic preparation process from start to finish, logically explained.

Metalog Master: a combined trouble-shooting guide and supply of in-depth information on the processes of mechanical preparation, including an expert system for the solving of preparation problems.

Consumables Specification: quick access to the relevant consumables for the chosen preparation methods.

A complete guide to materialographic specimen preparation. Contact your local dealer for a free copy of Metalog Guide™.

### **3. Maintenance**

#### **Daily**

##### *Suspension nozzles*

Shake the DP-suspensions lightly every morning. Blow the nozzles free from suspension every night.

##### *Cleaning of nozzles alone*

Should be carried out whenever the machine is not operated for more than 4-6 hours:

- Take a bottle of lukewarm water.
  - Replace the 1 $\mu$ m suspension bottle with the water bottle. If 1 $\mu$ m suspension is not used, choose the finest diamond suspension used. This should correspond to SUSP. 3.
1. Press Create in MODE.
  2. Remove the plastic suspension supply tube from suspension nozzle 3.
  3. Manually activate ABRASIVE for 5 seconds with a beaker under the nozzle.
  4. Refill tube.
  5. Repeat procedure for SUSP. 2 and SUSP. 1.

#### **Weekly**

##### *Coolant*

Check if the supply of coolant in the tank is sufficient and not too contaminated.

##### *Extended Cleaning of Nozzles and Tubes*

Must be performed every week or before breaks longer than 24 hours. The cleaning should be carried out in the following order:

1. Take a bottle of lukewarm water.
2. Replace the 1 $\mu$ m suspension bottle with the water bottle. If 1 $\mu$ m suspension is not used, choose the finest diamond suspension used. This should correspond to SUSP. 3.
3. Press Create in the MODE.
4. Press key for step no. 1 in CHOICE.
5. Press SUSP. 3 in DATA.
6. Press SUSP. in Manual Control and hold a beaker under the nozzle to avoid water splashing on the cloth. Press for 5-10 seconds.
7. Move the bottle with water from 1 $\mu$ m suspension to 3 $\mu$ m suspension. Press SUSP. 2 corresponding to 3 $\mu$ m suspension. Repeat point 6.
8. Repeat point 7 for 6 $\mu$ m, suspension, SUSP. 1.

**Yearly**  
*V-belts*

Check regularly and at least once a year that the V-belts are tight. Adjust if necessary.

**General Cleaning**

The machine and the recirculating cooling system should be kept as clean as possible in order to avoid contamination of the specimens.

*Cooling Tank*

Every time water is changed, the cooling tank should be cleaned carefully to remove all grinding waste. Cleaning will be easier if Struers plastic inserts are used. Remember to add Additive for Cooling Fluid.

**IMPORTANT**

Do not use benzine or petroleum for cleaning when the additive is used.

*Painted Surfaces*

Painted surfaces and keyboard should be cleaned with a moist cloth and common household detergents.

**IMPORTANT**

Never use alcohol, acetone or similar solvents.

*Air Filter*

The air filter should be cleaned regularly and water in the glass vessel around the filter let out. This is carried out by pressing the button at the bottom of the vessel. Access through hole in the rear plate (fig. 10.3).

*Lubricant Nozzles*

The lubricant nozzles can be blocked by impurities in the lubricant:

- Remove the nozzle and try to blow it clean with compressed air (blow from the point and inward).
- If this does not help, the nozzle has to be cleaned with a very fine needle under microscope.

*Manual Lubrication*

Before the nozzle is remounted, the whole system should be washed out by activating the Manual lubricant key.

*Abrapol-2*  
*Instruction Manual*

**Lubrication Chart**

<b>Point of Lubrication</b>	<b>Lubricant</b>	<b>Interval</b>
Quick coupling	Oil	200 h
Surfaces of columns of up-and-down unit	Oil	200 h
Grinding disc bearing	Grease	2,000 h
Gear of gear motor	Grease	10,000 h or 2 years

Grease: eg. Shell Albida Grease LX

## **4. Machine Adjustments**

The Adjustments listed below are only required after repairs or service.

- Adjustment of Counter-pressure** Counter-pressure is the pressure on the primary side of the compressed air cylinder, which counterbalances the weight of the up/down unit. The counter-pressure returns the unit to its initial top position. It is adjusted by a reduction valve (fig. 20b.8). The correct pressure is 2.7 bar (when using air to eliminate hysteresis). Measure with a precision manometer.
- Adjustment of Down-speed** The Down-speed is adjusted by an adjustable throttle valve (fig. 10.4), normal position of the valve being approx. 1 3/4 turn from home position. To obtain faster travel, turn the valve approx. 1/4 turn. Screwing the valve beyond the latter position may cause problems in the form of excessive initial grinding pressure (for approx. 1/4-1 sec), involving a risk of the specimen tearing the grinding paper or polishing cloth. When specimens are protruding more than 6 mm from the sample holder disc, the brake (see below) does not work, and in this case the throttle valve (fig. 10.4) is closed 1/2 to 1 turn to reduce the down-speed.
- Adjustment of brake action in bottom position of vertical movement** The brake incorporated in the compressed air cylinder comes into operation 5-10 mm above the polishing disc, depending on how far specimens protrude from the specimen holder. Access to the adjusting screw is gained through two holes at the left side plate. The upper screw adjusts the braking action of the up/down unit when it returns to its initial top position. Normal adjustment of the down-speed brake action is approx. 1/8 turn from the home position. Adjustment is very delicate, because the braking action should not prevent the specimen holder from contact with the grinding disc when working at a low process pressure. This may cause problems, particularly with thin specimens, therefore specimens should always protrude approx. 2.5 mm from the specimen holder.
- Adjustment of quick-coupling** After loosening the two screws (fig. 19.10) the position of the quick-coupling on the shaft can be vertically adjusted. Mounting the quick-coupling too close to the shaft end may transmit the grinding pressure to the specimen holder through the two lock-balls. This would damage the specimen holder stud and make removal of the specimen holder difficult. Incorrect mounting can be noticed either by measuring or by the fact that the lock slots of the specimen holder studs have become smooth from wear.



## 5. Trouble shooting

<b>Error</b>	<b>Cause</b>	<b>Action</b>
Display shows "ES"	The emergency stop button has been activated (Emergency Stop)	Release the emergency stop button by pulling the button.
	One of the motors has been overloaded.	Wait till "ES" disappears from the display. Reduce the grinding and polishing pressure. If "ES" immediately appears again, a phase is missing or the motor is defective.
No suspension escapes the nozzles	The nozzles are blocked.	Clean the nozzles.
	The suspensions do not have the right viscosity	You are using suspensions not manufactured by Struers: Change to the original Struers suspensions.  The suspensions are too old: Change the suspensions.  The suspension bottles have not been shaken for a long time: Shake the bottles.
Centre of specimen not polished (after previous grinding on Petrodisc-M or diamond grinding disc)	The specimens go too far beyond the outer edge of the grinding disc.	Adjust the position of the specimen holder disc in proportion to the grinding disc.
The polishing suddenly takes more time	The polishing cloth is filled up with removed material.	Change polishing cloth.
	The suspension ability has deteriorated and the diamonds sink to the bottom.	Shake the suspension bottle.
The supply of coolant is inadequate	The cooling pump rotates the wrong way.	Change two phases in the electric connection of the pump.
	There is no water in the tank.	Fill up with water and additive.

*Abrapol-2*  
*Instruction Manual*

<b>Error</b>	<b>Cause</b>	<b>Action</b>
The grinding paper is torn up by contact with the specimens	Excessive grinding pressure.	Reduce the pressure with the reduction valve (fig. 1.7).
	The specimens are not at the same level.	Use Struers levelling apparatus.
	Soft and brittle specimens.	Use adhesive grinding paper or run a dummy specimen over the whole surface of the grinding paper.
	Specimen holder lowered too fast (especially if the specimens project more than 6 mm over the edge).	The speed of the downward movement is controlled by a throttle valve (fig. 20.8). Adjust the speed by turning the screw at the top of the valve, clockwise for decreasing and counter-clockwise for increasing the speed.
	Incorrect adjustment of the shock reduction device in the compressed-air cylinder.	Adjust the screw on the cylinder. Access through the lower hole in the side plate (behind the bottles if doser unit is mounted). Normal setting 1/8 rotation from the starting point.
Display shows NO AIR	Missing or incomplete air supply.	Check that air hoses are tight and properly clamped.
Display defective	Defect in the circuit of the display or defect in the CPU circuit.	Replace the display.

## 6. Technical Data

### **Abrapol-2**

<i>Voltage/frequency</i>	3 x 200V, 50/60Hz, max. 6.8A 3 x 220V, 50/60Hz, max. 2 A 3 x 380V, 50Hz, max. 3.8 A 3 x 415V, 50Hz, max. 3.6 A 3 x 440V, 60Hz, max. 3.5 A 3 x 480V, 60Hz, max. 3.6 A
<i>Internal fuse</i>	1.6 AT
<i>Entry pressure</i>	6 bar $\pm$ 0.2 Maximum consumption: approx. 20 l/min.
<i>Main motor</i>	550W (0.75 hp) at disc rotation 150 rpm 1100W (1.5 hp) at disc rotation 300 rpm
<i>Specimen holder motor</i>	150 W (0.2 hp)
<i>Rotational speed</i>	Disc: 150 rpm or 300 rpm Specimen holder: 150 rpm
<i>Vertical working pressure</i>	50-700 N
<i>Timer intervals</i>	From 0:05 to 99:55 min at 0:05 min intervals
<i>Dimensions and Weight</i>	Height: 1420 mm Length: 850 mm Depth: 510 mm Weight: 250 kg
<i>Noise Level</i>	The noise level of the machine is 55 dB (A) measured at idle running at a distance of 1 m/39.4" from the machine.
<i>Safety standard</i>	IEC 204 / EN 60204-1 (VDE 0113)
<b>Recirculating Cooling Unit</b>	
<i>Motor</i>	110 W
<i>Pump</i>	64 l/min (1 m water column pressure at 0.1 atm)
<i>Dimensions and Weight</i>	Height: 410 mm Length: 400 mm Depth: 500 mm Weight: about 7.2 kg
<i>Standard</i>	IP 44

# Abrapol-2

## Spare Parts and Diagrams



Á

Manual No.: 14227001

Date of Release: 18.02.1999

Revised: 09.06.2000



## *Abrapol-2*

### *Spare Parts and Diagrams*

---

**Always state *Serial No* and *Voltage/frequency*  
if you have technical questions or when ordering spare parts.**

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

**Instruction Manuals:** Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

**Service Manuals:** Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2000.

**Struers A/S**  
Valhøjs Allé 176  
DK-2610 Rødovre/Copenhagen  
Denmark  
Telephone +45 36 70 35 00  
Telefax: +45 38 27 27 01

---

## 1. Spare parts and Diagrams

	<b>Table of contents</b>	<b>Drawing</b>
<b>Photos and drawings</b>	1 Front	
	2 Controls	
	3 Placing the program card	
	4 Inserting the specimen holder	
	5 Main switch	
	6 Non-return valve	
	7 Connection of recirculating pump	
	8 Valve for mains water	
	9 Dosing unit	
	10 Dosing unit, connection	
	11 Drip lubricator	
	12 Microstop	
	13 Contents of packing case	
	14 Recirculating Cooling Unit	
	15 Valve for mains water	
	16 Dosing unit	
	17 Drip lubricator	
	18 Microstop	
	19 Abrapol-2	
	20 Pneumatic parts	
	21 Pneumatic diagram	
	22 Contactor box	
	23 Circuit boards	
<b>Diagrams</b>	24 Circuit diagram for mains voltage.....	14221007G
	25 Wiring Diagram, Contactor Box .....	14221012H
	Control Voltage Wiring Diagram.....	14221034D
	26 Connecting Diagram Supply Voltage 100-230V....	14221013A
	27 Control Voltage Circuit Diagram.....	14221033D
	28 Circuit Diagram, Keyboard/Display Board .....	14221052
	29 Circuit Diagram, In/out, Supply Board .....	14221062B
	30 Pneumatic Doser Diagram.....	14220536A
	31 Pneumatic Mounting Diagram .....	14220231K
	32 Thermal relay settings.....	14221045C


The drawings are not to scale. Some of the drawings may contain position numbers not used in connection with this manual.

*Abrapol-2*  
*Spare Parts and Diagrams*

**1. Front**  
*Fig. 1*

**Pos. and Description**

**Spare Part No.**




*Abrapol-2*  
*Spare Parts and Diagrams*

**2. Controls**  
*Fig. 2*

**Pos. and Description**

**Spare Part No.**






*Abrapol-2*  
*Spare Parts and Diagrams*


**3. Placing the Program Card**  
*Fig. 3*

Pos. and Description	Spare Part No.
1 Switch.....	394MP040
2 Front plate with touch pad keyboard.....	14220087
Program set V02, 2 pcs. ....	422MP040
3 Pin for program card.....	422MP041
4 Ball cock .....	383MP020
5 Flex tube, complete .....	422MP003



**4. Inserting the Specimen Holder**


*Fig. 4*



*Abrapol-2*  
*Spare Parts and Diagrams*

**5. Main Switch**

*Fig. 5*




*Abrapol-2*  
*Spare Parts and Diagrams*

**6. Non-return valve**  
*Fig. 6*


**Pos. and Description**

**Spare Part No.**

- 1 Non-return valve, 3/4" for TRECI ..... 260MP082




**7. Connection of  
Recirculating Pump**



*Abrapol-2*  
*Spare Parts and Diagrams*


**8. Valve for mains water**

*Fig. 8*




*Abrapol-2*  
*Spare Parts and Diagrams*

**9. Dosing Unit**  
*Fig. 9*



*Abrapol-2*  
*Spare Parts and Diagrams*


**10. Dosing Unit, Connection**  
*Fig. 10*






*Abrapol-2*  
*Spare Parts and Diagrams*

**11. Drip lubricator**  
*Fig. 11*



**12. Microstop**  
*Fig. 12*




*Abrapol-2*  
*Spare Parts and Diagrams*

**13. Contents of Packing Case**  
*Fig. 13*

**Pos. and Description**

**Spare Part No.**

3 Water hose 1/2 " (10m).....RNU29316



② DP-disc or MD- disc to be ordered separately

**14. Recirculating Cooling Unit**

*Fig. 14*




Fig. 14

*Abrapol-2*  
*Spare Parts and Diagrams*

**15. Valve for Mains Water**  
*Fig. 15*

Pos. and Description	Spare Part No.
1 ABVEN Solenoid valve, complete 50Hz.....	04226902
1 ABVEN Solenoid valve, complete 60Hz.....	04226903
5a Solenoid 50Hz .....	375MP192
5a Solenoid 60Hz .....	375MP193
5b Solenoid valve for water, 50Hz.....	422MP004A
5b Solenoid valve for water, 60Hz.....	422MP004B
6 3-pole box, contact plug incl. ....	422MP005





Fig. 15

*Abrapol-2*  
*Spare Parts and Diagrams*

**16. Dosing Unit**  
*Fig. 16*

Pos. and Description	Spare Part No.
1 Bottle, 1L, cap with hole.....	376MP084
Cap plate for lubricant bottle (1000 ml).....	375MP028
2 9-pole box incl. contact plug .....	422MP008
3 Water trap.....	435MP017
4 Lubricant pump, complete.....	375MP025
5 Box without valves and pumps.....	422MP006
Solenoid valve .....	375MP187
O-ring $\varnothing 10.82-1.78$ (10 pcs) .....	130MP057
Air tube $\varnothing 5/\varnothing 3.2$ (2m), Blue.....	375MP178
6 Suspension bottle, 250 ml.....	375MP027
Non-return valve for cap plate .....	422MP018
Cap plate for suspension bottle (250 ml).....	422MP017
13 DP-Suspension nozzle, complete .....	375MP031
14 Lubricant nozzle, complete .....	375MP032
15 Nozzle holder plate .....	422MP013
Nozzle elbow.....	422MP014
Lubricant pump, complete.....	375MP025
Throttle gauge for air tubing, suspension nozzle.....	422MP016
16 Nozzle holder cover .....	422MP009



*Abrapol-2*  
*Spare Parts and Diagrams*

**16 A/B Dosing Unit**  
*Fig. 16 A/B*

<b>Pos. and Description</b>	<b>Spare Part No.</b>
<b>A</b>	
1 Lubricant nozzle head.....	375MP035
2 O-ring $\varnothing 3 \times 1$ for lubricant nozzle head .....	375MP036
3 Lubricant nozzle screw .....	375MP038
<b>A/B</b>	
4 O-ring $\varnothing 4,48 \times 1,78$ for nozzle screw.....	375MP037
5 Guide ball for nozzle .....	375MP270
6 O-ring $\varnothing 16 \times 2$ for nozzle ball.....	422MP015
<b>B</b>	
1 DP-Suspension nozzle head.....	375MP034
3 DP-Suspension nozzle screw .....	375MP039




Fig. 16a  
 Lubricant Nozzle





Fig. 16b  
 DP-Suspension Nozzle

*Abrapol-2*  
*Spare Parts and Diagrams*

**16C Dosing Unit**  
*16C (16.4)*

<b>Pos. and Description</b>	<b>Spare Part No.</b>
Lubricant pump, complete.....	375MP025
1 Diaphragm .....	375MP052
2 Piston .....	375MP053
3 Spring.....	375MP020
4 Nipple.....	375MP022
5 Banjo connection .....	375MP023
6 Check valve (non-return valve) .....	375MP024
Solenoid valve .....	375MP187
Filter.....	422MP007
Throttle gauge for air tubing, suspension nozzle.....	422MP016
Cap plate for suspension bottle .....	422MP017






*Abrapol-2*  
*Spare Parts and Diagrams*

**17. Drip Lubricator**  
*Fig. 17*


<b>Pos. and Description</b>	<b>Spare Part No.</b>
1 Bottle, 2 l.....	158MP016
2 Valve needle .....	158MP010
3 Valve outlet .....	158MP005
4 Solenoid valve .....	376MP044
5 Nozzle with fittings.....	422MP019
8 Lubricator cup.....	158MP014
9 Bottle fixture .....	158MP015
O-ring, ø11.91x2.62.....	249MP032
10 3-pole box incl. contact plug .....	422MP005
11 Solenoid Housing .....	14229048



*Abrapol-2*  
*Spare Parts and Diagrams*

**18. Microstop**  
*Fig. 18*

<b>Pos. and Description</b>	<b>Spare Part No.</b>
4 3-pole box incl. contact plug .....	422MP005
6 Microswitch.....	260MP202
13 Micrometer screw.....	260MP206




*Abrapol-2*  
*Spare Parts and Diagrams*

**19. Abrapol-2**  
*Fig. 19*

<b>Pos. and Description</b>	<b>Spare Part No.</b>
1 V-belt pulley (state frequency) .....	260MP068
Bushing for V-Belt pulley (260MP68) .....	260MP069
2 Nipple and connecting piece .....	260MP168
3 Emergency stop switch element .....	381MP440
4 Drain tray.....	260MP185
5 Emergency stop push-button, without switch element .....	381MP439
7 Cap for 260MP004 .....	260MP005
8 Push button, start.....	260MP004
9 Splash guard .....	260MP143
Lid for 260MP143 .....	260MP144
10 Quick-release coupling, complete .....	260MP027
11 Damper coupling, complete .....	260MP080
Damper coupling with lower flanger.....	12600330
12 Handle.....	260MP155
14 Specimen holder motor (state voltage and frequency) .....	260MP085
15 Sliding piece .....	422MP021
16 Dial rail .....	422MP022
17 Pressure spring.....	422MP023
18 Microswitch.....	260MP009
19 Double-acting compressed-air cylinder .....	260MP102
20 V-belt pulley for main motor .....	260MP071
21 Bushing for 260MP071 .....	260MP072
22 Driving belt, 50 Hz, 2 pcs. ....	260MP073
22 Driving belt, 60 Hz, 2 pcs. ....	260MP083
23 Main motor (state voltage and frequency) .....	260MP070
24 Hose, complete .....	422MP003
25 Water Cock 1/4" .....	383MP020
26 Column for Abra arm.....	12600312
27 Flexible water nozzle, complete.....	42MP003

*Abrapol-2*  
*Spare Parts and Diagrams*


**Fig. 19, Drawing**



*Abrapol-2*  
*Spare Parts and Diagrams*

**19A. Abrapol-2**  
*Fig. 19A*

Pos. and Description	Spare Part No.
Bearing unit complete, but without turntable, V-belt pulley and bearing flange .....	260MP065
1 Turntable w. quad ring and 4xM8 screw holes (300mm) .....	14590010
2 Axial bearing $\phi 100/\phi 135 \times 25$ .....	260MP017
Shell ALBIDA LX Grease .....	381MP415
3 Bearing flange .....	260MP197
4 Roller bearing, upper, 6007 .....	260MP018
5 V-belt pulley (please state frequency) .....	260MP068
6 Roller bearing, lower, 6007RS .....	260MP018
7 Bush for V-belt pulley .....	260MP069
8 Seeger ring I 62 x 2 .....	260MP019
9 Seeger ring U 35 x 1.5 .....	260MP067
10 Bearing protection ring .....	260MP183
11 V-ring VA 130 .....	260MP184
12 Thrust collar, set of 4 .....	315MP059
14 Quad-ring .....	2IQ04450
15 Quad-ring .....	2IQ04441
16 Quad-ring .....	2IQ04326



*Abrapol-2*  
*Spare Parts and Diagrams*



**20. Pneumatic Parts**

*Fig. 20*

**Pos. and Description**

**Spare Part No.**

1	Reducing valve .....	368MP330
2	Throttle sound absorber .....	422MP024
3	Manometer .....	14229027
4	Solenoid valve .....	260MP012
6	Quick exhaust valve.....	14222900
	Filter/reduction valve unit - SMC .....	375MP199
	Filter glass for 375MP199 .....	381MP435
8	Reduction valve.....	2YR00003
9	Non-return valve.....	260MP127
10	Non-return throttle valve .....	260MP128




*Abrapol-2*  
*Spare Parts and Diagrams*

**21. Pneumatic Diagram**

*Fig. 21 (Drawing 14220236)*

**Pos. and Description**

**Spare Part No.**




*Abrapol-2*  
*Spare Parts and Diagrams*

**22. Contactor Box**

*Fig. 22*

Pos. and Description	Spare Part No.
1 Thermo relay F1 (200-240V).....	381MP234
Thermo relay F1 (380-480V).....	394MP050B
2 Thermo relay F2 (200-240V).....	394MP049C
Thermo relay F2 (380-480V).....	381MP234
3 Mains noise filter .....	394MP045
4/5 Thermo relay F3/F4 (200-240V) .....	381MP233
Thermo relay F3/F4 (380-480V) .....	375MP263B
6 Transformer T1 (state voltage).....	394MP043
7 Main switch.....	2SE00149
8 Fuse F12, 2.5 AT (10 pcs.) .....	RFU14000
9 Fuse F11, 1.6 AT (10 pcs.) .....	422MP027
10 Auxiliary switch for K4.....	381MP237
11 Contactor (K1, K2, K3, K4) .....	376MP110







*Abrapol-2*  
*Spare Parts and Diagrams*

**23. Circuit Boards**

*Fig. 23*

<b>Pos. and Description</b>	<b>Spare Part No.</b>
1 Display Board .....	422MP010
2 CPU Board .....	422MP011
PROMP, latest version .....	422MP040
3 Interface/Power Supply Board.....	422MP012
4 3-pole cabled plug for solenoid valve.....	422MP029
5 9-pole cabled plug for solenoid valve.....	422MP030
Fuse on PCB F5 (2.5 AT) .....	422MP031
Fuse on PCB F6 (1.6 AT) .....	422MP032





Supply conductor and machine over-current protection.  
Sicherungen für Hauptkabel und Gerät.  
F = 20A  
Not supplied by Struers.  
Nicht von Struers geliefert.

F1-F4: See Fuse setting Table  
Sehe Sicherungstabelle

All wirings 1mm²

bk = black = schwarz  
bn = brown = braun  
ye/gn = yellow/green = gelb/grün

Connections for Struers recirc. pump  
3 x U/0.13 kW

Sample mover motor.  
Probenhalter Motor.

Polishing disc motor.  
Polierscheibe Motor.

G: lodning x2 S1 rd > bk 15.4.99  
F: 15 -> 75mm 10/11mm cv2/pt  
E: F=10A -> F=20A S.12.93 8H2/pt


D: 0.15 kW -> 0.18 kW 1.1/0.55 kW -> 0.15-0.18 kW 0.13/0.13	
C: 17AWG x4 + Power lif. bn -> rd; bk -> bn v. 4x1.5 1.5 7.96 83/8	ikke ang. rüheder: <input checked="" type="checkbox"/>
B: Spade term. -> gevinslag 2x 4 x 1.5 1.5 7.96 83/8	ikke ang. tolerancer: <input type="checkbox"/>
A: See Fuse ... lif. Control Voltage lif 1.1/0.55 83/8	Målestok: Overfladebehandl.
Matr: Abrapol - 2	Navn: BR
<b>STRUERS A/S</b>	Dato: 24.2.89
Kode: Anv. I type: Kfædsskema: Circuit Diagram for Mains Voltage. Schaltplan für Netzspannung.	Erst: 1006
	14221007G

Noise Filter.


Connections for Struers Rec. Pump.

See connecting diagram  
Sehe Verbindungs-schaltplan

Control voltage trafo  
Steuerspannungstrafo



H: F1-F3-F3-4 2- led. 00. 180400 BMA/JTV			
F: Lus 220V, v.T1.rd7, v.X2 0V/23/9-96BMA/JTV		G: Termovilla opdateret, tegnet på CAD. 19/11-98 CVE/JTV	
E: F12 : 1,6->2,5AT 8/7-93 BMA/JTV		D: Spødeterm -> gevind 24/8-91 BRY/VH	
C: Trafo delt i 2 1/9-90 HSB		B: F8->F12, F3->F11 30/4-90 LR/VH	
A: Div. tilretninger 16/8-98 BMA/VH			
Motr.:	Overf.bah.: Nej	Målforhold:	ikke ang. tol. efter DS/ISO 2768-
	Projektionsmetode:	1:1	Dato
			Tegn: 110498 BRY
Wiring Diagram, Contactor Box		Erat.:	
		14221012H	




All wiring red

*X10: ledning nr byttet om 3-5 og 4-2  
D: X10 til tekst og tem.nr. 28/1.97 BHJ/77J*


Matr.	Målestok: Overfladebehandl...
<b>STRUERS A/s</b>	
Kode: Anv. i type:	Navn: B.J. Kraft
Control Voltage Wiring Diagram	Dato 3.8.89
	Erst:
	14221034D

G.X10 pålænt 21/12.84 B7/6A


**Primary Supply**  
**Primär Versorgung**  
**100V AC**




**Primary Supply**  
**Primär Versorgung**  
**110/115V AC**





**Primary Supply**  
**Primär Versorgung**  
**200V AC**



**Primary Supply**  
**Primär Versorgung**  
**215/220V AC**



**Primary Supply**  
**Primär Versorgung**  
**230V AC**



Matr.

Målestok:



Ikke ang. tolerancer: ±

Ikke ang. ruheder:

Overfl.-behandl.

Ja   
 Nej

Dato Sign.

Tegn. 129.90 728  
 Kontr. 8/10-90 DM/PPY

Kode


Anv. i type: 394-419

Erst.:

**Connecting Diagram Supp. Volt. 100-230V**  
**Verbind.-schaltpl. Versorg. sp. 100-230V**

14221013A

Terminalnr. tilf.  
 150399 BMI/PPY



C1: 53: Not available  
 D: 09/04/1975  
 W: 14/04/1975  
 S: 14/04/1975  
 A: 14/04/1975

Kode	Art. type	Werk	Abt./opst.
		2	
<b>STRUERS A/S</b>			
Control Voltage Circuit Diagram			
14.221033D			

SIGN. FOR LUB UNIT  
 SIGN. FOR DOSER UNIT  
 All wirings 10

Art. No.	14.221033D
Line No.	1
Rev.	1
Date	4.7.75
Author	
Checked	
Material	
Quantity	
Unit	
Price	
Weight	
Volume	

W5 ye/gn From Contactor Box  
 rd1

Abrapol-2  
Spare Parts and Diagrams

28. Circuit Diagram,  
Keyboard/Display Board  
Diagram 14221052











FIG. II

Y0	Y1	Y2	Y3	Y4	Y5	Y6	Y7	
X0 P301a	P302a	P303a	P304a	HV301 SUSP1	HV309 WATER	HV310 GN STEP 1	HV314 GN STEP 5	
b	b	b	b	HV302 SUSP2	HV318 SINGLE STEP	HV310 RD STEP 1	HV314 RD STEP 5	
c	c	c	c	HV303 SUSP3	HV319 AUTO STEP	HV311 GN STEP 2	HV315 GN STEP 6	
d	d	d	d	HV304 RPM150	HV320 DISC ON	HV311 RD STEP 2	HV315 RD STEP 6	
e	e	e	e	HV305 RPM300	HV321 WATER ON	HV312 GN STEP 3	HV316 GN STEP 7	
f	f	f	f	HV306 LUB 1	HV306 EX 1	HV312 RD STEP 3	HV316 RD STEP 7	
g	g	g	g	HV307 LUB 2	HV307 EX 2	HV313 GN STEP 4	HV317 GN STEP 8	
dp	dp	dp	dp	HV308 µ -	HV308 µ -	HV313 RD STEP 4	HV317 RD STEP 4	
							LED'S	

FIG. I LED MATRIX






Rev.	Date	Author	Checked
1	20050322/1		

D: 230724 8:55:19  
 A: Full veld. R103 niet (2 x 100) 2/2  
 Matrik: Abropol - 2 2050322/1  
**STRUERS K/S**  
 Hout: 2005 / 2005  
 Chm: 13-10-05  
 E-mail:






Y: Solenoid Valves  
Magnet Ventile

P: Pumps  
Pumpen

bl: blue = Air  
blau = Luft

bk: black = Liquid  
schwarz = Flüssigkeit

		ikke ang. ruheder: <input checked="" type="checkbox"/>
A: Lgd + belegn. tilf. 20/9.96 BMJ/77N		ikke ang. tolerancer: ±
Matr. Abrapol-2, tilbehør: Dosersystem GB/D		Målestok: Overfladebehandl.:
<b>STRUERS A/s</b>		
Kode:	Anv. i type:	Navn: BR
Pneumatic Diagram Pneumatik Diagramm		Dato: 21.06.88
		Erst.:
		14220536 A



- A: PA  $\phi 5 \times 170$
- B: PA  $\phi 5 \times 220$
- C: PA  $\phi 8 \times 240$
- D: PA  $\phi 8 \times 310$
- E: PA  $\phi 8 \times 240$
- F: PA  $\phi 8 \times 125$
- G: PA  $\phi 8 \times 315$
- H: PA  $\phi 8 \times 300$



K. Pos 15 nykøbt fra Pos. 10 til pos. 35. 15/3.99 BHJ/PJ		Ikke ang. ruheder: <input checked="" type="checkbox"/>
J: Pos 220 fjernet pos 265 til F. 18/11.98 CVE/PJ		
I: Pos 130 + 135 v. manometer $\rightarrow$ Pos 125 + 126 29/10.98 BHJ/PJ		Ikke ang. tolerancer: $\pm$
H: Tekst $\rightarrow$ Pos. nr. Oriente pos 15 til F. Nippel 1/4-1/4 pos 325 til F. Slange led $\phi 8$ under. Pos skiver pos 76 til F. 25/9.98 BHJ/PJ		
Matr.		Målestok: Overfladebehandl.:
<b>STRUERS A/s</b>		
Navn: KK-BRH		Dato 18.4.88
Kode:	Anv. i type:	Erst: 14220231G
Pneumatik, monteringsstegn.		14220 231 K

Motor protection fuse setting:

Motor schutz Einstellung:

Rèlage du fusible protégeant le moteur:

U / f	F1/A	F2/A	F3/A	F4/A
200/50	0.4	0.96	4.8	5.1
200/60	0.4	0.96	4.0	4.8
220/50	0.4	0.96	4.5	4.7
220/60	0.4	0.96	4.5	4.7
380/50	0.23	0.55	2.6	2.7
415/50	0.23	0.55	2.5	2.6
440/50				
440/60	0.23	0.55	2.6	2.7
480/60	0.23	0.55	2.5	2.6
500/50	0.18	0.5	2	2.1

	Matr.:	Overfl.beh.:	Målforhold:	Ikke ang. tol. efter DS/ISO 2768-	
C: F3 og F4 opdateret for 200, 415 og 480 V 18/11-98 CVE/ B: F2/A alle værdier ændr. undt. 500/50 040196 BMJ/JTV A: 200V/50+60tilf. 090392 BRH/JTV		Projektionsmetode 			Sign.
				Tegn: 170889	BMJ
	<b>Indstilling af termosikring Abrapol-2</b>	Erst.:		14221045C	
		Kontr. 170889	v.H.		