# Abramin

## Instruction Manual

Manual No.: 13947001

Date of Release 01 .€G20F€

# Abramin

## Instruction Manual

---

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

**Instruction Manuals:** Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

**Service Manuals:** Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2010.

**Struers A/S**  
Pederstrupvej 84  
DK-2750 Ballerup  
Denmark  
Telephone +45 44 600 800  
Fax +45 44 600 801

---

## **Abramin Safety Precaution Sheet**

### **To be read carefully before use**

1. The operator should be fully aware of the use of the machine according to the Instruction Manual. The machine must be placed in an adequate working position.
2. Be sure that the actual voltage corresponds to the voltage stated on the back of the machine. The machine must be earthed.
3. Be sure that the water connections do not leak.
4. Keep clear of the rotational disc when operating. Be careful that you do not hurt your fingers/knuckles (the cooling water will prevent you feeling any pain while the injury is done!). Do not drop the sample while grinding.
5. Struers recommend that the mains water supply is shut off or disconnected if the machine is to be left unattended.
6. If you observe malfunctions or hear unusual noises - stop the apparatus and call technical service.

---

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

---

## **Disposal**

Equipment marked with a WEEE symbol  contain electrical and electronic components and must not be disposed of as general waste.

Please contact your local authorities for information on the correct method of disposal in accordance with national legislation.

# User's Guide

Table of Contents	Page
<b>1. Description</b> .....	3
<b>2. Technical Data</b> .....	4
<b>3. Accessories and Connected Equipment</b>	
Accessories.....	5
Grinding with Grinding Paper.....	5
Polishing on Al-Discs.....	5
<b>4. Preparation of Specimens</b>	
Preparation Methods.....	7
<b>PG</b> Plane Grinding.....	7
<b>FG</b> Fine Grinding.....	7
<b>DP</b> Diamond Polishing.....	7
<b>OP</b> Oxide Polishing.....	7
Struers Metalog Guide™.....	8
<b>5. Installation</b>	
Contents of Packing Box.....	9
Unpacking/ Removal of Transport Fittings.....	9
Connection to Compressed-air.....	9
Connection to Mains.....	10
Connection of Water.....	10
Lubricator.....	10
<b>6. Operation</b>	
Controls.....	11
Front Panel and Controls of Abramin.....	12
Mounting the Grinding and Polishing Disc.....	14
Inserting the Specimen Holder.....	14
Programming.....	14
Preparation and Tips.....	17
Manual Preparation.....	18

<b>7. Cleaning</b> .....	19
<b>8. Service and Repair</b>	
Technical Information.....	20
Mechanical Construction.....	20
Pneumatic System .....	21
Electric System .....	22
Trouble Shooting.....	23
Fault in the Power Circuit.....	23
Faults in the Operating Circuit .....	23
Faults in the $\mu$ P -circuit or Display .....	24
Index for Components.....	25
Maintenance.....	26
Adjustment .....	27
Adjustment of Down-speed.....	27
Adjustment of Up-speed .....	27
Adjustment of Slide Friction .....	27
Removal of the Quick-release Coupling.....	27

## **1. Description**

Abramin is a microprocessor controlled machine for the grinding and polishing of specimens fixed in a specimen holder disc. Grinding and polishing are carried out on a horizontal rotating disc and the rotating specimen holder with the specimen is pressed eccentrically against it.

Amount of pressure, lubricant or cooling water can be manually adjusted whereas the microprocessor controls time, disc speed and choice between water and lubricant.

Abramin has 8 program steps into which the above-mentioned data can be individually programmed. The disc and cooling water or lubricant can also be started manually.

With normal use there is, therefore, only the need to select the respective program step, adjust the pressure and start the process with the two push-buttons marked "Start".

## 2. Technical Data

Subject	Specifications				
<b>Electrical Data</b>	<b>Voltage/ frequency</b>	<b>Max. Power Consumption</b>	<b>Cable Specification</b>	<b>Fuse size</b>	
				<b>Min.</b>	<b>Max.</b>
	3 X 200V 50HZ	3.9 A	3x1,5mm <sup>2</sup> + PE	6	16
	3 x 200V 60Hz	3.9 A	3xAWG12 + PE	6	15
	3 X 220V 50HZ	3.8 A	3x1,5mm <sup>2</sup> + PE	6	16
	3 X 220V 60HZ	3.8 A	3xAWG12 + PE	6	15
	3 X 380V 50HZ	2.1 A	3x1,5mm <sup>2</sup> + PE	6	16
	3 X 415V 50HZ	1.9 A	3x1,5mm <sup>2</sup> + PE	6	16
	3 X 440V 60HZ	2.1 A	3xAWG12 + PE	6	15
	3 x 480V 60Hz	2.1 A	3xAWG12 + PE	6	15
<b>Residual Current Circuit Breaker</b>	type A, 30 mA (or higher) is recommended.				
<b>Compressed air supply</b>	6 bar (6 kp/cm <sup>2</sup> )				
<b>Compressed air consumption</b>	Approx. 15 l/min.				
<b>Effect transferred to the disc</b>	300 W at 150 rpm 550 W at 300 rpm				
<b>Effect, specimen holder motor</b>	50 W				
<b>Speed of rotation, turntable</b>	150/300 rpm				
<b>Speed of rotation, specimen holder motor</b>	150 rpm, 50 Hz 180 rpm, 60 Hz				
<b>Vertical specimen holder pressure:</b>	30-450 N				
<b>Polishing disc diameter</b>	Max. 320 mm				
<b>Max. height of specimen outside the specimen holder disc</b>	10 mm				
<b>Dimensions</b>	Height: 650 mm Width: 530 mm Depth: 780 mm				
<b>Weight</b>	99 kg (gross 121 kg)				
<b>Safety Standard</b>	Please refer to the Declaration of Conformity				
<b>Noise level</b>	Approx. 56 dB at a distance of 1 m from the machine				

### 3. Accessories and Connected Equipment

#### Accessories

All of the preparation steps necessary for the preparation of specimens can be carried out on Abramin, e.g. plane grinding, ordinary grinding, diamond disc grinding, polishing with diverse polishing cloths and finally, fine polishing with OP-S. Plane grinding can either be carried out with grinding paper or on the grindstone. In order to be able to use the grindstone on Abramin, however, it is necessary to mount a diamond truing device unit. If large specimens or many specimens are to be prepared, it will be more advantageous to carry out plane grinding on AbraPlan-20 as thereby much time and many pieces of grinding paper can be saved. Abramin can be used with specimen holder disc 160 mm dia. (see the Specimen Holders User's Guide) and in order to get full benefit from the consumables it is recommended that the various accessory discs are chosen in accordance herewith.

Specimen holder 160 mm dia. can be combined with the following preparation discs (Cat. No. in brackets):

#### *Grinding with Grinding Paper*

Ø 305 mm disc (02426909) for non-adhesive grinding paper in accordance with the Knuth-Rotor principle or Ø 230 mm disc (02426908) for quick-adhesion of grinding paper with Double-adhesive discs.

#### *Polishing on Al-Discs*

- Ø 250 mm (02426907) or Ø 300 mm (02426906)
- Microstop for Abramin (03946901).
  - Multidoser for 6 bottles. Automatic dosing, supplying constant amounts of suspension and lubricant.
  - Levelling apparatus Uniforce (04886102 / 04886101) which has its own levelling disc (please ask for separate leaflet).
  - Diamond truing device unit (04146903) for truing the grinding wheel (40800007) for plane grinding.
  - Metason 200-HT ultrasonic bath for cleaning the specimen holder with specimens between the various steps.
  - Drybox-2 drying apparatus for drying the specimen holder with specimens.

*Abramin*  
*Instruction Manual*

<b>Specification</b>	<b>Code</b>
<b>Accessories</b>	
<i>Wet Grinding Discs with Ring</i>	
250 mm dia., aluminium	02426915
305 mm dia., aluminium	02426909
<i>Polishing Discs</i>	
250 mm dia., aluminium	02426907
300 mm dia., aluminium	02426906
<i>Discs for Magnetic Fixation</i>	
MD-Disc, 250 mm dia.	02426919
MD-Disc, 300 mm dia.	02426918
<b>Connected Equipment</b>	
<i>Multidoser (6 Bottles)</i> for automatic dosing of 6 liquids (1-5 suspensions, 0-4 lubricants and 1 OP-Suspension)	04606204
<i>Multidoser Connection kit for Abramin</i>	04606907
<b>Spare Parts</b>	
Swivel	2NF10130

## 4. Preparation of Specimens

### **Preparation Methods**

#### ***PG Plane Grinding***

Plane Grinding (PG) may be carried out on PG-paper, SiC paper or, for hard materials, on Diamond Pad.

SiC paper is used on a wet grinding disc, according to the Knuth-Rotor principle. PG-paper or Diamond Pad is glued to a grinding/polishing disc.

#### ***FG Fine Grinding***

Fine Grinding (FG) is carried out on a cloth with low resilience such as DP-Plan, DP-Pan or DP-Dur, or on Petrodisc-M. SiC-Paper may also be used.

#### ***DP Diamond Polishing***

Diamond Polishing (DP) is carried out on a polishing cloth using diamond abrasive and a lubricant. The correct choice of cloth depends on the material.

#### ***OP Oxide Polishing***

Oxide Polishing (OP) is especially suited for soft and ductile material as a final polishing step.

## Struers Metalog Guide™

Abramin is designed for both grinding and polishing. In Struers Metalog Guide™ you will find a detailed description of grinding/polishing for automated mechanical specimen preparation. Struers Metalog Guide™ offers preparation methods for the most common materials, based on a simple analysis of two key properties: hardness and ductility. Finding the right method is easy, including choice of consumables. Always consult Struers Metalog Guide™ for the correct preparation method for the actual specimens.

Struers Metalog Guide™ contains 6 useful chapters:

Metalogram: a quick and safe guide to the right preparation method.

Metalog Methods: a complete catalogue of preparation methods, based on Struers' vast store of materialographic experience, and employing Struers' range of consumables.

Preparation Philosophy: the basics of modern specimen preparation, seen from a professional point of view.

Metalog Process: the materialographic preparation process from start to finish, logically explained.

Metalog Master: a combined trouble-shooting guide and supply of in-depth information on the processes of mechanical preparation, including an expert system for the solving of preparation problems.

Consumables specification: quick access to the relevant consumables for the chosen preparation methods.

### **Struers Metalog Guide**

A complete guide to materialographic specimen preparation.  
Contact your local dealer for a free copy of Metalog Guide.

## 5. Installation

### Contents of Packing Box

1	Abramin
1	Set of Instruction Manuals
100	Process cards
1	Protection ring of PVC
1	Lid for protection ring of PVC
1	Lubricant bottles, 250 ml
1	Lubricant stopper
1	Outlet tube 1" x 1.5 m
1	Elbow pipe for outlet (PVC)
3	Hose clips, 30 mm
1	Reinforced tube for compressed air, ¼" x 3 m
1	Hose connection for compressed air
2	Hose clips, 12 mm
1	Hose connection
1	Union for hose connection
1	Gasket for hose connection
1	Hose clip, 13 mm
1	Allen key, 6 mm
1	Allen key, 2.5 mm
1	Grease gun for spindle

### Unpacking/ Removal of Transport Fittings

Abramin is secured by 4 bolts to a base frame in the packing case. Abramin can be taken off the base frame, using the pipe wrench (13 mm), after which it can be placed where it will be used.

**Note:**

The Figures referred to in this manual can be found in the Spare Parts and Diagrams section.

### Connection to Compressed-air

Connect the reinforced compressed-air tube to the branch at the back of the apparatus (see fig. 3.© in the Spare Parts manual), use hose clamp (12 mm). Abramin is equipped with a water separator with filter and, therefore, no extra air purification is necessary.

Feed supply pressure should be 6 bar and can be supplied either from the central compressor, portable compressor with compressed air reservoir or compressed-air bottle. A capacity of 15 l/min at atmospheric pressure is sufficient.

### **Connection to Mains**

Check that the serial plate which is placed above the wire (fig. 3.⑦), is marked with the correct voltage and frequency.

Connect the wires (fig. 3.③) as follows:

Yellow/green: earth

Black: phase

Black: phase

Brown: phase

Check that the turntable rotates counter-clockwise when "Disc" is activated. If it rotates clockwise, reverse the direction by changing over two phases. The voltage must not deviate more than 10% from that indicated on the serial plate, otherwise damage can be caused to the electronics or the power can be reduced.

### **Connection of Water**

Connect the thin reinforced hose (fig. 3.④) to the tap and use a 13 mm hose clip to secure it. The supplied hose connection with union and gasket can also be mounted on the tap ( $\frac{1}{2}$ " pipe thread BS 2779:G  $\frac{1}{2}$ ).

Connect the outlet hose to the union on Abramin. The hose is led to the outlet in a gradual slope. This prevents the water from discharging too slowly as this may cause overflow or blocking of the outlet. The elbow pipe must be suitably placed in order to prevent elbows on the hose.

For the outlet hose use hose clip (30 mm). It must not be tightened too much (plastic tubing).

### **Lubricator**

If the built-in lubricator is to be used, the apparatus is ready for use when the 250 ml bottle is full of lubricant, the special stopper is set in and the bottle in position as shown on fig. 1 in the Spare Parts manual.

If Multidoser is to be used, connect it to the multiple socket on the right side of the apparatus.

## 6. Operation

### Controls

Main Switch

Please refer to the Spare Parts manual for figures.

The main switch is located on the left side of the apparatus (fig.3.①). The switch is only used in connection with servicing and can be locked when switched off during maintenance operations. The Stand-by/On switch  on the front panel is used during normal operating procedure.

Fig. 1.③ Start Buttons

The program step starts by pressing the two buttons simultaneously until the specimen holder begins to rotate. If the display shows 00.00 the sample holder is lowered without the motors being started. By activating the key Stop  the sample holder returns to its initial position. This is used in connection with Microstop (accessory).

① Emergency Stop

Interrupts all functions.

The specimen holder remains in lower position until the emergency button is released by pulling the button. When the emergency stop is activated the display shows intermittently E.S. or O.L., signifying Emergency Stop or Overload. (If the Emergency Stop has not been activated and the display nevertheless shows E.S or O.L., one of the two motors in the apparatus is overloaded - see chapter 8, Service and Repair).

⑫ Water Valve

Used to adjust the water flow.

The turning on/off of the water is effected by a built-in magnetic valve which is controlled by the microprocessor.

Needle Valve for Lubricant

For adjusting the flow of lubricant.

Opening and closing of the lubricant valve is controlled by the microprocessor.


Socket for Multidoser

For connection of automatic dosing with Multidoser.

⑤ Quick-coupling

For quick connection and disconnection to the specimen holder.

Front Panel and Controls of  
Abramin

- ① Stand-by/On Switch  The machine cannot operate in Stand-by and all buttons are blocked. Stand-by is indicated by a dot on the display.
- ② Program Step 8 keys for choice of preparation step. The 8 program steps can be programmed individually in regard to time, disc speed and water/lubricant.
- ③ Display Shows the set time in the relative step in minutes and seconds. It is also used to indicate various warnings.
- ④ Data When a program step has been selected the time can be entered with the two pointer keys to the left   and  . The two middle keys are for the selection of either 150 or 300 rpm for the polishing disc. The two keys Water  and Doser  are used for admission of either water or lubricant (lubricant and diamond suspension if Multidoser is connected).
- ⑤ Manual Keys The Stop key  interrupts the program step and the manual functions. The Disc key in the middle (150 rpm/300 rpm) starts the grinding or polishing disc and the Water key  opens the flow of water. By re-pressing Water or Disc, or by pressing Stop, the function is interrupted. If Doser has been selected in the program step the water will be cut off.
- ⑥ Process Card Holder Holder for fixing the process card, a key plan stating which program step should be used, which grinding/polishing backing should be put on the turntable and which pressure should be set (see examples of programming at the back of this manual).
- ⑦ Pressure Gauge For reading the grinding and polishing pressure.
- ⑧ Force (F) Adjusting knob for the setting of grinding or polishing pressure.
- ⑨ Eccentricity Markings By using the handwheel (fig. 3.②) the specimen holder can be moved radially in relation to the centre of the polishing disc so that the grinding paper and polishing cloths can be used to their full extent.

### **Mounting the Grinding and Polishing Disc**

Place the grinding and polishing disc on the turntable by inserting the three legs in the matching holes of the turntable. In order to prevent the polishing disc from wobbling the contact face of the polishing disc and the turntable should be clean.

### **Inserting the Specimen Holder**

Secure the specimen holder with the pre-levelled specimens (see chapter 3.) in the quick-coupling (fig. 1.⑤) in such a way that the stud on the specimen holder engages in the quick-coupling and while the specimen holder is being pressed up, press the quick-coupling's plastic flange down with the thumbs. In this way, the specimen holder moves further up in the coupling. While the specimen holder is still being pressed up, release the plastic flange so that it jumps back into place. The quick-coupling now grips the specimen holder stud and by turning the specimen holder the driving pins on the quick-coupling will engage the specimen holder so that it will rotate with it. (To remove the specimen holder, press the quick-coupling's plastic flange down with the thumbs while the other fingers are ready to 'catch' the specimen holder, which will fall down). Thereafter the quick-coupling with the specimen holder can be adjusted sideways in relation to the disc.

When grinding with grinding paper, adjustment shall be effected in such a way that the grinding paper is used to its full extent.

When grinding, the specimens should move out over the edge. When polishing, as much of the polishing cloth as possible should be used, but in such a way that the specimens do not move out over the edge.

With use of matching diameters for the various discs (see chapter 3), a fixed eccentricity can be used so that adjustment in between each step is avoided.

### **Programming**

Abramin is equipped with 8 program steps into which data such as time, number of revolutions and cooling medium (water or lubricant) can be entered. This presents the opportunity to build a specimen preparation with up to 8 different steps, for example 4 grinding steps and 4 polishing steps. Often the preparation data in the 4 grinding steps will be the same and, therefore, the same program step can be used for all 4 grinding steps.

Explanation		Key 1	Key 2	Key 3	Display: min/sec
<b>Programming of grinding step (example)</b>					
1	Turn on main switch.	<b>1</b>	<b>300 rpm</b>		
2	Change from "Stand-by" to "On". When grinding, operate with 300 rpm and cooling water. Only time needs to be adjusted.				00.00
3	Setting of time. Counts in steps of 5 seconds. At constant pressure change occurs quickly.				00.50
4	Reduction of time.				00.45
<b>Programming of polishing step (example)</b>					
1	Select the program step.	<b>2</b>	<b>300 rpm</b>		00.00
2	Change to 150 rpm		<b>150 rpm</b>		00.00
3	Turn off cooling water and change to Doser.				00.00
4	Set time (1 min 30 sec)				01.30

When the selected preparation step is programmed the various data can be entered on an empty process card (100 cards supplied with Abramin).

The left half of the card is for all necessary information needed when grinding or polishing, such as: which program step to select, which polishing disc or which cloth should be used, which grain size, paper or diamond type should be added and the grinding or polishing pressure to be set.

On the right half enter the data which is programmed in the steps. For example: time, speed and if water or lubricant is to be used.

The cards can be used as a file covering the different materials prepared in the laboratory. Fix the card in the holder on the front plate of Abramin so that the various steps in the preparation can be followed throughout the entire process. To make place for this, fold the card along the thick vertical line.

For grinding, use 300 rpm to get a quick grinding and to make sure that the grinding paper is held as firmly as possible when a grinding disc in accordance with the Knuth-Rotor principle (ROTAL) is used.

For polishing, 150 rpm - which is the same number of revolutions at which the specimen holder rotates - is generally used. This gives a perfect and uniform polishing all over the specimens without directional preferences which reduces the pull-out of inclusions. In the case of routine specimens, polishing can be carried out at 300 rpm to minimise the time. Grinding and polishing speeds are also dependent on the eccentricity between the axis of the turntable and the specimen holder. Large eccentricity gives relatively high speed and thereby a high rate of abrasion. By shifting the quick-coupling with the handwheel (fig. 3.②) the relative speed can be varied from 0 (when the turntable rotates at 150 rpm and by placing centre over centre) to max.

## Preparation and Tips

When the preparation is chosen and the individual programs entered, the following procedure should be carried out:

- Select the program step.
- Place the grinding or polishing disc.
- Secure the specimen holder in the quick-coupling.
- Use the handwheel (fig. 3.②) to adjust the eccentricity so that the grinding paper or polishing cloth can be used to full extent.

By using the disc diameter suggested in chapter 3, it is possible to find a permanent position which makes adjustment between the various steps unnecessary.

When grinding on Petrodisc-M and diamond grinding disc, it is important that the entire surface is used. This means that the specimens move out over the edge and in place over the recess in the centre. This prevents uneven wear.

- Adjust the pressure - see prep. examples.
- When polishing: dampen the polishing cloth

When plane grinding it is easiest to use adhesive paper. Double-adhesive Discs can be used for this purpose.

When grinding on ROTAL: pour water on the disc by activating the Manual Water key . The grinding paper and retaining ring can then be positioned.

When using ductile and soft materials which tend to tear the grinding paper in pieces, it is advantageous to place a grain 1000 paper - grains downwards - under the grinding paper to be used or to start the grinding disc and water manually before starting the program step.

- Start the step by simultaneously pressing the double start buttons until the specimen holder rotates.
- When polishing: adjust the lubricant flow on the adjustment screw and add diamonds for polishing.
- When grinding: regulate the water volume at the tap if necessary.
- The remaining time will be shown on the display which counts down to 0. When the time has elapsed the grinding and polishing disc will stop and the specimen holder will return to its initial position.
- The next program step can now be selected - please see the process card.

## **Manual Preparation**

Manual preparation is possible as the grinding- and polishing disc can rotate at both 300 and 150 rpm. The disc starts by pressing Disc (150 rpm/300 rpm) and will operate at 150 rpm or at the pre-set programmed speed as indicated with light diode. Water can be added by pressing Water . If the selected program step is programmed for Doser, the water cannot be turned on and a long beep-tone indicates this at the same time as the display shows Dose.

Disc and water are stopped by pressing Stop or by pressing Disc and Water  once again.

## **7. Cleaning**

Clean the painted surfaces and the plastic foil of the control panel with a damp cloth. The plastic foil must not be cleaned with alcohol.

The outlet tray and hose should be cleaned regularly to avoid blocking of the outlet.

After grinding, it is recommended to wipe-over the parts situated over the turntable before the polishing disc is placed on it. Contamination is thereby avoided.

## 8. Service and Repair

*(Service engineers only)*

### **Technical Information** *Mechanical Construction*

The turntable upon which the grinding/polishing discs are placed, is suspended by two precision ball bearings and a powerful two-speed electric motor drives it through a V-belt drive (fig. 5).

The motor and turntable are both mounted on individual bearing bridges which, in turn, are secured to a steel frame with rubber vibration dampened bolts (fig. 6).

The quick-coupling in which the specimen holder is fixed is built together with a vertically movable spindle with spindle casing to which a gear motor is secured. The spindle casing slides in a cast iron console and is prevented from rotating by a pin which glides in a track in the cast iron console (fig. 4).

The cast iron console with spindle, motor and quick-coupling is fixed to a horizontal traverse which is secured to the turntable's bearing bridge by two steel columns (fig. 6).

The cast-iron console can be moved horizontally along the beam by using a threaded spindle with handwheel (fig. 10).

The spindle casing with the quick-coupling and the motor is moved up and down by two compressed-air cylinders. The weight of the vertical, movable parts is counter balanced by the spiral spring.

*Pneumatic System*  
(see diagram 13940140)

Two compressed-air cylinders control the vertical movement of the spindle casing. The primary side is used for the up-movement and the secondary side for the down-movement. There is only pressure on one side at a time.

Pressure on the secondary side is set manually with the pressure regulating valve on the front of the apparatus and the connected manometer has a scale for compressive force (Newton).

A magnetic valve (which, when in normal position ventilates the cylinders through an adjustable throttle valve) admits air from the pressure regulating valve to the cylinders.

Pressure on the primary side comes direct from the compressed-air supply and a magnetic valve which - when in normal position - ventilates the cylinder's primary side through a throttle valve by which the downward movement speed can be regulated - admits air to the cylinder.

The upward movement speed is limited by a fixed throttle bush in the magnetic valve tube.

When the main switch on the apparatus is not activated, the cylinders - on account of their own weight - will enter middle position where the spring counter-balances the spindle housing with motor.

The compressed air is cleaned through a manually operated water separator with filter (5-8  $\mu\text{m}$ ).

*Electric System*

The electric system is divided into three parts: the microprocessor (CPU, 5 V DC), the interface (24 V AC, control voltage) for activating the relays, valves and for scanning the position-indicator and the high voltage unit (115-500 V) for supply to the main motor and the specimen holder gear motor. This is controlled by contactors with thermal cut-out so that no damage will be caused to the motors by overloading. If a thermal cut-out stops Abramin, the display panel will show "ES" or "OL" and restart cannot be effected before the thermal cut-out has cooled down 3-5 minutes.

The microprocessor circuit board is located vertically behind the membrane switch board. A 4-digit display with control circuit board is situated just behind the front plate. The interface circuit board and the 5 V DC voltage supply are located behind the CPU print. Abramin is equipped with main cut-off switch, which should only be used for servicing purposes as the micro-processors memory for the user's inset programmed data does not have battery back-up. In stead, Abramin has Stand-by which is controlled by a switch on the front plate. In case of power failure the user must re-program the individual steps, which is easy and quickly done on account of the simple programming procedure.

**Trouble Shooting**  
*Fault in the Power Circuit*

The power circuit is directly connected to the three-phased mains voltage.

If there are faults on the main motor or specimen holder motor refer to the following:

Use the "Power-Circuit-Diagram" and "Wiring Diagram". The wires for this circuit are easily found as they are all *black* and have white numerals - as also stated on the two above-mentioned diagrams. Manual activation of contact K1, K2 or K3 with an insulated object will immediately disclose whether the faults are in the operating circuit.

*Faults in the Operating Circuit*

The operating circuit is equipped with 24 V from transformer T1, the primary winding of which is connected to two of the mains phases.

If there are faults on:

Contacts K1, K2, K3 or magnetic valves Y1, Y2, Y3 or activating contacts S2, S3, S4, S5, S6 (See 8.2.3), refer to "24 V Control Voltage Diagram" and "Wiring Diagram".

The wires in this circuit are easily found as they are all *red* with black numerals - as stated on the above-mentioned diagrams.

The horizontal print, I/O and Power Supply Board are equipped with green diodes which all confirm the function printed just at the side of the relative diode on the circuit board.

3 diodes are placed in a separate group and mark the location of:

24 V control current

+5 V (1) supply voltage for  $\mu$ P

+5 V (2) supply voltage for  $\mu$ P

If the diode for 24 V control voltage does not light, change the 1A-fuse F7, which is situated on the right side of the machine, next to the socket for Multidoser. If the diode for the 24 V control voltage lights, check that the diode for activating the relative function lights. Assume, for example, that contact K2 does not function. The diode for K2 lights = the fault is in the contact or in the connection from the interface (with diode) to the contactor. If the diode does not light, the fault is on the  $\mu$ P -board or the I/O and Power Supply Board. Both circuit boards are connected with plug and can quickly be changed. First change I/O and the Power Supply Board.

*Faults in the  $\mu$ P -circuit or Display* The  $\mu$ P -circuit is equipped with 5 V from two parallel coupled stabilisers, 5 V (1) and 5 V (2), which are both protected by a common 1,5 A fuse, F6, which is on the I/O and the Power-Supply Board. Two green diodes indicate the presence of both voltages.

The mentioned circuit board cannot be repaired during servicing and, therefore, the trouble-shooting procedure given here must only be regarded as a guide line.

Fault	Procedure
The display is all black	Check V5 (1) & 5V (2)
The display shows unknown sign	Replace the :P-board
The display and front plate diodes have small defects	Use test routine (see below)
The activating contacts on the diodes on I/O and the supply board show positive, but the function is not carried out	Use test routine (see below)

**Test routine:**

Test all display and front plate diodes separately, without the possibility of any interference via contact activation from outside.

Next, test the RAM-memory without possible interference from outside. The result will appear on the display as GOOD or BAD.

When - after one second - the word GOOD disappears, the test routine is accessible for testing all contacts on the membrane control panel or the activating contacts.

Each activation will result in a singular showing on the display. This will disappear as soon as the contact is released. It must be stated, however, that the motor overload protector contacts F1, F2 and F3 and the emergency stop contact S3 are connected-up in series and normally closed (NC "Normally Closed"). It is, therefore, necessary to activate the emergency stop in order to be able to test all of the other contact functions which are all normally open (NO "Normally Open").

Activation of the test routine:

Press the STOP key on the control panel and *at the same time* set the change-over to Stand-by. After the RAM-test one can, if so required, go out of the test routine and back to the normal program by setting the switch on I (ON).

*Index for Components*

Main (Disc) motor .....	M1
Specimen holder motor .....	M2
Contactor, Disc motor 150 rpm .....	K1
Contactor, Disc motor 300 rpm .....	K2
Contactor, Specimen holder motor .....	K3
Magnet valve, "Air up" .....	Y1
Magnet valve "Air down" .....	Y2
Magnet valve, Water .....	Y3
Magnet valve, Lubricant .....	Y4
Main switch .....	S1
Stand-by/On switch .....	S2
Emergency stop (NC).....	S3
Micro switch, specimen holder motor .....	S4
Two-hand start contact, left.....	S5
Two-hand start contact, right.....	S6
(Micrometer stop) .....	S7

(The apparatus is prepared for connecting)

## **Maintenance**

### **Belt Tension**

Once a year, check the tension of the V-belt from the main motor to the turn table. Access to the V-belt is made by removing the large rear top sheet. When correctly adjusted, a pressure of 15 N in the middle of the V-belt causes the V-belt to move approx. 8 mm. The V-belt is tightened by pushing the motor plate backwards in its track after having loosened it.

### **Air Filter and Water Separator (fig. 6. ⑩)**

The compressed air filter must be cleaned approximately once a year - depending on the quality of the compressed air.

Further, the water separator - which is combined with the filter - should be emptied approximately once a month, depending on the water contents of the air, so that the water does not reach up to the filter. It is not necessary to remove the rear top sheet, just press carefully on the knob at the bottom of the water separator glass. The knob is accessible from the bottom side of the machine.

### **Lubrication**

All ball bearings are sealed and, therefore, require no servicing. The spindle which moves in the cast-iron console is the only part which needs grease.

Lubrication of the spindle is effected through a grease hole behind the front plate. Remove the front plate and press the grease in with a grease gun. Old grease should first be wiped off. Repeat the procedure with the spindle in bottom position.

## **Adjustment**

### *Adjustment of Down-speed*

Adjustment is carried out on throttle valve V3 which is screwed into the magnetic valve Y1's outlet gate. This can be done with a screwdriver through the hole in the left side of the machine (fig. 3.⑧). To reduce the speed, turn right.

### *Adjustment of Up-speed*

It is only possible to make a slight, fine adjustment of the up-speed as the speed is preadjusted by means of a throttle bush which is inserted in the compressed air tube leading to Y1.

The fine adjustment is made on the V4 throttle valve with a screwdriver through the hole in the left side of the machine (fig. 3.⑨).

### *Adjustment of Slide Friction*

In order to avoid vibration when grinding, the cast-iron console is equipped with four adjustment screws which ensure that the console is wobble-free on the horizontal traverse. These screws can be adjusted with a 2.5 mm Allen key without first having to remove the large front cover. Two of the screws are in the cast-iron, in front at the lowest part of the horizontal traverse and the other two are on the top on the reverse side. The screws should be tightened in such a way that it is still possible to move the spindle sideways with the handwheel.

## **Removal of the Quick-release Coupling**

The quick-coupling needs no adjustment. It is screwed on the shaft and is secured with loctite 222 (fig. 4). In order to remove the coupling do as follows:

- Block the shaft by placing a pin through the hole in the spindle pipe (4.⑩) and the shaft (4.⑪)
- Loosen the coupling with a wrench on the coupling-shaft (4.④)

The thread should be re-sealed with loctite when replacing the coupling. Worn parts should always be replaced.

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
*Abramin*  
*Spare Parts and Diagrams*

---

**Always state *Serial No* and *Voltage/frequency*  
if you have technical questions or when ordering spare parts.**

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
#### Diagrams

1.11	Pneumatic diagram .....	13940140
1.12	Survey of printed circuit boards.....	13940153
1.13	Survey of the complete $\mu$ processor circuit .....	13940154
1.14	Survey of inputs to 8279.....	13940165
1.15	Survey of output ports from 8155 .....	13940166
1.16	Keyboard/display board, circuit diagram .....	13940169
1.17	Keyboard/display board, component illustration.....	13940178
1.18	Component illustration.....	13940183
1.19	$\mu$ P-board circuit diagram .....	13940184
1.20	Power circuit diagram (black wiring) 3x220-500V/50-60Hz .....	13940221E
1.21	24V control voltage diagram (red wirings).....	13940224
1.22	Wiring diagram .....	13940227E
1.23	Diagram for cooling unit TRECI/connecting kit ABTRE.....	13940228
1.24	In/out and power supply, component illustration .....	13940245
1.25	In/out and supply board circuit diagram .....	13940246A
1.26	Connecting Diagram Supp. Volt. 100-230V .....	14221013A

<p>Some of the drawings may contain position numbers not used in connection with this manual.</p>
---

Fig. 1 Abramin Front

2	Front top sheet.....	394MP017
5	Quick-release coupling, complete.....	13940038
9	Rear top sheet.....	394MP091
12	Water cock.....	266MP002
3	Push button, Start (S5, S6).....	260MP004
3	Cap for push button Start.....	260MP005
1	Emergency stop push button without switch element.....	381MP439
1	Emergency stop switch element.....	381MP440


Description

- |   |  |   |                        |
|---|--|---|------------------------|
| ① | Emergency stop                                       | ⑪ | Preparation disc       |
| ② | Base unit  | ⑫ | Cooling water          |
| ③ | Double start buttons                                 | ⑬ | Protection ring        |
| ④ | Specimen holder                                      | ⑭ | Multidoser (accessory) |
| ⑤ | Quick-clamping device                                |   |                        |
| ⑥ | Screws for loosening the<br>specimen holder position |   |                        |
| ⑦ | Force selector                                       |   |                        |
| ⑧ | Microstop (accessory)                                |   |                        |
| ⑨ | Rear top sheet                                       |   |                        |
| ⑩ | Handle for adjusting<br>dosing position              |   |                        |


**Fig. 2 Abramin Control Panel**

- 1 Stand-by switch (S2) 394MP040
- 2 Spring for process card.....394MP027
- 3 Touch pad keyboard, complete ..... 13940235
- 4 Bottle, 250ml..... 158MP016A

**Fig. 3 Abramin,  
Seen from the Back**

- |   |                                   |          |
|---|-----------------------------------|----------|
| 1 | Main switch (S1) .....            | 368MP037 |
| 5 | Drain hose connection piece ..... | 394MP075 |
|   | Drain hose 1" x 2m .....          | 383MP027 |
| 6 | Hose connection for V1 .....      | 394MP054 |


**Fig. 4 Spindle Unit**  
Drawing 13940185A

	Quick coupling, complete.....	394MP083
1	Plastic flange .....	394MP002
2	Guide .....	260MP031
3	Driving pin with spring.....	260MP028
4	Bush with balls.....	394MP081
	Ball $\varnothing 7/32$ " (10 pcs.) .....	RBA00055
5	Coupling shaft.....	394MP082
5	Key, type A 6x6x30.....	394MP003
6	Seeger ring U 48 x 1.75.....	260MP032
7	Seeger ring U 25 x 1.2.....	260MP033
8	Screw M5 x 8 (20 pcs.) .....	394MP004
9	Spindle pipe with pin.....	394MP005
10	Pin.....	394MP006
11	Spindle shaft.....	394MP092
12	Flange.....	394MP094
13/14	Bearing (2 pcs.), 6006-2RS.....	394MP009
15	Seeger ring U 30 x 1.5.....	276MP005
16	Toothed coupling .....	276MP090
17	Spacer pipe.....	394MP011
18	Key.....	394MP012
	Spring .....	394MP014
	Clamp for spindle house .....	394MP015
19	Plate frame for gear motor/cylinder .....	394MP021
20	Gear motor.....	394MP039

Abramin  
Spare Parts and Diagrams

Fig. 4, Spindle Unit,  
Drawing 13940185A


**Fig. 5 Turntable Unit**  
*Drawing 13940186*

1	Shaft .....	401MP026
2	Turntable.....	401MP020
3	Hub .....	401MP021
4	Distance bushing, 51 x 2.5 .....	401MP022
5	V-belt pulley .....	401MP027
6/7	Ball bearing, 6205-27.....	401MP023
8	Pressure ring $\varnothing$ 41.6 /25.2 .....	2GF51016
9	Seeger ring U 25 x 1.2.....	2ZL10250
	Nut M20 x 1.5 .....	2TC61520

Abramin  
Spare Parts and Diagrams

Fig. 5, Turntable Unit,  
Drawing 13620002


**Fig. 6 Fixation Beam**  
Drawing 13940187

1	Bottom frame .....	394MP022
	Spill pan .....	394MP023
	Splash guard.....	394MP024
	Vent cover.....	394MP026
2	Rubber support .....	2GB00001
3	Turntable support bridge with screws and bushes.....	362MP006
4	Base plate for column .....	394MP031
5	Motor plate .....	394MP035
6	Motor support bridge with screws and bushes .....	394MP036
7	Spacer .....	394MP037
8	Motor V-belt pulley.....	394MP038
	V-belt .....	362MP038
9	Main motor (state voltage and frequency) .....	362MP018
10	Fixation beam .....	394MP032
11	Column .....	394MP033
12	V-ring .....	394MP034
13	Noise filter.....	394MP045
15	Bottom plate.....	394MP018
16	Filter unit-SMC (V1) .....	381MP189
	Filter glass for 381MP189.....	381MP435
	Quick coupling angle ø8, 1/4" .....	381MP432
	Quick coupling angle ø5, 1/4" .....	2NF10083

Abramin  
Spare Parts and Diagrams

Fig. 6, Fixation Beam,  
Drawing 13940187


**Fig. 7 Mounting Plate for  
Pneumatics/Electronics**  
*Drawing 13940188*

1	Transformer (T1).....	394MP043
2	Contactor (K1, K3).....	376MP110
3	Contactor (K2).....	376MP110
3	Over voltage protection circuit.....	375MP267
4	Interface circuit board.....	394MP046
5	CPU circuit board incl PROMs.....	394MP047
	PROMs, latest version.....	394MP084
6	Display circuit board.....	394MP048
7	Thermo relay (F1), (state voltage and frequency) ..	394MP049
8	Thermo relay (F3), (state voltage and frequency) ..	394MP050
10	Fuse (F6).....	422MP032
11	Microswitch (S4).....	394MP042
12	Reduction valve (V2).....	368MP330
	Quick coupling angle ø5, 1/4".....	13750406
13	Air cylinder (V5, V6).....	394MP056
14	Manometer.....	394MP057
15	Solenoid valve (Y1, Y2).....	394MP058
	Quick coupling, T ø5,(3 pcs).....	2NF10100
	Quick coupling angle ø5, 1/4".....	2NF10083
16	Throttle valve for Y1 and Y2 (V3, V4).....	394MP059
17	Air hose ø5/ø3.2, 2 m, blue.....	375MP178
18	Mounting plate for pneumatic and electronic components.....	394MP020
	Throttle bush.....	394MP060
19	Fuse Holder.....	381MP220
19	Fuse (F7).....	422MP027
20	Auxilliary switch K2.....	381MP237
21	Female plug 24V.....	260MP160


Abramin  
Spare Parts and Diagrams

Fig.7, Mounting Plate for  
Pneumatics/Electronics,  
Drawing 13940188


**Fig. 8 Water Cock**  
Drawing 13940189

- 1 Reinforced hose with union .....394MP061
- 2 Solenoid valve (Y3), (state frequency).....394MP062
- 3 Cock.....362MP013
- 4 Sprinkler pipe.....394MP063

**Fig. 9A Lubricator**  
Drawing 13940006B


1	Cork for lubricator bottle .....	158MP015
2	Lubricator cup .....	158MP014B
3	Valve needle .....	158MP005
4	Lubricant arm .....	158MP008
5A	Drip pipe .....	158MP010
5B	Nylon washer .....	158MP007A
6	Mounting bush .....	394MP065
7	Distance piece .....	394MP066
10	Silicone tube $\varnothing 2/\varnothing 4$ (10 m) .....	14602901
	Solenoid valve (Y4) .....	394MP076

Abramin  
Spare Parts and Diagrams


**Fig. 9B**  
Drawing 13940109A

- |   |                      |          |
|---|----------------------|----------|
| 8 | Lubricant pipe ..... | 394MP067 |
| 9 | Knob .....           | 394MP068 |


**Fig. 10 Handwheel**  
Drawing 13940191

1	Handwheel.....	394MP070
2	Spindle for handwheel .....	394MP088
3	Spindle bearing.....	394MP089
4	Holder plate for spindle.....	394MP073
5	Plastic bearings (2 pcs.) .....	394MP074
6	Spacing nut.....	394MP090
	Thumb screws (for retaining quick-coupling position) .....	276MP028
	Driving pin (for grinding/polishing discs).....	144MP018
	Ring (for 300 mm ø discs for grinding paper) .....	414MP068
	Micrometer for Microstop .....	2GR10101



1.11 Pneumatic Diagram  
13940140

1.12 Survey of  
Printed Circuit Boards  
13940153



1.13 Survey of the Complete  
 $\mu$  Processor Circuit  
13940154

1.14 Survey of Inputs to 8279

13940165


**1.15 Survey of  
Output Ports from 8155**  
13940166

24 V~ OUTPUTS	PA 0 : DISC MOTOR, 150 RPM, K1, CONTACTOR	0,16A
	PA 1 : " " , 300 RPM, K2, " "	"
	PA 2 : SAMPLE MOV. MOTOR , K3, " "	"
	PA 3 : AIR PRESSURE UP , Y1, SOL.VALVE	0,2A
	PA 4 : " " DOWN, Y2, " . "	"
	PA 5 : WATER (RECIRC.) , Y3 (K4), " . "	0,4A
	PA 6 : LUBRICANT (OPTIONAL Y4), " . "	0,2A
	PA 7 : NOT AVAILABLE	

DC OUTPUTS MAX 50V/10mA	PB 0 :	PROGRAM	PRG.1~000
	PB 1 :	NUMBER	PRG.2~001
	PB 2 :	CODE	⋮
			⋮
			⋮
			PRG.8~111
	PB 3 :	DOSER	
	PB 4 :	SPEED SIGNAL , PWM : PERIOD-100ms, RESOL.= 1 ms	
PB 5 :	BEEPER ( INTERNAL )		
PB 6 :	TORQUE SIGNAL		
PB 7 :	NOT AVAILABLE		


PC 0 - PC 7 : NOT AVAILABLE

1.16 Keyboard/Display Board,  
Circuit Diagram  
13940169


Abramin  
Spare Parts and Diagrams

1.17 Keyboard/Display Board,  
Component Illustration  
13940178

PCB 205-P-268/1

1.18 Component Illustration  
13940183


PCB 205 - P - 269/1


# Abramin Spare Parts and Diagrams

## 1.19 $\mu$ P-Board Circuit Diagram 13940184


Abramin  
Spare Parts and Diagrams

1.20 Power Circuit Diagram (Black Wiring)  
3x220-500V/50-60Hz  
13940221C


V / Hz	F1/A	F3/A	F4/A	(Optional)
200 / 50	2.8	0.35	0.4	
200 / 60	2.8	0.35	0.4	
220 / 50	2.7	0.35	0.4	
220 / 60	3.0	0.35	0.4	
380 / 50	1.5	0.22	0.23	
415 / 50	1.3	0.22	0.23	
440 / 60	1.5	0.22	0.2	
500 / 50	1.1	0.18	0.18	
480 / 60	1.3	0.22	0.2	
440 / 50	1.2	0.20	0.2	


1.21 24V Control Voltage  
Diagram (Red Wirings)  
13940224

all wirings is red 1° common lead without number.

V / Hz	F1/A	F3/A	F4/A	(Optional)
200 / 50	2,8	0,35	0,4	
200 / 60	2,8	0,35	0,4	
220 / 50	2,7	0,35	0,4	
220 / 60	3,0	0,35	0,4	
380 / 50	1,5	0,22	0,23	
415 / 50	1,3	0,22	0,23	
440 / 60	1,5	0,22	0,2	
500 / 50	1,1	0,18	0,18	
480 / 60	1,3	0,22	0,2	
440 / 50	1,2	0,20	0,2	

E: T1 sekunder 1.5V → 1.5A; 2.5V → 2.5A 22/3.99 BHJ/77V

D: F1+F3: Term.w. 98 rykket; 97 tilf. 081098 BHJ/77V

C: rd X9:3; rd K1:A2 → rd K3:A2 061196 BHJ/77V

B: F1:95 tilf. rd 2: F7; F7: rd 2 S1 → rd 2 F1:95 77V

A: udvalgt term. bet. v/st. type af br. 130396 BHJ/77V

Matr.: Abramim

**STRUERS K/S**

Anvendes i type: **Wiring Diagram**

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Ikke ang. tolerancer: ±

Målestok:






Navn: *ASB / G. H.*

Dato: *8.6.89*


Erstatter: *3940222 3/91 ASB*

**13940227 E**


Sample Mower Motor

W5 2 x 0.75  
W4 2 x 0.75







1.23 Diagram for Cooling Unit  
TREC/Connecting Kit ABTRE  
13940228

Abramin  
Spare Parts and Diagrams

1.24 Component Illustration,  
In/Out and Power Supply  
13940245

Rev. No.	1
Date	3.7.85
Author	ABRAMIN
Checker	2050322/1
Project No.	13940246 B
Sheet No.	1

**STRUERS K/S**  
IN/OUT AND SUPPLY BOARD  
CIRCUIT DIAGRAM

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



Herewith declares that

<i>Product Name:</i>	Abramin
<i>Type No:</i>	394
<i>Machine Type:</i>	Grinding and polishing machine

is in conformity with the provisions of the following directives:

**Safety of Machinery** 98/37/EC according to the following standard(s):  
 EN ISO 12100-1:2005, EN ISO 12100-2:2005, EN ISO 13850:2006, EN 60204-1:2006.

**EMC-Directive** 2004/108/EC according to the following standard(s):  
 EN 61000-6-1:2002, EN 61000-6-3:2002.

**Low Voltage Directive** 2006/95/EC according to the following standard(s):  
 EN 60204-1:2006.

**Supplementary Information** The equipment complies with the American standards:  
 FCC part 15, subpart B and UL508.

**The above has been declared according to the global method, module A**

Date: 21.06.2010

  
 Christian Skjold Heyde,  
 Vice President, R&D and Production, Struers A/S

**Hersteller** Struers A/S  
Pederstrupvej 84  
DK-2750 Ballerup, Danmark  
Telefon +45 44 600 800

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Christian Skjold Heyde,  
Vice- President, R& D et Production, Struers A/S

# Abramin

**Microprocessor controlled table top machine for automatic grinding, lapping and polishing of all materials**

- Easy to operate
- Eight microprocessor controlled steps
- Short preparation times
- Optimum reproducibility
- Quick-coupling device for specimen holder
- Microstop can be mounted

**Abramin is a universal grinding, lapping and polishing machine for preparation of specimens for inspection in microscope**

The built-in microprocessor facilitates the operation and the programmable methods make the Abramin a versatile tool for the preparation of all materials. Abramin is ideal for quality control facilities with a relatively high sample volume requiring a fast turn-around time. All preparation steps can be carried out on the machine and no additional manual work is needed.

The programming and change of preparation methods is easy and fast because the data selection of the programs is logically organized and the machine is easy to operate.

## **Eight Programmable Preparation Steps**

Abramin is easily programmed on a splash-proof keyboard, on which the functions are indicated by illumination and the preparation time is shown on a display. Eight preparation steps can be programmed by selecting time, rotational speed, water for grinding or lubricant for polishing. These eight steps can be combined in a number of ways to allow preparation of a wide variety of materials without reprogramming.

### **Grinding**

Grinding is carried out with SiC paper in 4 steps. As an alternative, MD-Allegro or MD-Largo can be used for one-step grinding which provides excellent planeness and is even more economical. For grinding sintered carbides and ceramics, MD-Piano is used. All discs are easily interchangeable.

### **Polishing**

Polishing disc and specimen holder of the Abramin rotate with 150 rpm and the

holder is eccentrically positioned in relation to the polishing disc. This prevents a preferential direction during preparation and ensures a fast and good polish without creating comet tails and pull-outs. For diamond polishing, 9-1/4  $\mu\text{m}$  is used on Struers MD- or DP-Cloths. Abramin can perform any type of preparation from a fast one-step polish to a perfect plane polish. For final polishing, without deformation or scratches, oxide polishing suspensions on MD/OP-Chem cloth are used.

### **The Abra System**

Abramin accepts the same specimen holders, grinding and polishing discs, and consumables as the other Abra machines and TegraPol/TegraForce. The accessories are all interchangeable on these machines, making the Abra System extremely flexible and allowing for an increase of capacity at any time, if required.

### **Microstop**

Abramin can be mounted with a microstop so that a certain removal from the specimen can be preset.

Insertion of the specimen holder is quick and uncomplicated

All preparation parameters are easily programmed on the touch-pad keyboard

## Design

Abramin has a pneumatic system for an automatic vertical movement of the specimen holder. The force applied to the specimens is continuously variable from 30 to 450 N, so that pressure sensitive as well as large specimens can be prepared.

A built-in water separator with filter ensures effective valve function. The 160 mm (6.3") dia. specimen holder discs are easy to handle and 3 to 12 specimens can be prepared simultaneously.

The specimen holder is easily clamped in the quick-coupling device, and preparation discs up to 305 mm (12") are placed directly on the turntable without using any tools.


The specimen holder can be displaced horizontally in relation to the grinding and polishing disc with a hand wheel to allow optimum utilisation of grinding discs and polishing cloths.

Abramin conforms with the strict safety standards laid down for machine tools (IEC 204, VDE 0113), and is easy to service.

## Accessories

The Multidoser applies DP-Suspensions, OP-Suspensions and Lubricants automatically, thus saving money as well as time. (See special leaflet)

The microstop ensures a preset stock removal.

## Technical data

	3-phased Supply fuse: Max. 12 A Compressed air supply: 6 bar
<b>Power:</b>	Main motor: 300/550 W at 50 Hz, 370/650 W at 60 Hz Specimen holder motor: 50 W Specimen holders: 160 mm (6.3") dia.
<b>Rotational speed:</b>	Turntable: 150/300 rpm Specimen holder: 150 rpm at 50 Hz 180 rpm at 60 Hz
<b>Dimensions:</b>	Height: 650 mm (25.6") Width: 555 mm (21.5") Depth: 780 mm (30.7") Weight: 99 kg (218 lb)

## Specifications

	Code
Abramin, Microprocessor controlled bench top machine for automatic grinding and polishing. Specimen holders and discs are ordered separately	ABMIN
Recirculation cooling unit	TRECI
MD-Disc, for magnetic fixation of MD-Consumables. 300 mm (12") dia.	DEMAL
Preparation Disc, for consumables with adhesive back. Aluminium. 300 mm (12") dia.	DEDAL
Wet-grinding disc, aluminium, for non-adhesive SiC paper, 300 mm (12") dia.	ROTAL
Microstop for Abramin	ABCIM
Please state voltage and frequency.	

Struers' equipment is in conformity with the provisions of the applicable International Directives and their appurtenant Standards. (Please contact your local supplier for details)

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice.

**Struers A/S**  
Pederstrupvej 84  
DK-2750 Ballerup, Denmark  
Phone +45 44 600 800  
Fax +45 44 600 801  
struers@struers.dk

**USA and CANADA**  
**Struers Inc.**  
24766 Detroit Road  
Westlake, OH 44145-1598  
Phone +1 440 871 0071  
Fax +1 440 871 8188  
info@struers.com

**SWEDEN**  
**Struers A/S**  
Smältvägen 1  
P.O. Box 11085  
SE-161 11 Bromma  
Telefon +46 (0)8 447 53 90  
Telefax +46 (0)8 447 53 99  
info@struers.dk

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**Struers GmbH**