Discotom-5

Instruction Manual

Manual No.: 15117001

Date of Release 01.06.2004

Discotom-5
Instruction Manual

Table of Contents	Page
User's Guide	1
Reference Guide.....	14
Quick Reference Guide.....	26

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2004.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Discotom-5

Safety Precaution Sheet

To be read carefully before use

1. The operator(s) should be fully instructed in the use of the machine and its cut-off wheels according to the Instruction Manual and the instructions for the cut-off wheels.
2. The machine must be placed on a safe and stable support table. All safety functions and guards of the machine must be in working order.
3. Use only intact cut-off wheels. The cut-off wheels must be approved for min. 50 m/s.
4. Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
5. The workpiece must be securely fixed in the quick-clamping device or similar. Large or sharp workpieces must be handled in a safe way.
6. The cutting handle should be lowered slowly and carefully, in order to avoid breaking the cut-off wheel.
7. Struers recommends the use of an exhaust system as the materials being cut may emit harmful gasses or dust.
8. The rubber curtains on the left hand side make it possible to cut long workpieces with a height of up to 80mm. Protruding workpieces up to Ø25mm can be cut.
9. Do not put your hand through the rubber curtain during operation of the machine.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.)

User's Guide

Table of Contents	Page
1. Getting Started	
Checking the Contents	2
Recirculation Cooling Unit.....	2
Unpacking Discotom-5.....	2
Placing Discotom-5	2
Getting Acquainted with Discotom-5.....	3
Discotom-5 Rear View	4
Noise Level.....	4
Power Supply	4
Direction of the Cut-off Wheel	4
Connecting the Recirculation Cooling Unit.....	5
2. Basic Operations	
Using the Controls.....	6
Front Panel Controls of Discotom-5	6
Front Panel Controls	7
Cooling Valve Positions	8
Cutting Position.....	8
Cleaning Positions.....	8
Light Diode and Beep Signals.....	9
Fitting or Changing the Cut-off Wheel.....	10
Clamping the Workpiece	10
Starting/Stopping the Cutting Process.....	11
Automatic Cutting (AUTO mode).....	11
Manual Cutting	12
Combining Manual and Auto Operation	13
Cleaning	13
Flushing the Cutting Chamber.....	13

1. Getting Started

Checking the Contents

In the packing box you should find the following parts:

- 1 Fork spanner (24mm), for changing the cut-off wheel
- 1 Connector pipe for water outlet
- 1 Set of Instruction Manuals

Recirculation Cooling Unit

- 1 Cooling water tank with lid
- 1 Recirculation pump
- 1 Recirculation sieve
- 1 Trolley
- 1 PVC hose ½", 2 m
- 1 Outlet hose, 2 m
- 2 Hose clamps
- 1 Large screw hose clamp
- 2 Disposable plastic inserts

Unpacking Discotom-5



Remove the screws underneath the pallet to free the Discotom-5.

Placing Discotom-5

Discotom-5 should be placed on a table that is able to carry a weight of min. 125 kg.

**Getting Acquainted
with Discotom-5**


Take a moment to familiarise yourself with the location and names of the Discotom-5 components.

(Right side view)

- ① Main power switch
- ② Instrument panel
- ③ Cut-off wheel locking handle
- ④ Cutting handle
- ⑤ Protection guard, with window
- ⑥ Cutting table
- ⑦ Quick clamping device and spring clamp
- ⑧ Cooling valve
- ⑨ Flush hose inlet

Discotom-5 Rear View

- ① Electrical cable for power supply
- ② Power connection for recirculation pump
- ③ Reset button for recirculation pump
- ④ Water inlet
- ⑤ Water outlet
- ⑥ Overflow drain
- ⑦ Exhaust hose
- ⑧ Cut-off wheel locking handle
- ⑨ Upper back plate

Noise Level

Approx. 67 dB (A) measured at idle running, at a distance of 1.0 m/39.4" from the machine.

Power Supply

- First check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.
- Discotom-5 is factory mounted with an electric cable. Mount a plug on the cable:
Black and brown = phase
Yellow/green = earth

Direction of the Cut-off Wheel

Check that when the power is turned on, the cut-off wheel rotates in the direction indicated by the moulded arrow in the guard of Discotom-5. If the direction of rotation is incorrect, switch the two phases.

Connecting the Recirculation Cooling Unit

Connect the Recirculation cooling unit as follows:

- To connect the electric cable for the recirculation pump, remove the upper back plate on Discotom-5 and follow the instructions on the wiring diagram found on the inside of the back plate.
- Check that when the power is turned on, the pump rotates in the direction indicated by an arrow on the top of the pump. If not, re-check the connection of the electric cable.
- Replace the back plate.
- Clamp the inlet hose to the pump and attach the other end to the water inlet on Discotom-5 ④.
- Attach the connector pipe to the water outlet ⑤, mount the large outlet hose and secure with the large hose clamp.
- Place the tank on the trolley.
- Place a disposable plastic insert in the tank and fold it over the rim.
- Fill the tank with water. The water level should be 8-10 cm below the upper rim of the tank.
- Add Struers Additive according to the instructions.
- Fit the tank lid and sieve.
- Roll the trolley into place and lead the outlet and overflow drain hoses ⑤ & ⑥ into the cooling water tank. Adjust the lengths if necessary.

Reset Button


If the recirculation pump is overloaded or blocked, a thermal cut-out will switch the pump off. This safety cut-out must be reset by pushing the red button on the back of the Discotom-5 ③ before the machine can be restarted.

IMPORTANT

- ① Too high a level of coolant in the tank might damage the pump. To avoid this, place the disposable insert so that the pre-punched hole is in front of the overflow aperture in the tank.
- ② Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the Additive container). Remember to top up with Struers Additive each time you refill with water.

2. Basic Operations








Using the Controls Front Panel Controls of Discotom-5

Control name	Symbol
① Main power switch	
② Power indicator light	Power
③ Start key	
④ Stop key	
⑤ Auto key	Auto
⑥ Auto mode indicator light	Auto
⑦ Forwards button.....	
⑧ Backwards button	
⑨ Feed Speed potentiometer	
⑩ Feed Speed Reduced indicator light.....	
⑪ Load bar	
⑫ Emergency stop button.....	
⑬ Emergency stop/overload indicator light	

Discotom-5
Instruction Manual


Front Panel Controls

Name	Key	Function	Name	Key	Function
① MAIN SWITCH		The main power switch is located on the right hand side of Discotom-5. Turn clockwise to switch on the power.	⑦ FORWARDS		Moves cutting table forward.
② POWER INDICATOR		Lights when main power is switched on.	⑧ BACKWARDS		Moves cutting table backwards.
③ START		Starts the machine. The cut-off wheel starts rotating and the cooling water is turned on.	⑨ FEED SPEED		Regulates the feed speed of the cutting table. The speed can be set from 0.1-2.5 mm/sec.
④ STOP		Stops the machine. The cut-off wheel stops rotating and the cooling water is turned off.	⑩ FEED SPEED REDUCED		Lights when the feed speed has been set too high and is reduced by the machine.
⑤ AUTO		Select or deselect automatic cutting mode	⑪ LOAD BAR		Reflects the present motor load status.
⑥ AUTO INDICATOR		Lights when automatic cutting mode selected.	⑫ EMERGENCY STOP		Push the red button to stop. Pull the red button to release.
			⑬ EMERG. STOP OVERLOAD		Lights when emergency stop has been pressed or the cutting motor is overloaded.

Cooling Valve Positions


In the back of the cutting chamber the flush hose, and a valve for cooling and cleaning, are located. During cutting the valve is placed in the vertical position to cool the cut-off wheel. During cleaning the valve is turned clockwise to regulate the water flow of the hose.

Cutting Position


Turn the knob counter-clockwise to vertical position to cool the cut-off wheel.

Cleaning Positions

Turn the knob to regulate the flow while flushing the cutting chamber.





Turn the knob clockwise to vertical position for flushing the cutting chamber with maximum flow.

Discotom-5
Instruction Manual

Light Diode and Beep Signals

Discotom-5 is provided with various light diodes indicating the status of the machine. Beep signals will sound if a key is pressed when the function cannot be activated.

Status	Indication	Comments
Green POWER indicator on	The main power switch has been turned on.	
White AUTO indicator on	AUTO mode is active.	
🔔 Beep signal	You pressed AUTO, START  , FORWARDS  or BACKWARDS  while the motor was running.	
	You pressed START  while the protection guard was open.	
	You pressed START  for cutting while cooling valve isn't in cutting position.	Turn cooling valve to cutting position.
Red  indicator on	The emergency stop has been pressed.	Release the emergency stop to shut off the  indicator.
Red  indicator on and 🔔 beep signal	The emergency stop has been pressed and you tried to press another key.	
Yellow  FEED SPEED REDUCED indicator on	In AUTO mode: the feed speed is too high.	The machine automatically turns the feed speed down.
	An AUTO cut has finished. Feed speed reduced.	Turn the feed speed down or press any key.
Red  and yellow FEED SPEED REDUCED indicators on	Cutting in AUTO mode not possible as the cutting table does not move.	Press any key to stop the indicators. Material may be too hard for the cut-off wheel chosen or cutting table may be jammed. See Trouble Shooting section at back of manual.
Red  indicator and 1 red cell on LOAD BAR	Pump motor overloaded.	Stop the machine. Push the red button on the back of the machine to reset the safety cut-out.
Red  indicator and 2 red cells on LOAD BAR	Main motor overloaded. Process halted.	The indicators will disappear when the motor has cooled down.
Red  indicator, 2 red cells on LOAD BAR and 🔔 beep signal		You tried to press a key before the motor had cooled down.

Fitting or Changing the Cut-off Wheel

- Note: The spindle for Discotom-5 is left-hand threaded. If Discotom-5 has been used for manual cutting, fasten the cutting handle with the black handle at the right hand side of the machine before changing the cut-off wheel.
- Press and hold down the black locking knob on the right-hand side of the cut-off wheel, turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner. Remove the flange and the old cut-off wheel.
- Mount the new cut-off wheel.

IMPORTANT

Conventional cut-off wheels based on Al_2O_3/SiC abrasives should be placed between two cardboard washers, to protect the cut-off wheel. For maximum precision with diamond or CBN cut-off wheels, do not use cardboard washers.

- Mount the flange and the nut. Tighten carefully and release the locking knob.

Clamping the Workpiece

- Place the workpiece between the quick clamping device and the spring clamp, and the back stop.
- Push the clamps towards the workpiece and lock the quick-clamping device with the locking handle.

IMPORTANT

The workpiece must be firmly held in the quick-clamping device.

- Carefully close the protection guard.

Starting/Stopping the Cutting Process

Automatic Cutting (AUTO mode)

- When the machine is used in AUTO mode the cut-off wheel is stationary and the cutting table moves. In MANUAL mode the cutting table is stationary and the operator moves the cut-off wheel.
- Clamp the workpiece.
- Loosen the black handle situated on the right hand side of the cutting chamber.
- Lower the cut-off wheel by drawing the cutting handle downwards until the cut-off wheel is positioned ready to cut the workpiece. Fasten the black handle again.
- Press ▲ FORWARDS to position the cutting table and the workpiece just in front of the cut-off wheel.
- Close the protection guard.
- Set the desired feed speed
- Press the AUTO mode key to select AUTO mode.
- Press START ◊ to start the machine. The cut-off wheel starts rotating and the cooling water is turned on. The cutting table moves forward toward the cut-off wheel at the pre-set feed speed.
- If the machine has been stopped with the STOP ⊙ key, the cutting table will stay in position. To return the cutting table to its start position press BACKWARDS ▼ once. Struers recommends manually lifting the cut-off wheel away before the cutting table is moved back. Please note that the protection guard has to be closed in order to carry out this operation. If the protection guard has been opened press and hold the BACKWARDS key ▼ to move the cutting table to its start position.
- If the cutting table reaches the back of the cutting chamber during the cutting process, the cut-off wheel automatically stops and the cutting table returns to its start position.

Note: the cutting table can be positioned with FORWARDS ▲ and BACKWARDS ▼.

Reduced Feed Speed

If the feed speed is set too high, the yellow FEED SPEED REDUCED indicator (△) above the feed speed control flashes and the feed speed is automatically reduced by 10%. Discotom 5 can repeat this reduction up to 5 times. Once the workpiece has been cut, the operator can reduce the feed speed using the feed speed control until the overload indicator stops flashing. Discotom-5 will then be able to cut another similar workpiece without feed speed reduction.

Manual Cutting

- Clamp the workpiece.
- Loosen the small black handle on the right side of the cutting chamber. The cutting handle is now easily moved up and down.
- Position the cutting table and the workpiece under the cut-off wheel with the FORWARDS ▲ and BACKWARDS ▼ keys
- Close the protection guard.
- Make sure that AUTO MODE is de-selected.
- Press START ◇ to start the machine. The cut-off wheel starts rotating and the cooling water is turned on.
- Pull down the cutting handle and let the cut-off wheel work itself into the workpiece. Increase the pressure and begin cutting.
- When the cut-off wheel has almost cut through the workpiece reduce the pressure.
- When the cut-off wheel has cut through the workpiece push the cutting handle back to its top position.
- Press the Ⓞ STOP key to stop the machine

Remember

When cutting manually, the motor load bar should be used to monitor the force on the workpiece. Try to keep the motor load in the green or the beginning of the yellow area.


Combining Manual and Auto Operation

Manual and Automatic cutting mode may be used in combination. You may start in manual mode and cut into the workpiece. If you lock this position you can continue in automatic mode, and the workpiece will move towards the cut-off wheel.

Cleaning

Flushing the Cutting Chamber

To ensure a longer lifetime for your Discotom-5, Struers strongly recommends that you clean the cutting chamber with the flush hose every day.

- Press STOP  to stop the cut-off wheel and cooling water.
- Open the protection guard and release the workpiece from the quick-clamping devices.
- Pull out the flush hose from the back of the cutting chamber and point it towards the bottom of the cutting chamber.
- Turn the cooling valve to cleaning position (see sketch at left) and adjust the flow as needed.
- Clean the cutting chamber thoroughly.

- **Note:** Do not clean the protection guard with the flush hose; use a damp cloth instead.

- Turn the cooling valve to cooling position.
- Leave the protection guard open to let the cutting chamber dry completely.

AVOID RUST!
Leave the protection guard open to let the cutting chamber dry completely.

Reference Guide

Table of Contents	Page
1. Advanced Operations	
Adjusting the Back Stop.....	15
Clamping Irregular Workpieces.....	15
Cutting Protruding Workpieces.....	15
Mounting an Adjustable Stop.....	15
Connection to an External Exhaust System.....	15
Mounting an Adjustable Cutting Table.....	16
Optimising the Cutting Results.....	17
Explanation of Safety Factors.....	17
2. Accessories.....	18
3. Consumables	
Cut-off Wheels.....	19
Other Consumables.....	19
4. Trouble-Shooting.....	20
5. Maintenance	
Daily Service.....	23
Weekly Service.....	23
Monthly service.....	23
Replacing the Cooling Water.....	23
Checking the Recirculation Cooling Unit.....	24
Maintenance of Cutting Tables.....	25
Maintenance of Cut-off Wheels.....	25
Storing Bakelite Bonded Al ₂ O ₃ Cut-off Wheels.....	25
Maintenance of Diamond and CBN Cut-off Wheels.....	25
Maintenance of Clamping Devices.....	25
6. Technical Data.....	26

1. Advanced Operations

Adjusting the Back Stop

- Use a 13 mm spanner to loosen the back stop.
- Place the workpiece. Adjust the back stop so that the centre of the workpiece is under the centre of the cut-off wheel:
- Tighten the back stop.
- If you also need to adjust the clamping devices, loosen them and position them parallel to the back stop. Tighten the clamping devices.

Clamping Irregular Workpieces

Irregular workpieces without plane clamping surfaces must be clamped using special clamping tools, as the workpieces must not move during cutting. This could result in damage to the cut-off wheel or to the workpiece itself. For this purpose a vertical clamping device with one flat clamping shoe is available. (See Accessories).

Cutting Protruding Workpieces

In the cutting chamber there is a hole in the right side allowing protruding workpieces to be cut. The workpiece can have a dia. of up to 25 mm (1").

Mounting an Adjustable Stop

If you need to repeatedly cut samples of the same length, an adjustable stop may be mounted. The adjustable stop is mounted to the right of the cutting table and allows samples of up to 57 mm length to be cut. (See Accessories).

Connection to an External Exhaust System

Struers recommends the use of an exhaust as workpieces may emit harmful gasses or dust when cut.

On the back of Discotom-5 you will find a joint for a 50 mm dia. exhaust hose.

Mount an exhaust hose from your local exhaust system onto the joint.

Mounting an Adjustable Cutting Table

If you need to cut plane parallel sections, an adjustable table may be mounted on top of the existing cutting table. The adjustable table positions the workpiece at a right angle to the cut-off wheel and is movable up to 60mm. (See Accessories).

- Dismount the clamping tools, if mounted.
- Move the cutting handle to the top position and the cutting table to the back in order to create working space.
- Move the centre part of the adjustable cutting table to the left by turning the handle until the holes for the Allen screws are visible.
- Loosely mount two T-nuts underneath the adjustable cutting table using Allen screws.
- Slide the adjustable cutting table onto the cutting table:
 - Position the adjustable cutting table on the left-hand side of the cutting table with the spindle on the left.
 - Open the hinged protection guard to let the spindle pass.
 - The guide pins and T-nuts on the bottom of the adjustable cutting table are slid into the third T-slot from the left; this automatically aligns the adjustable cutting table.
- Tighten the Allen screws for the T-nuts.
- Mount clamping tools using the T-slots in the adjustable cutting table.

Optimising the Cutting Results

The following table shows answers to a number of common questions about how to achieve better quality cuts:

Question	Answer
How can I avoid discoloration or burning of the sample?	Reduce the feed speed
	Change the cut-off wheel, as the hardness of the present cut-off wheel may be inappropriate for the hardness of the workpiece*
How can I avoid burrs?	Use a softer cut-off wheel*
	Clamp the workpiece securely with the left-hand quick-clamping device. If two clamping devices are fitted, tighten the right hand spring clamp just enough to prevent the workpiece from shifting on being cut.
How can I prevent the cut-off wheel wearing out too quickly?	Reduce the feed speed or use a harder cut-off wheel*
How can I cut more quickly?	Place the workpiece in a position which allows the cut-off wheel to cut the smallest possible cross-section. Increase the feed speed.

*) Please refer to the Selection Guide in the Struers [Cut-off Wheels brochure](#).

Explanation of Safety Factors

The protection guard has a safety switch to prevent the cut-off wheel from starting while the guard is open. Furthermore, a locking mechanism prevents the operator from opening the protection guard until the cut-off wheel stops spinning.

2. Accessories

Specification	Code
Adjustable cutting table for cutting parallel sections	TRETA
Adjustable stop	TREST
Quick clamping device, left	TRELQ
Spring clamp, right	TRERI
Vertical clamping system with one flat clamping shoe	TREVS
Riser block for vertical clamping system (TREVS)	TREKS
Set of 4 swivel shoes for vertical clamping system (TREVS)	TREVI
Recirculation cooling unit, 65 l	TRECA
Replacement stainless steel bands 1 pc. 40 mm and 1 pc. 20 mm 4 pcs. 40 mm 3 pcs. 20 mm	TRETT TREUT TREAT

3. Consumables

Cut-off Wheels

Please refer to the Selection Guide in the [Struers Cut-off Wheels brochure](#)

Other Consumables

Specification	Code
<i>Additive for Cooling Fluid</i> 1 l 5 l	ADDUN ADDFI
<i>Bacterial Cleaning Liquid</i> 0.5 l	BACTE
<i>Disposable Plastic Insert</i> for collecting waste material in the Recirculation cooling unit. Package with 10 pcs.	EXOSP

4. Trouble-Shooting

Error	Explanation	Action
Indications on the Front Panel		
FEED SPEED REDUCED \triangle and ∇ indicator lights come on.	Emergency stop active.	Solve the problem causing the emergency stop, for example, by reducing the feed speed, or using a different cut-off wheel, then release the emergency stop.
	Cutting motor overloaded due to heavy use.	Open the protection guard and let the cutting motor cool down for about 5 - 10 min.
Machine Problems		
The cut-off wheel does not rotate.	The protection guard may be open.	Press the safety guard completely down. If this does not help, call a Struers service technician.
Chamber light does not work.	Replace the lamp.	Remove the plastic cap to gain access to the fluorescent lamp. Pull out the lamp and replace it.
Water leaking.	Leak in a Recirculation Cooling Unit hose.	Check the hose and tighten the hose clamp.
	Water overflow in the cooling water tank.	Remove the excess water in the tank.
Workpieces or cutting chamber rusty.	Insufficient additive in cooling fluid.	Use Struers Additive for cooling fluid in the cooling water, at the correct concentration. Check with a refractometer. Follow the instructions in the 'Maintenance' Section.
	The protection guard is left closed after use.	Leave the protection guard open to let the cutting chamber dry.
Protection guard window blurred.	Insufficient cleaning.	Clean with mild soapy water (do not use the cleaning hose).
Quick-clamping device unable to hold the workpiece.	The quick-clamping device is not balanced.	Adjust the screw underneath the clamping column. Use a 3 mm Allen key.
	Clamp worn.	Call a Struers service technician.

Discotom-5
Instruction Manual

Error	Explanation	Action
Cutting Problems		
Discoloration or burning of the workpiece.	The hardness of the cut-off wheel is inappropriate for the hardness / dimensions of the workpiece.	See Consumables section, under heading 'Cut-off Wheels'.
	Inadequate cooling.	- Check that there is sufficient water in the Recirculation cooling unit. - Check the cooling water hoses.
Unwanted burrs	Wheel is too hard.	See Consumables section, 'Cut-off Wheels'.
	Feed speed too high at the end of the operation.	Reduce the feed speed near the end of the operation.
	Lack of support for the workpiece.	If possible support the workpiece on both sides.
Cutting quality differs.	Cooling water hose clogged.	Clean the cooling water hose and the cooling tube. Check the water flow by turning the cooling valve to cleaning position.
Cut-off wheel breaks.	Incorrect mounting of the cut-off wheel.	- Check that the hole has the correct diameter. - Check cardboard washer on both sides of the cut-off wheels and replace if worn. The nut must be tightened properly.
	Incorrect clamping of the workpiece.	Make sure that only the left quick-clamping device is tight. The spring clamp to the right should only press lightly. Use the vertical clamping system if the workpiece is an irregular shape.
	Wheel is too hard.	See Consumables section, 'Cut-off Wheels'.
	Feed speed is set too high.	Reduce the feed speed.
	Inadequate cooling.	Check that there is enough water in the cooling unit. Check the cooling water hoses.

Discotom-5
Instruction Manual

Error	Explanation	Action
The cut-off wheel wears down too quickly.	The feed speed is too high.	Reduce the feed speed.
	Insufficient cooling.	- Check that there is enough water in the Recirculation cooling unit. - Check the cooling water hoses.
	The cut-off wheel is too soft for the task.	See Consumables section, 'Cut-off Wheels'.
The cut-off wheel does not cut through the workpiece.	Incorrect choice of cut-off wheel.	See Consumables section, 'Cut-off Wheels'.
	Cut-off wheel worn.	Replace the cut-off wheel.
	The cut-off wheel gets caught in the workpiece.	Support the workpiece and clamp it on both sides of the cut-off wheel so that the cut stays open.
The workpiece breaks when clamped.	The workpiece is brittle.	Place the workpiece between two polystyrene plates. Note! Always cut brittle workpieces very carefully.
The sample is corroded	The sample has been left in the cutting chamber for too long.	Remove the sample directly after cutting. Leave the protection guard open when you leave the machine.
	Insufficient additive for cooling fluid.	Use Struers Additive for cooling fluid in the cooling water at the correct concentration. Check with a refractometer. See Maintenance section.

5. Maintenance

Daily Service

- Clean the cutting chamber, especially the cutting table with the T-slots.
- Clean the protection guard window with a damp cloth. Do not use the flush hose.
- Do **not** clean the lamp glass with alcohol. Use a damp cloth.
- Empty and clean the sieve in the Recirculation Cooling Unit.
- Refill the Recirculation Cooling Unit.

IMPORTANT

Always use the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to top up with Struers Additive each time you refill with water.

Weekly Service

- Clean the cutting chamber thoroughly.

Monthly service

Replacing the Cooling Water


- Replace the cooling water in the Recirculation Cooling Unit at least once a month.
- Place a disposable plastic insert in the tank and fold it over the rim.
- **Please Note:** Too high a level of coolant in the tank might damage the pump. To avoid this, place the disposable insert so that the pre-punched hole is in front of the overflow aperture in the tank.
- Fill the tank with 65 l of water. The water level should be 8-10 cm below the upper rim of the tank.
- Add Struers Additive according to the instructions.

Checking the Recirculation Cooling Unit

The cooling unit should be checked for cooling water every week and refilled when the flushing pump cannot reach the cooling water. Remember to add Struers Additive for Cooling Fluid: One part of Additive for 33 parts of water.

The concentration of additive should always lie between 2.7 and 3.3%; the pH should be around 9. Replace the cooling water if pH is ≤ 8.0 .

To check the concentration of additive, use a refractometer (0-10 Brix, see diagram.). Concentration = $2 \times \text{Brix}$. Add Struers Additive for Cooling Fluid if the concentration is lower than 2.7%. Add water if the concentration is higher than 3.3%.

Maintenance of Cutting Tables

The stainless steel bands are also available as spare parts and should be replaced if they become worn or damaged.

Maintenance of Cut-off Wheels

*Storing Bakelite Bonded Al_2O_3
Cut-off Wheels*

These cut-off wheels are sensitive to humidity. Therefore, do not mix new, dry cut-off wheels with used humid ones. Store the cut-off wheels in a dry place, horizontally on a plane support.

*Maintenance of Diamond and
CBN Cut-off Wheels*

The precision of diamond and CBN cut-off wheels (and thus the cut) depends on how carefully the following instructions are observed:

- Never expose the cut-off wheel to a heavy mechanical load, or heat.
- Store the cut-off wheel in a dry place, horizontally on a plane support, preferably under light pressure.
- A clean and dry cut-off wheel does not corrode. Therefore, clean and dry the cut-off wheel before storing. If possible, use ordinary detergents for the cleaning.
- Regular dressing of the cut-off wheel is also part of the general maintenance.

**Maintenance of
Clamping Devices**

Important

It is recommended to thoroughly clean and lubricate the Quick Clamping Device and Vertical Clamping Device at regular intervals.

6. Technical Data

Subject		Specifications	
		Metric/International	US
Cutting motor	Three-phased cutting motor Rotational speed (running idle)	3,7 kW 2850 rpm at 50Hz 3420 rpm at 60 Hz	5 HP 2850 rpm at 50Hz 3420 rpm at 60 Hz
Cut-off Wheel	Diameter x thickness	ø250x1.5x32mm	9.8 x 0.06 x1.25"
	Max. cutting diameter: 80mm: Max. length of cut: 140mm (at a workpiece height of 65mm and width of 205mm) The max. cutting capacity might be less, depending on the workpiece material and the cut-off wheel	140 mm	5.6"
Cutting Table Dimensions	Width Length	265mm 200mm	10.6" 8"
Recirculation Cooling Unit	Approx. flow: Tank volume	12 l/min 65 l	3.2 gallons/min 17 gallons
Dimensions and Weight	Discotom-5 Width Depth Height Weight	674mm 617mm 458mm (110kg)	27" 24.75" 18.33" 243lbs
Environment	Safety standard Dynamic Balancing Noise level	IEC 204-1/EN 60204-1(VDE 0113), 89/392/EEC, 89/339/EEC Max. admissible unbalance U_{per} according to ISO 1940/1, Balance Quality Grade G6.03: $U_{per} = 2 \text{ gmm/kg} = 110 \text{ gmm}$. Approx. 67 dB(A) during cutting, at a distance of 1.0 m/39.4" from the machine.	
Electrical Data	Voltage/frequency: 3 x 200-210V/50-60Hz 3 x 220-230V/50Hz 3 x 220-240V/60Hz 3 x 380-415V/50Hz 3 x 460-480V/60Hz	Max. Load: 21.3A 19.8A 19.8A 11.8A 11.8A	

Quick Reference Card

Clamping the Workpiece

- Place the workpiece between the clamps and the back stop.
- Push the clamps towards the workpiece and lock the quick-clamping device with the locking handle.

Auto Cutting

- Clamp the workpiece.
- Position the cut-off wheel so that it can cut the sample and lock it in place with the black cut-off wheel locking handle on the right of the machine.
- Position the cutting table with ▲ FORWARDS and ▼ BACKWARDS
- Carefully close the protecting guard.
- Set feed speed and Press AUTO.
- Press START ◊ to begin cutting.
- After cutting, the cutting table automatically returns to its start position.
- Open the protection guard and release the clamping device.
- Remove the workpiece.

Manual Cutting

- Loosen the black handle on the right hand side of Discotom-5 so that the cut-off wheel moves freely
- Clamp the workpiece.
- Position the cutting table under the cut-off wheel.
- If AUTO mode is selected, de-select it.
- Press START ◊. Lower the cut-off wheel until it contacts the workpiece and apply pressure until the sample is cut.
- Keep the Load Bar in the green, or first part of the yellow area. Relax pressure near the end of the cut.
- Open the protection guard and release the clamping device.
- Remove the workpiece.

Changing the Cut-off Wheel

- Fasten the black handle to stop the cutting handle from moving.
- Press the black locking knob on the right-hand side of the cut-off wheel, turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner. Remove the flange and the old cut-off wheel.
- Mount the new cut-off wheel.
- Mount the flange and the nut. Tighten carefully.

Flushing the Cutting Chamber

- Pull out the flush hose from the back of the cutting chamber. Point the flush hose towards the bottom of the cutting chamber.
- Turn the cooling valve clockwise to flushing position.
- Clean the cutting chamber.
- Turn the cooling valve counter-clockwise to cooling position.
- Always leave the protection guard open to let the cutting chamber dry.

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Date of Release 01.06.2004

Discotom-5
Spare Parts and Diagrams

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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Some of the drawings may contain position numbers
not used in connection with this manual.

Discotom-5
Spare Parts and Diagrams

Spare part list for Discotom-5

Drawing	Pos.	Spare Part	Cat no:
15110001	50	Discotom-5, complete PVC Hose transparent 13/32" (2 m)	RNU13313
15110010		Cover complete	
		Rubber tightening kit for front of cover	15112901
	20	Sealing list EPDM	
	30	Rubber plate	
	40	Fastening plate	
	50	Nuts M4 (6 pcs)	
		Rubber curtain kit	15112902
	90	Rubber curtain with slits	
	110	Rubber curtain	
	120	Fastening plate	
	125	Screws M4x12 A2	
	150	Lamp glass, complete	15110080
		O-ring (for lamp glass)	2IO18065
		Lamp 9W	2HG10827
15110020		Front plate complete	
		Adjusting button	13770101
	90	Button	
	100	Cover for button	
	110	Arrow	
15110005		Back-and bottom plate, complete	
	110	Buffer \varnothing 15X7/M5x5, Type D	2GS00090
	130	Handle for flushing hose	15110495
	140	O-ring \varnothing 9.3x2. (10 pcs.)	2IO24020
	150	Hose PVC transparent. 5/15" (800 mm)	2NU19310
	160	Hose clamp s11/9	2NS11109
	165	Bushing	2GK90129
	190	Hose PVC reinforced. \varnothing 12/12 (490 mm)	2NU21218
	255	0.25 AT fuse, glass type	2FU13500
	330	Knob for cooling valve	15332904

Discotom-5
Spare Parts and Diagrams

Spare part list for Discotom-5

Drawing	Pos.	Spare Part	Cat no:
15110045		Bottom plate, complete	
	600	Fuse holder	2XS12004
	605	Fuse 3,15AT	2FU14150
	680	Cardboard blotters ø33/85F.xxTRE (25 pcs)	368MP077
	690	Flange	14150030
	700	Nut, M16 Left-hand thread	2TA70160
	730	Washer for spindle ø40/ø16.5x5	14150012
15110030		Y-table, complete	
		Replacement stainless steel bands. 1x20 + 1x40	05116906
	100	Stainless steel band. 40 mm)	
	110	Screws M8x20 A4. Low head (6 pcs)	
	120	Stainless steel band. 20 mm	
	140	Guide pins 6m6x16 (4 pcs)	
		Replacement stainless steel bands. 4x40	05116907
	100	Stainless steel band. 40 mm (4 pcs)	
	110	Screws M8x20 A4. Low head (12 pcs)	
	140	Guide pins 6m6x16 (8 pcs)	
15110035		X-table complete	
		Replacement stainless steel bands. 1x20 + 1x40	05116906
	240	Stainless steel band. 40 mm	
	250	Stainless steel band. 20 mm	
	260	Screws M8x20 A4. Low head (6 pcs)	
	140	Guide pins 6m6x16 (4 pcs)	
		Replacement stainless steel bands. 3x20	05116908
	250	Stainless steel band. 20 mm (3 ps)	
	260	Screws M8x20 A4. Low head (9 pcs)	
	140	Guide pins 6m6x16 (6 pcs)	
		T-Bolt M8 for 10 mm T-slots	15112903
	210	Screw M8x20 A4. Low head	
	220	Nut M8 for 10 mm T-slots	
15110075		Quick clamping device. Left	
		Quick clamping device. Left. Complete	15110075
		Knob for handle	13810088
15110076		Quick clamping device. Right	
		Quick clamping device. Right. Complete	15110076
		Knob for handle	15110982

Discotom-5
Spare Parts and Diagrams

Spare part list for Discotom-5

Drawing	Pos.	Spare Part	Cat no:
14980068		Recirculation tank	
	90	Cooling Tank, 65 l	381MP140
	100	Lid	381MP142
	110	Clamp for sieve	381MP143
	110	Sieve holder	14980073
	110	Sieve bag	315MP053
	110	Sieve bag for tank, complete (Incl. clamp and holder)	315MP052
	115	Red Plastic Cover for tank (ø117)	2GT20117
14980076		Recirculation pump	
		Rec. pump, complete. 3x200-240V/50Hz	381MP408A
		Rec. pump, complete. 3x380-440V/50Hz	381MP408B
		Rec. pump, complete. 3x200-220V/60Hz	381MP409A
		Rec. pump, complete. 3x380-480V/60Hz	381MP409B
		Rec. pump, complete. 3x220-480V/60Hz CSA	13810958
	20	Sieve ø99 for recirculation pump	R4980074
30	Clamp for sieve ø99	2NS28012	


Discotom-5
Spare Parts and Diagrams


Spare part list for Discotom-5

Drawing

Pos.	Spare Part	Cat no:
	Various	
	Socket spanner (13 mm)	2GR00313
	Fork spanner (24 mm)	2GR00224
	Vertical clamping system for 10 mm T-slots	
	TREVS	
1	Operating key for TREVS	2GR27021
2	Spring for TREVS	2GR27023
3	Flat clamping shoe	2GR26150

**Vertical clamping system for
10 mm T-slots (TREVS)**

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

FT1: Thermal Over Relay: SETTING RANGES

MANS SUPPLY:	FT1: SET =
3x200 - 230V / 50Hz:	0.75 Amp.
3x200 - 240V / 60Hz:	0.46 Amp.
3x380V / 60Hz:	0.43 Amp.
3x380 - 415V / 50Hz:	0.43 Amp.
3x480 - 480V / 60Hz:	0.26 Amp.

- AUTO S1
- START S2
- STOP S3
- TABLE RETURN S5
- TABLE FORWARD S4
- MOVING TABLE POSITION RP1
- WATER VALVE IN COOLING POSITION. BQ1
- WATER VALVE IN SLUCING POSITION. BQ2
- OVERLOAD FROM M1 THERMOSWITCH FT1
- RECIRCULATION IN OVERLOAD. FT1
- SAFETY GUARD OPEN. S2
- EMERGENCY STOP ON. S1
- CUTTING MOTOR M1, CURRENT ALT. T2
- SAFETY GUARD RELEASE YM1

- CUTTING MOTOR
- RECIRCULATION MOTOR
- CUTTING MOTOR M1, NOT IN OVERLOAD.
- RECIRCULATING PUMP NOT IN OVERLOAD
- CUTTING MOTOR RUN
- RECIRCULATING PUMP ON.
- SAFETY GUARD LOCK
- STEPMOTOR, MOVING TABLE
- CUTTING CURRENT:
- AUTO FEED ON:
- EMERGENCY AND OVERLOAD:
- FEED INFORM:
- KEYBOARD BUZZER:

- GREEN Cut. Cur. 12.5%
- GREEN Cut. Cur. 25%
- GREEN Cut. Cur. 37.5%
- GREEN Cut. Cur. 50%
- GREEN Cut. Cur. 62.5%
- GREEN Cut. Cur. 75%
- GREEN Cut. Cur. 87.5%
- GREEN Cut. Cur. 100%
- YELLOW Cut. Cur. 125%
- YELLOW Cut. Cur. 150%
- RED Cut. Cur. 175%
- RED Cut. Cur. 200%
- YELLOW
- RED

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

Herewith declares that

<i>Product Name:</i>	Discotom-5
<i>Type No:</i>	511
<i>Machine Type:</i>	Cut-off machine

is in conformity with the provisions of the following directives:

Safety of Machinery 98/37/EEC according to the following standard(s):
EN292-1:1991, EN292-2:1991, EN418:1992, EN60204-1:1997.

EMC-Directive 89/336/EEC and 92/31/EEC according to the following standard(s):
EN50081-1:1992, EN50082-1:1997.

Low Voltage Directive 73/23/EEC according to the following standard(s):
EN60204-1:1997.

Supplementary Information The equipment complies with the American standards:
FCC part 15, subpart B and UL508.

The above has been declared according to the global method, module A

Date: 01.03.2006

Christian Skjold Heyde,
Vice President, R&D and Production, Struers A/S

Hersteller Struers A/S
Pederstrupvej 84
DK-2750 Ballerup, Danmark
Telefon +45 44 600 800

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<i>Produktname:</i>	Discotom-5
<i>Typennr.:</i>	511
<i>Maschinenart:</i>	Trennmaschine

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Mounting instructions

Contents

The following parts should be in your kit:

<u>Pcs</u>	<u>Name</u>	<u>Cat no</u>
1	Brake box 50 Hz:	15111135
	60Hz:	15111136
1	Box with PROMS (1.5&2.4) V1.5:	15113911
	V2.4:	15113912
2	Cable entries incl. nuts	293MP009
1	Auxiliary switch	2KH00122
4	PT screws for box	2TX04010
1	Wire connector	2XK71617
1	Mounting Instruction.....	15117003
1	Drilling jig drawing.....	15117675

Mechanical preparation

In order to mount the brake kit it is necessary to modify the back plate of the Discotom; it may be easier to drill the holes before travelling to site.

- Make two holes for the brake kit wiring. See fig 1.
- Drill according to drawing 15117675.
- Mount the 2 cable entries (293MP009) on the plate.

Fig. 1


Connecting wires

Important

Before working, disconnect the power to the Discotom.

- Disconnect the 2 existing cables going into the Discotom. Remember to mark their original position or reinsert the wires without mounting the back plate.
- Drill the holes in the back plate according to drawing 15117675 (if not already prepared).
- The brake box has two wires W21 and W22. Remove or lift the machine's cabinet to gain access to the top of the contactors. Cable with red wires (W21) through lower hole. Cable with black wires (W22) into top hole.
- Mount auxiliary contactor (cat no.: 2KH00122). See fig 2


Fig. 2

*Dynamic Brake Box
Mounting Instructions*

- Attach the wires according to diagram 15113475.
- Connect wire red 7 to wire red 17 (located on top of the contactor). See fig 3.

Fig. 3

- Connect wire red 1 to the wire terminal at pos 11. See fig 4

Fig.4

- Connect all other wires according to the diagram.
- Exchange PROM according to the version in the machine. Please return the old and the PROM not used to Struers A/S for recycling.

PROM

Adjusting the current and timer

- Open the brake box.
There are 2 controls on top of the black control module:
Current and time.

Current and time


⚡

Always make sure both dials are turned to below the 9 o'clock position before start of fine tuning.

- Fine tune the current and time so that the current will brake the disc in about 2 sec. and the timer will shut off at the same time as the cover lock activates.


Important
Take care when adjusting the current ! Turning the current dial to more than the 10 o'clock position may destroy the cutting motor!

- Mount the brake box on back of machine with 4 PT screws.

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Control Circuit for Discotom-5
with Dynamic Braking Box.

Tuesday, June 26, 2001

Discotom-5

The powerful and universal automatic cut-off machine for materialographic laboratories

- **Automatic and manual cut-off machine in one**
- **Very high cutting capacity**
- **Electronic feed speed reduction ensures excellent cutting results**
- **Effective cooling system preventing surface burns**
- **Cutting chamber made in corrosion resistant material**

Discotom-5 is a high performance table-top cut-off machine with a strong 3.7 kW (5.0 HP) motor. It is designed for cut-off wheels with a diameter of 250 mm (10"). Discotom-5 can be operated both automatically and manually, making it suitable for a wide range of applications.

Automatic Cutting made Easy

Automatic cutting mode is engaged by a simple press key. After pressing the Start key the cutting table with the clamped workpiece will move forward towards the cut-off wheel at the set feed speed. Before cutting, you can of course position the cut-off wheel to compensate for the wear of the cut-off wheel.

Manual mode for Fast and Uncomplicated Cutting

When operating the cut-off machine manually, the cutting table with the clamped workpiece is easily positioned under the cut-off wheel, again with simple press keys. You then press the Start key and pull the cutting handle to begin cutting.

Very High Cutting Capacity

You can cut workpieces with a diameter of up to 80 mm (3.1") on the Discotom-5. In automatic mode though, the cutting table is movable up to 200 mm (7.9") enabling cutting of very long samples. You can for instance cut a length of 140 mm (5.5") at a workpiece thickness of 65 mm (2.6").

Electronic Feed Speed Reduction

When cutting automatically, you choose a feed speed from 0.1 to 2.5 mm/sec. If the feed speed is set too high and the machine cannot cut through the workpiece, Discotom-5 will automatically reduce the feed speed and continue cutting. This may repeat itself up to 5 times and prevents damage or burns on the sample surface.

Effective Cooling of the Cutting Area

A very important design detail is that the water cooling tubes are integrated into the cut-off wheel guard. The cooling water is applied to the cut-off wheel under high pressure and at a right angle. The cooling water thereby "breaks" the thin layer of air surrounding the spinning cut-off wheel, and the result is excellent cooling of the cutting area. Discotom-5 is available with a large 65 l recirculating cooling unit, ensuring that the cooling water does not overheat.

Monitoring the Load on the Cutting Motor

The front panel includes a LED-bar indicating the load on the cutting motor. The bar is divided into three areas,

a green area showing that the load on the cutting motor is ok, a yellow area indicating that the load is growing towards the maximum and the red area indicating that the motor is overloaded.

In automatic cutting mode the LED-bar will show the load on the cutting motor and if the red area is reached the feed speed will be reduced. When cutting manually the LED-bar gives a good indication of how much force to put on the workpiece.

Corrosion resistant Cutting Chamber

Although Discotom-5 is a high performance table-top cutting machine, it takes up remarkably little room - only 650 mm (26") in width and 580 mm (23.2") in depth.

The cutting chamber, made of corrosion resistant materials, is provided with a 200 mm x 265 mm (7.9" x 10.4") cutting table, which is easily moved up to 200 mm (7.9") by control keys. The cutting table is made of corrosion resistant cast aluminium and is fitted with a top surface of stainless steel bands. These steel bands are easily exchanged if damaged or worn. A built-in, adjustable flush hose makes daily cleaning quick and efficient, and a powerful lamp lights up the cutting chamber allowing you to monitor the cutting process through the large window in the protection guard. An exhaust system may be connected.


Quick and Handy Clamping of the Workpiece

The integrated cutting table has three 10 mm T-slots, facilitating the use of

universal clamping tools. An optional adjustable stop can be used to set the length of cut-off samples when making repetitive cuts.

Adding an Adjustable Stop

If you need to cut samples of a specific length, an adjustable stop is available for Discotom-5. The adjustable stop allows samples of up to 57 mm (2.3") to be cut, and you only have to set the sample length once.

Adding an Adjustable Cutting Table for Parallel Sectioning

An adjustable cutting table is available for mounting on the Discotom-5. It positions the workpiece at a right angle to the cut-off wheel, allowing parallel sections to be cut. The table can move up to 60 mm (2.4").

High Safety Level

Discotom-5 cannot be started while the protection guard is open. During cutting the protection guard is locked, eliminating any risk of injury. A built-in electric, dynamic motor brake stops the cutting motor within 5 seconds, ensuring quick access to the cutting chamber after cutting.

The graph shows the cutting capacity in automatic cutting mode with a vew cut-off wheel and flexible clamping tools. The real cutting capacity depends on the sample material

Technical Data

Cutting Motor	Three-phased motor 3.7 kW (5.0 HP)
Cut-off Wheels	Cut-off wheel size: 250 mm x 1.5 mm x 32 mm (10" x 0.06" x 1.25") Rotational speed (running idle): 2850 rpm at 50 Hz, 3420 rpm at 60 Hz
Cutting Table	Cutting table dimensions: Width: 265 mm Length: 200 mm
Environment	Safety standards: IEC 204-1/EN 60204-1 (VDE 0113), 89/392/EEC, 89/339/EEC Noise level: Approx. 67 dB(A) during cutting, at a distance of 1.0 m (39.4") from the machine.
Dimensions and Weight	Discotom-5 without recirculating cooling unit Width: 674 mm (27") Depth: 617 mm (24.75") Height: 458 mm (18.33") Weight: 110 kg (243 lbs)
Specifications	Code
Discotom-5 Automatic cut-off machine. Automatic feeding with electronic control of feed speed. Ability for manual operation. With cutting table with 10 mm T-slots. Recirculation cooling unit (TRECA) and clamping tools are ordered separately	TREVA
Accessories	
Recirculating Cooling Unit For water cooling. Contents 65 l. With a rolling pallet.	TRECA
Adjustable Stop For setting the length of samples in repetitive cuts in a range up to 57 mm.	TREST
Adjustable Cutting Table For cutting plane parallel sections. It positions the workpiece at a right angle to the cut-off wheel. Movable up to 60 mm. With 10 mm T-slots and exchangeable stainless steel bands.	TRETA
Quick-clamping Device for 10 mm T-slots, Left. Complete with back stop. To be mounted at the left hand side of the cut-off wheel	TRELQ
Spring Clamp for 10 mm T-slots, Right. For retaining the cut-off specimen. Complete with back stop. To be mounted at the right hand side of the cut-off wheel	TRERI
Compact Quick-clamping Device for 10 mm T-slots, Left. To be mounted at the left hand side of the cut-off wheel. Adjustment range without repositioning: 7 mm. To be used together with back stop from (TRELQ)	TRELC
Compact Quick-clamping Device for 10 mm T-slots, Right. To be mounted at the right hand side of the cut-off wheel. Adjustment range without repositioning: 7 mm. To be used together with back stop from (TRERI)	TRERC
Vertical Clamping System for 10 mm T-slots. For clamping irregularly shaped workpieces on the cutting table. Clamping height adjustable up to 58 mm. Complete with operating key and one flat clamping shoe	TREVS
Riser Block for 10 mm T-slots. For elevating Vertical Clamping System (TREV) when clamping high workpieces. To increase the clamping height by 60 mm	TREKS
Swivel Shoes for Vertical Clamping System. Set of 4 multi-shaped swivel shoes for (TREV)	TREVI
Replacement Stainless Steel Bands For cutting table or adjustable cutting table (TRETA) on Discotom-5. Set of 1 pc. 40 mm and 1 pc. 20 mm steel bands.	TRETT

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Struers Ltd.
Struers Inc.