Discotom-50

Instruction Manual

Manual No.: 15137001

Date of Release 26.07.2006

Discotom-50
Instruction Manual

Table of Contents	Page
User's Guide	1
Reference Guide.....	21
Quick Reference Guide.....	33

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manuals may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manuals may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2006.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Discotom-50

Safety Precaution Sheet

To be read carefully before use

1. The operator(s) should be fully instructed in the use of the machine and its cut-off wheels according to the Instruction Manual and the instructions for the cut-off wheels.
2. The machine must be placed on a safe and stable support table. All safety functions and guards of the machine must be in working order.
3. The unit must be installed in compliance with local safety regulations.
4. Use only intact cut-off wheels. The cut-off wheels must be approved for min. 50 m/s.
5. The machine is not for use with saw-blade type cut-off wheels.
6. Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
7. The workpiece must be securely fixed in the quick-clamping device or similar. Large or sharp workpieces must be handled in a safe way.
8. Do not work on or around cutting table when the table is repositioned using the forwards/backwards key (▲ ▼).
9. To achieve maximum safety and lifetime of the machine, use only original Struers consumables.
10. The cutting handle should be lowered slowly and carefully, in order to avoid breaking the cut-off wheel.
11. Struers recommends the use of an exhaust system as the materials being cut may emit harmful gasses or dust.
12. The machine emits only moderate noise. However, the cutting process itself may emit noise, depending on the nature of the workpiece. In such cases, the use of hearing protection is recommended.
13. The machine must be disconnected from the mains prior to any service.

Discotom-50
Instruction Manual

- 14.** Do not put your hand through the rubber curtain during operation of the machine.
- 15.** Use of working gloves is recommended as workpieces may be both very hot and produce sharp edges.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.)

User's Guide

Table of Contents	Page
1. Getting Started	
Checking the Contents of the Packing Box	2
Unpacking Discotom-50.....	2
Placing Discotom-50	2
Getting Acquainted with Discotom-50.....	3
Front View	3
Rear View.....	4
Noise Level.....	4
Power Supply	4
Connecting the Recirculation Cooling Unit (optional).....	5
2. Basic Operation	
Using the Controls.....	6
Front Panel Controls of Discotom-50	6
Front Panel Controls	7
Cooling Valve Positions	8
Cutting position	8
Cleaning positions	8
Flush Hose	9
Cleaning of Flush Hose Nozzle	9
Light Diode and Beep Signals.....	10
Using the Load Bar.....	11
Overload Protection	11
Fitting or Changing the Cut-off Wheel.....	12
Clamping the Workpiece	12
Positioning the Cutting Table	12
Starting/Stopping the Cutting Process.....	13
Automatic Cutting (AUTO mode).....	13
Manual Cutting	14
Combining Manual and Automatic Operation.....	14
Programmable Stop.....	15
Variable Rotational Speed	16
Feed Speed.....	17
Material-specific Cut-off Wheels	18
Using CBN or Diamond Cut-off wheels	19
Cleaning	20
Flushing the Cutting Chamber.....	20

1. Getting Started

Checking the Contents of the Packing Box

In the packing box you should find the following parts:

- 1 Spanner (24 mm) for cut-off wheel
- 1 Spanner (13 mm) for transport bolts
- 1 Connector pipe for water outlet
- 1 Filter for drain
- 1 Set of Service Documentation for Siemens frequency converter
- 1 Set of Instruction Manuals

Unpacking Discotom-50

Remove the screws underneath the pallet to release the Discotom-50.

Placing Discotom-50



Discotom-50 should be placed on a table capable of carrying min. 150 kg.

Discotom-50
Instruction Manual

**Getting Acquainted
with Discotom-50**

Take a moment to familiarise yourself with the location and names of the Discotom-50 components.


Front View

- ① Cutting handle
- ② Instrument panel
- ③ Cooling valve
- ④ Flush hose inlet
- ⑤ Protection guard, with window
- ⑥ Main power switch
- ⑦ Secondary fuse, 3.15AT
- ⑧ Cut-off wheel locking handle

Discotom-50
Instruction Manual

Rear View

- ① Upper back plate
- ② Reset button for recirculation pump
- ③ Cable inlet for recirculation pump
- ④ Electrical cable for power supply
- ⑤ Cut-off wheel locking handle
- ⑥ Overflow drain
- ⑦ Water inlet
- ⑧ Exhaust hose
- ⑨ Water outlet

Noise Level

Approx. 67 dB (A) measured at idle running, at a distance of 1.0 m/39.4" from the machine.

Power Supply

- First check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.
- Discotom-50 is factory mounted with an electric cable. Mount a plug on the cable:
 - Black and brown = phase
 - Yellow/green = earth

Important
Residual Current Circuit-Breaker must be min. 20 mA.

Connecting the Recirculation Cooling Unit (optional)

Connect the Recirculation Cooling Unit (optional), product code TRECA as follows:

- To connect the electric cable for the recirculation pump, remove the upper back plate on Discotom-50 and follow the instructions on the wiring diagram found on the inside of the back plate.
- Check that when the power is turned on, the pump rotates in the direction indicated by an arrow on the top of the pump. If not, switch two of the phases to the pump motor and re-check.
- Replace the back plate.
- Clamp the inlet hose to the pump and attach the other end to the water inlet on Discotom-50 ⑦.
- Attach the connector pipe to the water outlet ⑨, mount the large outlet hose and secure with the large hose clamp.
- Place the tank on the trolley.
- Place a disposable plastic insert (EXOSP) in the tank and fold it over the rim.
- Fill the tank with 63 l water and 2 l Struers Additive. The water level should be 8-10 cm below the upper rim of the tank.
- Fit the tank lid and sieve.
- Roll the trolley into its proper position and place the outlet and overflow drain hoses ⑨ & ⑥ inside the cooling water tank. Adjust the hose lengths if necessary.


IMPORTANT

- ① Too high a level of coolant in the tank might damage the pump. To avoid this, place the disposable insert so that the pre-punched hole is in front of the overflow aperture in the tank.
- ② Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the Additive container). Remember to top up with Struers Additive each time you refill with water.

Re-set button

If the recirculation pump is overloaded or blocked, a thermal cut-out will switch the pump off. This safety cut-out must be reset by pushing the red button on the back of the Discotom-50 ② before the machine can be restarted.








Using the Controls
*Front Panel Controls
of Discotom-50*


Control name	Symbol
① Main power switch	
② Power indicator light	Power
③ Start key	
④ Stop key	
⑤ Auto key	Auto
⑥ Auto mode indicator light	Auto
⑦ Forwards key	
⑧ Backwards key	
⑨ Feed Speed knob	
⑩ Feed Speed Reduced indicator light.....	
⑪ Rotational speed knob	
⑫ Load bar	
⑬ Emergency stop.....	
⑭ Emergency stop/overload indicator light	

Discotom-50
Instruction Manual

Front Panel Controls


Name	Key	Function	Name	Key	Function
① MAIN SWITCH		The main power switch is located on the right hand side of Discotom-50. Turn clockwise to switch on the power.	⑦ FORWARDS		Moves cutting table forwards.
② POWER INDICATOR	Power	Lights when main power is switched on.	⑧ BACKWARDS		Moves cutting table backwards.
③ START		Starts the machine. The cut-off wheel starts rotating and the cooling water is turned on.	⑨ FEED SPEED		Regulates the feed speed of the cutting table. The speed can be set from 0.1-2.5 mm/sec.
④ STOP		Stops the machine. The cut-off wheel stops rotating and the cooling water is turned off.	⑩ FEED SPEED REDUCED		Lights when the feed speed has been set too high and is reduced by the machine.
⑤ AUTO	Auto	Select or deselect automatic cutting mode	⑪ ROTATIONAL SPEED		Regulates the rotational speed of the cut-off wheel. The speed can be set from 1000-3000 rpm.
⑥ AUTO INDICATOR	Auto	Lights when automatic cutting mode selected.	⑫ LOAD BAR		Reflects the present motor load status.
			⑬ EMERGENCY STOP		Push the red button to stop. Pull the red button to release.
			⑭ EMERG. STOP OVERLOAD		Lights when emergency stop has been pressed or the cutting motor is overloaded.

Cooling Valve Positions

Flush hose and
Cooling valve

The flush hose, and a valve for cooling and cleaning, are located on the front of the machine. During cutting, the valve is placed in the vertical position to cool the cut-off wheel. During cleaning the valve is turned clockwise to regulate the water flow of the hose.

Cutting position

Cleaning positions

Please Note

After cleaning, turn cooling valve back to vertical position. START  can not be activated unless the cooling valve is in the cutting position.


Flush Hose

The nozzle of the flush hose is adjustable for additional control of coolant flow. Turn the nozzle clockwise/counter-clockwise as appropriate.

IMPORTANT

Do not move the cooling valve to the cleaning position until the flush hose is pointing inside the cutting chamber.

Cleaning of Flush Hose Nozzle

The flush hose nozzle may collect swarf, inhibiting flow. To clean, unscrew nozzle head as shown, and rinse both parts under clean water.

Discotom-50
Instruction Manual

Light Diode and Beep Signals Discotom-50 is provided with various light diodes indicating the status of the machine. Beep signals will sound if a key is pressed when the function cannot be activated.

Status	Explanation	Comments
Green POWER indicator on	The main power switch has been turned on.	
Green AUTO indicator on	AUTO mode is active.	
🔔 Beep signal	You pressed AUTO, START ◊, FORWARDS ▲ or BACKWARDS ▼ while the motor was running.	
	You pressed START ◊ while the protection guard was open.	
	You pressed START ◊ for cutting while cooling valve isn't in cutting position.	Turn cooling valve to cutting position.
Red STOP Ⓢ indicator on	The emergency stop has been pressed.	Release the emergency stop to shut off the STOP Ⓢ indicator.
Red STOP Ⓢ indicator on and 🔔 beep signal	The emergency stop has been pressed and you tried to press another key.	
Yellow FEED SPEED REDUCED △ indicator on	In AUTO mode: the feed speed is too high.	The machine automatically turns the feed speed down.
	An AUTO cut has finished. Feed speed reduced.	Reduce feed speed or press any key.
Yellow FEED SPEED REDUCED △ and Red STOP Ⓢ indicator lights come on	Cutting motor overloaded due to heavy use: Material is too hard Feed speed setting too high Rotational speed setting too high Cutting table is jammed	Press any key to stop the indicators. Solve the problem causing the machine to stop, for example, by reducing the feed speed, or using a different cut-off wheel, then release the emergency stop.
Red STOP Ⓢ indicator and LOAD BAR in red area	Process halted. Main motor or frequency inverter overloaded or overheated.	Press the STOP key Ⓢ to re-set the warning indicator and wait until the motor or frequency inverter had cooled down.
Red STOP Ⓢ indicator and LOAD BAR in red area and 🔔 beep signal	START was pressed before re-setting the warning indicator.	

Using the Load Bar

Discotom-50 is equipped with an LED load bar, indicating the load on the cutting motor. The bar is divided into three areas; a green area showing that the load on the cutting motor is OK, a yellow area indicating that the load is near maximum and a red area indicating overload of the motor.

In automatic cutting mode, the feed speed is automatically reduced on overload, i.e. if the red area is reached.

When cutting manually, the load bar serves as an indicator of how much force is applied on the sample.

The machine may cut with the load bar in the yellow area for a period of time depending on how hot the motor is when cutting is started:

Starting with a hot motor, cutting will be stopped after a few minutes.

Starting with a cold motor, cutting will be stopped after up to 20 minutes.

For continuous load, the feed speed should be reduced to a green load level.

Overload Protection

Discotom-50 is equipped with an overload protection to prevent damage to the cutting motor and the recirculation pump motor.

Cutting Motor

If the cutting motor has been overloaded, a thermal cut-out will stop the motor.

- Open the protection hood and allow the motor to cool down before resuming cutting.

Recirculation Pump

If the recirculation pump has been overloaded, a thermal cut-out will stop the motor.

- To reset the pump, press the red re-set button on the rear of the machine. Allow the motor to cool down, before resuming cutting.

Fitting or Changing the Cut-off Wheel

IMPORTANT

The spindle of Discotom-50 is left-hand threaded.

- If Discotom-50 has been used for manual cutting, secure the cut-off wheel with the locking handle on the right hand side of the machine before changing the cut-off wheel.
- Press and hold down the black locking knob on the right-hand side of the cut-off wheel, turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner. Remove the flange, cardboard washers and the old cut-off wheel.
- Mount the new cut-off wheel.

IMPORTANT

Conventional cut-off wheels based on $\text{Al}_2\text{O}_3/\text{SiC}$ abrasives should be placed between two cardboard washers, to protect the cut-off wheel. For maximum precision with diamond or CBN cut-off wheels, do not use cardboard washers.

- Mount the flange and the nut. Tighten carefully and release the locking knob.

Clamping the Workpiece

- Clamp the workpiece with the clamping device of your choice e.g. a quick clamping device. Place the workpiece between the clamp and the back stop.
- Push the clamp towards the workpiece and lock the clamp on the quick-clamping device.

Generally, it is recommended to clamp the workpiece as far back in the cutting chamber as possible.

IMPORTANT

Prior to cutting, ensure that the workpiece is held firmly in place by a clamping device.

Positioning the Cutting Table

Prior to cutting, the cutting table is positioned using the FORWARDS ▲ and BACKWARDS ▼ keys.

Starting/Stopping the Cutting Process

Automatic Cutting (AUTO mode)

When the machine is used in AUTO mode the cut-off wheel is stationary and the cutting table moves. In MANUAL mode the cutting table is stationary and the operator moves the cut-off wheel.

- Clamp the workpiece.
- Loosen the cut-off wheel locking handle situated on the right-hand side of the cutting chamber.
- Lower the cut-off wheel by pulling the cutting handle downwards until the cut-off wheel is positioned ready to cut the workpiece. Re-fasten the locking handle.
- Press ▲ FORWARDS to position the cutting table and the workpiece just in front of the cut-off wheel.
- Close the protection guard.
- Set the desired feed speed and rotational speed.
- Select AUTO mode by pressing the AUTO mode key.
- Press START ◊; the cut-off wheel starts rotating and the cooling water starts to flow. The cutting table moves forward towards the cut-off wheel at the preset feed speed.
- When the cutting table reaches the back of the cutting chamber, the cut-off wheel automatically stops and the cutting table returns to its start position.
- If the machine is stopped with the STOP ⊕ key, the cutting table will stay at the present position. To return the cutting table to its start position, wait for the double beep and press BACKWARDS ▼ once.
- If the protection guard has been opened, press and hold the BACKWARDS key ▼ to move the cutting table to its start position.

Please Note

The cutting table will only move when AUTO mode is switched on.

Manual Cutting

Manual cutting may be the preferred cutting mode for simple cutting tasks.

For cutting of extremely hard materials (>1000 HV), manual cutting is recommended.

- Clamp the workpiece.
- Loosen the locking handle on the right side of the cutting chamber. The cutting handle is now easily moved up and down.
- Position the cutting table and the workpiece under the cut-off wheel with the FORWARDS ▲ and BACKWARDS ▼ keys.
- Close the protection guard.
- Ensure that AUTO mode is de-selected.
- Press START ◊ to start the machine. The cut-off wheel starts rotating and the cooling water is turned on.
- Pull down the cutting handle and let the cut-off wheel work itself into the workpiece. Increase the pressure and begin cutting.
- When the cut-off wheel has almost cut through the workpiece reduce the pressure.
- When the cut-off wheel has cut through the workpiece push the cutting handle back to its top position.
- Press the STOP ⊕ key to stop the machine.

Combining Manual and Automatic Operation

IMPORTANT

When cutting manually, the motor load bar should be used to monitor the force on the workpiece.

Try to keep the motor load in the green or the beginning of the yellow area.

Manual and Automatic cutting modes may be used in combination. You may start in manual mode and cut into the workpiece.

Secure the locking handle in this position, press AUTO and continue in automatic mode; the workpiece will move towards the cut-off wheel.

Programmable Stop

Discotom-50 is supplied with an automatic feature, allowing you to pre-program a specific stop position.

- Position the cutting table at the desired stop position using the FORWARDS ▲ and BACKWARDS ▼ keys.
- Set the stop position by pressing the STOP ⊕ key and holding it down whilst activating the FORWARDS ▲ key.
- Use BACKWARDS ▼ to move the cutting table clear of the cut-off wheel. Place and clamp the workpiece and press START ◇.
- The cutting table stops when it reaches the pre-set stop position and the cutting table moves to the position from where cutting was started.

Cancelling the Programmable Stop

- Cancel the stop position by pressing the STOP ⊕ key and holding it down whilst activating the BACKWARDS ▼ key.

Note

The stop position will also be cancelled when the machine is switched off.

Variable Rotational Speed

Discotom-50 is fitted with a variable rotational speed feature. This allows you to optimise cutting quality and cut many different materials using the same cut-off wheel. Thus materials as different as aluminium and hardened steel can be cut using the same cut-off wheel.

Recommended Settings, Rotational Speed

The table below is a guide to recommended rotational speed settings for popular materials using a single cut-off wheel (Struers 53UNI).

What happens when the rotational speed changes?
 Ordinary cut-off wheels, based on bakelite-bonded siliciumcarbide (SiC) or aluminiumoxide (Al₂O₃), change their properties with the rotational speed. The higher the speed, the 'harder' the properties of the wheel, and vice versa. The benefit of a higher speed is that the wheel wears at a slower rate. At the same time, the workpiece becomes more prone to heat damage, which is countered by reducing the feed speed. Using a slower rotational speed, the same wheel becomes softer, and capable of cutting through harder materials, with less risk of workpiece damage. The downside is increased wheel wear. Struers have optimised both rotational and feed speed, as per the table and graph shown.

Rotational Speed Settings using 53UNI cut-off wheel			
Metalog Guide Description	Material Tested	Hardness (HV)	Speed (RPM)
Non ferrous soft metal	Aluminium	50-110	3000
Very ductile metal	Stainless steel	220	2200
Medium soft ferrous metal	UHB IMPAX Cold work tool steel	300	2200
Medium hard metal	UHB IMPAX Cold work tool steel	480	2000
Hard ferrous metal	UHB ARNE Cold work tool steel	750	1750


Please note that the above settings are based on lab tests of the specific materials listed above. Any deviation in material will require a different setting. Interpolate the Vickers hardness to approximate the corresponding rotational speed.

Feed Speed
Recommended Settings,
Feed Speed

The corresponding feed speed setting depends on the sample size; please refer to the graph below.

To use the graph, find the maximum diameter of your workpiece on the x-axis, then move upwards to the curve relevant to your material. Then move left to the y-axis to the corresponding feed speed.

53UNI

Cut-off wheel recommendations

53UNI is recommended for all-round cutting of materials <800 HV.

51UNI is recommended for all-round cutting of materials >800 HV.

Note: Some materials cannot be cut using the 53UNI universal cut-off wheel; for example Titanium.

For extremely hard materials (>1000HV), manual cutting is recommended.

Important


Always adapt **both** rotational speed **and** feed speed to a different type for material.

Material-specific Cut-off Wheels

For most common materials, cutting with a single universal cut-off wheel is possible. If demands on cutting quality are particularly high we recommend the use of cut-off wheels designated for use with the specific material, eg 56UNI for aluminium (Please refer to the Selection Guide in the [Struers Cut-off Wheels Brochure](#)).

When using material-specific cut-off wheels the rotational speed settings remain the same. However the recommended feed speed settings change slightly, please refer to the graph below.

Each curve represents a series of cuts with a material-specific cut-off wheel (name of Struers cut-off wheel in brackets).

*Using CBN or Diamond
Cut-off wheels*

As previously mentioned in the information box on rotational speed, only bakelite-bonded cut-off wheels can be optimised by altering the rotational speed.

When CBN or diamonds are used as an abrasive, the rotational speed must be fixed. Compared to the SiC and Al₂O₃ wheels, CBN wheels should run at a significantly higher speed, whereas diamond wheels perform best at a lower speed.

Wheel Speed for Diamond Cut-off Wheels

Optimal wheel speed is 18 - 24 m/sec.

Set the rotational speed to n rpm, if possible and within safety limits.

$$n = (420,000 / D) \text{ rpm.}$$

D = the diameter of the cut-off wheel in mm.

Wheel Speed for CBN Cut-off Wheels

Optimal wheel speed is 25 - 35 m/sec.

Set the rotational speed to n rpm, if possible and within safety limits.

$$n = (600,000 / D) \text{ rpm.}$$

D = the diameter of the cut-off wheel in mm

Feed Speed

The feed speed may be adjusted within a limited range. Note that the feed speed depends on placement of the wheel as well as sample placement and size.

Automatic Feed Reduction


If the feed speed is set too high, the yellow FEED SPEED REDUCED indicator Δ (situated above the feed speed control) flashes and the feed speed is automatically reduced. Discotom-50 can repeat this reduction up to 5 times.

Once the workpiece has been cut, the operator can reduce the feed speed setting using the feed speed control until the FEED SPEED REDUCED indicator stops flashing. Discotom-50 will then be able to cut another similar workpiece without feed speed reduction.

Cleaning

Flushing the Cutting Chamber

To ensure a longer lifetime for your Discotom-50, Struers strongly recommends that you clean the cutting chamber with the flush hose every day.

- Press STOP ∇ to stop the cut-off wheel and cooling water.
- Open the protection guard and release the workpiece from the clamping device
- Remove clamping device if carrying out major cleaning.
- Pull out the flush hose from the holder in the front and point it towards the bottom of the cutting chamber.
- Turn the cooling valve to cleaning position (re drawing) and adjust the flow as needed.
- Clean the cutting chamber thoroughly.
- Turn the cooling valve back to cutting position.

AVOID CORROSION!

After cutting, leave the protection guard open to let the cutting chamber dry completely.

Reference Guide

Table of Contents	Page
1. Advanced Operations	
Clamping Irregular Workpieces.....	22
Mounting an Adjustable Cutting Table.....	22
Cutting Protruding Workpieces	22
Connection to an External Exhaust System	22
Optimising Cutting Results	23
Explanation of Safety Factors.....	23
2. Accessories.....	24
3. Consumables	
Cut-off Wheels.....	25
Other Consumables	25
4. Trouble-Shooting	26
5. Maintenance	
Daily Service	29
Weekly Service.....	29
Monthly Service	29
Replacing the Cooling Water	29
Maintenance of Cutting Table	30
Maintenance of Cut-off Wheels.....	30
Storing Bakelite Bonded Al ₂ O ₃ Cut-off Wheels	30
Maintenance of Diamond and CBN Cut-off Wheels	30
Maintenance of Clamping Devices.....	30
6. Technical Data	31

1. Advanced Operations

- Clamping Irregular Workpieces** Irregular workpieces without plane clamping surfaces must be clamped using special clamping tools, as the workpieces must not move during cutting, thereby damaging the cut-off wheel or the workpiece itself. For this purpose a vertical clamping device is available. Please refer to the section on Accessories.
- Mounting an Adjustable Cutting Table** If you need to cut plane parallel sections, an adjustable table may be mounted on top of the existing cutting table. The adjustable table positions the workpiece at a right angle to the cut-off wheel and has a sideways movement of 60mm. Please refer to the section on Accessories.
- Cutting Protruding Workpieces** Discotom-50 comes as standard with a rubber curtain and a hinged protection guard, allowing protruding workpieces to be cut.
- Connection to an External Exhaust System** Struers recommends the use of an exhaust as workpieces may emit harmful gasses or dust when cut. On the back of Discotom-50 you will find a joint for a 50 mm dia. exhaust hose. Mount an exhaust hose from your local exhaust system onto the joint.

Optimising Cutting Results

The following table offers answers to a number of common questions about how to achieve better quality cuts:

Question	Answer	Comments
How can I avoid discoloration or burning of the sample?	Reduce rotational speed.	Causes increased wheel wear.
	If reduced rotational speed fails to solve problem, change to softer cut-off wheel*	
How can I avoid unplane cut?	Reduce rotational speed.	Causes increased wheel wear.
How can I avoid smearing?	Reduce rotational speed.	Causes increased wheel wear.
How can I avoid burrs?	Use a softer cut-off wheel*	Causes increased wheel wear
	Secure the workpiece on both sides of the cut-off wheel	
How can I prevent the cut-off wheel wearing out too quickly?	Increase rotational speed.	May cause sample discoloration and unplane cut.
	Use a harder cut-off wheel*	
How can I cut faster?	Place the cut-off wheel as low as possible.	
	Place the workpiece in a position that allows the cut-off wheel to cut the smallest possible cross-section.	
	Increase feed speed.	May cause sample discoloration and unplane cut.
How can I avoid machine vibrations?	Minor vibrations: Increase feed speed in steps of 0.1 m/s.	May cause sample discoloration and unplane cut.
	Major vibrations: Increase rotational speed by 500 rpm.	May cause sample discoloration and unplane cut.

* Please refer to the Selection Guide in the Struers [Cut-off Wheels Brochure](#).

Explanation of Safety Factors

The protection guard has a safety switch to prevent the cut-off wheel from starting while the guard is open. Furthermore, a locking mechanism prevents the operator from opening the protection guard until the cut-off wheel stops spinning.

2. Accessories

Specification	Code
<p>Adjustable cutting table For cutting plane parallel sections. It positions the workpiece at a right angle to the cut-off wheel. Moveable up to 60 mm. With 10 mm T-slots and exchangeable steel bands.</p>	TRETA
<p>Quick clamping device for 10 mm T-slots, Left For securing the workpiece. Complete with backstop.</p>	TRELQ
<p>Spring clamp For retaining the cut-off specimen. Complete with backstop.</p>	TRESC
<p>Compact quick clamping device for 10 mm T-slots, Right For securing the workpiece. Occupies only 45 mm of the cutting table. To be used together with backstop from TRESC.</p>	TRERC
<p>Compact quick clamping device for 10 mm T-slots, Left For securing the workpiece. Occupies only 45 mm of the cutting table. To be used together with backstop from TRELQ.</p>	TRELC
<p>Vertical clamping system for 10 mm T-slots For clamping irregular shaped workpieces. Complete with operating key and one flat clamping shoe</p>	TREVS
<p>Riser block for 10 mm vertical clamping system For elevating Vertical Clamping System by 60 mm</p>	TREKS
<p>Swivel shoes for vertical clamping system Set of four multi-shaped swivel shoes for TREVS</p>	TREVI
<p>Recirculation Cooling Unit, 65 l For water cooling on Discotom-5/-50. With a rolling pallet</p>	TRECA
<p>Replacement stainless steel bands Set of 2 pcs. 225 x 40 mm steel bands</p>	TREUV
<p>Table Unit With compartment for recirculation cooling unit. Width: 900 mm, depth: 750 mm, height: 800 mm</p>	LABUL

3. Consumables

Cut-off Wheels

Please refer to the Selection Guide in the Struers [Cut-off Wheels Brochure](#).

Other Consumables

Specification	Code
<i>Additive for Cooling Fluid</i> 1 l 5 l	ADDUN ADDFI
<i>Bacterial Cleaning Liquid</i> 0.5 l	BACTE
<i>Disposable Plastic Insert</i> for collecting waste material in the Recirculation Cooling Unit. Package with 10 pcs.	EXOSP

4. Trouble-Shooting

Error	Explanation	Action
Indications on the Front Panel		
Red STOP  indicator on.	The emergency stop has been pressed.	Release the emergency stop to shut off the STOP  indicator.
Red STOP  indicator on and beep signal.	The emergency stop has been activated and you tried to press another key.	
Yellow FEED SPEED REDUCED  and Red STOP  indicator lights come on.	Cutting motor overloaded due to heavy use: Material is too hard Feed speed setting too high Rotational speed setting too high Cutting table is jammed	Press any key to stop the indicators. Solve the problem causing the machine to stop, for example, by reducing the feed speed, or using a different cut-off wheel.
Machine Problems		
Machine has stopped cutting.	The AutoStop function is active.	Switch unit off and on at Main Switch to reset AutoStop.
During cutting the cutting table stops when reaching a specific point.	The AutoStop function is active.	Switch unit off and on at Main Switch to reset AutoStop.
The cut-off wheel rotates and coolant runs, but cutting table does not move.	Cutting table is stationary because AUTO mode has been de-selected.	Cut manually by moving cutting handle or press AUTO to select automatic cutting.
The cut-off wheel does not rotate.	The protection guard may be open.	Press the safety guard completely down. If this does not help, call a Struers service technician.
Cut-off wheel rotates, but no cooling water.	The thermal overload protection of the recirculation pump has been activated.	Press red reset button on back of unit.
Chamber light does not work.	Replace the lamp.	Remove plastic cap to gain access to the fluorescent lamp. Pull out the lamp and replace it.
Water leaking.	Leak in a Recirculation Cooling Unit hose.	Check the hose and tighten the hose clamp.
	Water overflow in the cooling water tank.	Remove the excess water in the tank.
Workpieces or cutting chamber rusty.	Insufficient additive in cooling fluid.	Use Struers Additive for cooling fluid in the cooling water, at the correct concentration. Check with a refractometer. Follow the instructions in the 'Maintenance' Section.
	The protection guard is left closed after use.	Leave the protection guard open to let the cutting chamber dry.

Discotom-50
Instruction Manual

Error	Explanation	Action
Cutting Problems		
Discoloration or burning of the workpiece.	The hardness of the cut-off wheel is inappropriate for the hardness / dimensions of the workpiece.	Reduce rotational speed or change wheel. See Consumables section, under heading 'Cut-off Wheels'.
	Inadequate cooling.	Check that there is sufficient water in the Recirculation Cooling Unit.
Unwanted burrs.	Cut-off wheel too hard.	Reduce rotational speed or change wheel. See Consumables section, 'Cut-off Wheels'.
	Feed speed too high at the end of the operation.	Reduce the feed speed near the end of the operation.
	Lack of support for the workpiece.	If possible, support the workpiece on both sides.
Cutting quality differs.	Cooling water hose clogged.	Clean the cooling water hose and the cooling tube. Check the water flow by turning the cooling valve to cleaning position.
Cut-off wheel breaks.	Incorrect mounting of the cut-off wheel.	Check that the bore/centre hole has the correct diameter. Check cardboard washer on both sides of the cut-off wheels and replace if worn. The nut must be tightened properly.
	Incorrect clamping of the workpiece.	Make sure that only the left quick-clamping device is tight. The spring clamp to the right should only press lightly. Use the vertical clamping system if the workpiece is an irregular shape.
	Cut-off wheel is too hard.	Reduce rotational speed or change wheel. See Consumables section, 'Cut-off Wheels'.
	Feed speed is set too high.	Reduce the feed speed.
	Inadequate cooling.	Check that there is enough water in the cooling unit.

Discotom-50
Instruction Manual

Error	Explanation	Action
The cut-off wheel wears down too quickly.	The rotational speed is too low.	Increase rotational speed.
	The feed speed is too high.	Reduce feed speed.
	Insufficient cooling.	Check that there is enough water in the Recirculation Cooling Unit. Check the cooling water hoses.
The cut-off wheel does not cut through the workpiece.	The rotational speed is too low.	Increase rotational speed.
	The cut-off wheel is too soft for the task	See Consumables section, 'Cut-off Wheels'.
	Incorrect choice of cut-off wheel.	See Consumables section, 'Cut-off Wheels'.
	Cut-off wheel worn.	Replace the cut-off wheel.
The workpiece breaks when clamped.	The cut-off wheel gets caught in the workpiece.	Support the workpiece and clamp it on both sides of the cut-off wheel so that the cut stays open.
	The workpiece is brittle.	Place the workpiece between two polystyrene plates. Note! Always cut brittle workpieces very carefully.
The sample is corroded	The sample has been left in the cutting chamber for too long.	Remove the sample directly after cutting. Leave the protection guard open when you leave the machine.
	Insufficient additive for cooling fluid.	Use Struers Additive for cooling fluid in the cooling water at the correct concentration. Check with a refractometer. See Maintenance section.

5. Maintenance

Daily Service

- Clean the cutting chamber, especially the cutting table with the T-slots.
- Do **not** clean the lamp glass with alcohol. Use a damp cloth.
- Empty and clean the sieve in the Recirculation Cooling Unit.
- Refill the Recirculation Cooling Unit.

IMPORTANT

Always use the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to top up with Struers Additive each time you refill with water.

Weekly Service

- Remove clamping device
- Clean the cutting chamber thoroughly.

Checking the Recirculation Cooling Unit

The cooling unit should be checked for cooling water after 8 hours use or at least every week. The unit must be refilled if the flushing pump cannot reach the cooling water. Remember to add Struers Additive for Cooling Fluid: One part of Additive for 33 parts of water.

To check the concentration of additive, use a refractometer. Concentration = Brix value. The concentration of additive should always lie between 2.7 and 3.3 %. Add Struers Additive for Cooling Fluid if the concentration is too low.

Monthly Service

Replacing the Cooling Water

- Replace the cooling water in the Recirculation Cooling Unit (optional TRECA) at least once a month.
- Place a disposable plastic insert (EXOSP) in the tank and fold it over the rim.
- **Please Note:** Too high a level of coolant in the tank might damage the pump. To avoid this, place the disposable insert so that the pre-punched hole is in front of the overflow aperture in the tank.
- Fill the tank with 63 l of water and 2 l Struers Additive. The water level should be 8-10 cm below the upper rim of the tank.

Maintenance of Cutting Table

The stainless steel bands are available as spare parts and should be replaced if they become worn or damaged.

Maintenance of Cut-off Wheels

Storing Bakelite Bonded Al₂O₃ Cut-off Wheels

Bakelite bonded cut-off wheels are sensitive to humidity. Therefore, do not mix new, dry cut-off wheels with used humid ones. Store the cut-off wheels in a dry place, horizontally on a plane support.

Maintenance of Diamond and CBN Cut-off Wheels

The precision of diamond and CBN cut-off wheels (and thus the quality of the cut) depends on how carefully the following instructions are observed:

- Never expose the cut-off wheel to a heavy mechanical load, or heat.
- Store the cut-off wheel in a dry place, horizontally on a plane support, preferably under light pressure.
- A clean and dry cut-off wheel does not corrode. Therefore, clean and dry the cut-off wheel before storing. If possible, use ordinary detergents for the cleaning.
- Regular dressing of the cut-off wheel is also part of the general maintenance.

Maintenance of Clamping Devices

Important

It is recommended to thoroughly clean and lubricate the Quick Clamping Device and Vertical Clamping Device at regular intervals.

6. Technical Data

Subject		Specification	
		Metric/International	US
CUTTING SPECIFICATIONS			
Workpiece Dimensions (Max)	<i>Workpiece inside cutting chamber:</i>		
	Height	125 mm	4.92"
	Width	220 mm	8.66"
	Depth	195 mm	7.67"
	<i>Workpiece protruding cutting chamber:</i>		
	Height	125 mm	4.92"
	Depth	400 mm	15.75"
Cutting Capacity (Max) <i>Please refer to Cutting Diagram</i>	Max. workpiece diameter	105 mm	4.13"
	<i>Max. size of cut (example)</i>		
	Height	90 mm	3.54"
	Depth	140 mm	5.51"
PHYSICAL SPECIFICATIONS			
Cutting Motor	Cutting power	4 kW	5.4 HP
Cut-off Wheel	Diameter x Thickness x Centre-hole	300 x 2 x 32 mm	11.8 x 0.078 x 1.26"
	Rotational speed (running idle)	1000 – 3000 rpm	1000 – 3000 rpm
Positioning and Feed	Positioning range (of cut-off wheel)	0 – 120 mm	0 – 4.72"
	Positioning range (of cutting table)	0 – 195 mm	0 – 7.67"
	Max. positioning speed	15 mm/s	0.59"/s
	Feed Speed range	0.1 – 2.5 mm/s	0.04 – 0.1"/s
Cutting Table	Width	336 mm	13.23"
	Depth	225 mm	8.86"
	T-slots	10 mm	0.39"
Dimensions and Weight	Height (closed/open hood)	570/890 mm	22.44/35.04"
	Width	810 mm	31.89"
	Depth (closed/open hood)	785/816 mm	30.9/32.12"
	Weight	125 kg	276 lbs
	<i>Optional table unit (LABUL)</i>		
	Height	800 mm	31.5"
	Width	900 mm	35.4"
Depth	750 mm	29.5"	
Recirculation Cooling Unit	Tank volume (optional TRECA)	65 l	17 gallons
	Approx. flow	14 l/min	3.7 gallons/min

Discotom-50
Instruction Manual

Subject		Specification
Electrical Data	<i>Voltage / frequency:</i> 3 x 200 V / 50 Hz 3 x 220-230 V / 50 Hz 3 x 380-415 V / 50 Hz 3 x 200-210 V / 60 Hz 3 x 220-240 V / 60 Hz 3 x 380-415 V / 60 Hz 3 x 460-480 V / 60 Hz	<i>Max. Load</i> 19 A 19 A 11 A 18.6 A 18.6 A 10.8 A 10.8 A
Environment	Directives Noise level	Machinery Directive (MD) 98/37/EEC, 99/79/EEC EMC Directive (EMC) 89/336/EEC, 92/31/EEC Low Voltage Directory (LVD) 72/23/EEC Approx. 67dB(A) running idle, at a distance of 1.0 m / 39.4" from the machine.

Quick Reference Card

Clamping the Workpiece

- Secure the workpiece using the clamping device of your choice. The most common is the quick-clamping device. Place the workpiece between the clamp and the back stop.
- Push the clamp towards the workpiece and lock the quick-clamping device with the locking handle.

Automatic Cutting

- Clamp the workpiece.
- Position the cut-off wheel so that it can cut the sample and lock it in place with the black cut-off wheel locking handle on the right of the machine.
- Position the cutting table with FORWARDS ▲ and BACKWARDS ▼ keys.
- Carefully close the protection guard.
- Set feed speed and Press AUTO.
- Press START ◇ to begin cutting.
- After cutting, the cutting table automatically returns to its start position.
- Open the protection guard and release the clamping device.
- Remove the workpiece.

Manual cutting

- Loosen the locking handle on the right hand side of the machine, so that the cutting handle moves freely.
- Clamp the workpiece.
- Position the cutting table under the cut-off wheel.
- If AUTO mode is selected, de-select it.
- Press START ◇.
- Lower the cut-off wheel until it contacts the workpiece and apply pressure until the sample is cut.
- Keep the Load Bar in the green, or first part of the yellow area. Relax pressure near the end of the cut.
- Open the protection guard and release the clamping device.
- Remove the workpiece.

Changing the Cut-off Wheel

- Fasten the black handle to stop the cutting handle from moving.
- Press the black locking knob on the right-hand side of the cut-off wheel, turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner. Remove the flange and the old cut-off wheel.
- Mount the new cut-off wheel.
- Mount the flange and the nut. Tighten carefully.

Flushing the Cutting Chamber

- Pull out the flush hose from the front of the machine.
- Point the flush hose towards the bottom of the cutting chamber.
- Turn the cooling valve clockwise to flushing position.
- Clean the cutting chamber.
- Turn the cooling valve counter-clockwise to cooling position.
- Always leave the protection guard open to let the cutting chamber dry.

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Date of Release 26.07.2006

Discotom-50
Spare Parts and Diagrams

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

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Some of the drawings may contain position numbers
not used in connection with this manual.

Discotom-50
Spare Parts and Diagrams

Spare part list for Discotom-50

Drawing	Pos.	Spare Part	Cat no:	
15130001		Discotom-50, complete		
	50	PVC Hose transparent 13/32" (2 m)	2NU19310	
15130010		Cutting screen		
		Cutting screen complete	15130010	
	20	Seal for front of cover	2IP21038	
	90	Rubber curtain with slits	15130158	
	110	Fastening plate	15130150	
	140	Screws M4x12 A2	2TR80408	
	160	Lamp glass, complete O-ring (for lamp glass) Lamp 9W	15110080 2IO20060 2HG10827	
15130020		Front plate complete		
	140	Button	2GD14425	
	150	Cover for button	2GD54025	
	160	Arrow	2GD64025	
15130005		Back-and bottom plate, complete		
	110	Buffer $\varnothing 15 \times 7 / M5 \times 5$, Type D	2GS00090	
	130	Handle for flushing hose	15110495	
	140	Hose clamp	2NS11709	
	150	Hose for flushing inlet (Please specify length)	2NU21218	
	160	Flex cable	2NU3114	
	170	Cable holder	2GK31116	
	355	0.25 AT fuse, glass type Knob for flushing valve	2FU13500 15130065	
	15110045		Bottom plate, complete	
		600	Fuse holder	2XS12004
605		Fuse 4.0 AT	2FU14200	
660		Cardboard blotters $\varnothing 33 / 85 F . x x T R E$	50000028	
670		Flange	14150030	
680		Nut, M16 Left-hand thread	2TA70160	
700		Washer for spindle $\varnothing 40 / \varnothing 16 . 5 \times 5$	14150012	
730		Flushing head, complete O-ring for flushing head	15130090 2IO00005	
780		Hose for flushing head (1 m.)	2NU19310	

Discotom-50
Spare Parts and Diagrams

Spare part list for Discotom-50

Drawing	Pos.	Spare Part	Cat no:
15130030		Y-table, complete	
		Replacement stainless steel bands (2 pcs.)	05136904
	100	Stainless steel band (40 mm)	
	110	Screws M8x20 A4, low head (6 pcs)	
	130	Guide pins 6m6x16 (4 pcs)	
15110075		Quick clamping device. Left	
		Quick clamping device. Left. Complete	15110075
		Knob for handle	13810088
15130078		Quick clamping device. Right	
		Quick clamping device. Right. Complete	15130078
		Knob for handle	15110984
14980068		Recirculation tank	
	90	Cooling Tank, 65 l	381MP140
	100	Lid	381MP142
	110	Clamp for sieve	2NS20215
	110	Sieve holder	14980073
	110	Sieve bag	315MP053
	110	Sieve bag for tank, complete (Incl. clamp and holder)	315MP052
15110035		Adjustable cutting table	
		X-table, complete (Accessory)	
		Replacement stainless steel bands. 1x20 + 1x40	05116906
	240	Stainless steel band (40 mm)	
	250	Stainless steel band (20 mm)	
	260	Screws M8x20 A4 low head (6 pcs)	
	140	Guide pins 6m6x16 (4 pcs)	
		Replacement stainless steel bands, 3x20	05116908
	250	Stainless steel band, (20 mm, 3 ps)	
	260	Screws M8x20 A4, low head (9 pcs)	
140	Guide pins 6m6x16 (6 pcs)		
	T-Bolt M8 for 10 mm T-slots	15112903	
210	Screw M8x20 A4, low head		
220	Nut M8 for 10 mm T-slots		
15130080		Three-jaw chuck (Accessory)	
	160	Hand wheel	2GH00212
	50	Jaw incl. flange	2GR20801
14980076		Recirculation pump	
		Rec. pump, complete. 3x200-240V/50Hz	381MP408A

Discotom-50
Spare Parts and Diagrams

Spare part list for Discotom-50


Drawing


Pos.	Spare Part	Cat no:
	Rec. pump, complete. 3x380-440V/50Hz	381MP408B
	Rec. pump, complete. 3x200-220V/60Hz	381MP409A
	Rec. pump, complete. 3x380-480V/60Hz	381MP409B
	Rec. pump, complete. 3x220-480V/60Hz CSA	13810958
20	Sieve ø99 for recirculation pump	14980074
30	Clamp for sieve ø99	2NS28012
Various		
	Socket spanner (13 mm)	2GR00313
	Fork spanner (24 mm)	2GR00224
	Filter for drain in cutting chamber	15130583

Vertical clamping system for 10 mm T-slots
TREVS

1	Operating key for TREVS	2GR27021
2	Spring for TREVS	2GR27023
3	Flat clamping shoe	2GR26150

**Vertical clamping system for
 10 mm T-slots (TREVS)**

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

FT1: Thermal Overload Relay: SETTING RANGES

MANNS SUPPLY:	FT1: SET =
3x200 - 230V / 50Hz	0.60 Amp.
3x200 - 240V / 50Hz	0.53 Amp.
3x380V / 60Hz	0.30 Amp.
3x380 - 415V / 50Hz	0.33 Amp.
3x460 - 480V / 60Hz	0.30 Amp.

CUTTING MOTOR.

RECIRCULATION PUMP MOTOR

"ACCESSORIES" SAMPLE MOVER

SAMPLE MOVER SIGNAL START / STOP.

EMERGENCY STOP OFF AND SAFETY GUARD CLOSED.

RECIRCULATING PUMP MOTOR NOT IN OVERLOAD

INVERTER FOR CUTTING MOTOR NOT READY.

RECIRCULATING PUMP STOPPED.

SAFETY GUARD RELEASE

STEPMOTOR, MOVING TABLE

CUTTING CURRENT:






GREEN	Cut. Cur. 12.5 %
GREEN	Cut. Cur. 25%
GREEN	Cut. Cur. 37.5%
GREEN	Cut. Cur. 50%
GREEN	Cut. Cur. 62.5%
GREEN	Cut. Cur. 75%
GREEN	Cut. Cur. 87.5%
GREEN	Cut. Cur. 100%
YELLOW	Cut. Cur. 125%
YELLOW	Cut. Cur. 150%
RED	Cut. Cur. 175%
RED	Cut. Cur. 200%

AUTO FEED ON:

EMERGENCY OR OVERLOAD:

FEED INFORM:

KEYBOARD BUZZER:

COLOR CODES: BK = Black
 BN = Brown
 RD = Red
 OG = Orange
 YE = Yellow
 GN = Green
 BU = Blue
 VT = Violet
 GY = Grey
 WH = White
 PK = Pink

ALL WIRINGS 0.4 TO 0.75mm², EXCEPT OTHERWISE MARKED.

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












W13 TO BQ1
W14 TO BQ2
SEE PAGE 1

*Shield to Front Plate

W9 TO X4, SEE PAGE 1.
W10 TO X4, SEE PAGE 1.
W6 TO X2, K1, K2, FT1 SEE PG. 2.
W20 TO MMV400

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Herewith declares that

<i>Product Name:</i>	Discotom-50
<i>Type No:</i>	513
<i>Machine Type:</i>	Cut-off machine

is in conformity with the provisions of the following directives:

Safety of Machinery 98/37/EEC according to the following standard(s):
EN ISO 12100-1:2003, EN ISO 12100-2:2003, EN294:1992, EN349:1994, EN954-1:1997,
EN1037:1995, EN1050:1996, EN60204-1:1997.

EMC-Directive 2004/108/EEC according to the following standard(s):
EN61000-6-1:2002, EN61000-6-3:2002.

Low Voltage Directive 2006/95/EEC according to the following standard(s):
EN60204-1:1997.

Supplementary Information The equipment complies with the American standards:
FCC part 15, subpart B and UL508, NFPA 79, JIC EGP-1.

The above has been declared according to the global method, module A

Date: 04.12.2007

Christian Skjold Heyde,
Vice President, R&D and Production, Struers A/S

Hersteller Struers A/S
Pederstrupvej 84
DK-2750 Ballerup, Danmark
Telefon +45 44 600 800

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<i>Produktname:</i>	Discotom-50
<i>Typennr.:</i>	513
<i>Maschinenart:</i>	Trennmaschine

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Discotom-50

 **Struers**

**Combined automatic
and manual cut-off
machine for 300 mm
cut-off wheels**

The moveable cutting table permits cutting of large, irregular samples

Total flexibility in cutting regardless of material type

Use only one cut-off wheel for a wide range of materials

Optimise cutting with CBN and diamond cut-off wheels

Highly reproducible results

Flexible and user-friendly

Ideal for large and deep samples

Optional: Three-jaw chuck for cutting virtually all materials

Discotom-50 successfully combines the newest cutting technologies with unprecedented flexibility and user-friendliness. It drastically widens the range of your cutting applications to very large and hard materials without additional time or cost. Additionally Struers has done extensive research to unveil the optimum combination of spindle speed and feed speed for a range of popular materials thereby bringing you to a much higher level of control and flexibility.

Deal with a wide range of materials using only one cut-off wheel

The Discotom-50 spindle speed ranges from 1000 - 3000 rpm and makes it possible for customers dealing with many different materials to use one single wheel. By adjusting spindle speed, Discotom-50 is able to cut through a wide range of materials as diverse as aluminium and hardened steel to give you full flexibility in cutting. No time is wasted changing wheels and no money is wasted by storing many different wheel types.

Using Struers' three-jaw chuck, Discotom-50 is capable of dealing with even the most difficult materials (>800 HV). The sample is constantly rotated during cutting in order to reduce the contact area between the cut-off wheel and the sample, ensuring optimum cooling at the same time.

Optimise cutting quality with material-specific cut-off wheels

For the highest quality demands and cutting of difficult or very hard materials, it is recommended to use material specific cut-off wheel and optimise feed speed accordingly. High quality cutting allows shorter grinding and polishing steps reduce your total preparation times.

Optimise cutting with CBN and diamond cut-off wheels

Use of CBN and/or diamond cut-off wheels requires fixed spindle speed. Compared to SiC and Al₂O₃ wheels, CBN and diamond wheels deliver their best at a lower speed, which is made possible on Discotom-50. The spindle speed may be adjusted within a range of 1000 – 3000 rpm.

Highly flexible and reproducible operations

If the feed speed is set too high and the machine cannot cut through the workpiece, Discotom-50's Electronic Feed Speed Reduction automatically reduces the feed speed and continues cutting at a lower speed. By automatically reducing feed speed, Discotom-50 prevents samples from burning and cut-off wheels from being damaged. This may be repeated up to 5 times.

Using the set-and-store function, the optimum feed speed is stored and can be used for successive cutting on similar workpieces making your cutting jobs easy and time saving.

Adjusting the speed of the moveable cutting table can be done prior to and during automatic cutting.

Three-jaw chuck for cutting of very hard materials (Option)


Programmable stop

Discotom-50's programmable stop ends cutting at a specific, pre-programmed position and returns to the start position after having completed the cut. This feature becomes particularly time-saving when cutting many similar, smaller samples. No time is wasted by

Rotational Speed Settings using 53UNI Cut-off wheel Discotom-50

Metalog Guide Description	Tested Material	Hardness HV	Speed RPM
Non ferrous soft metals	Aluminium	50-110	3000
Very ductile metals	Stainless steel	220	2200
Medium soft ferrous metals	UHB IMPAX Cold work tool steel	300	2200
Medium hard ferrous metals	UHB IMPAX Cold work tool steel	480	2000
Hard ferrous metals	UHB ARNE Cold work tool steel	750	1750

Feed Speed Settings using 53UNI Cut-off wheel Discotom-50

- Non ferrous soft metals
- Medium soft ferrous metals
- Medium hard ferrous metals
- - - Hard ferrous metals
- Very ductile metals

Moveable cutting table

Front-mounted flush hose and cooling valve

Adjustable cutting table for cutting many similar specimens (optional)

waiting for the table to move all the way to the back. Likewise, the fast-forward function allows a fast positioning of the cutting table where it is needed when cutting new samples.


User-friendly and multi-functional

Discotom-50 is easy to operate and allows quick positioning and easy cutting. Combining manual positioning of the sample with automatic cutting, ensures high sample throughput.

A wide range of cutting applications are available by choosing either automatic or manual cutting with Discotom-50's Automatic On /Off function. From quality demanding to high throughput cutting jobs, from very soft to very hard materials or from small to large samples, Discotom-50 grants extreme flexibility, can be operated quickly, easily and at a minimum cost.

Cut large and deep samples

Discotom-50 has a large cutting chamber and a movable cutting table. By moving the large table instead of the wheel, Discotom-50 cuts larger and deeper samples, up to 105 mm (4.1") in diameter or 90 x 140 mm (3.5 x 5.5") square. The cutting capacity is further increased to 125 mm (4.9") diameter when using Struers' optional three-jaw chuck.

Reliable and user-friendly design

The cutting motor is located outside the cutting chamber to prevent any contact with the coolant. This ensures maximum reliability of the motor and provides additional space inside the cutting chamber.

Discotom-50's Load meter on the cutting motor is accurately monitored on the 12-level LED-bar on the control panel. The bar is divided into three areas, a green area showing that the load is ok, a yellow area showing levels close to overload and the red area indicating that the motor is overloaded.

Cutting certain materials may produce both harmful and unpleasant fumes. In line with our general recommendation to install an exhaust from the cutting chamber, the machine is prepared for direct connection to an external exhaust system. The flush hose and the cooling

valve are placed on the front for maximum convenience.

High safety standards

Discotom-50 complies with international safety standards and its safety features include an emergency stop and cut-off switches for all machine functions and motor, ensuring a high level of safety for the operator.

The cutting process cannot be started until the impact resistant protection guard has been closed completely. After the start button has been pressed, the guard is locked instantly. Should the guard be forced open during the cutting process, the main motor will be switched off immediately.

Specification

Discotom-50

Automatic cut-off machine for 300 mm cut-off wheels. With variable spindle speed and variable, automatic feed. 4 kW (5.4 HP) motor. Recirculation cooling unit (TRECA) and clamping tools are ordered separately

Accessories

Three-jaw chuck

Adjustable cutting table with motorised three-jaw chuck. For extra large and difficult workpieces up to 125 mm (4.9") diameter

Adjustable cutting table

100 x 225 mm (3.9" x 8.8"). Supplied for cutting plane parallel sections. It positions the workpiece at a right angle to the cut-off wheel. Moveable up to 60 mm (2.4"). With 10 mm T-slots and exchangeable steel band. 2 pcs

Quick-clamping Device for 10 mm T-slots, Left

For securing the workpiece. Complete with back stop

Spring clamp

For securing the workpiece. Complete with back stop

Vertical Clamping System for 10 mm T-slots

For clamping irregularly shaped workpieces. Complete with operating key and one flat clamping shoe

Code

TREVL

TRECH

TRETA

TRELQ

TRESC

TREVS

Riser Block for 10 mm T-slots

For elevating Vertical Clamping System by 60 mm.

Code

TREKS

Swivel Shoes for Vertical Clamping System

Set of 4 multi-shaped swivel shoes.

TREVI

Recirculation Cooling Unit

For water cooling on Discotom-5/-50. Capacity 65 l (17 US Gal.). With a rolling pallet.

TRECA


Struers A/S

 DK-2750 Ballerup, Denmark
 Phone +45 44 600 800
 Fax +45 44 600 801
 struers@struers.dk

Cutting capacity

Maximum cutting diameter: 105 mm (4.1")

with three jaw chuck option: 125 mm (4.9")

Physical specifications

Cutting Motor	Cutting power	5.5 kW	7.4 HP
Cut-off Wheel	Diameter x Thickness x Center-hole Rotational speed (running idle)	300 x 2 x 32 mm 1000 – 3000 rpm	11.8 x 0.08 x 1.3" 1000 – 3000 rpm
Positioning & Feed	Positioning range (of cut-off wheel) Positioning range (of cutting table) Max. positioning speed Feed Speed range	0 – 120 mm 0 – 195 mm 15 mm/s 0.1 – 2.5 mm/s	0 – 4.7" 0 – 7.7" 0.8"/s 0.04 – 0.1"/s
Cutting table	Width Depth T-slots	336 mm 225 mm 10 mm	13.2" 8.9" 0.4"
Dimensions and Weight	Height (closed/open hood) Width Depth (closed/open hood) Weight Optional table unit (LABUN) Height Width Depth	570/890 mm 810 mm 785/816 mm 125 kg 795 mm 700 mm 750 mm	22.4/35.2" 31.9" 30.9/32.1" 276 lbs 31.3" 27.5" 29.5"
Environment	Noise level	Approx. 67dB(A) running idle, at a distance of 1.0 m / 39.4" from the machine	

Struers' equipment is in conformity with the provisions of the applicable International Directives and their ap-
 purtenant Standards. (Please contact you local supplier for details)

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Struers' products are subject to constant product development. Therefore, we reserve the right to introduce
 changes in our products without notice

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Struers Inc.
Struers Ltd.
Struers (Shanghai) Ltd.