Discotom-60/-65

Instruction Manual

Manual No.: 15907001

Date of Release: 09.09.2013

Table of Contents	Page
User's Guide	1
Reference Guide.....	35
Quick Reference Guide.....	61

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the *Serial No.* and *Voltage* on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manuals may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manuals may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

Original instructions. The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2013.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Discotom-60/ -65

Safety Precaution Sheet

To be read carefully, before use

1. The operator(s) should be fully instructed in the use of the machine and its cut-off wheels according to the Instruction Manual and the instructions for the use of cut-off wheels.
2. The machine must be installed in compliance with local safety regulations.
3. The machine must be placed on a safe and stable support table. All safety functions and guards of the machine must be in working order.
4. Use only intact cut-off wheels. The cut-off wheels must be approved for use with rotational speeds between 1000 and 3460 rpm.
5. The machine is not for use with saw-blade type cut-off wheels.
6. Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
7. The workpiece must be securely fixed in the quick-clamping device or similar. Large or sharp workpieces must be handled in a safe way.
8. Do not work on or around cutting table when the table is repositioned using the Y-table positioning joystick.
9. To achieve maximum safety and lifetime of the machine, use only original Struers consumables.
10. The cutting handle should be lowered slowly and carefully, in order to avoid breaking the cut-off wheel.

Discotom-60 /-65
Instruction Manual

- 11.** Struers recommends the use of an exhaust system as the materials being cut may emit harmful gasses or dust.
- 12.** The machine emits only moderate noise. However, the cutting process itself may emit noise, depending on the nature of the workpiece. In such cases, the use of hearing protection is recommended.
- 13.** The machine must be disconnected from the mains prior to any service.
- 14.** Use of working gloves is recommended as workpieces may be both very hot and produce sharp edges. Wearing of gloves is also recommended when flushing and cleaning the machine.
- 15.** Struers recommends the use of safety shoes when working with heavy samples.
- 16.** Use of safety goggles is recommended when using the flushing hose.
- 17.** If any of the cutting-chamber hood springs are damaged (at the rear of machine), they must be replaced before the machine is used again.
- 18.** When a recirculation cooling unit is used, observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
Do not use flammable cooling fluid.
- 19.** When lifting the machine using a forklift, lift from front or rear - never lift the machine from the side.
- 20.** When lifting the machine using lifting straps, ensure that the straps are crossed and do not press on the sides of the machine.

Discotom-65

- 1.** Prior to any service, disconnect the machine then wait until residual potential on the capacitors is discharged.
- 2.** Do not cycle mains power more than once every three minutes.
Damage to the frequency converter will result.

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.)

User's Guide

Table of Contents	Page
1. Getting Started	3
Checking the Contents of the Packing Box.....	3
Getting Acquainted with Discotom	4
Noise Level	5
Power Supply	6
Connection to an External Exhaust System.....	8
Connecting a Recirculation Cooling Unit.....	9
2. Basic Operation.....	10
Using the Controls.....	10
Front Panel Controls of Discotom	10
Front Panel Controls	11
Cooling Valve Positions	12
Flush Hose	12
Cleaning of Flush Hose Nozzle	12
Setting the Language	13
Moveable Table.....	13
Y-Table	13
Reading the Cutting Display.....	14
Changing Cutting Mode and Cutting Parameters.....	15
Changing Cutting Mode	15
Changing Cutting Parameters.....	16
Using Motor Load and Temperature Display	22
Automatic Feed Speed Control.....	22
Overload Protection	23
Fitting or Changing the Cut-off Wheel	24
Clamping the Workpiece	24
Positioning the Cutting Table	24
Starting/Stopping the Cutting Process	25
Automatic Cutting.....	25
Manual Cutting.....	26
Combining Manual and Automatic Operation	26
Stop Modes	27
Relative Stop Position	27
Absolute Stop Position	28
Automatic Stop	29
Variable Motor RPM (D-65 only)	30
Recommended Settings:.....	30
Cleaning.....	32
Flushing the Cutting Chamber	32

3. Routine Maintenance	33
Daily Service	33
Weekly Maintenance	33
Monthly Maintenance	33
Replacing the Cooling Water	33
Yearly Service	34
Inspection of Guard.....	34

1. Getting Started

Checking the Contents of the Packing Box

In the packing box you should find the following parts:

- 1 Fork spanner (24 mm) for cut-off wheel
- 1 Connector pipe for water outlet
- 3 Elbow pipes for water outlet
- 1 Outlet hose 2 m,
for connection to external cooling unit
- 3 Hose clamps, 35-60 mm
- 1 Grease for maintenance/lubrication of the spindle
- 1 Oil for maintenance of cutting table
- 1 Set of Instruction Manuals

Unpacking Discotom

- Remove the bolts from all of the transport brackets that secure the Discotom to its transport pallet.
- Remove the brackets.

Placing Discotom


Discotom should be placed on a table strong enough to support a minimum of 171 kg/ 377 lbs.

Struers recommends the use of the Table Unit, which is designed for use with Discotom machines, see "[Accessories](#)".

Getting Acquainted with Discotom


Take a moment to familiarise yourself with the location and names of the Discotom's components.

Front View


- ① Cutting handle
- ② Control panel
- ③ Y-table joystick
- ④ Emergency stop
- ⑤ Cooling valve

Side view, right

- ⑥ Main power switch
- ⑦ Cut-off wheel locking handle

Rear View

- ① Cooli unit connection
- ② Service socket (RS232)
- ③ Water inlet
- ④ Drainage hose
- ⑤ Water outlet flange
- ⑥ Flushing hose
- ⑦ Cutting chamber hood springs
- ⑧ Exhaust connection

Noise Level

Approx. 67 dB (A) measured at idle running, at a distance of 1.0 m/39.4" from the machine.

Power Supply

- First check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.

Discotom-60

- Discotom-60 is factory mounted with an electric cable. Mount a plug on the cable according to the following:

EU cable		UL cable	
L ₁	Brown	L ₁	Black
L ₂	Black	L ₂	Red
L ₃	Black or grey	L ₃	Orange/ turquoise
Earth	Yellow/green	Earth	Green (or Yellow/green)
Neutral	Blue (Not used)	Neutral	White (Not used)

Direction of the Cut-off Wheel

Check that the cut-off wheel rotates in the direction indicated by the arrow on the cut-off wheel guard. If the direction of rotation is incorrect:


EU cable	UL cable
Switch two of the phases	Switch phases L ₁ and L ₂

Discotom-65

Discotom-65 is factory mounted with an electric cable and an additional wire for earth connection. Both MUST be hard-wired into the mains (ref. EN 50178 / 5.2.11.1):

EU cable		UL cable	
L ₁	Brown	L ₁	Black
L ₂	Black	L ₂	Red
L ₃	Black or grey	L ₃	Orange/ turquoise
Earth	Yellow/green	Earth	Green (or Yellow/green)
Neutral	Blue (Not used)	Neutral	White (Not used)

- The second Protective Earth must be earthed by connection to a terminal marked ⊕ or PE.
Alternatively, use a system whereby there is an automatic disconnection of the supply in case of loss continuity of the protective conductor.

Important:

For Electrical Installations with Residual Current Circuit Breakers

For Discotom-65 connected to electrical installations with residual current circuit breakers, a residual current circuit breaker, type B time delayed, 30 mA is REQUIRED (ref. EN 50178 / 5.2.11.1).


For Electrical Installations without Residual Current Circuit Breakers

The equipment must be protected by an insulation transformer (double-wound transformer)

Please contact a qualified electrician to verify which option is suitable for the local installation setup.

Both requirements refer to the European standard EN 50178 / 5.2.11.1. Similar standards apply in North America.

Connection to an External Exhaust System

Struers recommends the use of an exhaust system, as workpieces may emit harmful gases when cut. The exhaust system will also reduce the level of water condensation on the sides of the protection guard.

To connect the Discotom to an exhaust system:

- Remove the red cap from the flange.
- Mount an exhaust hose from your local exhaust system onto the flange (50 mm (approx. 2") dia.).

Important

If no exhaust is connected, damp air (produced by the cutting process) may escape from the cutting chamber and penetrate into other areas of the cabinet. This may cause damage to components and shorten the lifetime of the machine.

Connecting a Recirculation Cooling Unit

To ensure optimal cooling, Discotom can be fitted with a Cooli unit. Cooling System 5 is a Cooli configuration designed for use with Discotom.

Note:

Cooling System 5 includes a static filter. For intensive use, and for materials generating a lot of swarf, a bandfilter such as Coolimat-200 is recommended.


Note

Before connecting the cooling unit to the Discotom, follow the instructions in the Cooling Units Instruction Manual to prepare it for use.


To connect the Discotom to a Recirculation Cooling Unit:

- Plug the Cooli control unit's communication cable into the Discotom's control socket.
- Connect the water inlet hose to the Cooli pump using the quick coupling (A).
- Connect the other end of the hose to the Discotom water inlet.

- Insert the connector pipe ① in the water outlet on the back of Discotom and mount the 90° elbow pipe ②. Lubricate the sealing ring with grease or soap to facilitate insertion. (Use the other elbow pipe if more suitable).
- Mount the outlet hose ③ onto the elbow pipe and clamp using a hose clamp.
- Check that the outlet hose slopes downwards when connected. If necessary adjust the length of the hose.
- Insert the other end of the hose into the mounting hole in the bracket on top of the Cooli filter unit (B).
- Connect the cooling unit to the mains power supply.

IMPORTANT


Before connecting, check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.

Drainage Hose

- Place the drainage hose in the tank of the Cooling Unit or lead to drain.




2. Basic Operation

Using the Controls Front Panel Controls of Discotom

Control name	Symbol
① DISPLAY	
② MENU KEY	
③ ESCAPE	ESC
④ MULTI-FUNCTION KNOB	
⑤ JOYSTICK	
⑥ START.....	
⑦ STOP.....	
⑧ EMERGENCY STOP.....	


Front Panel Controls

Name	Item	Function	Name	Item	Function
① DISPLAY		The Discotom's display.	⑤ JOYSTICK		Move up-or down to position the Y- table.
② MENU KEY		Menu dependent multi-function key. See the bottom line of the individual screens.	⑥ START		Starts the machine and recirculation unit and/or band filter.
③ ESCAPE	Esc	Moves one step backward in menus. If modified parameters have not been stored, changes are lost.	⑦ STOP		Stops the machine and recirculation unit and/or band filter.
④ MULTI-FUNCTION KNOB		Push knob to select function. Turn knob to adjust settings. Push knob to store modified settings.	⑧ EMERGENCY STOP		Push the red button to stop. Pull the red button to release.

Cooling Valve Positions

The cooling valve switches the flow of cooling fluid between the cutting chamber (cutting) and the flushing hose (cleaning). During cutting the valve is placed in the vertical position to cool the cut-off wheel. During cleaning the valve is turned clockwise to divert cooling fluid to the flushing hose.

Please Note

After cleaning, turn cooling valve back to vertical position. START  can not be activated unless the cooling valve is in the cutting position.

Flush Hose

■ Press down the lever, on the top of the flush hose, to open the hose's on/off valve and to adjust the flow of water. For instructions on how to use the flushing hose when cleaning the Discotom, see "Cleaning".

IMPORTANT

Do not move the cooling valve to the cleaning position until the flush hose is pointing inside the cutting chamber.

Cleaning of Flush Hose Nozzle

The flush hose nozzle may collect swarf, inhibiting the flow of cooling fluid. To clean, unscrew nozzle head and rinse under clean water.

Setting the Language

When the Discotom is started for the first time, a screen display prompts selection of a preferred language. To change the language:

From the CUTTING menu, press MENU  button once to select CONFIGURATION Menu.

Turn knob to navigate between parameters in the CONFIGURATION Menu.

Push knob to select LANGUAGE. A pop-up menu appears.


From the pop-up menu, select your preferred language and enter the setting by pushing the knob.

Press ESC button to move from CONFIGURATION Menu to Cutting menu.

Moveable Table Y-Table

The Y-table is a motor driven, movable table. Using the joystick (see “Control Panel Functions”) the table can be moved backwards and forwards. The Y-table is especially useful for cutting wide workpieces.


Note:

When Discotom is switched on, the cutting table will move to the front reference position (the latest position of the cutting table). However, when Discotom has been switched on and off 10 times, the cutting table will firstly move to the front of the cutting chamber and then to the back. This automatic movement helps keep the two guide shafts free of cutting swarf.

Reading the Cutting Display

The Cutting Display displays three types of information:

- A:** Cutting Mode
- B:** Cutting Parameters, and
- C:** Motor Information

Cutting Mode

The top area of the display (**A**) displays the selected Cutting Mode: Manual Cutting or Automatic Cutting mode.

Note:

In Manual Cutting mode, Cutting Parameters are not displayed.
This area of the display will be blank.

Cutting Parameters

In Automatic Cutting mode, the middle area of the display (**B**) displays information about the Cutting Parameters: Feed speed, Stop position and Motor RPM (D-65 only).

The Cutting Parameters can be set both before and during cutting. The set value is displayed to the left of the bar graph. The actual value (during cutting) is displayed to the right of the bar graph.

Motor Information

The bottom area of the display (**C**) displays Motor information: Motor load and Motor temperature. The values displayed are in relative (%) values.


Changing Cutting Mode and Cutting Parameters

Changing Cutting Mode


The Discotom has two cutting modes: Manual and Automatic. To toggle between these two modes:

Turn knob to highlight cutting mode.

Push knob to toggle between the two cutting modes.


Changing Cutting Parameters


In automatic cutting mode the Discotom applies the selected cutting parameter values for: Feed speed, Stop mode and Motor RPM (D-65 only).

To adjust the values of these parameters:

- Turn the knob to highlight a selected cutting parameter.
- Push the knob to allow editing of the highlighted parameter. Turn the knob to change the value of the parameter.
- Push the knob to store the new value.

Turn knob to select Cutting Parameter.
 Push knob to allow editing of Cutting Parameter.
 Turn knob to adjust value.
 Push knob to store new value.

Actual values of the cutting parameters: Feed speed, Stop position and Motor RPM (D-65 only) are displayed to the right of the columns (A).


Turn the knob to change the value of the selected cutting parameter. The arrow above the column will move to reflect the new setting (B).

Setting the Feed Speed

To set the Feed speed:


Turn knob to highlight the current Feed speed.

Push knob to select the current Feed speed.

Turn knob to change the Feed speed.


Push knob to store the new Feed speed value.

Selecting a Stop Mode

To select Stop mode:


Turn knob to highlight the current Stop mode.


Push knob to display available modes.

Turn knob to highlight the mode to be used.

Push knob to select highlighted mode.


Changing Motor RPM


To set the Motor RPM:


Turn knob to highlight the current Motor RPM.

Push knob to select the current Motor RPM.


Turn knob to change the Motor RPM.

Push knob to select the Motor RPM.

Note:

Motor RPM can be varied during cutting.
After selecting the current RPM, the speed can be varied
by turning the knob whilst cutting is in progress.
The optimal RPM value can then be stored by pressing the knob.

*Using Motor Load and
Temperature Display*

The motor load and temperature values displayed are relative percentage (%) values.

Manual Cutting Mode

Motor load and temperature displays indicate how much force is being applied to the sample and how hot the motor is. As more force is applied, motor load is increased resulting in a higher temperature. A heavy force applied for an extended period of time, can result in the motor's temperature becoming greater than the safe working limit. Discotom will then automatically stop cutting to prevent damage to the motor.

Applying heavy force for a long period of time may also shorten the life of the cut-off disc.

Automatic Cutting Mode

Feed speed is automatically reduced on overload by the Automatic Feed Speed Control feature.

Automatic Feed Speed Control

This feature protects the Discotom from damage due to continuous overloading of the motor (150% overload). It can also help quickly find the optimal feed speed for cutting samples.


If the feed speed is set too high and the motor is overloaded:

- The cutting table automatically stops.
- The table retracts 2 mm.
- The table returns to the position where the cutting stopped.
- Cutting resumes at the original preset speed for approximately 0.2mm (to remove damaged/ burnt material on both cut-off wheel and sample).
This is called the dressing phase.
- If, during the dressing phase, the motor load stays within operating limits (the motor is not overloaded), cutting resumes using the preset feed speed.
- If, however, during the dressing phase, the motor is overloaded, the feed speed is reduced by 20%, after which cutting is resumed.
- If the motor continues to be overloaded, this cycle is repeated until the feed speed has been reduced to 10% of the original, preset value.

When cutting of the sample is completed, without overloading the motor, the optimum feed speed for this particular workpiece has been found. The operator can then use this feed speed to cut similar workpieces.

Feed Speed Reduced

If during a cutting operation the feed speed has been reduced one or more times, a warning triangle will be displayed beside the Feed speed bar graph (see, below).

Overload Protection

Discotom is equipped with overload protection to prevent damage to the cutting motor.

If the cutting motor has been continuously overloaded, a thermal cut-out will stop the motor.

- Open the protection hood and allow the motor to cool down before resuming cutting.

Fitting or Changing the Cut-off Wheel

IMPORTANT

The spindle of the Discotom is a left-hand thread.

- If Discotom has been used for manual cutting, secure the cut-off wheel with the locking handle on the right hand side of the machine before changing the cut-off wheel.
- Press and hold down the spindle locking knob on the right-hand side of the cut-off wheel, turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner. Remove the flange, cardboard washers and the old cut-off wheel.
- Mount the new cut-off wheel.

IMPORTANT

Conventional cut-off wheels based on Al_2O_3/SiC abrasives should be placed between two cardboard washers, to protect the cut-off wheel. For maximum precision with diamond or CBN cut-off wheels, do not use cardboard washers.

- Mount the flange and the nut. Tighten carefully and release the locking knob.

Clamping the Workpiece

- Clamp the workpiece with the clamping device of your choice, for example, a quick clamping device. Place the workpiece between the clamp and the back stop.
- Push the clamp towards the workpiece and lock the clamp on the quick-clamping device.

Generally, it is recommended to clamp the workpiece as far back in the cutting chamber as possible.

IMPORTANT

Prior to cutting, ensure that the workpiece is held firmly in place by a clamping device.

Positioning the Cutting Table

Prior to cutting, the cutting table is positioned using the Y-table joystick.

Starting/Stopping the Cutting Process

Automatic Cutting

When the machine is used in Automatic cutting mode the cut-off wheel is stationary and the cutting table moves. In Manual cutting mode the cutting table is stationary and the operator moves the cut-off wheel.

- Clamp the workpiece.
- Loosen the cut-off wheel locking handle situated on the right-hand side of the cutting chamber.
- Lower the cut-off wheel by pulling the cutting handle downwards until the cut-off wheel is positioned ready to cut the workpiece. Re-fasten the locking handle.
- Press the Y-table joystick upwards to position the cutting table and the workpiece just in front of the cut-off wheel.
- Close the protection guard.
- Set the desired feed speed.
- If necessary, select Automatic cutting mode (see, "Changing Cutting Mode").
- Press the start button ; the cut-off wheel starts rotating and the cooling water starts to flow. The cutting table moves towards the cut-off wheel at the preset feed speed.
- When the cutting table reaches the back of the cutting chamber, the cut-off wheel automatically stops.
- Dependent on the return position setting, the Y-table will stay where it is, return to where cutting started or return to its zero position (at the front of the cutting chamber).

Note

The cutting table will only move when Automatic Cutting mode is selected.

Note

During cutting, it is possible to move the sample away from the cut-off wheel by pressing the Y-table joystick downwards.

Manual Cutting

- Clamp the workpiece.
- Loosen the locking handle on the right side of the cutting chamber.
- Position the cutting table and the workpiece under the cut-off wheel with the Y-table joystick.
- Close the protection guard.
- Ensure that Manual Cutting mode is selected .
- Press the start button  to start the machine. The cut-off wheel starts rotating and the cooling water is turned on.
- Pull down the cutting handle and let the cut-off wheel work itself into the workpiece. Increase the pressure and begin cutting.
- When the cut-off wheel has almost cut through the workpiece reduce the pressure.
- When the cut-off wheel has cut through the workpiece push the cutting handle back to its top position.
- Press the stop button  to stop the machine.

Combining Manual and Automatic Operation

IMPORTANT

When cutting manually, the motor load information should be used to monitor the force on the workpiece.

Manual and Automatic cutting modes may be used in combination.

- Cut into the workpiece in Manual mode.
- Secure the locking handle in this position, select Automatic Cutting mode and continue in Automatic mode; the workpiece will move towards the cut-off wheel.

Stop Modes


Discotom has 3 different stop modes:
Relative stop position
Absolute stop position
Automatic stop

Relative Stop Position

The Y-table is set to stop in a position relative to where cutting starts. After entering the approximate sample size, the cutting process will stop immediately after the sample has been cut through. The positioning range is 0-200 mm.


Insufficient Y-Table Travel

If, when the relative stop position is being adjusted, there is insufficient Y-table travel left for the selected stop setting, a warning triangle will be displayed beside the Stop (Rel.) bar graph (see, below). For example, if the Y-table is in its mid-position (100 mm) and the relative stop position is set to 150 mm. Moving the Y-table towards the front of the cutting chamber increases the amount of travel and this automatically updates the relative stop data. When the Y-table reaches the point where there is sufficient travel for the selected stop setting, the warning triangle will no longer be displayed.

Absolute Stop Position

The cut-off wheel is set to stop in a position measured from the zero position where the cutting table is at the very front of the cutting chamber. The positioning range is 0-200 mm.

Automatic Stop

When Automatic stop mode is selected the machine automatically stops when the workpiece has been cut through. For normal cutting, Automatic stop function is recommended.

Note:

When cutting tubes or other non solid samples, Automatic stop is **not** recommended. Use either **Relative** or **Absolute** stop instead.


Automatic Stop Explained

The automatic stop function is based on monitoring the electrical current used by the cut-off wheel motor. For automatic stop to function properly, two criteria have to be met:

- The current used at the **start** of a cutting operation has to exceed a minimum value. If it does not, **Automatic** stop will not be enabled.
- The current used **during** the cutting operation has to remain above a minimum value. When it falls beneath this value, cutting will stop.

In some samples, for example tubes, the current may fall beneath the minimum value when the wall of the tube has been cut through and the centre (hollow part) of the tube is reached during the cutting operation. If this happens, cutting will stop even though the sample has not been fully cut through.

When automatic stop is enabled during a cutting operation, this is indicated by a graph icon beside the Auto stop bar graph (see, below).

Variable Motor RPM (D-65 only)

Discotom- 65 has a variable rotational speed (Motor RPM) feature. This allows the operator to optimise cutting quality and cut many different materials using the same cut-off wheel. Materials as different as aluminium and hardened steel can then be cut using the same cut-off wheel.


Recommended Settings:
 Rotational Speed
 (Discotom-65 only)

The table below is a guide to the recommended RPM settings for popular materials using a single cut-off wheel (Struers 40A30).

Rotational Speed Settings using 40A30cut-off wheel			
Description	Material Tested	Hardness (HV)	Speed (RPM)
Non ferrous soft metal	Aluminium	50-110	3000
Very ductile metal	Stainless steel	220	2200
Medium soft ferrous metal	UHB IMPAX Cold work tool steel	300	2200
Medium hard metal	UHB IMPAX Cold work tool steel	480	2000
Hard ferrous metal	UHB ARNE Cold work tool steel	750	1750

Please note that the above settings are based on laboratory tests of the specific materials listed above. Any deviation in material will require a different setting. Interpolate the Vickers hardness to approximate the corresponding rotational speed.

Feed Speed: The corresponding feed speed setting depends on the sample size; Recommended Settings for 40A30 please refer to the graph below.

To use the graph, find the maximum diameter of your workpiece on the x-axis, then move upwards to the curve relevant to your material. Then move left to the y-axis to find the corresponding feed speed.

Cut-off wheel recommendations

40A30 is recommended for all-round cutting of materials <800 HV.
50A30 is recommended for all-round cutting of materials >800 HV.

Note: Some materials cannot be cut using the 40A30 universal cut-off wheel; for example Titanium.

Important

Always adapt **both** rotational speed (D-65 only) **and** feed speed to a different type for material.



Cleaning

Flushing the Cutting Chamber

To maximise the Discotom's lifetime Struers strongly recommends daily cleaning of the cutting chamber with the flush hose.

- Press the stop button  to stop the cut-off wheel and cooling water.
- Open the cutting chamber hood and release the workpiece from the clamping device. Remove the clamping device if carrying out a more thorough cleaning.
- Pull out the flush hose from its holder and point it towards the bottom of the cutting chamber.
- Turn the cooling valve clockwise to cleaning position and adjust the flow as needed.

The display will show the following:

- Clean the cutting chamber thoroughly.
- Turn the cooling valve back to cutting position.

Note:

When flushing, do not spray cooling fluid on the unsealed areas indicated in the illustration to the left as this may cause leakage.

AVOID CORROSION!

After cutting, leave the cutting chamber hood open to let the cutting chamber dry completely.

3. Routine Maintenance

Accumulated dirt and swarf can restrict or cause damage to the movement of the cutting table.

To ensure a longer lifetime for your Discotom Struers strongly recommends daily cleaning of the cutting chamber.

Clean the cutting chamber thoroughly if the Discotom is not to be used for a longer period of time.

Note:

For more detailed service and maintenance instructions relevant to the cooling unit, read the Cooling Units Instruction Manual.

Daily Service

- Clean the cutting chamber, especially the T-slot cutting table.
- Do **not** clean the lamp glass with alcohol. Use a damp cloth only.

Weekly Maintenance

- Remove clamping device
- Clean the cutting chamber thoroughly:
Move the cutting table forwards and backwards to access the whole of the cutting chamber.
 - Clean along the length of the guide shafts with the flushing gun and a brush to remove accumulated swarf.
 - Clean under the cutting table with the flushing gun and a brush to remove accumulated swarf.
- Wipe the gutter on the inside of the cutting-chamber hood and remove accumulated swarf.

Cooling Unit

- Check the level of the cooling water after 8 hours use or at least every week.

Monthly Maintenance

Replacing the Cooling Water

- Replace the cooling water in the Recirculation Cooling Unit at least once a month.

Yearly Service
Inspection of Guard

- Visually inspect the guard and the glass for signs of wear or damage.


Important

Carry out inspection at more regular intervals if Discotom is used for more than one 7 hour shift a day.

Struers recommends that the PETG glass is replaced after 5 years of routine operation.

The guard should be replaced immediately if it has been weakened by collision with projectile objects or if there are visible signs of deterioration as a result of using a cooling fluid other than those produced by Struers.

A label on the guard indicates when the glass is due to be replaced.


Reference Guide

Table of Contents	Page
1. Advanced Operations	36
Configuration Menu	36
Selecting and Using the Configuration Menu	36
Configuration Parameters	37
New Pass Code	38
Changing Operation Mode	39
Clamping Irregular Workpieces	40
Clamping and Cutting Long, Protruding Workpieces	40
Mounting an Adjustable Cutting Table	40
Optimising Cutting Results	41
Explanation of Safety Factors	41
2. Accessories	42
3. Consumables	43
4. Trouble-Shooting	44
Error Messages	47
5. Maintenance	54
Service Information	54
Maintenance of Cutting Table	55
Lubricating the Cutting Table	55
Maintenance of Cut-off Wheels	57
Storing Bakelite Bonded Al₂O₃ Cut-off Wheels	57
Maintenance of Diamond and CBN Cut-off Wheels	57
Maintenance of Clamping Devices	57
6. Technical Data	58
Cutting Capacity	60

1. Advanced Operations

Configuration Menu

Selecting and Using the Configuration Menu

From the CUTTING menu, press the  button once to select CONFIGURATION menu.


Turn knob to highlight parameters in the CONFIGURATION menu.
Push knob to allow editing of the highlighted parameter.
Turn knob to adjust the setting.
Push knob to enter the new setting.

Press Esc button to move out of the CONFIGURATION menu back to the CUTTING menu.

Configuration Parameters

- Display contrast: The contrast settings of the display can be adjusted to suit individual preferences (default value: 30, range 0-50).
- Units: The Feed and Stop values in the display panel can be set to be displayed in either mm (default) or inches.
- Return position: After cutting or after pressing stop button , the return movement of the cut-off wheel can be set to one of three different modes:
Start Position: Discotom automatically retracts the Y-table to its original position at the time Start  was pressed.
Zero Position: Discotom automatically fully retracts the Y-table to the front of the machine.
Stay: Discotom does not move the Y-table after cutting.
- Important**
Use the Stay function for bakelite bonded diamond or CBN cut-off wheels, as retraction might destroy the rim of the cut-off wheel.
- Language: The language can be set to English (default), German, French, Spanish, Japanese or Chinese.
- Operation mode: It is possible to select three different operation modes. Different operation modes provide different levels of access to parameters as follows:
Configuration: Full functionality, access to all parameters.
Development: No access to parameters in CONFIGURATION menu, except Display contrast:
Production No access to parameters except Display contrast and Stop parameters for Absolute and Relative stop positions if they are being used.
- Cooli type: Indicates the type of Cooli unit that is connected to the Discotom.

New Pass Code

Please Note


When a pass code is set the operator has 5 attempts to enter the correct pass code after which the Discotom will be locked.
Re-start Discotom using the Main Switch
then enter the correct pass code.

To set a pass code:

- From the CONFIGURATION menu, select Operation mode to enter the OPERATION MODE PARAMETERS menu.

Push knob to select Pass code.

Use the menu key  to move the cursor (digit selected). Pressing the key once moves the cursor one place to the left. (The default pass code is '2750'). Turn knob to change the digit then press the menu key to move to the next digit.

Push knob to enter the pass code.
If correct, the pass code will be accepted and the display will change.

Turn knob to move the cursor to New pass code.


Push knob to select New pass code. Enter your new pass code as above.

Changing Operation Mode

To change the operation mode:

- From the CONFIGURATION menu, select Operation mode to enter the OPERATION MODE PARAMETERS menu.


 Push knob to select Pass code.

OPERATION MODE PARAMETERS
Operation mode: Configuration
Pass code: [REDACTED]

 Use the menu key  to move the cursor (digit selected). Pressing the key once moves the cursor one place to the left. (The default pass code is '2750'). Turn knob to change the digit then press the menu key to move to the next digit. Push the knob to enter the code.

OPERATION MODE PARAMETERS
Operation mode: Configuration
Pass code: ****
New pass code: ----

 Push knob to select Operation mode.

SELECT OPERATION MODE

Production	▲
Development	
Configuration	▼

 Select the desired operation mode and push knob to confirm.

- Clamping Irregular Workpieces** Irregular workpieces (without flat clamping surfaces) must be clamped using special clamping tools, as the workpieces must not move during cutting. This could result in damage to the cut-off wheel or to the sample itself. Use the T-slots to mount the special clamping tools. Struers offers a large selection of Clamping Tools (See "[Accessories](#)").
To achieve faster cutting, position the workpiece so that the wheel will cut the smallest possible cross-section.
- Clamping and Cutting Long, Protruding Workpieces** The left-hand side cover of the protection hood can be replaced with a specially designed safety box extension. Please refer to the section on Accessories for details.
- Mounting an Adjustable Cutting Table** To cut flat parallel sections, an adjustable table may be mounted on top of the existing cutting table. The adjustable table positions the workpiece at a right angle to the cut-off wheel and has a sideways movement of 60mm.
Please refer to the section on [Accessories](#).

Optimising Cutting Results

The following table offers answers to a number of common questions about how to achieve better quality cuts:

Question	Answer	Comments
How can I avoid discoloration or burning of the sample?	Reduce rotational speed. (D-65 only).	Causes increased wheel wear.
	If reduced rotational speed fails to solve problem, change to a softer cut-off wheel*	
How can I avoid unplane cuts?	Reduce rotational speed. (D-65 only).	Causes increased wheel wear.
How can I avoid smearing?	Reduce rotational speed. (D-65 only).	Causes increased wheel wear.
How can I avoid burrs?	Use a softer cut-off wheel*	Causes increased wheel wear
	Secure the workpiece on both sides of the cut-off wheel	
How can I prevent the cut-off wheel wearing out too quickly?	Increase rotational speed.	May cause sample discoloration and unplane cut.
	Use a harder cut-off wheel*	
How can I cut faster?	Place the cut-off wheel as low as possible.	
	Place the workpiece in a position that allows the cut-off wheel to cut the smallest possible cross-section.	
	Increase feed speed.	May cause sample discoloration and unplane cuts.
How can I avoid machine vibrations?	Minor vibrations: Increase feed speed in steps of 0.1 m/s.	May cause sample discoloration and unplane cuts.
	Major vibrations: Increase rotational speed by 500 rpm. (D-65 only).	May cause sample discoloration and unplane cuts.

* Please refer to the Selection Guide in the Struers [Struers Cut-off Wheels brochure](#).

Explanation of Safety Factors

The cutting chamber hood has a safety switch to prevent the cut-off wheel from starting while the hood is open. Furthermore, a locking mechanism prevents the operator from opening the hood until the cut-off wheel stops spinning.

2. Accessories

Please refer to the [Discotom-60/-65 brochure](#) for details of the range available.

Specification	Cat. No.
Safety Extension Box For extending the cutting chamber when cutting long samples. To be mounted on the left hand side of Discotom-60/-65. For samples with a maximum size of h x w: 80 x 140 mm / 3.15" x 5.5". Complete with rubber curtain. Box dimensions (h x d x w): 180 x 460 x 600 mm.	05906911

Clamping Tools

Please refer to the [Struers Clamping Tools brochure](#) for details of the range available.

3. Consumables

Cut-off Wheels

Please refer to the Selection Guide in the Struers
[Struers Cut-off Wheels brochure](#).

Other Consumables

Specification	Cat. No:
<i>Corrozip</i> <i>Additive for Cooling Fluid</i> Environment friendly. To protect the machine from corrosion and to improve cutting and cooling qualities. 1 l 5 l	 49900045 49900046
Grease for maintenance/lubrication of the spindle	16080802
Oil for maintenance of cutting table	16080845

4. Trouble-Shooting

Error	Explanation	Action
Machine Problems		
Machine has stopped cutting.	The AutoStop function is active.	Switch unit off and on at Main Switch to reset AutoStop.
During cutting the cutting table stops when reaching a specific point.	The AutoStop function is active.	Switch unit off and on at Main Switch to reset AutoStop.
Chamber light does not work.	Replace the lamp.	Remove plastic cap to gain access to the fluorescent lamp. Pull out the lamp and replace it.
Water leaking.	Leak in a Recirculation Cooling Unit hose.	Check the hose and tighten the hose clamp.
	Water overflow in the cooling water tank.	Remove the excess water in the tank.
	During cleaning, water has been sprayed onto unsealed parts of the cutting chamber	Avoid spraying unsealed parts of the cutting chamber during cleaning, see " Cleaning ".
Workpieces or cutting chamber rusty.	Insufficient additive in cooling fluid.	Use Struers Additive for cooling fluid in the cooling water, at the correct concentration. Check with a refractometer. Follow the instructions in the ' Maintenance ' Section.
	The protection guard is left closed after use.	Leave the protection guard open to let the cutting chamber dry.

*Discotom-60 /-65
Instruction Manual*

Error	Explanation	Action
Cutting Problems		
Discoloration or burning of the workpiece.	The hardness of the cut-off wheel is inappropriate for the hardness / dimensions of the workpiece.	Reduce rotational speed (D-65 only) or change wheel. See "Consumables" section, under heading 'Cut-off Wheels'.
	Inadequate cooling.	Check that there is sufficient water in the Recirculation Cooling Unit.
Unwanted burrs.	Cut-off wheel too hard.	Reduce rotational speed (D-65 only) or change wheel. See "Consumables" section, 'Cut-off Wheels'.
	Feed speed too high at the end of the operation.	Reduce the feed speed near the end of the operation.
	Lack of support for the workpiece.	If possible, support the workpiece on both sides.
Cutting quality differs.	Cooling water hose clogged.	Clean the cooling water hose and the cooling tube. Check the water flow by turning the cooling valve to cleaning position.
Cut-off wheel breaks.	Incorrect mounting of the cut-off wheel.	Check that the bore/centre hole has the correct diameter. Check cardboard washer on both sides of the cut-off wheels and replace if worn. The nut must be tightened properly.
	Incorrect clamping of the workpiece.	Make sure that only the left quick-clamping device is tight. The spring clamp to the right should only press lightly. Use the vertical clamping system if the workpiece is an irregular shape.
	Cut-off wheel is too hard.	Reduce rotational speed (D-65 only) or change wheel. See "Consumables" section, 'Cut-off Wheels'.
	Feed speed is set too high.	Reduce the feed speed.
	Inadequate cooling.	Check that there is enough water in the cooling unit.

*Discotom-60 /-65
Instruction Manual*

Error	Explanation	Action
The cut-off wheel wears down too quickly.	The rotational speed is too low.	Increase rotational speed (D-65 only).
	The feed speed is too high.	Reduce feed speed.
	Insufficient cooling.	Check that there is enough water in the Recirculation Cooling Unit. Check the cooling water hoses.
The cut-off wheel does not cut through the workpiece.	The rotational speed is too low.	Increase rotational speed (D-65 only).
	The cut-off wheel is too soft for the task	See "Consumables" section, 'Cut-off Wheels'.
	Incorrect choice of cut-off wheel.	See "Consumables" section, 'Cut-off Wheels'.
	Cut-off wheel worn.	Replace the cut-off wheel.
The workpiece breaks when clamped.	The cut-off wheel gets caught in the workpiece.	Support the workpiece and clamp it on both sides of the cut-off wheel so that the cut stays open.
	The workpiece is brittle.	Place the workpiece between two polystyrene plates. Note! Always cut brittle workpieces very carefully.
The sample is corroded	The sample has been left in the cutting chamber for too long.	Remove the sample directly after cutting. Leave the cutting chamber hood open when you leave the machine.
	Insufficient additive for cooling fluid.	Use Struers Additive for cooling fluid in the cooling water at the correct concentration. Check with a refractometer. See "Maintenance" section.

Error Messages

Error messages are divided into two classes:
Messages
Errors

Messages

Messages are intended to inform the operator of the machine's progress and advise about minor operational errors.

Errors

Errors must be rectified before cutting can be continued. In some cases, cutting cannot continue before an authorised technician has rectified the error. Turn off the unit at the Main Switch immediately. Do not attempt to operate unit before a technician has rectified problem.

Message/Error	Type	Explanation	Action
Start denied, process menu not selected	MESSAGE (#0)	The start button  was pressed in a menu where it's not possible to start cutting	Return to the Cutting menu, ensure that all cutting parameters are set to an appropriate value, press start button  again to begin cutting.
Are you sure you want to change pass code?	MESSAGE (#2)		Press the menu button  to accept the new pass code. Press Esc button to reject the new pass code and continue to use the current one.
Protection hood not closed! Close the hood and press the menu button  .	MESSAGE (#3)	The cutting chamber hood was open when the Discotom was switched on.	Close the hood and press the menu button  .
Discotom is searching for reference position, please wait	MESSAGE (#4)		Wait until the reference position has been found and the display changes. On every 10th start up, Discotom will perform an extended reference routine during which the Y-table will be moved fully forward and backward. This is to ensure that the table's tracks are kept clean.
The reference position is found, Discotom is ready.	MESSAGE (#5)		
Cutting finished, stop position reached	MESSAGE (#6)		

*Discotom-60 /-65
Instruction Manual*

Message/Error	Type	Explanation	Action
Cutting finished, end position reached	MESSAGE (#7)		
Process in progress!	MESSAGE (#8)	Operator tried to change certain parameters that it is not possible to change during a cutting operation, for example, "Units".	Wait until the cutting operation has finished before carrying out this change.
Process stopped	MESSAGE (#9)	Cutting operation terminated because the STOP button was been pressed.	
Do not operate the valve during cutting!	MESSAGE (#10)	The cooling valve was turned to the cleaning position during cutting.	Turn the cooling valve back to the cutting position. Wait until after the cutting operation has finished before using the cooling valve.
Too many feed reductions, process is aborted!	MESSAGE (#11)	It was not possible to complete a cutting operation. The automatic feed speed control feature reduces in steps of 20% a maximum of 4 times after which, if motor load is still too high, the cutting operation is terminated.	The three most likely causes are: The feed speed was much too high, before cutting started. The cut-off wheel being used is not suitable for cutting the sample. The motor RPM is not appropriate (D-65 only). Consult Struers or your local dealer for information on how to cut this type of material.
Editing restricted by operating mode	MESSAGE (#12)	Operator tried to change certain parameters that are locked when Discotom is in the current operating mode.	If the parameter has to be changed, either the operation mode has to be changed to a menu level where the parameter can be changed, or the operation mode can be changed to "Configuration" mode (in which all parameters can be changed). Changing the operation requires that the operator enters the correct pass code.

*Discotom-60 /-65
Instruction Manual*

Message/Error	Type	Explanation	Action
Feed speed reduced: From x.xx mm/s to y.yy mm/s	MESSAGE (#14)	An automatic feed reduction occurred during cutting.	For information only.
Do you want to continue with the reduced feed speed?	MESSAGE (#15)	During cutting, feed speed was automatically reduced.	Press "Yes" to use the adjusted (reduced) feed speed for the next cutting operation. Press "No", to use the original feed speed.
Cutting motor brake error!	ERROR (#16)	The braking operation to reduce the cut-off wheel's rotation took too long.	The braking circuit has developed a fault. Contact a Struers service technician.
LIN-bus error during power on, please call service technician	ERROR (#17)	During power up, 1 or more LIN-bus modules were defective. The LIN-bus modules are small control units connected to a single cable. They perform control functions, such as switching the cutting chamber light on and off, and moving the Y-table.	Switch off and on at Main Switch. If the error message is still displayed, contact a Struers service technician. In certain circumstances (dependent on which module has failed) it may still be possible to operate the Discotom in "Manual" cutting mode.
Cutting table blocked!	MESSAGE (#18)	Operator used the Y-table joystick and the table did not move.	Check that the Y-table's runners are not obstructed. If they are, clear them and try to operate the table again. If they are not, or clearing them does not help, contact a Struers service technician.
Protection hood not closed!	MESSAGE (#19)	The start button  was pressed whilst the cutting chamber hood was open.	Close the hood and press the start button  again to start cutting.

*Discotom-60 /-65
Instruction Manual*

Message/Error	Type	Explanation	Action
Frequency inverter communication error!	ERROR (#20)	Data communication with the frequency inverter failed (D-65 only).	Switch off and on at Main Switch. If the error message is still displayed, contact a Struers service technician.
Contactors K1-K3 not activated, Please check the hood	MESSAGE (#21)	The start button  was pressed, but none of the 3 motor contacts were activated (D-60 only). The cutting chamber hood may not be fully closed.	Check that the hood is closed properly. If the error message is still displayed, contact a Struers service technician.
Cutting motor blocked! Please reduce the motor load	MESSAGE (#22)	During cutting, motor load exceeded its safe limit.	Reduce feed speed or cutting force. Check also that the sample is properly secured as a sample that is not correctly secured can trap the cut-off wheel with enough force to stop it rotating.
Cutting motor will not stop! Please call service technician	ERROR (#23)	A motor terminal failure occurred. After pressing the stop button  , the motor did not stop. (D-60 only).	Switch off the Discotom at the Main Switch and call a Struers service technician.
Frequency inverter fault! Lenze error code: XX	ERROR (#25)	Frequency inverter failure (D-65 only).	Switch off the Discotom at the Main Switch and call a Struers service technician.
Feed motor not stopped!	ERROR (#26)	After pressing stop button  , the motor did not stop.	Switch off the Discotom at the Main Switch and call a Struers service technician.
Cutting table position not found!	ERROR (#27)	During the movement of Y-table, a position was not found within a set amount of time. The feed motor module is most likely defective.	Switch off the Discotom at the Main Switch and call a Struers service technician.

*Discotom-60 /-65
Instruction Manual*

Message/Error	Type	Explanation	Action
Start ignored by frequency inverter. Please check the hood.	MESSAGE (#28)	The start button \diamond was pressed and cutting did not start. The most likely cause is that the cutting chamber hood has not been properly closed (D-65 only).	Check that the hood is closed properly. If the error message is still displayed, Switch off the Discotom at the Main Switch and call a Struers service technician.
Auto stop function not calibrated !	ERROR (#29)	During power up, a motor control function was not activated. (D-65 only).	The Discotom can be used, but the Auto stop function will not work. Call a Struers service technician.
Cutting motor supervision error, contactor K1 not activated	ERROR (#30)	The start button \diamond was pushed, but the K1 motor contact was not activated (similar messages can be displayed showing contacts K2, K3, K4 and K5). (D-60 only).	Switch off and on at Main Switch. If the error message is still displayed, contact a Struers service technician.
Cutting motor supervision error, contactor K1 not deactivated	ERROR (#31)	The start button \diamond was pushed, but the K1 motor contact was not deactivated (similar messages can be displayed showing contacts K2, K3, K4 and K5). (D-60 only).	Switch off and on at Main Switch. If the error message is still displayed, contact a Struers service technician.
Emergency stop activated	MESSAGE (#32)		Pull the button out to re-set the Emergency stop and allow Discotom to be re-started.
Cutting finished, stopped by auto stop	MESSAGE (#33)	Cutting stopped because the amount of current used by the motor fell beneath the Auto stop value.	If the sample has not been cut through, this could be because cutting the sample has not required much motor power. This can be the case with thin walled tube samples. Change the stop position to Relative or Absolute instead of Auto.
No cutting motor current!	ERROR (#34)	At the start of a cutting operation, the motor current is less than 1.0 Ampere. (D-60 only)	Check that the Discotom is correctly connected to the mains power supply and that all 3 phases are functional. Switch off and on at Main Switch. If the error message is still displayed, contact a Struers service technician.

*Discotom-60 /-65
Instruction Manual*

Message/Error	Type	Explanation	Action
Don't activate joystick during power on. Please restart Discotom	ERROR (#35)	Joystick was operated during power up. The joystick must not be operated during start up since the Y-table has to be stationary while a reference position is measured.	Switch off and on at Main Switch (make sure that the joystick is not touched during power up). If the error message is still displayed, contact a Struers service technician.
No cutting motor rotation	ERROR (#36)	Motor rotation of less than 1,000 RPM detected at the start of cutting (D-60 only).	Clear any obstruction away from the cut-off wheel. If the error message is still displayed, and the cut-off wheel is actually turning, contact a Struers service technician.
Cutting table movement blocked!	MESSAGE (#37)	The Y-table stopped moving during a cutting operation.	Clear any obstruction away from the Y-table's runners. If the error message is still displayed, and the cut-off wheel is actually turning, contact a Struers service technician.
Manual cleaning mode. Turn valve back in cooling position as shown when finished. (pictogram)	MESSAGE (#100)	The Discotom is stopped and the cooling valve has been turned to the flushing position.	Follow the on-screen instructions.
Turn valve in cooling position as shown. (pictogram)	MESSAGE (#101)	The cooling valve was not in the cooling position when the cutting operation started.	Follow the on-screen instructions and restart the cutting operation.
Reference pos. for cutting table not found! Discotom can only be used in manual mode - try to restart the machine. If you get this message again, please contact a Struers service technician.	ERROR (#102)	The Y-table was unable to find a reference position.	Switch unit off and on at Main Switch. If the error message is still displayed, contact a Struers service technician. The Discotom can still be used in Manual cutting mode.
The selected stop position is temporary reduced, because the mechanical stop will be reached before the selected position.	MESSAGE (#103)	The relative stop position entered is beyond the remaining Y-table travel.	Move the Y-table toward the front of the machine until there is sufficient Y-table travel for the cutting operation to be performed.

*Discotom-60 /-65
Instruction Manual*


Message/Error	Type	Explanation	Action
<p>It is now time to service your Discotom, please call for a service visit. SERVICE INFO: Total operation time: 3100h Time since last service:1600h Service exceeded by: 100h</p>	MESSAGE (#105)	The time between services has exceeded the recommended number of hours.	<p>Press "F1:OK". Contact a Struers Service Technician to service the machine.</p>
<p>The cutting motor is overheated! You can not use Discotom until the cutting motor has cooled down.</p>	MESSAGE (#106)	Motor temperature has exceeded its safe working limit. The motor has been stopped.	<p>Restart Discotom in manual mode and run unloaded. The motor's built-in fan will cool the motor. Monitor the motor's temperature on the display.</p>

5. Maintenance

Service Information

Struers recommends that a regular service check be carried out after every 1500 hours of use.

Information on total operation time and servicing of the machine is displayed on the screen at start-up:

SMM connected: No/ Yes

Information for Service Technicians on the status of the Struers Memory Module (SMM) that is used to save database parameters, default is "No".

A pop-up message will appear after 1,000 hours operation time to remind the user that a service check should be scheduled.

After the 1,500 hours operation time has been exceeded the pop-up message will change to alert the user that the recommended service interval has been exceeded.

- Contact a Struers Service Technician to service the machine.

Maintenance of Cutting Table

The stainless steel bands are available as accessories and should be replaced if they become worn or damaged. The bands are available as spare parts.

To allow humidity to escape from the cutting table and chamber, it is recommended to leave the hood open when the machine is not in use.

Lubricating the Cutting Table


To maintain the optimum performance of Discotom-60/ -65, lubricate the cutting table at regular intervals (approx every 50 hours).

Check the service counter on the right hand side of the machine regularly to monitor the actual number of usage hours.

After lubrication, note the date and number of service hours on the Maintenance Log table.

Spindle

- Move the cutting table to its rear position using the joystick.


- Place the grease gun on the nipple in front of the spindle for the cutting table and push twice to grease the spindle.

A grease gun with grease for lubricating the Y-table spindle is supplied with the Discotom-60/ -65. When all the grease has been used, refill the grease gun.

Slide Bearings

- Remove the stainless steel bands located over the guide shafts (see illustration).
- Remove the 4 M6 set screws.

- Fill the holes with oil eg. Shell TELLUS oil S100.
Use approx. 25ml oil for each side or until the reservoirs are full.
A bottle with oil for lubricating the slide bearings is supplied with the Discotom-60/ -65. When all the oil has been used, refill with eg. Shell TELLUS oil S100.

Important

If the oil in the reservoir appears milky;
this indicates that water has entered the oil reservoir.
Call a Struers Service Technician to clean the reservoir.

- Re-fit the screws and remount the stainless steel bands.

Maintenance of Cut-off Wheels

Storing Bakelite Bonded Al_2O_3 Cut-off Wheels

Bakelite bonded cut-off wheels are sensitive to humidity. Therefore, do not mix new, dry cut-off wheels with used damp ones. Store the cut-off wheels in a dry place, horizontally on a flat support.

Maintenance of Diamond and CBN Cut-off Wheels

The precision of diamond and CBN cut-off wheels (and thus the quality of the cut) depends on how carefully the following instructions are observed:

- Never expose the cut-off wheel to a heavy mechanical load, or heat.
- Store the cut-off wheel in a dry place, horizontally on a flat support, preferably under light pressure.
- A clean and dry cut-off wheel does not corrode. Therefore, clean and dry the cut-off wheel before storing. If possible, use ordinary detergents for the cleaning.
- Regular dressing of the cut-off wheel is also part of the general maintenance.

Maintenance of Clamping Devices

Important

It is recommended to thoroughly clean and lubricate the Quick Clamping Device and Vertical Clamping Device at regular intervals.

6. Technical Data


Subject		Specification	
		Metric/International	US
CUTTING SPECIFICATIONS			
Workpiece Dimensions (Max)	<i>Workpiece inside cutting chamber:</i>		
	Height	125 mm	4.92"
	Width	220 mm	8.66"
	Depth	195 mm	7.67"
	<i>Workpiece protruding cutting chamber:</i>		
	Height	125 mm	4.92"
	Depth	400 mm	15.75"
Cutting Capacity (Max) <i>Please refer to Cutting Diagram</i>	Max. workpiece diameter	105 mm	4.13"
	<i>Max. size of cut (example)</i>	160 mm at 65 mm height	6.30" at 2.56" height
		120 mm at 90 mm height	4.72" at 3.54" height
PHYSICAL SPECIFICATIONS			
Cutting Motor	Cutting power, constant [S1]	4 kW	5.4 HP
	Cutting power, intermittent [S3]	4.7 kW	6.3 HP
	Maximum power	8 kW	10.8 HP
Cut-off Wheel	Diameter x Thickness x Centre-hole	300 x 2 x 32 mm	11.8 x 0.078 x 1.26"
Rotational Speed	Discotom-60 RPM	2860 rpm at 50 Hz 3460 rpm at 60 Hz	
	Discotom-65 RPM	1000 - 3000	
Positioning and Feed	Positioning range (of cut-off wheel)	0 – 120 mm	0 – 4.72"
	Positioning range (of cutting table)	0 – 200 mm	0 – 7.87"
	Max. positioning speed	20 mm/s	0.79"/s
	Feed Speed range	0.01 – 2.5 mm/s	0.001 – 0.1"/s
Cutting Table	Width	435 mm	17.1"
	Depth	265 mm	10.4"
	T-slots	10 mm	0.39"
Dimensions and Weight	Height (closed/open hood)	570/890 mm	22.44/35.04"
	Width	910 mm	35.8"
	Depth (closed/open hood)	785/816 mm	30.9/32.12"
	Weight	171 kg	377 lbs
	<i>Optional table unit</i>		
	Width	900 mm	35.4"
	Depth	750 mm	29.5"
	Height	800 mm	31.5"
Fume Extraction	Dia. for connecting tube	50 mm	2"
	Recommended capacity at 0mm/0" water gauge	350 m ³ /h	12360 ft ³ /h
Recirculation Cooling Unit	Tank volume (optional)	100 l	26.4 gallons
	Approx. flow	42 l/min	11 gallons/min

*Discotom-60 /-65
Instruction Manual*

Subject	Specification				
Dynamic Balancing	Max. admissible unbalance Upper according to ISO 1940/1, Balance Quality Grade G6.03		Upper = 2 gmm/kg = 110 gmm.		
ELECTRICAL DATA					
Discotom-60					
<i>Mains Cable Specification*</i>	Voltage/ frequency	Cable size	Min. Fuse size	Max. Fuse size	
	3x200-210/50-60	3xAWG12 + PE	25 A	40 A	
	3x220-230/50	3x2,5mm ² + PE	25 A	40 A	
	3x220-240/60	3xAWG12 + PE	25 A	40 A	
	3x380-415/50	3x2,5mm ² + PE	20 A	40 A	
	3x460-480/60	3xAWG12 + PE	20 A	40 A	
	<p><i>*Important:</i> Local standards may overrule the recommendations for the main supply cable. If necessary, please contact a qualified electrician to verify which option is suitable for the local installation setup.</p>				
<i>Residual Current Circuit Breaker</i>	type A, 30 mA (or better) is recommended.				
Discotom-65					
<i>Mains Cable Specification*</i>	Voltage/ frequency	Cable size	Min. Fuse size	Max. Fuse size	Protective Earth Cable
	3x200-210/50-60	3x4mm ² + PE	35	40	4 mm ²
	3x220-230/50	3xAWG8 + PE	35	40	AWG 8
	3x220-240/60	3x4mm ² + PE	25	40	4 mm ²
	3x380-415/50	3xAWG8 + PE	25	40	AWG 8
	3x460-480/60	3x4mm ² + PE	35	40	4 mm ²
	<p><i>*Important:</i> Local standards may overrule the recommendations for the main supply cable. If necessary, please contact a qualified electrician to verify which option is suitable for the local installation setup.</p>				
<i>Residual Current Circuit Breaker</i>	type B time-delayed, 30 mA (or better) is REQUIRED.				
EU Directives		Please refer to the Declaration of Conformity			
Environmental	Noise level	Approx. 67dB(A) running idle, at a distance of 1.0 m / 39.4" from the machine.			

Cutting Capacity

The graph shows the projected cutting capacity under the following conditions: A new cut-off wheel.
The workpiece is laid directly on the cutting table, with overhang where appropriate. Vertical clamping is used. The actual cutting capacity depends on the sample material, cut-off wheel and clamping technique.

Quick Reference Guide

Clamping the Workpiece

- Secure the workpiece using the clamping device of your choice, usually the quick-clamping device. Place the workpiece between the clamp and the backstop.
- Push the clamp towards the workpiece and lock the quick-clamping device with the locking handle.

Automatic Cutting

- Clamp the workpiece.
- Position the cut-off wheel so that it can cut the sample and lock it in place with the cut-off wheel locking handle on the right of the machine.
- Position the cutting table with the Y-table joystick.
- Carefully close the protection hood.
- Select Automatic Cutting mode.
- If necessary, set the feed speed.
- Press the start button  to begin cutting.
- After cutting, the cutting table automatically returns to its start position.
- Open the protection hood and release the clamping device.
- Remove the workpiece.

Manual cutting

- Loosen the locking handle on the right hand side of the machine, so that the cutting handle moves freely.
- Clamp the workpiece.
- Position the cutting table under the cut-off wheel.
- Carefully close the protection hood.
- Select Manual Cutting mode.
- Press the start button .
- Lower the cut-off wheel until it contacts the workpiece and apply pressure until the sample is cut.
- Reduce pressure near the end of the cut.
- Open the protection hood and release the clamping device.
- Remove the workpiece.

Changing the Cut-off Wheel

- Fasten the cut-off wheel locking handle on the right of the machine to stop the cutting handle from moving.
- Press the knob for the spindle lock on the right-hand side of the cut-off wheel, whilst turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner. Remove the flange and the old cut-off wheel.
- Mount the new cut-off wheel.
- Mount the flange and the nut. Tighten carefully.

Flushing the Cutting Chamber

- Point the flush hose towards the bottom of the cutting chamber.
- Turn the cooling valve clockwise to flushing position.
- Press down the lever, on the top of the flush hose, to open the hose's on/off valve and adjust the flow of water.
- Clean the cutting chamber.
- Turn the cooling valve counter-clockwise to the cutting position.
- Always leave the protection hood open to let the cutting chamber dry.

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
Date of Release 09.09.2013

Discotom-60/ -65
Spare Parts and Diagrams

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

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The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2013.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Telefax: +45 44 600 801

Spare Parts and Diagrams

Table of contents

	Drawings	Drawing
	Stepper Motor assembly	15480018J
	Communication plate, assembled	15730070C
	Discotom-60, complete	15900001P
	Base and back plates, assembly.....	15900005V
	Cover, assembly	15900010D
	Cabinet, assembly.....	15900015B
	Front Plate, assembly	15900020J
	Base Assembly	15900025J
	Y-table assembly.....	15900030G
	Spindle Assembly.....	15900040E
	Brake circuit, assembly	15900047C
	Trapezoid spindle assembly.....	15900060G
	Gas Spring Assembly.....	15900065B
	Locking System, assembly.....	15900070A
	Spindle lock, assembled	15900080B
	Valve housing, assembly	15900090C
	Valve knob, assembly	15900095A
	Discotom-65, complete	15910001P
	Base and Back plates, assembly	15910005U
	Front Plate, assembly	15910020J
	Diagrams	Diagram
<i>Discotom-60</i>	Main Supply Circuit page (3 pages)	15903100G
	LIN Lamp Driver (2 pages).....	15903201F
	Wiring Diagram for Trafo.....	15903470A
	Wiring Diagram	15903450H
<i>Discotom-65</i>	Main Supply Circuit page (3 pages)	15913100G
	LIN Lamp Driver (2 pages).....	15903201F
	Wiring Diagram for Trafo.....	15903470A
	Wiring Diagram	15913450H

<p>Some of the drawings may contain position numbers not used in connection with this manual.</p>

Discotom-60/ -65
Spare Parts and Diagrams

The following is a list of the spare parts that may need replacement during the lifetime of the equipment.

To check the availability of other replacement parts, please contact your local Struers Service Technician. It may help identify the part by referral to its position number on the assembly drawings included in this manual.

Spare part list for Discotom-60/ -65

Drawing	Pos.	Spare Part	Cat no:
15900001		Discotom-60, complete	
	15	Motor 3x240VD/50Hz 4kW 2p CSA	2ME59013
		Motor 3x210VD/50Hz 4kW 2p CSA	2ME59012
		Motor 3x420VD/50Hz 4kW 2p CSA	2ME59011
	60	Sealing plate	15900168
	90	Sealing plate	15900162
	110	Bracket for switch	15900181
	120	Switch AZM 170-02ZRKA 24V	2SS00007
	160	Cutting handle	15900177
	190	O-RING 6.30-2.40 72 NBR 872	2IO24008
	200	Front plate, assembl.	15900020
	220	Cover, assembl.	15900010
	260	Spring 3,6x40x146, 9934, 2pcs	2GF29934
	320	Lamp, DULUX S/E 7W I, 2pcs	2HG11782
	330	Reflector, 2pcs	15900144
	405	4.0 AT Cl.-CC fuse 10x38 CSA, 2pcs	2FC11040
		4.0 AT aM-fuse 10x38, 2pcs	2FC10040
		2.0 AT Cl.-CC fuse 10x38 CSA, 2pcs	2FC11020
		2.0 AT aM-fuse 10x38, 2pcs	2FC10020
	510	W1 NETKABEL 60Hz UL/CSA	15113500
		W1 NETKABEL 50Hz	15113501
		W1 NETKABEL 60Hz UL/CSA	15113506
	520	A3 PCB lamp driver, tested	15903001
	560	Cover for PCB	15900208

Discotom-60/ -65
Spare Parts and Diagrams

Spare part list for Discotom-60/ -65

Drawing	Pos.	Spare Part	Cat no:
15900005 Discotom-60 / 15910005 Discotom-65		Base and Back plates assembly	
	10	Base, assembly	15900025
	60	Spindle housing, assembl.	15900040
	70	Key A 8x7x32	2ZF10832
	75	Distance bushing	15900206
	100	Sealing plate	15900163
	130	Guard f cut-off wheel, painted	15909002
	150	Fixed flange	15900160
	160	Dowel pin 6m6x18 D6325	2ZS03458
	170	Flange	14150030
	180	Washer ø40/ø16.5x5. Stainl.st.	14150012
	210	Nut MV. M16 A2 DIN934	2TA70160
	230	Spindle lock, assembl.	15900080
	240	Valve housing, assembly	15900090
	250	Valve knob, assembly	15900095
	280	Hose PVC w. steel coil ø12/ø18, 1.6m	2NU21218
	300	Hose PVC w. steel coil ø12/ø18, 1.6m	2NU21218
	305	Flex hose PMA PIST-23S, 0.5m	2NU31220
	310	Spray gun, assembl.	15320007
	340	Hose PVC, transp. ø8/ø11, 1m	2NU19310
	440	Buffer ø20x15/M6x8. Type D, 2 pcs	2GS00093
	450	Neoprene grommit ø18½/ø25/ø32½-1.5, 2 pcs	2GK90143
	460	Rubber grommit TET 14-20 (IP67)	2GK91420
	470	Shock Absorber OEM 25MB adjust.	2YS10025
	530	Sealing band black. U=14/32mm, 0.446m	2IP28459
	540	Flec.rubb.profile 36510. 5x10, 0.12m	2IP30510
	545	Sealing band black. 31018, 0.2m	2IP30001
	550	Sealing profile 209-0203, 0.245m	2IP00203
	570	Transformer 10-100A	2MT80100
Discotom-60	580	Main switch KG32 K300E	2SE20317
Discotom-65	580	Switch 3x32A, 2x CAM AUX NC	2SE20318
Discotom-60	660	Wiring set for Discotom-60	15903590
Discotom-65	660	Wiring set for Discotom-65	15913590
	800	GEKA hose nipple 1/2in	2NF60000

Spare part list for Discotom-60/ -65

Drawing	Pos.	Spare Part	Cat no:
15900010		Cover, assembly	
	20	Window	15900157
	30	Side plate	15900141
	40	Interlock key	15670378
	70	Sealing profile	15900179
	150	Window, 2pcs	15900180
	190	Sealing band 1038-02 EPDM Sh.A (0.46m)	2IP21038
15900015		Cabinet, assembly	
	60	Brush, 2pcs	15900183
	100	Sealing profile 209-0205 (1 x 208mm + 2 x 367mm)	2IP00205
	110	Sealing strip Grey 1.5x10 (2 x 425mm))	2IP10151
15900020 Discotom-60/ 15910020 Discotom-65		Front plate	
	20	Touch pad, Discotom-60	R5900189
	30	Display 320x240 LED. WH	2HD32024
	40	Window for display	15480465
	70	Main PCB, Discotom 6X programm.	15903000
	75	O-ring 72 NBR 872. 4.30-2.40, 2pcs	2IO24003
	115	2 chann.opt.encod. w.print 24p	2HR12411
	150	O-ring 72 NBR 872. 21.20-3.00	2IO30008
	190	Sealing strip Grey 1.5x10, 0.36m	2IP10151
15900025		Base, assembly	
	20	Rubb.foot grey No.281 ø30x15, 6 pcs	2GB00001
	60	Y-table, assembly	15900030
	80	O-RING 25.02-2.62 72 NBR 872, 2 pcs	2IO26213
	90	O-ring 72 NBR 872. 9.30-2.40, 2 pcs	2IO24021
	105	SHELL TELLUS OIL S 100, 50ml	10162905
	130	Trapezoid spindle, assembly	15900060
	250	Gas spring, assembly, 2 pcs	15900065
	260	Locking system, assembly	15900070
	270	Plate for motor, assembly	15900045
	440	O-ring NBR 70. 50.00-4.00	2IO40080
	450	Hose nipple 2601-9-1/4	2NF40085

Spare part list for Discotom-60/ -65

Drawing

Pos. Spare Part Cat no:

15900030

Y-table, assembly

10	Y-table, assembly	R5900103
20	TRELB, Replacement steel band	05906901
30	Spindle nut w. screws	R5900104
40	Spindle protection	R5900187
50	Bellow protection	15900108
	Y-table guide rail bearing kit consisting of:	15902902
60	Radial bearing SMS777 30/40x30, 2 pcs	
70	Bearing machined, 2 pcs	
80	O-ring 34,50-3,00 70 NBR, 4 pcs	
90	Simmer ring SD 30x40x4, 4 pcs	
100	Oil ring AS 30-40-7-10, 4 pcs	
110	Hex.sock.scr. M5x6 A2 DIN912, 2 pcs	
120	O-ring NBR 70. 16.0-1.50	2IO15045
130	O-ring NBR 70. 32.00-2.00, 2 pcs	2IO20059
180	Sock.setScr.M6x10. 70A2 DIN916 , 6 pcs	2TI10610
190	Lubr.nipple H2 M6x1, stainl.st.	2GN20009
210	Buffer	15900203
230	Plastic cap, Grease nipple	2GT70001

15900040

Spindle housing, assembl.

	Kit, Sealing, Motor consisting of:	15112922
30	Ball bearing 6007 2RS1 ø35/ø62, 2pcs	
60	Key A 8x7x28	
80	Sealing ring	
90	O-ring NBR 70. 31.42-1.78	
110	Sim.ring40-62-7BAUM3SLX7 75FPM	
120	Plate for V-Ring	
130	V-ring V-40A	
140	Namur A01G81.5 5m M8x30	
150	Rotex-28 coupling, machined	
180	O-ring NBR 70	


Spare part list for Discotom-60/ -65

Drawing	Pos.	Spare Part	Cat no:
15900047		Brake circuit, assembly	
	20	35A 800V KBPC3508 BRIGDE RECT.	2VB30750
	40	10R BREAK RESISTOR HSD40A	2RK05601
15900060		Trapezoid spindle, assembly	
	30	Ball bearing 66201-A-2RSR	2BK00150
	40	Stepper motor, assembl.	15480018
	50	Coupling w.element CPS 15-1 12	2JH70016
	60	Seeger-H-half moon ring H 12	2ZJ30120
	70	Internal circlip J32 DIN472	2ZL20320
15900065	90	O-RING 46.00-2.00 88 NBR 156	2IO20067
		Gas spring, assembly	
15900070	10	Gas spring 350N Lift-O-Mat S80	15900065
		Locking system, assembly	
15900080	10	Thread bracket, locking handle	R5670421
	20	Bearing housing,locking handle	R5110320
	40	Locking handle, assembl.	15900055
	90	O-ring 72 NBR 872. 8.10-1.60, 2 pcs	2IO16010
	100	BOLT M10x80 KV 8.8 FZB DIN931, 4 pcs	2TY11080
15900090		Spindle lock, assembl.	
	10	Cartridge	15900080
15900095	60	Oil ring SA ø10/16x4.5	2IT00010
		Valve housing, assembly	
	10	Valve housing	15900090
15900095	20	Inductive sensor A01G142, 2pcs	2HQ00023
	30	O-ring 72 NBR 872. 25.00-2.00	2IO20055
15900095		Valve knob, assembly	
	10	Valve knob	15900095
	20	O-ring NBR 70. 16.0-1.50, 2pcs	2IO15045





Discotom-60/ -65
Spare Parts and Diagrams

Spare part list for Discotom-60/ -65

Drawing	Pos.	Spare Part	Cat no:
15480018		Stepper motor, assembl	
	10	Stepper motor assembl. w. plug	15483532
	70	SMU PCB and assy. kit	15482905
15730070		Communication plate, assembled	
	10	Plate for connetors	15730070
	20	Chassis Con 4-pol FEMALE IP68	2XM00107
	30	CAP for chassis connector IP68	2XM90000
		Wiring Set for Discotom-60	15903590
	660	Contact block 1 NO 3/4. MTI, 2 pcs	2SB10072
661	Contact block 1 NC 1/2. MTO, 2 pcs	2SB10071	
662	Power Supply 230-500V. 24V/7.5A	2PA90300	
663	Brake circuit, assembly	15900047	
664	Contactora 25A 3SL 0BR 24VDC, 5 pcs	2KM40025	
666	Aux.contact gold 3RH1921-1FE22, 5 pcs	2KH40302	
667	Reverserings kit 3RA1924-2B, 2 pcs	2KH40001	
668	5.00A T Fuse glass 6.3x32 250V	2FU14250	
669	8.00A T FUSE GLASS 6,3x32 250V	2FU14350	
670	Communication plate, assembled	15730070	
671	W9 Joystick Cable A1:X111 - RQ1	15903509	
		Wiring set for Discotom-65	15913590
661	Contact block 1 NO 3/4. MTI, 2 pcs	2SB10072	
662	Contact block 1 NC 1/2. MTO, 2 pcs	2SB10071	
663	Power Supply 230-500V. 24V/7.5A	2PA90300	
664	8.00A T FUSE GLASS 6,3x32 250V	2FU14350	
660	Contactora 32A 3SL 0BR 24VDC	2KM40133	
		Various	
		Pipe for drain	15110997
		Spanner for cutting wheel 24 mm	2GR00224
		Grease gun (GadusS2 V22000)	16080802
		Lubrication oil (Shell Tellus S100)	16080845


Kabel udtag.

50 is assembled in 60 with OmniFit 100 - thread must be **maximum** 3 mm.

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




















Power supply, motor brake, contactors, terminals and fuses*:
See wire set

*Fuses F1 + F2 are variant parts
See Discotom-60, complete


Pos 60: Joint between this part and pos 130 to be lubricated using ALBIDA EP2 grease
Pos 120, 230: To be secured using Omnitilt 50M
Pos 140: To be sealed using Omnitilt 46M
Pos 160: To be secured using Omnitilt 230M
Pos 350: Inside joints to be sealed using Silicone TRV 109 Alu. Joints between this part and pos 10 to be sealed using Elasticon 33N, transparent
Pos 540, 545: To be glued using Loctite 409 Gel
Manufacturing instructions: 15907305 and 15907306


V	2012.09.06	2NS11169-v2NS14139 (2PAR1240) A10 added, see wire set	SPE	2012.09.06	
A	07-11-08		SJE		
Revision	Draw date	Revision description	Draw. Int.	Appr. date	Appr. Init.
			1:5		
			Scale	Format	Tolerance: ISO/ISO 2768 - mK
			A1		Surface treat: None
Description: 15900005 Base and back plates, assembly					Rev:
					V



SECTION D-D
SCALE 1:1





SECTION C-C
SCALE 1:1


SECTION E-E
SCALE 1:1


DETAIL G
SCALE 1:1

DETAIL F
SCALE 1:1

SECTION A-A



SECTION B-B




L = 463mm
To be mounted using Loctite 409

- 10 Inside joints to be sealed using grey Wetroomsilicone 512.
- 20 To be mounted using grey Wetroomsilicone 512.
Gap in top and bottom of pos 20 must be filled with black Elasticon 33N
Curing time 24 H
- 30 Joints to be sealed using grey Wetroomsilicone 512.
- 60 130 To be secured after adjustment using OmniFit 50M

ID	2009.12.08	Pos.160 20Z44208-20Z44206; Pos.180 2ZA20004-2	A1000@E	2009.12.08		
Revision	A	11-05-06		DJE	Drawn by	Appr. by
Revision	Cruc. date	Revision description		Form:	Scale:	Surface treat.
				A1	1:5	None
						Rev: D
Description:						
15900010 Cover, assembly						


SECTION A-A


C	30-07-07	BB: XPE9 -> XPE11		RJE		
A	11-05-06			RJE		
Revision	Crea. date	Revision description		Draw. Init	Appr. date	Appr. Init
	dd-mm-yy				dd-mm-yy	
		Material:	Scale: 1:5	Format: A2	Tolerance: DS/ISO 2768: mK	Surface treat: None
ID:	Description:					Rev:
	15900015 Cabinet, assembly					B

Struers
 Roskildevej 69
 DK-2300 Lyngby
 Denmark
 Phone: +45 44 65 80 00
 Fax: +45 44 65 80 04

DETAIL A
SCALE 2 : 1
General distance
around front plate




SECTION D-D
SCALE 2 : 1
Place black metal frame
behind rubber sleeve


SECTION C-C


Joystick, assembly
See 15903509

* = See wire set 15903590











DETAIL B
SCALE 1 : 1









Pos 76: Remove 2x M5 nuts from joystick before placing this component

J	2012-01-18	(140) Nut for BR1 changed	BRY	2012-01-18	JTV
A	10-03-08		RJE		
Revision	Crea. date	Revision description	Draw. Init	Appr. date	Appr. Init
	yyyy-mm-dd			yyyy-mm-dd	
		Material:	Scale: 1:2	Format: A2	Tolerance: DS/ISO 2768-
ID: 15900020 Front plate, assembly					Rev: J
<small> Høvedvej 14 DK-7270 Ballerup-Centerpark Denmark Phone: +45 44 600 000 Fax: +45 44 600 804 </small>					





Manufacturing instruction: 15907303

J	2010-02-18	to K14. Test about first edition	BRV	2010-02-18	JTV
A	09-05-09				
Revision	Chng. data	Revision description	Draw. Int.	Appr. date	Appr. Init.
			Scale	Format	Draw. Code
			1:5	A1	None
15900025 Base, assembly					File: J


Pos 20: Joint between this part and pos 10 to be lubricated using ALBIDA EP2 grease
 Pos 140: Thread to be lubricated using Anti Size grease
 Pos 160, 200: To be secured using Omnifit 50M
 Pos 190: To be secured using Omnifit 50H
 Manufacturing instructions: 15907302

G	2011-02-18	Pos.170 long again	JTV	2011-02-18	JTV
A	11-04-05		RJE		
Revision	Creo. date	Revision description	Draw. Int	Appr. date	Appr. Int
			Scale: 1:2	Form: A1	Material: 6061-T6
Description: 15900030 Y-table, assembly					
					G

DETAIL B
SCALE 2 : 1

DETAIL C
SCALE 2 : 1

Manufacturing instruction: 15677300



E	2013-08-16	BQ3 note added	JTV	2013-08-16	JTV
A	05-01-06		RJE		
Revision	Crea. date	Revision description	Draw. Init	Appr. date	Appr. Init
	2013-08-16			2013-08-16	
	Material:		Scale:	Format:	Tolerance: DS/ISO 2768- Weight : 6077,076 g
			1:1	A2	mk
ID: Description: 15900040 Spindle housing, assembly			Rev:		
					E

1

2

3

4

C 2012.04.27

Thermal paste added (Kølepasta)

SPE

2012.04.27

A 09-01-07

RJE

Rev

Crea. date
yyyy-mm-dd

Revision description

Draw. Init

Appr. date
yyyy-mm-dd

Appr. Init

Pederstrupvej 84
DK-2750 Ballerup
Copenhagen
Denmark
Phone: +45 44600 800
Fax: +45 44600 804

Material:

Scale:

1:1

Format:

A4

Tolerance: DS/ISO 2768 - mK

Weight : g




ID:

Description:


15900047 Brake circuit, assembly

Rev:

C



SECTION A-A


(20)(30): Omnifit 230L
 (50): Coupling to be fitted on stepper motor axle (40) before fitting on trapezoid spindle (20)

G	2011.10.17	15480015->15990029	SPE	2011.10.17	
A	13-03-06		RJE		
Revision	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
		Material:	Scale: 1:2	Format: A3	Tolerance: DS/ISO 2768- mK Surface treat.: None
ID:	Description: 15900060 Trapezoid spindle, assembly			Rev:	G

Paderbursvej 84
 DK-2750 Ballerup/Copenhagen
 Denmark
 Phone: +45 44 600 800
 Fax: +45 44 600 804

B	2012-02-28	Pos 30 ny model	JTV	2012-02-28	JTV
A	24-10-05		RJE		
Rev	Crea. date yyyy-mm-dd	Revision description	Draw. Init	Appr. date yyyy-mm-dd	Appr. Init
F	 Pederstrupvej 64 DK-2750 Ballerup Copenhagen Denmark Phone: +45 44600 800 Fax: +45 44600 804	Material:	Scale: 1:2	Format: A4	Tolerance: DS/ISO 2768 - mK
					Weight : 67.49 g
ID:		Description: 15900065 Gas spring, assembly			Rev: B

SECTION A-A

- 40 Threaded end of shaft to be greased with ALBIDA EP2
- 60 To be greased with ALBIDA EP2

C					
B					
A	31-10-05			RJE	
Revision	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
		Material:	Scale: 1:2	Format: A3	Tolerance: DS/ISO 2768- mK
<small>Pederstrupvej 84 DK-2750 Ballerup/Copenhagen Denmark Phone: +45 44 600 800 Fax: +45 44 600 804</small>		ID:	Description:		Surface treat.: Non
15900070 Locking system, assembly					Rev: A


1

2

3

4

A



B

C

D

E


SECTION A-A


B	2011-11-15	D,1 D,3 and C,3: Text added	JJO		
A	05-01-06		RJE		
Rev	Crea. date yyyy-mm-dd	Revision description	Draw. Init	Appr. date yyyy-mm-dd	Appr. Init
F	 Pederstrupvej 84 DK-2750 Ballerup Copenhagen Denmark Phone : +45 44600 800 Fax : +45 44600 804	Material:	Scale: 2:1	Format: A4	Tolerance: DS/ISO 2768 - mK
					Weight : 141.01 g
ID:		Description:	Rev:		
		15900080 Spindle lock, assembly	B		

BQ2, length to be shortened to 240 mm

BQ1, length to be shortened to 90 mm

SECTION A-A



Manufacturing instruction: 15907306

C	27-11-2007	Pos.no. changed	SPE	27-11-2007	JTV
B	14-12-06	B2, B3: Designations BQ1 and BQ2 exchanged	RJE		
A	19-06-06		RJE		
Revision	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
		Material:	Scale: 1:1	Format: A3	Tolerance: DS/ISO 2768- mK Surface treat.: None
ID:	Description: 15900090 Valve housing, assembly				Rev: C









Pedersstrupvej 84
DK-2720 Ballerup/Copenhagen
Denmark
Phone: +45 44 600 800
Fax: +45 44 600 804






SECTION A-A

Manufacturing instruction: 15907306

C					
B					
A	19-06-06			RJE	
Rev	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
 Frederstrupvej 84 DK-2750 Ballerup Copenhagen Denmark Phone: +45 44600 800 Fax: +45 44600 804	Material:	Scale: 1:1	Format: A4	Tolerance: DS/ISO 2768 - mK Surface treat.: Non	ID: 15900095 Valve knob, assembly
					Rev: A





P	2012-09-13	2P480300_400 ->2P481240 moved to wire set	SDE	2012-09-13	
A	20-06-06		BJE		
Revision	Chng. date	Revision description	Draw. Init.	Appr. date	Appr. Init.
CD	Description	Scale	Form	Tolerance: DIN/ISO 2768 - msp	Surface treatm: None
		1:5	A1		
15910001 Discotom-65, complete					P

Power supply, contactor, terminals and fuses*:
See wire set


*Fuses F1, F2, F6, F7 and F8 are variants:
See Discotom-65, complete


Pos 60: Joint between this part and pos 130 to be lubricated using ALBIDA EP2 grease
 Pos 120, 230: To be secured using Omnitilt 50M
 Pos 140: To be sealed using Omnitilt 45M
 Pos 160: To be secured using Omnitilt 230L
 Pos 350: Mating joints to be sealed using Silicone RTV 109 Abu. Joints between this part and pos 10 to be sealed using Elasticon 33N, transparent
 Pos 540, 545: To be glued using Loctite 409 Gel
 Manufacturing instructions: 15907306 and 15907306

U	2012.09.06	2NS14109-2NS14139 (ZPA91240) A10 added; see list sup	2012.09.06		
A	19-11-08				
Revision	Draw date	Revision description	Draw Inst	Appr date	Appr Inst
		Scale: 1:5	Form: A1	Distance: 0,0/ISO 2286	Unit: mm
		Description:		None	Rev:
15910005 Base and backplates, assembly					U


DETAIL A
SCALE 2 : 1
General distance
around front plate


SECTION D-D
SCALE 2 : 1
Place black metal frame
behind rubber sleeve

SECTION C-C

- * SB1 110
- 105

BR1

- 115 120 130 140 150

* SB2

- 110

* SB3


- 105

* SB3

- 100


- 60
- 50
- A5
- 30

DETAIL B
SCALE 1 : 1

Pos 76: Remove 2x M5 nuts from joystick before placing this component

* = See wire set 15913590


J	2012-01-18	(140) Nut for BR1 changed	SPE	2009-05-06	
A	10-03-09		RJE		
Revision	Cre. date	Revision description	Draw. Init	Appr. date	Appr. Init
	03.03.09			03.03.09	
		Material:	Scale: 1:2	Format: A2	Tolerance: DS/ISO 2768- mK
ID:		Description:		Surface treat.: None	Rev:
		15910020 Front plate, assembly			J

SEE REQUIREMENTS FOR MAINS CABLE
- note #1

MUST BE EXTERNALLY FUSED
Fmax: See note #2

All wires AWG 12

VOTAGE / FREQ.	note #1 W1 - mains cable	note #2 max. ext. fuse	F1+F2 (fuse size) See note #3	F3 (fuse size)	T2 Number of wires
3 x 200-210V / 50-60Hz	AWG = 12	3 x 40AT	2 x 4AT (CC)	5AT	2
3 x 220-230V / 50Hz	2.5mm²	3 x 40AT	2 x 4AT (aM)	5AT	2
3 x 220-240V / 60Hz	AWG = 12	3 x 40AT	2 x 4AT (CC)	5AT	2
3 x 380-415V / 50Hz	2.5mm²	3 x 40AT	2 x 2AT (aM)	5AT	2
3 x 460-480V 60Hz	AWG = 12	3 x 40AT	2 x 2AT (CC)	5AT	2

Note #3:
F1, F2, F3, F5 are time delay fuses
CC...Class-CC characteristic
aM...aM characteristic

REV E: SPE 2009-05-07
- See page 2
REV F: AKR 2009-11-26
- X1 added
- A10 changed to variant part.
- LED lighting added
- Embedded SMU added
REV G: AKR 2012-08-17
- A10 changed so it's no longer variant part










Rev B: 27/07/07 TKO:Trafa T1 replaced by SMPS A10, correction on diagram, variant table changed, X3 corrected to X1, GNYE-4 changed to GNYE5
Rev C: FTH 13-05-2008
Connections on M1 changed
Rev D: FTH 21-05-2008
New aux.contacts on K1 to K5

Svensk A/S
Pødenhusvej 84
DK-2750 Ballerup
Denmark
Telephone: +45 44 600 800






DISCOTOM-60, MAIN SUPPLY CIRCUIT

Size A2	CAGE Code	DWG NO 15903100	Rev G
Scale	Sheet	1 of 3	

Machine Control Board


W17 12x0,75mm2
W18 12x0,75mm2


ALL WRINGS 0,75mm², EXCEPT OTHERWISE MARKED.
 COLOR CODES: BK = Black, BN = Brown, RD = Red, OG = Orange, YE = Yellow, GN = Green, BU = Blue, VT = Violet, GY = Grey, WH = White, PK = Pink








REV: 0: 2008-08-07 - XPE23 added REV: 1: 2008-11-26 - X114 added - A10 changed to support part - LED lighting added - 24V protected SMU added REV: 2: 2008-12-18 REV: 3: 2009-02-28 - A10 changed to 10A longer version part	Rev B: 27/07/07 TNO/Trafo T1 replaced by SMP5 A10, grounding changed Rev C: FTH 13-05-2008 Connections on M1 changed Rev D: FTH 21-05-2008 New aux.cables on K1 to K5	Drawn AS Checked by DC:2008/06/06 Date: 2008-06-06 Reference: 45144000-000
DISCOTOM-60, MAIN SUPPLY CIRCUIT		
Rev: A1	Drawn Code	DWG NO
Scale		15903100
File: August 11, 2002		2 of 3


Machine Control Board


REV E: SPS 2009.05.07 - 3 pin shield REV F: AFR 2009.11.26 - 3 pin shield - A10 changed to variant part - LED display added - Emulation SW added REV G: AFR 2012.08.17 - A10 changed to 18 pin shield variant part Rev D: FTH 21-05-2008 Connections on M1 changed	Rev B: 270707 TKO Trfo T1 replaced by SMPS A10, X1 corrected to X3 See other pages Rev C: FTH 13-05-2008 Rev D: AFR 2012.08.17 Connections on M1 changed Rev D: FTH 21-05-2008	Simon AS Prepared by Dec 2008/2009 Checked Approved Telephone: +44 (0)20 800
DISCOTOM-60, MAIN SUPPLY CIRCUIT		
New sub.coresets on RT to KG A1	Code Code Scale	DWG NO 15903100
Friday, August 11, 2012	Scale	3 of 3

MICROCONTROLLER

GATE DRIVER


A3
UL
94-V0
PCB No. 12050430-4

Rev. C: 20-04-2007 TK
Correction for layout made, transil added.
Rev. D: 28-06-2007 TK
Heart beat indicator connection changed.


STRUERS A/S
PEDERSTRUPVEJ 84
DK-2750 BALLERUP
DENMARK
PHONE: +45 44 600 800

LIN Lamp driver			
See other pages	Size A3	CAGE Code	DWG NO 15903201
Thursday, August 23, 2007	Scale	TK	Rev F
Sheet		1 of 2	


CURRENT SENSE A


CURRENT SENSE B


LAMP DRIVER A


LAMP DRIVER B

ASTABLE MULTIVIBRATOR A


ASTABLE MULTIVIBRATOR B


Rev. C: 20-04-2007 TK Changed starting of lamp principle, changed frequency, transils added, V3, V4, V7, V8 transistors changed.		STRUERS A/S PEDERSTRUPVEJ 84 DK-2750 BALLERUP DENMARK PHONE: +45 44 600 800	
Rev. E: 02-08-2007 TK: C1,C8 changed to 560pF		LIN Lamp driver	
Rev. F: 21-08-2007 TK: C1,C8 changed to 470pF			
See other pages	Size A3	CAGE Code <Cage Code>	DWG NO
Thursday, August 23, 2007	Scale	TK	Sheet 2 of 2
			15903201 Rev F

TRAFO NO. AA-72329 (200-240V)

200 - 210V

max. 220V


220 - 240V

max. 253V


- COLOR CODES:**
- BK = BLACK
 - BN = BROWN
 - RD = RED
 - OG = ORANGE
 - YE = YELLOW
 - GN = GREEN
 - BU = BLUE
 - VT = VIOLET
 - GY = GREY
 - WH = WHITE

TRAFO NO. AA-72329 (400-460V)


380 - 415V

max.440V


408 - 450V

max.473V

430 - 460V

max. 483V


460 - 480V

max. 506V

* Jumper
Max. 10% overvoltage pr. coupling.








STRUERS A/S PEDERSTRUPVEJ 84 DK-2750 BALLERUP DENMARK PHONE.: +45 44 600 800		WIRING DIAGRAM: WIRING AT MAINS TRANSFORMERS 2MT72329: POWER SUPPLY 200V - 480V 50/60Hz.			
		FILE NAME: 15903470.DSN	Size A4	CAGE CODE	DWG NO 15903470
		SCALE	13.07.2006/KFC	Sheet 1	of 1

PL A10 and/ of anden type, X51&2 bytt. plac.
 Kabel ved A1: X118: W5→W41. 2012-09-12
 BRV/AJ05620
 G. (09-10): Color on wires by X118 changed.
 2010-08-23 BRV/JTV

Mod:	Overlsh:	Udsendt:	Side ang. tol. efter DS/SD 2768
			Line
			Item 180607 BRV
Wiring diagram Discotom-65			15913450 H

Machine Control Board

REV D: CHD (2009/10/24)
- A12 recommended
- Q31 changed
- BK03, BK01, BK02 added
- K1 added
- BK03, BK01, BK02 added
- CHD E3 changed to 200/10
- V11 US, CA updated
- See also changes on other pages
REV E: SFE 2009/05/07
- XPE2 added

REV F: JKR 2009/11/25
- LED type added
- Enclosed SMU added
- A13 added on page 1
- A10 changed to safety part
REV G: JKR 2012/06/17
- A10 changed to fit no longer available

REV D: CHD (2009/10/24)
- A12 recommended
- Q31 changed
- BK03, BK01, BK02 added
- K1 added
- BK03, BK01, BK02 added
- CHD E3 changed to 200/10
- V11 US, CA updated
- See also changes on other pages
REV E: SFE 2009/05/07
- XPE2 added

REV F: JKR 2009/11/25
- LED type added
- Enclosed SMU added
- A13 added on page 1
- A10 changed to safety part
REV G: JKR 2012/06/17
- A10 changed to fit no longer available

Rev B: 27/07/07 TKO A10 added, T1 removed, X15 corrected to X105, grounding changed
See other pages

Rev C: CHD (2008/07/09)
- all GNYE wires changed to AWG10
- see also changes on other pages

Rev A: CHD (2008/07/09)
- see also changes on other pages


ALL WRINGS 0.75mm², EXCEPT OTHERWISE MARKED.
COLOR CODES: BK = Black, BN = Brown, RD = Red, OG = Orange, YE = Yellow, GN = Green, BU = Blue, VT = Violet, GY = Grey, WH = White, PK = Pink


Discotom-65
Photomicro 65
Discotom Ballroom
Chermark
Telephone: +49 44 600 800

Rev	DATE	BY	NO	15913100	G
A1					
Scale					

Friday, August 17, 2012

Sheet 2 of 3


Color codes:
 BK = BLACK
 BN = BROWN
 RD = RED
 OG = ORANGE
 YE = YELLOW
 GN = GREEN
 BL = BLUE
 W = WHITE
 GR = GREY
 PI = PINK
 GNYE = GREEN/YELLOW

G. (G8-10): Color on wires by X118 changed.
 2010-08-23 BRV/JV
 F. (5.5) K3 BK33->BK36; (3,6) New SMU and LED
 light added. 2009.12.04 SPE

		Date: 21/02/07 BRV Issue: 2/168
Wiring diagram Discotom-60		Part no. 15903450 H

Color code:
 BK = BLACK
 BRN = BROWN
 RD = RED
 ORG = ORANGE
 YEL = YELLOW
 BLU = BLUE
 VLT = VIOLET
 GR = GREEN
 GRN = GREY
 WH = WHITE
 PNK = PINK
 GRN/Y = GREEN/YELLOW

5. (08-10): Color on wires by X118 changed.
 2010-08-23 BRV/AV
 F1 (6,5) X3 BK33-2BK36; (3,6) New SIMU and LED
 light added. 2006.12.9 SPK

Metu	Overblik	Uitvoerder	Site eng. tot. efter D5/ISO 2768-
			Rev. 1.00
			11/0007/01
Wiring diagram Discotom-65			15913450 G

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is in conformity with the provisions of the following directives:

Safety of Machinery 2006/42/EC according to the following standard(s):
EN ISO 12100:2010, EN ISO 13850:2006, EN ISO 13849-1:2008/AC:2009, EN ISO 13849-2:2008,
EN 60204:2006.

EMC-Directive 2004/108/EC according to the following standard(s):
EN 61000-6-1:2007, EN 61000-6-3:2007+A1:2011, EN 61000-6-4:2007/A1:2011.

Low Voltage Directive 2006/95/EC according to the following standard(s):
EN 60204:2006, EN 50178:1997.

RoHS 2011/65/UE according to the following standard(s):
EN 50581:2012.

Supplementary Information The equipment complies with the following standards:
NFPA70:2011, NFPA79:2007
FCC 47 CFR Part 15 Class A

The above has been declared according to the global method, module A

Date: 21.01.2013

Christian Skjold Heyde,
Vice President, R & D and Production, Struers A/S

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Account no.: 2191-0126305060

VAT no.: 15683309

Letter of conformity: Safety hood/window for the following equipment

AbraPlan-20	PETG thickness: 8 mm
AbraPol-20	PETG thickness: 8 mm
Discotom-5/-6/-10/-60/-65/-100	PETG thickness: 4 mm
Labotom-3	PETG thickness: 4 mm
Labotom-5	PETG thickness: 3 mm
Labotom-15	PETG thickness: 5 mm
Axitom/-5	PETG thickness: 4 mm
Exotom-100/-150	PETG thickness: 4 mm
Magnutom-400/-500	PETG thickness: Front: 10mm, Side: 8mm

The safety hoods/window primary function is to protect against samples not clamped properly.

The safety hood is made from PETG, which is generally used in all Struers machines as "safety glass".

PETG is resistant to scratches and chemicals and provides the highest degree of safety for the operator.

Commodity Supplier: RIAS A/S, Industrivej 11, 7000 Roskilde, Denmark

Producer: NEOPLAST ApS, Tingbjergvej 4, 4632 Bjæverskov, Denmark

Supplier: Struers A/S, Pederstrupvej 84, 2750 Ballerup, Denmark

This is to certify that the product is produced according to regulations from the Commodity Supplier.

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Yours faithfully

Flemming Perret-Gentil
R&D Manager
R&D Department, Equipment

+45 44 600 913 (direct)

+45 27 130 890 (mobile)

flemming.perret.gentil@struers.dk

www.struers.com

Designed for 300 mm cut-off wheels and equipped with a large cutting table of 435 x 265 mm (wxd)

Versatile automatic and manual cutting machine in one – offering exceptional flexibility

Large display and extremely simple operation by joystick and turn/push knob

Automatic Feed Speed Control for fast cutting of hard materials without risk of damaging wheel or sample

Programmable stop offering three different Stop Modes

Variable rotational speed option for optimum cutting, and cutting of hard and soft materials without changing cut-off wheel

Powerful illumination of cutting chamber for easy clamping and cutting

Corrosion resistant cutting chamber for a long product lifetime

Discotom-60/-65 offers you a spacious cutting table and maximum user-friendliness, allowing for easy clamping and sectioning of your large and irregular workpieces

The Discotom-60/-65 successfully combines the latest cutting technologies with unprecedented flexibility and user-friendliness

The cutting table is provided with 10 mm T-slots allowing various clamping tools to be used

Extended cutting capacity

The Discotom-60/-65 allows for easy clamping and sectioning of all metallographic samples with a diameter of up to 105 mm (4.1"). Using automatic cutting, samples of e.g. 160 mm depth x 65 mm height or 178 mm depth x 35 mm height can be cut. It is possible to cut samples as long as 435 mm in automatic mode. Please refer to the cutting capacity diagram, which shows the capacity for automatic cutting. If manual and automatic cutting are used in combination, the capacity increases significantly.

Thus, when operated automatically as well as manually, Discotom is offering exceptional flexibility.

Save time by Automatic Feed Speed Control

The unique Automatic Feed Speed Control does not just reduce feed speed, but is also capable of testing a preset feed speed after a motor overload, and resume cutting using this speed. In other words, feed reduction is done only as and when absolutely necessary, meaning that no time is wasted.

If, during testing, the feed speed does not cause a second motor overload, cutting will continue at the pre-set feed speed. If, instead, cutting fails a third time, the feed speed will be reduced by 20% at which

cutting will resume. This process will repeat itself until the feed speed is reduced to 10% of the pre-set value.

The optimum feed speed

When the operator has found the optimum feed speed for a particular workpiece he or she can cut similar workpieces, without any problems. The process can be tracked on the Motor Load and Feed Speed bars on the display.

Programmable stop

Discotom-60/-65 features three different Stop Modes for fast and easy cutting: AutoStop, Absolute Stop Position and Relative Stop Position.

AutoStop

Cutting stops when the sample has been cut through.

Absolute Stop Position

Cutting stops at a pre-set position counted from zero position where the table is at the very front of the cutting chamber. The positioning range is 0-200 mm.


Relative Stop Position


Cutting stops at a position relative to the position where cutting starts. Thus, after entering the approximate sample size, the cutting process will stop immediately after the sample has been cut through. The positioning range is 0-200 mm.

Maximum user-friendliness

Discotom-60/-65 is equipped with a large display and a joystick-controlled cutting-table. Cutting parameters and cutting mode are fast and easily selected via the turn/push-knob.

The Struers Cooling System extends the lifetime of the cooling water, and effectively removes swarf (option).

Discotom-60/-65 is equipped with a large easy-to-read display and user-friendly turn/push knob and joy-stick.

Optional table unit

Parallel cutting table (option) for cutting of multiple slices.

Graphic display

The large LCD display permits the graphic display of a lot of information in an easy-to-read fashion:

Cutting mode

The top bar lets the user select between Automatic Cutting (with moving cutting-table) and Manual Cutting (manually moving the cut-off wheel into the sample).

Feed Speed bar

The Feed Speed bar presets the Feed Speed. The Automatic Feed Speed Control may interfere on this preset speed as and when necessary.

Rotational speed (rpm)

On Discotom-65 the user can select a spindle speed between 1000 and 3000 rpm. By changing the spindle speed, the properties of the cut-off wheel change. In effect, one cut-off wheel can be used for a wide range of materials, wasting no time changing wheels.

Ordinary SiC and Al₂O₃ wheels get 'harder' with higher spindle speeds, and vice versa. The benefit of the higher speed is that the wheel wears at a slower rate. At the same time, the workpiece becomes more prone to heat damage that is countered by reducing feed speed. Using a slower spindle speed, the same wheel becomes softer, and capable of cutting through harder materials, with less risk of workpiece damage. The downside is increased wheel wear.

Thus, by adjusting the spindle speed, the user is in control of the cutting properties of the wheel as well as the wheel wear.

Stop position

Pop-up menu for choice of 3 different ways of programming the stop. Please refer to separate section, called Programmable Stop.

Motor Load bar

The load bar indicates the load on the cutting motor. The load bar range is 0-200%.

Effective cooling during cutting

Most Discotom-60/-65's are configured with a Cooling System 5, consisting of a 100 litre tank, a static filter and a small pump. This recirculation unit will fit into the practical table unit, 05116916, (specifications overleaf).

For intensive use, and for materials generating a lot of swarf, we recommend using the Coolimat-200. This is a larger, external unit comprising a 200 litre tank and a band filter. This machine offers the choice of small and large pump and the possibility to connect 2 pumps in order to service 2 machines.

Built-in flush hose for cleaning

Discotom comes with an external flushing hose. The waterflow is adjusted directly on the watergun.

High level of safety

The Discotom-60/-65 has an electric, dynamic motor which allows to stop the cutting motor within a few seconds. The safety lock, which ensures that the operator cannot be injured by a rotating cut-off wheel, will open for access to the cutting chamber as soon as a sensor detects that the cut-off wheel is not rotating. During cutting, the protection guard is locked. If the protection guard is forcibly opened, the cut-off wheel will immediately stop spinning.

For cutting of very long workpieces which cannot fit into the cutting chamber, an opening is made in the left side wall, which is accessed via a 'tunnel' on the outside. This effectively prevents anybody from being harmed on the cut-off wheel. The tunnel is offered as an option.

A long product lifetime

The top of the cutting table is fitted with stainless steel bands which are very resistant to corrosion. If damaged, the steel bands are easily exchanged.

Struers cut-off wheels are especially designed for Struers machines, taking into consideration the most recent developments in wet cutting techniques.

Rotational Speed Settings using 40A30 Cut-off wheel

Metalog Guide Description	Tested Material	Hardness HV	Speed RPM
Non-ferrous soft metals	Aluminium	50-110	3000
Very ductile metals	Stainless steel	220	2200
Medium soft ferrous metals	UHB IMPAX Cold work tool steel	300	2200
Medium hard ferrous metals	UHB IMPAX Cold work tool steel	480	2000
Hard ferrous metals	UHB ARNE Cold work tool steel	750	1750

SPECIFICATIONS

Discotom-60/-65

Automatic cutting machines. Automatic feed with electronic control of feed speed.
Manual operation possible.
For 300 mm (12") cut-off wheels.
Cutting table with 10 mm T-slots. 4.0 kW (5.4 HP) motor.
Recirculation cooling unit, clamping devices and table unit are ordered separately.

Discotom-60

Cat.no:	Voltage:
05906129	3 x 200-210 V / 50-60 Hz
05906135	3 x 220-230 V / 50 Hz
05906136	3 x 220-240 V / 60 Hz
05906146	3 x 380-415 V / 50 Hz
05906154	3 x 460-480 V / 60 Hz

Discotom-65

Cat.no:	Voltage:
05916135	3 x 220-230 V / 50 Hz, with EU cable
05916136	3 x 200-240 V / 50-60 Hz, with UL cable
05916146	3 x 380-415 V / 50-60 Hz, with EU cable
05916154	3 x 460-480 V / 60 Hz, with UL cable

Accessories

Recirculation System-5

With 100 l water tank (05766905), Small pump (05766xxx), 100 l static filter* (05766907) and Cooli-1* control unit (05761116)

Cat.no:	Voltage:
05766816	1 x 100 V / 50 Hz
05766822	1 x 220 V / 50 Hz
05766823	1 x 120 V / 60 Hz
05766824	1 x 240 V / 60 Hz
05766916	1 x 100-120 V / 50/60 Hz CSA

*) For intensive use, Coolimat-200 (061611xx) is recommended.

Table unit

For Labotom, Discotom and Unitom.
With compartment for recirculation cooling unit.
Width: 900 mm, height: 800 mm, depth: 750 mm.
Cat.no: 05116916

Parallel Cutting Table

for 10 mm T-slots
For cutting uniform slices without repositioning and re-clamping the workpiece.
140 x 265 mm cutting table.
Positioning range: 60 mm. With detachable handle.
Table height 58 mm.
Cat.no: 05876919

Lifting Brackets

For lifting of Discotom. Set of 4 brackets.
No tools required.
Cat.no: 05906902

T-slot cleaner

T-slot cleaner for 10 and 12 mm T-slots
Cat.no: 05486910

Quick-clamping tool for 10 mm T-slots, Left

For securing the workpiece. Complete with back stop.
To be mounted on the left-hand side of the cut-off wheel.
Cat.no: 05876909

Quick-clamping tool for 10 mm T-slots, Right

For securing the work piece. Complete with back stop.
To be mounted on the right-hand side of the cut-off wheel.
Cat.no: 05876910

Replacement Stainless Steel Band

For cutting table on Discotom-60/-65.
1 pc. of 40 mm x 265 mm.
Cat.no: 05906901

TECHNICAL DATA

Subject	Specifications
Cut-off wheel size	300 mm / 12"
Cutting motor	Constant cutting power (S1): 4.0 kW / 5.4 HP Intermittent cutting power (S3): 4.7 kW / 6.3 HP
Rotational speed	Discotom-60: 2860 rpm at 50 Hz, 3415 rpm at 60 Hz Discotom-65: 1000-3000 rpm at 50/60 Hz
Positioning and Feed	Positioning range (of cut-off wheel) 0 – 120 mm 0 – 4.72" Positioning range (of cutting table) 0 – 200 mm 0 – 7.87" Max. positioning speed 20 mm/s 0.79"/s Feed Speed range 0.01 – 2.5 mm/s 0.001 – 0.1"/s
Cutting Table w x d Cutting table movement	435 mm x 265 mm (17.1" x 10.4") 200 mm (7.8")
Cutting capacity	Ø 105 mm (4.1") 160 mm at 65 mm height 120 mm at 90 mm height
Recirculation cooling unit, Cooling System 5, 05766xxx, (option) Capacity of:	100 l (26.4 gallons) tank volume Large pump: 125 l/min (33 gal/min) flow at 1 bar Small pump: 60 l/min (15.9 gal/min) flow at 1 bar
Dimensions and weight	Height (closed/open hood) 570/890 mm / 22.44/35.04" Width 910 mm / 35.8" Depth (closed/open hood) 785/816 mm / 30.9/32.12" Weight 171 kg (377 lbs)
Table unit, 05116916 (option)	Dimensions w x h x d 900 mm x 800 mm x 750 mm 35.4" x 29.5" x 31.5"
EU-directives	98/37/EEC - Safety of Machinery 73/23/EEC - Low Voltage Directive 89/336EEC, 92/31/EEC, 93/68/EEC - EMC Directive
Dynamic Balancing	Max. admissible unbalance Upper according to ISO 1940/1, Balance Quality Grade G6.03: Upper = 2 gmm/kg = 110 gmm.
Noise Level	Approx. 67 dB(A) during cutting, at a distance of 1.0 m/39.4" from the machine.

Struers' equipment is in conformity with the provisions of the applicable International Directives and their appurtenant Standards. (Please contact your local supplier for details)

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup, Denmark
Phone +45 44 600 800
Fax +45 44 600 801
struers@struers.dk
www.struers.com

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CANADA

Struers Ltd.
7275 West Credit Avenue
Mississauga, Ontario L5N 5M9
Phone +1 905-814-8855
Fax +1 905-814-1440
info@struers.com

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Struers Inc.