# Exotom-100

## Instruction Manual

Manual No.: 15047001

Date of Release 30.07.2003

*Exotom-100*  
*Instruction Manual*

<b>Table of Contents</b>	<b>Page</b>
User's Guide .....	1
Reference Guide.....	19
Quick Reference Card .....	48

---

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the *Serial No.* and *Voltage* on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

**Instruction Manuals:** Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

**Service Manuals:** Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2003.

**Struers A/S**  
Pederstrupvej 84  
DK-2750 Ballerup  
Denmark  
Telephone +45 44 600 800  
Fax +45 44 600 801

---

## **Exotom-100**

### **Safety Precaution Sheet**

#### **To be read carefully before use**

1. The operator(s) should be fully instructed in the use of the machine and its cut-off wheels according to the Instruction Manual and the instructions on the cut-off wheels.
2. Use the supplied Allen key to remove the two red transport screws at the rear of the machine. Do not attempt to open the protection guard before the screws have been removed.
3. The machine must be placed on a safe and stable support. All safety functions of the machine must be in working order. The machine must be levelled by means of the adjustable legs provided.
4. The unit must be installed in compliance with local safety regulations.
5. Prior to lifting the unit by the built-in truck lifting point, ensure that the boom is properly secured with the locking pins provided. Before transport, secure the cutting arm with the locking system provided.
6. To achieve maximum safety and lifetime of the machine, use only original Struers consumables.
7. Use only intact cut-off wheels. The cut-off wheels must be approved for a spindle speed of min. 1950 rpm / 42 m/s. Do not use sawblades.
8. Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
9. The workpiece must be securely clamped in the quick-clamping device or the like. Large or sharp workpieces must be handled in a safe way.
10. Struers recommend the use of an exhaust system as the cutting materials may emit harmful gasses or dust. See the instructions of the cut-off wheel.
11. The machine emits only moderate noise. However, the cutting process itself may emit noise, depending on the nature of the workpiece. In such cases, use of hearing protection is recommended.
12. The machine must be disconnected from the mains prior to any service.
13. Ensure that the cut-off wheel is secured before working on or around the cutting table.
14. Do not put hands through the rubber curtain during operation of the machine.

*Exotom-100*  
*Instruction Manual*

- 15.** Protruding workpieces should be shielded or marked.
- 16.** Use of working gloves is recommended as workpieces may be both very hot and produce sharp edges.
- 17.** If any unusual noise is heard when the protection hood is operated, refrain from further use of the machine, and contact Struers service technician.

---

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

---

# User's Guide

Table of Contents	Page
<b>1. Getting Started</b>	
Checking the Contents of Packing.....	3
Recirculation Cooling Unit.....	3
Unpacking and Placing Exotom-100.....	3
Mounting the Control Panel.....	4
Getting Acquainted with Exotom-100 .....	5
Side view, left.....	6
Side view, right.....	6
Cooling unit compartment .....	7
Power Supply .....	8
Mounting the Cut-off wheel .....	8
Direction of the Cut-off Wheel .....	8
Connection to an External Exhaust System .....	8
Setting Up the Recirculation Cooling Unit .....	9
Direct Water Supply for the Cooling Unit.....	9
<b>2. Operation</b>	
Using the Controls.....	10
Control Panel of Exotom-100 .....	10
Control Panel Functions.....	11
Display Types.....	11
Setting the Language .....	12
Reading the Cutting Display.....	13
Changing Cutting Mode and Cutting Parameters.....	14
Changing Cutting Mode .....	14
Changing Cutting Parameters.....	14
Reading the Motor Information .....	14
Sleep Mode .....	14
Changing the Cut-off Wheel.....	15
Clamping the Workpiece .....	15
Positioning the Cut-off Wheel .....	15
QuickPosition .....	16
Cutting Parameters .....	16
Feed .....	16
Force .....	16
Stop.....	16
AutoStop.....	16

*Exotom-100*  
*Instruction Manual*

Setting the Cutting Parameters .....	17
Cutting on Exotom-100.....	17
Starting the Cutting.....	17
Fast advance .....	17
Stopping the Cutting (Manual Stop) .....	17
Re-starting cutting .....	17
Direct Cut.....	18
ExciCut.....	18
AxioCut (option).....	18
Additional cooling.....	18

## 1. Getting Started

### Checking the Contents of Packing

The packing box contains the following items:

- 1 Exotom-100
- 1 Control Panel for Exotom-100
- 1 Allen key (5 mm) for mounting Control Panel
- 3 Screws (5 mm)
- 1 Fork spanner (30 mm) for cut-off wheel
- 1 Spanner (comb. 12/17 mm) for adjustment of the Control Panel
- 1 Drain hose, 2 m
- 1 Fitting for drain hose
- 2 Cabinet doors
- 1 Set of Instruction Manuals


### *Recirculation Cooling Unit*

- 1 Recirculation cooling unit
- 1 Trolley for cooling unit


### Unpacking and Placing Exotom-100

- Unscrew the nuts from the four transport bolts fixing the machine to the pallet.
- Lift the machine from the pallet by means of a forklift truck from the front, and place in a suitable location.
- Remove the safety-springs from the front crossbar, and remove bar.
- Use the supplied Allen key to remove the two red transport screws at the rear of the machine. Do not attempt to open the protection hood before the screws have been removed.
- Take out the loose parts (trolley, tank, drain hose etc.).
- Mount the cabinet doors.
- Turn the adjustable feet so that the machine stands firmly and is level.

**NB!**  
Please remove the transportation support before use.


### Mounting the Control Panel

- Unpack the Control Panel and mount on right side of unit using the two 5 mm screws and Allen key supplied (use holes ① and ② as indicated on drawing).
  - Mount the earth wire below the mounting plate using the third 5 mm screw ④.
  - Press cover towards cabinet wall until it snaps onto the mounting plate.
  - Connect Control Panel plug ⑥ to socket ⑤.
- The position of the Control Panel is controlled by a friction mounting. After mounting the control panel it is possible to adjust the mobility of the Control Panel this way:
- Adjust the main joint with the nut ③ using the combined 12/17 mm spanner supplied.
  - The joint just behind the Control Panel may be adjusted too, using the same spanner.

**Getting Acquainted with  
Exotom-100**


Take a moment to familiarise yourself with the location and names of the Exotom-100 components.

- ① Cutting chamber
- ② Control Panel
- ③ Flexible water jets
- ④ Protection guards
- ⑤ Emergency stops
- ⑥ Flushing hose with flush nozzle
- ⑦ Cut-off wheel
- ⑧ Spindle lock button
- ⑨ Guard for cut-off wheel
- ⑩ Protection hood


*Exotom-100*  
*Instruction Manual*

*Side view, left*


- ① Exhaust flange
- ② Water inlet
- ③ Water outlet

*Side view, right*

- ④ Name Plate
- ⑤ Type plate
- ⑥ Main Switch
- ⑦ Electrical connection box
- ⑧ Connection for electrical cable for power supply
- ⑨ Connection for external warning light
- ⑩ Fuse socket for external warning light
- ⑪ Connection for Control Panel

*Cooling unit compartment*

- ① Cooling unit drain pump on/off
- ② Filter drawer
- ③ Recirculation cooling unit
- ④ Trolley for cooling unit

## Power Supply

Before connecting the machine, check that the mains voltage is correct by referring to the type plate.


Open the electric connection box and connect a 4-lead cable in the following way:

PE: earth  
L1: phase  
L2: phase  
L3: phase

### **IMPORTANT**

Check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.

## Mounting the Cut-off wheel

- ① Nut
- ② Flange
- ③ Washer

- Press the knob for the spindle lock on the left-hand side of the cut-off wheel, while turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner (30mm).
- Remove the flange and the cut-off wheel.
- Mount the new cut-off wheel with a cardboard flange on each side.
- Mount the flanges and nut. Tighten carefully.

## Direction of the Cut-off Wheel

To check that the cut-off wheel rotates in the direction indicated on the wheel-guard, do the following:

- Close the protection hood.
- Turn Main Switch to On.
- Press START  $\blacklozenge$ .
- Check the rotational direction of the wheel.
- Press STOP  $\circledcirc$ .

If the direction of rotation is incorrect, turn mains off and switch two of the phases.

## Connection to an External Exhaust System

Struers recommends the use of an exhaust system as workpieces may emit harmful gases when cut. The unit is prepared for connection to an exhaust system via an 80 mm fitting on the left hand side of the cabinet.

Recommended capacity for exhaust system: 350m<sup>3</sup>/h at 0mm water gauge.

### **Setting Up the Recirculation Cooling Unit**

- Check that the nylon sieve in the filter drawer is properly placed and close the drawer.
- Mount the drain hose on the back of the machine.
- Place the recirculation cooling unit in the compartment.

***IMPORTANT***

The cooling unit must be positioned inside the compartment as indicated by the arrows on the edge of the cooling unit.

- 
- Fill the cooling unit with 150 l of water (5-10 cm from the upper edge).
- Use additive enclosed with the machine and add to cooling water as specified on the bottle label and stir.

### **Direct Water Supply for the Cooling Unit**

Exotom-100 may be directly connected to the water mains.


- Connect the water inlet fitting ( $\frac{1}{2}$ " internal,  $\frac{3}{4}$ " external thread).
- Open the external water supply.
- Fill the cooling unit with 150 l of water (5-10 cm from the upper edge).
- Shut off the water.
- Add additive in the quantity stated on the bottle label and stir.

***IMPORTANT***

When the cooling unit is pushed into the compartment, care should be taken that the flushing hose is not trapped behind it. To counter this, pull the flushing hose all the way out of its holder, then replace. The flushing hose, should then rest inside the cooling unit.





## 2. Operation

### Using the Controls Control Panel of Exotom-100

## Control Panel Functions

Name	Key	Function	Name	Key	Function
1 POWER	<b>Power</b>	Lights when the Main Switch has been turned to ON.	5 MENU		Push button to toggle between Cutting menu (daily use) and Configuration menu (basic settings).
2 START		Starts the machine and recirculation unit.	6 JOYSTICK		Move up- or downwards to position cut-off wheel.
3 STOP		Stops the machine and recirculation unit.	7 MULTI-FUNCTION KNOB		Push knob to select function. Turn knob to adjust settings.
4 CUTTING MODE		Push button to select desired cutting mode: Direct cutting, ExciCut, AxioCut/Step, AxioCut/Sweep.			

## Display Types

The Cutting Display will be shown on the Control Panel when turning Exotom-100 on. This display is for everyday use. By pressing the MENU  button once, the CONFIGURATION menu will appear on the display. This menu will normally only be accessed during installation.

## Setting the Language

Press MENU  button once to select CONFIGURATION Menu.

Turn knob to toggle between parameters in the CONFIGURATION Menu.


Push knob to select LANGUAGE. A pop-up menu appears.


From the pop-up menu, select your preferred language by pushing the knob.

Press MENU  button to move from CONFIGURATION Menu to Cutting Display.

### Reading the Cutting Display

The Cutting Display offers three types of information:  
Cutting Mode  
Cutting Parameters, and  
Motor Information

- A Cutting Mode.
- B Cutting Parameters.
- C Motor Information.

### Cutting Mode

The top bar displays the selected Cutting Mode:  
Direct Cutting, ExciCut, AxioCut/Step or AxioCut/Sweep  
(AxioCut/Step and -/Sweep are optional).

### Cutting Parameters and Motor Information

The large window on the Cutting Display, displays information about the Cutting Parameters (FEED, FORCE, STOP) as well as Motor information (TEMP, LOAD).

The Cutting Parameters (FEED, FORCE and STOP) can be set both prior to and during cutting. The set value is displayed below the column. Actual value is displayed above the bar graph.

The Motor Information columns TEMP and LOAD, inform about the condition of the motor during cutting. Read-out in relative value (%).


## Changing Cutting Mode and Cutting Parameters

### Changing Cutting Mode


Toggle between the four Cutting Modes, by pressing the CUTTING MODE button. The selected mode, Direct Cutting, ExciCut, AxioCut/Step or AxioCut/Sweep, appears on the top bar of the Cutting Display.

### Changing Cutting Parameters

Toggle between the three Cutting Parameters by pushing the knob. The highlighted Cutting Parameter is selected.

Push knob to select Cutting Parameter.  
Turn knob to adjust setting.

Actual values of the Cutting Parameters FEED, FORCE and STOP are displayed on top of the columns (A).

The actual position of the cut-off wheel (relative its starting position) is displayed graphically by the small icon  to the right of the STOP column.

Turn the knob to change the setting of the selected Cutting Parameter. The arrow on the left of the column will move to reflect the new setting (B).

By moving the pointer of the STOP column below the column, AUTO mode is selected.

### Reading the Motor Information

The highlighted area of the bar graphs TEMP and LOAD, shows the status of the cutting motor:

**TEMP.** Temperature indicator of the cutting motor.

**LOAD.** Load indicator of the cutting motor.

### Sleep Mode

To increase the lifetime of the display, the backlight is automatically switched off if Exotom-100 has not been used for 15 min. Push any key on the Control Panel to re-activate the backlight.

## Changing the Cut-off Wheel

- Press the knob for the spindle lock on the left-hand side of the cut-off wheel, while turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner (30 mm).
- Remove the flange and the cut-off wheel.
- Mount the new cut-off wheel with a cardboard washer on each side.
- Mount the flange and nut. Tighten carefully and close the guard.

### **IMPORTANT**

Conventional cut-off wheels based on  $Al_2O_3/SiC$  abrasives should be placed between two cardboard washers, to protect the cut-off wheel.  
For maximum precision with diamond or CBN cut-off wheels, do not use cardboard washers.

## Clamping the Workpiece

- Clamp the workpiece with the clamping device of your choice e.g. a quick clamping device. Place the workpiece between the clamp and the back stop.
- Push the clamp towards the workpiece, and lock the quick clamping device with the locking handle.

Ensure that only one of the quick clamping devices is tight, the other device should only press lightly. Use support tools if the geometry of the workpiece makes support necessary.

### **IMPORTANT**

When using the AxioCut cutting mode, the cut-off wheel cover may hit the jaw of the Quick Clamping Device, if same is mounted in the shown forward position; ie using the cross-directional T-slot closest to the front of the machine.

To avoid this situation, place the quick-clamping device as far back as possible, using the front-most T-slot.

Test possible infringement by starting the machine with the cut-off wheel completely clear of the clamping device.

Should the cut-off wheel hit the clamping device anyway, no danger is involved. Exotom-100 will stop automatically and display the following message: "AxioCut motor blocked".

## Positioning the Cut-off Wheel


After clamping, test the correct position of the cut on the workpiece, by lowering the cut-off wheel into position, min. 1 mm above the workpiece:

- Push the joystick on the Control Panel downwards, to lower the cut-off wheel. The cut-off wheel stops when the joystick is released.

**QuickPosition**

Positioning of the cut-off wheel may be done, by simply pushing the joystick downwards until the cut-off wheel has contact with the workpiece. The force is automatically reduced on contact, preventing damage to the cut-off wheel and the workpiece. After contact with the workpiece, the cut-off wheel is automatically retracted 2 mm, to be ready for cutting.

**Cutting Parameters**

**Feed**

The feed speed can be set to values between 0.05-5.00 mm/sec (0.002-0.2 "/s).

**Force**

The maximum permitted force between cut-off wheel and workpiece can be set to values between 20-700 N (4-150 lbs). A built-in measuring cell constantly computes the force. If the force limit is reached, the feed speed will automatically be reduced to a value that allows the force to stay just below the set limit. As soon as the force drops below the set limit, the speed will be increased to the original setting.

**NOTE**

When cutting with a force below 50 N, beware of hysteresis (actual force applied may differ from selected cutting force).

**Stop**

There are two ways to set the stop position: AutoStop and Fixed Stop, explained under Advanced Operations.

*AutoStop*

When the AutoStop function is selected, the machine automatically stops when the workpiece has been cut through. For normal cutting, the AutoStop function is recommended.

## Setting the Cutting Parameters

The Cutting Parameters, FEED, FORCE and STOP, are set in the Cutting Menu.

If in CONFIGURATION Menu, press MENU  once to select Cutting Menu.

Push knob to toggle between parameters in the Cutting Menu.

Turn knob to edit the value of the selected parameter.


Push knob to move to next parameter.

## Cutting on Exotom-100

### *Starting the Cutting*

- Position the cut-off wheel.
- Close the protection hood.
- Press START . The cut-off wheel starts rotating, the cooling water starts running and the cut-off wheel will slowly move down into the workpiece at the pre-set Feed Rate.

### *Fast advance*

After START  has been pressed, the cut-off wheel may also be advanced towards the workpiece manually. To advance the cut-off wheel, push the joystick downwards. The cut-off wheel will advance towards the workpiece at maximum speed (5mm/s).

### *Stopping the Cutting (Manual Stop)*

Exotom-100 automatically stops cutting when the workpiece is cut through.


- Press STOP  to manually interrupt the cutting process. The machine can be stopped at any time during operation by pressing STOP .

### *Re-starting cutting*

Having interrupted cutting, the cut-off wheel has been automatically moved out of the workpiece.


- Press START  to resume cutting.

### Direct Cut

Direct Cutting is the normal cutting mode. The cut-off wheel is moved into the workpiece in a slightly curved, vertical movement. Direct Cutting is intended for ordinary materials.

### ExciCut

ExciCut cutting mode is ideal for cutting of very hard materials (HV >400). The oscillating movement of the cut-off wheel has three main advantages: less wear on the cut-off wheel, less risk of damage to workpiece and less risk of motor overheating. To select ExciCut, press the CUTTING MODE button, until ExciCut appears in the top bar of the display.


### AxioCut (option)

The AxioCut option permits cutting of extra large workpieces (adding 150mm to max. depth). Two different cutting modes are provided, Step and Sweep. The Sweep-mode offers better cut-off wheel economy than the Step-mode. The latter method, however, is faster.

To select AxioCut/Step or AxioCut/Sweep, press the CUTTING MODE button, until the desired mode appears in the top bar of the display.

### Additional cooling

Two flexible water jets are provided for additional cooling. Cooling water from the integrated water jets will flow away from the upper part of the cutting area, particularly when cutting hollow workpieces. To counter this, use the flexible water jets provided.

- Position the cooling jets to the left and right of the cutting area.
- Turn the valve on the jet to a position parallel to the hose to activate the cooling jet. The cooling water will flow as soon as the cutting starts.

**Note**

When using the flexible water jets, cooling water is diverted from the integrated water jets positioned over the cut-off wheel.

# Reference Guide

Table of Contents	Page
<b>1. Advanced Operations</b>	
Configuration Menu .....	21
Display contrast.....	22
Language.....	22
Return position .....	22
Units.....	22
Operation mode .....	22
New Pass Code.....	23
Changing Operation Mode .....	24
Cutting Offset .....	25
Cutting Display.....	26
Cutting Modes.....	27
Direct Cut.....	27
ExciCut.....	27
AxioCut (option).....	27
AxioCut/Step.....	28
AxioCut/Sweep.....	28
Stop Settings .....	29
AutoStop.....	29
Fixed Stop .....	29
QuickPosition.....	30
OptiFeed.....	30
Cutting Long Workpieces .....	31
Clamping Irregular Workpieces.....	31
Connection to External Exhaust System.....	31
Other Safety Features .....	32
Optimising Cutting Results .....	33
<b>2. Accessories</b> .....	34
<b>3. Consumables</b>	
Cut-off Wheels.....	35
Other Consumables .....	35

## **4. Trouble-Shooting**

Error Messages .....	39
Messages .....	39
Errors .....	39
Fatal Errors .....	39

## **5. Maintenance**


Daily Service .....	42
Cleaning the Cutting Chamber .....	42
Weekly Service .....	42
Monthly Service .....	42
Replacing the Cooling Water .....	42
Maintenance of Recirculation Unit .....	43
Low Level Warning .....	43
Refilling the Cooling Unit .....	43
Emptying the Cooling Unit .....	43
Filter Drawer .....	44
Additive for Cooling Water .....	44
Maintenance of Cutting Table .....	45
Maintenance of Cut-off Wheels .....	45
Storage of Bakelite Bonded Al <sub>2</sub> O <sub>3</sub> Cut-off Wheels .....	45
Maintenance of Diamond and CBN Cut-off Wheels .....	45
Maintenance of Clamping Devices .....	45

## **6. Technical Data** .....

Cutting Capacity .....	47
------------------------	----

## 1. Advanced Operations

### Configuration Menu


Press MENU  button once to select CONFIGURATION Menu.

Turn knob to highlight different parameters in the CONFIGURATION Menu.

Push knob to edit the chosen parameter.

Press MENU  button to move from CONFIGURATION Menu to Cutting Display.

*Display contrast*

The contrast settings of the display can be adjusted to suit individual preferences (default value: 22, adjustment interval: 0-50).


*Language*

The language can be set to English (default), German, French, Spanish or Japanese.

*Return position*

After cutting or after pressing STOP , the return movement of the cut-off wheel can be set to three different functions:

*Top:* Exotom-100 automatically retracts the cut-off wheel to the top position.

*Start:* Exotom-100 automatically retracts the cut-off wheel to the original position of the cut-off wheel, at the time you pressed START  (default).

*Stay:* The cut-off wheel stays down.

***Important***

Use the Stay function for bakelite bonded diamond or CBN cut-off wheels, as retraction might destroy the rim of the cut-off wheel.

*Units*

The Feed, Force and Stop values in the display panel can be set to displayed in either mm/Newton (default) or inches/pounds force.

*Operation mode*

It is possible to select three different operation modes:

*Configuration:* Full functionality

*Development:* No access to parameters in CONFIGURATION menu, except Display contrast


*Production:* Access to START, STOP, Fixed Stop and movement of cut-off wheel, and to Display contrast in the CONFIGURATION menu

To select your own pass code, go to the CONFIGURATION menu. Select Operation mode to get access to CHANGE OPERATION MODE menu.

## New Pass Code

To select your own pass code, go to the CONFIGURATION menu. Select Operation mode to get access to CHANGE OPERATION MODE menu.

 Push knob to select Pass code


 Turn knob until default pass code '176' appears.

 Push knob to select default pass code, and a new line (New pass code) appears in the CHANGE OPERATION MODE menu.


 Turn knob to move the cursor to New Pass Code.

 Push knob to select and turn knob to set your own three-digit pass code.


 Press Menu button to confirm your unique pass code and exit the CHANGE OPERATION MODE menu.

## Changing Operation Mode

To change the operation mode, highlight the parameter in the CONFIGURATION menu.

 Push knob to move to Change Operation Mode menu.

 Enter your unique pass code by turning the knob.  
Confirm by pushing knob.

 Push knob  
A pop-up will appear.


 Select desired operation mode and push knob to confirm.

Press Cutting Mode button to move to CONFIGURATION menu.

**Cutting Offset**

**Note:**  
 Cutting Offset requires installation of the AxioCut option.


The horizontal position (cutting offset) of the cut-off wheel can be set manually, prior to the cutting. In the ExciCut and Direct Cut cutting modes, the cut-off wheel may be moved 10 mm towards the front and 60 mm towards the back of the cutting chamber ('+10/-60', default value is 0). This allows for greater flexibility in terms of cutting difficult work pieces and clamping.


To change the position, select the parameter by pushing the knob. Turn knob to change the value. Push knob to confirm the desired value.

Once the new value has been confirmed, the cut-off wheel moves to the desired position - even if the protection hood is open.

If the cutting offset has been changed (i.e. not the default value of 0), the Cutting Display heading will remind the operator which offset has been selected.

## Cutting Display

The Cutting Display appears on the Control Panel when Exotom-100 is turned on. The Cutting Display is designed as a single-layer user interface, i.e. without sub-menus. Whereas this menu is for everyday use, the CONFIGURATION menu (as described above) will normally only be used at installation.

If in CONFIGURATION menu, press MENU  to select Cutting Display.

Push knob to toggle between parameters in the Cutting Display. (FEED, FORCE, STOP).

Turn knob to edit the value of the selected parameter.

Push knob to move to the next parameter.

## **Cutting Modes**

The Cutting menu has four cutting modes:

Direct Cut


ExciCut

AxioCut/Step (optional)

AxioCut/Sweep (optional)


To select Cutting Mode, press the CUTTING MODE button, until the desired mode appears in the top bar of the Cutting menu. The latest Cutting Mode used is saved, and will appear in the Cutting Display the next time Exotom-100 is switched on.

### **Direct Cut**

Direct Cut is the normal cutting mode. The cut-off wheel is moved into the workpiece in a slightly curved, vertical movement, without any separate movement on the horizontal axis. Direct Cutting is intended for ordinary materials.

### **ExciCut**


ExciCut cutting mode is ideal for fast cutting of very hard materials (HV >400). The oscillating movement of the cut-off wheel has three main advantages: less wear on cut-off wheel, less risk of damage to workpiece and less risk of motor overheating. To select ExciCut, press the CUTTING MODE button, until ExciCut appears in the top bar of the display.

### **AxioCut (option)**


The AxioCut option permits cutting of extra large workpieces (adding 150 mm to max depth). Two different modes, Step and Sweep, are provided; Step cutting offers faster cutting, whereas Sweep cutting offers better wheel economy. Both methods are based on the Direct Cutting method, i.e. without oscillation of the cut-off wheel.

*AxioCut/Step*

In Step-mode, the cut-off wheel enters the workpiece in three alternating, pre-programmed steps of 10mm. This method offers fast cutting of even very hard materials. To select AxioCut/Step, press the CUTTING MODE button, until the desired mode appears in the top bar of the display. Note that the two initial steps of a cycle, 1 & 2, are 5mm only. The cutting depth on steps 3 to 5 is 10mm. After completing step 5, steps 3 to step 5 are repeated until the workpiece has been cut through. The arrows on the drawing point towards the position of the cut-off wheel *after* cutting.

*AxioCut/Sweep*

Sweep-mode offers better cut-off wheel economy than Step-mode, as a harder cut-off wheel can be used which in turn is more resistant to wear. The latter method, however, is faster.

AxioCut/Step is intended for cutting of large items. Attempting to cut small items in this mode may prove inefficient. Select alternative Cutting Mode.

In Sweep-mode, the cut-off wheel enters the workpiece at a pre-programmed feed speed. A Step column replaces the far-left column in the Cutting Display (Feed) when in Sweep-mode. FORCE is automatically set to maximum (700 N), to ensure utilisation of all available force, no matter which STEP size has been selected.

To select AxioCut/Sweep:

- Press the CUTTING MODE button, until the desired mode appears in the top bar of the display.
- Push knob to select STEP parameter.
- Turn knob to set STEP to the desired value (0.1-4.0 mm in steps of 0.1 mm).
- Push knob to select STOP parameter.
- Turn knob to set the desired STOP position.
- Press START ◊.

Note that the first step of a cycle moves the cut-off wheel from the back to the front of the workpiece without touching it.

The arrows on the drawing point towards the position of the cut-off wheel *after* cutting/re-positioning. After completing step 6, steps 1 to 6 are repeated until the workpiece has been cut through.

## Stop Settings

### AutoStop

There are two ways to set the stop position: AutoStop or Fixed Stop.

When the AutoStop function is selected, the machine automatically stops when the workpiece has been cut through. For normal cutting, the AutoStop function is recommended. In the Cutting Menu, toggle between FORCE/FEED/STOP by pushing the knob, until STOP is highlighted. Turn knob to move pointer to setting AUTO at the bottom of the STOP column.

When the pointer in the STOP column is not on AUTO, Exotom-100 will only stop when it reaches the preset stop position or the STOP  key is pressed.

Note: AutoStop acts on sudden load changes. In some cases, eg using very low feed rate and/or low force, subtle load changes may not be detected by the machine, preventing AutoStop from working properly. This may particularly be the case when cutting soft materials, pipes or workpieces with changing cross-section. If you are facing problems repeatedly, the sensitivity of the AutoStop-function may be adjusted to fit your specific application. Please contact your Struers service engineer for details.

### Fixed Stop


Fixed stop is used when a specific stop position is desired. When cutting pipes or other workpieces with changing cross sections, the cut-off wheel may retract before the workpiece has been cut. To overcome this, use the Fixed Stop function. Clamp the workpiece and position the cut-off wheel just above the workpiece.

This position is automatically set to 0 (zero) and indicated by the cut-off wheel icon  to the right of the Stop column. Accordingly, as soon as Start has been pressed, the actual position of the cut-off wheel, becomes a relative starting point (zero), from where the cutting depth is calculated.

Select the STOP parameter and set the desired stop position by turning the knob.

Exotom-100 will now stop when it reaches the pre-set stop position.

Remember to take the wear of the cut-off wheel into consideration.


*QuickPosition*

Positioning of the cut-off wheel may be done automatically, by simply pushing the joystick downwards until the cut-off wheel is in contact with the workpiece. The force is automatically reduced on contact, preventing damage to the cut-off wheel and the workpiece.

After contact with the workpiece, the cut-off wheel is automatically retracted 2mm, to be ready for cutting.

**OptiFeed**


During cutting, Exotom-100 continuously measures the load on the cutting arm.

The pre-set Feed and Force values are interpreted as maximum values. Throughout cutting, Exotom-100 will try to stay as close to these values as possible.

The factors that determine the load, are the shape and properties of the workpiece.

Whenever the set maximum Force limit is reached, Exotom-100 will reduce the Feed rate to be able to maintain the present force.

The figure below illustrates the increase in force as the cut-off wheel gets closer to the centre of a round workpiece.

### **Cutting Long Workpieces**

Workpieces that are wider than the cutting chamber can be positioned so that one or both ends can protrude through the rubber curtains on both sides of the cutting chamber.

### **Clamping Irregular Workpieces**

Irregular workpieces without plane clamping surfaces must be clamped using special clamping tools, as the workpieces must not move during cutting. This could result in damage to the cut-off wheel or to the sample itself. Use the T-slots to mount the special clamping tools. Struers offers a kit of Clamping Tools (See Accessories).

To achieve faster cutting, position the workpiece so that the wheel will cut the smallest possible cross-section.

### **Connection to External Exhaust System**

Struers recommends the use of an exhaust system as the workpieces may emit harmful gasses or dust during cutting. A joint for connection of a 80 mm dia. exhaust hose is located on the left side of Exotom-100.

- Remove the red cap.
- Mount an exhaust hose from your local exhaust system onto the flange.

## **Other Safety Features**

To prevent accidental access to the cutting chamber via the cutting curtains, each opening is covered by a spring-loaded guard, which only opens outwards.

The protection hood has a safety switch to prevent the cut-off wheel from starting while the hood is open. Furthermore, a locking mechanism prevents the opening of the protection hood before the cut-off wheel has come to a complete stop.

The motors of Exotom-100 are protected against overload. Should the motors overheat and/or overload, the motors will disengage until a normal temperature has been obtained.

Exotom-100 is prepared for the connection of an external warning light. The connection socket is situated on the right hand side of the unit.

## Optimising Cutting Results

The following table shows possible answers to a number of common questions:

<b>Optimising the Cutting Results</b>	
<b>Question</b>	<b>Answer</b>
How can I avoid discoloration or burning of the sample?	Use a lower Feed Speed.
	Change the cut-off wheel as the hardness of the present cut-off wheel may be inappropriate for the hardness of the sample. *)
How can I avoid burrs?	Use a softer cut-off wheel. *)
	Clamp the workpiece securely at the right hand clamping device. Tighten the left hand clamping device just enough to prevent the workpiece from moving when being cut.
How can I avoid the cut-off wheels wearing too quickly?	Use a lower Feed Speed, a different cutting mode or a harder cut-off wheel. *)
How can I achieve faster cutting?	Position the workpiece so as to cut the smallest possible cross-section. Use a high Feed Speed. If workpiece shape and properties permit, change to ExciCut or AxioCut/Step (optional) Cutting Modes

\*) See Consumables Section, Cut-off Wheels

## 2. Accessories

Specification	Code
<b>AxioCut Kit for Exotom-100</b> Increases cutting capacity (depth) by 150 mm. Requires installation by Struers Service Technician. Gross weight 9 kg	EXOMA
<b>Replacement Stainless Steel Bands</b> For Cutting Table with 12 mm T-slots. Available in three sizes, each supplied in sets of two: 100 x 225 mm 100 x 300 mm 60 x 548 mm	EXOBS EXOBM EXOBL
<b>Quick-clamping Device for 12 mm T-slots, Left</b> For securing the workpiece. Complete with back stop. To be mounted on a cutting table and placed to the left of the cut-off wheel.	EXOLE
<b>Quick-clamping Device for 12 mm T-slots, Right</b> For securing the workpiece. Complete with back stop. To be mounted on a cutting table and placed to the right of the cut-off wheel.	EXORI
<b>Clamping Tools for 12 mm T-slots</b> For clamping irregularly shaped workpieces. Complete with clamps, supports and bolts.	MAGOF
<b>Vertical Clamping System for 12 mm T-slots</b> For clamping irregularly shaped workpieces. Complete with operating key and one flat clamping shoe.	EXOVS
<b>Chain Spanner for 12 mm T-slots</b> For clamping cylindrical or irregularly shaped workpieces. Complete with anchor block and operating key.	EXOCS
<b>Swivel Shoes for Vertical Clamping System</b> For clamping irregularly shaped workpieces. Set of four multi-shaped swivel shoes.	TREVI
<b>Riser Block for 12 mm Vertical Clamping System</b> For elevating Vertical Clamping System when clamping high workpieces	EXOKS
<b>Arm Extension for 12 mm Vertical Clamping System</b> For extending arm of Vertical Clamping System	EXOAR
<b>Spare Nylon Filter for Recirculation Cooling System</b> For insertion into filter-drawer in cooling unit compartment.	EXOFI

### 3. Consumables

#### Cut-off Wheels

Application	Dimensions	Grinding agent	Code
For ExciCut and direct cutting of extremely hard steels (350-800 HV)	432x3x32mm	Al <sub>2</sub> O <sub>3</sub>	101MA
For ExciCut and direct cutting of very hard and ductile ferrous metals. Strong, fibre-reinforced wheel (250-700 HV)	432x3x32mm	Al <sub>2</sub> O <sub>3</sub>	202MA
For ExciCut and direct cutting of steels (< 500 HV). Extremely hard (white cast iron) or ductile (18/8 or St60) steel types are cut with ExciCut. Gives good wheel economy	432x3x32mm	Al <sub>2</sub> O <sub>3</sub>	102MA
For ExciCut and direct cutting of steels (< 500 HV). Gives very good wheel economy.	432x3x32mm	Al <sub>2</sub> O <sub>3</sub>	104MA
For ExciCut and direct cutting of non-ferrous metals (< 300 HV). Primarily for ExciCut.	432x3x32mm	SiC	106MA
For direct cutting of titanium and its alloys, as well as other non-ferrous metals (< 400 HV)	432x3x32mm	SiC	106MA
For direct cutting of sintered carbides and ceramics	305x1.8x32mm	Diamond	25EXO
	350x1.8x32mm	Diamond	26EXO

#### Other Consumables

Specification	Code
<i>Additive for Cooling Fluid</i> 1 l 5 l	ADDUN ADDFI
<i>Bacterial Cleaning Liquid</i> 0.5 l	BACTE

## 4. Trouble-Shooting

<b>Error</b>	<b>Explanation</b>	<b>Action</b>
<b>Machine Problems</b>		
Water leaking.	Leak in recirculation water hose.	Check the hose and tighten the hose clamp.
	Water overflow in the recirculation water tank.	Remove the excess water in the tank.
Samples or cutting chamber corroded.	Insufficient additive for cooling fluid.	Add Struers Additive for cooling fluid to the cooling water, using the correct concentration. Check with a refractometer. Follow the instructions in the Maintenance section.
	The machine is left with closed protection hood.	Leave the protection hood open to let the cutting chamber dry.
Quick-clamping device unable to hold the workpiece.	The quick-clamping device is not in balance.	Adjust the screw underneath the clamping column. Use a 3mm Allen key.
	Clamping heart worn.	Call Struers service technician.

*Exotom-100*  
*Instruction Manual*

<b>Error</b>	<b>Explanation</b>	<b>Action</b>
<b>Cutting Problems</b>		
Discoloration or burning of the sample.	The hardness of the cut-off wheel is inappropriate for the hardness / dimensions of the sample.	See Consumables section, Cut-off Wheels.
	Inadequate cooling.	Check that there is enough water in the recirculation cooling unit. Check the condition of the cooling tray.
	Feed Speed too high.	Reduce the Feed Speed.
Unwanted burrs.	Wheel too hard.	See Consumables section, Cut-off Wheels.
	Insufficient support of workpiece.	Add further support to the workpiece.
The cutting quality differs.	Cooling water hose clogged.	Clean the cooling water hose and the cooling tube. Check the water flow by turning the cooling valve to cleaning position.
	Insufficient cooling water.	Refill tank with water. Remember to add Struers Additive.
The cut bends to one side.	Feed Speed is too high.	Reduce the Feed Speed.
The cut-off wheel breaks.	Incorrect mounting of the cut-off wheel.	Check that the centre-hole has the correct diameter. Check cardboard washer on both sides of the cut-off wheel. The nut must be tightened properly.
	Incorrect clamping of the workpiece.	Make sure that only one of the quick-clamping devices is tight. The other device should only press lightly. Use support tools if the geometry of the workpiece makes support necessary.
	Wheel too hard.	See Consumables section, Cut-off Wheels.
	Feed Speed too high.	Reduce the Feed Speed.
	Inadequate cooling.	Check that there is enough water in the recirculation cooling unit. Check the cooling water hoses.
The cut-off wheel wears down too quickly.	Feed Speed too high.	Reduce the Feed Speed.
	Insufficient cooling.	Check that there is enough water in the recirculation cooling unit. Check the cooling water hoses.
	The cut-off wheel is too soft for the task.	See Consumables section, Cut-off Wheels.
	Exotom-100 vibrates (worn bearings).	Call Struers service technician.

*Exotom-100*  
*Instruction Manual*

<b>Error</b>	<b>Explanation</b>	<b>Action</b>
The cut-off wheel does not cut through the sample.	Incorrect choice of cut-off wheel.	See Consumables section, Cut-off Wheels.
	Cut-off wheel worn.	Replace the cut-off wheel.
	The cut-off wheel gets caught in the workpiece.	Support the workpiece and clamp it on both sides of the cut-off wheel in such a way to allow the cut to stay open.
	Incorrect choice of Cutting Mode. AxioCut/Step (option) is intended for large workpieces.	See section on Operation, Cutting Mode.
The workpiece breaks when clamped.	The workpiece is brittle.	Place the workpiece between two polystyrene plates. NB! Always cut brittle workpieces very carefully.
The sample is corroded.	The sample is not resistant to water.	Use a neutral liquid as cooling fluid or cut without using cooling fluid at all. <b>DO NOT USE AN INFLAMMABLE LIQUID</b>
	The sample has been left in the cutting chamber for too long.	Leave the protection hood open, when you leave the machine.
	Insufficient additive for cooling fluid.	Add Struers Additive for cooling fluid to the cooling water in the correct concentration. Check with a refractometer. See Maintenance section.
AutoStop does not stop the cutting action.	Feed Speed and/or Force is set too low to allow load change to be detected.	If problem persists, contact local Struers service engineer, for adjustment of AutoStop sensitivity.

**Error Messages**

Error messages are divided into three classes:  
Messages  
Errors  
Fatal Errors

*Messages*

Messages are intended to inform the operator of the machine's progress and advise about minor operational errors.

*Errors*

Errors must be rectified before cutting can be continued.

*Fatal Errors*

In case of Fatal Errors, cutting cannot continue before an authorised technician has rectified the error. Turn off the unit at the Main Switch immediately. Do not attempt to operate unit before technician has rectified problem.

<b>Message</b>	<b>Explanation</b>	<b>Action</b>
Process in Progress.	Start is not possible as other process in progress e.g. emptying of recirculation tank.	
Process already stopping.	Appears if Stop is pressed repeatedly.	
Water level in tank too low; or pump filter clogged.	Water level in recirculation tank is too low. Alternatively, recirculation jumper filter clogged.	Refill circulation water if needed. Check condition of pump filter. Clean if needed.
Water level low; emptying finished.	Emptying of recirculation cooling unit completed.	Switch off recirculation pump to complete emptying process.
Flush nozzle not parked.	Start not possible. Start was pressed while flush nozzle was out of holder.	Park flush nozzle to continue.
Flushing was interrupted to avoid pump damage. Please park the flush nozzle to continue.	The pump automatically stops after 3 minutes to avoid pump damage.	Park flush nozzle to continue.
Emptying of tank was stopped to avoid pump damage.	The recirculation pump automatically stops after 8 minutes to avoid pump damage.	To re-start the pump, turn the pump handle horizontally, then re-activate the pump by turning the pump handle to the vertical position.
Valve handle in emptying tank pos.	Drain pump handle has not been moved to off position after use.	Re-position drain pump handle.
Do not operate handle now.	During cutting, an attempt is made to activate tank emptying function manually.	
Reference position not found.	Stop was pressed during search for reference point.	Press Start again to do new search. If message remains, switch unit off and on at Main Switch to start new search for reference point.
Cutting stopped by flange sensor.	Flange sensors next to the cut-off wheel flange will stop the movement of the cutting arm if it meets an obstacle.	Replace cut-off wheel or re-position workpiece.

*Exotom-100*  
*Instruction Manual*

<b>Message</b>	<b>Explanation</b>	<b>Action</b>
Searching for reference position.	When unit is switched on, a reference position process is performed.	Wait until process has stopped.
Protection hood not closed.		Close hood properly.
Emergency stop is active.		Release emergency stop(s).
Out of Cutting Range.	Wheel at bottom position when Start is pressed.	Move cut-off wheel to higher start position.
Cutting finished.	Message after each cutting which has been completed normally. Protection hood can be opened.	
Process stopping.	Appears after Stop has been pressed once.	
<b>Error Message</b>	<b>Explanation</b>	<b>Action</b>
Main supply voltage too low.	Supply voltage insufficient.	Await normalisation of power grid.
Cut-off motor overloaded.	Wait until the motor has cooled down.	Press Start.
ExciCut motor overloaded.	Wait until the motor has cooled down.	Press Start.
Protection hood not locked.	Safety lock has failed to lock hood.	Re-close hood.
Cut-off motor blocked.	Cut-off motor may fail to move the wheel if caught in the workpiece.	Move cut-off wheel clear of workpiece before pressing Start.
Contact K3 not activated.	Contact fails to operate.	Re-start. If error remains, contact Struers service technician.
Contact K4 not activated.	Contact fails to operate.	Re-start machine using Start and Stop buttons. Alternatively, re-start using Main Switch If error remains, contact Struers service technician.
Step motor or encoder error.	Exotom-100 cannot complete reference point search.	Re-start using Main Switch. If error remains, contact Struers service technician.
Reference sensor not activated.	Search for reference point cannot be conducted.	Check whether cutting arm is blocked. Re-start using Main Switch. If error remains, contact Struers service technician.
Step motor may be blocked.	Cutting arm cannot move, most likely due to blockage of arm.	Remove blockage. If error remains, contact Struers service technician.
AxioCut motor blocked.	Cutting arm cannot move, most likely due to blockage of arm or AxioCut motor may fail to move the wheel if gets caught in the workpiece.	Move obstacle or move cut-off wheel clear of workpiece before pressing Start.
Recirculation pump overloaded.	Water hose blocked or too much swarf in coolant.	Replace coolant. Inspect pump.

*Exotom-100*  
*Instruction Manual*

<b>Fatal Error</b>	<b>Explanation</b>	<b>Action</b>
Protection hood lock not released.		Re-start machine using Start and Stop buttons. Alternatively, re-start using Main Switch.
Cut-off motor will not stop.		Turn off unit at Main Switch. Contact Struers service technician.
Contactora K3 not deactivated.	Error in control system.	
Contactora K4 not deactivated.	Error in control system.	
Safety relay error during cutting.		
Force system uncalibrated.	Force-measuring system is not calibrated.	
No serial communication.	No contact between machine system and operating panel.	
Nominal motor current undefined.	Motor current has not been defined during set-up of unit.	
Program version mismatch.	Conflict between software in machine system and operating panel.	
AxioCut motor driver error.	When unit is switched on, AxioCut motor starts operating.	

## **5. Maintenance**

### **Daily Service**

- Clean the cutting chamber, especially the cutting table and the T-slots.

### **Cleaning the Cutting Chamber**

To ensure a longer lifetime for your Exotom-100, Struers strongly recommends daily cleaning of the cutting chamber with the flushing hose:

- Lift the flushing hose from the holder on the rear wall. The recirculation pump is activated.
- Point the flushing hose towards the bottom of the cutting chamber.
- Press the lever on the flush nozzle to turn the water on. Adjust the flow using the lever.
- Clean the cutting chamber thoroughly.
- Turn off the water by releasing the lever on the flush nozzle.

***IMPORTANT***

Do **not** clean the lamp glass with alcohol.  
Use a soft, damp cloth.

### **Weekly Service**

- Refill the Recirculation Cooling Unit. The water level should be 5-10 cm below the upper edge of the tank.
- Clean the cutting chamber thoroughly.
- Clean the filter in the filter drawer.
- Check and if necessary empty and clean the sieve in the drawer of the Recirculation Cooling Unit.

***IMPORTANT***

Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to add Struers Additive each time you refill with water.

### **Monthly Service**

#### *Replacing the Cooling Water*

- Replace the cooling water in the Recirculation Cooling Unit at least once a month.
- Fill the tank with 150 l of water. The water level should be 5-10 cm below the upper edge of the tank.
- Add Struers Additive according to the instructions.

## Maintenance of Recirculation Unit

### Low Level Warning

The Exotom-100 will automatically display a warning if the level of cooling water gets too low. However, it is recommended that the cooling water level is checked every week.

### Refilling the Cooling Unit

- Refill the cooling unit with water using an external water hose or via direct connection to the water mains.
- Stop refilling when the water level is 5-10 cm from the upper edge.

**NB!**

Never refill by pouring the water into the cutting chamber, as it would be impossible to check the increasing level in the tank.

**IMPORTANT**


Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to add Struers Additive each time you refill with water.

Use only Struers Additive ADDUN in the Recirculation Cooling Unit.

ADDUN is specially selected for use with Struers cut-off machines.

Other additives may not be compatible with some of the components of the cut-off machine.

### Emptying the Cooling Unit

When replacing the cooling water, the cooling unit is emptied using the recirculation pump. Turn the black pump handle in the cooling unit compartment to the vertical position to activate the pump.

By design, the pump cannot empty the tank completely, preventing cutting waste from entering the cooling system. Tilt the tank, to release the remaining water.


**IMPORTANT**

Remember to turn the pump handle horizontally (off position) after use.

**Note**

The recirculation pump automatically stops after 8 minutes to avoid pump damage. To re-start the pump, turn the pump handle horizontally, then re-activate the pump by turning the pump handle to the vertical position.

Filter Drawer

The filter drawer is situated on the left in the cooling unit compartment. A nylon sieve filters larger cutting waste particles from the cooling water. It is recommended to always have the filter drawer mounted during cutting.

Depending on use, the filter must be checked and emptied regularly, to ensure its effectiveness. Should the filter get blocked, an overflow function will prevent flooding, and direct the unfiltered cooling water into the tank.


Additive for Cooling Water

Remember to add Struers Additive for Cooling Fluid: One part of Additive for 33 parts of water.

The concentration of additive should always lie between 2.7 and 3.3%; the pH should be around 9. Replace the cooling water if pH is  $\leq 8.0$ .

To check the concentration of additive, use a refractometer (0-10 Brix, see diagram.). Concentration =  $2 \times \text{Brix}$ .

Add Struers Additive for Cooling Fluid if the concentration is lower than 2.7%. Add water if the concentration is higher than 3.3%.

### **Maintenance of Cutting Table**

The stainless steel bands forming the cutting table should be replaced if they have become worn or damaged. The bands are available as spare parts.

To allow humidity to escape from the cutting table and chamber, it is recommended to leave the hood open when the machine is not in use.

### **Maintenance of Cut-off Wheels**

*Storage of Bakelite Bonded  
Al<sub>2</sub>O<sub>3</sub> Cut-off Wheels*

This type of cut-off wheel is sensitive to humidity. Therefore, do not mix new, dry cut-off wheels with used humid ones. Store the cut-off wheels in a dry place, horizontally on a plane support.

*Maintenance of Diamond and  
CBN Cut-off Wheels*

The precision of diamond and CBN cut-off wheels (and thus the cut) depends on how carefully the following instructions are observed:

- Never expose the cut-off wheel to overload, such as heavy mechanical load, or heat.
- Store the cut-off wheel in a dry place, horizontally on a plane support, preferably under light pressure.
- A clean and dry cut-off wheel does not corrode. Therefore, clean and dry the cut-off wheel before storing. If possible, use ordinary detergents for the cleaning.

### **Maintenance of Clamping Devices**

***Important***

It is recommended to thoroughly clean and lubricate the Quick Clamping Device and Vertical Clamping Device at regular intervals.


## 6. Technical Data

Subject		Specification	
		Metric/International	US
<b>CUTTING SPECIFICATIONS</b>			
<b>Workpiece Dimensions (Max)</b>	<i>Workpiece inside cutting chamber:</i>		
	Height	250 mm	10"
	Width	650 mm	25.6"
	Depth	550 mm	20"
	<i>Workpiece protruding cutting chamber:</i>		
	Height	140 mm	5.5"
	Depth	210 mm	8.3"
<b>Cutting Capacity (Max)</b> <i>Please refer to Cutting Diagram</i>	Max. workpiece diameter	160 mm	6.3"
	<i>Max. size of cut. In brackets, with AxioCut option</i>		
	Height	100 mm (100 mm)	3.9" (3.9")
	Depth	300 mm (450 mm)	11.8" (17.7")
<b>PHYSICAL SPECIFICATIONS</b>			
<b>Cutting Motor</b>	Cutting power*	10.5 kW	14 HP
	Maximum power*	15 kW	20 HP
	*) Except 460-480 V / 60 Hz		
<b>Cut-off Wheel</b>	Diameter x Thickness x Centre-hole	432 x 3 x 32 mm	17 x 0.12 x 1.26"
	Rotational speed (running idle)	1950 rpm	1950 rpm
<b>Positioning &amp; Feed</b>	Positioning range (of cut-off wheel)	0 – 250 mm	0 – 9.8"
	Max. positioning speed	60 mm/s	2.4"/s
	Feed Speed range (adjustable in steps of)	0.05 – 5 mm/s (0.05mm/s)	0.002 – 0.2"/s (0.002"/s)
	Cutting force	50-700 N	10-150 lbf
<b>Cutting table</b>	Width	855 mm	33.6"
	Depth	550 mm	21.6"
	T-slots	12 mm	0.48"
<b>Dimensions and Weight</b>	Height (closed/open hood)	1756/2300 mm	69/90.5"
	Width (ex-/including control panel)	1050/1350 mm	41.3/53.1"
	Depth	1500 mm	59.1"
	Weight	820 kg	1810 lbs
<b>Recirculation Cooling Unit</b>	Tank volume	150 l	39.6 gallons
	Approx. flow	25 l/min	6.6 gallons/min

*Exotom-100*  
*Instruction Manual*

Subject		Specification			
<b>Electrical Data</b>	Voltage / frequency:	Cutting Power*	Max. Power	Nom. Load	Max. Load
	3 x 200 V / 50 Hz	10.5 kW	15 kW	36 A	69 A
	3 x 200-210 V / 60 Hz	10.5 kW	15 kW	35 A	67 A
	3 x 220-230 V / 50 Hz	10.5 kW	15 kW	31 A	60 A
	3 x 220-240 V / 60 Hz	10.5 kW	15 kW	29 A	55 A
	3 x 380-415 V / 50 Hz	10.5 kW	15 kW	18 A	34 A
	3 x 380-415 V / 60 Hz	10.5 kW	15 kW	17 A	32 A
3 x 460-480 V / 60 Hz	12.6 kW	18 kW	18 A	34 A	
	*) at intermittent duty, S3 15%				
<b>Environment</b>	Noise level	Approx. 78 dB(A) running idle, at a distance of 1.0 m / 39.4" from the machine.			

**Cutting Capacity**

The graph shows the projected cutting capacity under the following conditions: A new cut-off wheel. The workpiece is laid directly on the cutting table, with overhang where appropriate. Vertical clamping is used. The actual cutting capacity depends on the sample material, cut-off wheel and clamping technique.

# Quick Reference Card

## Clamping the Workpiece

- Place the workpiece between the clamp and the back stop on the right hand cutting table.
- Push the clamp towards the workpiece and lock the quick-clamping device with the locking handle.

## Starting the Cutting

- Position the cut-off wheel.
- Close the protection hood.
- Press START  $\blacklozenge$ . The cut-off wheel starts rotating and the cooling water starts running.

## Stopping the Cutting

### *Auto Stop*

- In the Cutting Menu, toggle between FORCE/FEED/STOP pushing the knob, until STOP is highlighted.
- Turn knob to move pointer to setting AUTO at the bottom of the STOP column.

### *Manual Stop*

- Press STOP  $\ominus$  and the cutting process stops. The cut-off wheel stops rotating and the cooling water stops.

## Setting the Cutting Parameters

- Press MENU  $\square$  once if CONFIGURATION MENU is selected.
- Push knob to toggle between parameters in the Cutting Menu.
- Turn knob to edit the value of the selected parameter.
- Push knob to move to next parameter.

## Changing the Cut-off Wheel

- Press the knob for the spindle lock and turn the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner (30 mm).
- Remove the flange and the cut-off wheel.
- Mount the new cut-off wheel.
- Mount the flange and nut. Tighten carefully and close the guard.

## Cleaning the Cutting Chamber

- Point the flushing hose towards the bottom of the cutting chamber.
- Press the lever on the flush nozzle to turn the water on. Adjust the flow using the lever.
- Clean the cutting chamber thoroughly.
- Turn off the water by releasing the lever on the flush nozzle.

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*Exotom-100*  
*Spare Parts and Diagrams*

---

**Always state *Serial No* and *Voltage/frequency*  
if you have technical questions or when ordering spare parts.**

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Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2003.

**Struers A/S**  
Pederstrupvej 84  
DK-2750 Ballerup  
Denmark  
Telephone +45 44 600 800  
Telefax: +45 44 600 801

---

## **Spare Parts and Diagrams**

### **Table of contents**

### **Drawing**

## **Exotom-100**

### *Drawings*

Exotom-100, complete.....	15040001E
Chassis, complete.....	15040003L
Pump, assembly.....	15040009C
Chassis, assembly.....	15040010K
Wire drum, assembly .....	15040013B
Counter weight, assembly.....	15040014
Cutting wheel guard.....	15040022D
Flushing Gun, complete .....	15040060
Shock absorber, complete.....	15040076A
Cutting table, complete .....	15040005G
Long-cut kit, complete.....	15040040C
Fixed rod .....	15040041A
Long-cut motor .....	15040044A
ExiCut motor .....	15040046D
High power electrical assembly .....	15040011D
Protection Guard, complete.....	15040016I
Lamp, complete.....	15040017
Top hat, assembly.....	15040036B
Cutting arm, assembly .....	15040021E
Tilting frame, assembly.....	15040025B
Cutting motor, complete.....	15040028D
Door for recirculation, RHS.....	15040031A
Door for recirculation, LHS.....	15040032A
Control panel, complete.....	15040050
Control panel assembly, front.....	15040055C
Water Panel, complete.....	15040065D
Feed Actuator, complete.....	15040070F
Stepper Motor, assembly .....	15040074B
Snap lock tube, complete.....	15040067
Quick clamping device for 12 mm T-slots, RHS compl. ....	15040081
Quick clamping device for 12 mm T-slots, LHS compl. ....	15040098
Back Stop for 12 mm T-slots, RHS complete.....	15040084
Back Stop for 12 mm T-slots, LHS complete.....	15040093

*Exotom-100*  
*Spare Parts and Diagrams*

*Diagrams*

Exotom-100 Block Diagram .....	15043050
Exotom-100, Control voltage diagram (5 pages) .....	15043110A
Exotom-100, Main supply circuit .....	15043100A
Exotom-100, Control panel.....	15043200B
Exotom-100, Machine control board, A1 Power supply (8 pages) .....	15043210A
Exotom-100, LED circuit board .....	15043220
Exotom-100, Wiring diagram (5 pages) .....	15043450B
Primary connections on transformer AA-70951B .....	15043470A

Some of the drawings may contain position numbers  
not used in connection with this manual.

*Exotom-100*  
*Spare Parts and Diagrams*

<b>Drawing</b>	<b>Spare Part list for Exotom-100</b>		<b>Cat no.</b>
	<b>Pos.</b>		
		Rough Filter	15040397
		101MA Cut-off wheel ø432mm (5 pcs)	40009001
		104MA Cut-off wheel ø432mm (5 pcs)	40009004
		202MA Cut-off wheel ø432mm (5 pcs)	40009003
		ADDFI Additive 5L	49900003
		Fork Spanner 30mm	2GR00230
		Fork Spanner 13mm	2GR00213
		Pressure hose, water, white	2NU93020
15040001		<b>Exotom-100, complete</b>	
	150	Curtain, back wall	15040315
15040003		<b>Chassis, complete</b>	
	20	T1: Transformer 200-460V/ 600VA	2MT70951
	140	PCB A1, Machine control	15043393
	390	Fuse 4.00AT Glass (10 pcs)	RFU14200
	590	PD Eagle Belt W1792	2JT61792
	1090	Cardboard blotters 25 Pcs, ø32.5/120	381MP431
	1100	Loose Flange	R5040160
	1340	Shock Absorber with screws, complete	R5040076
15040009		<b>Pump, assembly</b>	
	20	Kit for impeller pump	15042901
15040010		<b>Chassis, assembly</b>	
	326	Snap Lock hose, complete	15040067
	390	Filter for impeller pump, complete	15040096
	400	Flushing gun, complete	15040060
	520	Perlator 3/8"	2YI20038
15040013		<b>Wire drum, assembly</b>	
	50-60	Bearing house, complete ø28	15042905
	100	Bearing Flange w. House PFT 25 FM	RBK40025
15040014		<b>Counter weight, assembly</b>	
	20	Glide guides for weight (4 pcs)	RBG00121
15040022		<b>Cutting wheel guard, assembly</b>	
	30	Sensor strap with Screws	R5040166
	50	Proximity Switch DJ2 (Namur)	2HQ00022


*Exotom-100*  
*Spare Parts and Diagrams*


<b>Drawing</b>	<b>Pos.</b>		<b>Cat no.</b>	
15040060		<b>Flushing gun, complete</b>		
	10	Flushing gun, complete	15040060	
	30	Perlator for flushing gun	2YI26211	
		Shower Brush	15042938	
15040076	40	Shower Hose	2NU99010	
		<b>Shock Absorber, complete</b>		
	40	Update shock absorber cartridge kit. Up to s/n 5040020.	RYS10016	
	40	Replacement shock absorber cartridges. Use only from s/n 5040020.	15040319	
15040005	70-80	Bearing kit for Shock Absorber	15042906	
		<b>Cutting table, complete</b>		
	110	Timing Belt GT-36-5MR-09	2JT20450	
	120	EXOBL Replacement steel bands for 12 mm T-slots 60x548 (Set of 2 pcs)	05046909	
	150	EXOBM Replacement steel bands for 12 mm T-slots 100x300 (Set of 2 pcs)	05046908	
15040040	160	EXOBS Replacement steel bands for 12 mm T-slots 100x225 (Set of 2 pcs)	05046907	
	340	Spring kit for side flap hinge	15042904	
		<b>Long-cut kit, complete</b>		
15040041	100	M5, Long cut motor assembly	15040044	
	120	Timing belt, GT 290 3MR09	2JT60270	
15040044		<b>Fixed rod</b>		
	20	Vibration damper	2GS41260	
15040046	40	Joint bearing GA20	2BL20020	
		<b>Long-cut Motor, M5</b>		
15040044	10	M5, Long cut motor assembly	15040044	
		<b>Excicut Motor, M3</b>		
15040046	10	M3-Motor 1x230V/24W	15040167	
		<b>High power electrical assembly</b>		
15040011				
	F1 (3x200-240V/50-60Hz):	51	6.00 A Midget fuse 10x38 (3 pcs)	RFC10060
	F1 (3x380-480V/50-60Hz):	51	4.00 A Midget fuse 10x38 (3 pcs)	RFC10040
	F2 (3x200-240V/50-60Hz):	52	6.00 A Midget fuse 10x38 (3 pcs)	RFC10060
	F2 (3x380-480V/50-60Hz):	52	4.00 A Midget fuse 10x38 (3 pcs)	RFC10040
	F3 (3x200-240V/50-60Hz):	53	6.00 A Midget fuse 10x38 (3 pcs)	RFC10060
F3 (3x380-480V/50-60Hz):	53	4.00 A Midget fuse 10x38 (3 pcs)	RFC10040	

*Exotom-100*  
*Spare Parts and Diagrams*

<b>Drawing</b>	<b>Pos.</b>	<b>Chain Spanner for 12 mm T-slots EXOCS</b>	<b>Cat no.</b>
	1	Chain Spanner without key	2GR26910
	2	Spring for EXOCS, w. screw	2GR26023
	3	Chain, 500 mm	2GR26925
	4	Operating key for EXOCS (12 mm T-slots)	2GR26031
	5	Anchor block	2GR26930

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



EMERGENCY STOP

- Outputs:
- 1x 24Vac contactor for main motor (K1)
  - 3x 24Vac contactors for star/delta motorstarter (K2-K4)
  - 1x 24Vac contactor for rec. water pump (K5)
  - 1x 24Vac contactor for excit motor (K6)
  - 1x 24Vac solenoid valve for flushing
  - 1x 24Vac solenoid valve for flushing
  - 1x 24Vac stepmotor driver for feed
  - 1x 5A dc-motor driver for y-cut

Rev. A: (FTH 04-10-99) XS7 added	STRUBER A/S Vejlevej 176 DK-8260, Hvidebø +45-3670-3900	Scale	Sheet	of	5
EXOTOM-100. CONTROL VOLTAGE DIAGRAM		Size	15043110A		
24 VAC OUTPUTS, STEP MOTOR AND DC MOTOR OUTPUT		A2			
FTH / FTH		Scale			
Tuesday, October 05, 1999					

COLOR CODES:

- BK = BLACK
- BN = BROWN
- RD = RED
- OG = ORANGE
- YE = YELLOW
- GN = GREEN
- BU = BLUE
- VT = VIOLET
- GY = GREY
- WH = WHITE

S1, EMERGENCYSTOP


S2, EMERGENCYSTOP

To A1:X12-2 (page 3)










To A1:X12-1 (page 3)

To A1:X10-1 (page X)

To A1:X10-2 (page X)

Rev. (FTH 04-10-99) XS7 added		Source A/B Vejle, Aale 1.76 Vejle, Aale 1.76 Dinmark Halskøbenhavn - 46.3000.9650
Size	Page Code	Page No
A2		15043110A
Scale		Sheet 4 of 5
Tuesday, October 05, 1999		FTH / FTH


COLOR CODES: BK = BLACK  
 BN = BROWN  
 RD = RED  
 OG = ORANGE  
 YE = YELLOW  
 GN = GREEN  
 BU = BLUE  
 VT = VIOLET  
 GY = GREY  
 WH = WHITE

Rev.A: (FTH 04-10-99) XS7 added		Source: A/S Vejlevej 179 DK-8000 Ålborg Danmark Illustration: L&S 9903369
Size	CAGE Code	EWG NO
A2		
Scale		Sheet
Tuesday, October 06, 1999		5 of 5


EXOTOM-100, CONTROL VOLTAGE DIAGRAM		
CONTROL PANEL CIRCUIT		
Rev	1	
Part No	15043110A	

All wires AWG 8

AWG16


For connection see other diagram

Must be short circuit protected with external fuses  
Fmax: See note #1.

All wires AWG 10


VOLTAGE / FREQ. (from nameplate)	Note #1: (max. ext. fuse)	F1, F2, F3 (fuse size)	M1 (nom. power)	M1 S3, 10 min. (cutting power)	M2 (power)	Note #2: (M2-connection)
3x200V / 50Hz (CSA)	3x 80 AT	3x 6 AT	7.5 kW	10.5 kW	470 W	Delta
3x200V / 60Hz (CSA)	3x 80 AT	3x 6 AT	7.5 kW	13.5 kW	590 W	Delta
3x220-240V / 50Hz	3x 80 AT	3x 6 AT	7.5 kW	10.5 kW	480 W	Delta
3x220-240V / 60Hz	3x 80 AT	3x 6 AT	7.5 kW	10.5 kW	480 W	Delta
3x380V-415V / 50Hz	3x 50 AT	3x 4 AT	7.5 kW	10.5 kW	480 W	Star
3x380V-415V / 60Hz	3x 50 AT	3x 4 AT	7.5 kW	10.5 kW	480 W	Star
3x480V-480V / 60Hz	3x 50 AT	3x 4 AT	9.0 kW	12.6 kW	660 W	Star
3x480V / 60Hz (CSA)	3x 50 AT	3x 4 AT	9.0 kW	12.6 kW	660 W	Star

System A/E: 4.76  
Rev: A: FTH-01-10-09  
Wire size corrected, CSA models added  
Date: 2010-10-09  
Drawing No.: 443-3070-3300


EXOTOM-100, MAIN SUPPLY CIRCUIT  
DWG. NO. 15043100A  
Scale: 1:1  
Date: Monday, October 11, 1999  
Sheet: 1 of 1

SUPPLY IN AND RS485 COMMUNICATION


POWER IN AND TOPFROM MACHB CONTROL BOARD

INCREMENTAL DIAL KNOB

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TOUCH PAD AND PUSH BUTTONS INPUT


TOUCHPAD INPUT

START STOP BUTTONS INPUT


OUTPUT TO POWER SIGNALS (AND DECOUPLING)


POWER SIGNALS


CPU AND MEMORY CIRCUIT

RESET CIRCUIT


JOYSTICK INPUT


JOYSTICK INPUT


ADDRESS SELECT & RW LOGIC

DISPLAY INTERFACE


CSFT INVERTER


REV. A	PROLOGIC CONTROL BOARD FOR EXOTOM 100	REV. B
REV. C	REV. D	REV. E
REV. F	REV. G	REV. H
REV. I	REV. J	REV. K
REV. L	REV. M	REV. N
REV. O	REV. P	REV. Q
REV. R	REV. S	REV. T
REV. U	REV. V	REV. W
REV. X	REV. Y	REV. Z
REV. AA	REV. AB	REV. AC
REV. AD	REV. AE	REV. AF
REV. AG	REV. AH	REV. AI
REV. AJ	REV. AK	REV. AL
REV. AM	REV. AN	REV. AO
REV. AP	REV. AQ	REV. AR
REV. AS	REV. AT	REV. AU
REV. AV	REV. AW	REV. AX
REV. AY	REV. AZ	REV. BA
REV. BB	REV. BC	REV. BD
REV. BE	REV. BF	REV. BG
REV. BH	REV. BI	REV. BJ
REV. BK	REV. BL	REV. BM
REV. BN	REV. BO	REV. BP
REV. BQ	REV. BR	REV. BS
REV. BT	REV. BU	REV. BV
REV. BW	REV. BX	REV. BY
REV. BZ	REV. CA	REV. CB
REV. CC	REV. CD	REV. CE
REV. CF	REV. CG	REV. CH
REV. CI	REV. CJ	REV. CK
REV. CL	REV. CM	REV. CN
REV. CO	REV. CP	REV. CQ
REV. CR	REV. CS	REV. CT
REV. CU	REV. CV	REV. CW
REV. CX	REV. CY	REV. CZ
REV. DA	REV. DB	REV. DC
REV. DD	REV. DE	REV. DF
REV. DG	REV. DH	REV. DI
REV. DJ	REV. DK	REV. DL
REV. DM	REV. DN	REV. DO
REV. DP	REV. DQ	REV. DR
REV. DS	REV. DT	REV. DU
REV. DV	REV. DW	REV. DX
REV. DY	REV. DZ	REV. EA
REV. EB	REV. EC	REV. ED
REV. EE	REV. EF	REV. EG
REV. EH	REV. EI	REV. EJ
REV. EK	REV. EL	REV. EM
REV. EN	REV. EO	REV. EP
REV. EQ	REV. ER	REV. ES
REV. ET	REV. EU	REV. EV
REV. EW	REV. EX	REV. EY
REV. EZ	REV. FA	REV. FB
REV. FC	REV. FD	REV. FE
REV. FF	REV. FG	REV. FH
REV. FI	REV. FJ	REV. FK
REV. FL	REV. FM	REV. FN
REV. FO	REV. FP	REV. FQ
REV. FR	REV. FS	REV. FT
REV. FU	REV. FV	REV. FW
REV. FX	REV. FY	REV. FZ
REV. GA	REV. GB	REV. GC
REV. GD	REV. GE	REV. GF
REV. GG	REV. GH	REV. GI
REV. GJ	REV. GK	REV. GL
REV. GM	REV. GN	REV. GO
REV. GP	REV. GQ	REV. GR
REV. GS	REV. GT	REV. GU
REV. GV	REV. GW	REV. GX
REV. GY	REV. GZ	REV. HA
REV. HB	REV. HC	REV. HD
REV. HE	REV. HF	REV. HG
REV. HH	REV. HI	REV. HJ
REV. HK	REV. HL	REV. HM
REV. HN	REV. HO	REV. HP
REV. HQ	REV. HR	REV. HS
REV. HT	REV. HU	REV. HV
REV. HW	REV. HX	REV. HY
REV. HZ	REV. IA	REV. IB
REV. IC	REV. ID	REV. IE
REV. IF	REV. IG	REV. IH
REV. II	REV. IJ	REV. IK
REV. IL	REV. IM	REV. IN
REV. IO	REV. IP	REV. IQ
REV. IR	REV. IS	REV. IT
REV. IU	REV. IV	REV. IW
REV. IX	REV. IY	REV. IZ
REV. JA	REV. JB	REV. JC
REV. JD	REV. JE	REV. JF
REV. JG	REV. JH	REV. JI
REV. JJ	REV. JK	REV. JL
REV. JM	REV. JN	REV. JO
REV. JP	REV. JQ	REV. JR
REV. JS	REV. JT	REV. JU
REV. JV	REV. JW	REV. JX
REV. JY	REV. JZ	REV. KA
REV. KB	REV. KC	REV. KD
REV. KE	REV. KF	REV. KG
REV. KH	REV. KI	REV. KJ
REV. KK	REV. KL	REV. KM
REV. KN	REV. KO	REV. KP
REV. KQ	REV. KR	REV. KS
REV. KT	REV. KU	REV. KV
REV. KW	REV. KX	REV. KY
REV. KZ	REV. LA	REV. LB
REV. LC	REV. LD	REV. LE
REV. LF	REV. LG	REV. LH
REV. LI	REV. LJ	REV. LK
REV. LM	REV. LN	REV. LO
REV. LP	REV. LQ	REV. LR
REV. LS	REV. LT	REV. LU
REV. LV	REV. LW	REV. LX
REV. LY	REV. LZ	REV. MA
REV. MB	REV. MC	REV. MD
REV. ME	REV. MF	REV. MG
REV. MH	REV. MI	REV. MJ
REV. MK	REV. ML	REV. MM
REV. MN	REV. MO	REV. MP
REV. MQ	REV. MR	REV. MS
REV. MT	REV. MU	REV. MV
REV. MW	REV. MX	REV. MY
REV. MZ	REV. NA	REV. NB
REV. NC	REV. ND	REV. NE
REV. NF	REV. NG	REV. NH
REV. NI	REV. NJ	REV. NK
REV. NL	REV. NM	REV. NN
REV. NO	REV. NP	REV. NQ
REV. NR	REV. NS	REV. NT
REV. NU	REV. NV	REV. NW
REV. NX	REV. NY	REV. NZ
REV. OA	REV. OB	REV. OC
REV. OD	REV. OE	REV. OF
REV. OG	REV. OH	REV. OI
REV. OJ	REV. OK	REV. OL
REV. OM	REV. ON	REV. OO
REV. OP	REV. OQ	REV. OR
REV. OS	REV. OT	REV. OU
REV. OV	REV. OW	REV. OX
REV. OY	REV. OZ	REV. PA
REV. PB	REV. PC	REV. PD
REV. PE	REV. PF	REV. PG
REV. PH	REV. PI	REV. PJ
REV. PK	REV. PL	REV. PM
REV. PN	REV. PO	REV. PP
REV. PQ	REV. PR	REV. PS
REV. PT	REV. PU	REV. PV
REV. PW	REV. PX	REV. PY
REV. PZ	REV. QA	REV. QB
REV. QC	REV. QD	REV. QE
REV. QF	REV. QG	REV. QH
REV. QI	REV. QJ	REV. QK
REV. QL	REV. QM	REV. QN
REV. QO	REV. QP	REV. QQ
REV. QR	REV. QS	REV. QT
REV. QU	REV. QV	REV. QW
REV. QX	REV. QY	REV. QZ
REV. RA	REV. RB	REV. RC
REV. RD	REV. RE	REV. RF
REV. RG	REV. RH	REV. RI
REV. RJ	REV. RK	REV. RL
REV. RM	REV. RN	REV. RO
REV. RP	REV. RQ	REV. RR
REV. RS	REV. RT	REV. RU
REV. RV	REV. RW	REV. RX
REV. RY	REV. RZ	REV. SA
REV. SB	REV. SC	REV. SD
REV. SE	REV. SF	REV. SG
REV. SH	REV. SI	REV. SJ
REV. SK	REV. SL	REV. SM
REV. SN	REV. SO	REV. SP
REV. SQ	REV. SR	REV. SS
REV. ST	REV. SU	REV. SV
REV. SW	REV. SX	REV. SY
REV. SZ	REV. TA	REV. TB
REV. TC	REV. TD	REV. TE
REV. TF	REV. TG	REV. TH
REV. TH	REV. TI	REV. TJ
REV. TK	REV. TL	REV. TM
REV. TN	REV. TO	REV. TP
REV. TQ	REV. TR	REV. TS
REV. TU	REV. TV	REV. TW
REV. TX	REV. TY	REV. TZ
REV. UA	REV. UB	REV. UC
REV. UD	REV. UE	REV. UF
REV. UG	REV. UH	REV. UI
REV. UJ	REV. UK	REV. UL
REV. UM	REV. UN	REV. UO
REV. UP	REV. UQ	REV. UR
REV. US	REV. UT	REV. UV
REV. UW	REV. UX	REV. UY
REV. UZ	REV. VA	REV. VB
REV. VC	REV. VD	REV. VE
REV. VF	REV. VG	REV. VH
REV. VI	REV. VJ	REV. VK
REV. VL	REV. VM	REV. VN
REV. VO	REV. VP	REV. VQ
REV. VR	REV. VS	REV. VT
REV. VU	REV. VV	REV. VW
REV. VX	REV. VY	REV. VZ
REV. WA	REV. WB	REV. WC
REV. WD	REV. WE	REV. WF
REV. WG	REV. WH	REV. WI
REV. WJ	REV. WK	REV. WL
REV. WM	REV. WN	REV. WO
REV. WP	REV. WQ	REV. WR
REV. WS	REV. WT	REV. WU
REV. WV	REV. WX	REV. WY
REV. WZ	REV. XA	REV. XB
REV. XC	REV. XD	REV. XE
REV. XF	REV. XG	REV. XH
REV. XI	REV. XI	REV. XI
REV. XJ	REV. XJ	REV. XJ
REV. XK	REV. XK	REV. XK
REV. XL	REV. XL	REV. XL
REV. XM	REV. XM	REV. XM
REV. XN	REV. XN	REV. XN
REV. XO	REV. XO	REV. XO
REV. XP	REV. XP	REV. XP
REV. XQ	REV. XQ	REV. XQ
REV. XR	REV. XR	REV. XR
REV. XS	REV. XS	REV. XS
REV. XT	REV. XT	REV. XT
REV. XU	REV. XU	REV. XU
REV. XV	REV. XV	REV. XV
REV. XW	REV. XW	REV. XW
REV. XX	REV. XX	REV. XX
REV. XY	REV. XY	REV. XY
REV. XZ	REV. XZ	REV. XZ
REV. YA	REV. YA	REV. YA
REV. YB	REV. YB	REV. YB
REV. YC	REV. YC	REV. YC
REV. YD	REV. YD	REV. YD
REV. YE	REV. YE	REV. YE
REV. YF	REV. YF	REV. YF
REV. YG	REV. YG	REV. YG
REV. YH	REV. YH	REV. YH
REV. YI	REV. YI	REV. YI
REV. YJ	REV. YJ	REV. YJ
REV. YK	REV. YK	REV. YK
REV. YL	REV. YL	REV. YL
REV. YM	REV. YM	REV. YM
REV. YN	REV. YN	REV. YN
REV. YO	REV. YO	REV. YO
REV. YP	REV. YP	REV. YP
REV. YQ	REV. YQ	REV. YQ
REV. YR	REV. YR	REV. YR
REV. YS	REV. YS	REV. YS
REV. YT	REV. YT	REV. YT
REV. YU	REV. YU	REV. YU
REV. YV	REV. YV	REV. YV
REV. YW	REV. YW	REV. YW
REV. YX	REV. YX	REV. YX
REV. YY	REV. YY	REV. YY
REV. YZ	REV. YZ	REV. YZ
REV. ZA	REV. ZA	REV. ZA
REV. ZB	REV. ZB	REV. ZB
REV. ZC	REV. ZC	REV. ZC
REV. ZD	REV. ZD	REV. ZD
REV. ZE	REV. ZE	REV. ZE
REV. ZF	REV. ZF	REV. ZF
REV. ZG	REV. ZG	REV. ZG
REV. ZH	REV. ZH	REV. ZH
REV. ZI	REV. ZI	REV. ZI
REV. ZJ	REV. ZJ	REV. ZJ
REV. ZK	REV. ZK	REV. ZK
REV. ZL	REV. ZL	REV. ZL
REV. ZM	REV. ZM	REV. ZM
REV. ZN	REV. ZN	REV. ZN
REV. ZO	REV. ZO	REV. ZO
REV. ZP	REV. ZP	REV. ZP
REV. ZQ	REV. ZQ	REV. ZQ
REV. ZR	REV. ZR	REV. ZR
REV. ZS	REV. ZS	REV. ZS
REV. ZT	REV. ZT	REV. ZT
REV. ZU	REV. ZU	REV. ZU
REV. ZV	REV. ZV	REV. ZV
REV. ZW	REV. ZW	REV. ZW
REV. ZX	REV. ZX	REV. ZX
REV. ZY	REV. ZY	REV. ZY
REV. ZZ	REV. ZZ	REV. ZZ

CPU AND MEMORY CIRCUIT


ANALOG/DIGITAL AND DIGITAL/ANALOG CONVERTERS


LINKS TO OTHER PAGES


RS-485 INTERFACE

RESET CIRCUIT





DECOUPLING


Sheet	5	of	8
Scale			
DWG NO	15043210		
Rev	I A		
<b>EXOTOM-100 MACHINE CONTROL BOARD A1</b> <b>STEP MOTOR CIRCUIT</b>			
Author	Suresh A/S - Rev A, Prototype corr. and 06/07/09 PTH		
Designer	DK-2110 Redwood		
Estimator	44530203800		

SUPERVISION OF K4

SUPERVISION OF K3

NOT USED

DOOR CLOSED AND LOCKED SUCCESSFULLY

COMMON


FLUSH PISTOL OFF HOOK (NAMUR SENSOR)

TANK EMPTYING (NAMUR SENSOR)

FLANGE COLLISION DETECTOR (NAMUR SENSOR)

Z-TOP POSITION (NAMUR SENSOR)

Rev A: Prototype contr. board (05.07.00 FTH)		Sheet	6
Drawn by: S. B. 1.76	Checked by: J. A.	Scale	1:1
DK-2010 (Rev. 0)	EXCITE Code: <C>Capp Code: <C>	DWG NO	15043210
1488183000 - 4151_367038000	EXOTOM-100: MACHINE CONTROL BOARD AI	Rev	1 A
Transfery, August 01, 2000	FTH / FTH	Scale	1:1

Low level switch  
in rec. bank

Overload switch  
in rec.pump motor

Overload switch  
in rec.pump motor

Emergency stop  
switch (left side)


Emergency stop  
switch (right side)

Status from  
door lock switch


page 2

REV: 01/05		DWG NO		REV	
Rev A: Prototype come out (05/07/99 FTH)		743		15043210	
VIA: 1/16" Dia		PAGE COPIES		7	
Dwg: 3/10 Release		SHEET		7 of 8	
MATERIAL: 4141_3670_3800		FTH / FTH			
EXOTOM-100, MACHINE CONTROL BOARD, A1					
Tuesdays, August 01, 2000					


ANALOG INPUTS


INPUT FOR WHEEL SPEED DETECTION


VREF BUFFER

DECOUPLING

Rev A: Prototype come out 06.07.08 (PTH)		REVISED BY	
Date of this revision		DATE	
Manufacture: +46.3670.3800		MANUFACTURE	
<b>EXOTOM-100, MACHINE CONTROL BOARD A1</b>			
Size	CAGE Code	DWG NO	Rev
742		15043210	1
Scab		PTH/PTH	A
Tuesday, August 01, 2000	Sheet	8	of 8


STRUERS A/S  
 Valbøjs Allé 176  
 DK-2610 Rødovre  
 Denmark  
 +45 3670 3500

### Exotom-100, LED Circuit board

Size	CAGE Code	DWG NO	Rev
A4	<Cage Code>	15043220	
Scale		FTH / FTH	Sheet 1 of 1
Tuesday, August 01, 2000			

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








Herewith declares that

<i>Product Name:</i>	Exotom-100
<i>Type No:</i>	504
<i>Machine Type:</i>	Cut-off machine

is in conformity with the provisions of the following directives:

**Safety of Machinery** 98/37/EEC according to the following standard(s):  
EN292-1:1991, EN292-2:1991, EN418:1992, EN954-1:1996, EN1050:1996, EN60204-1:1997.

**EMC-Directive** 89/336/EEC, 92/31/EEC and 93/68/EEC according to the following standard(s):  
EN50081-1:1994, EN50082-1:1997.

**Low Voltage Directive** 73/23/EEC and 93/68/EEC according to the following standard(s):  
EN60204-1:1997.

**Supplementary Information** The equipment also complies with the following standards:  
IEC204-1, UL508:1988, UL508A:1993, NFPA70:1996, NFPA79:1997  
and FCC part 15, subpart B.

The above has been declared according to the global method, module A

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# Exotom-100

**Automatic cut-off  
machine for cutting of  
workpieces up to 160 mm  
in diameter.  
Designed for the  
production environment**

## **A revolution in cutting technology**

- **Extremely simple operation**  
*Just two operating keys*  
*No programming required*
- **Unprecedented capacity and accessibility**
- **Up to four different cutting modes**

### **Unique technologies and features:**


- ExciCut** - easy through the hardest workpieces
- OptiFeed** - intelligent feed control
- AxioCut** - cutting mode module
- Variable cutting offset** - flexible wheel positioning
- Three operation modes** - block non-essential controls

Effortless operation even wearing working gloves

### ExciCut -for cutting the hardest materials

The Exotom-100 comes standard with the unique ExciCut feature. The technology was first introduced by Struers a decade ago, and appears on the

Exotom-100 in a highly refined version. Through an oscillating movement of the cut-off wheel, the contact area between the wheel and the workpiece is kept at a minimum throughout the cutting process.

ExciCut is ideal for hardened materials and means faster cutting of standard materials. Since the contact area between sample and cut-off wheel constantly changes, cooling is greatly improved. In turn, this means reduced risk of damage to the workpiece, and reduced work in the processes following the cutting.



### Direct Cut -the classic cutting mode

For smaller samples and some materials it might be an advantage to use the direct cutting mode where the cut-off wheel is fed through the workpiece in a linear movement.

### AxioCut -for extra large workpieces

For cutting of larger workpieces, additional 150 mm cutting capacity is available with AxioCut. This new feature allows cutting of workpieces with a dimension of e.g. 100 x 450 mm (for


other dimensions, please refer to Cutting Diagram on the back of this brochure). The field-upgradable option includes two different cutting modes: Step and Sweep.

Step cutting offers very fast cutting, whereas Sweep cutting offers fast cutting combined with better wheel economy.

### OptiFeed -Protection and speed

Cutting of workpieces with uneven cross-section or composition, represents a great challenge to the operator. The risk of damage to the workpiece or the cut-off wheel is high. Bearing this in mind, Struers has developed the new OptiFeed principle, based on a preset maximum force limit.


As the contact area (between workpiece and wheel, cf. below diagram) is getting larger, the


force increases to maintain the preset feed speed. If the preset force limit is reached, the feed speed is automatically reduced to avoid damage.

When the cut-off wheel has passed the centre of the workpiece, the contact area is getting smaller again, and with it the resulting force. In turn, the feed speed is automatically increased to the preset maximum.



Thanks to OptiFeed you will always cut through the workpiece at the fastest possible speed, within the set parameters - and without overloading the workpiece and the cut-off wheel.

Two flexible water jets for additional cooling. Ideal for hollow workpieces

*Single-level display:  
Full control of cutting  
process + motor information  
in same view*

### Controls of Unique Simplicity

Two simple keys combined with a large graphic display are incorporated on a swivelling control panel. With the multi-function knob all cutting parameters are easily controlled even wearing working gloves. The joystick is used for fast and precise positioning of the cut-off wheel. The side-mounted control panel is kept away from the machine for maximum safety and to avoid soiling of display and controls. Also the risk of slipping samples damaging the display is eliminated. The control panel can easily be moved if very long workpieces have to be placed in in the cutting chamber.

### Large Cutting Chamber

The corrosion resistant cutting chamber is exceptionally spacious with its 855 x 550 mm cutting table, offering unrestricted access. Bulky workpieces are easily placed manually or using a forklift truck, since the cutting table is the highest point in the chamber. Maximum mounting flexibility is provided by multiple, bi-directional 12 mm T-slots. The powerful, built-in flush hose provides for easy cleaning of the cutting chamber. Extra large windows and fully, indirect lightning, give an excellent, non-glare view of the cutting chamber. A balanced mechanism ensures effortless opening and closing of the hood.

### Variable Cutting Offset

The horizontal position (cutting offset) of the cut-off wheel can be set manually. In the ExciCut and Direct Cut cutting modes, the cut-off wheel may be moved 10 mm towards the front and 60 mm towards the

back of the cutting chamber. This allows for maximum flexibility when cutting difficult workpieces (e.g. turbine blades) and for clamping generally.

### Three Operation Modes

Particularly in an industrial environment with many users, operator errors may cause damage to machine and samples. To overcome this, non-essential controls may be blocked partly or fully by using the pass code-protected operation modes.

### High-capacity Cooling

The workpiece and cut-off wheel are effectively cooled via a large capacity (150 litres) cooling unit. For extreme cases, additional cooling may be added, using the two flexible water jets provided. This is particularly useful when cutting hollow workpieces. The cooling unit comes standard with an automatic emptying function and a

warning on low coolant level. A special filter tray is provided to collect cutting debris, greatly reducing the need for cleaning of the recirculation tank.

### AutoStop

After the workpiece is securely clamped and Start has been pressed you can leave the machine unattended. OptiFeed makes sure that cutting is completely under control. With AutoStop, Exotom-100 automatically detects when the workpiece has been cut through, returns the wheel to its start position and stops all motors.

### QuickPosition

Positioning the cut-off wheel can be sped up dramatically. Simply push the joystick downwards and the cut-off wheel will speed towards the workpiece. On contact with the workpiece, the cut-off wheel stops automatically without harming workpiece or wheel.

*Exceptionally spacious and accessible cutting chamber:  
Even the largest workpieces are easily placed on the cutting table*

Rubber curtains mean trouble-free cutting of very long workpieces. Move the swivelling control panel for optimum access

After contact, the wheel is automatically retracted a little from the workpiece, to be ready for cutting. QuickPosition can be used both before and during cutting.

### Cutting Long Workpieces

If the user wants to cut workpieces protruding the cutting chamber, it is possible to swivel the control panel aside, for easy access. The cutting chamber comes standard with side-mounted rubber curtains, for protruding workpieces.

### Programmable Return Position

The return position of the cut-off wheel can be set to three different positions: Exotom-100 retracts the cut-off wheel to the top position, to the position where the cutting was started or it lets the cut-off wheel stay in the cut. In this way, it is very easy to get convenient access to the workpiece or an extremely quick exchange when cutting uniform workpieces.

### Electronic Monitoring

More than a dozen electronic sensors monitor the complete cutting process and all machine functions to ensure flawless cutting. Never has sectioning on a large cut-off machine been this safe.

### Safety Features

Exotom-100 complies with international safety standards, and safety features include emergency stops and safety sensors for all machine functions and motors. The cutting motor cannot be started before the protection hood is closed, and after Start, the hood cannot be opened before the machine has stopped. The windows of the protection hood are made of impact resistant material.



### Cut-off wheels

A range of cut-off wheels for all ferrous and non-ferrous metals is available for Exotom-100.

Diamond cut-off wheels are available for cutting of ceramic and mineralogical materials and sintered carbides.

Application	Dimensions	Abrasive	Code
For ExciCut and direct cutting of extremely hard steels (350-800 HV)	432 x 3.0 x 32 mm	Al <sub>2</sub> O <sub>3</sub>	101MA
For ExciCut and direct cutting of steels (< 500 HV). Extremely hard (white cast iron) or ductile (18/8 or St60) steel types are cut with oscillating cutting. Gives good wheel economy	432 x 3.0 x 32 mm	Al <sub>2</sub> O <sub>3</sub>	102MA
For ExciCut and direct cutting of very hard and ductile ferrous metals. Strong, fibre-reinforced wheel (250-700 HV)	432 x 3.0 x 32 mm	Al <sub>2</sub> O <sub>3</sub>	202MA
For ExciCut and direct cutting of steels (< 500 HV). Gives very good wheel economy	432 x 3.0 x 32 mm	Al <sub>2</sub> O <sub>3</sub>	104MA
For ExciCut and direct cutting of non-ferrous metals (< 300 HV). Primarily for ExciCut	432 x 3.0 x 32 mm	SiC	106MA
For direct cutting of hardened steels	406 x 1.8 x 32 mm	CBN	108MA
For direct cutting of sintered carbides and ceramics	350 x 1.8 x 32 mm 305 x 1.8 x 32 mm	Diamond Diamond	26EXO 25EXO

## Technical data

Subject		Metric/International	US
<b>CUTTING SPECIFICATIONS</b>			
<b>Workpiece Dimensions (max.)</b>	<i>Workpiece inside cutting chamber:</i>		
	Height	250 mm	10"
	Width	650 mm	25.6"
	Depth	550 mm	20"
	<i>Workpiece protruding cutting chamber:</i>		
	Height	140 mm	5.5"
Depth	360 mm	14.2"	
<b>Cutting Capacity*</b>	Max. cutting diameter	160 mm	6.3"
	Rectangular sample:** (In brackets with AxioCut option)		
	Height	100 mm (100 mm)	3.9" (3.9")
	Depth	300 mm (450 mm)	11.8" (17.7")
<b>PHYSICAL SPECIFICATIONS</b>			
<b>Cutting Motor</b>	Cutting power	10.5 kW	14 HP
	Maximum power	15-18 kW	20-24 HP
<b>Cut-off Wheel</b>	Diameter x Thickness x Centre-hole	432 x 3 x 32 mm	17 x 0.12 x 1.26"
	Rotational speed (running idle)	1950 rpm	1950 rpm
<b>Positioning &amp; Feed</b>	Positioning range (of cut-off wheel)	0 – 250 mm	0 – 9.8"
	Max. positioning speed	50 mm/s	2.0"/s
	Feed Speed range (adjustable in steps of)	0.05 – 5 mm/s (0,05 mm/s)	0.002 – 0.2"/s (0.002"/s)
	Cutting force	Max. 700 N	Max. 150 lbf
<b>Cutting table</b>	Width	855 mm	33.6"
	Depth	550 mm	21.6"
	T-slots	12 mm	0.48"
<b>Dimensions and Weight</b>	Height (closed/open hood)	1756/2300 mm	69/90.5"
	Width (ex/including control panel)	1050/1350 mm	41.3"/53.1"
	Depth	1500 mm	59.1"
	Weight	820 kg	1810 lbs
<b>Recirculation Cooling Unit</b>	Tank volume	150 l	39.6 gallons
	Approx. flow	25 l/min	6.6 gallons/min
<b>ENVIRONMENT</b>			
<b>Noise level</b>	78 dB(A) measured at idle running, at a distance of 1.0 m/39.4" from the machine.		
<b>ELECTRICAL DATA</b>			
<b>Voltage/Frequency</b>	200-240 V / Max. power 15 kW / Max. load 69 A		
	380-480 V / Max. power 18 kW / Max. load 34 A		

\*Projected cutting capacity with new cut-off wheel. Actual cutting capacity depends on sample material, cut-off wheel and clamping technique

\*\*Example of cutting capacity. Please refer to Cutting Diagram above

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice.

Specifications	Code
<b>Exotom-100, AxioCut</b> Automatic cut-off machine with OptiFeed, AxioCut, ExciCut and Direct Cut. 10.5 kW (14 HP) motor. For 432 mm cut-off wheels. With recirculation cooling unit (capacity 150 l) and cutting table with 12 mm T-slots. Quick-clamping devices (EXOLE + EXORI) and vertical clamping system (EXOVS) respectively are ordered separately.	EXOCM
<b>Exotom-100</b> Automatic cut-off machine with OptiFeed, ExciCut and Direct Cut. 10.5 kW (14 HP) motor. For 432 mm cut-off wheels. With recirculation cooling unit (capacity 150 l) and cutting table with 12 mm T-slots. Quick-clamping devices (EXOLE + EXORI) and vertical clamping system (EXOVS) respectively are ordered separately. Field-upgradable to AxioCut by optional kit (EXOMA)	EXOCD
Accessories	Code
<b>AxioCut Kit for Exotom-100</b> Increases cutting capacity (depth) by 150 mm. Requires installation by Struers Service Technician.	EXOMA
<b>Replacement Stainless Steel Bands</b> For Cutting Table with 12 mm T-slots. Available in three sizes, each supplied in sets of two: 100 x 225 mm 100 x 300 mm 60 x 548 mm	EXOBS EXOBM EXOBL
<b>Quick-clamping Device for 12 mm T-slots, Left</b> For securing the workpiece. Complete with backstop. To be mounted on a cutting table and placed to the left of the cut-off wheel.	EXOLE
<b>Quick-clamping Device for 12 mm T-slots, Right</b> For securing the workpiece. Complete with backstop. To be mounted on a cutting table and placed to the right of the cut-off wheel.	EXORI
<b>Clamping Tools for 12 mm T-slots.</b> For clamping irregularly shaped workpieces. Complete with clamps, supports and bolts.	MAGOF
<b>Vertical Clamping System for 12 mm T-slots</b> For clamping irregularly shaped workpieces. Clamping height up to 92 mm. Complete with operating key and one flat clamping shoe.	EXOVS
<b>Swivel Shoes for Vertical Clamping System.</b> For clamping irregularly shaped workpieces. Set of four multi-shaped swivel shoes.	TREVI
<b>Riser Block for 12 mm Vertical Clamping System</b> For elevating Vertical Clamping System when clamping high workpieces. To increase clamping height by 74 mm.	EXOKS
<b>Arm Extension for 12 mm Vertical Clamping System</b> For extending arm of Vertical Clamping System.	EXOAR
<b>Chain Spanner</b> For 12 mm T-slots. Complete with anchor block and operating key.	EXOCS
<b>Spare Nylon Filter for Recirculation Cooling Unit</b> For insertion into filter-drawer in cooling unit compartment.	EXOFI
<b>Bandfilter connection kit</b>	EXOCO

## Clamping Tools

Struers offer a wide range of clamping tools for Exotom-100. Clamping requirements differ greatly, which is why Struers have chosen to supply Exotom-100 without any 'standard' clamping devices. Rather, the user can order exactly the clamping system matching his requirements. Please refer to separate brochure.

*Large cylindrical and uneven samples (e.g. crankshafts) are easily clamped with the Struers chain spanner*

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#### Struers Ltd.
#### Struers Inc.