LaboPress-1

Instruction Manual

Manual No.: 15057001

Date of Release 11.01.2008

LaboPress-1
Instruction Manual

Table of Contents	Page
User's Guide	1
Reference Guide.....	12
Quick Reference Card	22

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2008.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

LaboPress-1

Safety Precaution Sheet

To be read carefully before use

1. The operator should be fully aware of the use of the machine according to the Instruction Manual.
2. The machine must be placed in a well ventilated room on a working table with adequate height for convenient operation.
3. Be sure that the actual voltage corresponds to the voltage stated on the back of the machine and on the heating/cooling unit. The machine must be earthed.
4. Be sure that the water connections are mounted correctly and without leaks. The main water supply should be turned on when the machine is in use.
Struers recommend that the mains water supply is shut off or disconnected if the machine is to be left unattended
5. Be sure that the outlet hose is safely attached to the water outlet system.
6. Be sure that the mounting unit is correctly assembled on the press before starting the process.
7. Be sure that the top closure with upper ram is either correctly mounted on the mounting cylinder or completely removed from the mounting cylinder before starting the press.
8. Do not operate the mounting press with a higher force/pressure than recommended for the actual cylinder diameter and resin in Struers Application Guide for Hot Mounting.
9. Following a heating cycle, ensure the mounting cylinder is cooled for a minimum of two minutes before opening.
10. Do not operate the machine whilst assembling or disassembling the mounting unit.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

User's Guide

Table of Contents	Page
1. Getting Started	
Checking the Contents	2
LaboPress-1.....	2
Mounting Unit	2
Unpacking LaboPress-1.....	2
Placing LaboPress-1	2
Transport Screw Cap	2
Getting Acquainted with LaboPress-1.....	3
Noise Level.....	4
Supplying Power	4
Supplying Water	4
Water Inlet.....	4
Water Outlet.....	4
Recirculating Cooling Unit.....	4
Assembling the Mounting Unit.....	5
Removing the Cover	5
Installing the Lower Ram	5
Installing the Mounting Unit.....	5
Installing the Cover.....	6
Installing the Dust Protection Ring	6
Disassembling the Mounting Unit.....	6
Dropping the Lower Ram	6
Removing the Dust Protection Ring	6
Removing the Cover	6
Removing the Mounting Unit	7
Removing the Lower Ram	7
Changing the Mounting Unit.....	7
2. Basic Operations	
Using the Controls.....	8
Front Panel Controls of LaboPress-1	8
Main Switch	8
Placing the Specimen	10
Pouring Resin over the Specimen.....	10
Placing Two Specimens	10
Installing the Top Closure.....	11
The Mounting Process	11
Heating On.....	11
Heating Off	11
Removing the Top Closure	11

1. Getting Started

Checking the Contents

LaboPress-1

In the packing box you should find the following parts:

- 1 LaboPress-1 machine
- 2 Measuring spoons for mounting resin
- 1 Funnel
- 1 Air filter
- 1 Set of Instruction Manuals

Mounting Unit

- 1 Mounting unit
- 1 Top closure with upper ram
- 1 Lower ram
- 1 Dust protection ring
- 1 Piston pin
- 1 Mould release agent (FASTI)
- 1 Scraper (PROAN)


Unpacking LaboPress-1

LaboPress-1 is detached from the bottom of the packing case by removing the four screws from below.

Placing LaboPress-1

LaboPress-1 should be placed on a steady table with an adequate working height. The machine must be placed close to the power supply, water mains and water outlet facilities.

Transport Screw Cap


Carry out the following procedure before using the machine for the first time. Exchange the transport screw cap ① with the enclosed air filter ②, to equalise the pressure in the hydraulic system.

The transport screw cap ① is situated underneath the cover for the mounting unit.

**Getting Acquainted with
LaboPress-1**

Take a moment to familiarise yourself with the location and names of the LaboPress-1 components.

- ① Front panel
- ② Top closure
- ③ Main switch
- ④ Mounting unit
- ⑤ Dust protection ring
- ⑥ Cover for mounting unit
- ⑦ Top closure holder

Noise Level

The noise level of the machine is 60 dB (A) measured when the pump is running, at a distance of 1.0 m/39.4" from the machine.

Supplying Power

LaboPress-1 is factory mounted with an electric cable. Mount a plug on the cable:

Brown: phase

Blue: neutral

Yellow/green: earth


IMPORTANT

Check that the mains voltage corresponds to the voltage stated on the type plate on the machine. Do not operate the machine before the transit screw cap has been exchanged with the enclosed air filter.

Supplying Water

Cooling water can be supplied from the water mains or a Recirculating Cooling Unit (available separately).

Water Inlet

Connect the pressure hose to the water mains tap for cold water:

- Mount the reduction ring with gasket on the water mains tap, if necessary.
- Insert the gasket and tighten the coupling nut completely.

- ① Water inlet
- ② Water outlet tube

IMPORTANT

Only connect to cold water

Water Outlet

- Lead the water outlet tube (pos. 2) to the drain and be absolutely sure the hose has a steady downward slope and is without any obstructions.
- Attach the water outlet tube to the water outlet system.

Recirculating Cooling Unit

An optional cooling water recirculation unit is available for users with special environmental or other needs that demand cooling water be recirculated. (See Accessories section.)


When connected to LaboPress-1 it must be operated manually, using the mains power switch to start and stop the cooling process.

Separate directions for installing and using the cooling water recirculation unit are provided with the unit.

Assembling the Mounting Unit

Removing the Cover

- Press the catch (pos. 4).
- Press gently on both sides of the cover to release it, and lift it up.

④ Catch


Installing the Lower Ram

- Place the lower ram on the top of the piston rod.
- Turn the lower ram so that the holes in the axle journal on top of the rod and ram line up, and insert the piston pin.
- Make sure that the ends of the pin do not protrude.

IMPORTANT

The axle journal is part of the machine's safety mechanism to protect against damage to the machine. Please see: Maintenance. If it breaks, it must only be replaced with a Struers replacement part or the safety mechanism may not function.


Installing the Mounting Unit

- ① Mounting cylinder
- ② Fixation screw
- ③ Upper quick coupling for cooling water
- ④ Lower quick coupling for cooling water
- ⑤ Plug

IMPORTANT


Do not operate the machine whilst installing the mounting unit

- Unscrew the fixation screw about 10 mm.
- Place the mounting cylinder over the lower ram with the fixation screw in the position shown.
- Turn the mounting unit in a clockwise direction until it stops.
- Tighten the fixation screw completely.
- Mount the tube with the straight quick coupling on the lower quick coupling (position 4) of the mounting unit and push to connect. Make sure that the ring comes to stop in the extreme end of the quick coupling.
- Mount the tube with the elbow quick coupling on the upper quick coupling of the mounting unit (pos. 3) and push to connect. Make sure that the ring comes to stop in the extreme end of the quick coupling.
- Mount the plug (pos. 5) in the socket. Tighten the coupling nut.

Installing the Cover

- Insert the barb into the front of the opening for the cover.
- Insert the barbs on both sides of the cover into the cabinet. Press the sides of the cabinet gently.
- Press the rear end of the cover in so that the catch engages.

④ Catch

Installing the Dust Protection Ring

Place the dust protection ring around the mounting cylinder. The concave side should face upwards.

Disassembling the Mounting Unit

Dropping the Lower Ram

- Switch on the main power of the machine.
- Press the key RAM DOWN ▼ to lower the ram to its lower limit.


Removing the Dust Protection Ring

Lift out the dust protection ring.

Removing the Cover

- Press the catch (position 4).
- Press gently on both sides of the cover, to release it, and lift it up.

Removing the Mounting Unit

- ① Mounting cylinder
- ② Fixation screw
- ③ Upper quick coupling for cooling water
- ④ Lower quick coupling for cooling water
- ⑤ Plug

IMPORTANT

Do not operate the machine whilst removing the mounting unit

- Disconnect the plug (pos. 5) from the socket. Loosen the coupling nut and pull.
- Disconnect the tube with elbow quick coupling on the upper quick coupling on the mounting unit (pos. 3). Pull the ring. Wait 5 seconds to allow the water to flow out of the cooling coil.
- Disconnect the tube with straight quick coupling on the lower quick coupling on the mounting unit (pos. 4). Pull the ring.
- Unscrew the fixation screw about 10 mm.
- Turn the mounting unit in an anti-clockwise direction until it stops.
- Lift the mounting unit.

Removing the Lower Ram

- Push the piston pin out of the lower ram and lift it off.

IMPORTANT

The axle journal is part of the machine's safety mechanism to protect against damage to the machine. Please see: Maintenance. If it breaks, it must only be replaced with a Struers replacement part or the safety mechanism may not function.

Changing the Mounting Unit

Follow the instructions for Removing the Mounting Unit and Installing the Mounting Unit.


2. Basic Operations

Using the Controls

Front Panel Controls of
LaboPress-1







Main Switch

- ① Front panel
- ③ Main switch

LaboPress-1
Instruction Manual

Pos. No.	Key/Knob	Function	Pos. No.	Key	Function
① MAIN SWITCH		Turns the main power of the machine on/off. The main switch is located to the right of the cabinet.	⑧ HEATING OFF		Switches off the heating unit.
② POWER		Lights up when the main power is on.	⑨ HEATING		Selects a temperature of 180°C. The heating unit is <i>not</i> switched on.
③ RAM UP		Starts the upwards movement of lower ram. Creates and maintains the pressure in the mounting cylinder.	⑩ HEATING INDICATOR		Lights up when the temperature reaches 180°C.
④ RAM DOWN		Starts the downward movement of the lower ram. Relieves the pressure in the mounting cylinder.	⑪ HEATING		Selects a temperature of 150°C. The heating unit is <i>not</i> switched on.
⑤ FORCE		Gauge indicating the actual force (pressure).	⑫ HEATING INDICATOR		Lights up when the temperature reaches 150°C.
⑥ HEATING ON		Switches on the heating unit.			There is no specific cooling indicator on LaboPress-1. Cooling is controlled by the water mains tap.
⑦ HEATING INDICATOR		Lights up when heating unit is on.			

Placing the Specimen

- Press and hold RAM UP ▲ to raise the lower ram to its upper limit. NOTE: The ram has reached its upper limit when the force needle rises sharply.
- Apply "*Mould Release Agent*" to the surface of the lower ram.
- Place the specimen on the ram. The specimen must be clean, dry and free from grease. The distance between the specimen and the cylinder wall must be minimum 3 mm to avoid cracks in the resin.

Pouring Resin over the Specimen

- Press and hold RAM DOWN ▼ to lower the ram to its lower limit. NOTE: The ram has reached its lower limit when the force needle falls to zero.
- Pour a suitable amount of resin into the cylinder by means of the enclosed funnel.
- **Always** make sure there is sufficient resin to cover the sample after compression. Please note that the volume of the resin reduces when the granulate becomes compressed. If insufficient resin is used, the rams may come in contact with the sample, and the cylinder may be damaged.


IMPORTANT

The "*Mould Release Agent*" must always be applied to the mounting rams as a thin layer to prevent the resins sticking to the surface. By means of Struers AntiStick a thin layer of stearate powder can easily be dabbed on the rams.

Placing Two Specimens

- Follow the instructions in Placing the Specimen.
- Apply "*Mould Release Agent*" to all surfaces of the intermediate ram.
- Place the intermediate ram (available separately) on top of the resin.
- Place the second specimen on the intermediate ram and fill the cylinder with an adequate quantity of resin.
- Again, make sure you use adequate resin!

Installing the Top Closure

- Remove resin dust from the upper part of the mounting cylinder.
- Clean the cylindrical surface of the upper ram. Cured resin can easily be removed without damage to the surface of the ram using the scraper supplied.
- Apply "Mould Release Agent" to all accessible surfaces of the upper ram.
- Place the top closure with the upper ram on the mounting cylinder.
- Press the top closure straight down, turning it counter-clockwise until you hear a click.
- Press the top closure down and turn it clockwise until its lower limit.
- Turn the top closure a quarter of a turn back.

IMPORTANT

If the ram does not fit in the cylinder check ram and cylinder for cured resin. The tolerance between the cylinder and the ram is very fine and even small amounts of resin from previous mounting may cause problems.

The Mounting Process

Heating On

- Press HEATING **180°C** or **150°C** to select the temperature.
- Press HEATING ON ☼ to switch on the heating unit.
- Press RAM UP ▲ until the force, indicated on the gauge, is as desired. If the force should become too high: Press RAM DOWN ▼ briefly.
- Maintain the force (pressure) during the heating period by pressing the appropriate RAM UP/DOWN ▲▼ key.

Heating Off

- When the heating period is finished, press HEATING OFF ☼ to switch off the heating unit.
- Open the cooling water tap on the water mains.
- Maintain the force (pressure) during the cooling period.
- When the cooling period is finished, close the cooling water tap on the water mains.

If you have stopped the machine during the mounting process:
Cool the mounting cylinder a minimum 2 min before opening, after a heating period. Please note that the mount might be destroyed.

Removing the Top Closure

When the mounting process is finished:

- Press RAM DOWN ▼ briefly to relieve the pressure.
- Turn the top closure counter-clockwise until released from the thread.
- Press RAM UP ▲ to raise the raise ram to its upper limit.
- Place the top closure in the top closure holder.

Reference Guide

Table of Contents	Page
1. Struers Metalog Guide™	13
2. Application Guide for Hot Mounting	13
3. Accessories	14
4. Consumables	15
5. Trouble-Shooting	16
6. Maintenance	
Daily Service	18
Monthly Service	18
Cleaning under the Lower Ram	18
Cleaning the Rams	18
Removing the Upper Ram	19
Lubricating the Threads for the Top Closure	19
Decalcifying the Cooling Coil in the Mounting Unit	19
7. Technical Data	20

1. Struers Metalog Guide™

Struers Metalog Guide™ offers preparation methods for the most common materials, based on a simple analysis of two key properties: hardness and ductility. Finding the right method, including choice of consumables is easy. Always consult Struers Metalog Guide™ for the correct preparation method for the actual specimens.

Struers Metalog Guide™ contains 6 useful chapters:

Metalogram

A quick and safe guide to the right preparation method.

Metalog Method:

A complete catalogue of preparation methods, based on Struers' extensive materialographic experience, and employing Struers' range of consumables.

Preparation Philosophy

The basics of modern specimen preparation, seen from a professional point of view.

Metalog Process

The materialographic preparation process from start to finish, logically explained.

Metalog Master

A combined trouble-shooting guide and supply of in-depth information on the processes of mechanical preparation, including an expert system for the solving of preparation problems.

Metalog Code

Quick access to the relevant consumables for the chosen preparation methods.

Struers Metalog Guide™

A complete guide to materialographic specimen preparation. Contact your local dealer for a free copy of Metalog Guide™.

2. Application Guide for Hot Mounting

Useful mounting data and hints can be found in the Struers [Application Guide for Hot Mounting](#) or visit the Struers website on <http://www.struers.com> and see the section on **Knowledge**.

3. Accessories

Type	Description	Code
Mounting Unit consisting of cylinder, upper and lower ram, heating/cooling unit and top closure	25 mm dia. 30 mm dia. 1 ¼" dia. 1 ½" dia. 40 mm dia. 50 mm dia.	CYNMI CYNTE CYNQU CYNHA CYNIF CYNAM
Recirculation Cooling Unit (must be operated manually on LaboPress-1)		TRECI
Intermediate Ram for making 2 simultaneous mounts	25 mm dia. 30 mm dia. 1 ¼" dia. 1 ½" dia. 40 mm dia. 50 mm dia.	CYLDO CYLDO CYLDO CYLDO CYLDO CYLDO

4. Consumables

Resin	Application	Specific property	Material	Contents	Code
ConduFast (Resin 1)	Electropolishing	Electrically conductive Very low shrinkage	Acrylic resin with iron powder filler Thermoplastic, Grey	1 kg	RESFE
SpeciFast (Resin 3)	Glass clear mounts Porous samples Surface electrical insulator for ConduFast	Transparent Medium shrinkage	Acrylic Thermoplastic	1.5 kg 7.5 kg 25 kg	RESTH RESYV RESFA
PolyFast	Edge retention	Very low shrinkage Medium removal rate	Phenolic resin with carbon filler Thermosetting, Black	1 kg 7.5 kg	FAPSA FAPME
	Examination in scanning electron microscope	Electrically conductive Low emission in the vacuum chamber			
IsoFast (Resin 4)	Edge retention Planeness	Low shrinkage Low removal rate	Dialylphthalate with glass fibre filler Thermosetting, Green	1 kg 7.5 kg	RESDI RETTY
DuroFast (Resin 5)	Edge retention Planeness	Very low shrinkage Good adhesion Very low removal rate	Epoxy with mineral filler Thermosetting, Black	1 kg 7.5 kg	RESIF RESEN
MultiFast	Routine examination, Backup resin	Low shrinkage Medium removal rate	Phenolic resin with wood flour filler Thermosetting, Black	2.5 kg 7.5 kg 25 kg	RESRU RESUK RESLA
MultiFast Green	Routine examination Colour marking	Low shrinkage Medium removal rate	Phenolic resin with wood flour filler Thermosetting, Green	2.5 kg 7.5 kg 25 kg	FAGSA FAGME FAGLA
MultiFast Brown	Routine examination Colour marking	Low shrinkage Medium removal rate	Phenolic resin with wood flour filler Thermosetting, Brown	2.5 kg 7.5 kg 25 kg	FABSA FABME FABLA

Product	Application	Specific property	Material	Contents	Code
Pre-Mount	Serial mounting of uncomplicated shapes	Easy to handle Low shrinkage Medium removal rate	Phenolic resin with wood flour filler In tablet form Thermosetting, Black	25 mm/ 1" dia. approx. 1100 pcs. 10 kg	RESUN
				30 mm/ 1.1/4" dia. approx. 760 pcs. 10 kg	RESGI
				40 mm/ 1 1/2" dia. approx. 460 pcs. 10 kg	RESAM

Product	Application	Material	Contents	Code
Mount Release Wax	To stop mounts sticking to the rams	Wax	24 g	PREWA
AntiStick	To stop mounts sticking to the rams	Very fine stearate powder	8 g	FASTI

5. Trouble-Shooting

Error	Cause	Action
Insufficient compression.	Incorrect force/pressure.	Set the correct parameter.
	Defect in the pressure system.	Call a Struers service technician.
Insufficient heating.	Incorrect temperature setting. Heating time too short.	Set the correct temperature. Increase the heating time.
	Defect in the heating system.	Call a Struers service technician.
Insufficient cooling	Cooling time too short.	Increase the cooling time.
	Mains tap supplying cooling water is either closed or not opened sufficiently. Filter at the water inlet is blocked. The water is too hot.	Open the tap. Clean the filter. Use only cold water.
	Insufficient water within the Recirculation Cooling Unit. The water is too hot.	Fill to the correct water level. See the section: 'Maintenance'.
	Defect in the cooling system.	Call a Struers service technician.
Cooling water drips underneath the machine.	The quick coupling is not mounted correctly.	Remove the cover from the mounting unit and check the quick coupling connections.
The top closure will not start threading on the mounting cylinder.	The top closure is not mounted correctly.	Press the top closure straight down, turning it counter-clockwise until you hear a click. Then turn the top closure clockwise.
	The upper ram is too hot.	Let the top closure and the upper ram cool down. Reduce the mounting temperature
	Cured resin inside the mounting cylinder.	Clean the mounting cylinder with a brass wire brush.
	Cured resin on the cylindrical surface of the upper ram.	Clean the ram with the scraper supplied.
	The top closure has been dropped on the floor/table, causing a dent in the edge of the ram.	Replace the upper ram.
	Damage to the threaded connection in the top closure or to the upper ram.	Remove the upper ram from the top closure (see instructions in chapter Maintenance). Try both the following procedures to identify the problem: <ul style="list-style-type: none"> ▪ Try to mount the top closure, without the upper ram. If this is not possible, call a Struers service technician. ▪ Try to insert the disconnected upper ram into the mounting cylinder. If this is not possible, call a Struers service technician.

LaboPress-1
Instruction Manual

Error	Cause	Action
The top closure cannot be screwed down completely	Dirt in the threads on the top closure and the mounting cylinder.	Clean the threads. Use only a dry lubricant powder.
	The disc for thermal insulation, placed on the top of the upper ram, has a larger diameter than the upper ram.	Call a Struers service technician.
The top closure cannot be loosened.	<p>The top closure has not been turned ¼ turn back before starting the process.</p> <p>Cured resin on the cylindrical surface of the top ram.</p> <p>Dirt in the threads of the top closure.</p>	<p>To release the top closure: Move the lower ram up and down several times.</p> <p>If this does not help: -Put heating on for 1 min.</p> <p>If this does not help: - Set the force or pressure to zero. - Set the heating time and cooling time to 15 min. - Complete a mounting process.</p> <p>If this does not help: - Remove the two handles on the top closure. - Remove the plastic cover from the top closure (See Maintenance). -Loosen the top closure with a fork spanner.</p>

6. Maintenance

Daily Service

- Clean all accessible surfaces with a moist cloth.


Monthly Service

Cleaning under the Lower Ram

When mounting, some old cured resin will fall down from the lower ram and build up underneath. The cured resin must be cleaned away, to stop it harming the machine. The axle journal securing the lower ram has been designed to break in this case, thus preventing serious damage to the machine.)

- Remove the mounting unit.
- Remove the resin under the lower ram with a cloth or a soft brush.
- Install the mounting unit.

Cleaning the Rams

Before each mounting, check that there is no resin left on the flat surface of the rams from previous operations. Clean the cylindrical surface of the upper ram and optional intermediate ram, as this may cause problems with travel in the cylinder and mounting the top closure. Cured resin can easily be removed without any damage to the surface of the rams using the scraper supplied.


If the rams have been scratched on the polished surfaces, disassemble and polish them with 6 or 3 μm diamonds on a polishing cloth.

IMPORTANT

Do not polish the surface on the sides of the rams. If the sides of the rams are heavily scratched, they must be scrapped.

Removing the Upper Ram

- Take the cover plugs out of the ends of the handles and unscrew them with an allen key. (pos. 1).
- Remove the plastic cover (pos.2).
- Lift the upper ram and remove the lock washer (pos.3).
- Pull out the upper ram. Do not remove the rod (pos. 4) from the upper ram unless absolutely necessary.
- If it is necessary to remove the rod, hold the upper ram in a vice or a similar device. The ram surface must be protected with plastic or soft metal.

- ① Handle
- ② Plastic cover
- ③ Lock washer
- ④ Rod

Lubricating the Threads for the Top Closure

In case there is any friction in the threads, they should be lubricated with a dry lubricant such as MoS₂ powder or graphite.

IMPORTANT

Please observe that the threads in the joint between the top closure and the mounting cylinder must not be lubricated with oil or grease. As the mounting cylinder operates at high temperatures, only a dry lubricant, resistant to high temperatures, can be used.

Decalcifying the Cooling Coil in the Mounting Unit

When using cooling water from the mains tap in areas with a high chalk or mineral content, deposits can build up in the cooling coil. This may reduce the cooling effect, so once a year the cooling coil should be decalcified.

- Remove the mounting unit. (See this section.)
- Flush the coil with a mild decalcifying fluid, such as that used for coffee machines.
- Flush the cooling coil with clean water.
- Re-install the mounting unit.

7. Technical Data

Subject		Specifications	
		Metric/International	US
Mounting Units (Optional)	Diameter	25, 30, 40, 50 mm	1 ¼", 1½"
Compression	Force on piston rod	Variable between 0 and 50 kN	Variable between 0 and 12000 Lbf
Heating	Temperature	150/180 °C	302/356°F
Water Supply	Tap Water or Recirculation Cooling Unit		
	Pressure for tap water	1 - 10 bar	14.5 - 145 psi
	Inlet	1/2" or 3/4"	1/2" or 3/4"
Supply Voltage	Power phases	1	
	Power consumption:		
	Max. Total Mounting Unit	1500W 750W	
	Voltage/frequency:	Max. continuous load	
	100-120V / 50Hz 100-120V / 60Hz 220-240V / 50Hz 220-240V / 60Hz	10A (100V) - 12A (120V) 10A (100V) - 12A (120V) 4.5A (220V)- 5.0A (240V) 4.5A (220V)- 5.0A (240V)	
Dimensions and Weight	Width	420 mm	
	Depth	580 mm	
	Height (Mounting Unit and Top Closure installed)	420 mm	
	Weight	38 kg	

Quick Reference Card

Placing the Specimen

- Press and hold RAM UP ▲ to raise the lower ram to its upper limit.
- Apply "Mould Release Agent" to the surface of the lower ram.
- Place the specimen on the ram.

Pouring Resin over the Specimen

- Press and hold RAM DOWN ▼ to lower the ram to its lower limit.
- Pour a suitable amount of resin into the cylinder using a funnel. Use sufficient resin to allow for shrinkage.

Installing the Top Closure

- Remove resin dust from the upper part of the mounting cylinder.
- Clean the cylindrical surface of the upper ram. Cured resin can easily be removed without damage to the surface of the ram using the scraper supplied.
- Apply "Mould Release Agent" to all accessible surfaces of the upper ram.
- Place the top closure with the upper ram on the mounting cylinder.
- Press the top closure straight down, turning it counter-clockwise until you hear a click.
- Press the top closure down and turn it clockwise until it reaches its lower limit.
- Turn the top closure a quarter of a turn back.

The Mounting Process

Heating On

- Press HEATING **180°C** or **150°C** to select the temperature.
- Press HEATING ON ☀ to switch on the heating unit.
- Press RAM UP ▲ until the force (pressure), indicated on the gauge, is as desired.
- Maintain the force during the heating period.

Heating Off

- Press HEATING OFF ☀ to switch off the heating unit.
- Open the cooling water tap on the water mains.
- Maintain the force (pressure) during the cooling period.
- When the cooling period is finished, close the cooling water tap on the water mains.

Removing the Top Closure

When the mounting process is finished:

- Press RAM DOWN ▼ shortly to relieve the pressure.
- Turn the top closure counter-clockwise until released from the thread.
- Press RAM UP ▲ to raise the lower ram to its upper limit.
- Place the top closure in the top closure holder.

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Spare Parts and Diagrams

LaboPress-1/-3
Spare Parts and Diagrams

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2005.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Telefax: +45 44 600 801

Spare Parts and Diagrams

Table of contents

Drawing

Common Parts for LaboPress-1 and LaboPress-3

Drawings

Bottom, Complete	15050050R
LaboPress-1, cabinet, complete.....	15050001C
Cylinder Unit	15360036C
Top Closure	15360026A

LaboPress-1

Diagrams

Hydraulic Diagram	15051001D
Wiring Diagram	15053450E

LaboPress-3

Diagrams

Hydraulic Diagram	15081004D
Wiring Diagram	15083450F

Please note

Some of the drawings may contain position numbers
not used in connection with this manual.

LaboPress-1/-3
Spare Parts and Diagrams




Drawing	Pos.	Spare Part	Cat No.
15050050		Bottom, complete	
	30+40	Rubber foot with screw and spacer	R4910140
	145	Fuse 3.15A T, glass tube	359MP109A
	350	Axle journal	15360228
	450	Coupling 45°	2NH11445
	510	Straight coupling with hose nipple	2NH10014
		Air Filter (lid for hydraulic oil)	2YF90005
15050001		LaboPress-1, cabinet, complete	
	50	Top cabinet 505	15050115
	70	Top closure holder	15050385
		Various parts	
		Funnel ø80/15	2GP38508
		Measuring spoon	50300000
		Allan Key 5	2GR00050
		Pipe Wrench 13mm	2GR00006
		Timer clock	62000025
		LaboPress-3 Only	
		Pressure water hose	2NU93020
		Gasket with filter ¾"	2IX20410
		Gasket ø11/ø24x1.5	13590359
		Reducing piece ½" ¾"	2NG30013
	15360036		Cylinder Unit ø25
155		Cover plate ø25	15368005
190		Fixation screw	15050040
		Lower ram, ø25	15360200
		Piston pin, for ø25 ram	15360215
		Dust protection ring	15360313
		Scraper (PROAN)	05366901

LaboPress-1/-3
Spare Parts and Diagrams

Drawing	Pos.	Spare Part	Cat No.
15360026		Top closure ø25	
		Top closure ø25, complete	15360020
	30	Insulation, disc for ø25 top closure	15360300
	40	Pressure spring ø12.5x1.0	2GF10125
	50	Top part of ø25 top closure	15360170
	60	Lock washer	15050241
	70	Insulation shell for ø25 top closure	15360180
	80	Handle for top closure	15360027
15360036		Cylinder unit ø30	
	155	Cover plate ø30/1¼"	15368006
		Lower ram ø30 for ø7	15360201
		Piston pin, ø7	15360216
15360026		Top Closure ø30	
		Top Closure ø30, complete	15360021
	30	Insulation disc for ø30/ø1¼" top closure	15360301
	50	Top part of ø30/1¼" top closure	15360171
	70	Insulation shell for ø30/ø1¼" top closure	15360181
15360036		Cylinder unit ø 1 ¼"	
		Lower ram, ø7	15360202
		Piston pin, ø7	15360216
15360026		Top closure ø 1 ¼"	
		Top closure ø1¼", complete	15360022
15360036		Cylinder unit ø 1 ½"	
		Lower ram ø 1 ½"	15360203
	155	Cover plate ø40/1½"	15368007
		Piston pin, ø7	15360217
15360026		Top closure ø 1 ½"	
		Top closure ø1½", complete	15360023
	30	Insulation disc for ø40/ø1½"/ø50 top closure	15360302
	50	Top part of ø40/1½" top closure	15360172
	70	Insulation shell for ø40/ø1½" top closure	15360182

LaboPress-1/-3
Spare Parts and Diagrams



Drawing	Pos.	Spare Part	Cat No.
15360036		Cylinder unit ø40	
		Lower ram ø40	15360204
		Piston pin, ø7	15360217
15360026		Top closure ø40	
		Top closure ø40, complete	15350024
15360036		Cylinder unit ø50	
	155	Cover plate ø50	15368008
		Lower ram ø50	15360205
		Piston pin, ø7	15360218
15360026		Top closure ø50	
		Top closure ø50, complete	15360025
	50	Top part of ø50 top closure	15360173
	70	Insulation shell for ø50 top closure	15360183

- 450 510 Smøres med Unisiikon TK M 1011
- 60 340 100M Omnifit
- 80 90 Limes med 409 GEL
- 90 High vacuum grease på 2 udv. kanter.
- 150 Opslides. Flap placeres inde i hydraulik-rum.
- 480 Placeres 360mm fra kobling.

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



1. Hydraulic Powerpack
2. Double-acting Hydraulic Cylinder
3. $\varnothing 5$ Polyamid Tube
4. Pressure Gauge
5. Valve Manifold
6. Relief Valve
7. Accumulator
8. Throttle Plug

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



ALL WIRINGS 0.75mm²
EXCEPT OTHERWISE MARKED.

COLOR CODES
BK = BLACK
BN = BROWN
RD = RED
GN = GREEN
BU = BLUE
GY = GREY
WH = WHITE

X4:8	W3	BK3	N110220V	POWER IN
X4:7	W3	BK4	M110220V	POWER IN
X4:6	W4	BK6	N110220V	HEATING ON
X4:5	W4	BK6	OUT	HEATING ON
X4:4	W5	BK1	N110220V	PUMP/MOTOR ON
X4:3	W5	BK2	IN	PUMP/MOTOR
X4:2	W5	BK2	OUT	PUMP/MOTOR
X4:1	W5	BK4	M110220V	THERMISTOR/FCH
X5:1	W4	RE3	0ND	
X5:1	W4	RE4	SIGNAL IN	

REV. A: X4:1 → X4:8 ... X4:8 → X4:1	3 PHASES AC VARS RES ALLS 170 DIC 25 TO FRODOVRE DRUMMERS P/N: 12050378#
REV. B NEW SYMBOL, Y1+Y2	
REV. C NEW T1 → NEW CONNECTION REV. D: W1 140220V → 120V AC REV. E: W6 IN X1:1 → X1:16 @ 19-02-99/CEM	
Thursday, June 08, 2000	Scale: 021095/CEM
Monday, October 02, 1995	Page Code: 15053450

1. Hydraulic Powerpack
2. Double-acting Hydraulic Cylinder
3. ø5 Polyamid Tube
4. Pressure Transmitter
5. Valve Manifold
6. Relief Valve
7. Accumulator
8. Throttle Plug

D: 1 stk. Throttle Plug
03-08-2004 MPE/JTV

Matr.: LaboPress-3

Overfl.beh.:

Målforhold:

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C: Y1,Y2 ændret.
270197 BMJ/

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A: symbol ved Pos 5
rettet. 200995 BMJ/

Struers

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Dato

Sign.


Tegn: 070895 BMJ

Kontr.

Hydraulic Diagram

Erst.:

15081004 D

MAINS CABLE
(North America)
100V-240V/50Hz
3xAWG 12
(WH) (BK) (GN)

100V-240V/50Hz
3x2.5mm²
YE/GN

ALL WIRINGS 0.75mm² EXCEPT OTHERWISE MARKED.
COLOR CODES:
BK = BLACK
BN = BROWN
RD = RED
YE = YELLOW
GN = GREEN
BL = BLUE
GY = GREY
WH = WHITE

SEE CONNECTIONS
DIAGRAM 15083150
WIRING FOR 120V AND 240V

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Herewith declares that

<i>Product Name:</i>	LaboPress-1
<i>Type No:</i>	505
<i>Machine Type:</i>	Mounting Apparatus

is in conformity with the provisions of the following directives:

Safety of Machinery 98/37/EEC according to the following standard(s):
EN ISO 12100-1:2003, EN ISO 12100-2:2003, EN60204-1:1997.

EMC-Directive 2004/108/EEC according to the following standard(s):
EN61000-6-1:2002, EN61000-6-3:2002.

Low Voltage Directive 2006/95/EEC according to the following standard(s):
EN60204-1:1997.

Supplementary Information The equipment complies with the American standards:
FCC part 15, subpart B and UL508.

The above has been declared according to the global method, module A

Date: 04.12.2007

Christian Skjold Heyde,
Vice President, R&D and Production, Struers A/S

Hersteller Struers A/S
Pederstrupvej 84
DK-2750 Ballerup, Danmark
Telefon +45 44 600 800

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Christian Skjold Heyde,
Vice- President, R & D et Production, Struers A/S

LaboPress-1

Semi-automatic mounting press for hot mounting of all materials. Equipped with low energy consuming heating/cooling unit

Semi-automatic press with electro-hydraulic pump

Instead of the tedious process of using a hand operated pump, simply press a button. LaboPress-1 is equipped with an electro-hydraulic unit. The force is displayed on a gauge and adjusted easily. When the pressure is stabilized after the resin has been compressed, the press will maintain this pressure due to a built in reservoir.

Contrary to traditional, hand operated mounting presses it is not necessary to stay at the press.

Saving energy

New designed mounting units produce more heat using less energy

Innovative design has merged the heating/cooling coils and the mounting cylinder into one compact and lightweight unit. This allows for a very rapid heat transfer and thus for a very fast heating and cooling.

The power consumption of the heating unit could be reduced to 750 W while the heating rate compared to conventional presses could be increased by about 100 %. This means that less than half the energy is sufficient to heat up the mounting unit in a far shorter time than a conventional press could do. After the temperature has been reached, the energy consumption will be reduced even more as the heating will be switched off and only be used to keep the temperature constant.


Two heating levels

On LaboPress-1, the temperature can be set to different levels. In all normal mounting situations 180°C will be used to get the fastest possible result. However, for heat sensitive materials or resins the temperature setting can be reduced to 150°C.

Easy exchange and operation

Mounting units up to 50 mm diameter

The new mounting unit is equipped with quick couplings for the water supply and a plug for the electrical connection, making it possible to exchange the complete unit in seconds. Six different cylinder diameters in metric and American units are available, accommodating for all possible applications.

Intermediate ram

The new mounting units are designed also to be used with an intermediate ram. The capacity of the press can then be doubled as two mounts can be produced at the same time.

New top closure

LaboPress-1 is equipped with a new top closure that combines the advantages of both bayonet and threaded closure.

Equipped with a triple external thread, the new closure will always "catch" like a bayonet and the threads will force the upper ram down into the cylinder, highly facilitating the work for the operator. As the threads are external, contamination with resin is rather unlikely.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup, Denmark
Phone +45 44 600 800
Fax +45 44 600 801
struers@struers.dk
www.struers.com

Technical Data		
	LaboPress-1	
Mounting unit	Mounting dia. Interchangeable	25, 30, 40, 50 mm and 1¼", 1½"
	Dual Mounts	Optional intermediate ram
Mounting parameters	Force range, kN	0-50
	Heating temperature °C	150 / 180
Electronics	Controls	Touch pads
	Read out during mounting	Actual force process Heating on Preset heating temperature
Safety standards conforming to: CE and FCC	General	89/392/EEC & 91/386/EEC & 93/44/EEC EN292-1, EN292-2, EN60204-1 (IEC 204-1)
	EMC	89/336/EEC & 92/231/EEC EN50081-2, EN50082-2 FCC part 15, subpart B
Supply	Electrical, no of phases	1
	Power consumption, W	
	Max. total Mounting unit	1500 750
	Water	Tap water Min 1 bar, max 10 bar
Dimensions and weight	Width x depth x height, mm (Mounting Unit with Top Closure installed)	420 x 580 x 420
	Weight, kg	38

Specification

LaboPress-1

Semi-automatic, electro-hydraulic mounting press. Mounting units in a range from 25-50 mm dia. are accommodated. The mounting unit is ordered separately

Mounting Unit for ProntoPress-20, ProntoPress-10, LaboPress-3 and LaboPress-1

Consisting of cylinder, upper and lower ram, heating/cooling unit and top closure.

	Code
Cylinder diameter 25 mm	CYNMI
Cylinder diameter 30 mm	CYNTE
Cylinder diameter 1¼"	CYNQU
Cylinder diameter 1½"	CYNHA
Cylinder diameter 40 mm	CYNIF
Cylinder diameter 50 mm	CYNAM

Intermediate ram

For making 2 mounts in ProntoPress-20, ProntoPress-10, LaboPress-3, LaboPress-1, Primopress, Prontopress-2, Predopress and Prestopress-3. Please state cylinder diameter

CYLD0

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice.

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Struers Ltd.
Struers Inc.