Labotom-3

Instruction Manual

Manual No.: 15327001

Date of Release 04.04.2011

Labotom-3
Instruction Manual

Table of Contents	Page
User's Guide	1
Reference Guide.....	15
Quick Reference	27

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

Original instructions. The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2011.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Labotom-3 Safety Precaution Sheet

To be read carefully before use

1. The operator(s) should be fully instructed in the use of the machine and its cut-off wheels according to the Instruction Manual and the instructions for the cut-off wheels.
2. The machine must be placed on a safe and stable support table. All the machines safety functions and guards must be in working order.
3. Use only intact cut-off wheels. The cut-off wheels must be approved for min. 50 m/s.
4. Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
5. The workpiece must be securely fixed in the quick-clamping device or the like. Large or sharp workpieces must be handled in a safe way.
6. Protruding workpieces should be shielded or marked.
7. The cutting handle should be lowered slowly and carefully, in order to avoid breaking the cut-off wheel.
8. Struers recommends the use of an exhaust system as the materials to be cut may emit harmful gasses or dust.
9. Do not open the protection guard before the cut-off wheel has stopped completely.
10. Do not put your hands through the rubber curtain during operation of the machine.

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

User's Guide

Table of Contents	Page
1. Getting Started	
Checking the Contents of the Crate	2
Unpacking Labotom-3	2
Placing Labotom-3	2
Getting Acquainted with Labotom-3	3
Cooling Valve	4
Back view	5
Noise Level	6
Supplying Power	6
Direction of the Cut-off Wheel	6
Connecting a Cooli-1/-3/-5 Cooling Unit.....	7
Mounting the Quick-Clamping Device and Spring Clamp	8
2. Basic Operations	
Using the Controls.....	9
Front Panel Controls of Labotom-3.....	9
Controls	10
Mounting a Cut-off Wheel	11
Dismounting a Cut-off Wheel	11
Clamping the Workpiece	12
Starting/Stopping the Cutting Process	13
3. Regular Maintenance	
Recirculation Cooling Unit.....	14
Daily	14
Weekly	14

1. Getting Started

Checking the Contents of the Crate

In the packing crate you should find the following parts:

- 1 Fork spanner 24mm, for changing the cut-off wheel
- 1 Inlet hose, 13mm dia., 2m
- 1 Outlet hose, 50mm dia., 2m
- 4 Hose clamps
- 1 Elbow pipe, 90°
- 1 Elbow pipe, 45°
- 1 Drain pipe
- 1 Set of Instruction Manuals

Unpacking Labotom-3

- Use a 13 mm socket spanner to remove the screws underneath the pallet.


Placing Labotom-3

Labotom-3 should be placed on an optional Struers table unit, or other suitable table. Ensure that the table is level. The table should be able to carry a weight of min. 100 kg. The depth should be at least 750 mm. Generally, the best working position will be achieved using a table 800 mm high or lower. Make sure there is enough room behind the table for the inlet and outlet hoses and for the protection guard to be opened fully.


Note: Placing Labotom-3 on a table unit
Lead the electric cable, inlet and outlet hose through the opening in the back of the table.

**Getting Acquainted
with Labotom-3**

Take a moment to familiarise yourself with the location and names of the Labotom-3 components.

- ① Front panel controls
- ② Lifting handle for protection guard
- ③ Protection guard
- ④ Cutting handle
- ⑤ Opening for protruding workpieces
- ⑥ Main power switch


- ① Spindle lock
- ② Tilting cutting unit
- ③ Connection for external fume extraction system (optional kit required)
- ④ Cut-off wheel
- ⑤ Cutting table
- ⑥ Quick-clamping device (optional)
- ⑦ Spring clamp (optional)

Cooling Valve

The cooling valve lever is mounted behind the guard for the cut-off wheel.


Cutting the lever is turned clockwise and the water is directed to cool the cut-off wheel.

Cleaning the lever is turned counter-clockwise to flush the cutting chamber.

*Labotom-3
Instruction Manual*

Back view

- ① Electrical cable for power supply
- ② Back plate
- ③ Name plate
- ④ Flush hose
- ⑤ Water outlet
- ⑥ Water inlet
- ⑦ Cooli Unit connection

Noise Level

Approx. 64 dB (A) measured at idle running, at a distance of 1.0 m/39.4" from the machine.

Supplying Power

<p><i>IMPORTANT</i> Check that the mains voltage corresponds to the voltage stated on the type plate on the back of the machine.</p>

- Labotom-3 is factory mounted with an electric cable. Mount a plug on the cable:
 - Black and brown: phase
 - Yellow/green: earth

Direction of the Cut-off Wheel



Check that the cut-off wheel rotates in the direction indicated by the arrow on the guard for the cut-off wheel. If the direction of rotation is incorrect, switch two of the phases.

Connecting a Cooli-1/-3/-5 Cooling Unit

To ensure optimal cooling, Labotom-3 can be fitted with a Cooli unit. Cooling System 3 is a Cooli configuration designed for use with Labotom.

Note

Before connecting the cooling unit to the Labotom, follow the instructions in the Cooli-1/-3/-5 Instruction Manual to prepare it for use.


To connect the Labotom-3 to a Cooli unit:

- Plug the Cooli-1/-3/-5 control unit's communication cable into the Labotom's control socket.
 - Connect the Cooli water inlet hose to the Cooli pump using the quick coupling (A).
 - Connect the other end of the hose to the quick coupling of the Labotom water inlet.
-
- Insert the drain pipe ① in the water outlet on the back of Labotom-3 and mount the 90° elbow pipe ②. Lubricate the sealing ring with grease or soap to facilitate insertion. (Use the other elbow pipe if more suitable).
 - Strip the steel spring from approx. 3 cm of the outlet hose and cut. Bend the cut end towards the centre of the hose. Mount the outlet hose ③ onto the elbow pipe and clamp the stripped section using a hose clamp.
 - Check that the outlet hose slopes downwards when connected. If necessary adjust the length of the hose.
 - Insert the open end of the hose into mounting hole in the bracket on top of the Cooli filter unit (B).
 - Connect the cooling unit to the mains power supply.

IMPORTANT

Before connecting, check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.

Mounting the Quick-Clamping Device and Spring Clamp


- ① Quick clamping device
- ② Spring clamp

The quick clamping device and spring clamp are ordered separately.

- Mount the back stop for the quick clamping device on the left hand side of the cutting table. The cut out corner should be on the right. Do not tighten the screws.
- Mount the back stop for the spring clamp on the right hand side of the cutting table, with the flat surface facing the front. Do not tighten the screws.
- Place the back stops as needed. The usual position is at the back of the cutting table parallel to the front. (Use a ruler). Tighten the screws using the socket spanner.
- Mount the quick clamping device on the left hand side of the cutting table and the spring clamp on the right. Their positions can be adjusted for the dimensions of the workpiece. Tighten the screws using the socket spanner.

2. Basic Operations

Using the Controls Front Panel Controls of Labotom-3

Location of the Main
Switch




Please see drawing 1 in the Section *Getting Acquainted with Labotom-3*.

Location of the Flush Hose


Please see drawing 2 in the Section *Getting Acquainted with Labotom-3*.

*Labotom-3
Instruction Manual*

Controls

Name	Key	Function	Name	Key	Function
<p>① MAIN SWITCH</p>		<p>The main switch is located on the right hand side of Labotom-3. Turn clockwise to switch on the power. The main switch also serves as an emergency stop.</p>	<p>④ FLUSH</p>		<p>Starts the pump for flushing. The flush nozzle is located in the rear of the cutting chamber, to the right. The flushing hose can be pulled out and pushed back. Press nozzle to flush.</p>
<p>② START</p>		<p>Starts the machine. The cut-off wheel starts rotating and the cooling water is turned on. Cannot be activated with the protection guard open or overload of the cutting motor.</p>			
<p>③ STOP</p>		<p>Stops the machine. The cut-off wheel stops rotating and the cooling water is turned off. Stops the flushing, when activated.</p>			

Mounting a Cut-off Wheel

- ① Nut
- ② Washer
- ③ Flange

- Put the rim into the groove in the top of the guard for the cut-off wheel and slot the cut-off wheel into position.
- Mount the flange, washer and nut.
- Press the knob for the spindle lock on the right-hand side of the cut-off wheel, turning the cut-off wheel until the spindle lock clicks.
- Tighten the nut moderately with the fork spanner and release the spindle lock.

Note: The spindle on Labotom-3 is left-hand threaded.

IMPORTANT

Conventional cut-off wheels like $\text{Al}_2\text{O}_3/\text{SiC}$ should be placed between two cardboard discs, to protect the cut-off wheel.
For maximum precision with diamond or CBN cut-off wheels, do not use cardboard discs.

Dismounting a Cut-off Wheel

- Push the cutting handle backwards to place the tilting cutting unit in the back position.
- Press the knob for the spindle lock on the right-hand side of the cut-off wheel, turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with the fork spanner. Remove the washer, flange and old cut-off wheel.

Note: The spindle on Labotom-3 is left-hand threaded.

Clamping the Workpiece

The following covers the use of the quick-clamping device and spring clamp, which can be ordered separately.

- Clean the cutting table with the flush hose.
- Place the workpiece between the clamp in the quick-clamping device and the back stop on the left hand side of the cutting table.
- Adjust the position of the clamps, if needed. Use the socket spanner.
- Pull the handle for the spring clamp and place the workpiece in the desired position for the cut. Lower the cut-off wheel to check the position.
- Turn the handle on the quick-clamping device to vertical position. Push the clamp device towards the workpiece and lock it firmly by pulling the locking handle back.

IMPORTANT

It is very important that the workpiece is firmly and securely fixed in the quick-clamping device.

Note: Placing the Back Stop

Generally the back stop should be placed at the rear of the cutting table in order to leave maximum room for the workpiece.
When cutting ductile materials like aluminium, a smoother cutting action can be achieved by placing the workpiece further to the front of the machine.

Starting/Stopping the Cutting Process

- The protection guard is fitted with a hydraulic brake action. Carefully close the protection guard, release and allow to close automatically.
- Switch on the main power.
- Press the START  button to start the machine. The cut-off wheel starts rotating and the cooling water is turned on (check that the cooling valve lever is in cutting position -lever horizontal).
- Carefully move the cut-off wheel towards the workpiece by pulling the cutting handle until it is in contact with the workpiece.
- Let the cut-off wheel make a small notch in the workpiece. Then increase the force and continue cutting. Adapt the speed at which the cut-off wheel is fed through the workpiece to suit the material and the wheel.
- When the cut-off wheel has almost cut through the workpiece, reduce the cutting force.
- Return the cutting handle to its back position.
- Press the STOP  button to stop the cut-off wheel and cooling water. For safety reasons, wait until the cut-off wheel has stopped completely before opening the protection guard.

IMPORTANT

The protection guard on Labotom-3 is equipped with a safety switch. The cutting motor will not start with the protection guard open.

REMEMBER...

Leave the protection guard open when the machine is not used to let the cutting chamber dry completely.

3. Regular Maintenance

Recirculation Cooling Unit

For Maintenance of the Recirculation Cooling Unit please refer to the Cooli Instruction Manual.

Daily

Cleaning the Cutting Chamber

Keep the working area around the cutting table clean at all times. To ensure a longer lifetime for your Labotom-3, Struers strongly recommends that you clean the whole cutting chamber with the cleaning hose whenever the machine is left unused for a period.

- Turn the cooling valve lever counter-clockwise (lever vertical).
- Press the FLUSH  button to start the pump for flushing.
- Pull the flush hose out from the back of the cutting chamber.
- Clean the cutting chamber thoroughly. Adjust the flow as needed. Avoid flushing around the connection to the fume extraction system (if mounted).
- Push the STOP  button to stop the pump.
- Tuck the cleaning hose back into place.
- Leave the protection guard open to let the cutting chamber dry completely.

REMEMBER...

Leave the protection guard open when the machine is not in use to let the cutting chamber dry completely.

Protection Guard

- Clean the protection guard window with a damp cloth.

Weekly

- Clean the cutting chamber thoroughly.

Reference Guide

Table of Contents	Page
1. Advanced Operations	
Clamping Irregular Workpieces	16
Long Workpieces	16
Mounting the Rubber Curtain	16
Cutting Long Workpieces	16
Adjusting the Cutting Handle.....	17
Connection to an External Exhaust System	17
Optimising the Cutting Results	18
2. Accessories	19
3. Consumables	19
4. Trouble-Shooting	20
5. Maintenance	
Adjustments	23
Maintenance of Cutting Tables	23
Maintenance of Cut-off Wheels	23
Storing of Conventional Cut-off Wheels.....	23
Maintenance of Diamond and CBN Cut-off Wheels.....	23
Maintenance of Clamping Devices.....	23
6. Technical Data	24

1. Advanced Operations

Clamping Irregular Workpieces Irregular workpieces without plane clamping surfaces must be clamped using special clamping tools, as the workpieces must not move during the cutting. This could result in damage to the cut-off wheel or to the workpiece itself. Use the T-slots to mount the special clamping tools.

Struers offers a selection of Clamping Tools (Please refer to the [Struers Clamping Tools brochure](#) for details of the range available).

To achieve faster cutting, position the workpiece so that the wheel will cut the smallest possible cross-section.

Long Workpieces

To cut workpieces exceeding the width of the cutting chamber, a plate in the cutting chamber wall may be removed and replaced by rubber curtains.

A set of rubber curtains can be ordered separately as part of a spare part kit (Cat. No. 15322902).

Please note

In some countries it is a requirement to have a fixed side wall in the cutting chamber and rubber curtains may not be used.
Please check local regulations before ordering.


Mounting the Rubber Curtain

- Remove the left-hand side cover of the protection hood.
- Mount the rubber curtains and the attachment rail inside the hood and secure with the 3 M5 screws supplied.


Cutting Long Workpieces

- Place the workpiece without clamping it.
- Position the workpiece so that it can pass through the rubber curtain when the protection guard is closed.
- Check that the tilting cutting unit can pass, when cutting.
- Clamp the workpiece.
- If needed, put a cloth around the protruding workpiece to absorb any cooling water that may leak out.

Adjusting the Cutting Handle

Connection to an External Exhaust System

The cutting handle is factory-mounted in a position which will suit most people when Labotom-3 is placed on a table 800 mm high. However it is possible to change the position of the cutting handle:

- Loosen the 2 screws. Take care not to displace the cutting motor axially.
- Turn the handle to the desired position.
- Tighten the screws.
- Check that the cut-off wheel will pass through the middle of the slot in the cutting table.

Struers recommends the use of an exhaust system as workpieces may emit harmful gasses when cut.

A passage for exhaustion is located on the left hand side of Labotom-3. An optional fume extractor kit enables the connection of a 50mm dia. exhaust hose.

- Remove the cover and seal.
- Mount the connecting branch ① outside the protection guard using the seal ② and the long screws supplied with the fume extractor kit.
- Assemble the spacers ③, splash guard ④ and cover ⑤ inside the protection guard. Secure with the short screws originally used for mounting the cover. Make sure that the splash guard lies directly under the cover.
- Strip the steel spring from approx. 2 cm of the exhaust hose and cut. Bend the cut end towards the centre of the hose. Mount the exhaust hose onto the connecting branch ① and clamp the stripped section using a hose clamp.
- Connect the exhaust hose to your local exhaust system.

Optimising the Cutting Results The following table shows possible answers to a number of common questions:

Optimising the Cutting Results	
Question	Answer
How can I avoid discoloration or burning of the sample?	Apply a lower cutting force
	Change the cut-off wheel as the hardness of the present cut-off wheel may be inappropriate for the hardness of the workpiece*)
How can I avoid burrs?	Use a softer cut-off wheel*)
	Where a quick-clamping device (optional) and spring clamp (optional) are mounted: Clamp the workpiece securely using the quick-clamping device. Position the spring clamp to provide a force just sufficient to prevent the cut-off piece from shifting at the end of the cutting action.
How can I avoid the cut-off wheels wearing too quickly?	Apply a lower cutting force or use a harder cut-off wheel*)
How can I achieve faster cutting?	Position the workpiece so as to cut the smallest possible cross-section. Apply a higher cutting force.

*) Please refer to the Selection Guide in the [Struers Cut-off Wheels brochure](#).

2. Accessories

Clamping Tools

Please refer to the [Struers Clamping Tools brochure](#) for details of the range available.

Specification	Cat. No.
Fume extractor kit. Connecting tube 50 mm dia.	05326902
Table unit. With compartment for recirculation cooling unit. 800 mm high	05116916
T-slot cleaner for 10 and 12 mm T-slots	05486910
Cooli System 3 For water cooling of Labotom-3 With 50 l tank, small pump, Cooli-1 and filter bag	
1 x 100 V / 50 Hz	05766516
1 x 220-240 V / 50 Hz	05766522
1 x 100-120 V / 60 Hz	05766523
1 x 220-240 V / 60 Hz	05766524
1 x 100-120 V / 50-60 Hz CSA	05766616

3. Consumables

Cut-off Wheels

Please refer to the Selection Guide in the [Struers Cut-off Wheels brochure](#).

Other Consumables

Specification	Cat. No.
<i>Corrozip-LF</i> Additive for Cooling Fluid	
1 l	49900038
5 l	49900039

4. Trouble-Shooting

Error	Explanation	Action
Machine Problems		
The cut-off wheel does not rotate.	The protection guard is open.	Press the protection guard completely down. If this does not help: Call a Struers service technician.
	Cutting motor overloaded due to heavy use.	Open the protection guard and let the cutting motor cool down for about 5 - 10 min.
No cooling water.	Cooling valve lever is in cleaning position.	Turn the cooling valve lever clockwise (lever horizontal).
No water in the cleaning hose.	The nozzle is clogged.	Unscrew the jet and clean. If this does not help: Unscrew the housing on the nozzle side and clean inside. If possible use compressed air.
	Cooling valve lever is in cutting position.	Turn the cooling valve lever counter-clockwise (lever vertical).
Water leaking.	Leak in a hose of Recirculation Cooling Unit.	Check the hose and tighten the hose clamp.
	Water overflow in the cooling water tank.	Remove the excess water from the tank.
Workpieces or cutting chamber rusty.	Insufficient additive for cooling fluid.	Add Struers Additive for cooling fluid to the cooling water, using the correct concentration. Please follow the instructions in the Recirculation Cooling Unit Instruction Manual.
	The machine is left with closed protection guard.	Leave the protection guard open to let the cutting chamber dry.
Protection guard window blurred.	Insufficient cleaning.	Clean with mild soapy water.
Quick-clamping device unable to hold the workpiece.	The quick-clamping device is not balanced.	Adjust the screw underneath the clamping column. Use a 3mm Allen key.
	Clamping heart worn.	Call Struers service technician.

*Labotom-3
Instruction Manual*


Error	Explanation	Action
Cutting Problems		
Discoloration or burning of the workpiece.	The hardness of the cut-off wheel is inappropriate for the hardness / dimensions of the workpiece.	See Consumables Section, Cut-off Wheels.
	The force on the cut-off-wheel is too high.	Apply a lower force on the cut-off wheel.
	Inadequate cooling.	-Check that there is enough water in the recirculation cooling unit. -Check the cooling water hoses.
Unwanted burrs.	Wheel too hard.	See Section Consumables, Cut-off Wheels.
	Too high force on the cut-off wheel near the end of the operation.	Reduce the cutting force near the end of the operation.
	Lack of support.	If possible, support the workpiece at both sides.
The cutting quality differs.	Cooling water hose clogged.	Check the cooling water flow by pressing the button for cleaning. Clean the cooling water hose and the cooling tube. If possible use compressed air.
	Insufficient cooling water.	Refill tank with water. Remember Struers Additive.
The cut bends to a side.	Initial cutting rate too fast.	Let the cut-off wheel make a small notch in the workpiece before the actual cutting is carried out.
	The force on the cut-off-wheel is too high.	Apply a lower force on the cut-off wheel.
The cut-off wheel breaks.	Incorrect mounting of the cut-off wheel.	- Check that the hole has the correct diameter. - Check there is a cardboard washer on both sides of the cut-off wheel (Conventional cut-off wheels only). - The nut must be tightened moderately
	Incorrect clamping of the workpiece.	Make sure that only one side of the workpiece is clamped firmly. The other side should only be fixed lightly. Use support tools (optional) if the geometry of the workpiece makes support necessary.
	Not enough support of the workpiece.	Support the free end of the workpiece.
The cut-off wheel breaks.	Cut-off wheel too hard.	See Consumables Section, Cut-off Wheels.

*Labotom-3
Instruction Manual*

Error	Explanation	Action
	The force on the cut-off-wheel is too high.	Apply a lower force on the cut-off wheel.
	Inadequate cooling.	- Check that there is enough water in the recirculation cooling unit. - Check the cooling water hoses.
The cut-off wheel wears down too quickly.	The force on the cut-off-wheel is too high.	Apply a lower force on the cut-off wheel.
	The cut-off wheel is too soft for the task.	See Consumables Section, Cut-off Wheels.
	Labotom-3 vibrates (worn bearings).	Call Struers service technician.
The cut-off wheel does not cut through the workpiece.	Incorrect choice of cut-off wheel.	See Consumables Section, Cut-off Wheels.
	Cut-off wheel worn.	Replace the cut-off wheel.
	The cut-off wheel gets caught in the workpiece because of internal stress in the workpiece.	Support the workpiece and clamp it on both sides of the cut-off wheel in such a way to allow the cut to stay open.
		Make a relief cut: Cut about halfway through the workpiece. Turn the workpiece 180° and position the cut approx. 1 - 2 mm off-centre.
The workpiece breaks when clamped.	The workpiece is brittle.	Place the workpiece between two polystyrene plates. NB! Always cut brittle workpieces very carefully.
The sample is corroded.	The sample has been left in the cutting chamber for too long.	Leave the protection guard open, when you leave the machine.
	Insufficient additive for cooling fluid.	Add Struers Additive for cooling fluid to the cooling water in the correct concentration. See Maintenance Section.

5. Maintenance

Adjustments

The tilting cutting unit should offer a slight resistance when being moved. Too much friction will be tiring during operation whilst too little might cause excessive vibrations during cutting. The friction is adjusted from the factory, but it might change during use. To adjust the friction:

- Use the socket spanner, also used for positioning the quick-clamping device and back stop.
- Adjust the friction in the bearing for the tilting cutting unit.

Maintenance of Cutting Tables

The stainless steel bands (available as spare parts) should be replaced if damaged or worn.

Maintenance of Cut-off Wheels

Storing of Conventional Cut-off Wheels

These cut-off wheels are sensitive to humidity. Therefore, do not mix new, dry cut-off wheels with used humid ones. Store the cut-off wheels in a dry place, horizontally on a plane support.

Maintenance of Diamond and CBN Cut-off Wheels

The precision of diamond and CBN cut-off wheels (and thus the cut) depends on how carefully the following instructions are observed:

- Never expose the cut-off wheel to a heavy mechanical load, or heat.
- Store the cut-off wheel in a dry place, horizontally on a plane support, preferably under light pressure.
- A clean and dry cut-off wheel does not corrode. Therefore, clean and dry the cut-off wheel before storing. If possible, use ordinary detergents for the cleaning.
- Regular dressing of the cut-off wheel is also part of the general maintenance.

Maintenance of Clamping Devices

Important

It is recommended to thoroughly clean and lubricate the Quick Clamping Device and Vertical Clamping Device at regular intervals.


6. Technical Data

Subject		Specifications	
		Metric/International	US
Cutting	Cutting power constant [S1] Cutting power intermittent [S3] Maximum power constant Rotational speed (running idle) Wheel speed (running idle), cut-off wheel 250 mm dia.	2.2 - 3.0 kW 3.2 - 4.0 kW 4.4 - 6.0 kW 2845 rpm (50 Hz) 3450 rpm (60 Hz) 37.2 m/s (50 Hz) 45.1 m/s (60 Hz)	3.0 - 4.0 HP 4.3 - 5.4 HP 6.0 - 8.0 HP 3450 rpm (60 Hz) 101 mph (60 Hz)
Cut-off Wheel	Diameter x thickness x hole	250 x 1 - 1.5 x 32 mm	9.8 x 0.04 - 0.06 x 1.26"
Flange for Cut-off Wheel	Diameter	62.5 mm	2.5"
Cutting Chamber Geometry	<i>Max. workpiece dimensions:</i> To the left of the cut-off wheel Width Depth Thickness To the right of the cut-off wheel Width Depth <i>Passage for protruding workpieces:</i> At left hand side Width Height At right hand side Width/height	145 mm Up to 500 mm Up to 250 mm 75 mm Up to 210 mm 103 mm 95 mm 25 mm	5.7" Up to 19.7" Up to 9.8" 3.0" Up to 8.3" 4.1" 3.7" 1.0"
Cutting Table	Width Depth T-slots	200 mm 215 mm 10 mm, 3 pcs	7.9" 8.5"
Cutting Capacity	Max. cutting diameter For other dimensions: Please see graph. (The actual cutting capacity depends on the sample material, cut-off wheel and clamping technique).	90 mm	3.5"
Fume Extraction	<i>Optional fume extractor kit:</i> Dia. for connecting tube	50 mm	2.0"
Dimensions and Weight	<i>Labotom-3</i> Width Depth Height Weight <i>Optional Struers Table Unit</i> Width: Depth: Height:	651 mm 645 mm 410 mm 80 kg 900 mm 750 mm 800 mm	25.6" 25.4" 16.1" 180 lbs 35.5" 29.5" 31.5"

*Labotom-3
Instruction Manual*

Subject	Specifications			
Environment	Directives	Please refer to the Declaration of Conformity		
	Noise level	Approx. 64 dB(A) running idle, at a distance of 1.0 m / 39.4" from the machine.		
Electrical Data	Voltage/frequency:	Max. Load:		
	3 x 200 V / 50 Hz 3 x 200-210 V / 60 Hz 3 x 220-230 V / 50 Hz 3 x 220-240 V / 60 Hz 3 x 380-420 V / 50 Hz 3 x 460-480 V / 60 Hz	14.5 A 11.5 A 12.5 A 11.5 A 7.0 A 6.0 A		
Mains Cable Specification	Voltage/ frequency	Cable size	Min. Fuse size	Max. Fuse size
	3x200V/50HZ	3x4mm ² + PE	20	25
	3x200-210/60CSA	3xAWG12 + PE	20	25
	3x220-230V/50	3x4mm ² + PE	20	25
	3x220-240/60CSA	3xAWG12 + PE	20	25
	3x380-415V/50	3x4mm ² + PE	10	20
	3x380-415/60CSA	3xAWG12 + PE	10	20
	3x460-480/60CSA	3xAWG12 + PE	10	20
<i>Important:</i> Local standards may overrule the recommendations for the main supply cable. If necessary, please contact a qualified electrician to verify which option is suitable for the local installation setup.				
Residual Current Circuit Breaker	type A, 30 mA (or better) is recommended.			

Cutting Capacity

The graph shows the projected cutting capacity under the following conditions: A new cut-off wheel. The workpiece is laid directly on the cutting table, with overhang where appropriate. Vertical clamping is used. The actual cutting capacity depends on the sample material, cut-off wheel and clamping technique.

Quick Reference

Clamping the Workpiece

- Place the workpiece between the clamp and the back stop on the left hand side of the cutting table.
- Push the clamp towards the workpiece and lock the quick-clamping device with the locking handle.

Starting/Stopping the Cutting

- Carefully close the protection guard.
- Press START . The cut-off wheel starts rotating and the cooling water starts running.
- Lower the cut-off wheel onto the workpiece by pulling the cutting handle slowly until contact with the workpiece. Do not apply too heavy a cutting force.
- Let the cut-off wheel make a small notch in the workpiece. Increase the pressure and keep it constant while cutting.
- When the cut-off wheel is almost through the workpiece the cutting pressure should be reduced.
- After cutting, return the cut-off wheel to top position. Stop the cut-off wheel and the cooling water by pressing STOP . For safety reasons, wait about 5 seconds before opening the protection guard.
- Open the protection guard and release the clamping device.
- Remove the workpiece.

Dismounting a Cut-off Wheel

- Push the cutting handle backwards to place the tilting cutting unit in the back position.
- Press the knob for the spindle lock on the right-hand side of the cut-off wheel, turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with the fork spanner. Remove the washer, flange and old cut-off wheel.

Mounting a Cut-off Wheel

- Put the rim into the groove in the top of the guard for the cut-off wheel and then catch the hole of the cut-off wheel with the spindle.
- Mount the flange, washer and nut.
- Press the knob for the spindle lock on the right-hand side of the cut-off wheel, turning the cut-off wheel until the spindle lock clicks.
- Tighten the nut moderately with the fork spanner and release the spindle lock.

Cleaning the cutting chamber

- Press the Flush  button to start the pump for flushing.
- Pull the flush hose out from the back of the cutting chamber.
- Clean the cutting chamber thoroughly. Adjust the flow as needed.
- Push the STOP  button to stop the pump.
- Tuck the cleaning hose back into place.
- Leave the protection guard open to let the cutting chamber dry completely.

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Date of Release 04.04.2011

Labotom-3
Spare Parts and Diagrams

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2010.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Telefax: +45 44 600 801

Spare Parts and Diagrams

Table of contents

Drawing

Labotom-3

Drawings

Labotom-3, complete	15320001D
Bottom, complete	15320010L
Protection guard, complete	15320040E
Cutting screen	15320050B
Quick Clamping Tool for 10 mm T-slots, Left	15870040H
Quick Clamping Tool for 10 mm T-slots, Right.....	15870041H
Stopper for 10mm T-slots.....	15870045C

Diagrams

Circuit diagram	15323115D
Wiring diagram	15323450C
Transformer connections (200-240V)	14923150A
Transformer connections (380-480V)	14923160B

Some of the drawings may contain position numbers
not used in connection with this manual.

Labotom-3
Spare Parts and Diagrams

Spare part list for Labotom-3


Drawing	Pos.	Spare Part	Cat no:
15320001		Labotom-3, complete	
	110	Linear damper	2YS20090
	120	Flushing hose complete with tap	15320007
	130	Bushing	2GK90129
	140	Label with warning of hand injury	10160401
15320010		Bottom, complete	
	180	Buffer	2GS10315
	190	Vibration damper (Stop in front)	2GS00200
	210	Steel band 20 mm	15110271
	210 + 220	Replacement steel bands 1x20+1x40 mm	05116906
	220	Replacement steel bands 4x40 mm	05116907
	430	Spindle extender	15320220
	450	Guide pin ø6 (10 pcs)	315MP037
	460	Fixed flange	15320226
	470	Cardboard blotters ø33/85F.xxTRE (25pcs)	368MP077
	480	Loose flange	15320225
	490	Washer for spindle	14150012
	500	Nut for spindle	2TA70160
	630	Ball knob	2GH00135
650-660	Rubber feet with screw (4 pcs)	10162901	
680	O-ring 50.00-4.00 NBR70 (10 pcs)	RIO40080	
15320040		Protection guard complete	
	90 -110	Rubber curtain kit, left	15322902
	110 -130	Rubber curtain kit, right	15322903
	150 - 160	Sealing band, complete	15322904
15320050		Cutting screen	
	20 - 80	Locking pin, complete	15322901

Drawing

Spare part list for Labotom-3

Pos.	Spare Part	Cat no:
	Various	
	Angle pipe wrench, 13 mm	RIP30458
	Fork spanner, 24 mm	RGR00224
	Pipe piece for drain	15110997
	Elbow 45 ⁰	381MP360
	Elbow 87 ⁰	2NG20587
	Hose clamp 40-60 (5 pcs)	RNS24060
	Hose clamp S19/9Zy (5pcs)	RNS11909
	Drain hose ø51 (6 metre)	RNU30251
	Water hose ½" (10 metre)	RNU29316
	<i>Replacement stainless steel bands</i>	
	1 pc. 40 mm and 1 pc. 20 mm	05116906
	4 pcs. 40 mm	05116907
	Vertical clamping system for 10 mm T-slots	
	TREVS	
1	Operating key for TREVS	2GR27021
2	Spring for TREVS	2GR27023
3	Flat clamping shoe	2GR26150


Vertical clamping system for 10 mm T-slots (TREVS)

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













Omnifit 230M between house and insert, both parts must be degreased first!

Special tool are used for assembling "15871505".

Omnifit 230M

Special tool are used for assembling "15871506".

Omnifit 230M

Left quick clamp.tool:
Shorten the grease nipple 2 mm until base plate is change.

Filled with Grease SHELL Alvania EP2 (50g)

DETAIL C

(2x) Omnifit 50M

(4x) Omnifit 50M

End of pin must be flush or below top of base plate.

- 30 Use mounting tool to align rod and jaw.
- 50 30 Use mounting tool to align rod and jaw.
- 30 The house is to be given "Stainless steel Polish", after its assembled.
- 230 Adjust for easy rod - movement.

Revision E: Pos 190 and 210 removed ; Pos 250 and 270 change ; Pos 290 and 300 added ; Text for Omnifit 230M between Pos 40 and 70 added ; Omnifit 230M for pos 200 removed. ; Text for pos 270 added (Special tool). See the difference between drawing revision D and E for more detail.




H	2010-01-15	Pos.200. Material changed: Stainl.St.->Hard.St.	BRY	2010-01-15	JTV
A	18-04-06		MD	18-04-06	FPG
Revision	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
		Material:	Scale: 1:1	Format: A3	Tolerance: DS/ISO 2768- mK Surface treat.: None
		ID:	Description:		Rev:
		15870040 Quick Clamping Tool 10mm Left, Assy			H

Struers
Pederstrupvej 84
DK-2750 Ballerup/Copenhagen
Denmark
Phone: +45 44 600 800
Fax: +45 44 600 804

Omnifit 230M between house and insert, both parts must be degreased first!

Special tool are used for assembling "15871505".

Special tool are used for assembling "158715062".

- 250 Tighten with torque wrench 7 Nm.
- 50 30 Use mounting tool to align rod and jaw "15871501".
- 40 The house is to be given "Stainless steel Polish", after its assembled.
- 230 Adjust for easy rod - movement.

Revision E: Pos 190 and 210 removed ; Pos 250 and 270 change ; Pos 290 and 300 added ; Text for Omnifit 230M between Pos 40 and 70 added ; Omnifit 230M for pos 200 removed. ; Text for pos 250 added (Special tool). See the difference between drawing revision D and E for more detail.

H	2010-01-15	Pos.300. Material changed: Stainl.St.->Hard.St.	BRY	2010-01-05	JTV
A	18-04-06		MD	18-04-06	FPG
Revision	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
		Material:	Scale: 1:1	Format: A3	Tolerance: DS/ISO 2768- mK
			ID:		Description:
15870041 Quick Clamping Tool 10mm Right, Assy					Rev: H

Struers
 Pederstrupvej 84
 DK-2750 Ballerup/Copenhagen
 Denmark
 Phone: +45 44 600 800
 Fax: +45 44 600 804

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


1

MAINS SWITCH

ALL WIRINGS 1st (AWG16), EXCEPT OTHERWISE MARKED.

COLOR CODES:
 BK = BLACK
 BN = BROWN
 RD = RED
 OG = ORANGE
 YE = YELLOW
 GN = GREEN
 BU = BLUE
 VT = VIOLET
 GY = GREY
 WH = WHITE

**FT1: Thermal Overload Relay:
 With Short Circuit Trip.
 SETTING RANGES**

MAINS SUPPLY:	FT1: SET =
3x200 - 230V / 50Hz.	0,4 Amp.
3x200 - 240V / 60Hz.	0,4 Amp.
3x380V / 60Hz.	0,23 Amp.
3x380 - 415V / 50Hz.	0,23 Amp.
3x460 - 480V / 60Hz.	0,20 Amp.

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















Herewith declares that

<i>Product Name:</i>	Labotom-3
<i>Type No:</i>	532
<i>Machine Type:</i>	Cut-off machine

is in conformity with the provisions of the following directives:

Safety of Machinery 2006/42/EC according to the following standard(s):
EN ISO 12100-1:2005, EN ISO 12100-2:2005, EN 60204-1:2006, EN ISO 13857:2008,
EN 954-1:1999, EN ISO 14121-1:2007, EN 614-1:2007, EN 1088:1998.

EMC-Directive 2004/108/EC according to the following standard(s):
EN 61000-6-1:2002, EN 61000-6-3:2002.

Low Voltage Directive 2006/95/EC according to the following standard(s):
EN 60204-1:2006.

Supplementary Information The equipment complies with the American standards:
FCC part 15, subpart B, UL508, ED14 22-2-88 and NFPA 79:1997.

The above has been declared according to the global method, module A

Date: 21.06.2010

Christian Skjold Heyde,
Vice President, R&D and Production, Struers A/S

Hersteller
Struers A/S
Pederstrupvej 84
DK-2750 Ballerup, Danmark
Telefon +45 44 600 800

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<i>Produktname:</i>	Labotom-3
<i>Typennr.:</i>	532
<i>Maschinenart:</i>	Trennmaschine

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Christian Skjold Heyde,
Vice- President, R& D et Production, Struers A/S

Labotom-3

**Powerful manual
cut-off machine with high
cutting capacity**

- Labotom-3 has a powerful 3.0 kW / 3.2 kW cutting motor and a cutting capacity of up to 90 mm (3.5") diameter.
- Cutting table with top surface of stainless steel and integrated T-slots.
- Bright LED light in the cutting chamber.
- Adjustable cutting handle to suit the height of the operator and the table unit.
- Flush hose for cleaning the cutting chamber.

Labotom-3 is a desirable option for materialographic laboratories preparing a moderate number of specimens or for sampling stations situated in production units. It is simple to operate and practically no training is needed.

Easy and comfortable Operation

Labotom-3 is a compact table-top model for abrasive wet cutting.

When opening the large protection guard, there is free access for clamping the workpiece and exchange of the cut-off wheel. The closing of the protection guard activates a hydraulic brake action.

A powerful LED lamp lights up the cutting chamber and facilitates clamping and monitoring the cutting operation through the large window in the protection guard. Cooling water is applied directly onto the 250 mm (10") dia. cut-off wheel resulting in an excellent cooling of the cutting

area. The spring balanced cutting unit gives a relaxed advance at the start of cutting and a relaxed return when cutting is finished. After cutting, the cutting chamber can be cleaned with the built-in flush hose.

The smooth surface of the base and a large drain (50 mm dia.) make it easy to flush away the swarf produced during cutting.

The position of the cutting handle can be adjusted to suit the height of the operator and laboratory table. The simple design of Labotom-3 makes it a very reliable machine with low maintenance.

Cutting Capacity

The left side wall is prepared for the installation* of an optional rubber curtain. The projected cutting capacity for a new cut-off wheel is up to 90 mm dia. As workpieces can be placed in various positions on the cutting table, the best utilisation of the cut-off wheel can be achieved, resulting in optimum cutting capacity. Please see graph.

Cutting Capability

Struers' entire range of 250 mm dia. cut-off wheels, including diamond and CBN cut-off wheels, can be used, therefore practically all types of materials can be cut.

*) Only in countries where permitted by local legislation.

The ample motor power and efficient cooling system enable fast cutting and perfect, clean cuts.

Flexible clamping of the Workpiece

The integrated cutting table has three 10 mm T-slots, facilitating the use of universal clamping tools.

Safety Features

The cutting motor will not start with the protection guard open. Any attempt to open the guard during cutting results in a safety switch stopping the cutting motor.

A thermal fuse stops the cutting motor in the event of overload. The cut-off machine can be connected to a fume extraction system to vent any harmful fumes produced during cutting.

Technical Data

Cutting Motor:	Three-phased motor:	3.2 kW (4.3 HP) at 50 Hz 3.0 kW (4.1 HP) at 60 Hz
Cut-off Wheels:	Rotational speed (running idle):	2845 rpm (50 Hz) 3450 rpm (60 Hz)
Cut-off Wheel size:	250 mm x 1-1.5 mm x 32 mm (9.8" x 0.04" - 0.06" x 1.26")	
Cutting Chamber:	Maximum workpiece dimensions: To the left of the cut-off wheel: Width: 145 mm (5.7"), depth: up to 500 mm (19.7"), height up to 250 mm (9.8") To the right of the cut-off wheel: Width: 75 mm (3.0"), depth: up to 210 mm (8.3") Side opening for protruding workpieces: Left hand side (optional): 103 x 95 mm (4.1" x 3.7") Right hand side: 25 x 25 mm (1.0" x 1.0")	
Cutting Table:	Table size: 200 x 215 mm (7.9" x 8.5") with three 10 mm T-slots	
Cutting Capacity*	Max. cutting diameter: 90 mm (3.5"). Other dimensions: Please see graph. <i>*The cutting capacity stated is the projected cutting capacity with a new cut-off wheel. Actual cutting capacity depends on the sample material, cut-off wheel and clamping technique.</i>	
Environment:	Noise level: Approx. 64 dB(A) at a distance of 1.0 m (39.4") from the machine, running idle.	
Dimensions and weight:	Width: 651 mm (25.6") Depth: 645 mm (25.4")	Height: 410 mm (16.1") Weight: 80 kg (180 lbs.)

Specifications

Labotom-3

Manual cut-off machine. Complete with cutting table with 10 mm T-slots. Motor power 3.0 kW / 3.2 kW. Recirculation cooling unit (05766xxx) and clamping tools are ordered separately.

Cat. no:

053261xx

Accessories

Recirculation Cooling Unit

With 50 l tank, small pump, Cooli-1 control unit and 50 l static filter with filter paper

05766xxx

Quick-clamping Device, left

For securing the workpiece. Complete with backstop. Max clamping width: 50 mm.

05876909

Spring Clamp, right

For retaining the cut-off specimen. Complete with backstop.

05116904

Vertical Clamping System

For clamping irregularly shaped workpieces on the cutting table. Clamping height adjustable up to 58 mm.

05116905

Fume Extractor Kit

For connection of Labotom-3 to an external fume extraction system. (connecting tube 50 mm dia.)

05326902

Table Unit

With Compartment for recirculation cooling unit. Width: 900 mm, depth: 750 mm, height 800 mm.

05116916

Safety kit for Labotom-3

Factory mounted emergency stop and safety lock. Safety kit is ordered separately.

05326905


Struers' equipment is in conformity with the provisions of the applicable International Directives and their appurtenant Standards. (Please contact your local supplier for details).

Struers' products are subject to constant product development. Therefore, we reserve ourselves the right to introduce changes in our products without notice.

Quick Clamping device for Labotom-3

Vertical clamping system with riser block and extension arm

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup, Denmark
Phone +45 44 600 800
Fax +45 44 600 801
struers@struers.dk
www.struers.com

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CANADA

Struers Ltd.
7275 West Credit Avenue
Mississauga, Ontario L5N 5M9
Phone +1 905-814-8855
Fax +1 905-814-1440
info@struers.com

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Struers Inc.