Movipol-3

Instruction Manual

Manual No.: 13927003

Date of Release F1 01 00J

Movipol-3
Instruction Manual

Table of Contents	Page
User's Guide	1
Reference Guide	16
Quick Reference Guide	21

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the *Serial No.* and *Voltage* on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2009.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Movipol-3

Safety Precaution Sheet

To be read carefully before use

1. The operator must be fully instructed in the use of the apparatus, according to the User's Guide.
2. The operator must be fully instructed in the use of the electrolytes used with Movipol-3.
3. Observe the current safety precautions regarding the handling, mixing, filling and disposal of the electrolyte (see section Safety Precautions). The Material Safety Data Sheet for each of the Struers electrolytes can be supplied upon request.
4. Be sure that the actual voltage corresponds to the voltage stated on the type plate of the apparatus.
5. Check that the apparatus is standing firmly on the support. If the apparatus is carried or hung by the shoulder strap, check that the strap is intact and that the buckle is correctly fastened.
6. Always check that the pump motor is working before taking the pencil out of its holder.
7. Always place the pencil in its holder when it is not used.
8. Do not leave the apparatus unattended when filled with electrolyte.
9. Make sure that the apparatus is safely secured during transportation and that it does not contain electrolyte.
10. The maximum polishing voltage allowed at the work site must be observed.

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

Disposal

Equipment marked with a WEEE symbol  contain electrical and electronic components and must not be disposed of as general waste.

Please contact your local authorities for information on the correct method of disposal in accordance with national legislation.

User's Guide

Table of Contents	Page
1. Getting Started	
Checking the Contents of Packing	3
Placing Movipol-3.....	3
Getting Acquainted with Movipol-3	4
Supplying Power	5
Connecting the Transformer	5
Charging the Batteries	5
Connecting the Pump Motor and Cathode	5
Connecting the Anode.....	5
External Etching Operations.....	5
2. Operation	
Using the Controls.....	6
Front Panel Controls	6
Safety Functions	7
Cut-out of Pump Motor.....	7
Current Limiter Function of the Polish/Etch Circuit	7
Cut-out of Polish/Etch Circuit	7
Cut-out in case of under-voltage	7
Mounting the Container.....	7
Dismounting the Container.....	7
Mounting the Pump-Container Unit	8
Dismounting the Pump-Container Unit	8
Filling Electrolyte into the Container	8
Emptying the Container of Electrolyte	8
Mounting the Polishing Chamber	8
Starting the Polishing/Etching Process.....	9
Stopping the Polishing/Etching Process.....	9
Cleaning the Container and the Pump	11
Voltage Examples	11
Cast Iron. 3.4% C, 0.29% Si, 0.6% Mn, 0.3% P	11
Duplex Steel.....	11
Commercial Pure Aluminium 99.5% Al.....	11
Stainless Steel	12
Low Carbon Steel 0.1% C.....	12
High Carbon Steel.....	12
Titanium Alloys.....	12

3. Safety precautions

Electricity.....	13
Electrolytes in General	13
Perchloric Acid in Particular	13
Training of Operators	14
Mixing the Solution	14
Storage of Perchloric Acid or Solution.....	14
Fire and Explosion Hazards	14
Disposal	15
Working with Movipol-3 in the field.....	15
Dangerous Localities.....	15
Safety Measures	15
Dangerous situations	15

1. Getting Started

Checking the Contents of Packing

In the packing box you should find the following parts:

- 1 Movipol-3 with control unit, pump, pencil
- 1 Transformer
- 1 Container for electrolyte
- 1 Anode wire
- 1 Anode clamp
- 1 Magnet
- 1 Shoulder strap
- 10 Polishing Chambers
- 1 Funnel
- 1 Set of Instruction Manuals


Placing Movipol-3

Movipol-3 is a portable apparatus for field electropolishing. If possible the apparatus should be placed on a stable support. If not it can be hung by the shoulder strap in a suitable place.

During transportation the apparatus must not contain electrolyte and should be secured in order not to turn upside down.

**Getting Acquainted with
Movipol-3**

Take a moment to familiarise yourself with the location and names of the Movipol-3 components.

- | | |
|--|------------------------------------|
| ① Charge/pump switch | ⑧ Charge indicator |
| ② Start button | ⑨ Voltage indicator for pump motor |
| ③ Adjusting knob for polishing voltage | ⑩ Polishing indicator |
| ④ Adjusting knob for etching voltage | ⑪ Etching indicator |
| ⑤ Polishing time | ⑫ Socket pump/cathode |
| ⑥ Etching time | ⑬ Socket anode |
| ⑦ Ammeter | ⑭ Socket for external etching |
| | ⑮ Pencil |

Supplying Power

IMPORTANT

Check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.

Connecting the Transformer

- Connect the mains cable to the mains.
- Connect the 42V multi plug to the left side of Movipol-3.

Charging the Batteries

- Switch O/Charge - I/Pump to top position, the batteries of Movipol-3 are now being charged. When the batteries are fully charged, the Charge Indicator lights green. If the batteries need charging the Charge Indicator will light red.

Even though the batteries are factory charged, Struers recommends the charging of Movipol-3 batteries 1-2 hours before operating the apparatus.

WARNING

A fully charged set of batteries will self-discharge in 6 months, and a complete discharging may reduce the lifetime. Therefore the batteries should be charged at regular intervals.

IMPORTANT

When the transformer is being used Movipol-3 should only be operated when the Charge Indicator lights green or yellow.

Connecting the Pump Motor and Cathode

- Connect the multi plug to the socket in the control unit with the screw cap. Both the electrolyte pump motor and the cathode in the pencil are now connected to the control unit.

Connecting the Anode

- Place the banana plug of the anode wire in the socket, Anode. The connection to the sample is established with the enclosed magnet or clamp.

External Etching Operations

- When external etching operations should be carried out, i.e. when using a cathode connection other than the built-in cathode in the pencil, the black socket, Output Volt 2-90V DC should be used as cathode connection.

2. Operation

Using the Controls

Front Panel Controls

Name	Function	Name	Function												
<p>①</p> <p>O/CHARGE I/PUMP SWITCH</p>	When switched to O/Charge the batteries are being charged if the apparatus is connected to mains. When switched to I/Pump the pump motor starts operating.	<p>⑦</p> <p>AMMETER</p>	For reading the amperage during the polishing/etching process. The range of 0-2 is used for polishing. The range of 0-0.2 A is used for etching. Changing between the ranges is done automatically when the timers switch from polishing to etching.												
<p>②</p> <p>START BUTTON</p> 	Push button for starting the polishing/etching process.	<p>⑧</p> <p>CHARGE INDICATOR</p>	<p>Indicates the battery charging level. When connected to 42V AC possible readings are as follows:</p> <table border="1"> <thead> <tr> <th>Battery Charging Level</th> <th>O/ Charge</th> <th>I/ Pump</th> </tr> </thead> <tbody> <tr> <td>Medium charged</td> <td>Yellow</td> <td>Green</td> </tr> <tr> <td>Fully charged</td> <td>Green</td> <td>Green</td> </tr> <tr> <td>Not Charged</td> <td>Red</td> <td>Red</td> </tr> </tbody> </table>	Battery Charging Level	O/ Charge	I/ Pump	Medium charged	Yellow	Green	Fully charged	Green	Green	Not Charged	Red	Red
Battery Charging Level	O/ Charge	I/ Pump													
Medium charged	Yellow	Green													
Fully charged	Green	Green													
Not Charged	Red	Red													
<p>③</p> <p>ADJUSTING KNOB FOR POLISHING VOLTAGE</p>	Adjusting knob for setting the voltage for polishing. It can be set in the range of 10-90 V. The voltage range of 60-90 volt DC is marked with a voltage danger warning.	<p>⑨</p> <p>VOLTAGE INDICATOR FOR PUMP MOTOR</p>	<p>Indicates whether the pump is supplied with sufficient voltage:</p> <p><i>0</i>: The apparatus is switched off.</p> <p><i>Yellow</i>: Pump motor in operation at sufficient voltage.</p> <p><i>Red</i>: Pump motor blocked due to overload of pump motor or too low voltage. Electrical circuit switched off due to overload.</p>												
<p>④</p> <p>ADJUSTING KNOB FOR ETCHING VOLTAGE</p>	Adjusting knob for setting the voltage for etching.	<p>⑩</p> <p>POLISHING INDICATOR</p>	The Polishing Indicator lights red when the polishing process is running and is switched off when the process is finished.												
<p>⑤</p> <p>POLISHING TIME</p>	Adjusting knob for setting the time of polishing. When set to ∞ the process continues until the knob is set to zero.	<p>⑪</p> <p>ETCHING INDICATOR</p>	The Etching Indicator lights red when the etching process is running and is switched off when the process is finished.												
<p>⑥</p> <p>ETCHING TIME</p>	Adjusting knob for setting the time of etching. When set to ∞ the process continues until the knob is set to zero.	<p>⑫</p> <p>SOCKET PUMP/ CATHODE</p>	For connecting the pump motor and cathode to the control unit.												

Name	Function	Name	Function
<p>⑬</p> <p>SOCKET ANODE</p>	For connecting the anode wire to the control unit.	<p>⑭</p> <p>SOCKET FOR EXTERNAL ETCHING</p>	For connecting the cathode when external etching operations should be carried out.

Safety Functions

Cut-out of Pump Motor

In case of higher load than 1 A the cut-out function is activated. The Voltage Indicator for Pump Motor will light red. The apparatus can be restarted when the error has been eliminated.

Current Limiter Function of the Polish/Etch Circuit

In case of a voltage range of 70-90 V the current is limited corresponding to the maximum effect of the transistors (140 VA).

Cut-out of Polish/Etch Circuit

In case of a voltage range of 2-70 V at 2 A the cut-out function is activated and an alarming signal is heard. When the error has been corrected the apparatus can be restarted. In case of overload the circuit is cut out and the Voltage indicator for Pump Motor lights red. The apparatus can be restarted after cooling when the indicator lights yellow.

Cut-out in case of under-voltage

Under-voltage may occur using insufficiently charged batteries. The cut out function for under-voltage will disconnect the process and an alarm signal is heard.

Mounting the Container

- Place the pump in the container.
- Press down the black knobs and turn them until the markings turn towards the pump.
- Check that the container is firmly locked.

Dismounting the Container

- Press the two black knobs down and turn them so that the markings turn away from the pump, a click is heard and the pump is released from the container.

Mounting the Pump-Container Unit

- Lower the pump-container unit into the cabinet.
- Press the unit down until a click is heard. The pump-container unit is now locked in the cabinet.

Dismounting the Pump-Container Unit

- Press the two black locks to the sides from the pump in direction of the cabinet wall, a click is heard and the pump-container unit is dismounted and can be removed from the cabinet.

Filling Electrolyte into the Container

- Before filling electrolyte into the container, the section, Safety Precautions should be carefully observed.
- Remove the pump and container from the cabinet.
- Remove the pump from the electrolyte container.
- Pour the electrolyte into the electrolyte container. The level must be between the minimum (350 ml) and maximum (750 ml) marks.
- Mount the container on the pump and mount the unit in the cabinet of Movipol-3.

Emptying the Container of Electrolyte

- Before emptying the container of electrolyte, the section, Safety Precautions should be carefully observed.
- Remove the pump-container unit from the cabinet and remove the container.
- If the electrolyte is to be used again pour the electrolyte carefully back into the electrolyte bottle. Use the funnel supplied with the apparatus.
- If the electrolyte is not to be used again it should be filled into a container suitable for disposal. Please observe the local regulations and instructions for disposal.

IMPORTANT

Funnel, gloves, ventilation as well as all other prescribed equipment must be used during filling/emptying.

Mounting the Polishing Chamber

- Press the polishing chamber onto the pencil until it fits tightly.

A polishing chamber has a limited lifetime as the plastic hardens because of contact with the electrolyte. The chamber should be renewed regularly.

Starting the Polishing/Etching Process

- Make sure that the correct type and amount of electrolyte has been filled in the electrolyte container.
- Connect the anode to the sample by means of the clamp or magnet supplied with the apparatus.
- Set the polishing voltage and polishing time.
- Set the etching voltage and etching time.
- Activate the pump motor by switching O/Charge - I/Pump to bottom position.
- Take the pencil from the pencil holder and place it on the spot which is to be prepared.
- Press the pencil against the surface ensuring that the polishing chamber is in total contact with the surface.
- Keep the pencil perpendicular (at a angle of 90°) to the surface. The electrolyte will now begin to circulate.

IMPORTANT

Take care that no air is drawn around the periphery of the polishing chamber.

- Press the button START \diamond . A "sensor circuit" records whether there is electric connection from the cathode in the pencil through the electrolyte to the anode.
- If the electric circuit is closed the timer starts the process.

IMPORTANT

Only press the start button when the noise from the pump indicates that the flow through the pencil is constant.

- Read the ammeter to check if the correct amperage has been obtained.

Stopping the Polishing/Etching Process

- When both the Polishing and Etching Indicators are switched off the process has been completed and the pencil should be placed in the pencil holder.
- When the pencil has been placed in the holder, switch off the pump motor by switching the O/charge - I/Pump to O/Charge.
- Clean the polished spot with water/alcohol and dry carefully.
- The polishing/etching process can be stopped at any time during operation by setting the adjusting knob for polishing/etching time to zero.

External Etching

After polishing with one electrolyte or after mechanical polishing, the sample can be etched using a different electrolyte and the external etching connection. (External etching accessory required)

- Connect the anode to the sample using the clamp or magnet supplied with Movipol-3.
- Connect the tongs for external etching to the black socket on the Movipol-3.
- Set the etching voltage and etching time.
- Pick up a ball of cotton wool with the tongs.
- Dip the cotton wool ball in a suitable electrolyte (e.g. 10% oxalic acid for stainless steel).
- Start the pump and then the polishing/etching process.
- Swab the specimen surface with the cotton wool, and observe the movement on the Ammeter.
- Clean the polished spot with water/alcohol and dry carefully.

Cleaning the Container and the Pump

The electrolyte container and the pump system must be cleaned with water after use.

- Empty the container of electrolyte.
- Fill the container with water and start the pump by switching O/Charge - I/Pump to bottom position to circulate water from the electrolyte container through the pump system.

Voltage Examples

Due to graphite inclusions, electrolytic polishing is not ideal, but it is satisfactory for "testing on the spot". The metallic structure develops excellently and only the graphite appearance will be less successful. Due to the softness of the graphite this may be helped by rubbing the polished surface with a piece of wet cotton wool, perhaps with AP-Alumina Suspension FF.

Cast Iron. 3.4% C, 0.29% Si, 0.6% Mn, 0.3% P

Pregrinding with Transpol	Up to grit 400
Electrolyte	AC-2
Polishing voltage	40V
Polishing current	0.5A
Polishing time	15 sec

Duplex Steel

Pregrinding with Transpol	Up to grit 400
Electrolyte	A-2
Polishing	
Voltage	35V
Current	0.7A
Time	15 sec
Etching	
Voltage	8V
Current	0.2A
Time	12 sec

*Commercial Pure Aluminium
99.5% Al*

Pregrinding with Transpol	Up to grit 400
Electrolyte	A-2
Polishing voltage	85V
Polishing current	0.3A
Polishing time	20 sec

*Movipol-3
Instruction Manual*

Stainless Steel

Pregrinding with Transpol	Up to grit 400
Electrolyte	E-5
Polishing voltage	80V
Polishing current	0.6A
Polishing time	10 sec

External etch in 10% Oxalic acid or 10% (NH)₂S₂O₃

Low Carbon Steel 0.1% C

Pregrinding with Transpol	Up to grit 400
Electrolyte	A-2
Polishing	
Voltage	60V
Current	0.8A
Time	10 sec
Etching	
Voltage	5V
Current	0.18A
Time	5 sec

High Carbon Steel

Pregrinding with Transpol	Up to grit 220
Electrolyte	A-2
Polishing	
Voltage	60V
Current	1.0A
Time	15 sec
Etching	
Voltage	3-6V
Current	0.013A
Time	3 sec

Titanium Alloys

Pregrinding with Transpol	Up to grit 220
Electrolyte	A-3
Polishing voltage	50V
Polishing current	1.0A
Polishing time	15 sec

3. Safety precautions

Electricity

In some working areas a maximum safety voltage is specified. Such a maximum range would typically be 42V AC (60 DC). Movipol-3 is used with a 42V AC transformer through a 54V DC re-chargeable battery set. Movipol-3 can be used according to these restrictions. Alternatively Movipol-3 may work on batteries only.

WARNING

The polishing circuit may generate voltages up to 90V DC. The Adjusting knob for Polishing Volt indicates the critical range from 60 to 90V in yellow.

Electrolytes in General

When using electrolytes all necessary safety measures should be taken. Apart from the safety precautions listed in the following, a material Safety Data Sheet for each of the electrolytes delivered by Struers can be supplied upon request.

- It is essential that the user(s) is/are fully instructed in the work procedure of these electrolytes.
- Movipol-3 is designed for use with electrolytes recommended by Struers. Other electrolytes, e.g. electrolytes containing strong bases or acids, may harm the construction or endanger the safety of the user.
NOTE that many electrolytes contain alcohol or other inflammable solvents. Make sure that all safety precautions are followed when using such electrolytes.
- When cleaning the apparatus after use, make sure that no electrolyte is allowed to dry and/or crystallise inside the apparatus or on the polished material.
- It is essential that the user is fully trained in the use of Movipol-3 and electrolytes involved.

Perchloric Acid in Particular

Electrolytes from Struers with the prefix A consist of approx. 1 l stock solution to which 15 to 90 ml perchloric acid (60%) should be added. Before mixing the perchloric acid to the stock solution it is of great importance that the following safety precautions are carefully observed. Furthermore a Material Safety Data Sheet for the perchloric acid in question can be supplied upon request.

Training of Operators

- All personnel involved in the mixing, use, storage, transportation, and disposal of the electrolytes or its components, should be thoroughly trained in the precautions for handling perchloric acid.
- Great importance must be attached to precautions against inhalation of vapours of the solution or its components, against skin contact, mixing and overheating, and concerning storage and disposal methods.

Mixing the Solution

- Place the solvent/water mixture in a water bath with temperature control. Carefully add the perchloric acid to the solvent/water mixture, stir continuously.
- The mixing must take place in a ventilated chemical fume hood designed for perchloric acid use.
- The operator(s) must use the listed protective clothing or devices: full face shield or splash goggles, rubber gloves and lab coat or coveralls.
- Avoid the use of any combustible or carbonaceous containers, reaction vessels, spill pans, storage shelves, or materials of this type when dealing with the acid.

Storage of Perchloric Acid or Solution

- Do not permit any acid to crystallise on bottle necks, caps, or anywhere else.
- Store in secure, cool, and ventilated area with metal, glass, or ceramic spill catch pan.
- Store away from other chemicals, combustible and organic materials.
- Do not, under any circumstances, permit solutions to dry out.

Fire and Explosion Hazards

- 60% perchloric acid is a strongly corrosive and oxidising product. Heating may cause an explosion and contact with combustible material may cause fire.
- Fire fighting should be done from a protected location. Extinguish with water spray only. Do not use dry chemicals or carbon dioxide.
- Do not produce anhydrous perchloric acid, either from its salts or from aqueous solutions, e.g. by heating with high boiling acids or dehydrating agents, such as sulphuric acid or phosphorous pentoxide. In addition to spontaneous explosion, the anhydrous acid explodes instantaneously on contact with oxidizable organic materials.
- The use or storage of perchloric acid should be limited to quantities less than 500 g per hood.

Disposal

Follow local regulations for disposal of spillage and waste. Dilution and/or neutralisation are the normally recommended methods of disposal of the electrolyte.

**Working with Movipol-3
in the field**

Dangerous Localities

Certain safety measures must be taken when working with Movipol-3 in the field as some working locations can be dangerous. This is the case when working in:

- Closed containers, e.g. boilers, tanks and fractionating columns.
- Narrow areas such as shafts, pits, pipes and tunnels.

Safety Measures

Working in these above mentioned locations can be dangerous as poisonous substances can be accumulated and as there can be too little oxygen in the air.

Consequently it is recommended that tanks or other narrow areas are not opened or entered without written permission.

Necessary preparations and safety precautions should be mentioned in the permission and it must be issued by the manager in charge of safety for the actual process line or plant.

No working process should begin without this permission. The safety precautions and work processes stated for the place of work should be followed. In case of doubt regarding the contents of the authorisation, the issuer should be consulted.

Dangerous situations

Dangerous situations can occur when working in locations with choking or poisonous gases. In order to make sure that the air can be breathed without risk it must be analysed. Electric insulation defects may also involve a safety risk.

However, safe working procedure may be carried out if adequate safety precautions are taken before and during the work process.

Reference Guide

Table of Contents	Page
1. Consumables and Accessories	
Accessories for Movipol-3	17
Supplementary Equipment	17
2. Trouble-shooting	
Electrolyte	18
Flow	18
Pencil and Cathode	18
Polishing Chamber	18
Power Supply	18
3. Maintenance	
Pump and Electrolyte Container	19
Every day	19
Cabinet	19
Every day	19
Every week	19
Once a year	19
Exchange of Batteries	19
4. Technical Data	20

1. Consumables and Accessories

Accessories for Movipol-3

Specification	Cat. No.
<p><i>42 Volt Charging/supply Transformer for Continuous Use</i> for Movipol-3. Insulating transformer with outdoor casing. 1 x 115 V / 50 - 60 Hz 1 x 230 V / 50 - 60 Hz</p>	03926521 03926533
<p><i>Polishing Chambers</i> flexible type, for highly curved surfaces Internal diameter approx. 9 mm, external diameter approx. 13.5 mm, 10 pcs.</p>	03926904
<p><i>Batteries</i> set of 9 pcs.</p>	03926901
<p><i>Tongs with wire and plug</i> for external etching.</p>	325MP062

Supplementary Equipment

Specification	Cat. No.
<p><i>Microscope PSM-2</i> portable microscope. Complete with 10 x eyepiece, 10 x objective, lamp house and carrying case with battery. Batteries are not included</p>	04286101
<p><i>Transpol-2</i> portable apparatus for non-destructive grinding and polishing. Complete with flexible cable, DC-motor, straight handle, right angle handle. 1 x 120 V / 50-60 Hz 1 x 230 V / 50-60 Hz</p>	05956122 05956137
<p><i>Transcopy</i> reflective metallographic replica, set with 50 pcs. Complete with Transcopy fluid, replica foils 20 x 30 mm and 25 microscope slides.</p>	40200008

2. Trouble-shooting

Note:

The Figures referred to in this manual can be found in the Spare Parts and Diagrams section.

Polishing defects can normally be corrected using the following trouble-shooting schedule:

Electrolyte

- Check the age of the mixed electrolyte. It should not be more than 2-3 months old. The A-2 electrolyte has a particularly short lifetime (approximately 2 months) in a mixed condition.
- Check the number of polishings made with the electrolyte. It should not be used for several hundred polishings.
- Check that the right electrolyte - material combination is used.

Flow

- The pump suction is adjusted to 300 "20 mmhg at the factory (some materials, such as Nimonic and Titanium alloys need a flow corresponding to 400-450 mmhg). The suction should be measured with a vacuum meter using the suction hose nipple (fig. 7.3).
- Check that all connections are tight and that air is not sucked into the stream of electrolyte.
- The pump suction (flow) can be adjusted by means of the screw (fig. 7.2). Screw inward increases the flow. Screw outward decreases the flow.

Pencil and Cathode

- Check that the cathode (fig. 9.100) is even with, or retracted approximately 0.5 mm from the edge of the pencil tip without polishing chamber.
- Check that the preset voltage can be established between the cathode and the specimen (the anode).
- The resistance can be very high in an electrolyte which is incorrectly mixed.

Polishing Chamber

- Check that the polishing chamber is properly mounted and still flexible. If not flexible, exchange the chamber; see section INSTALLATION.

Power Supply

- Check that the batteries are charged.
- Check that the transformer voltage is 42V "10%.

3. Maintenance

Pump and Electrolyte Container

Every day

The electrolyte container must be emptied and rinsed with water after use. The pump system must also be rinsed out. To do this, circulate water from the electrolyte container through the pump system.

Cabinet

Every day

- Electrolyte spilled on the front plate or other parts of the cabinet should be wiped off.

Every week

- Stains of metal oxide can be removed by means of the cleansers and sealants used for cars and fiberglass boats.

Once a year

- Treat the cabinet surface with a surface sealant once a year in order to preserve it.

For general cleaning purposes use a cloth and a weak detergent.

Exchange of Batteries

The batteries supplied will last for several years provided they are used in a correct way. In case of a full discharge, however, the batteries may be damaged and will have to be exchanged.

A battery kit (9 pcs.) can be ordered from your Struers representative; see section CONSUMABLES AND ACCESSORIES.

IMPORTANT

Connections must be done correctly (+ pole to ÷pole); see fig. 11.
Never leave damaged batteries in the apparatus for an extended period.

Follow the procedure for changing the battery fuse; see fig. 11.

- Remove the 3 rails (fig. 11.3, 6 screws (fig. 11.4)).
- Disconnect the old batteries and remove them.
- Place the new batteries in exactly the same position as the old ones.
- Replace the 3 rails and the cabinet.
- Charge the batteries; see section INSTALLATION.

4. Technical Data

Power supply to transformer	220V/110V (AC), 50/60Hz
Power consumption	Max. 20 VA
Internal circuit fuse	5A
Internal battery fuse	10A
Internal transformer fuse Automatic internal thermal protection of electric circuit	2.5A
Noise Level	Approx. 52 dB (A), measured at a distance of 1 m/39.4" from the machine.
Polishing power	Max. 140W
Polishing area	Approx. 9 mm dia.
Dimensions	Width: 410 mm Depth: 200 mm Height: 280 mm
Weight	8.8 kg (including batteries)

Quick Reference Guide

Filling of the electrolyte

- The surface must be plane ground, e.g. with a Transpol apparatus.
- Choose the correct electrolyte for the test material.
- Pour the electrolyte into the electrolyte container. The liquid level must be between the minimum and maximum marks.
- Place the pump part in the container.

Mounting the electrolyte container

- Lower the entire pump-container assembly into the cabinet part.
- Press on the top part until the assembly locks.
- Insert the multi plug in the socket on the front plate and tighten the nut.

Activating the process

- Adjust the two voltages on the scale to a voltage corresponding to the desired amperage.
- Adjust polishing and etching time.
- Switch on the pump.
- Take the pencil from the pencil holder.
- Press the pencil against the prepared surface of the test piece, ensuring that the polishing chamber is in total contact with the surface.
- Keep the pencil perpendicular (at a 90° angle) to the surface during the whole process. The electrolyte starts to circulate. Take care that no air is drawn in around the periphery of the polishing chamber.
- Press the start button and the process (polishing/etching) will begin.

IMPORTANT

Only press the start button, when the noise from the pump indicates that the flow through the pencil is constant.

Surveillance of the process

- Read the ammeter to check if the right amperage has been obtained.
- When both LED-indicators for Polishing and Etching switch off, remove the pencil from the test piece.

Stop the process

- Place the pencil in the pencil holder. Press it down carefully.

Do not switch off the pump before the pencil is placed in the pencil holder

- When the pencil has been securely placed in its holder, switch off the pump.

Cleaning the apparatus

- Clean the polished spot with water/alcohol and dry carefully.
- Pour the electrolyte back into the storage bottle. Use a funnel.

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Date of Release 16.11.2005
updated F6.FG2009

Movipol-3
Spare Parts and Diagrams

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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
<p>Some of the drawings may contain position numbers not used in connection with this manual.</p>

Movipol-3
Spare Parts and Diagrams

Pos. and Description

Cat. No.


Fig. 1 Movipol-3, Complete

Movipol-3
Spare Parts and Diagrams

Fig. 2 Movipol-3, Details
Drawing 13920003

Pos. and Description	Cat. No.
5 3-pole socket 99-0211-30-04.....	392MP065
5 Socket cover 08-0201-28-00	392MP066
6 Hub for handle	392MP007
6 Handle bushing.....	392MP008
6 Collar bushing, 1 pc. Disc 15 x 21 x 2.5 mm, 1 pc. Screw MC 5 x 8 SS, 1 pc. Disc no. 8, 1 pc. Disc CW 10, 1 pc.	
Korrex button SW4, 1 pc.	392MP011
7 Fixture for shoulder strap	392MP009
8 Handle, complete	392MP003
8 Stanchion	392MP010
9 Lock-pin, 1 pc. Disc M8, 1 pc. Top nut M8, 1 pc.	392MP012
10 Screw cap EFV 3/4".....	392MP016
11 Lock washer	392MP017
12 Cabinet, upper part	392MP002
13 Base	392MP001
14 Stub for flex hose	392MP018
Container holder	392MP019
Spring suspension.....	392MP020
Spring.....	392MP021
Clamp, complete	130MP038
Magnet, complete.....	130MP039
Anode wire, complete.....	392MP070


Pos. and Description

Cat. No.

Fig. 3 Supply transformer
Drawing 14440002

3-pole transformer plug 692-2-99-0210-29-04392MP0064

Pos. and Description

Cat. No.

Fig. 5 Battery fuse exchange

ENGLISH

Battery Fuse Exchange

- Remove the pump assembly including the electrolyte container
- Lay down Movipol-3 on the side.
- Unscrew the 4 assembly screws ① in the bottom frame.
- Slide the bottom frame to the normal horizontal position.
- Exchange the 10A fuse ② in the fuse holder next to the batteries.
- Reassemble the cabinet.

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

Pos. and Description


Cat. No

2	Fuse 10AT (10 pcs.).....	392MP075
3	MOVAB, Set of batteries	03926901

Movipol-3
Spare Parts and Diagrams

Fig. 6 Pump and Electrolyte Container
Drawing 13920202 (1 of 2)

Pos. and Description	Cat. No.
1 ELMA wing button 023-3325 ELMA lid 040-3625	392MP006
3 Eccentric, 1 pc. Axle, 1 pc. Comp. spring, 1 pc. Disc M4, 1 pc. Nut M4, 1 pc.....	392MP013
4 Locking pawl Comp. spring.....	392MP014
5 Pencil holder	392MP026
6 Pump assembly platform.....	392MP027
8 Locking plate.....	392MP015

Movipol-3
Spare Parts and Diagrams

**Fig. 7 Pump and Container,
 Cross Section**
Drawing 13920202 (2 of 2)

Pos. and Description	Cat. No.
Pump, complete (with pencil).....	R3920202
Pump, complete (without pencil)	13920202
3 Return pipe socket.....	392MP029
4 Hose connection	13920155
5A Pump house with return pipe	
Hose 1/8" x 20 mm	
Screw M3 x 8 + nut M3	
Screw M4 x 20.....	392MP036
5B Return pipe	13920251
6 Glass disc	392MP037
7 Glass disc spring.....	392MP038
Vane wheel	
Vane wheel roller, 3 pcs.	392MP039
8 Electrolyte container	392MP040
9 Handle 1307/20B-M5.....	392MP041
10 Motor flange.....	392MP044
11 Spacer bushing	
O-ring 5.3 x 2.4	392MP045
12 Sling disc.....	392MP046
Gasket flange.....	392MP047
13 Pump axle	392MP048
13 Cyl. Pin 2x10 A4 (30 pcs)	RZS11210
15 Set of seals	
14 a V-ring VS 008, 1 pcs.	
15 b O-ring 108 x 4, 1 pcs.	
9.20 c O-ring 15 x 2, Viton, 1 pcs.	
9.40 d O-ring 4.48 x 1.78, Viton, 4 pcs.	392MP78
17 PVC hose 5/32" x 20 mm, 1 pcs.....	RNU19305
20 Air tube, complete	13920258
21 Axle house.....	392MP035
22 Motor screen lid.....	392MP050
23 Motor screen	13929005
24 Motor	392MP042
25 Feed pipe, complete	392MP034
26 Spring.....	13140108

Movipol-3
Spare Parts and Diagrams

*(Fig. 7, Drawing 13920202 (2 of 2),
Pump and Container, Cross Section)*


Fig. 8 Control panel
Drawing 13920382

Pos. and Description	Cat. No.
1 Switch with cover U 851.....	213MP002
2 Push button 18235,	392MP071
Cover U 224.....	391MP006
5 Ammeter KM6	392MP053
Ammeter bracket	392MP069
10 4-pole socket 09-0112-00-04.....	2XM00112
Socket cover 08-0044-00-12	2XM90044
Control PCB 205-P-395/1 A2	392MP060
Tape cable	392MP061
SMPS-board 205-P-3041 A1	392MP062
11 Banana socket, red	2XG08751
Tongs with wire and plug	325MP062
12 Banana socket, black	2XG08750
Terminal strip 25.600.1253	392MP057
13 Front plate, Movipol-3	13922901
14 ELMA button 020-4425, 1 pcs. ELMA lid 040-4025, 1 pcs. ELMA wing disc 041-4025, 1 pcs. V-ring VS-0007NBR60, 1 pcs.	392MP005




Fig. 9 Pencil, Cross Section
Drawing 13920152

Pos. and Description	Cat. No.
Pencil, complete	392MP022
2 Polishing Chamber, short 14mmxr8	13920323
10 Pencil casing	13920156
20 O-ring Ø15.00x2.00, part of Set of Seals.....	392MP078
30 Pencil pipe stub.....	13920153
40 O-ring Ø4.48x1.78 (10 pcs.)	375MP037
40 O-ring Ø4.48x1.78, part of set of sealing	392MP078
70 Pencil hoses, complete with Cathode wire, 1 pcs. Hose 1.8" x 1200, 2 pcs. Flex hose, 1 pcs.	392MP033
Red cathode wire, complete with banana plugs.....	392MP068
90 Pencil handle.....	13920158
100 Cathode pipe with O-rings.....	392MP028
110 Cover A ø14.5 ELMA	2GD53025

Pos. and Description

Cat. No.


Fig. 10 Shoulder Strap

4 Strap.....392MP074

Fig. 11 Batteries
Drawing 13920018

Pos. and Description	Cat. No.
2 Fuse holder	392MP067
2 Battery Fuse, 10 A (10 pcs)	392MP075
3 Rail (3 pcs).....	392MP072
4 Phillips screw & Disc M4, 7 pcs.	392MP073

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




Herewith declares that

<i>Product Name:</i>	Movipol-3
<i>Type No:</i>	392
<i>Machine Type:</i>	Electrolytic polishing apparatus

is in conformity with the provisions of the following directives:

EMC-Directive 2004/108/EEC according to the following standard(s):
 EN61000-6-1:2002, EN61000-6-3:2002.

Supplementary Information The equipment complies with the American standards:
 FCC part 15, subpart B

The above has been declared according to the global method, module A

Date: 04.12.2007

 Christian Skjold Heyde,
 Vice President, R& D and Production, Struers A/S

Hersteller Struers A/S
Pederstrupvej 84
DK-2750 Ballerup, Danmark
Telefon +45 44 600 800

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Date: 04.12.2007

Christian Skjold Heyde,
Vice- President, R& D et Production, Struers A/S

On-site materialographic preparation and examination

Non-destructive preparation

Conventional materialographic preparation entails cutting a sample from a material. This method is not always appropriate as it destroys the component which is being examined.

Many applications require the ability to perform a materialographic preparation and examination without ensuing damage to the subject. The non-destructive method is used in these cases.

The usual preparation steps such as grinding, polishing and etching are performed on the part to be examined in a limited area where no damage can be done, whereupon the surface can be analysed.

On-site preparation is particularly suitable for quality checks in power stations, oil platforms, bridges, aircraft etc. The non-destructive preparation method enables the necessary inspection to be carried out on site.

Struers manufactures a complete range of portable equipment for non-destructive preparation, from basic grinding through mechanical and electrolytic polishing to electrolytic etching. The range also includes a portable microscope for on-site use and methods for producing replicas of the prepared surface.

Using replica methods makes it possible to perform the microscopic examination in the laboratory - under perfect working conditions.

The entire Struers non-destructive preparation range has been designed for field work. The equipment is thus as compact and as light as possible. It is easy to transport and can be used literally under any conditions.

Transpol-2

All the necessary accessories are contained in Transpol-2's carrying case.

Transpol-2 - Mechanical grinding and polishing under all working conditions

Transpol-2 is a portable grinding and polishing apparatus. It is extremely compact and is designed for field work under even the most difficult conditions. Transpol-2 fits into a practical shoulder bag with all the necessary accessories and weighs an approximate total of only 5 kg. Transpol-2 consists primarily of 2 parts: a control unit and a grinding and polishing unit. The control unit comprises the power supply and speed regulator unit (100-7000 rpm). The grinding and polishing unit is connected to the control unit by a soft, flexible cable. It consists of a direct current motor with a holder with an easily exchangeable rubber disc, upon which the grinding paper and polishing cloth are mounted.

Transpol-2 is available with two types of holder: one for straight mounting and one for right-angle mounting, thus making it possible to use Transpol-2 in even the most inaccessible places.

Grinding

Transpol-2 uses 32 mm diameter SiC grinding paper discs in grit sizes 60, 120, 240 and 500#. Transpol-2 also includes a special flapper wheel for grinding of extremely uneven surfaces.

Polishing

Transpol-2 uses DP-Dac, DP-Dur, DP-Mol, DP-Nap and OP-Felt polishing cloths, all 32 mm diameter, and DP-Paste, DP-Stick or DP-Spray. Transpol-2 is furthermore delivered with a special felt cone for polishing of uneven and curved surfaces.

Movipol-3

Due to the special design of the polishing pencil, electrolytic polishing with Movipol-3 is quick and uncomplicated

Movipol-3 - Electrolytic polishing and etching in the field

Movipol-3 is a portable electrolytic metal polishing and etching apparatus. It is extremely compact and robust and can be used anywhere. Electrolytic preparation is a particularly fast and efficient method of non-destructive metallographic preparation. It can be applied directly to critical surfaces and is widely used for metal safety inspection, especially for critical parts of larger units, e.g. for welds and other joints. The method effectively reveals changes in the microstructure of the metal and allows to prevent possible damage because of cracks and leaks.

Automatic polishing and etching

Movipol-3 consists of a compact unit comprising the power supply, control unit, pump motor, electrolyte container and polishing unit. When the polishing unit is pressed against the metal surface, an electrolytic reaction is started. After 5-10 seconds of use, the prepared area is as reflective as a mirror and etching can commence.

PSM-2 Transcopy

The entire process is automatic: after completion of the polishing, Movipol-3 proceeds to the etching process. When this is finished, the unit emits an audio signal. In less than one minute a surface is obtained which is ready for analysis, either on-site or by means of replica for laboratory examination.

Complete Freedom

Movipol-3 is equipped with rechargeable batteries and supply transformer. In other words, the apparatus is as suitable for stationary laboratory work as it is for field work, depending on the requirements.

PSM-2 and Transcopy constitute an indispensable tool for non-destructive materialography.

PSM-2 - Fast and easy microscopy outside the laboratory

PSM-2 is a small, portable microscope for use in the field. It is battery-powered and provides magnifications of 100 x, 200 x and 400 x. PSM-2 comes in a shoulder bag which also contains the batteries.

PSM-2 is particularly suitable for preliminary examination of the prepared surface. A replica is recommended if a more detailed analysis is required. PSM-2 is highly suitable for preliminary replica examinations when used together with Transpol-2 and Movipol-3.

Transcopy Replica Foil


Transcopy Replica Foil is used to make a replica of any polished and etched surface. It is especially used if on-site analysis using a portable microscope is not sufficient or if micrographs of the prepared structure are required. Transcopy Replica Foil creates a permanent replication of microstructures, cracks and defects for future analysis and comparison in the laboratory.

Function

Transcopy Replica Foil consists of reflecting plastic film with a self-adhesive back. A replica is made by applying Transcopy fluid to the film and pressing it firmly onto the prepared and etched surface. A slightly over etched structure results in a higher contrast of the replica. After 4-5 minutes the foil is removed from the surface. By removing the cover paper, the replica can be adhered to a plain plate and then examined under an optical microscope.

Transcopy

With Transcopy, Replica foil, replicas can be taken on the spot for subsequent examination in the laboratory.

RepliSet

The **RepliSet** system is designed to produce an exact 3D copy of a surface. It is used for non-destructive testing and field applications allowing the structure or irregularities on critical components to be examined and measured under laboratory conditions.

RepliSet is a specially formulated fast curing two-part silicone rubber with a good releasing ability for flexible high-resolution 3D replicas, which behave like a metallic surface when examined in an optical microscope. RepliSet compounds are supplied in cartridges and are dispensed using a hand-operated dispensing gun. The cartridges contain both polymer and curing agent, which are automatically mixed in a disposable static-mixing nozzle during application to the surface. The application system offers superior and fast results, regardless of the conditions.

RepliFix is a less advanced parallel to RepliSet. RepliFix and RepliSet are designed to bond together. The two components are mixed and applied by hand. RepliFix is used as support for RepliSet or as a stand alone product for moulding of surface shape for low tech applications.

A backing slide bonds to the RepliSet or RepliFix replica. The backing slide serves to maintain the original profile and ensures a flat back to the replica.


A replica can be taken from all metallic materials and most other solid materials like ceramics, plastics and glass. There are no size, shape or thickness limitations on the replicas that can be made and it is even possible to take replicas from otherwise inaccessible surfaces. The RepliSet system can produce replicas, which at the same time are dimensionally correct, with ultra fine detail reproduction and with a flat back.

An operator dispenses RepliSet onto a sample, and then peels off the cured replica.

*Pure copper.
Sample etched with cupric chloride and ammonia.
Magnification 100 x*

Sample Replica

Further to examination by optical microscopes, 3D examination can be carried out using non-contact measuring instruments such as laser measuring equipment or measuring projectors with 2D or 3D facilities. A replica of a cavity, for instance inner holes, can be examined using measurement and profile projectors. The replicas are suitable for 3D examination by SEM.

Metallographic applications

Typical applications are on-site non-destructive testing in connection with quality control, inspection and maintenance of power plants, oil platforms, bridges, aircraft, etc.

Water or high temperatures do not affect RepliSet so the application range is very broad.

Typical tasks are:

- Quality control of edges, corners, heights, angles, surface finish and other dimensions.
- Inspection of internal surfaces such as bolt hole threads.
- Detection and monitoring of pitting, corrosion, cracking, creep and wear.
- Metrology.
- Forensic investigation.

RepliSet 50 ml system.

Case for RepliSet

The RepliSet Case is designed for transportation and use of the 50 ml RepliSet system. It is made of aluminium and is at the same time elegant and sturdy. It can be carried as hand luggage on flights and is compact and sturdy to such an extent, that it can be taken to locations with narrow or difficult access.

The contents is either fixed by straps in the lid or placed in compartments in the two detachable foam rubber inserts. Each item has its fixed position. The user has access to all that is needed to perform a regular replication by just opening the lid of the RepliSet Case. The lower insert carries a small stock of consumables.

Transpol-2

Technical Data

Voltages	230 V, 50/60 Hz; 120 V, 50/60 Hz;
Speed	100-7000 rpm
Dimensions (without shoulder bag)	Width 180 mm
	Height 75 mm
	Depth 150 mm
Total weight incl. accessories	5 kg

Specifications

Transpol-2, complete with soft cable, DC motor 230 or 120 V, straight handle, right angle handle, flapper wheel, felt cone, 4 rubber discs	059561xx
Rubber disc, 30 mm dia.	40800055
SiC grinding papers, adhesive, bundle of 100	
Grit 60	40400049
Grit 120	40400129
Grit 240	40400130
Grit 500	40400131
DP cloths, adhesive, dia. 32 mm, bundle of 25	
DP-Dur	40500040
DP-Dac	40500145
DP-Mol	40500041
DP-Nap	40500042
OP-Felt	40500043
Flapper wheel, 80 grit, 40 mm dia. x 10 mm	40800053

Movipol-3

Technical Data

Voltage supply to transformer	220 V 50/60 Hz, 115V 50/60 Hz
Voltage supply from transformer	42 V AC
Input fuse	5 A
Battery fuse	10 A
Thermal overload protection	
Polishing power	Max. 140 VA
Polishing area	About 9 mm dia
Dimensions	Width 410 mm
	Height 280 mm
	Depth 200 mm
Weight with batteries, but without electrolyte	Net 8.80 kg

Specifications

Movipol-3, complete with batteries and supply transformer, 220 V or 115 V	039265xx
42 V charging/supply transformer for continuous use for Movipol 3, 230 V or 115 V	044462xx
Set of extra batteries (9 pcs.)	03926901
Polishing Chambers for Movipol-3, Flexible type, 10 pcs.	03926904

PSM-2

Technical Data

Batteries	4 x 1.5 V, type LR20	Length	210 mm (8.3")
Lamp (threaded base)	5 V, 6 W	Diameters	15-60 mm (0.6"-2.4")
		Weight	600 g

Specifications

Portable microscope PSM-2, complete with 10 x eyepiece, 10 x objective, lamp housing and carrying case, without batteries	04286101
Objective, magnification 20 x	04286901
Objective, magnification 40 x	04286902
Measuring set consisting of measuring eye piece and plate micrometer	04286903

Transcopy

Transcopy Kit

Set consisting of 40 ml Transcopy Liquid, 50 replica foils 20 x 30 mm (0.8" x 1.2"), 1 pipette, spray nozzle and 50 microscope slides 25 x 75 mm (1" x 3")	40900090
Transcopy Replica Foils, 50 pcs.	40900091
Transcopy Liquid, 40 ml	40900092

RepliSet

Technical Data

Resolution of cu red replica	Down to 0.1 micron
Shrinkage	Negligible
Tear Strength	15-20 kN/m ²
Temperature range for the surface to be examined	-10°C to +180°C
Life span of the finished replicas is practically indefinite provided they are stored according to the instructions.	
Content in static-mixing nozzle	1.1 ml in nozzle for 50 ml cartridge 9.3 ml in nozzle for 265 ml cartridge

Specifications

Replication system for non-destructive testing of a microstructure or a 3D structure. Fast curing two-part silicon rubber compound for flexible high-resolution 3D replicas. For the 50 ml system, the hand-operated dispensing gun (40900066) and the static mixing nozzles (40900088) are used in combination with the 50 ml cartridges. For the 265 ml system, the hand-operated dispensing gun (40900065) and the static-mixing nozzles (40900056) are used in combination with the 265 ml cartridges.

Struers A/S
 Pederstrupvej 84
 DK-2750 Ballerup, Denmark
 Phone +45 44 600 800
 Fax +45 44 600 801
 struers@struers.dk

Specifications

Cat. no.

RepliSet-F1

Particularly useful for replicating horizontal or sloping surfaces in low temperature conditions or where rapid results are required. Fluid rapid curing compound with working life of 0.5-1 min. and curing time of 4 min. at 25°C
 1 cartridge of 50 ml
 5 cartridges of 50 ml
 2 cartridges of 265 ml

40900069
 40900047
 40900051

RepliSet-F5

General purpose material. Particularly useful for replicating horizontal or sloping surfaces in normal or high temperature conditions. Fluid fast curing compound with working life of 5 min. and curing time of 18 min. at 25°C
 1 cartridge of 50 ml
 5 cartridges of 50 ml
 2 cartridges of 265 ml

40900068
 40900046
 40900050

RepliSet-T1

Particularly useful for replicating vertical or overhead surfaces in low temperature conditions or where rapid results are required. Thixotropic rapid curing compound with working life of 0.5-1 min. and curing time of 4 min. at 25°C
 1 cartridge of 50 ml
 5 cartridges of 50 ml
 2 cartridges of 265 ml

40900071
 40900049
 40900053

RepliSet-T3

General purpose material. Particularly useful for replicating vertical or overhead surfaces in normal or high temperature conditions. Thixotropic fast curing compound with working life of 3 min. and curing time of 10 min. at 25°C
 1 cartridge of 50 ml
 5 cartridges of 50 ml
 2 cartridges of 265 ml

40900070
 40900048
 40900052

RepliSet-GF1

Replication system especially for comparative microscopy and metrology. Particularly useful for replicating horizontal or sloping surfaces and filling holes. Fluid rapid curing compound with working life of 0.5 - 1 min. and curing time of 4 min. at 25°C.
 1 cartridge of 50 ml
 5 cartridges of 50 ml

40900078
 40900076

RepliSet-GT1

Replication system especially for comparator microscopy and metrology. Particularly useful for replicating vertical or overhead surfaces. Thixotropic rapid curing compound with working life of 0.5 - 1 min. and curing time of 4 min. at 25°C.
 1 cartridge of 50 ml
 5 cartridges of 50 ml

40900079
 40900077

RepliFix

Specially formulated hand mixed fast curing two-part silicone rubber. Bonds to RepliSet. Particularly useful in combination with RepliSet for producing a rigid backing. It can be used directly for moulding of surface shape for profile measurement

RepliFix-2

For low temperature conditions or where rapid results are required. Working life of 2-3 min. and curing time of 10 min. at 25°C. Net 500 g

40900084

RepliFix-20

For high temperature conditions or for taking replicas of complicated geometry or large areas. Working life of 20 min. and curing time of 60 min. at 25°C. Net 500 g

40900086

Accessories

Hand-operated dispensing gun

For 50 ml cartridges of RepliSet
 For 265 ml cartridges of RepliSet

40900066
 40900065

Static-mixing nozzles

For 50 ml cartridges. 35 pcs.
 For 265 ml cartridges. 10 pcs.

40900088
 40900056

Nozzle Tips

For replicating flat surfaces. Fishtail spreaders, 10 mm width.
 To be mounted on 50 ml static mixing nozzle (40900088). 30 pcs.

40900089

For replicating small holes. Luer needle, 1 mm dia., 30 mm long. To be mounted on 50 ml static-mixing nozzle (40900088). 10 pcs.

40900060

For replicating larger holes. Flexible hose, 6 mm dia., 100 mm long. To be mounted on 50 ml static-mixing nozzle (40900088) 10 pcs.

40900061

Backing Slides

A flexible plastic slide, which bonds to the replica and ensures a flat back to the replica. For levelling of replicas to assist microscopic examination, as dimensional support for metrology and for well-ordered labelling, transport and storage of RepliSet replicas 26 x 76 x 1 mm, 50 pcs.

40900087

Backing paper

For RepliSet replication system. Bonds to the replica and facilitates labelling, handling and the levelling of replicas to assist microscopic examination. 60 x 70 mm. 100 pcs.
 A4 (210 x 297 mm), for cutting up to the required size. 10 pcs.

40900062
 40900063

Case for RepliSet 50 ml system

Aluminium case with room for all necessities for field applications.
 The contents of the RepliSet Case is ordered separately

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CANADA

Struers Ltd.
 7275 West Credit Avenue
 Mississauga, Ontario L5N 5M9
 Phone +1 905-814-8855
 Fax +1 905-814-1440
 info@struers.com

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All rights reserved. © Struers 200J.
Struers Limited
Struers Inc.
Struers GmbH