PSM-2

Portable Microscope
Instruction Manual

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Portable Microscope

Instruction Manual

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the *Serial No.* and *Voltage* on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2009.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801


User's Guide

Table of Contents	Page
1. Getting Started	
Description	2
2. Operation	
Mounting the base ring.....	3
Transcopy Replica	3
3. Measuring with PSM-2	4
4. Technical Data	4
5. Accessories	4
6. Spare Parts	15

1. Getting Started

Description

PSM-2 is a monocular, non-inverted microscope with a battery driven lamp for primary structural analyses in reflected light. It is suited for on the spot examination of large objects during both the building and the operation processes. It can be used both on external surfaces and in corners, cavities and tubes. Further, PSM-2 is intended for control of structural impression, carried out with Transcopy Replica.

- ① 10x objective
- ② Locking ring to prevent the objective from touching the sample
- ③ Detachable base ring for examination on plane surfaces
- ④ Detachable lamp with cable
- ⑤ Carrying case with battery cassette
- ⑥ Switch and socket
- ⑦ 10x ocular

2. Operation

- Place the batteries in the battery cassette according to the diagram.
- Place the lamp in the socket of the microscope and the plug on the lead into the socket of the battery cassette. Set the switch to "On".
- Check that the base ring (4) has been mounted on the microscope, when using this for examination of plane or slightly curved surfaces and when examining Transcopy Replica.

Mounting the base ring

- The base ring has a bayonet fitting closure. Move the base ring, with the arrow pointing to the lamp house, into the threaded part and turn clockwise until the base ring is fixed.
- To remove the base ring, turn anti-clockwise until the arrow points to the lamp house and take the ring out.
- Fix the locking ring (3) in a position preventing the objective from touching the sample. Then the microscope is ready for use. The focusing is done by adjusting the focusing screws in the middle of the microscope.
- When changing the objective the base ring must be removed and the objective set to the lowest position. The objective is then unscrewed and removed to allow a new one to be fitted (Only tighten with fingers). The objective is set to a higher position and the base ring remounted.

TIPS

- When using the microscope for internal examinations of corners, cavities, etc., the base ring must be removed.
- Take care that the objective is not damaged by touching the sample.

Transcopy Replica

If the microscope is to be used for examination of structural impressions (Transcopy Replica), this must be glued onto a slide and then placed on a table with the structural surface uppermost. The microscope is positioned above the replica, which can be moved for viewing.

Remember...

Do not forget to switch off the lamp after use.

3. Measuring with PSM-2

Remove the eye piece by loosening the small screw on the side of the tube. Instead mount the measuring eye-piece and tighten the screw lightly. Focus on the micrometer plate and bring the scale in the measuring eye piece above the scale on the micrometer plate by moving and turning PSM-2. A possible correction factor is calculated according to proportion between, for instance, the full scale of the eye piece and the indication to which it corresponds on the measuring plate.

4. Technical Data

Ocular: 10x magnification
Objective: 10x magnification
Batteries: 4 pcs. 1.5 V, LR20, AM1, Size D
Lamp: 5 V, 6 W
Diameter: 15-60 mm
Length: 210 mm
Weight: 600 g

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
6. Spare Parts

Pos.	Description	Spare Parts No.
1	Objective, 10x	428MP001
2	Locking ring with screw.....	428MP003
3	Base ring	428MP004
4	Lamp, complete, but without bulb	428MP005
5	Carrying case, complete	428MP007
6	Switch (in carrying case).....	428MP008
7	Eye piece (10x).....	428MP009
8	Optional eye piece (see accessories)	
9	Bulb	428MP006
10	Gear shaft with thumb wheels	428MP010
10A	Gear wheels (2 pcs).....	428MP020
10B	Gear shaft / pinion	428MP030
10C	Rack.....	428MP012
10	Gear shaft with thumb wheels	428MP010
11	Cord with socket	428MP011
12	Main body	428MP002

On-site materialographic preparation and examination

Non-destructive preparation

Conventional materialographic preparation entails cutting a sample from a material. This method is not always appropriate as it destroys the component which is being examined.

Many applications require the ability to perform a materialographic preparation and examination without ensuing damage to the subject. The non-destructive method is used in these cases.

The usual preparation steps such as grinding, polishing and etching are performed on the part to be examined in a limited area where no damage can be done, whereupon the surface can be analysed.

On-site preparation is particularly suitable for quality checks in power stations, oil platforms, bridges, aircraft etc. The non-destructive preparation method enables the necessary inspection to be carried out on site.

Struers manufactures a complete range of portable equipment for non-destructive preparation, from basic grinding through mechanical and electrolytic polishing to electrolytic etching. The range also includes a portable microscope for on-site use and methods for producing replicas of the prepared surface.

Using replica methods makes it possible to perform the microscopic examination in the laboratory - under perfect working conditions.

The entire Struers non-destructive preparation range has been designed for field work. The equipment is thus as compact and as light as possible. It is easy to transport and can be used literally under any conditions.

Transpol-2

All the necessary accessories are contained in Transpol-2's carrying case.

Transpol-2 - Mechanical grinding and polishing under all working conditions

Transpol-2 is a portable grinding and polishing apparatus. It is extremely compact and is designed for field work under even the most difficult conditions. Transpol-2 fits into a practical shoulder bag with all the necessary accessories and weighs an approximate total of only 5 kg. Transpol-2 consists primarily of 2 parts: a control unit and a grinding and polishing unit. The control unit comprises the power supply and speed regulator unit (100-7000 rpm). The grinding and polishing unit is connected to the control unit by a soft, flexible cable. It consists of a direct current motor with a holder with an easily exchangeable rubber disc, upon which the grinding paper and polishing cloth are mounted.

Transpol-2 is available with two types of holder: one for straight mounting and one for right-angle mounting, thus making it possible to use Transpol-2 in even the most inaccessible places.

Grinding

Transpol-2 uses 32 mm diameter SiC grinding paper discs in grit sizes 60, 120, 240 and 500#. Transpol-2 also includes a special flapper wheel for grinding of extremely uneven surfaces.

Polishing

Transpol-2 uses DP-Dac, DP-Dur, DP-Mol, DP-Nap and OP-Felt polishing cloths, all 32 mm diameter, and DP-Paste, DP-Stick or DP-Spray. Transpol-2 is furthermore delivered with a special felt cone for polishing of uneven and curved surfaces.

Movipol-3

Due to the special design of the polishing pencil, electrolytic polishing with Movipol-3 is quick and uncomplicated

Movipol-3 - Electrolytic polishing and etching in the field

Movipol-3 is a portable electrolytic metal polishing and etching apparatus. It is extremely compact and robust and can be used anywhere. Electrolytic preparation is a particularly fast and efficient method of non-destructive metallographic preparation. It can be applied directly to critical surfaces and is widely used for metal safety inspection, especially for critical parts of larger units, e.g. for welds and other joints. The method effectively reveals changes in the microstructure of the metal and allows to prevent possible damage because of cracks and leaks.

Automatic polishing and etching

Movipol-3 consists of a compact unit comprising the power supply, control unit, pump motor, electrolyte container and polishing unit. When the polishing unit is pressed against the metal surface, an electrolytic reaction is started. After 5-10 seconds of use, the prepared area is as reflective as a mirror and etching can commence.

PSM-2 Transcopy

The entire process is automatic: after completion of the polishing, Movipol-3 proceeds to the etching process. When this is finished, the unit emits an audio signal. In less than one minute a surface is obtained which is ready for analysis, either on-site or by means of replica for laboratory examination.

Complete Freedom

Movipol-3 is equipped with rechargeable batteries and supply transformer. In other words, the apparatus is as suitable for stationary laboratory work as it is for field work, depending on the requirements.

PSM-2 and Transcopy constitute an indispensable tool for non-destructive materialography.

PSM-2 - Fast and easy microscopy outside the laboratory

PSM-2 is a small, portable microscope for use in the field. It is battery-powered and provides magnifications of 100 x, 200 x and 400 x. PSM-2 comes in a shoulder bag which also contains the batteries.

PSM-2 is particularly suitable for preliminary examination of the prepared surface. A replica is recommended if a more detailed analysis is required. PSM-2 is highly suitable for preliminary replica examinations when used together with Transpol-2 and Movipol-3.

Transcopy Replica Foil


Transcopy Replica Foil is used to make a replica of any polished and etched surface. It is especially used if on-site analysis using a portable microscope is not sufficient or if micrographs of the prepared structure are required. Transcopy Replica Foil creates a permanent replication of microstructures, cracks and defects for future analysis and comparison in the laboratory.

Function

Transcopy Replica Foil consists of reflecting plastic film with a self-adhesive back. A replica is made by applying Transcopy fluid to the film and pressing it firmly onto the prepared and etched surface. A slightly over etched structure results in a higher contrast of the replica. After 4-5 minutes the foil is removed from the surface. By removing the cover paper, the replica can be adhered to a plain plate and then examined under an optical microscope.

Transcopy

With Transcopy, Replica foil, replicas can be taken on the spot for subsequent examination in the laboratory.

RepliSet

The **RepliSet** system is designed to produce an exact 3D copy of a surface. It is used for non-destructive testing and field applications allowing the structure or irregularities on critical components to be examined and measured under laboratory conditions.

RepliSet is a specially formulated fast curing two-part silicone rubber with a good releasing ability for flexible high-resolution 3D replicas, which behave like a metallic surface when examined in an optical microscope. RepliSet compounds are supplied in cartridges and are dispensed using a hand-operated dispensing gun. The cartridges contain both polymer and curing agent, which are automatically mixed in a disposable static-mixing nozzle during application to the surface. The application system offers superior and fast results, regardless of the conditions.

RepliFix is a less advanced parallel to RepliSet. RepliFix and RepliSet are designed to bond together. The two components are mixed and applied by hand. RepliFix is used as support for RepliSet or as a stand alone product for moulding of surface shape for low tech applications.

A backing slide bonds to the RepliSet or RepliFix replica. The backing slide serves to maintain the original profile and ensures a flat back to the replica.


A replica can be taken from all metallic materials and most other solid materials like ceramics, plastics and glass. There are no size, shape or thickness limitations on the replicas that can be made and it is even possible to take replicas from otherwise inaccessible surfaces. The RepliSet system can produce replicas, which at the same time are dimensionally correct, with ultra fine detail reproduction and with a flat back.

An operator dispenses RepliSet onto a sample, and then peels off the cured replica.

*Pure copper.
Sample etched with cupric chloride and ammonia.
Magnification 100 x*

Sample Replica

Further to examination by optical microscopes, 3D examination can be carried out using non-contact measuring instruments such as laser measuring equipment or measuring projectors with 2D or 3D facilities. A replica of a cavity, for instance inner holes, can be examined using measurement and profile projectors. The replicas are suitable for 3D examination by SEM.

Metallographic applications

Typical applications are on-site non-destructive testing in connection with quality control, inspection and maintenance of power plants, oil platforms, bridges, aircraft, etc.

Water or high temperatures do not affect RepliSet so the application range is very broad.

Typical tasks are:

- Quality control of edges, corners, heights, angles, surface finish and other dimensions.
- Inspection of internal surfaces such as bolt hole threads.
- Detection and monitoring of pitting, corrosion, cracking, creep and wear.
- Metrology.
- Forensic investigation.

RepliSet 50 ml system.

Case for RepliSet

The RepliSet Case is designed for transportation and use of the 50 ml RepliSet system. It is made of aluminium and is at the same time elegant and sturdy. It can be carried as hand luggage on flights and is compact and sturdy to such an extent, that it can be taken to locations with narrow or difficult access.

The contents is either fixed by straps in the lid or placed in compartments in the two detachable foam rubber inserts. Each item has its fixed position. The user has access to all that is needed to perform a regular replication by just opening the lid of the RepliSet Case. The lower insert carries a small stock of consumables.

Transpol-2

Technical Data

Voltages	230 V, 50/60 Hz; 120 V, 50/60 Hz;
Speed	100-7000 rpm
Dimensions (without shoulder bag)	Width 180 mm
	Height 75 mm
	Depth 150 mm
Total weight incl. accessories	5 kg

Specifications

Transpol-2, complete with soft cable, DC motor 230 or 120 V, straight handle, right angle handle, flapper wheel, felt cone, 4 rubber discs	059561xx
Rubber disc, 30 mm dia.	40800055
SiC grinding papers, adhesive, bundle of 100	
Grit 60	40400049
Grit 120	40400129
Grit 240	40400130
Grit 500	40400131
DP cloths, adhesive, dia. 32 mm, bundle of 25	
DP-Dur	40500040
DP-Dac	40500145
DP-Mol	40500041
DP-Nap	40500042
OP-Felt	40500043
Flapper wheel, 80 grit, 40 mm dia. x 10 mm	40800053

Movipol-3

Technical Data

Voltage supply to transformer	220 V 50/60 Hz, 115V 50/60 Hz
Voltage supply from transformer	42 V AC
Input fuse	5 A
Battery fuse	10 A
Thermal overload protection	
Polishing power	Max. 140 VA
Polishing area	About 9 mm dia
Dimensions	Width 410 mm
	Height 280 mm
	Depth 200 mm
Weight with batteries, but without electrolyte	Net 8.80 kg

Specifications

Movipol-3, complete with batteries and supply transformer, 220 V or 115 V	039265xx
42 V charging/supply transformer for continuous use for Movipol 3, 230 V or 115 V	044462xx
Set of extra batteries (9 pcs.)	03926901
Polishing Chambers for Movipol-3, Flexible type, 10 pcs.	03926904

PSM-2

Technical Data

Batteries	4 x 1.5 V, type LR20	Length	210 mm (8.3")
Lamp (threaded base)	5 V, 6 W	Diameters	15-60 mm (0.6"-2.4")
		Weight	600 g

Specifications

Portable microscope PSM-2, complete with 10 x eyepiece, 10 x objective, lamp housing and carrying case, without batteries	04286101
Objective, magnification 20 x	04286901
Objective, magnification 40 x	04286902
Measuring set consisting of measuring eye piece and plate micrometer	04286903

Transcopy

Transcopy Kit

Set consisting of 40 ml Transcopy Liquid, 50 replica foils 20 x 30 mm (0.8" x 1.2"), 1 pipette, spray nozzle and 50 microscope slides 25 x 75 mm (1" x 3")	40900090
Transcopy Replica Foils, 50 pcs.	40900091
Transcopy Liquid, 40 ml	40900092

RepliSet

Technical Data

Resolution of cu red replica	Down to 0.1 micron
Shrinkage	Negligible
Tear Strength	15-20 kN/m ²
Temperature range for the surface to be examined	-10°C to +180°C
Life span of the finished replicas is practically indefinite provided they are stored according to the instructions.	
Content in static-mixing nozzle	1.1 ml in nozzle for 50 ml cartridge 9.3 ml in nozzle for 265 ml cartridge

Specifications

Replication system for non-destructive testing of a microstructure or a 3D structure. Fast curing two-part silicon rubber compound for flexible high-resolution 3D replicas. For the 50 ml system, the hand-operated dispensing gun (40900066) and the static mixing nozzles (40900088) are used in combination with the 50 ml cartridges. For the 265 ml system, the hand-operated dispensing gun (40900065) and the static-mixing nozzles (40900056) are used in combination with the 265 ml cartridges.