Instructions for Use
Gebrauchsanweisung
Mode d'emploi

Planopol-2

App. No.	
V/Hz	

Please refer to the above when inquiring about the apparatus and when ordering spare parts.

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1. DESCRIPTION
2. TECHNICAL DATA
3. INSTALLATION
 - 3.1 Contents of packing case
 - 3.1.1 Bench model
 - 3.1.2 Built-in model
 - 3.2 Mounting
 - 3.2.1 Bench model
 - 3.2.2 Built-in model
 - 3.3 Connection to power
 - 3.4 Connection to water
4. CONVERSION KIT
 - 4.1 Contents of packing case
 - 4.2 Mounting of Conversion Kit
5. OPERATION
6. ACCESSORIES
7. MAINTENANCE
8. TROUBLE-SHOOTING
9. SPARE PARTS LIST

Struers Planopol-2

All the information and dates contained in this Manual were valid when it was printed. The policy of Struers A/S being, however, to make current improvements of our products, we reserve ourselves the right to modify our products without notice. Besides, we ask you to note that this Manual concerns apparatus and all accessories. Consequently, it may happen that the Manual mentions equipment which is not of relevance for your purpose at the present time.

3. INSTALLATION

3.1 Contents of packing case

3.1.1 Bench model

- 1 Planopol-2
- 1 protection ring
- 1 drain elbow pipe
- 2 hose clips 20-32 mm
- 1 hose connection
- 1 union nut for hose connection
- 1 sealing ring for hose connection
- 1 hose clip 13 mm
- 1 hexagon spanner 6 mm for turntable
- 1 fuse 0.5 AT 220 V

3.1.2 Built-in model

- 1 Planopol-2
- 1 protection ring
- 1 spill pan with drain pipe
- 1 sealing ring 25 mm dia for speed selector
- 1 packing, rubber, for sealing ring
- 1 pipe key NV13 and miscellaneous screws and nuts for assembly
- 1 hose connection
- 1 union nut for hose connection
- 1 packing for hose connection
- 1 hose clip 13 mm
- 1 hexagon spanner 6 mm for turntable

3.2 Mounting

3.2.1 Bench model

Planopol-2 is detached from the bottom of the packing case by removing the four screws from below, after which Planopol-2 can be placed on the bench.

3.2.2 Built-in model

The built-in model may be mounted in a working bench of suitable height (about 850-900 mm) with a width under the bench of min. 678 mm (see fig. 5) and a solid bench top 20-30 mm thick.

A number of machines may be placed so close to each other that their spill pans meet. But for convenience of operation a somewhat wider spacing is recommended.

Make a cut-out in the bench top and drill 4 holes according to fig. 5, observing the minimum distances to the front and rear edges of the bench.

Mount the bottom frame with motor etc. under the bench by means of the four fittings (fig. 12.51). Pass the screws through the four holes in the bench top. Adjust the positioning in the fittings in such a way that a polishing disc placed on the turntable will be parallel to the bench top at a height of 2 mm above it.

Dismantle the turntable (fig. 12.52) by removing the 3 screws. Mount the speed selector according to fig. 6 and fit the water cock. Stick the spill pan to the bench top with contact cement, carefully centering its holes with the components of the machine. Finally, mount the turntable and the bottom guard.

3.3 Connection to power

Mount a plug on the electric cable, connecting the leads as follows:

3-phase machine:

Yellow/green lead to earth

Brown lead to phase

Black lead to phase

Black lead to phase

The turntable should rotate anti-clockwise. If this is not the case after the connection, change the direction by changing over two of the phase leads.

Single-phase machine:

Yellow/green lead to earth

Brown lead to phase

Blue lead to neutral

3.4 Connection to water

Connect the thin reinforced hose to a water outlet and fit the 13 mm hose clip. The accompanying hose nipple with packing and union nut may, where not already provided, be mounted on the water outlet (1/2" pipe thread, BS 2779: G 1/2).

If a recirculating cooling unit is used, connect the hose to the discharge branch of the pump and fit a hose clip.

Lead the thick hose to the drain and be very careful to place it with a steady downward slope to prevent the water from being discharged too slowly which may cause overflow or blocking of the water outlet. In order to avoid elbows on the tubing, the enclosed elbow pipe can be mounted on the tubing in an appropriate place.

1. DESCRIPTION

Planopol-2 is a standard machine for lapping, wet grinding and aluminium and diamond polishing of metallographic, ceramic and mineralogical specimens.

It has a circular turntable accepting grinding, lapping and/or polishing discs of diameters up to 300 mm.

Supported by precision ball bearings, the turntable is driven via a V-belt drive by a powerful two-speed electric motor, which ensures smooth and noiseless running.

Motor and turntable are mounted on a steady steel frame also supporting the cabinet. The latter is provided with a water cock, a flexible sprinkler tube and a spill pan with drain.

On the control panel (bench model) are arranged, to the left, a main switch, and to the right, three push buttons with light indication for setting the speed (high/low) and for stop.

2. TECHNICAL DATA**Voltage/frequency**

3-phase machines:
 3 x 220 V/50 Hz
 3 x 380 V/50 Hz
 3 x 415 V/50 Hz
 3 x 500 V/50 Hz
 3 x 220 V/60 Hz
 3 x 440 V/60 Hz

1-phase machines:
 115 V/60 Hz
 220 V/50 Hz

Motor output

3-phase:
 50 Hz 550/300 W
 60 Hz 650/370 W

1-phase:
 50 Hz 370/220 W
 60 Hz 450/300 W

Disc speeds

Standard: 150/300 rpm
 With Conversion Kit: 75/150 rpm
 or: 300/600 rpm
 or: 600/1200 rpm

Max. diameter of polishing/grinding disc

300 mm

Dimensions

Bench model:
 Height: 340 mm
 Width: 462 mm
 Depth: 743 mm
 Weight: 53 kg

Built-in model:
 Height: 350 mm (under the bench)
 Width: 453 mm
 Depth: 644 mm (over the bench)
 Depth: 678 mm (under the bench)
 Weight: 43 kg

5. OPERATION

Planopol-2 is started by first operating the mains switch (fig. 1.1). This actuates the red lamp in the uppermost of the three push buttons at the right side (fig. 1.2). The button in the middle of the three (fig. 1.3) starts the principal motor at low speed, whereas the lowest button (fig. 1.4) switches to the higher speed. Activation of these buttons causes the lamp built into each button to glow. The machine is stopped by pressing the uppermost, red push button. The push buttons may be activated in any order. When stopped by the mains switch or by a power failure, the motor does not start on reconnection of power, but has to be started by operation of the push buttons. This prevents unintentional motor start during exchange of grinding/polishing disc.

In the case of prolonged overloading of the motor a built-in thermal relay cuts out the motor, re-start being possible only after a period of cooling.

The grinding, lapping or polishing disc is mounted on the turntable by inserting the three pins of the disc in the corresponding holes in the turntable. Make sure that the contact surfaces are clean. The water flow is controlled by a cock (fig. 1.5), and the distribution of water is controlled by means of the flexible sprinkler tube (fig. 1.6).

Below are stated the appropriate speeds for the various preparation processes:

- 75 rpm: Lapping on cast iron discs and PAD-K pellen cloth. Diamond polishing on Petrodisc and PA-W, PAN-W pellen cloth.
- 150 rpm: Diamond polishing (grinding/lapping on Petrodisc-M).
- 300 rpm: Wet grinding on grinding paper/wheel. Grinding on Petrodisc-M or DP-Net.
- 600 rpm: Alumina polishing (automatic preparation).
- 1200 rpm: Alumina polishing (manual preparation).

The grinding paper may be affixed to a plane polishing disc by means of Spray Adhesive, or it may be placed on a special disc according to the Knuth-Rotor principle. By the latter method the disc is filled with water while at rest, after which paper and retaining ring are put on. Then the disc is made to rotate t min. 250 rpm while the entire surface of the grinding paper is pressed gently against the disc with a specimen, or the like.

6. ACCESSORIES

A complete range of different discs of diameters between 200 and 300 mm is available, including discs for grinding with grinding paper, lapping on cast iron with a variety of surfaces, fine lapping on plastics or special lapping discs, as well as polishing discs of PVC and aluminium. In addition, we supply a complete programme of consumables, such as polishing cloths, diamond paste or spray, SiC paper or powder, alumina, etc. Please refer to special leaflets.

Planopol-2 has been prepared for mounting of PdM-Force, a motor-driven sample mover for simultaneous lapping and/or polishing of 1-9 specimens. For automatic preparation from 3 to 12 specimens simultaneously also Pedemax-2 may be mounted on Planopol-2 (3-phase model) thus making it possible to grind, lap and polish a number of specimens at a time. Both PdM-Force and Pedemax-2 reduce the preparation time and provide reproducible results. For details you are referred to special leaflets.

4. CONVERSION KIT

Planopol-2 is supplied as standard with 150/300 rpm. by mounting of a Conversion Kit (codeword: PLATT) the following rotary speeds may be obtained: 75/150, 300/600, and 600/1200 rpm.

4.1 Contents of packing case

- 1 V-belt pulley with two grooves (fig. 8.3)
- 1 fitted bolt (fig. 8.2)
- 1 washer (fig. 8.4)
- 1 adjusting plate (fig. 8.1)
- 1 V-belt Z 17 (10 x 425 mm) (60 Hz)
- 1 V-belt Z 18 (10 x 450 mm) (50 Hz)
- 1 V-belt HY-T SPZ (10 x 1312 mm) (60 Hz)
- 1 hexagon spanner (8 mm) for fitted bolt
- 1 hexagon spanner (2.5 mm) for V-belt pulleys
- 3 plastic plates stating the speeds 75, 600 and 1200 rpm

4.2 Mounting of Conversion Kit

The rotary speed is changed by changing the transmission ratio of motor and turntable. The V-belt pulley of the motor (fig. 3.5) is provided with two grooves of different diameters. In the standard version of 150/300 rpm the drive is effected directly from the small diameter groove to the V-belt pulley of the turntable. If the motor V-belt pulley is inverted so that the drive is established from the large diameter groove, the speed of the turntable is doubled. Like the motor V-belt pulley, the extra V-belt pulley (fig. 3.6) has two grooves of different diameters. By inserting this extra pulley as an intermediate wheel a speed either twice or half the one mentioned above may be obtained, depending on which way the intermediate wheel and the motor V-belt pulley are fitted.

The various combinations possible will appear from fig. 7, which also shows which V-belts to use. Always make sure that the V-belt pulleys are mounted in such a way that the belt pulls between two grooves running in the same plane.

Change of rotary speed is done as follows (after disconnection of power plug): Take out the 6 screws securing the rear section of the cabinet and lift this section off. Remove the nuts securing the

motor bridge, after which the V-belt pulley may be mounted according to fig. 8. Fit the V-belt pulley as shown in fig. 7, which indicates the correct positioning of the motor V-belt pulley. The motor V-belt pulley is fastened by means of two pointed screws in the bottom of the V-belt grooves (fig. 11.7 and fig. 12.27). When these screws are slackened, the V-belt pulley can be pulled off.

for exchange of the long V-belt the turntable must be taken out. This is done after removal of three screws at the top of the turntable.

Position the belts as shown in fig. 7. First tighten the short V-belt by means of the adjusting plate (fig. 8.1), then clamp the latter. Remount the motor bridge and tighten the long belt by pushing the motor rearwards and clamp it into position.

It is important to ensure adequate belt tension since undertensioning will cause the belt to slip on the pulleys, impairing the drive transmitted to the turntable, whereas overtensioning produces noise.

the long belt is correctly tensioned when a force of 5 N applied halfway between the pulleys causes a deflection of about 6 mm. When applied to the short belt, the same force should produce a deflection of 2-3 mm.

the push button plates indicating the new speeds are exchange by first unscrewing the metal ring round the buttons and then pressing off the green plastic plate. Fit the new indications and assemble the buttons in reverse order. Finally, remount the rear section of the cabinet.

7. MAINTENANCE

Planopol-2 is designed and constructed for many years of troublefree running. It will not need much attention, as its bearings are sealed and require no further lubrication. So maintenance will generally be limited to cleaning the exterior of the machine, the spill pan and the drain hose. Where the machine is used, for instance, for lumina polishing of SiC lapping, the spill pan should be cleaned sufficiently often to prevent drying out of the suspension and thus any risk of blocking of the hose. In order to avoid erratic disc movements, the contact surfaces between disc and turntable

should be kept clean.

Now and again the V-belt may need tightening, the transmission being less efficient with inadequate belt tension. To tighten the belt, remove the rear part of the cabinet, slacken the nuts securing the motor plate and push plate and motor towards the rear until the belt tension is correct (see section 4.2).

Exchange of th belt is effected by additionally removing the turntable and loosening the nuts securing the motor plate so as to slacken the belt, after which the belt can be taken off the turntable pulley from above. Mounting of the new belt is done in reverse order.

8. TROUBLESHOOTING

Trouble	Cause	Remedy
Motor does not run. Red light in stop button.	Motor overload. Thermal relay has cut out power.	Wait 5 minutes.
Motor does not run.	Internal transformer fuse blown.	Disconnect power supply. Remove rear cabinet section and replace transformer fuse.
Motor does not run, but buzzes.	One phase lacking (3-phase machine)	Check all fuses of local wiring.
Motor will only run at one speed.		Call service man.
No light in push button, otherwise normal machine function.	Bulb burnt out.	Disconnect power supply. Exchange bulb by unscrewing ring round push button. Insert new bulb and turn left.
Poor drive.	Slipping V-belt(s).	Tighten, see chapt. 7.
Water not drained.	Hose pressed flat. Dirt trapped in hose. Presence of air stops discharge.	Adjust hose position. Clean hose. Discharge hose should have a steady downward slope.
Machine produces noise.	Grinding/polishing disc vibrates. Bearings of motor or turntable defective.	Clean contact surfaces between disc and turntable. Exchange bearings/motor.

In the case of other problems, or if abovementioned faults cannot be remedied by following the instructions, call a service man.

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9. SPARE PARTS LIST

Apparatus: Planopol-2

Fig. and Pos.No.	Description	Spare part No.
11.1	Rear part of cabinet	362 MP 1
11.2	Front part of cabinet	362 MP 39
11.3	Control panel	362 MP 40
11.4	Bottom plate	362 MP 4
11.5	Steel frame	362 MP 5
11.6	Turntable support bridge with screws and bushes	362 MP 6
11.7	Motor support bridge with screws and bushes	362 MP 7
11.8	Motorplate with distance pieces and screws	362 MP 8
11.9	Spill pan	362 MP 9
11.10	Protection ring	362 MP 10
11.11	Drainhose 1" x 2000 mm with hose clip	362 MP 11
11.12,1.6	Flexible sprinkler tube	362 MP 12
11.13,1.5	Water cock	362 MP 13
11.14	Admission hose with union	362 MP 14
11.15	Short hose 3/8" x 300 mm with hose clip	362 MP 15
11.16	Cover plug	362 MP 16
11.17	Rubber knob with screw	362 MP 17
11.18	Motor, state voltage and frequency	362 MP 18
11.19	V-belt pulley for motor for 50 Hz version	362 MP 19
11.21	Vent cover	362 MP 21
11.22	Turntable	362 MP 22
11.23	Ball bearing, 6205-2Z	362 MP 23
11.24	Axle journal	362 MP 24
11.25	Hub, 80/52	362 MP 25
11.26	Distance bush, 51 x 2,5	362 MP 26
11.27	V-belt pulley	362 MP 27
11.29	Main switch	362 MP 41
1.2/3/4	Push button without switch element	362 MP 29
11.30	Switch element for push button	362 MP 30
	Bulb 24V/0,08 A	362 MP 31
	Coloured plastic plate, state colour	362 MP 32
11.31	Fuse, 0,5 AT	362 MP 33
11.32	Transformer	362 MP 34
11.33	Contacto, state frequency	362 MP 35
11.34	Thermo relay, state voltage and frequency	362 MP 36
	Capacitor, state voltage and frequency	362 MP 42
	V-belt (10 x 1237)	362 MP 38
	V-belt pulley for motor 60 Hz version	362 MP 37

Built-in model

12.1	Sprinkler pipe	266 MP 1
12.2	Tap	266 MP 2
	O-ring for tap	266 MP 2
	Admission hose with unions	266 MP 4
	Tap complete with hose, sprinkler pipe and outlet hose	266 MP 5
	Outlet hose with hose clip	266 MP 6
12.7	Protection ring	266 MP 7
	Clamping device for specimen holder, complete	266 MP 8
12.9	Cover plug	266 MP 9
	Cable lead through	266 MP 10
12.19	Ball bearing 6205-2Z	311 MP 19
12.21	Spring clip	311 MP 59
12.22	Pressure ring	311 MP 58
12.23	Switch with knob	266 MP 23
12.24	Motor, state voltage and frequency	266 MP 24
	Capacitor, state voltage and frequency	266 MP 25
12.26	V-belt pulley for motor	266 MP 26
12.27	Key for motor V-belt pulley	266 MP 27
12.29	V-belt Z44	266 MP 29
12.42	Spill pan	266 MP 42
12.43,6	Nameplate with sealings	266 MP 43
	Mains lead, three wires	266 MP 45
	Mains lead, four wires	266 MP 46
	Bottom frame	266 MP 47
	Nut for axle journal	311 MP 20
12.51	Mounting fitting with screws	266 MP 51
12.52	Turntable	266 MP 52
	Bottom plate	266 MP 53
	Motor bridge	266 MP 54
	Two V-belts for 75/150 rpm (Z44, Z19)	266 MP 55
	Two V-belts for 600/1200 rpm (Z46, Z19)	266 MP 56
Conversion Kit for built-in model (codeword: PLAKI)		
	Intermediate pulley, complete with bronze bearing, tight-fitting bolt, washer, and adjusting plate	266 MP 30
	Intermediate pulley without bronze bearing	266 MP 31
	Bronze bearing and tight-fitting bolt for intermediate pulley	266 MP 32
	V-belt (Z19)	266 MP 33
	V-belt (Z46)	266 MP 34
Conversion Kit for Planopol-2 (codeword: PLATT)		
8.3/4	Intermediate pulley with bronze bearing and tight-fitting bolt	362 MP 43
	V-belt (Z17) (60 Hz)	362 MP 44
	V-belt (Z18) (50 Hz)	362 MP 45
	V-belt HY-T SPZ (60 Hz)	362 MP 46

9. SPARE PARTS LIST

Apparatus: Planopol-2

Fig.and Pos.No.	Description	Spare part No.
11.1	Rear part of cabinet	362 MP 1
11.2	Front part of cabinet	362 MP 39
11.3	Control panel	362 MP 40
11.4	Bottom plate	362 MP 4
11.5	Steel frame	362 MP 5
11.6	Turntable support bridge with screws and bushes	362 MP 6
11.7	Motor support bridge with screws and bushes	362 MP 7
11.8	Motorplate with distance pieces and screws	362 MP 8
11.9	Spill pan	362 MP 9
11.10	Protection ring	362 MP 10
11.11	Drainhose 1" x 2000 mm with hose clip	362 MP 11
11.12,1.6	Flexible sprinkler tube	362 MP 12
11.13,1.5	Water cock	362 MP 13
11.14	Admission hose with union	362 MP 14
11.15	Short hose 3/8" x 300 mm with hose clip	362 MP 15
11.16	Cover plug	362 MP 16
11.17	Rubber knob with screw	362 MP 17
11.18	Motor, state voltage and frequency	362 MP 18
11.19	V-belt pulley for motor for 50 Hz version	362 MP 19
11.21	Vent cover	362 MP 21
11.22	Turntable	362 MP 22
11.23	Ball bearing, 6205-2Z	362 MP 23
11.24	Axle journal	362 MP 24
11.25	Hub, 80/52	362 MP 25
11.26	Distance bush, 51 x 2,5	362 MP 26
11.27	V-belt pulley	362 MP 27
11.29	Main switch	362 MP 41
1.2/3/4	Push button without switch element	362 MP 29
11.30	Switch element for push button	362 MP 30
	Bulb 24V/0,08 A	362 MP 31
	Coloured plastic plate, state colour	362 MP 32
11.31	Fuse, 0,5 AT	362 MP 33
11.32	Transformer	362 MP 34
11.33	Contacto, state frequency	362 MP 35
11.34	Thermo relay, state voltage and frequency	362 MP 36
	Capacitor, state voltage and frequency	362 MP 42
	V-belt (10 x 1237)	362 MP 38
	V-belt pulley for motor 60 Hz version	362 MP 37

Built-in model


12.1	Sprinkler pipe	266 MP 1
12.2	Tap	266 MP 2
	O-ring for tap	266 MP 2
	Admission hose with unions	266 MP 4
	Tap complete with hose, sprinkler pipe and outlet hose	266 MP 5
	Outlet hose with hose clip	266 MP 6
12.7	Protection ring	266 MP 7
	Clamping device for specimen holder, complete	266 MP 8
12.9	Cover plug	266 MP 9
	Cable lead through	266 MP 10
12.19	Ball bearing 6205-2Z	311 MP 19
12.21	Spring clip	311 MP 59
12.22	Pressure ring	311 MP 58
12.23	Switch with knob	266 MP 23
12.24	Motor, state voltage and frequency	266 MP 24
	Capacitor, state voltage and frequency	266 MP 25
12.26	V-belt pulley for motor	266 MP 26
12.27	Key for motor V-belt pulley	266 MP 27
12.29	V-belt Z44	266 MP 29
12.42	Spill pan	266 MP 42
12.43,6	Nameplate with sealings	266 MP 43
	Mains lead, three wires	266 MP 45
	Mains lead, four wires	266 MP 46
	Bottom frame	266 MP 47
	Nut for axle journal	311 MP 20
12.51	Mounting fitting with screws	266 MP 51
12.52	Turntable	266 MP 52
	Bottom plate	266 MP 53
	Motor bridge	266 MP 54
	Two V-belts for 75/150 rpm (Z44, Z19)	266 MP 55
	Two V-belts for 600/1200 rpm (Z46, Z19)	266 MP 56

Conversion Kit for built-in model (codeword: PLAKI)

	Intermediate pulley, complete with bronze bearing, tight-fitting bolt, washer, and adjusting plate	266 MP 30
	Intermediate pulley without bronze bearing	266 MP 31
	Bronze bearing and tight-fitting bolt for intermediate pulley	266 MP 32
	V-belt (Z19)	266 MP 33
	V-belt (Z46)	266 MP 34

Conversion Kit for Planopol-2 (codeword: PLATT)

8.3/4	Intermediate pulley with bronze bearing and tight-fitting bolt	362 MP 43
	V-belt (Z17) (60 Hz)	362 MP 44
	V-belt (Z18) (50 Hz)	362 MP 45
	V-belt HY-T SPZ (60 Hz)	362 MP 46

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
Fig. 4


Fig. 7

Note: Length of V-belts refers to internal length Foto

STRUERS K/S

Navn Kaa / AP

Dato 11.10.84

Anvendes i type

Erstatter 362- M-100R

Planopol-2, seen from below
change of rotational speed

362-M-102 R



Fig. 8

Matr.:

STRUERS K/S

Anvendes i type:

Mounting of cone polley

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±

Målestok:


Navn:

Dato: 21 - 4 - 81

Erstatter:

362-M-101 R

PRAC ORIGINAL

CAPACITOR

- 220V, 50Hz 25 μ F
- 220V, 60Hz 25 μ F
- 110V, 50Hz 100 μ F (4x 25)
- 115V, 60Hz 100 μ F (4x 25)

12 x 0,75[□] Cable

SWITCH

3 x 0,75[□]

MAINS

- bk = black
- bl = blue
- bn = brown
- vt = violet
- rd = red
- wh = white
- or = orange
- tr = transparent
- ye = yellow
- pk = pink
- gn = green
- be = beige

COLOR CODE

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Ikke ang. tolerancer:

Fig. 9

Matr.:

Foto

Målestok:

STRUERS K/S

Navn: HC/ *AP*


Dato: 16-3-82

Anvendes i type:

Erstatter:

Circuit Diagram
Single - Phase - Motor

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