Planopol-V

Instruction Manual
Gebrauchsanweisung
Mode d'emploi

Manual No.: 14017001

Date of Release 22.09.1994

Planopol-V
Instruction Manual

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Safety Precaution Sheet

To be read carefully before use

1. The operator should be fully aware of the use of the machine according to the Instruction Manual.
2. The machine must be placed in an adequate working position.
3. Be sure that the actual voltage corresponds to the voltage stated on the back of the machine. The machine must be earthed.
4. Make sure that the water connections are without leaks.
5. Keep clear of the rotational disc when operating. Do not drop the specimen while hand grinding.
6. The main water supply should be turned off if you leave the machine unattended.
7. If you observe malfunctions or hear unusual noises - stop the machine and call technical service.

The machine is designed for use with consumables and accessories supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the machine, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

IMPORTANT

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the front cover of this manual or on an attached label below. If in doubt consult the rating plate of the machine itself.

We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:
Instruction Manuals

A Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals

A Struers Service Manual may only be used by a trained technician authorized by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

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Table of Contents

Page

1.	Installation	5
2.	Operation	7
3.	Consumables and Accessories	9
4.	Maintenance	10
5.	Trouble-Shooting	11
6.	Technical Data	12

1. Installation

1.1 Packing List

- 1 Planopol-V
- 1 DP-disc of pvc, 300 mm dia.
- 1 Splash guard
- 1 Hexagon spanner, 6 mm
- 1 Hexagon spanner, 2.5 mm
- 1 Outlet hose, Ø 25 x 1500 mm
- 1 Angular outlet hose
- 1 Straight outlet hose
- 2 Clips, 20-32 mm
- 1 Hose connection
- 1 Nut
- 1 Gasket
- 1 Clip, 13 mm

1.2 Connection to Mains

Before the apparatus is connected it must be checked that the voltage which is indicated on the label of the machine, is in accordance with the local voltage.

Planopol-V is equipped with a three-core cable, which is connected to mains as follows:

1-phase machine

Yellow/green: earth
Brown: phase
Blue lead: neutral

3-phase machine

Yellow/green: earth
Brown: phase
Black: phase
Black: phase

The plug must be connected to a socket with switch.

1.3 Connection to Water

The thin reinforced hose is connected to a water cock by means of a 13 mm clip, if necessary the hose connection complete with gasket and union should be used. The amount of water is regulated with the needle valve (fig. 1.8). The thick hose, Ø 25 mm, is led to outlet. It is important to place and possibly shorten the hose so that it slopes evenly towards the outlet, otherwise there will be a very slow discharge, which will eventually result in clogging and flooding from the discharge tray (fig. 1.9). Use the elbow pipe and the straight tube when connecting and placing the discharge hose.

2.2 Operation

Place the chosen grinding wheel, lapping or polishing disc on the turntable, so that the pins fit in the corresponding holes.

NB: before the placing, please check that the contact faces are clean and smooth.

The rate of revolution is selected on the dial knob (fig. 1.1), and the apparatus is started with the ON/OFF push button (fig. 1.2). Then the speed can be read on the tachometer (fig. 1.4) and, if necessary, it can be fine adjusted.

At start and stop it is unnecessary to operate the speed selector, if the speed remains unchanged.

The motor-load is read on the ammeter. Exceeding for a shorter time of max. amperage, indicated with a red mark, will not damage the motor.

NB: at constant overload of the motor the built-in thermal cut-out will be released. After a couple of minutes, the thermal cut-out can be switched on again by pressing the button (fig. 1.6).

The power supply to the apparatus must be switched off (Power I/O) while the thermal cut-out is reset. If not there is a risk that the internal fuse will melt.

If specimen holders, Pedemax-2, or PdM- Force are used, the ammeter can be used for constant check of the grinding, lapping or polishing process.

3. Consumables and Accessories

A large number of discs from 200 to 300 mm Ø are available for Planopol-V for a.o. grinding with self-adhesive grinding paper or mounted according to the Knuth-Rotor system, lapping on cast iron with various kinds of surfaces, fine lapping on special lapping discs, as well as polishing discs of both pvc and aluminium. Moreover, a complete program of consumables is available, e.g. polishing cloths, diamond paste/spray, SiC-paper/powder, aluminium oxide etc. See separate leaflets.

Planopol-V is prepared for mounting of PdM-Force (fig. 2) and Pedemax-2 (fig. 3), which are motor operated specimen movers for simultaneous lapping, grinding, or polishing of 1-12 specimens. Thus, the process time can be reduced, and it is possible to work reproducibly. PdM-Force and Pedemax-2 are described in separate leaflets.

Furthermore, the reproducibility can be improved when using Multidoser for automatic dosing of lubricant and diamond suspensions.

4. Maintenance

Planopol-V is designed to work for many years without major service, as all the bearings are sealed, lubricated precision bearings, which do not require further lubrication. Therefore, the normal maintenance will be limited to cleaning the machine on the outside, as well as in the spillpan and the outlet tube. If working with e.g. aluminium oxide polishing or SiC lapping the spillpan should be cleaned so often that the suspension is not allowed to dry, which would clog the outlet hose. In order to avoid play of the disc the contact faces between the disc and the turntable should be kept clean.

At intervals it may be necessary to tighten the V-belt, as the torque will be reduced if the belt is not tight enough. This is done by removing the rear part of the cabinet, loosening the nuts which hold the motor plate and displacing the same and the motor backwards until the belt is tight enough, i.e. that a pressure on the V-belt of 100N between the pulleys may displace it 10-20 mm.

Replacement of the belt is carried out by also removing the turntable and, having loosened the nuts which hold the motor plate, so that the belt is slack, it is possible from above to remove the belt from the V-belt pulley.

Inserting a new belt is carried out in the reverse order.

5. Trouble-Shooting

Error	Cause	Action
Motor does not run, lamp (fig. 1.3) lights	Motor overloaded, thermal cut-out has switched off	Let motor cool for approx. 5 min, thermal cut-out to be switched on again (fig. 1.6)
Reduced torque	V-belt slips	Adjust belt
Water does not run out	Outlet hose squeezed	Straighten the hose
	Outlet hose clogged	Clean the hose
	Outlet hose does not slope downwards	Adjust the hose to an even slope
Machine makes noise	Grinding wheel, lapping or polishing disc vibrates	Clean the contact faces between disc and turntable

If there are other faults, or if the above mentioned error(s) cannot be remedied by the user, please call a service technician.

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Pos.	Description	Spare Part No.
<i>Electric parts</i>		
5.30	PC board, assembled, 115 V, incl. 20 kOhm speed potentiometer	401 MP 30
5.31	PC board, assembled, 220 V, incl. 20 kOhm speed potentiometer	401 MP 31
5.32	Mounting plate for PC board	401 MP 32
5.33	Speed control PC board, assembled	401 MP 33
5.34	RFI suppression filter Z1	401 MP 34
5.35	Smoothing choke L3	401 MP 35
5.36	RFI suppression choke L1	401 MP 36
5.37	Terminal strip	401 MP 37
5.38	Brake resistor (115 V)	401 MP 38
5.39	Brake resistor (220 V)	401 MP 39
7.40	RPM Instrument P1	401 MP 40
7.41	Ammeter P2 (115 V)	401 MP 41
7.42	Ammeter P2 (220 V)	401 MP 42
7.43	Thermal fuse FR1, 8A (115 V)	401 MP 43
7.44	Thermal fuse FR1, 4A (220 V)	401 MP 44
7.45	Pilot lamp H1 (115 V)	401 MP 45
7.46	Pilot lamp H1 (220 V)	401 MP 46
7.47	Potentiometer RP1, 20 kOhm	401 MP 47
7.48	Mounting plate for instruments	401 MP 48
7.49	Motor, 115 V 50-60 Hz	401 MP 49
5.50	Motor, 220 V 50-60 Hz	401 MP 50
7.61	On/Off switch element S1, complete	401 MP 61
7.62	Soft start printed circuit board	401 MP 62
<i>Water in- and outlet</i>		
4.51	Reinforced tubing 5/16 inch x 2000 mm, with 1/4" adapter	401 MP 51
4.52	Needle valve	401 MP 52
4.53	Reinforced tubing 3/8" x 300 mm	401 MP 53
4.54	Adjustable metal hose	401 MP 54
4.55	Angular hose connection	401 MP 55
4.56	Reduction socket 1/4"-1/8" RG	401 MP 56
5.57	Outlet tubing Ø 25 x 350 mm, transparent PVC	401 MP 57
5.58	PVC tube Ø 25 x 90 mm, with Socket MIFV 25 mm x 3/4" and adapter NIFV 25 mm x 3/4"	401 MP 58

2. Drawings and Diagrams

The drawings are not to scale. Some of the drawings may contain position numbers not used in connection with this manual.

Fig. 1


Fig. 2



Fig. 3


Fig. 4

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

2 SIX DIGIT ORDER #25-11



RS1 10K0791 M/S 15
 C18 0u6/x
 C19 0u6/x
 C20 0u6/x
 C21 100K
 C22 1u7
 C23 1u7
 C24 1u7
 R23 10K
 R24 10K
 R25 10K
 R26 10K
 R27 10K
 R28 10K
 R29 10K
 R30 10K
 R31 10K
 R32 10K
 R33 10K
 R34 10K
 R35 10K
 R36 10K
 R37 10K
 R38 10K
 R39 10K
 R40 10K
 R41 10K
 R42 10K
 R43 10K
 R44 10K
 R45 10K
 R46 10K
 R47 10K
 R48 10K
 R49 10K
 R50 10K
 R51 10K
 R52 10K
 R53 10K
 R54 10K
 R55 10K
 R56 10K
 R57 10K
 R58 10K
 R59 10K
 R60 10K
 R61 10K
 R62 10K
 R63 10K
 R64 10K
 R65 10K
 R66 10K
 R67 10K
 R68 10K
 R69 10K
 R70 10K
 R71 10K
 R72 10K
 R73 10K
 R74 10K
 R75 10K
 R76 10K
 R77 10K
 R78 10K
 R79 10K
 R80 10K
 R81 10K
 R82 10K
 R83 10K
 R84 10K
 R85 10K
 R86 10K
 R87 10K
 R88 10K
 R89 10K
 R90 10K
 R91 10K
 R92 10K
 R93 10K
 R94 10K
 R95 10K
 R96 10K
 R97 10K
 R98 10K
 R99 10K
 R100 10K

LAYER 5 3105
 800888 8105

X
 0
 Y


Struers Planopol-V · 401-M-41A · 02.89
 Component illustration
 Motor speed regulation

See adjustment procedure.

All measurements are made at 220V/50Hz mains voltage
 All voltages in the power circuit (the upper half of the diagram) are measured with respect to V.ref. (term. 16-18) as zero. Use TP1 "0V" for Osc. etc.
 All voltages in the control circuit (the lower half of the diagram) are measured with respect to V.ref. (term. 16-18) as zero. Use TP1 "0V" for Osc. etc.

Struers Planopol-V • 401-M-66 • 02.89
Wiring diagram 115/220 V, 50-60 Hz

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All rights reserved. © Struers 1994.