Prepomatic-2

Instruction Manual

Manual No.: 15347001

Date of Release 30.07.2009

*Prepomatic-2
Instruction Manual*

| Table of Contents | Page |
|--------------------------|-------------|
| User's Guide | 1 |
| Reference Guide..... | 21 |

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2009.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Prepamatic-2 Safety Precaution Sheet

To be read carefully before use

- 1.** The operator should be fully aware of the use of the machine according to the Instruction Manual.
- 2.** The machine must be installed in compliance with local safety regulations.
- 3.** The actual voltage must correspond to the voltage stated on the back of the machine. The machine must be earthed.
- 4.** Water connections must be leak-proof. Struers recommend that the mains water supply is shut off or disconnected if the machine is to be left unattended.
- 5.** Ensure that the emergency stop and other safety devices are in correct working order.
- 6.** In the event of malfunction or unusual noise, the machine should be stopped immediately and technical service called.
- 7.** Prepamatic-2 should be connected to an exhaust to remove the fumes from lubricants containing alcohol and from the alcohol used during cleaning.

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

User's Guide

| Table of Contents | Page |
|---|------|
| 1. Getting Started | |
| Checking the Contents of Packing | 2 |
| Placing Prepamatic-2 | 2 |
| Removal of Transport Fittings | 2 |
| Connecting Prepamatic-2 | 3 |
| Connection to Water Outlet..... | 3 |
| Connection to Water Supply | 3 |
| Connection to Compressed Air | 3 |
| Connection to Power Supply..... | 4 |
| Direction of Rotation | 4 |
| Filling with Consumables | 5 |
| Cooling Water | 5 |
| Alcohol and Soap..... | 6 |
| Suspension and Lubricant | 6 |
| Getting Acquainted with Prepamatic-2 | 7 |
| Software Settings | 8 |
| Configuration Menu..... | 8 |
| Setting the Language..... | 9 |
| 2. Basic Operations | |
| Using the Controls..... | 10 |
| Groups of Keys | 10 |
| Acoustic Signals..... | 10 |
| Location of Main Switch..... | 10 |
| Front Panel Controls | 11 |
| Reading the Display | 12 |
| Sleep Mode | 12 |
| Changing/Editing Values | 13 |
| Numeric Values..... | 13 |
| Alphanumeric Values | 14 |
| Preparing specimens | 15 |
| Selecting a Preparation Method..... | 15 |
| Selecting between Specimen Mover Plates and Holders | 17 |
| Using the Specimen Mover Disc for Single Specimens | 17 |
| Using the Specimen Holder | 18 |
| Starting the Process..... | 18 |
| Changing the Grinding Stone | 19 |

1. Getting Started

Checking the Contents of Packing

The crate in which Prepamatic-2 is delivered contains the following items:

- 1 Prepamatic-2
 - with 5 grinding/polishing stations, or
 - with 1 grinding station and 4 grinding/polishing stations
- 1 box with:
- 5 Splash rings
 - (4 pcs for Prepamatic-2 with grinding station)
- 1 Splash ring for OP-S
- 1 Water inlet hose 3/4" x 2 m
- 2 Filter gaskets
- 1 Reduction ring with gasket
- 1 Water outlet hose 32 mm x 3 m
- 2 Hose clamps 25-40 mm
- 1 Pressure hose 1/4"
- 1 Hose nipple for pressure hose
- 2 Hose clamps 12 mm
- 1 Soap,
- 4 Rubber blocks
- 1 Allen key 5 mm,
 - (only Prepamatic-2 with grinding station)
- 3 Specimens, prepared
- 1 Set of Instruction Manuals

Placing Prepamatic-2

Prepamatic-2 should be placed near water supply, drain, 3 -phase current and compressed air.

Position the machine leaving a space of at least 0.5m on the right, for servicing electronics, and 0.3m on the left, for access to the bottle unit.


Do not bolt Prepamatic-2 to the floor. Support each corner with the enclosed rubber blocks until the machine is level.

Removal of Transport Fittings

After installing Prepamatic-2 in its final position, the support block which is placed under the coupling for the specimen holder head can be removed and Prepamatic-2 is ready for connection.

Connecting Prepamatic-2

All the connections for Prepamatic-2 are located on the back of the machine.

- ① Water outlet
- ② Water inlet
- ③ Compressed air
- ④ Electrical connection box

Connection to Water Outlet

Connect the large outlet tube to the pipe marked "Water Outlet" on the back of Prepamatic-2. Lead the tube to the drain. Make sure that the tube slopes evenly downwards to prevent dirt from settling and blocking the tube.

Connection to Water Supply

Connect the water inlet hose to the threaded pipe marked "Water Inlet" on the back of the machine. Connect the other end of the tube to the water tap.

IMPORTANT

Do not forget to insert the two filters in the connecting tube. If the water installation is new, flush thoroughly before connecting Prepamatic-2.
Water pressure: min. 1 bar, max. 10 bar

Connection to Compressed Air

- Mount the hose nipple on the pressure hose and secure with the 12 mm hose clamp.
- Connect the hose nipple to the inlet marked "Compressed Air".
- Connect the other end of the pressure hose to the compressed air supply. Air pressure: min. 5.0 bar, max. 9.9 bar.

Connection to Power Supply

IMPORTANT

Check that the mains voltage corresponds to the voltage stated on the type plate on the back of the machine.

Open the electric connection box and connect a 4-lead cable* in the following way:

- PE: earth
- L1: phase
- L2: phase
- L3: phase

*Please see the section on *Technical Data* for recommended cable specifications.

Direction of Rotation

The direction of rotation for the grinding and polishing discs must be counter-clockwise. This is checked as following:

- Check that there is nothing to prevent the carrousel from turning, e.g. transport protection.
- Make sure that no specimen holder or specimen mover plate is mounted.
- Switch on the main switch (on the right hand side). If the safety guards are closed, the carrousel will turn, and the reference position of the carrousel can be found.
- Select a preparation program in the STRUERS METHODS.
- Press "Start" and check the direction of rotation of the preparation discs.
- Press "Stop" after the check is completed.

The discs should rotate counter-clockwise, see also the indication arrows on the supporting discs.

If the direction of rotation is clockwise, switch off Prepamatic-2 and unplug the machine.

- Change two of the phases.
- Repeat the rotation check.

Filling with Consumables

Cooling Water

Before Prepamatic-2 can be used, the different consumables required for the preparation have to be loaded.

Water from the Recirculation Cooling Unit is used for cooling and lubrication when grinding on the first position.

To fill the Recirculation Cooling Unit:

- Open the door at the front of Prepamatic-2.
- Pull out the “drawer” with the Recirculation Cooling Unit.
- Place a disposable plastic insert (see Accessories, section 5) in the tank and fold it over the rim.
- Fill the tank with 29.1 water and 900 ml Struers Additive. The water level should be about 7 cm below the upper rim of the tank.

IMPORTANT

Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to top up with Struers Additive each time you refill with water.

*Prepamatic-2
Instruction Manual*

Alcohol and Soap

Both soap and alcohol are used for cleaning the specimens. The two containers for these liquids are located in front of the Recirculation Cooling Unit and accessible through the door in front of Prepamatic-2. The container with the blue mark on the lid is to be filled with alcohol of min. 93% and connected to the blue tube. The container with the black mark on the lid is to be filled with a mixture of 100 ml concentrated soap solution, delivered together with the machine, and 4.9 l of water. Connect to the black tube.

IMPORTANT

Connect the correct tube to the proper container.

Suspension and Lubricant


Seven bottles for suspensions and lubricants are located in the bottle unit on the left hand side of Prepamatic-2. Bottles 1 - 6 can be used for both diamond suspension and lubricant, bottle no. 7 is for oxide polishing suspensions only. We recommend using bottles 1 - 4 for DP-Suspensions, e.g. 9µm in bottle 1, 6µm in bottle 2, 3µm in bottle 3 and 1µm in bottle 4. Bottles 5 and 6 can be used for DP-Lubricant Blue and Red. Bottle 7 should be filled with either OP-S or OP-U, depending on the materials to be prepared. Struers range of DP-Suspensions, lubricants or OP-Suspension are delivered in bottles that fit conveniently in the bottle unit.

IMPORTANT

After placing the different consumables in the bottle unit they have to be entered into the *Bottle Configuration* in the *Configuration Menu*. Please see the section *Configuration Menu* in the chapter *Advanced Operations*.

Getting Acquainted with Prepamatic-2

Take a moment to familiarise yourself with the location and names of the Prepamatic-2 components.


- ① Control unit
- ② Main switch (on the right hand side of the machine)
- ③ Bottle unit
- ④ Door to Recirculation Cooling Unit and cleaning solutions
- ⑤ Access to water and oil filter

Software Settings


Configuration Menu

Switch on the power at the main switch located at the right hand side of the machine. The following display will appear briefly:

The display will then change to the last screen shown before Prepamatic-2 was switched off, usually a preparation method. When switching on Prepamatic-2 for the first time, the MAIN MENU display will appear. If the heading in the display is different, press Esc, until the MAIN MENU appears (a long beep can be heard).


The MAIN MENU is the highest level in the menu structure. (See also the *Menu Overview*, Section 6). Pre-defined Struers methods, your own methods, manual functions or the configuration menu can then be selected.

↓
Press MENU DOWN ↓ to select CONFIGURATION

↓
Press ENTER ← to activate the CONFIGURATION menu

*Prepomatic-2
Instruction Manual*

Setting the Language


- Press MENU UP/DOWN to select Language.
- ↓
Press ENTER to activate the Language pop-up menu.
- ↓

LANGUAGE

 - English
 - Deutsch
 - Français
 - 日本語
- ↓
Press MENU UP/DOWN to select the language you prefer.
- ↓
Press ENTER to accept the language.
- The CONFIGURATION Menu now appears in the language you have chosen.
- Esc** Press Esc to return to the Main Menu.

2. Basic Operations

Using the Controls Front Panel

Groups of Keys

- A** Programming functions.
- B** Movement of the carrousel and dosing of suspension and lubricant
- C** START \diamond /STOP \ominus .

Acoustic Signals









- Bell *Short Beep*: when a key is pressed, a short beep indicates that the command has been accepted.
- Bell *Long Beep*: a long beep indicates that the key is currently inactive.

Location of Main Switch

The main switch is located at the right hand side of the machine.


*Prepomatic-2
Instruction Manual*

Front Panel Controls

| Name | Key | Function | Name | Key | Function |
|----------------|---|--|-----------------|--|---|
| MAIN SWITCH | | The main switch is located on the right side of the machine. | Esc |  | Leaves the present menu or aborts functions/changes. |
| EMERGENCY STOP |  | - Push the red button to stop. - Pull the red button to release. | ENTER |  | Accepts a selected parameter value or chooses a menu. |
| FUNCTION KEY |  | Controls for various purposes. See the bottom line of the individual screens. | Carrousel Left |  | Moves the carrousel left, to the next preparation disc. |
| FUNCTION KEY |  | Controls for various purposes. See the bottom line of the individual screens. | LUBRICANT |  | Manual dosing of lubricant |
| FUNCTION KEY |  | Controls for various purposes. See the bottom line of the individual screens. | Carrousel Right |  | Moves the carrousel right to the next preparation disc. |
| FUNCTION KEY |  | Controls for various purposes. See the bottom line of the individual screens. | ABRASIVE |  | Manual dosing of abrasive |
| MENU |  | Scrolls up in the menu tree structure. When setting a parameter the value increases. | START |  | Starts the preparation process according to the pre-set method. |
| MENU |  | Scrolls down in the menu tree structure. When setting a parameter the value decreases. | STOP |  | Stops the preparation process. |


Reading the Display

The display can show various kinds of information, for example parameters of a preparation method, or the ongoing process. An example of the screen for a preparation method is shown below:

- A Heading.
- B Inverted text: cursor position
- C Function key options.

During the process the screen could look as follows:

- A Name of the method used.
- B Preparation step currently running.
- C Remaining time in this step.

Sleep Mode

To increase the lifetime of the display, the backlight is automatically switched off if Prepamatic-2 has not been used for 30 min. Press any key to re-activate the backlight.

Please Note

The sample screens in this Instruction Manual show a number of possible texts. The actual display screen may differ from the samples in the Instruction Manual.


Changing/Editing Values

Numeric Values

Depending on the type of value, there are two different ways of editing.


STRUERS METHODS can not be edited or changed in any way, this is only valid for USER METHODS and CONFIGURATION.


 Press MENU UP/DOWN  to select the value to be changed, e.g. Force:


 Press ENTER  to edit the value.

 Two square brackets [] appear around the value.

 Press MENU UP/DOWN  to increase or decrease the numeric value.

 Press ENTER  to accept the new value. (Pressing Esc, aborts the changes, preserving the original value.)

*Prepomatic-2
Instruction Manual*


Alphanumeric Values

Press MENU UP/DOWN to select the alphanumeric value to be changed, e.g. Surface

↓
 Press ENTER ↲ to edit the value.

↓ A pop-up menu appears

↓
 Press MENU UP/DOWN to select the correct surface.

↓
 Press ENTER ↲ to accept the new preparation surface and to return to the previous screen. (Pressing Esc, aborts the changes, preserving the original setting.)

Preparing specimens

Selecting a Preparation Method

In the Main Menu, select STRUERS METHODS and press ENTER ↴ . The following screen will appear:

Depending on your material, select one of the 10 methods.
Consult your Metalog Guide for details on the selection of methods
for various materials.

The following table will give a rough guideline.

| Method | Materials to be prepared |
|---|---|
| Metalog Method A | Very soft, non ferrous materials, Al, Cu |
| Metalog Method B | Soft, ductile materials, Cu, Ti |
| Metalog Method C | Soft, less ductile materials, Cu-alloys, soft steels |
| Metalog Method D | Medium hard materials, steel, grey cast iron |
| Metalog Method E | Harder materials, tool steels, white cast iron |
| Metalog Method F | Very hard materials, ductile, sintered carbides |
| Metalog Method G | Very hard materials, brittle, ceramics |
| Metalog Method X | short method for soft materials, Al, Cu |
| Metalog Method Y | short method for harder materials, steel |
| Metalog Method Z | short method for hard materials, ceramics, sintered carbides |
| SiC-Paper (Prepamatic-2 with grinding station only) | As soft materials can not be ground on a grinding stone, this method is designed for the plane grinding of materials normally prepared according to methods A, B and C. Please select one step for the plane grinding of your material, continue with the appropriate method afterwards. |

*Prepomatic-2
Instruction Manual*

After determining the correct method, e.g. Metalog Guide Method D, move the cursor down to that method and press Enter to show the individual preparation steps.

| 4 Metalog Method D | | VIEW METHOD | | |
|--------------------|------------|-------------|-----------|---------|
| Pos. | Surface | Grit | Lubricant | Removal |
| | | Grain | | Time |
| 1. | MD-Piano | #220 | Water | 3m 0s |
| 2. | MD-Allegro | 9P | DP-Blue | 6m 0s |
| 3. | MD/DP-Dac | 6P | DP-Blue | 5m 0s |
| 4. | MD/DP-Dac | 3P | DP-Blue | 4m 0s |
| 5. | MD/OP-Chem | OP-U | None | 1m 30s |
| No of specimens | | 3 | | |

F1 **F2** Start at step **F3** Manu. Clean **F4**

Place the required grinding discs and polishing cloth on Prepomatic-2 and the machine is ready for the preparation of your specimens. If not all of the suspensions or lubricants necessary for this preparation method are in place, a message will appear on the bottom of the screen.

| 4 Metalog Method D | | VIEW METHOD | | |
|--------------------|------------|-------------|-----------|---------|
| Pos. | Surface | Grit | Lubricant | Removal |
| | | Grain | | Time |
| 1. | MD-Piano | #220 | Water | 3m 0s |
| 2. | MD-Allegro | 9P | DP-Blue | 6m 0s |
| 3. | MD/DP-Dac | 6P | DP-Blue | 5m 0s |
| 4. | MD/DP-Dac | 3P | DP-Blue | 4m 0s |
| 5. | MD/OP-Chem | OP-U | None | 1m 30s |
| No of specimens | | 3 | | |

Consumable missing: OP-U

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Insert the required consumables and change the BOTTLE CONFIGURATION in the CONFIGURATION MENU.

*Prepamatic-2
Instruction Manual*

*Selecting between Specimen
Mover Plates and Holders*

Prepamatic-2 can be used to prepare both single specimens and those clamped into specimen holders.

When using a specimen holder, 3 specimens need to be prepared at a time. As they are clamped in a holder these specimens do not have to have a regular shape, but should have similar surface areas.


For the preparation of 1 or 2 specimens, single specimen mode must be used. Here regular, round specimens are required.

*Using the Specimen Mover Disc
for Single Specimens*

Specimen mover discs are available for specimens in 25, 30, 40 and 50 mm dia. and in 1", 1¼" and 1½" dia.

The specimen mover discs for 25mm /1", 30mm /1¼" and 40 mm /1½" diameter specimens are identical. By using the correct plastic inserts (white for metric, black for imperial) they are set to the correct size.

First, the specimens have to be pre-levelled in the pre-levelling device, see Fig.

- Place the specimen in the opening with the correct diameter.
- Place a stainless steel ring with a corresponding diameter over the specimen.
- Use the Allen key to clamp the specimen securely. Make sure that both screws are tightened equally.
- Place the head with the specimen mover plate in the levelling device.
- Insert the specimens into the specimen mover plate.
- Move the pressure feet down onto the specimens.
- Secure the pressure feet by tightening the screws.
- Insert the head with the specimen mover plate in Prepamatic-2.

Using the Specimen Holder

When using a specimen holder, 3 specimens need to be prepared at a time. These specimens should have approximately the same size and hardness to avoid uneven grinding.

- Place the specimen holder in the levelling device and tighten the screw on the side to fix the holder in position.
- Insert the specimens into the openings and clamp them by tightening the screw.
- Insert the head with the specimen mover plate in Prepamatic-2.

IMPORTANT

Select screws of the correct length. Screw that are too long may stick out from the specimen holder and result in collisions with parts of Prepamatic-2 during operation.

Starting the Process

- Press START  and the process will be carried out automatically.

Changing the Grinding Stone

This chapter is only valid for Prepamatic-2 with a grinding stone.

For safety reasons the grinding stone is completely encapsulated. For access to the stone the cover of the grinding chamber has to be removed.

- Remove the 4 finger screws on top of the lid covering the stone and the cleaning chamber.
- Remove the cover.
- Swing the tube for the cooling water aside.
- Use the 5 mm Allen key to remove the screw from the centre of the grinding stone.
- Remove the flange.
- Remove the stone.
- Make sure that the rubber mat is clean and lying flat on the steel disc.
- Insert a new stone or a diamond grinding disc.
- Insert the flange.
- Insert the screw and tighten.
- Swing the cooling tube back into position.
- Remount the lid and secure it with the 4 finger screws.

If a new stone has been mounted see the section *Dresser Adjustment* on how to readjust the dresser.

*Prepomatic-2
Instruction Manual*

Reference Guide


| Table of Contents | Page |
|---|------|
| 1. Advanced Operations | |
| Configuration Menu | 23 |
| Setup Cleaning Programs | 24 |
| Setup Bottle Configuration | 24 |
| Setup User Consumables | 25 |
| Display contrast | 25 |
| Language | 25 |
| F1, Default value | 25 |
| Method Options | 26 |
| Copying a Method | 26 |
| Pasting a Method | 27 |
| Resetting a Method | 28 |
| Renaming a Method | 29 |
| Name Editing Principles | 30 |
| Saving a Method | 31 |
| Step Options | 32 |
| Copying a Step | 32 |
| Pasting a Step | 32 |
| Method Parameters | 33 |
| Surface | 33 |
| Grit/grain | 33 |
| Lubricant | 34 |
| Pre-dosing | 34 |
| Mode | 34 |
| Direction | 34 |
| Disc speed | 35 |
| Cleaning | 35 |
| Process Time; Stock Removal; Removal/Time | 35 |
| Force | 35 |
| Abrasive dosing | 35 |
| Lubricant dosing | 35 |
| Manual Functions | 36 |
| Cleaning of all tubes | 37 |
| Cleaning of OP-X tube | 37 |
| Manual dressing | 37 |
| Dresser Adjustment | 38 |
| Error sources when using stock removal | 41 |
| Calibration curves | 42 |

| | |
|----------------------------------|----|
| 2. Accessories | 43 |
| 3. Trouble-Shooting | 44 |
| 4. Maintenance | |
| Daily Service | 46 |
| Weekly Service | 46 |
| Checking the Cooling Water | 46 |
| 5. Technical Data | 47 |
| 6. Menu Overview..... | 49 |

1. Advanced Operations


Configuration Menu

Press Esc until the Main Menu appears (a long beep can be heard).

↓
Press MENU DOWN ↓ to select CONFIGURATION

↓
Press ENTER ← to activate the CONFIGURATION menu

*Prepomatic-2
Instruction Manual*

Setup Cleaning Programs


User defined cleaning programs can be edited here. After selecting SETUP CLEANING PROGRAMS and pressing ENTER ↲, chose one of the Empty Programs and press ENTER ↲. Select the step to be changed and press ENTER ↲.

From the pop-up menu select the required cleaning media then specify the length of time that this media should be applied. Up to 10 different steps can be combined in a Cleaning Program. The name of the Cleaning Program can be edited afterwards.

Setup Bottle Configuration

The contents of the bottle unit can be defined in this sub-menu.


Pressing ENTER ↲ will show a pop-up menu from which all available suspensions and lubricants can be selected. To add new consumables, please see the following section, SETUP USER CONSUMABLES.

The different suspensions and lubricants used must be defined here so that Prepomatic-2 can locate the correct suspension or lubricant.

*Prepomatic-2
Instruction Manual*


Setup User Consumables

Additional, user defined consumables can be added. When this option is selected, pressing ENTER ↴ will bring up a new screen:

Now you can add new preparation surfaces, suspensions and lubricants.

If a new preparation surface is added, the SURFACE STRUCTURE has to be defined. Select one of the user surfaces in USER SURFACES STRUCTURE and press ENTER ↴. There are two possibilities to chose between:

If the new preparation surface already contains abrasive, select the first option. If abrasive has to be added, chose the second.

Depending on the selection, different options will be available when editing the method parameters.

If ABRASIVE was selected, abrasive dosing will be set to OFF and with the lubricant dosing only water can be used and set to ON or OFF.

If NON-ABRASIVE was selected, the suspension to be used can be selected under Grit/Grain, a lubricant can be chosen under Lubricant, and the level for both Abrasive and Lubricant dosing can be adjusted.

Display contrast

The contrast settings of the display can be adjusted for maximum readability.

Language

The language can be set to either English, German, French or Japanese.

F1, Default value


The factory setting of some parameters can be restored by pressing the function key F1 when the appropriate value is highlighted.

Method Options


Methods can be adjusted to suit the particular requirements of the user. As the Struers methods can not be changed, the first step will usually be to copy a Struers method into the User Methods database.

Copying a Method

In the Main Menu select STRUERS METHODS.


↓
Press ENTER ↲ to activate the STRUERS METHODS menu.


↓
F1 Select the correct method and press F1:COPY.
A pop-up menu will appear.
↓

*Prepamatic-2
Instruction Manual*


Pasting a Method

- ↓ When asked for confirmation,
ENTER ↲ Press ENTER ↲ .
- ↓ The method is now resident in the buffer.
↓ To save the method in a different location:
Esc Press Esc to return to the Main Menu.
- ↓ In the Main Menu select USER METHODS.
ENTER ↲ Press ENTER ↲ to activate the USER METHODS menu.

- ↓ Select a method group for the new method.
ENTER ↲ Press ENTER ↲ to open this group.

*Prepomatic-2
Instruction Manual*


↓ Select where you want to paste the method from the buffer.

F2 Press F2:PASTE.

↓ When you are asked for confirmation,

← Press ENTER ↲

Resetting a Method

If a method is no longer used it can be reset to the default values.

↓ Select the method you want to reset.

F3 Press F3:RESET

↓ **←** Press ENTER ↲ to confirm resetting of the method.

The names of the methods in the USER METHODS database can be edited and changed. After copying a method from the Struers Methods the name can be changed to a name of your choice.


*Prepomatic-2
Instruction Manual*

Renaming a Method

Select the method you want to rename.

↓
F4 Press F4:RENAME.

The following screen will appear:

Use F1 or F3 to select the character to be changed.

↓
▼ Press MENU DOWN ▼ to move to the character set.
↓


Name Editing Principles

- Place the main cursor on the character you want to select, using F1:LEFT or F3:RIGHT. Use MENU UP/DOWN to move to the character set in the next line. An auxiliary cursor in the text line shows the position in the method name.
 - Write the new name using the following keys:
 - F1** Moves the cursor to the left
 - F2** Deletes one character in the text
 - F3** Moves the cursor to the right
 - F4** Inserts a space in the text
- ← ENTER places the new character in the method name and moves the auxiliary cursor in the name to the right. Repeat the procedure for each character.

Write the new name using the name editing procedures.

↓
Esc Press Esc to leave the editor.

↓
← Press ENTER ← to save the new name.

Saving a Method

While working with a user method, changes made in the database can be saved.

- As soon as a parameter is changed in one of the method steps, F4:SAVE will be shown on the bottom line of the display when you return to the method overview.

| 2.1 My own method | | EDIT METHOD | | |
|-------------------|------------|-------------|-----------|---------|
| Pos. | Surface | Grit | Lubricant | Removal |
| | | Grain | | Time |
| 1. | MD-Piano | #220 | Water | 3m 0s |
| 2. | MD-Allegro | 8P | DP-Blue | 6m 0s |
| 3. | MD/DP-Dac | 6P | DP-Blue | 5m 0s |
| 4. | MD/DP-Dac | 3P | DP-Blue | 4m 0s |
| 5. | MD/OP-Chem | OP-U | None | 1m 30s |
| No. of specimens | | 3 | | |

F1 **F2** Start at step **F3** Manu. Clean **F4** **ESC** **SAVE**

F4 Press F4: SAVE to save the changes after all necessary parameters have been altered.

ENTER ↲ Press ENTER ↲ to confirm saving the modified method.

If you started to generate a method from the very beginning, instead of copying one of the STRUERS METHODS, the name will automatically change from Empty method to Unnamed method. This will show that the method is no longer empty.

Important

When saving changes, the original method will be overwritten. To preserve the original method, a copy should be made and renamed. The new method can then be changed as required.

For details, see *Copying a Method*.

Step Options

Copying a Step

Single steps can also be copied to customise preparation methods.

Select the method and the preparation step to be copied.

↓
F3 Press F3 Copy Step, and the step will automatically be copied into memory.

Pasting a Step

Select the method and step where the copied step is to be inserted.


↓
F4 Press F4 Paste Step, and the step will automatically be inserted, overwriting the existing parameters.

*Prepamatic-2
Instruction Manual*


Method Parameters

For every method step a series of parameters can be changed:


Surface

Pressing ENTER ↲ allows you to select between all the surfaces in the Prepamatic-2 database:

Grit/grain


Pressing ENTER ↲ allows you to select between the different grain sizes available:

*Prepomatic-2
Instruction Manual*

Lubricant

Pressing ENTER ↲ allows you to select between the different lubricants available

It is only possible to select lubricant if the surface chosen is one with a non-abrasive surface. Otherwise water will be selected automatically.

Pre-dosing

Pressing ENTER ↲ allows you to adjust the pre-dosing between 0 and 30 seconds.

Mode

Pressing ENTER ↲ allows you to select between the 3 different modes.

Time

The preparation time can be specified in both a main and a final phase. Usually the main phase is used to remove a large amount of material, by using high forces and longer times, whereas the final phase is used to improve the scratch pattern by applying lower forces for a short while at the end of the preparation step.

Removal

When using a defined removal, only the main phase is active. The force is kept constant to ensure that the specified amount of material is actually removed.

Removal/Time

A combination of the above. First a specified amount of material is removed using a relatively high force. Thereafter the specimens are run for a short while at a low force to improve the scratch pattern.

Direction


The direction of rotation of the specimen mover can be set to either co- (>>) or counter (><) rotation. Co-rotation will give the best preparation result and the most uniform polishing. Counter-rotation can be used when polishing using oxide polishing suspensions to keep the suspension on the polishing cloth.

Disc speed

Two speeds can be selected for every position.
On Prepamatic-2 with the grinding stone the speed on pos. 1 can be set to either 500 or 1000 rpm.
On the next 3 positions and on the first four positions of the Prepamatic-2 with 5 grinding/polishing positions the speed can be set to either 150 or 300 rpm. On the last position of both machines the speed can be set to 75 or 150 rpm.

Cleaning

Pressing ENTER ↲ allows you to select between the different cleaning programs.

Apart from the 6 pre-defined cleaning programs, additional 4 programs can be user defined in CONFIGURATION.

*Process Time; Stock Removal;
Removal/Time*

Depending on the mode selected either the process time, the stock removal or a combination of both can be adjusted.
The maximum preparation time is 20 minutes in each phase.
The maximum stock removal is 1000 µm.

Force

Pressing ENTER ↲ allows you to adjust the force between 5 and 50 Newton. This value is applied per specimen, i.e. if the force is set to 50 N and 3 specimens are prepared, a total force of 150 N will be applied.
When a specimen holder is used the no. of specimens has always to be set to 3 to obtain the correct force.

Abrasive dosing

Pressing ENTER ↲ allows you to adjust the abrasive dosing level between 0 and 20.

IMPORTANT

If a grinding or polishing disc with a surface containing abrasive is selected, dosing is automatically set to OFF.

Lubricant dosing

Pressing ENTER ↲ allows you to adjust the lubricant dosing level between 0 and 20.

IMPORTANT


If a grinding or polishing disc with a surface containing abrasive is selected, dosing can only be set to OPEN or OFF.

*Prepomatic-2
Instruction Manual*

Manual Functions


Several manual functions are included in the software of Prepomatic-2.

Press Esc until you reach the Main Menu (a long beep can be heard).

↓
Press MENU UP/DOWN to select MANUAL FUNCTIONS.

↓
 Press ENTER to activate the MANUAL FUNCTIONS Menu.

*Prepamatic-2
Instruction Manual*

Cleaning of all tubes

The tubes leading suspension and lubricant to the dosing nozzle should be cleaned once a week to avoid blockage.

Select Cleaning of all tubes and press ENTER ↵. Then follow the instructions on the display.

For changing a single suspension or lubricant, remove all the tubes from the bottles when told, but only insert the tube you want to clean in the container with water. Then follow the rest of instructions.

Cleaning of OP-X tube

If an oxide polishing suspension has been used, the tube should be cleaned every day before leaving the machine overnight.

Select Cleaning of OP-X tube and follow the instructions on the display.

Manual dressing

This function is used for the dressing/opening of MD-Piano, diamond grinding discs, diamond pads and the grinding stone.

All Consumables except Grinding Stone

- Select Manual Dressing and press ENTER ↵.
- Move the disc to be dressed to the front and press ENTER ↵.
- The discs will start rotating with the low speed, water is applied and the protection screens can be opened.
- Use the dressing stick to dress the disc, and press ENTER ↵ when finished.

Dressing of Grinding Stone

(Prepamatic-2 with grinding station only.)

- Select Manual Dressing and press ENTER ↵.
- Move the grinding stone to the front and press ENTER ↵.
- Press ENTER ↵ to start the dressing of the grinding stone.
- Please wait while the stone is being dressed.

*Prepamatic-2
Instruction Manual*


Dresser Adjustment

This chapter is only valid for the Prepamatic-2 with a grinding stone. If a new grinding stone is placed on Prepamatic-2 or the grinding stone is reinserted after a diamond grinding disc has been used, the dresser has to be adjusted.

A screen like the following can appear:

A message will appear if the machine is started without adjusting the dresser:


Open the dresser door and follow the instructions on the screen:

Press ENTER ↲, and follow the instructions.

*Prepomatic-2
Instruction Manual*

If you are turning the handle in the wrong direction this screen will appear:


Turn the handle in the opposite direction until this screen appears:

Press ENTER ← to start the rotation of the grinding stone.

*Prepomatic-2
Instruction Manual*

Follow the instructions on the screen and press ENTER ↵ when the dresser just touches the grinding stone.

Close the dresser door and the adjustment of the dresser is completed.

*Error sources when using
stock removal*

Prepomatic cannot distinguish whether the stock removal takes place on the sample or on the grinding stone. Therefore, a smaller stock removal than the set one can take place. In order to compensate for this a higher removal can be set - how much higher only tests can show, as the removal is dependent on the sample material.

Therefore, where a very accurate stock removal is required, it is recommended to make a series of tests and then measure the actual removal. Where an accurate removal is required, it is not recommended to prepare more than one sample at a time.

Another source of inaccurate stock removal is the zero setting of the measuring system, as this does not take place until after one rotation of the sample holder head (13 rotations, when grit sizes finer than or equal to 2400 are specified). If the sample is very soft some removal will take place during the rotation. This inaccuracy can be compensated for if the above-mentioned measuring series is carried out.

In general it is recommended to use grinding forces and grinding bases/grain sizes which give relatively slow removals, i.e. more than 20 seconds when grinding on stone or more than 180 seconds on MD-Allegro.

When grinding on the grinding stone, the set stock removal may be higher than what the stone can grind off without dressing (opening). Therefore, an automatic dressing of the grinding stone takes place when the stock removal has been reduced to less than 15 µm per 10 sec. If many dressings are necessary a reduction in accuracy must be anticipated.

If a set stock removal e.g. on MD-Allegro is too big and therefore will take a very long time (and is maybe even impossible to obtain) the grinding stops automatically after 20 minutes.

Calibration curves

The error sources and their contribution to inaccuracy in the stock removal cannot be sorted out individually. Therefore, it is necessary to collect some statistical material. From this a calibration curve can then be made for the actual sample. When this curve is established the value to be set in the method to achieve a certain removal can then easily be found.

In order to get a reasonably accurate curve it is recommended to perform about 20 test measurements. Be sure that the measuring points are chosen in the requested work area.

The enclosed curve sheet can be copied and used to plot a calibration curve.

The following variables must be maintained during the measuring and the subsequent application of the curves:

- Sample diameter
- Sample material
- Sample test material area
- Grinding base
- Grinding force
- Mounting material
- Grain size/dosing
- Lubricant dosing

2. Accessories

| Specification | Cat. No |
|---|----------------|
| MD-Disc for magnetic fixation. | 02426920 |
| DP-Polishing disc. | 03756902 |
| Head for specimen holders | 05346902 |
| Specimen holder for 3 specimens, 24 - 40 mm dia. | 05346901 |
| Specimen holder for 6 specimens, 25 mm dia. | 05346908 |
| Specimen holder for 6 specimens, 30 mm dia. | 05346909 |
| Head for specimen mover plates | 05346903 |
| Set of specimen mover plates | 05346905 |
| Levelling device | 05346904 |
| <i>Disposable Inserts</i> for 30l Recirculation Cooling Unit for collecting waste material. Package with 20 pcs. | 49900001 |

3. Trouble-Shooting

| Error | Explanation | Action |
|--|---|---|
| Error Message | | |
| Dresser movement error. | Check air. | Check the compressed air supply and the free movement of the dresser. |
| Stone at minimum, new stone required. | | Replace the grinding stone. |
| Specimen holder moving while carrousel is moving. | | Call technical service |
| Specimen holder moving error. | The Specimen holder cannot move up or down. | Check that the specimen holder head moves freely upwards. If necessary, lubricate spindle through the lubricating nipple in the dresser room. If that does not help, call technical service. |
| Low water pressure. Check water supply. | | Check that the water inlet and all hoses are connected correctly and that the water pressure is sufficient. If the water supply has been interrupted, there might be air in the system. Re-start 2 - 3 times by pressing F1:Continue. |
| Timeout in removal mode. Reduce stock removal or increase force. | | Start the process with a higher force or reduce the total amount of material to be removed. |
| Removal error, force too low. Reduce stock removal or increase force. | | Start the process with a higher force or reduce the total amount of material to be removed. |
| Carrousel blocked. Remove a possible hindrance. | | Check that the carrousel moves freely and remove a possible hindrance. |
| Water is low in recirculation tank. | | Refill the Recirculation Cooling Unit. |
| Stone cover open. | | Replace the stone cover. |

*Prepamatic-2
Instruction Manual*

| Error | Explanation | Action |
|--|--|--|
| Error Message | | |
| Overload | Overload in low speed disc motor. Overload in high speed disc motor. | Wait until the motor has cooled down. When the display returns to the previous message, the machine can be re-started. |
| | Overload in water high pressure pump. | Wait until the pump has cooled down. When the display returns to the previous message, the machine can be re-started. |
| | Overload in recirculation pump. | Wait until the pump has cooled down. When the display returns to the previous message, the machine can be re-started. |
| | Overload in sample mover. | Wait until the motor has cooled down. When the display returns to the previous message, the machine can be re-started. |
| Acoustic Signal | | |
| Long beep. | The selected function is protected. It cannot be activated. | Check the setting. |
| Machine Problems | | |
| Prepamatic-2 is switched on but the display is blank | The backlight of the display has been switched off; Prepamatic-2 is in Sleep Mode. | Press any key to re-activate the backlight. |
| Disc rotates in the wrong direction. | Wrong connection to the power supply. | Switch two of the phases. |
| Motor buzzes. | One phase lacking. | Check all fuses of the local power supply. |
| Water not draining away. | Drain hose squeezed. | Straighten the hose. |
| | Drain hose clogged. | Clean the hose. |
| | Drain hose does not slope downward. | Adjust the hose to an even slope. |
| Water dripping underneath the machine. | The water does not drain away. | See above. |
| Cooling water stops. | Water tap on mains closed. | Open for water. |
| | Filter at the water inlet blocked. | Clean filter. |
| | Wrong software setting. | Check software setting. |

4. Maintenance

Daily Service

Clean the OP-X tube every day when you have finished working with Prepamatic-2.

- Go to MANUAL FUNCTIONS in the MAIN MENU.
 - Select CLEANING OP-X TUBE.
 - Follow the instructions on the screen.
-
- Clean Prepamatic-2 with a damp cloth.
-
- Clean all suspension and lubricant tubes.
 - Go to MANUAL FUNCTIONS in the MAIN MENU.
 - Select CLEANING OF ALL TUBES.
 - Follow the instructions on the screen.
-
- Clear the drainage ducts under the carrousel.
 - Remove the plastic drain covers on the carrousel.
 - Insert a brush or flat scraper into the hole and turn the carrousel using the keypad controls **◀▶**.
 - Carefully brush away any debris along the ducts and down the drain.
 - Finally, carefully flush the ducts with water.

WARNING!

Do not use alcohol, acetone or similar solvents
for cleaning Prepamatic-2.

Checking the Cooling Water

The cooling water should be checked after 8 hours use or at least every week and refilled when necessary.

Remember to add Struers Additive for Cooling Fluid: One part of Additive for 33 parts of water.

To check the concentration of additive, use a refractometer. Concentration = Brix value. The concentration of additive should always lie between 2.7 and 3.3 %. Add Struers Additive for Cooling Fluid if the concentration is too low.


5. Technical Data

| Subject | Specifications | | | | |
|---|--|--|--------------------------------|----------------|--------------------------------|
| Software and Electronics | Display: | 128 x 240 dots (16 x 40 characters) | | | |
| | Controls: | touch pad | | | |
| | Database: | 10 Struers methods 70 user definable methods | | | |
| Power Supply | 3 x 200 V / 50 Hz 3 x 200-210V / 60 Hz 3 x 220-230 V / 50 Hz 3 x 220-240 V / 60 Hz 3 x 380-415 V / 50 Hz 3 x 460-480 V / 60 Hz | max. cont. load 7.3 A max. cont. load 7.3 A max. cont. load 7.3 A max. cont. load 7.3 A max. cont. load 3.9 A max. cont. load 3.9 A | | | |
| Mains Cable Recommendation | | Min. Fuse size | Minimum cable size @ Min. fuse | Max. Fuse size | Minimum cable size @ Max. fuse |
| | 3 x 200 V / 50 Hz | 40 A | 3x4.0mm ² + PE | 50 A | 3x4.0mm ² + PE |
| | 3 x 200-210 V / 60 Hz | 40 A | 3xAWG8 + PE | 50 A | 3xAWG8 + PE |
| | 3 x 220-230 V / 50 Hz | 40 A | 3x2.5mm ² + PE | 50 A | 3x4mm ² + PE |
| | 3 x 220-240 V / 60 Hz | 40 A | 3xAWG10 + PE | 50 A | 3xAWG10 + PE |
| | 3 x 380-415 V / 50 Hz | 32 A | 3x1.5mm ² + PE | 50 A | 3x4mm ² + PE |
| | 3 x 380-415 V / 60 Hz | 30 A | 3xAWG14 + PE | 50 A | 3xAWG10 + PE |
| | 3 x 460-480 V / 60 Hz | 30 A | 3xAWG14 + PE | 50 A | 3xAWG10 + PE |
| | <p><i>Important:</i> Local standards may overrule the recommendations for the main supply cable. If necessary, please contact a qualified electrician to verify which option is suitable for the local installation setup.</p> | | | | |
| <i>Residual Current Circuit Breaker</i> | type A, 30 mA (or higher) is recommended. | | | | |
| Electrical Data | Main motor: | 750 / 120 W | | | |
| | Sample mover motor: | 30 W | | | |
| | Carrousel motor: | 20 W | | | |
| | High pressure water pump: | 180 W | | | |
| | Recirculation pump for cooling water: | 140 W | | | |
| Water Supply | min 1.0 bar, max. 10.0 bar | Consumption: approx 1.5 l/cleaning | | | |
| Air Supply | min. 5.0 bar, max. 9.9 bar | max. 60 l/min, average 15 l/min | | | |

*Prepomatic-2
Instruction Manual*

| Subject | | Specifications | |
|-----------------------------------|---|---|----------------------|
| Specimen-holder head | Revolutions: | 125 rpm | |
| | Force: | 5 - 150 N | |
| Stock removal | Smallest step: | 5 µm | |
| | 0-position error: | 5 - 20 µm | |
| | Grinding error: | ± 3 µm | |
| | Max. stock removal: | 1000 µm | |
| | Repetition accuracy: | ± 5 µm | |
| Griding polishing disc | Revolutions: | 75, 150, 300 | |
| | | 75, 150, 300, 500, 1000 (Prepomatic-2 with grinding station) | |
| Dresser | Removal: | 0,06 mm Prepomatic-2 with grinding station only | |
| Dosing system | DP-Lubricant: | DP-Suspension | OP-Suspension |
| | steps | 20 | 20 |
| | quantity | 0.2 – 12.25 ml/min | 0.2 - 4.08 ml/min |
| Dimensions and Weight | Width: | 1100 mm / 43" | |
| | Depth: | 900 mm / 35" | |
| | Height: | 1710 mm / 67" | |
| | Weight: | 430 kg / 900 lbs | |
| Recirculation cooling unit | Capacity: | 35 l | |
| Safety standards | Please refer to the Declaration of Conformity | | |

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
Manual No.: 15347001

Date of Release GE01 .200J

Prepomatic-2
Spare Parts and Diagrams

| Table of Contents | Page |
|--------------------------|-------------|
| Spare Parts | 1 |
| Diagrams..... | 4+ |

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 200J.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Telefax: +45 44 600 801

1. Spare Parts

| Table of contents | Page /Drawing |
|--|---------------|
| Prepamatic-2, 5 Polishing units (5xP) | 3 |
| <i>Drawings</i> | |
| Prepamatic-2 (5xP), complete | 15340004A |
| Carousel wheel (5xP), complete | 15340569B |
| Polishing disc (5xP), pos 1, complete | 15340635 |
| Top section (5xP)..... | 15340031C |
| Front plate (5xP), complete | 15340041E |
| Carousel top (5xP), complete..... | 15340429C |
| | |
| Prepamatic-2, 4 Polishing units and Grinding stone (4xP+1G) | 11 |
| <i>Drawings</i> | |
| Prepamatic-2 (4xP+1xG), complete..... | 15340005A |
| Carousel wheel, complete (4xP+1xG) | 15340568B |
| Polishing disc, pos 5 (4xP+1xG), complete | 15340639 |
| Grinding disc, pos 1, complete | 13750660 |
| Top section (4xP+1xG)..... | 15340030D |
| Diamond dresser, complete | 15340362A |
| Throttle valve, complete | 13750400C |
| Front plate, complete (4xP+1xG) | 15340040E |
| Carousel top (4xP+1xG), complete..... | 13750429C |
| Lid f. cleaning chamber (4xP+1xG)..... | 13750495A |

Some of the drawings may contain position numbers
not used in connection with this manual.

*Prepamatic-2
Spare Parts and Diagrams*

| Table of contents | Page /Drawing |
|--|----------------------|
| Prepamatic-2, common spare parts | 23 |
| <i>Drawings</i> | |
| Bottom cabinet, complete..... | 15340007H |
| Lift for spindle, complete | 15340010B |
| Main spindle, complete | 15340012F |
| Left hand side cupboard, complete (air & water section) | 15340017G |
| M1, high pressure pump, complete | 15340018B |
| Peristaltic pumps module | 15340022D |
| Pump, complete | 14600031E |
| Bottle section, complete | 15340025F |
| Cleaning nozzle block, complete..... | 13750460D |
| Carousel drain, complete | 13750529B |
| Recirculation cooling unit, complete..... | 13750721C |
| Coupling (female), complete | 14290506A |
| Coupling (male) for sample holder | 15340050D |
| Coupling (male) for single sample holder..... | 15340051A |
| Levelling device, complete | 15340052 |
| Rack with PCB's, complete | 15340736D |
| High power module, complete..... | 15340738C |

Some of the drawings may contain position numbers
not used in connection with this manual.


Prepamatic-2
Spare Parts and Diagrams


**Spare part list for Prepamatic-2,
5 Polishing units (5xP)**



| Drawing | Pos. | Spare Part | Cat no: |
|----------|------|--|----------|
| 15340004 | | Prepamatic-2 (5xP), complete | |
| | 230 | Sealing v-profile, 2101 (1 m) | 2IP30002 |
| 15340569 | | Carousel wheel (5xP), complete | |
| | 80 | Bearing, spherical 1210 | 375MP107 |
| | 150 | Timing belt 16 T10/810 | 2JT30810 |
| | 160 | M4, disc motor 3x200-240v/50Hz ø24 | 13751024 |
| | 160 | M4, disc motor 3x200-220v/60Hz ø24 | 13751025 |
| | 160 | M4, disc motor 3x200-240v/50Hz ø24 | 13751026 |
| | 160 | M4, disc motor 3x380-480v/50-60Hz ø24 | 15340780 |
| | 330 | Brake cylinder | 375MP179 |
| | 400 | Brake pad for cylinder | 13750598 |
| | 600 | Carousel motor, complete | 13750624 |
| | 700 | Intermediate station, complete | 375MP119 |
| | 710 | Polishing disc, pos 2+3+5 (5xP), complete | 375MP120 |
| | 720 | Polishing disc, pos 4, complete | 375MP122 |
| | 740 | Polishing disc, pos 1 (5xP), complete | 15340635 |
| | 750 | Polishing disc, pos 2+3+5 (5xP), complete | 375MP120 |
| | 746- | Set of round belts (5xP - 50 Hz) | 15342908 |
| | 770 | | |
| | 746- | Set of round belts (5xP - 60Hz) | 15342909 |
| | 770 | | |
| 15340635 | | Polishing disc (5xP), pos 1, complete | |
| | 30 | Ball bearing 6004-2RSR | 2BK00045 |
| | 80 | Belt pulley (small) | 375MP118 |
| | 140 | Turntable with screws | 300MP026 |
| 15340031 | | Top section (5xP) | |
| | 50 | Cap-screw M4x8 A2 (10 pcs) | RTR80408 |
| | 60 | Crinkle washer M4 (10 pcs) | RZI20403 |
| | 90 | Cap-screw M4x8 A2 (10 pcs) | RTR80408 |
| | 100 | Crinkle washer M4 (10 pcs) | RZI20403 |
| | 500 | Micro switch, complete | 13750948 |
| | 630 | Nozzle block, complete | 15340175 |

Prepematic-2
Spare Parts and Diagrams

| Drawing | Pos. | Spare part | Cat no: |
|----------------|-------------|--|----------------|
| 15340041 | | Front plate (5xP), complete | |
| | 20-30 | Touch pad keyboard, complete | R5340201 |
| | 60 | Display PCB w. con. & back light, pos 502 | 14930045 |
| | 65 | Display protection glass | 14930385 |
| | 70 | Prep-2 CPU-PCB, pos 501 | 15343260 |
| | 71 | Set of PROM's, latest version (V.2.x) | R5343901 |
| | 90 | Prep-2 pressure con. PCB, pos 503 | 15343280 |
| 15340429 | | Carousel top (5xP), complete | |
| | 20 | Drain tray, pos 1-5(5xP), 2-5 (4xP+1xG) | 375MP094 |
| | 70 | Lid for cleaning compartment | R3750441 |
| | 320 | Flexible water nozzle, complete | 422MP003 |
| | 340 | Gasket PVC 1/4" (10 pcs) | RIF00012 |
| | 370 | Lid f. cleaning chamber (5xP) + thumb screws | R3758100 |
| | 380 | Thumb screw M4 | 375MP096 |

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




**Spare part list for Prepomatic-2,
4 Polishing units and Grinding stone (4xP+1xG)**


| Drawing | Pos. | Spare Part | Cat no: |
|----------|------|--|----------|
| 15340005 | | Prepomatic-2 (4xP+1xG), complete | |
| | 310 | Cardboard blotters ø35/58 | 13752902 |
| | 320 | Grinding stone ø200mm #150 | 40800005 |
| | 330 | Flange for grinding stone | 375MP126 |
| 15340568 | | Carousel wheel, complete (4xP+1xG) | |
| | 80 | Bearing, spherical 1210 | 375MP107 |
| | 150 | Timing belt 16 T10/810 | 2JT30810 |
| | 330 | Brake cylinder | 375MP179 |
| | 400 | Brake pad for cylinder | 13750598 |
| | 600 | Carousel motor, complete | 13750624 |
| | 700 | Intermediate station, complete | 375MP119 |
| | 710 | Polishing disc, pos 2+3+5(5xP), complete | 375MP120 |
| | 720 | Polishing disc, pos 4, complete | 375MP122 |
| | 730 | Polishing disc, pos 5 (4xP+1xG), complete | 375MP123 |
| | 740 | Grinding disc, pos 1 (4xP+1xG), complete | 375MP124 |
| | 750 | Polishing disc, pos 2+3+5 (5xP), complete | 375MP120 |
| | 746- | Set of round belts (50 Hz) | 375MP097 |
| | 770 | Set of round belts (60 Hz) | 375MP101 |
| 15340639 | | Polishing disc, pos 5 (4xP+1xG), complete | |
| | 30 | Ball bearing 6004-2RSR | 2BK00045 |
| | 80 | Belt pulley (large) | 375MP117 |
| | 140 | Turntable with screws | 300MP026 |
| 13750660 | | Grinding disc, pos 1, complete | |
| | 30 | Ball bearing 6004-2RSR | 2BK00045 |
| | 80 | Belt pulley (large) | 375MP117 |
| | 140 | Turntable for grinding stone with screws | 375MP121 |
| 15340030 | | Top section (4xP+1xG) | |
| | 50 | Cap-screw M4x8 A2 (10 pcs) | RTR80408 |
| | 60 | Crinkle washer M4 (10 pcs) | RZI20403 |
| | 90 | Cap-screw M4x8 A2 (10 pcs) | RTR80408 |
| | 100 | Crinkle washer M4 (10 pcs) | RZI20403 |
| | 290 | Air hose ø5/ø3.2mm (2m) blue | 2NU12445 |
| | 500 | Micro switch, complete | 13750948 |
| | 630 | Nozzle block, complete | 15340175 |
| 15340362 | | Diamond dresser, complete | |
| | 10 | Cylinder | 375MP072 |
| | 150 | Friction element (3 pcs) | 382MP028 |
| | 300 | Slotted optical switch (PCB B8-B13) | 375MP239 |
| | 500 | PLAMT, diamond point for dressing tool | 40800045 |



Prepomatic-2
Spare Parts and Diagrams


| Drawing | Pos. | Spare part | Cat no: |
|----------|-------|--|----------|
| 13750400 | 80 | Throttle valve, complete Air hose PU-3 blue ø2.9/ø4.3 (2m) | 375MP172 |
| 15340040 | 20-30 | Front plate, complete (4xP+1xG) Touch pad keyboard, complete | R5340201 |
| | 60 | Display PCB w. con. & back light, pos 502 | 14930045 |
| | 65 | Display protection glass | 14930385 |
| | 70 | Prep-2 CPU-PCB, pos 501 | 15343260 |
| | 71 | Set of PROMs, latest version (v.2.x) | R5343901 |
| | 90 | Prep-2 pressure con. PCB, pos 503 | 15343280 |
| | 120- | Lock for door | 375MP013 |
| | 150 | | |
| | 160 | Proximity sensor | 2HQ00023 |
| 13750429 | 20 | Carousel top (4xP+1xG), complete Drain tray, pos 1-5(5xP), 2-5 (4xP+1xG) | 375MP094 |
| | 50 | Marker washer | 13750515 |
| | 70 | Lid for cleaning compartment | R3750441 |
| | 185 | Drain tray, pos 1 (4xP+1xG) | 375MP095 |
| | 280 | Air hose PU-6 ø5.9/ ø8.3 (5m) | RNU14061 |
| | 250 | Proximity sensor | 2HQ00023 |
| | 320 | Flexible water nozzle, complete | 422MP003 |
| | 340 | Gasket PVC 1/4" (10 pcs) | RIF00012 |
| | 370 | Lid f. cleaning chamber (4xP+1xG) | 375MP093 |
| | 380 | Thumb screw M4 | 375MP096 |
| 13750495 | 30 | Lid f. Cleaning chamber (4xP+1xG) Rubber curtain for stone cover | 13750500 |
| | 60 | Rubber curtain for dresser track | 375MP275 |
| | 110 | Drip catcher for diamond point | 375MP268 |

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








Prepomatic-2
Spare Parts and Diagrams

Spare part list for Prepomatic-2, common

| Drawing | Pos. | Spare Part | Cat no: |
|----------|-------|--|----------|
| 15340007 | | Bottom cabinet, complete | |
| | 130 | Black cap & tube for 5l container (detergent) | 375MP131 |
| | 135 | Cap w. blue tube for 5l container (alcohol) | 13750729 |
| | 140 | 5l container with labels | 375MP129 |
| | 460 | Double tubing (blue 6 black) for 5l containers | 375MP180 |
| 15340010 | | Lift for spindle, complete | |
| | 80 | Ball bearing 625-2Z | 2BK00005 |
| | 130 | Nylonit slide bearing | 13750307 |
| | 140 | Gear wheel with pin & Omnifit 200m | R3750322 |
| | 160 | Lift motor M6, 24VDC w. gear wheel | 375MP203 |
| 15340012 | | Main spindle, complete | |
| | 130 | Shell EPOO grease 20ml syringe | 14790651 |
| | 150 | Air hose ø5/ø3.2 (2m) | 2NU12445 |
| | 160- | Bellows w. housing | 375MP045 |
| | 170 | | |
| | 210 | Ball bearing 6205-2ZR | 2BK00065 |
| | 270 | Ball bearing 6004-2RSR | 2BK00045 |
| | 350 | Timing belt HTD 375x15 metric | 2JT20375 |
| | 360 | Spe. holder motor M3 3x200-240V /50-60Hz | 375MP201 |
| | 360 | Spe. holder motor M3 3x380-500V /50-60Hz | 375MP202 |
| | 490 | Proximity sensor | 2HQ00023 |
| | 550 | Ball bearing LR 200 NPPU | 2BK10010 |
| | 710 | Slotted optical PCB B8-B13 | 375MP239 |
| 15340017 | | Left hand side cupboard, complete (air & water section) | |
| | 20 | Y3-Y8 solenoid valve | 2YM10004 |
| | 40-50 | Throttle valve V3,V4 | 394MP059 |
| | 70 | Gasket PVC 1/8" (10 pcs) | RIF00011 |
| | 80 | O-ring ø10.82x1.78 NBR 872 (10 pcs) | 2IO17815 |
| | 90 | Gasket PVC 1/8" (10 pcs) | RIF00011 |
| | 130 | Gasket PVC 1/8" (10 pcs) | RIF00011 |
| | 160 | Filter reduction valve unit, SMC | 2YF00002 |
| | 170 | Fine filter unit, SMC | 2YF00003 |
| | 190 | Gasket PVC 1/4" (10 pcs) | RIF00012 |
| | 195 | Air hose ø8/ø5 (5 m) | RNU12508 |
| | 220 | Gasket PVC 1/4" (10 pcs) | RIF00012 |
| | 390 | Gasket PVC 1/4" (10 pcs) | RIF00012 |
| | 410 | Water hose 1/2" (10 m) | RNU29316 |
| | 440 | Quick angle swivel fitting ø8-1/4" | 2NF10087 |
| | 460 | Gasket PVC 1/4" (10 pcs) | RIF00012 |
| | 470 | Air hose ø8/ø5 (5 m) | RNU12508 |
| | 480 | Vacuum hose 5/16" (3m) | 2NU29310 |

Prepamatic-2
Spare Parts and Diagrams

| Drawing | Pos. | Spare part | Cat no: |
|----------|-------|--|-----------|
| 15340018 | | M1, high pressure pump, complete | |
| | 10-50 | Motor 3x200-440v/50-60Hz | 375MP210 |
| | 10-50 | Motor 3x415-480v/50-60Hz | 375MP211 |
| | 90 | High pressure pump PROCON 1309 | 2yP01309 |
| | 120 | Pressure switch 1-10 bar for Prep. + Maps | 2SP10027 |
| | 130 | Gasket, PVC 1/4" (10 pcs) | RIF00012 |
| | 160 | Gasket, PVC 1/4" (10 pcs) | RIF00012 |
| | 175 | Gasket, PVC 3/8" (10 pcs) | RIF00013 |
| | 200 | Gasket, PVC 3/8" (10 pcs) | RIF00013 |
| | 230 | Gasket, PVC 1/4" (10 pcs) | RIF00012 |
| | 250- | Danfoss water solenoid valve 50Hz, complete | 2YM3002 |
| | 270 | | |
| | 250- | Danfoss water solenoid valve 60Hz, complete | 394MP062B |
| | 270 | | |
| | 300 | Gasket, PVC 1/4" (10 pcs) | RIF00012 |
| 15340022 | | Peristaltic pumps module | |
| | 30 | Bosch motor, complete with shaft | R4600550 |
| | 140 | Solenoid valve 1/8"(l) | 2YM10140 |
| | 170 | Air hose, blue Pu-4 (2m) | 2NU14041 |
| | 190 | Air hose, blue Pu -4 (2m) | 375MP177 |
| | 200 | Silicone tube ø2/ø4 (10 m) | 14602901 |
| | 210 | Silicone tube ø4/ø6 (10 m) | 14602902 |
| | 220 | Silicone tube ø4/ø6 (10 m) | 14602902 |
| 14600031 | | Pump, complete | |
| | 10-60 | Pump roller, complete | 14600601 |
| | 90 | Flange bearing (7 pcs) | RBG20040 |
| | 100- | Lid for pump with bearing | 14600602 |
| | 110 | | |
| | 120 | Locking spring | R4600024 |
| 15340025 | | Bottle section, complete | |
| | 50 | Lid for 300 ml bottle, complete | R4600048 |
| | 50 | Bottle with lid, 300 ml for DP-suspensions | 375MP027 |
| | 55 | Bottle, 1 l (Red, Blue & Green lub.), complete | 15340027 |
| | 60 | Bottle, 1 l (OP-suspens.), complete | 14600046 |
| | 70-80 | Tilting bottle screen | 375MP016 |
| | 110- | Tygon tube, 2m ENFT 21 ø2.06 | 2NU91221 |
| | 140 | | |
| | 150 | Silicone tube ø2/ø4 (10 m) | 14602901 |
| | 160 | Tygon tube, 2m ENFT 21 ø2.06 | 2NU91221 |
| | 170 | Silicone tube ø2/ø4 (10 m) | 14602901 |
| | 180- | Tygon tube, 2m ENFT 21 ø2.06 | 2NU91221 |
| | 190 | | |
| | 200 | Silicone tube ø4/ø6 (10 m) | 14602902 |

Prepomatic-2
Spare Parts and Diagrams


| Drawing | Pos. | Spare part | Cat no: |
|----------|-------|---|-----------|
| 13750460 | | Cleaning nozzle block, complete | |
| | 20 | Air nozzle 6508, complete | 375MP099 |
| | 30 | Water nozzle 6502, complete | 375MP100 |
| | 70-80 | Soap/alcohol nozzle, complete | 375MP098 |
| | 90 | Air hose, blue PU-4 (2m) | 2NU14041 |
| | 100 | Air hose, black PU-4 (2m) | 2NU14042 |
| | 110- | Air hose, blue PU-4 (2m) | 2NU14041 |
| | 120 | | |
| | 130- | Air hose, black PU-4 (2m) | 2NU14042 |
| | 140 | | |
| | 150 | Silicone tube ø2/ø4 (10 m) | 14602901 |
| 13750529 | | Carousel drain, complete | |
| | 70-80 | Drain hose ø32 (3m) | 2NU29401 |
| | 100- | Drain hose ø32 (3m) | 2NU29401 |
| | 110 | | |
| | 260- | Emergency stop push button w. yellow disc & | 381MP439 |
| | 280 | M22x1.5 brass nut | |
| | 300 | Proximity sensor DJ 2G (inductive) | 2HQ00023 |
| 13750721 | | Recirculation cooling unit, complete | |
| | 10 | Lid for tank | 12530021 |
| | 20 | Pump well w. shield | 253MP005 |
| | 30 | Plastic tank, 35 l | 253MP003 |
| | 50 | Dampener | 12530031 |
| | 90 | M2, recir. Motor 3x200-415v/50-60Hz | 253MP001A |
| | 90 | M2, recir. Motor 3x420-480v/50-60Hz | 12530068 |
| | 200 | O-ring ø21.20-3.00 (10 pcs) | 21030008 |
| | 210 | Water lever sensor | 03756912 |
| 14290506 | | Coupling (female), complete | |
| | 20 | Compression spring (3 pcs) | 14290510 |
| | 30 | Spring taps (3 pcs) | 14290511 |
| | 40-50 | Disc w. screws | 14292900 |
| | 60 | Balls KU9 (3 pcs) | 429MP007 |
| | 70 | Key | 14290513 |
| | 80 | Special screw | 14290514 |
| 15340050 | | Coupling (male) for sample holder | |
| | | Head for specimen holder | 05346902 |
| | 20 | Coupling stud | 14290580 |
| | 40 | Thrust pad | 375MP062 |
| | 60 | V-ring VA-130 | 2IV10310 |
| | 110 | O-ring ø8.0x2.00 VITON (3 pcs) 2x required | RIO20012 |
| | 120 | O-ring ø7.0x3.00 VITON (3 pcs) | RIO30004 |

Prepomatic-2
Spare Parts and Diagrams

| Drawing | Pos. | Spare part | Cat no: |
|----------|------|---|----------|
| 15340051 | | Coupling (male) for single sample holder | |
| | | Head for specimen mover plate | 05346903 |
| | 20 | Coupling stud | 14290580 |
| | 50 | V-ring VA-130 | 2IV10130 |
| | 100 | Trust pad | 375MP062 |
| | 110- | Knobs & screws (3 pcs) | 429MP018 |
| | 130 | | |
| | 140- | Spacer (3 pcs) | 429MP019 |
| | 150 | | |
| 15340052 | | Levelling device, complete | |
| | | Levelling device for Prepomatic | 05346904 |
| | 30 | Friction block | 15250234 |
| 15340736 | | Rack with PCB's, complete | |
| | 30 | Main switch | 2SE00149 |
| | 110 | SMP-bus (pos 101-110) | 375MP230 |
| | 120 | CPU PCB w. PROM, Prep-2 pos 101 | 15343254 |
| | 120 | Prep-2, set of PROM's, ver. 2.x | R5343901 |
| | 130 | Output PCB, pos 102 (see 15343170) | R4593342 |
| | 140 | Output PCB, pos 103 (see 15343170) | R4593342 |
| | 150 | Input PCB, pos 105 (see 15343170) | R4593345 |
| | 160 | Output PCB, pos 107 (see 15343170) | R4593342 |
| | 180 | Power supply PCB, 24VAC/5.5vDC | 2PI05017 |
| | 200 | Pump output PCB, 24VDC, pos 201 | 15343229 |
| | 210 | Output PCB DC to 2x9 AC, pos 202 | R5022366 |
| | 220 | Input PCB for inductive and optical switches, pos 204 | 375MP233 |
| | 230 | PCB, motor control 24vDC | 429MP075 |
| | 240 | Power supply Prep-2, pos 210 | 15343221 |

Prepematic-2
Spare Parts and Diagrams

| Drawing | Pos. | Spare part | Cat no: |
|----------------|-------------|--|----------------|
| 15340738 | | High power module, complete | |
| | 20 | Transformer T1, 200 VA | 375MP219 |
| | 70 | Contactor CI-9, 24v AC | 2KM10232 |
| | 80 | Auxiliary switch | 2KH00111 |
| | 90 | Barring device for contactor | 2KH10100 |
| | 100 | FT1 thermo relay (1.2 - 1.9A) 200-240v/50-60Hz | 2FS00125 |
| | 100 | FT1 thermo relay (0.19-0.29A) 380-415v/50Hz & 460-480v/60Hz | 2FS00123 |
| | 110 | FT2 thermo relay (0.40-0.63A) 200-240v/50-60Hz | 2FS00123 |
| | 110 | FT2 thermo relay (0.19-0.29A) 380-415v/50Hz & 460-480v/60Hz | 2FS00121 |
| | 120 | FT3 thermo relay (1.8-2.8A) 200-240v/50-60Hz | 2FS00127 |
| | 120 | FT3 thermo relay (1.2 - 1.9A) 380-415v/50Hz & 460-480v/60Hz | 2FS00125 |
| | 130 | FT4 thermo relay (4.0-6.2A) 200-240v/50-60Hz | 2FS00129 |
| | 130 | FT4 thermo relay (1.8-2.8A) 380-415v/50Hz & 460-480v/60Hz | 2FS00127 |
| | 141 | 2.00A fuse DIAZED E16 (5 pcs) | RFC01101 |
| | 141 | 4.00A fuse DIAZED E16 (5 pcs) | RFC01375 |

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











* Only for PAMTA (4x polishing. + 1x grinding)

Rev.A: (16-01-98 FTH)
PAMIF & PAMTA on sam
diagram

PREPAMATIC-2, CIRCUIT DIAGRAM - TOP ENCLOSURE

| | | | | | |
|-------------------------|------------|-----------|--------|-----------------|--------------------|
| | Size A3 | CAGE Code | DWG NO | 15343180 | Rev I I A |
| Thursday, June 08, 2000 | Scale | | | Sheet 1 of 1 | |

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
3065

Thursday, June 08, 2009

Y1, SOLONOID VALVE FOR WATER CLEANING STATION

Y3, SOLONOID VALVE FOR AIR CLEANING STATION

Y4, SOLONOID VALVE FOR ALKOHOL CLEANING STATION

Y5, SOLONOID VALVE FOR SOAP CLEANING STATION

Y6, SOLONOID VALVE FOR CARUSEL UNLOCK

Y7, SOLONOID VALVE FOR CARUSEL LOCK

Y8, SOLONOID VALVE FOR DRESSER DOWN

(Only for PAMTA (4x polishing + 1x grinding))

K1, CONTACTOR FOR HIGH PRESSURE PUMP MOTOR

K2, CONTACTOR FOR RECIRCULATION PUMP MOTOR

K3, CONTACTOR FOR SAMPLE HOLDER ROTATION MOTOR (C.W.)

K4, CONTACTOR FOR SAMPLE HOLDER ROTATION MOTOR (C.C.W.)

K5, CONTACTOR FOR DISK MOTOR LOW SPEED


K6, CONTACTOR FOR DISK MOTOR HIGH SPEED

Rev.A: (16-01-08 FTH)
PAMTA & PAMTA on same diagram

STRATEGIS AG
VALCH-SALLE 126
DK-2610 RODEBOE
DENMARK
PHONE +45 3670 8000


PREPMATIC-2
CIRCUIT DIAGRAM
POS 103-202

| Size | CAGE Code | WG NO | Page |
|-------------------------|---------------|----------|--------|
| D | < Cage Code > | 15343066 | 1 A |
| Thursday, June 08, 2000 | Scale | Cheat | 0 of 0 |

* Only for PAMTA (4x polishing + 1x grinding)

| | |
|---|---|
| Rev A (15-01-08 FT4) PAMTA & PAMTF on same diagram | STRUKTUR TECH A/S VALDEØDE ALLE 176 DK-3400 FREDSVÆRDE DENMARK PHONE: +45 3570 8000 |
| PREPAMAT IC-2 | |
| CIRCUIT DIAGRAM | |
| POS 105-2 04 | |
| Sheet D CAGE Code < Cage Code> | DWG NO 15343 068 |
| Thursday, June 08, 2000 | Date FTH/FTH |
| Code | Sheet 1 of 1 |

| | | | |
|---|---------------|--|-----|
| | | STRATEGIS AG VALCHER STRASSE 176 DK-2610 RODEBODIE DE-2610 RODEBODIE PHONE +49 45 30 800 | |
| PREPAMATIC-2 CIRCUIT DIAGRAM POS 107-207 | | | |
| Size | CAGE Code | WG NO | Rev |
| D | < Cage Code > | 15343070 | I |

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






















PREPAMATIC-2, POSITION 102-201, SIGNAL DESCRIPTION

| POSITION 102 (2x9 output type 205P261-1) | | W102 | | POSITION 201 (8x +/-24Vdc, 2x +24Vdc output) | | Power: c30:0Vac, c28:24Vac | | |
|--|------------------------|-----------|---------|--|------------------------|----------------------------|--------------------|-----------------|
| Port adr. | Function | X2 | | Connector | Type | X1 | W21 Lead colour | Actuator |
| | | collector | emitter | input | GND | + | - | |
| 01:B0 | Dosing pump1 - forward | a03 | c03 | X3-02 | X3-01 +/- 24Vdc driver | a2 | c2 | RD1 M11+ |
| 01:B1 | Dosing pump1 - back | a04 | c04 | X3-04 | X3-03 | | | RD2 M11- |
| 01:B2 | Dosing pump2 - forward | a05 | c05 | X3-06 | X3-05 +/- 24Vdc driver | a4 | c4 | RD3 M12+ |
| 01:B3 | Dosing pump2 - back | a06 | c06 | X3-08 | X3-07 | | | RD4 M12- |
| 01:B4 | Dosing pump3 - forward | a07 | c07 | X3-10 | X3-09 +/- 24Vdc driver | a6 | c6 | RD5 M13+ |
| 01:B5 | Dosing pump3 - back | a08 | c08 | X3-12 | X3-11 | | | RD6 M13- |
| 01:B6 | Dosing pump4 - forward | a09 | c09 | X3-14 | X3-13 +/- 24Vdc driver | a8 | c8 | RD7 M14+ |
| 01:B7 | Dosing pump4 - back | a10 | c10 | X3-16 | X3-15 | | | RD8 M14- |
| 02:C1 | Water valve for OPS | a12 | c12 | X3-20 | X3-19 + 24Vdc driver | a18 | a20 | RD15 Y2 RD16 |
| 02:C7 | not used | a21 | c21 | X2-02 | X2-01 + 24Vdc driver | c18 | c20 | |
| 00:A0 | Dosing pump5 - forward | a23 | c23 | X2-06 | X2-05 +/- 24Vdc driver | a10 | | RD9 M15+ |
| 00:A1 | Dosing pump5 - back | a24 | c24 | X2-08 | X2-07 | c10 | RD10 | M15- |
| 00:A2 | Dosing pump6 - forward | a25 | c25 | X2-10 | X2-09 +/- 24Vdc driver | a12 | c12 | RD11 M16+ |
| 00:A3 | Dosing pump6 - back | a26 | c26 | X2-12 | X2-11 | | RD12 | M16- |
| 00:A4 | Dosing pump7 - forward | a27 | c27 | X2-14 | X2-13 +/- 24Vdc driver | a14 | c14 | RD13 M17+ |
| 00:A5 | Dosing pump7 - back | a28 | c28 | X2-16 | X2-15 | c17 | RD14 | M17- |
| 00:A6 | not used | a29 | c29 | X2-18 | X2-17 +/- 24Vdc driver | a16 | c16 | |
| 00:A7 | not used | a30 | c30 | X2-20 | X2-19 | | | |

PREPAMATIC-2, POSITION 103-202, SIGNAL DESCRIPTION

| POSITION 103 (2x9 output type 205P261-1) | | X2 | | W103 | | POSITION 202 (2x9 24Vac triac output) | | Power: c26:0Vac, c28:24Vac, c32:12Vac | | Actuator |
|--|--|-----------|---------|-------|-------|---------------------------------------|------|---------------------------------------|-------------|----------|
| Port adr. | Funtion | collector | emitter | gate | -12V | Connector | Type | X1 | Lead colour | Actuator |
| 05:B0 | Solenoid valve f. cleaning station water | a03 | c03 | X3-02 | X3-01 | Triac output | a02 | 24Vac | 0Vac | |
| 05:B1 | Solenoid valve f. cleaning station air | a04 | c04 | X3-04 | X3-03 | Triac output | a04 | 24Vac | 0Vac | |
| 05:B2 | Solenoid valve f. cleaning station alcohol | a05 | c05 | X3-06 | X3-05 | Triac output | a06 | W22 | W22:RD1 | Y1 |
| 05:B3 | Solenoid valve f. cleaning station soap | a06 | c06 | X3-08 | X3-07 | Triac output | a08 | RD2 | W22:RD3 | Y3 |
| 05:B4 | Solenoid valve f. carousel unlock | a07 | c07 | X3-10 | X3-09 | Triac output | a10 | c | W22:RD6 | Y6 |
| 05:B5 | Solenoid valve f. carousel lock | a08 | c08 | X3-12 | X3-11 | Triac output | a12 | 02-12 | W22:RD7 | Y7 |
| 05:B6 | Solenoid valve f. dresser (PAMTA only) | a09 | c09 | X3-14 | X3-13 | Triac output | a14 | W22 | W22:RD8 | Y8 |
| 05:B7 | | a10 | c10 | X3-16 | X3-15 | Triac output | a16 | | | |
| 06:C1 | | a12 | c12 | X3-20 | X3-19 | Triac output | c14 | c | | |
| 06:C7 | | a21 | c21 | X2-02 | X2-01 | Triac output | c20 | 16-18 | | |
| 04:A0 | Contactor f. high pressure pump | a23 | c23 | X2-06 | X2-05 | Triac output | a18 | | RD41-42 | K1 |
| 04:A1 | Contactor f. recirculation pump | a24 | c24 | X2-08 | X2-07 | Triac output | a20 | | RD43-44 | K2 |
| 04:A2 | Contactor f. sample holder motor, clockwise | a25 | c25 | X2-10 | X2-09 | Triac output | a22 | RD40 | RD45 | K3 |
| 04:A3 | Contactor sample holder motor medløb | a26 | c26 | X2-12 | X2-11 | Triac output | a24 | c | RD46 | K4 |
| 04:A4 | Contactor f. disc motor low speed | a27 | c27 | X2-14 | X2-13 | Triac output | a26 | 22-26 | RD47-48 | K5 |
| 04:A5 | Contactor f. disc motor high speed | a28 | c28 | X2-16 | X2-15 | Triac output | a28 | | RD49-50 | K6 |
| 04:A6 | Solid state relay f. hot air drying (option) | a29 | c29 | X2-18 | X2-17 | Triac output | a30 | | | |
| 04:A7 | Solid state relay f. hot air drying (option) | a30 | c30 | X2-20 | X2-19 | Triac output | a32 | | | K10 |

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Herewith declares that

| | |
|---------------|--------------------------------|
| Product Name: | Prepomatic-2 |
| Type No.: | 534 |
| Machine Type: | Grinding and polishing machine |

is in conformity with the provisions of the following directives:

Safety of Machinery 98/37/EEC according to the following standard(s):
EN ISO 12100-1:2005, EN ISO 12100-2:2005, EN ISO 13850:2006, EN 60204-1:2006,
EN 954-1:1999.

EMC-Directive 2004/108/EEC according to the following standard(s):
EN 61000-6-1:2002, EN 61000-6-3:2002.

Low Voltage Directive 2006/95/EEC according to the following standard(s):
EN 60204-1:2006.

Supplementary Information The equipment complies with the American standards:
FCC part 15, subpart B and UL508.

The above has been declared according to the global method, module A

Date: 27.10.2008

Christian Skjold Heyde,
Vice President, R& D and Production, Struers A/S

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Christian Skjold Heyde,
Vice- President, R & D et Production, Struers A/S

Prepomatic-2

Fully automatic
preparation of
materialographic
specimens

- from plane grinding to
cleaning and drying



- Fully automatic - allowing unattended operation
- Sturdy construction - vibration-free preparation
- Use of grinding stone - high material removal
- Unsurpassed reproducibility
- Preparation of single specimens and specimens clamped in holders
- Peristaltic pumps for accurate dosing of suspensions and lubricants
- Controlled material removal
- Metalog Guide preparation methods included
- Low cost per sample

Prepomatic-2 is a fully automatic, microprocessor-controlled machine for the preparation of all kinds of materials. Prepomatic-2 is especially well suited for laboratories where high reproducibility is essential, and laboratories handling difficult specimens of a variety of materials. Prepomatic-2 will prepare from 1 to 3 single specimens at a time (or 3-9 small specimens in special fixtures). Alternatively, 3 or 6 specimens can be prepared clamped in specimen holders.

Automatic operation means identically prepared specimens and, consequently, outstanding reproducibility. Another advantage of this principle is precisely controlled consumption of diamond suspensions and lubricants. Therefore, Prepomatic-2 offers very good value for money when it comes to daily costs.

As Prepomatic-2 works completely unattended, the operator can spend his precious time analysing the specimens instead of preparing them.

Prepomatic-2 is available in two versions, one with five grinding/polishing stations and one with a grinding stone in the first position followed by 4 grinding/polishing stations.

The automatic dosing system with peristaltic pumps guarantees aerosol-free, precise dosing of both suspension and lubricant

Fully automatic operation

Prepomatic-2 takes care of the entire preparation process, from plane grinding to final polishing - including cleaning and drying - without operator attendance.

A graphic L.C. display on the operating panel shows all relevant information and allows the operator to select a method. The operator then mounts the grinding and polishing discs specified on the display, inserts the specimens - and presses Start. When this is done, he may attend to other jobs, e.g. microscopy - returning to the machine only to remove the finished, cleaned and dried specimens.

Sturdy construction for vibration-free preparation

Prepomatic-2 is a floor model, built on a strong, welded steel frame. On the carousel five grinding and/or polishing stations and one, completely encapsulated, cleaning station are mounted.

Depending on the actual preparation step, the carousel turns and the respective disc or cleaning station is positioned underneath the specimen mover. This sturdy construction allows for the preparation of all materials

without unwanted vibrations, giving the best possible preparation results.

Use of grinding stone - high material removal

For the preparation of most ferrous metals a grinding stone will allow for the fastest possible plane grinding. The grinding stone is completely encapsulated, thus allowing for a speed of 1,000 rpm. This high speed results in high material removal rates and short grinding times.

The necessary dressing of the stone is carried out automatically prior to the grinding step.

Unsurpassed reproducibility

As all preparation parameters are controlled by the microprocessor the same results will be obtained time after time. Thus changes in the structure can be related to actual changes during the production of the material instead of possible variations in the preparation method.

Preparation of single specimens and specimens clamped in holders

With easy exchangeable heads both single specimens and specimen holders can be prepared.


For the occasional specimen or to repeat just the final polishing step, single specimen mode is the most advantageous. There is no reason to waste consumables by preparing "dummy" specimens because a specimen holder has to be filled.

On the other hand, if 3 or 6 specimens have to be prepared, the use of a specimen holder provides several advantages. Handling is facilitated, as all specimens are treated as a single one.


Additionally, the accuracy when measuring the removal rate is improved and overall flatness of the specimens is better.

Peristaltic pumps for dosing of suspensions and lubricant

For the accurate dosing of both diamond suspensions, oxide polishing suspensions and lubricants, peristaltic pumps are used. They not only supply the correct amount exactly where it is needed, they also completely avoid the creation of aerosols. Thus the choice of peristaltic pumps is also very environmentally friendly.

The control panel and its large display makes Prepomatic-2 very easy to operate.

Prepomatic-2 can handle both single specimens and specimens clamped in holders.

Controlled material removal

For the examination of electronic components or parts for failure analysis the defined removal of material is crucial.

On Prepomatic-2 an advanced strain gauge is used to very accurately control the material removal. This allows for target preparation in the cases mentioned above, but it can also be used to ensure that enough material has been removed, i.e. scratches and deformed material introduced in the preceding step.

This system saves time and consumables, and it prevents an entire preparation cycle from being performed on an inadequately ground specimen.

Metalog Guide Methods included

The Metalog Guide Methods, covering the preparation of most materials are included in the software of Prepomatic-2.

When new materials have to be prepared a suitable method can easily be found, and if necessary changed and saved for future use. Thus, even operators with limited preparation experience can prepare all kinds of mate-

rials immediately. Modified programs can be saved in one of seven groups, each containing ten methods.

Low cost per sample

Because of the fully automatic preparation it is possible to optimise each method and to get the necessary quality in the shortest possible time. This, together with the automatic dosing of consumables results in very low, consistent costs per prepared specimen.

There are a number of other important features built into Prepomatic-2:

Easy programming

To establish a new preparation method the user simply transfers one of the Metalog Guide methods to an "open method" and changes the necessary parameters. All methods can be used for the preparation of 1, 2 or 3 specimens. The necessary parameters are recalculated automatically if a different no. of specimens have to be prepared.

Process control

Once a preparation program has been selected and Start has been pressed, Prepomatic-2 carries out all necessary tasks. However, the user has several additional possibilities to monitor the preparation of his specimens. When developing new methods, the preparation can automatically be stopped after a pre-selected step to check the quality of the specimen, different cleaning programs can be selected and


steps can be omitted or repeated. A control program prevents incorrect operation.

Wide range of specimen mover plates and holders

A wide range of different specimen mover plates for specimens 25, 30, 40 and 50 mm, 1", 1 1/4" and 1 1/2" diameter is available. In addition, specimen holders for 3 or 6 mounted or unmounted specimens and five small fixtures for unmounted specimens are available, making it possible to prepare up to nine specimens at a time.

Safety

Prepomatic-2 is provided with all necessary safety guards and complies with the requirements for machine tools (IEC 204, EN60204 -1).

Technical Data

| | | |
|---------------------------------------|--|---|
| Three-phase | Fuse: Compressed-air supply: Consumption: Water supply: Grinding/polishing discs: | max. 16 A min. 5 bar, max. 9.9 bar max. 60 l/min min. 1.0 bar, max. 10.0 bar 175-200 mm dia. |
| Speed: | Grinding/polishing discs: Prepomatic-2 with grinding stone: Prepomatic-2 with five grinding/polishing steps: | 75, 150, 300, 500, 1000 rpm (two rates of speed per station) 75, 150, 300 rpm (two rates of speed per station) |
| Specimen holder: | 125 rpm, bi-directional rotation | |
| Specimen height: | Min. 12 mm, max. 32 mm | |
| Grinding/ polishing force: | 5-150 N | |
| Dimensions: | Height: Width: Depth: Weight: | 1710 mm 1130 mm 940 mm 430 kg |

Specifications

Code

Prepomatic-2

Microprocessor controlled, fully automatic grinding and polishing machine with 5 grinding/polishing stations. Complete with recirculation cooling unit (capacity 30 l) and coupling unit (PAMUN). Heads, specimen mover plates, specimen holders, levelling device, MD-Discs for magnetic fixation or polishing discs are ordered separately.

PAMIF

Prepomatic-2

Microprocessor controlled, fully automatic grinding and polishing machine with 1 plane grinding station and 4 grinding/polishing stations. Complete with grinding stone and recirculation cooling unit (capacity 30 l) and coupling unit (PAMUN). Heads, specimen mover plates, specimen holders, levelling device, MD-Discs for magnetic fixation or polishing discs are ordered separately.

PAMTA

Head for Specimen Holders

For mounting of specimen holders on Prepomatic-2 or Prepomatic.

PAMHO

Specimen Holder

For clamping of 3 mounted or unmounted specimens 24-40 mm dia., aluminium.
To be mounted in head for specimen holder (PAMHO).

PAMLU

Specimen Holder

For clamping of 6 mounted or unmounted specimens 20-25 mm dia., aluminium.
To be mounted in head for specimen holder (PAMHO).

PAMLO

Specimen Holder

For clamping of 6 hot or cold-mounted specimens 23-25.4 mm / 1" dia., aluminium.
To be mounted in head for specimen holder (PAMHO).

PAMBA

Specimen Holder

For clamping of 6 hot or cold-mounted specimens 28 - 30 mm dia., aluminium.
To be mounted in head for specimen holder (PAMHO).

PAMBE

Head for Specimen Mover Plates

For mounting of specimen mover plates on Prepomatic-2 or Prepomatic.

PAMMO

Specimen Mover Plates

Set with specimen mover plates for the preparation of 3 x 25 mm / 1", 3 x 30 mm / 1¼", 3 x 40 mm, 3 x 1½", and 1 x 50 mm single specimens. Are mounted in head for specimen mover plates (PAMMO).

PAMDI

Levelling Device

For levelling of specimens in specimen holders or specimen mover plates.

PAMFU

Specimen Fixtures

For inserting into Prepomatic-2 or Prepomatic specimen mover plates.

PAMFI

For 3 specimens 7-15 mm dia., for specimen mover plate with holes 40 mm dia.

PAMME

For 3 specimens 12-20 mm dia., for specimen mover plate with hole 50 mm dia.

PAMIK

For 3 specimens max. 9 x 19 mm, for specimen mover plate with holes 40 mm dia.

PAMAX

For 3 specimens max. 12 x 23 mm, for specimen mover plate with hole 50 mm dia.

PAMON

For 1 specimen max. 23 x 26 mm, for specimen mover plate with holes 40 mm dia.

PAMST

Grinding Stone

200 mm dia., grit 150 for Prepomatic-2 (PAMTA)

PAMST

MD-Disc

For magnetic fixation 200 mm dia.,

DEMLA

DP-Polishing Disc

200 mm dia.

PAMLA

Struers' equipment is in conformity with the provisions of the applicable International Directives and their appurtenant Standards. (Please contact your local supplier for details)

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup, Denmark
Phone +45 44 600 800
Fax +45 44 600 801
struers@struers.dk

USA
Struers Inc.
24766 Detroit Road
Westlake, OH 44145-1598
Phone +1 440 871 0071
Fax +1 440 871 8188
info@struers.com

CHINA
Struers Ltd.
Office 702 Hi-Shanghai
No. 970 Dalian Road
Shanghai 200092, P.R. China
Phone +86 (21) 5228 8811
Fax +86 (21) 5228 8821
struers.cn@struers.de

CANADA
Struers Ltd.
7275 West Credit Avenue
Mississauga, Ontario L5N 5M9
Phone +1 905-814-8855
Fax +1 905-814-1440
info@struers.com

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