Prontopress-2

Instruction Manual
Gebrauchsanweisung
Mode d'emploi

Valid from Serial No.: 3772001

Manual No.:13777003

Date of Release: 09.05.1994

Safety Precaution Sheet

To be read carefully before use

1. The operator should be fully aware of the use of the apparatus according to the Instruction Manual.
2. Be sure that the actual voltage corresponds to the voltage stated on the back of the apparatus and the heating/cooling unit.
3. Be sure that the water tubings/connections are without leaks.
4. Be sure that the heating/cooling unit is correctly mounted on the press with the mounting cylinder before starting the process.
5. Be sure that the top closure is correctly mounted, before starting the process.
6. Keep the mounting pressure below the red area of the manometer.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

IMPORTANT

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the *Serial No.* and *Voltage* on the front cover of this manual or on an attached label below. If in doubt consult the rating plate of the machine itself.

We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals

A Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals

A Struers Service Manual may only be used by a trained technician authorized by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 1994.

Struers A/S

Valhøjs Allé 176
DK-2610 Rødovre/Copenhagen
Denmark
☎ Telephone +45 36 70 80 90
☎ Fax +45 31 41 65 44

Table of Contents

Page

1.	Installation	5
2.	Operation	7
3.	Accessories and Consumables	12
4.	Maintenance	15
5.	Trouble-Shooting	16
6.	Technical Data	17

1. Installation

1.1 Packing List

- 1 Prontopress-2
- 1 Tube spanner for fixing bolts, width 13 mm
- 1 Allen key 5 mm for oil screw, fig. 2.9, and draining screw.
- 2 Spanners for securing the lower ram (19 and 22 mm)
- 1 Bayonet closure (fig. 1.12)
- 1 Funnel
- 1 Measuring spoon
- 1 Water hose
- 1 Filter fitting
- 1 Fitting
- 1 Transmission fitting
- 1 Hose for water outlet (1,5 m)
- 1 Hose for connecting water from base unit to embedding unit (0.65 m)
- 5 Hose clips
- 1 Fuse
- 3 Resin samples (resin 3, resin 5 and resin 6)
- 1 Mould Release Wax

1.1.1 Packing List for Mounting Cylinder (Accessory)

The easily exchangeable cylinders come in the following sizes: 25, 30 and 40 mm and 1", 1¼" and 1½" inside diameter. Apart from the cylinders the case further contains:

- 1 Lower ram
- 1 Upper ram
- 1 Syringe with heat conductive paste

1.2 Before First Time Use of the Press

The transit screw on the rear plate of the press mentioned in the text under "See manual" must be unscrewed before using the press. Under the screw in the fitting situated in the rear plate is placed a piece of felt (meant) for purifying the air intake in the press. When the press is working the hydraulic oil gets hot and it should be possible for it to expand without causing an excess pressure in the press.

Keep the screw for later user when moving the press from the premisses. Do not mount another screw with head or the like.

1.3 Connection to Power, Removal of Oil Screw

Check the voltage and frequency on the identification plate on the rear side of the press.

Connect the base unit to the mains voltage with a single-phase plug to be connected as follows:

Yellow/green: earth
Brown: phase
Blue: neutral

Remove oil screw (fig. 2.9) to allow the press to "breathe". A felt disc under the screw prevents dirt from entering. Keep the screw for mounting in case of transport of press.

1.4 Connection of Cooling Water

Mount the hoses as follows:

- The short hose between the branches fig. 2.4 and 2.5.
- The outlet hose on the branch fig. 2.6.
- The inlet hose on the branch fig. 2.3.
- The inlet hose is connected to a water tap.
- The outlet hose is taken to a drain.

When the hoses are connected to branch fig. 2.5 and fig. 2.6, make sure that they are *not* pushed all the way up towards the back plate of the mounting unit. The hose clamps should be firmly tightened to avoid water leakage.

1.5 Assembly and Exchange of Mounting Cylinder

Screw the lower ram onto the piston rod using the two keys 19 and 22 mm. Pull out the locking pin on the back of the flange (fig. 2.8). Place the cylinder over the lower ram and screw it into the threaded ring until the locking pin engages with the hole in the cylinder. Smear a layer of 'Heat conductive paste' on the outer surface of the cylinder. Lower the embedding unit onto the cylinder, after which the press is ready for operation.

To exchange the cylinder lift off the embedding unit and pull out the locking pin. Unscrew the cylinder and demount the lower ram. Mount the new cylinder and ram as described above. If it proves difficult to pull the unit off the cylinder, cool the cylinder down and then heat it for 10-20 seconds. This causes the unit to expand, and it may now easily be pulled off.

2. Operation

2.1 Introduction

Prontopress-2 is an automatic hydraulic mounting press for metallographic laboratories.

2.2 Embedding Unit

The press is supplied with an embedding unit incorporating a heating/cooling element. Moulding cylinders and rams are delivered separately according to requirement.

The embedding unit will accept standard moulding cylinders of inside diameters between 25 and 40 mm including 1", 1 ¼" and 1 ½". A thermistor ensures accurate control of the temperature of the heating/cooling element.

The pressure in the moulding cylinder is established between the upper and lower rams. The upper ram is kept down in the cylinder by means of a bayonet closure on the upper rim of the cylinder, and the lower end of the cylinder is held in place on the base unit by a threaded ring. The lower ram is placed on the piston rod, see below.

2.3 Base Unit

The base unit consists of a motor driven hydraulic pump, a high pressure cylinder and automatic controls for compressive force, heating- and cooling time and heating temperature.

The compressive force is produced by a hydraulic system comprising a pump, a high-pressure cylinder and a pressure regulator.

The pump is driven by a reversible electric motor. When the motor runs in one direction, the oil is pumped from the reservoir into the primary side of the high-pressure cylinder, and the piston rod ascends. When the motor reverses, oil is pumped into the secondary side of the high-pressure cylinder, whereby the piston rod descends.

2.4 Controls

Power Switch (fig. 2.10)

Switch on the press by pressing the power switch (fig. 2.10). The main pilot light (fig. 1.8) turns on.

Three-way Toggle Switch (fig. 1.2)

Controls the position of the piston rod. Only when the switch is set to the neutral position in the middle the press can perform the embedding cycle.

Auto Start (fig. 1.1)

Starts the embedding cycle, with the preset heating/cooling time.

Timers (fig. 1.4-5)

Heating and cooling times can be preset, but cannot be adjusted during the operation of the press. During heating the (fig. 1.9) light comes on. During cooling the (fig. 1.10) light comes on.

"Force" (fig. 1.6)

The pressure can be preset by turning this knob and it can be adjusted while the press is operating. The pressure can be read on the manometer and should not exceed 40 kN (fig. 1.7).

Knob °C (fig. 1.3)

The process temperature can be adjusted by turning the knob to the desired temperature. Normal setting: 180°C.


Thermo Fuse (fig. 2.7)

The thermo fuse protects the pump motor against overload. Press for reset.

Melting Fuse (fig. 6.86)

Protects the heating/cooling element against overheating.

2.5 Mounting

Take the lower ram to its top position by setting the three-way switch (fig. 1.2) to its upper position. Apply Mould Release Wax to the lower ram and place the specimen on the lower ram. Then take the latter to its bottom position by setting the three-way switch to the lower position. Pour resin into the cylinder and apply Mould Release Powder to the upper part of the cylinder and to the upper ram. Pass the upper ram so far down into the cylinder that the bayonet closure can engage with the cylinder and then turn the bayonet closure 90° clockwise by turning the handle (green point on the bayonet closure should be in position opposite the corresponding green point on the heating/cooling unit). Set the process parameters (see section 6.2). Press  (fig. 1.1), after which the entire operation will proceed automatically. When the cooling has stopped, the specimen can be removed. Set the three-way switch (fig. 1.2) briefly to "down". Release the bayonet closure by turning it 90° back using the handle (green point on the bayonet closure should be in position opposite the corresponding green point on the heating/cooling unit). Set the three-way switch to "up" for ejection of bayonet closure and specimen from the cylinder. When the specimen has been ejected, set the three-way switch to "neutral", and the press is ready for the next mounting.

2.6 Mounting Data

The data stated below should be considered as a guide for obtaining a perfect, completely cooled mount.

Depending on the individual demands to the mount, the specimen size and the type and quantity of resin, the listed times can be optimized. See also the chapter "Hints".

Suggested data for all resins:

Cylinder Dia.	Quantity	Pressure	Heating Time	Cooling Time
25 mm/1"	20 ml	15 kN	6 min	4 min
30 mm/1 1/4"	30 ml	20 kN	8 min	4 min
40 mm/1 1/2"	40 ml	35 kN	10 min	4 min

NB Thermostat on maximum.

When using Resin 3, the supply of cooling water must be reduced to 0.1 l/min.

2.7 Hints

Resin 1 and 3

Resin 3 is not electrically conductive, can be used in conjunction with Resin 1 to provide an insulating surface with Resin 1 as the conductive backing.

Resin 4, 5 and 6

When using the rather expensive resins 4 and 5, these can be used in the bottom of the cylinder, which is then filled up with the less expensive bakelite resin.

Pre-mount

Pre-mount is best suited for routine mounting of large numbers of specimens where, due to their compact form and the addition of an intermediate ram, two mounts can be produced in one process (max. \varnothing 40 mm).

Because of its compact form, Pre-mount already exerts a certain pressure on the specimen from the start of the mounting process, making it suitable for use only with pressure-resistant specimens.

Mould Release Wax

Must always be applied to the mounting rams in a thin layer before the mounting process begins. This prevents the adhesive qualities of the resins from making it difficult to remove the mounts afterwards.

Distance to cylinder wall

The distance between the specimen and the cylinder wall must be minimum 3 mm (1/8") to avoid cracks in the resin. This goes especially for specimens with sharp edges.

Small specimens

Small thin specimens can be stabilized by means of Struers mounting clips, available in metal and in plastic (see the chapter "Consumables").

Clean specimens

The specimens must be clean, dry and free from grease. If necessary clean with alcohol or another suitable degreasing fluid.

Pressure

The recommended maximum pressures vary from 15 kN for \varnothing 25 mm to 35 kN for \varnothing 40 mm specimens. The pressure must never exceed 55 kN and be less than 15 kN.

Preheating

For porous and/or pressure sensitive specimens, such as minerals, electronic parts etc., it is an advantage to soften the resin by heating, before applying pressure.

This is mostly advantageous when using thermoplastic resins and especially Resin 3.

Example of time parameters to be used together with Resin 3:

Preheating with low pressure: 12 min

Heating with pressure: ½ min

Cooling: 5 min

<i>Heating time</i>	The heating time can be optimized depending on the specimen size and the quantity of resin used.
<i>Cooling time</i>	A cooling time of 4-5 min is mostly suitable. Note that when Resin 3 is used, the supply of cooling water should be reduced to 0.1 l/min. It is also sometimes necessary to increase the cooling time for Resin 3.
<i>Temperature</i>	The mentioned heating/cooling times are based on the thermostat being set on maximum. If the temperature is reduced the recommended heating time should be longer. The temperature used during the mounting process can be reduced to a minimum of 110°C for all resins.
<i>Porous specimens</i>	Thermoplastic resins (especially Resin 3) penetrate well into porous specimens. The best results are obtained by preheating the resin (see "Preheating"). <i>NB</i> For very porous specimens, hot compression mounting should be avoided. Cold mounting should be used instead (Epofix with vacuum impregnation).
<i>Pressure sensitive specimens</i>	Thermoplastic resins (especially Resin 3) should be used. The best results are obtained by preheating the resin (see "Preheating"). <i>NB</i> For very pressure sensitive specimens, hot compression mounting should be avoided. Cold mounting should be used instead.
<i>Temperature sensitive specimens</i>	The temperature used during the mounting process can be reduced to a minimum of 110°C for all resins. <i>NB</i> For very temperature sensitive specimens, hot compression mounting should be avoided. Cold mounting (Epofix) should be used instead.

3. Accessories and Consumables

3.1 Accessories

Soecification	Code
<i>Top closures</i> Standard bayonet closure	PREBA
Advanced bayonet closure. Makes it possible to place the upper ram individually, according to the height of the specimen. Built-in device secures easy opening and closing	PRENY
<i>Mounting cylinder set</i> Mounting cylinder set, including upper ram and lower ram, interchangeable: Diameter 25 mm Diameter 1" Diameter 30 mm Diameter 1 ¼" Diameter 1 ½" Diameter 40 mm	PROMI PROON PROTE PROQU PROHA PROIF
<i>Intermediate ram</i> Makes it possible to mount two specimens simultaneously. Please state cylinder diameter	PRODO
<i>Struers Engraver, complete</i> With 1 hard point, 220 V, 50-60 Hz As above, but 115 V, 50-60 Hz	ENGRA ENGRA

3.2 Consumables

Hot compression mounting resins

Resin 1 and 3 are thermoplastic. Resin 4, 5, 6 and Pre-mount are thermosetting.

Resin	Characteristics	Application	Code (Quantity)
Resin 1	Material: acrylic Classification: thermoplastic Powder color: dark grey Filler: Fe-powder Specific property: electrically conductive	Mounting of ferrous specimens and aluminium alloys for e.g. electrolytic polishing	RESFE (net 0.75 kg)
Resin 3 (transparent)	Material: acrylic Classification: thermoplastic Powder color: white Specific properties: transparent, low viscosity	Mounting of irregularly shaped specimens and of specimens with a rough surface. Penetrates well into cracks and pores. Provides crystal-clear mounts - a major benefit if the objective is to grind to a certain point in the specimen	RESNO (net 0.75 kg) RESYV (net 7.5 kg)
Resin 4	Material: diallylphtalate Classification: thermosetting Powder color: green Filler: glass fibers Specific properties: hard, wear resistant, highly resistant to chemicals	The glass fibers result in a very hard mount. This resin is a good choice when specimens with high edge retention and planeness are required	RESDI (net 0.75 kg) RESTDY (net 7.5 kg)
Resin 5	Material: epoxy Classification: thermosetting Powder color: black Filler: minerals, glass fibers Specific properties: hard, wear resistant, highly resistant to chemicals	This epoxy resin, containing minerals and glass fibers, gives a very hard mount. Due to its outstanding adhesion to the specimen it provides an extremely good edge retention. Suitable for all purposes requiring optimal edge retention and planeness	RESIF (net 0.75 kg) RESEN (net 7.5 kg)
Resin 6	Material: phenolic (bakelite) Classification: thermosetting Powder color: brown Filler: wood Specific property: inexpensive	Serial mounting. Especially suitable for routine mounts	RESBA (net 0.70 kg) RESUK (net 7.5 kg)
Pre-mount	Material: phenolic (bakelite) Classification: thermosetting Color: black Filler: wood Specific properties: easy to handle, extremely uniform mounts Shape: cylindrical pills for cyl.dia. 25 mm/1", height 16 mm (RESUN) for cyl.dia. 30 mm/1¼", height 16 mm (RESGI) for cyl.dia. 40 mm/1½", height 16 mm (RESAM)	Serial mounting of specimens of uncomplicated shapes. The pills can be immediately placed in the mounting press. Quick, easy and clean to use	RESUN (net 10 kg/ about 1100 pcs.) RESGI (net 10 kg/ about 760 pcs.) RESAM (net 10 kg/ about 460 pcs.)

		Code
<i>Mould Release Wax</i>	Prevents the specimens from sticking to the rams, net 53 g	PREWA
<i>Heat Conductive Paste</i>	Paste for improving the heat transmission between the mounting cylinder and the heating/cooling unit, net 100 g	PREPA
<i>Fixation Clips</i>	Metal clips for retaining small specimens, 9 mm, 100 pcs.	CLIPS
	Same as above, but 6 mm, 100 pcs.	CLIMO
<i>Multiclips</i>	Plastic clips for retaining small specimens, 50 pcs.	CLIFE
<i>Sample Identification Discs</i>	Struers Sample Identification Discs, pack of 20 sheets, each containing 15 discs	IDENT

4. Maintenance

Apart from cleaning and exchange of oil the base unit requires no regular maintenance. However, moving parts may be greased.

Fuse (fig. 2.11)

For replacement of the fuse, unscrew the black plastic button and remove the fuse. Replace with fuse of the right amperage.

4.1 Cleaning

Clean with a dry cloth (do not use compressed air).

Moulding cylinders and rams must be kept clean. The cylinder may be cleaned by passing the lower ram up through the cylinder without producing a mount. For a more thorough cleaning remove the cylinder. Cylinder and rams are cleaned in alcohol or acetone or possibly in an ultrasonic cleaner. Rams may be polished with 6 µm or 3 µm diamond on a soft polishing cloth (DP-Nap).

4.2 Exchanging the Oil

The hydraulic oil must be exchanged at intervals. The first oil exchange must take place after 1000 hours of operation and, after that, at every 3000 hours of operation, or at intervals of max. 10 years.

The press is drained of oil through the drain plug in the bottom of the press. When the press has been drained of oil, the embedding unit (fig. 1.11) is removed, and the press is turned upside down. Now 3 l of oil are filled through the drain hole and the plug is mounted.

Recommended type of oil Mobil DTE 18 or similar hydraulic oil having a viscosity of 95c St (12°E) at 40°C.

5. Trouble-Shooting

Error	Cause	Action
Resin 1 or 3		
Lack of conductivity	No contact to specimen due to use of an excessive amount of Resin 3	Make a new mount. Possibly place the entire mount with specimen in a slightly larger cylinder
Lack of conductivity	Insufficient curing of resin	Reinsert mount in the press for re-processing at prolonged heating times
The metal particles are removed in the electrolytic process	Excessive grinding time	Make a new mount
	Too little amount of Resin 3	Use more Resin 3
Resin 3		
Whiteness on specimen surface or in the middle of specimen	Insufficient curing (moist resin)	Prolong heating time possibly reprocess the faulty mount. The resin may be dried by exposing the open container to 30°-70°C
Mount white or cotton balled, especially if the mount is high	Too rapid cooling insufficient curing (moist resin)	Reduce height of mount. Lower the temperature and prolonged the heating time. Reduce the cooling water supply and prolong cooling time.
Resin 4		
Bulging mount	Insufficient cooling	Increase the cooling time
Individual grains visible in the resin	Insufficient curing	Increase the heating time and/or temperature
Specimen sticks to the ram	No Mould Release Wax applied/heating too brief	Apply Mould Release Wax. Heat longer
Resin 6 and Pre-moulds		
Bulging mounts	Insufficient cooling	Increase cooling time
Parts of mount become light or colorless on contact with alcohol	Insufficient curing	Increase heating time and/or temperature

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Table of Contents

Fig. 1
Fig. 2
Fig. 3
Fig. 4
Fig. 5
Fig. 5A
Fig. 6

1.1	Front	
1.2	Back	
1.3	Pressure regulator	377-M-45R
1.4	Piston rod.	377-M-46R
1.5	Inside view	13770047A
	Inside view, photo	
1.6	Heating element	377-M-48R

The drawings are not to scale. Some of the drawings may contain position numbers not used in connection with this manual.

Drawing


Fig. 1

Pos. Description

Spare Part No.

1.1 Front

1	Process starting button377MP005
2	Control button377MP004
3, 4, 5	Adjusting button13770101
6	Pressure adjusting button377MP010
7	Pressure gauge377MP011
8	Lamp, complete377MP002
9	Heating indicator lamp377MP008
10	Cooling indicator lamp377MP009
11	Embedding unit, complete, see fig. 5377MP014
12	Advanced bayonet closure, complete (PRENY)0776910
	Spring for advanced bayonet closure (PRENY)377MP111
13	Handle for advanced bayonet closure (PRENY), complete377MP112
14	Clamping screw with knob (PRENY)377MP114
15	Black Release Knob377MP119


Drawing

Fig. 2

Pos. Description Spare Part No.

1.2 Back

1	Mains cable.377MP018
2	Number plate with indication of voltage377MP019
7	Thermo fuse for motor377MP023
8	Locking pawl for cylinder377MP039
9	Transit screw for oil reservoir377MP025
10	Fuse holder.377MP003
	Fuse (state amperage).377MP015
11	Main switch.377MP001

Drawing

Fig. 3
 Drawing 377-M-45R


Pos. Description

Spare Part No.

1.3 Pressure regulator

28	Ball377MP028
30	Spring for pressure adjusting valve377MP030
31	Stop screw with copper sealing377MP031
32	Pressure regulator with copper sealing377MP032
33	Housing for pressure regulator, consisting of housing, spindle, catch, stop ring and 2 O-rings.377MP033

A-A


Drawing

Fig. 4
 Drawing 377-M-46R

Pos.	Description	Spare Part No.
1.4 Piston rod		
37	Screw thread377MP037
38	Flange for screw thread377MP038
39	Locking pawl for cylinder, complete377MP039
40	Allen screw for piston rod377MP040
	<u>Set of sealings no. 1</u>377MP109
41	Gasket plate	
42	Holder for Teflon ring	
43	Teflon ring	
44	Brass bearing	
45	Retaining ring for sealing ring	
46	Sealing ring	
47	O-ring 74.5 x 3.0 mm	
57	Locking ring	
46	Sealing ring, 5 pcs..377MP046
48	Cylinder.377MP048
49	Piston rod.377MP049
50	Piston.377MP050
	<u>Set of sealings no. 2</u>377MP110
47	O-ring 74.5 x 3.0 mm	
51	O-ring 69.2 x 5.7 mm	
52	U-sealing ring	
53	Valve block377MP053
54	O-ring377MP054
55	Filter377MP055
56	Top flange377MP056

(Fig. 4, Piston rod, Drawing 377-M-46R)


Drawing

Fig. 5
 Drawing 13770047A

Pos.	Description	Spare Part No.
1.5 Inside view		
61	Terminal strip for tool.377MP061
63	Terminal strip for mains connection377MP063
67	Rubber feet.377MP067
72	Lower ram, please state diameter.377MP072
73	Aluminium jacket for 25 mm, 1", 30 mm and 1¼" cylinder diameter.377MP073
	Spring ring for aluminium jacket.377MP117
74	Upper ram (please state diameter)377MP074
75	Retaining ring377MP075
80	Embedding cylinder (please state cylinder diameter)377MP080
97	Motor, state voltage and frequency377MP097
98	Coupling complete (motor coupling, pump coupling and coupling house)377MP098
99	O-Ring377MP099
100	Pump (up to serial no. 3771329)377MP100
100	Pump (from serial no. 3771330 to 3771419 please state serial no.)377MP115
100	Pump (from serial no. 3771420).377MP116
101	O-Ring377MP101
104	Bayonet body.377MP104
105	Spindle handle377MP105
106	Stop screw377MP106
117	Spring.377MP117

(Fig. 5, Inside view, Drawing 13770047A)


Drawing

Fig. 5A

Pos.	Description	Spare Part No.
Inside view, photo		
1	Print board, complete with potentiometer377MP093
2	Relay, 3-pole, 110V	2KL44311
2	Relay, 3-pole, 220V	2KL44322
3	Relay, 2-pole, 110V377MP121
3	Relay, 2-pole, 220V377MP120
4	Timer potentiometer, shortened	13590621
5	Triac, V1359MP024
6	Timer, socket.	2KT90511
7	Holder spring.	2KT90077
8	Cooling timer, 110V	2KT24110
8	Cooling timer, 220V	2KT24220
9	Heating timer, 110V	2KT25110
9	Heating timer, 220V	2KT25220
10	Water valve.377MP118
11	Cable entry377MP063
12	Cable entry for mains cable377MP062

(Fig. 5A, Inside view, photo)



Drawing

Fig. 6
 Drawing 377-M-48R

Pos. Description Spare Part No.

1.6 Heating element

83	Heating element377MP083
84	Clamp for melting fuse and clamp for cable.377MP084
86	Melting fuse359MP052
87	Thermistor359MP050
88	Bell-mouthed screwed connection377MP088
89	Terminal strip.377MP089
90	Division plate of plastic.377MP090
91	Lower insulating ring.377MP091
92	Upper insulating ring.377MP092
	Measuring spoon.377MP094
	Hydraulic oil, 4 liters377MP095
	Funnel377MP096
	Female socket for Timer Kit260MP159

2. Diagrams

Table of Contents


2.1	Hydraulic diagram	377-M-28
2.2	Heat control circuit diagram	377-M-33B
2.3	Main circuit diagram, 115 V/50-60 Hz, 1400 W	13773142A
2.4	Wiring diagram, heating unit	377-M-55A
2.5	Wiring diagram, general	13773144

Diagram

Description

377-M-28

2.1 Hydraulic diagram


1. Double-acting hydraulic cylinder.
2. Pressure regulator, adjustable.
3. Pressure gauge.
4. Oil pump.
5. Electro-motor.
6. Non-return valve.
7. Non-return valve.
8. Oil tank.
9. Oil filter.
10. Pressure regulator, fixed.

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



Herewith declares that

<i>Product Name:</i>	Prontopress-2
<i>Type No:</i>	377
<i>Machine Type:</i>	Mounting Apparatus

is in conformity with the provisions of the following directives:


Safety of Machinery 89/392/EEC and 91/368/EEC and 93/44/EEC according to the following standard(s):
EN292-1:1991, EN292-2:1991, EN60204-1:1992.

EMC-Directive 89/336/EEC and 92/231/EEC according to the following standard(s):
EN50081-1:1991, EN50082-1:1992.

Low Voltage Directive 73/23/EEC and 93/68/EEC according to the following standard(s):
EN60204-1:1992.

Supplementary Information The equipment complies with the American standards:
UL508.

The above has been declared according to the global method, module A

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Klaus Kisbøll, Entwicklungsleiter, Struers A/S
Fabricant Struers A/S
Klaus Kisbøll, Chef du service developpement, Struers A/S