TegraPol

-21/ -25/ -31/ -35

Instruction Manual

Manual No.: 15557001


Date of Release: €J.0J.20F€

TegraForce-3/ -5

Instruction Manual

TegraDoser-1

Instruction Manual

FCC Notice

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the Instruction Manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Pursuant to Part 15.21 of the FCC Rules, any changes or modifications to this product not expressly approved by Struers A/S could cause harmful radio interference and void the user's authority to operate the equipment.

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manuals may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manuals may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 20F€.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

TegraPol-21/ -25/ 31/ -35 Safety Precaution Sheet

To be read carefully before use

1. The operator(s) should be fully instructed in the use of the machine, any connected equipment and accessories and the applied consumables according to the relevant Instruction Manuals.
2. The machine must be placed on a safe and stable table with an adequate working height. All functions on the machine and any connected equipment must be in working order.
3. Operators should ensure that the actual voltage corresponds to the voltage on the back of the machine. The machine must be earthed. Follow the local regulations. Always turn the power off and remove the plug or the cable before opening the machine or installing additional components like TegraForce or TegraDoser-5.
4. Connect only to cold water. Make sure that the water connections are leak-proof and that the water outlet is working.
5. Struers recommend that the mains water supply is shut off or disconnected if the machine is to be left unattended.
6. Alcohol based consumables: follow the current safety rules for handling, mixing, filling, emptying and disposal of the alcohol-based liquids.
7. Make sure that the grinding or polishing discs are correctly mounted, the driving pins engaged with the holes in the turntable.
8. Keep clear of the rotating disc and the specimen mover during operation. While grinding manually, be careful not to touch the grinding disc. Do not attempt to collect a specimen from the tray while the disc is running.
9. If you observe malfunctions or hear unusual noises - stop the machine and call technical service.
10. The machine must be disconnected from the mains prior to any service. Wait until residual potential on the capacitors is discharged.
11. Do not cycle mains power more than once every three minutes. Damage to the drive will result.

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

TegraForce-3/ -5 Safety Precaution Sheet

To be read carefully before use

1. The operator should be fully aware of the use of the machine according to the Instruction Manual.
2. Keep your hands clear of the specimen holder or specimen mover plate when lowering TegraForce-5.
3. Do not touch the rotating parts during operation.
4. When working at machines with rotating parts care has to be taken that clothes and/or hair cannot be caught by the rotating parts. Appropriate safety clothing must be used.
5. If you observe malfunctions or hear unusual noises - stop the machine and call technical service.
6. Alcohol based consumables: follow the current safety rules for handling, mixing, filling, emptying and disposal of the alcohol-based liquids.
7. If the TegraPol with TegraForce has to be moved to a new location the TegraForce has to be detached from the TegraPol. (See section: *Installing TegraForce*)

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

Disposal

Equipment marked with a WEEE symbol  contain electrical and electronic components and must not be disposed of as general waste.

Please contact your local authorities for information on the correct method of disposal in accordance with national legislation.

Table of Contents	Page
1. Getting Started	
Checking the Contents.....	4
Getting Acquainted with TegraPol.....	5
Rear of TegraPol.....	6
Unpacking TegraPol.....	7
Placing TegraPol.....	7
Supplying Water.....	7
Connection to Water Mains.....	7
Connection to Water Outlet.....	7
Connecting a Recirculation Cooling Unit.....	8
Connecting the Water Inlet.....	8
Connecting the Water Outlet.....	8
Connecting the Communication Cable.....	8
Supplying Power.....	9
Single-phase Supply.....	9
2-phase Supply.....	9
Connection to the Machine.....	9
2. Basic Operation of TegraPol	
Front Panel.....	10
Front Panel Controls.....	11
Reading the Display.....	12
Acoustic Signals.....	12
Software Settings.....	13
Changing the Language and Display Contrast.....	13
Editing Numeric Values.....	15
Editing Alphanumeric Values.....	16
Preparation Process Setup.....	17
Mounting a Disc.....	18
Setting the Preparation Parameters.....	18
Starting the Process.....	20
Stopping the Process.....	20
Spin Function.....	20
3. Installing TegraForce-3 on TegraPol	
Checking the Packing Contents.....	21
Getting Acquainted with TegraForce-3.....	21
Setting up TegraForce-3 on TegraPol.....	22
Electrical Connections.....	23
Compressed Air Connections.....	23
Emptying the Water / Oil Filter.....	24
Removing TegraForce-3 from TegraPol.....	25

4. Operating TegraPol with TegraForce-3

Inserting Specimen Holders	26
Switching on the TegraPol with a TegraForce-3 installed	27
Setting the Preparation Parameters	27
Force Reduction	28
Pause Between Steps	29
Preparation with TegraForce-3	30

5. Installing TegraForce-5 on TegraPol

Checking the Packing Contents	31
Getting Acquainted with TegraForce-5	31
Setting up TegraForce-5 on TegraPol	32
Adjusting the Friction	33
Electrical Connections	33
Compressed Air Connections	34
Emptying the water / oil filter	34
Removing TegraForce-5 from TegraPol	35

6. Operating TegraPol with TegraForce-5

Inserting Specimen Holders	36
Inserting a Specimen Holder	36
Inserting a Specimen Mover Plate	37
Aligning the TegraForce-5	37
Switching on the TegraPol with a TegraForce-5 installed	38
Setting the Preparation Parameters	38
Manual Preparation	39
Force Reduction	40
Preparation with TegraForce-5	41
Recommendations for Grinding Single Specimens	42

7. Operating TegraPol with TegraForce and TegraDoser-1

Checking the Packing Contents	43
Getting Acquainted with TegraDoser-1	43
Filling TegraDoser-1	43
TegraDoser-1 Mounting	44
Fitting TegraDoser-1	44
Removing TegraDoser-1	44
Configuring TegraDoser-1	45
Setting Preparation Parameters	46
Running an Entire Preparation Method by Changing TegraDoser-1	47
Cleaning TegraDoser-1	49

8. Accessories	50
9. Trouble-shooting	52
10. Maintenance	
Daily Service	55
Weekly Service	55
Service Information	56
11. Technical Data	57

1. Getting Started

Checking the Contents

In the packing box you should find the following parts:

TegraPol-21 & -25


- 1 TegraPol-21 or TegraPol-25
- 1 Splash ring
- 2 Mains cables
- 1 Inlet hose 3/4" (2 m)
- 1 Filter gasket
- 1 Gasket $\varnothing 11/\varnothing 24$ (1.5 mm)
- 1 Reduction ring with gasket
- 1 Outlet hose $\varnothing 32$ (1.5 m)
- 1 88° bend outlet pipe $\varnothing 32$
- 2 Hose clamps
- 1 Set of Instruction Manuals

TegraPol-31 & -35

- 1 TegraPol-31 or TegraPol-35
- 1 Splash ring
- 2 Mains cables
- 1 Inlet hose 3/4" (2 m)
- 1 Filter gasket
- 1 Gasket $\varnothing 11/\varnothing 24$ (1.5 mm)
- 1 Reduction ring with gasket
- 1 Outlet hose $\varnothing 38$ (1.5 m)
- 1 88.5° bend outlet pipe $\varnothing 40$
- 2 Hose clamps
- 1 Set of Instruction Manuals

Getting Acquainted with TegraPol


Take a moment to familiarise yourself with the location and names of all the TegraPol components:

- A Water tap with sprinkler tube
- B Turntable
- C Front panel control(s)

TegraPol-21/ -25/ -31/ -35
TegraForce-3/ -5
Instruction Manual

Rear of TegraPol

- ① Water inlet tube
- ② Water outlet for TegraDoser-5
- ③ Water outlet tube
- ④ Service socket
- ⑤ Sample mover (TegraForce) socket
- ⑥ Doser unit (TegraDoser-5) socket
- ⑦ Auxiliary socket (not currently in use)
- ⑧ Main switch
- ⑨ Fuses
- ⑩ Mains connection
- ⑪ Communication socket (Cooli Control Box)
- Emergency Stop

Unpacking TegraPol

- Remove the TegraPol from the packing case.

Placing TegraPol

- Place TegraPol on a steady table with an adequate working height. The machine must be close to the power supply, water mains and water outlet facilities.

Supplying Water

Water for wet grinding may be supplied from the water mains or a Recirculation Cooling Unit (instructions on the next page).

Connection to Water Mains

IMPORTANT

The cold water supply must have a head pressure in the range 1 – 10 bar (14.5 – 145 psi).

- Mount the 90° end of the inlet hose onto the water inlet tube on the back of TegraPol (see [Getting Acquainted with TegraPol](#)):
 - Insert the filter gasket in the coupling nut with the flat side against the pressure hose.
 - Tighten the coupling nut completely.
- Mount the straight end of the inlet hose on the water mains tap for cold water:
 - Mount the reduction ring with gasket on the water mains tap, if necessary.
 - Introduce the gasket and tighten the coupling nut completely.

Connection to Water Outlet

- Mount the outlet hose onto the water outlet tube. (Lubricate with grease or soap to facilitate insertion.) Use a hose clamp for fastening.
If there is only limited space behind the TegraPol, use the pipe bend and attach it directly to the outlet tube on the back of the TegraPol. Then attach the drain hose to the other end of this bend.
- Lead the other end of the drain hose to the water outlet. Arrange the tube so that it slopes downward towards the drain throughout its length. Shorten the hose, if necessary.

REMEMBER

Make sure that the drain hose slopes downward towards the drain throughout its length and avoid sharp bends in the drain hose.

Connecting a Recirculation Cooling Unit

To ensure optimal cooling, TegraPol can be fitted with a Struers Recirculation Cooling Unit.

Please note:

The Recirculation Cooling Unit cannot be used with TegraPol units running methods that include the dosing of oxide polishing suspensions.

Note


Before connecting the cooling unit to the TegraPol, follow the instructions in the Struers Cooling Units Instruction Manual to prepare it for use.

Connecting the Water Inlet

- Take the inlet hose delivered with TegraPol and cut off the coupling nut on the straight end.
- Slide the hose clamp onto the hose and insert the quick coupling delivered with the pump for the recirculation cooling unit. Tighten the hose clamp.
- Connect the quick coupling on the inlet hose directly to the cooling unit's pump outlet (**A**).
If necessary, use one of the hoses delivered with the pump to extend the inlet hose.
- Mount the other end to the water inlet tube.

Connecting the Water Outlet

- Mount the outlet hose onto the water outlet tube. Use a hose clamp to secure the hose.
If there is only limited space behind the TegraPol, use the pipe bend and attach it directly to the outlet tube on the back of the TegraPol. Then attach the drain hose to the other end of this bend.
- Lead the other end into the mounting hole in the bracket on top of the static filter unit (**B**). Arrange the tube so that it slopes downward towards the drain throughout its length. Shorten the hose, if necessary.

Connecting the Communication Cable

- Connect the communication cable to the Cooli Control Box and then connect to the socket ⑪ at the rear of the TegraPol.

IMPORTANT

Before connecting, check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.

Supplying Power

Always remember to switch the power off when installing electrical equipment.

DANGER!


The machine must be earthed

IMPORTANT

Check that the mains voltage corresponds to the voltage stated on the type plate on the back of the machine.

The TegraPol is shipped with 2 types of Mains cables:

Single-phase Supply

The 2-pin (European Schuko) plug is for use on single-phase connections.

If the plug supplied on this cable is not approved in your country, then the plug must be replaced with an approved plug. The leads must be connected as follows:

Yellow/green: earth
Brown: line (live)
Blue: neutral

2-phase Supply

The 3-pin (North American NEMA) plug is for use on 2-phase power connections.

If the plug supplied on this cable is not approved in your country, then the plug must be replaced with an approved plug. The leads must be connected as follows:

Green: earth
Black: line (live)
White: line (live)

Connection to the Machine


Both cables are on the other end equipped with an IEC 320 cable connector that has to be connected to the TegraPol.

WARNING!








The output voltage from this cable is 200 – 240V and not 110V. DO NOT use this cable to connect equipment that use a 110V power supply. Failure to adhere to this may result in material damage.

2. Basic Operation of TegraPol

Front Panel

Front Panel Controls

Key	Function	Key	Function
	Moves cursor in display upwards or increases the value when editing parameter values.		Only active when TegraDoser-1 or TegraDoser-5 is fitted. Manual override – push button to apply lubricant from the doser bottle.
	Moves cursor in display downwards or decreases the value when editing parameter values.		Only active when TegraDoser-1 or TegraDoser-5 is fitted. Manual override – push button to apply diamond suspension from the doser bottle.
Esc	Leaves the present menu or aborts functions/changes.		Push button for starting the preparation process or to start Spin function for drying discs and consumables.
	Enables selected parameter values to be activated for editing. Saves the edited parameter values. Toggles when only 2 parameters available.		Push button for stopping the preparation process.
	Manual override - push button & hold to apply water (applies water during preparation processes and when no process is running).		For use with TegraForce-5. Lowers & raises TegraForce-5.
	EMERGENCY STOP - Push the red button to stop. - Turn the red button clockwise to release.		


Reading the Display

The display on the front panel provides different levels of status information. For example, when the machine is switched on using the Mains switch located at the rear of the machine, the display informs you about the physical configuration of the TegraPol and the version of software that is installed:

When operating the TegraPol, this display is the user-interface to TegraPol's software.

The display is primarily divided into 2 areas. The position of these areas and the information they contain are explained in the illustration below, which uses the CONFIGURATION Menu as an example:

- A** Heading: this is a navigational aid, telling you where you are in the software's hierarchy.
- B** Information fields: these will be either numerical values or text fields, providing information associated with the process shown in the heading. The inverted text shows the cursor position.


Acoustic Signals

When pressing a key, a short beep indicates that the command has been accepted, whereas a long beep indicates that the key cannot be activated at the moment.

Software Settings





When the machine is started for the first time, a screen display prompts you to select your preferred language. If you want to change the language after this, refer to “*Changing the Language and Display Contrast*”.


During normal operation, immediately after start up, where the splash screen is displayed, the software goes to the MAIN MENU. This is the highest level in the menu structure. From this menu, you can enter the Manual Preparation and Configuration menus.



Changing the Language and Display Contrast

Before continuing, we recommend that you select the language that best suits you (if you have not done so already).


↓
 Use the Arrow buttons  ▲▼ to select the language you prefer.

↓
 Press ENTER ↵ to accept the language.
↓
The CONFIGURATION Menu now appears in the language you have chosen.


If the Display contrast setting is OK, press Esc to return to the MAIN MENU.
To change the Display contrast setting, continue from here.

 Use the Arrow buttons  ▲▼ to select Display contrast.

↓
 Press ENTER ↵ to activate the Display contrast pop-up menu.

↓
 Use the Arrow buttons  ▲▼ to select the setting you prefer.

Contrast settings are: Default: 25
Contrast range: 0 – 50
Change increment: 1


↓
 Press ENTER ↵ to accept the setting.


↓
 Press Esc to return to the MAIN MENU.

Editing Numeric Values


Use the Arrow buttons ▲▼ to select the value to be changed, e.g. Display contrast:

Press ENTER  to edit the value.

A popup box showing the current value and up/down arrows is displayed.

Note:

If there are only two options, the popup box is not displayed. Use ENTER  to toggle between the two options and ignore the following steps.


Use the Arrow buttons ▲▼ to increase or decrease the numeric value.



Press ENTER  to accept the new value. (Pressing **Esc**, aborts the changes, preserving the original value.)

Editing Alphanumeric Values


Use the Arrow buttons   to select the text value to be changed, e.g. Language

Press ENTER  to edit the value.

A popup box showing all the options available is displayed:

Note:

If there are only two options, the popup box is not displayed. Use ENTER  to toggle between the two options and ignore the following steps.

Use the Arrow buttons   to select the correct option.

Press ENTER  to accept the new value. (Pressing **Esc**, aborts the changes, preserving the original value.)

Preparation Process Setup

TegraPol-21/ -25/ -31/ -35 are designed for the preparation of most materials. In Struers Metalog Guide™ you will find a detailed description of grinding/polishing for automated mechanical specimen preparation.

Struers Metalog Guide™ offers preparation methods for the most common materials, based on a simple analysis of two key properties: hardness and ductility. Finding the right method is easy, including choice of consumables.

Always consult Struers Metalog Guide™ on the Struers website for the correct preparation method for the actual specimens.

Metalog Guide™

A complete guide to materialographic specimen preparation.

[struers.com/KNOWLEDGE/Metalog Guide](https://www.struers.com/KNOWLEDGE/Metalog%20Guide)

Mounting a Disc

Important!
 Ensure that the back of the disc is clean and has nothing attached to it that may prevent it lying flat on the turntable.

- Place the disc on the turntable and rotate it until the 3 driving pins engage with the holes in the turntable.

Setting the Preparation Parameters


Before starting the preparation process, the disc speed and water application need to be specified. This is done using the control panel software:

↓ If the MAIN MENU is not currently displayed, press **Esc** until the MAIN MENU appears.

Use the Arrow buttons ▲ ▼ to select Manual Preparation:

↓


↓

Press ENTER ↵ to activate the PREPARATION Menu.

↓

↓ Within this menu, there are 3 preparation options.

Options	Speed(s) [rpm]	Water	Change Increment [rpm]
Grinding	300	ON	
Polishing	150	OFF	
User def for: TegraPol-21 & -31 TegraPol-25 & -35	50-300 40 – 600	ON or OFF ON or OFF	10

TegraPol-21/ -25/ -31/ -35
TegraForce-3/ -5
Instruction Manual

Use the Arrow buttons  ▲ ▼ to select the desired preparation option.

If User def is selected, use the standard editing procedures to specify the desired preparation conditions.

The TegraPol software is now ready to begin the preparation process.

Starting the Process

Once the desired process has been selected, press  to start the turntable rotating.

If water is required to an otherwise dry process, press  to override the process settings and apply water. This occurs as long as  is held.

Stopping the Process

- To stop the turntable (and water application if using a wet process), press .
- Or,
- Press the emergency stop button to stop the process. The emergency stop button must be turned to release it before the machine can be restarted. The sample mover will return to its initial position before the process can be re-started.

Spin Function

After a grinding process has been finished, the built-in Spin function can be used to rotate the preparation disc at high speed to remove water from the surface of the disc.

This function can be used to remove water from a MD-grinding disc or a SiC-paper before removing it, or to dry an MD-Disc or an MD-Chem cloth.

Note:

The Spin function can only be used when the specimen mover is moved into the top position.


- Press and hold the Start key  for at least one second to start the Spin function. (Two short beeps can be heard). If the specimen mover is in the down position the first short beep will be followed by one long beep.
- Release the Start key to stop the Spin function.

3. Installing TegraForce-3 on TegraPol

Checking the Packing Contents

- 1 TegraForce-3
- 1 Pressure hose, 1/4" (3 m)
- 1 Air hose $\varnothing 4/\varnothing 6$ (blue) (1 m)
- 1 Hose clamp
- 1 Quick coupling
- 1 Hose nipple
- 1 Allen key 2.5 mm
- 1 Allen key 4 mm
- 1 Allen key, 6 mm with handle


Getting Acquainted with TegraForce-3


- ① Handle to unlock/lock sideward movement on TegraForce-3
- ② D-Sub connector for TegraDoser-1 or nozzle block of TegraDoser-5
- ③ LED-light to illuminate preparation disc
- ④ Release outlet valve from water/oil filter for compressed air

Setting up TegraForce-3 on TegraPol

- Remove the cap covering the support hole placed on the top of the TegraPol cabinet.

- Remove the two caps covering the fastening screws on the left-hand side of the TegraPol cabinet by lifting them off with a knife.
- Lead the supporting column of TegraForce-3 down into the support hole in TegraPol.
- Turn the column until the cross pin at the end of the column engages in the slot in the casting of the TegraPol while the specimen mover is placed over the preparation disc.

- Tighten the two fastening screws on the left-hand side of the TegraPol cabinet and mount the two caps again.
- If it is not possible to swing the TegraForce-3 over the left hand side of the machine when a sample holder is mounted, loosen the two fastening screws, lift the specimen mover slightly up to allow the TegraForce to swing freely and tighten the fastening screws again.

Electrical Connections

There is only a single electrical connection between the TegraPol and the TegraForce. This provides both the 24V power supply to the TegraForce and a data bus, which enables the TegraPol and TegraForce to communicate.

- Switch off the TegraPol.
- Use the single cable attached to the TegraForce and connect it to the Sample mover connector on the rear of the TegraPol. To locate this connector, refer to the diagram for the back of the TegraPol in the section Getting Acquainted with TegraPol.

Compressed Air Connections

The quality of the air in the compressed air system used with this machine must meet the standards stated in ISO 8573-1. To connect the TegraSystem:

- Mount the quick coupling on the pressure hose and secure it with the hose clamp supplied.
- Connect the air hose to the quick coupling and fit the other end into the compressed air inlet on TegraForce-3.


IMPORTANT

The air pressure must be between 5 bar (72 psi) and 10 bar (145 psi) and have a quality specified in ISO 8573-1.

Emptying the Water / Oil Filter

TegraForce is fitted with a water / oil filter that removes excessive amounts of these substances from the compressed air supply. As a result of this, it is necessary to empty the filter periodically:

- Remove the preparation disc from the work area. Failure to do this may result in the preparation disc being contaminated by water and oil particles landing on it.
- Locate the release outlet valve on the bottom of the filter glass.

Release Outlet Valve

- Hold a container under the filter to retain any water released and press the release valve.
- Replace the preparation disc.

Removing TegraForce-3 from TegraPol

To remove TegraForce from TegraPol follow these steps:

- Remove the preparation disc.
- Switch off TegraPol with the switch at the back of the machine.
- Disconnect TegraForce from the compressed air supply.
- Remove the TegraDoser-1 or the dosing nozzle from TegraDoser-5 from the TegraForce.
- Disconnect the TegraForce from the TegraPol by removing the plug.
- Lock the handle to fix TegraForce on the column.
- Remove the two caps covering the fastening screws on the left-hand side of the TegraPol cabinet by lifting them off with a knife.
- Loosen the two fastening screws while holding the TegraForce in position.

WARNING!

TegraForce is relatively heavy. When loosening the two fastening screws it might slide downwards. Keep your hands off the column and away from between specimen mover disc and turntable.


- Lift the TegraForce straight up and place it on a safe and stable support.

4. Operating TegraPol with TegraForce-3

Inserting Specimen Holders

The TegraForce-3 can be operated with specimen holders only.

- If there isn't sufficient space between TegraForce-3 and the turntable disc to insert a specimen holder, release the black handle on the side of TegraForce-3 and swing it to the side until it is clear of the TegraPol.
- Push the black button on the TegraForce-3.

- Insert the specimen holder and rotate it until the three pins are aligned and then push the holder upwards until it locks in position.
- Release the black button on the TegraForce-3 and ensure that the specimen holder is securely fixed.
- If TegraForce-3 was swung to the side, swing it back and align it correctly over the preparation disc. The specimens should come close to the edge of the preparation disc. They may extend over the edge by 2 - 3 mm.

WARNING!

When working with specimen holders make sure that the screws clamping the specimens do not stick out of the specimen holder.
Use different length of screws for samples with different diameters.

Switching on the TegraPol with a TegraForce-3 installed

When a TegraPol has been fitted with a TegraForce-3 and switched on, the software inside the TegraPol is able to detect the TegraForce-3, automatically, and informs you of its presence in the start-up display.

Setting the Preparation Parameters

Before starting the preparation process, the disc speed, preparation time, force and water application need to be specified. This is done using the control panel software:


Note!
If a TegraDoser-5 is mounted on the TegraForce and connected to the TegraPol, then refer to the TegraDoser-5 Instruction Manual for the correct procedure for setting the preparation parameters.

↓
If the MAIN MENU is not currently displayed, press **Esc** until the MAIN MENU appears.

Use the Arrow buttons  ▲ ▼ to select Preparation:

↓

After a TegraForce has been mounted and connected to the TegraPol there is one more topic on the MAIN MENU, Preparation. Preparation is selected when the TegraForce is to be used together with the TegraPol.

If a sample has to be prepared manually, Manual preparation is selected and the same possibilities as described earlier are available.

Force Reduction

As soon as the TegraForce is connected, some additional topics are available in the CONFIGURATION Menu. Therefore the Configuration Menu is also divided into two parts, General Configuration and Machine Configuration. Force Reduction is used to decrease the force gradually during the last 30 seconds of the preparation time. The force is reduced in steps of 5 N until the final force of 10 N is reached. Thus the surface finish is finer and the preparation time in the following step can be reduced. Force Reduction can be set to either On or Off and this setting is valid for all the preparation steps carried out.

↓
 Press ENTER ↵ to activate the CONFIGURATION Menu.

```
CONFIGURATION
General Configuration
Machine Configuration
```

 Use the Arrow buttons ▲▼ to select Machine Configuration.

 Press ENTER ↵ to activate the Machine Configuration pop-up menu.

```
MACHINE CONFIG.
Pause between steps      No
Force Reduction:        ON
Sample holder cleaning  Press
                        ⏏
```

 Use the Arrow buttons ▲▼ to select Force Reduction and press ENTER ↵ to toggle between ON and OFF.

Pause Between Steps

This function can only be used when a TegraDoser-5 is connected. When two following preparation steps are based on the same surface, suspension and lubricant, these two steps will be joined so that the preparation continues directly from one step to the next. The advantage is that force and dosing level can be changed from one step to the next without having to restart the machine. However, it may be desirable to stop the machine, e.g. when using two subsequent steps on SiC with the same grit size. Therefore it is possible to set the *Pause between steps* function to Yes.

↓
 Press ENTER ↵ to activate the CONFIGURATION Menu.


↓

```
CONFIGURATION
General Configuration
Machine Configuration
```

↓

 Use the Arrow buttons   to select Machine Configuration.

↓

 Press ENTER ↵ to activate the Machine Configuration pop-up menu.

↓

```
MACHINE CONFIG.
Pause between steps      No
Force Reduction:        ON
Sample holder cleaning  Press
                        
```

↓

 Use the Arrow buttons   to select Pause between steps and press ENTER ↵ to toggle between Yes and No

Preparation with TegraForce-3

Press ENTER ↵ to activate the SAMPLE HOLDER Menu.

Within this menu, there are 4 preparation parameters.

Options	Settings for		Change Increment
	TegraPol-21 & 31	TegraPol-25 & 35	
Speed Disc Holder (TegraForce)	50 – 300	40 - 600	10 rpm
Time	50 – 150		
	60m 00s 05m 00s 01m 00s 00m 05s		30s 10s 5s
Force	30 - 400		10 N
Water	ON or OFF		

Use the Arrow buttons ▲▼ to select the desired preparation option.

The TegraPol software is now ready to begin the preparation process.

Important!


The TegraForce must be locked in position before the preparation process can be started. An error beep is sounded and an error message is displayed if it is not locked.

5. Installing TegraForce-5 on TegraPol

Checking the Packing Contents

- 1 TegraForce-5
- 1 Pressure hose, 1/4" (3 m)
- 1 Air hose $\varnothing 4/\varnothing 6$ (blue) (1 m)
- 1 Hose clamp
- 1 Quick coupling
- 1 Hose nipple
- 1 Allen key 2.5 mm
- 1 Allen key 4 mm
- 1 Allen key, 6 mm with handle


Getting Acquainted with TegraForce-5


- ① D-Sub connector for TegraDoser-1 or nozzle block of TegraDoser-5
- ② Release outlet valve from water/oil filter for compressed air
- ③ LED-light to illuminate preparation disc

Setting up TegraForce-5 on TegraPol

- Remove the cap covering the support hole placed on the top of the TegraPol cabinet.


- Remove the two caps covering the fastening screws on the left-hand side of the TegraPol cabinet by lifting them off with a knife.
- Lead the supporting column of TegraForce-5 down into the support hole in TegraPol.

- Turn the column until the cross pin at the end of the column engages in the slot in the casting of the TegraPol while the specimen mover is placed over the preparation disc.
- Tighten the two fastening screws on the left-hand side of the TegraPol cabinet and mount the two caps again.
- Connect the TegraForce-5 to the back of the TegraPol and to compressed air.
- Switch the TegraPol on and open for the compressed air.

- Mount a preparation disc on the TegraPol and insert a specimen mover plate (for single specimens) into TegraForce-5.
- Press the up/down key on the TegraPol to move the TegraForce in the down position.
- Adjust the distance between the disc and specimen mover plate to about 1 - 1.5 mm.

- To adjust the distance loosen the screw in the friction ring with an Allen key and move the ring up or down, depending on the direction of adjustment.
- When the distance is correct tighten the screw again to clamp the friction ring in position.

Adjusting the Friction

To adjust the friction of the sideways movement of TegraForce-5 loosen the screw in the friction ring and turn the lower part of the ring counter clockwise to increase the friction and clockwise to reduce the friction. Make sure that you do not move the friction ring up or down during that operation to avoid changing the distance between disc and specimen mover plate.

Electrical Connections

There is a single electrical connection between the TegraPol and the TegraForce-5. This provides both the 24V power supply to the TegraForce-5 and a data bus, which enables the TegraPol and TegraForce to communicate.

- Switch off the TegraPol.
- Use the single cable attached to the TegraForce-5 and connect it to the Sample mover connector on the rear of the TegraPol. To locate this connector, refer to the diagram for the back of the TegraPol in the section Getting Acquainted with TegraPol.

Compressed Air Connections

The quality of the air in the compressed air system used with this machine must meet the standards stated in ISO 8573-1. To connect the TegraForce:

- Mount the quick coupling on the pressure hose and secure it with the hose clamp supplied.
- Connect the air hose to the quick coupling and fit the other end into the compressed air inlet on TegraForce-5.


IMPORTANT

The air pressure must be between 5 bar (72 psi) and 10 bar (145 psi) and have a quality specified in ISO 8573-1.

Emptying the water / oil filter

TegraForce is fitted with a water / oil filter that removes excessive amounts of these substances from the compressed air supply. As a result of this, it is necessary to empty the filter periodically:

- Remove the preparation disc from the work area. Failure to do this may result in the preparation disc being contaminated by water and oil particles landing on it.
- Locate the release outlet valve on the bottom of the filter glass.

Release Outlet Valve

- Hold a container under the filter to retain any water released and press the release valve.
- Replace the preparation disc.

Removing TegraForce-5 from TegraPol

To remove TegraForce from TegraPol follow these steps:

- Remove the preparation disc and specimen mover disc if mounted.
- Switch off TegraPol with the switch at the back of the machine.
- Disconnect TegraForce from the compressed air supply.
- Remove the TegraDoser-1 or the dosing nozzle from TegraDoser-5 from the TegraForce.
- Disconnect the TegraForce from the TegraPol by removing the plug.
- Remove the two caps covering the fastening screws on the left-hand side of the TegraPol cabinet by lifting them off with a knife.
- Loosen the two fastening screws while holding the TegraForce in position.

WARNING!

TegraForce is relatively heavy. When loosening the two fastening screws it might slide downwards. Keep your hands off the column and away from between specimen mover disc and turntable.

- Lift the TegraForce straight up and place it on a safe and stable support.


6. Operating TegraPol with TegraForce-5

Inserting Specimen Holders

The TegraForce-5 can be operated with either specimen mover plates for single samples or specimen holders.

Inserting a Specimen Holder

- Press the Raise/Lower ▲▼ button to ensure that TegraForce-5 is fully raised.
- Push the black button on the TegraForce-5.


- Insert the specimen holder and rotate it until the three pins are aligned and then push the holder upwards until it locks in position.
- Release the black button on the TegraForce-5 and ensure that the specimen holder is securely fixed.

WARNING!

When working with specimen holders make sure that the screws clamping the specimens do not stick out of the specimen holder.
Use different length of screws for samples with different diameters.

Inserting a Specimen Mover Plate

- Use the Raise/Lower ▲▼ button to ensure that TegraForce-5 is fully raised.
- Push the black button on the TegraForce-5.

- Insert the specimen mover plate and rotate it until the three pins are aligned and then push the specimen mover plate upwards until it locks in position.
- Release the black button on the TegraForce-5 and ensure that the specimen mover plate is securely fixed.

Aligning the TegraForce-5

TegraForce-5 is equipped with an automatic alignment mechanism taking care of the fine alignment of the head relative to the preparation disc. It ensures that when the TegraForce-5 head is lowered, it returns to its previous position. To change the alignment of the head:

- With a specimen mover plate or specimen holder fitted, press the Raise/Lower button to lower the head.
- Push the entire head sideways (left or right) until it is aligned correctly. The specimens should be positioned as close to the edge of the preparation disc as possible, in fact they may extend over the edge by 2-3 mm.
- To adjust the friction of the sideways movement see: Adjusting the friction in the section: Setting up TegraForce-5 on TegraPol

IMPORTANT!

The automatic alignment is only a fine alignment of the TegraForce-5. It has to be positioned in approximately the correct position, max. $\pm 20^\circ$ from the correct position.

Switching on the TegraPol with a TegraForce-5 installed

When a TegraPol has been fitted with a TegraForce-5 and switched on, the software inside the TegraPol is able to detect the TegraForce-5, automatically, and informs you of its presence in the start-up display.


Setting the Preparation Parameters

Before starting the preparation process, the disc speed, preparation time, force, sample holder/mover direction, and water application need to be specified. This is done using the control panel software:


Note!

If a TegraDoser-5 is mounted on the TegraForce and connected to the TegraPol, then refer to the TegraDoser-5 Instruction Manual for the correct procedure for setting the preparation parameters.

↓ If the MAIN MENU is not currently displayed, press **Esc** until the MAIN MENU appears.

 Use the Arrow buttons  ▲ ▼ to select Single Sample Preparation or Sample Holder Preparation:

↓


After a TegraForce-5 has been mounted and connected to the TegraPol there are two more topics on the Main Menu, *Single Samples Preparation* and *Sample Holder Preparation*. One of these two choices is selected when the TegraForce is to be used together with the TegraPol.

Manual Preparation

If a specimen cannot be prepared using a standard specimen mover plate or specimen holder, it can be prepared manually, with the automatic dosing.

- Select the required preparation step and make sure that the specimen mover is in the upright position.
- Press Start.


- Press Enter to confirm selection of a Manual preparation step


- Press Start again to start the preparation step.
Dosing of suspension and lubricant runs automatically according to the pre-set parameters and the step is stopped automatically when the set preparation time has expired.

Force Reduction

As soon as the TegraForce is connected, an additional topic is available in the Configuration Menu. Force Reduction is used to decrease the force gradually during the last 30 seconds of the preparation time. The force is reduced in steps of 5 N until the final force of 10 N is reached. Thus the surface finish is finer and the preparation time in the following step can be reduced. Force Reduction can be set to either On or Off and this setting is valid for all the preparation steps carried out.

Press ENTER ↵ to activate the CONFIGURATION Menu.

Select Force Reduction and press ENTER ↵ to toggle between ON and OFF.

Preparation with TegraForce-5

Press ENTER ↵ to activate the Preparation selection:
SINGLE SAMPLES or SAMPLE HOLDER:

SINGLE SAMPLES	SAMPLE HOLDER
RPM Disc/Holder: 300/150	RPM Disc/Holder: 300/150
Time: 1 m 00s	Time: 1 m 00s
Force/Dir.: 30N /	Force/Dir.: 180N /
Water: <input checked="" type="checkbox"/> ON	Water: <input checked="" type="checkbox"/> ON

or

Within this menu, there are 4 preparation parameters.

Settings for SINGLE SAMPLES

Options	TegraPol-21 & 31	TegraPol-25 & 35	Change Increment
Speed * Disc	50 - 300	40 - 600	10 rpm
Holder (TegraForce)	50 - 150		10 rpm
Time	60m 00s 05m 00s 01m 00s 00m 05s		30s 10s 5s
Force/Dir.:	10 – 100/cw or ccw		5 N
Water	ON or OFF		

* Please see [Recommendations for Grinding Single Specimens](#) on page 42.

Settings for SAMPLE HOLDER

Options	TegraPol-21 & 31	TegraPol-25 & 35	Change Increment
Speed Disc	50 - 300	40 - 600	10 rpm
Holder (TegraForce)	50 - 150		10 rpm
Time	60m 00s 05m 00s 01m 00s 00m 05s		30s 10s 5s
Force/Dir.:	30 - 400/cw or ccw		10 N
Water	ON or OFF		

Use the Arrow buttons   to select the desired preparation option.

The TegraPol software is now ready to begin the preparation process.

Important!

Press the Raise/Lower button to lower the TegraForce-5 head. Otherwise manual preparation can be started or the Spin function can be activated.

***Recommendations for Grinding
Single Specimens***

Do not use plane grinding with coarse abrasives when preparing single specimens. It is normally not necessary, and the use of coarse abrasives can result in unplane specimens.

If, for whatever reason, it is necessary to grind using coarse abrasive, the planeness may be improved by carrying out the following:


- The height of the specimen should be between 8 – 40 mm and not exceed 0.7 x specimen diameter.
Example: A specimen with a diameter of 30 mm should not be higher than $30 \times 0.7 = 21$ mm.
- Use as small a grain size as possible.
- Use a mounting resin with a wear resistance similar to the specimens wear resistance.
- Use 150 rpm for both grinding disc and specimen mover.
(when using lower speeds decrease the speed on both the disc and the specimen mover).
- Use counter rotation on TegraForce-5.
- Use low forces.
- Position TegraForce-5 so that the specimens do *not* go over the centre of the preparation disc.

7. Operating TegraPol with TegraForce and TegraDoser-1

Checking the Packing Contents

- 1 TegraDoser-1
- 1 Funnel
- 1 Flat blade screwdriver
- 1 Sheet with labels

Getting Acquainted with TegraDoser-1

Filling TegraDoser-1

The TegraDoser bottle has a snap tight lid, which provides easy access when filling the bottle.

- Ensure that no preparation process is active.
- Open the snap lid on the top of the bottle and insert the funnel supplied with the TegraDoser.
- Fill/refill the bottle with lubricant / suspension and snap the top closed.

Important!

TegraDoser-1 can only be used with lubricants and diamond suspensions. It is perfectly suited for use with the new all-in-one DiaPro Suspensions allowing unattended preparation. It cannot be used with oxide polishing suspensions. Oxide polishing suspensions are to be used on TegraDoser-5.


TegraDoser-1 Mounting

TegraDoser is mounted on the side of the TegraForce. The interface between the TegraForce and TegraDoser is achieved through the D-Sub connector located in the recess on the side of the TegraDoser.

Fitting TegraDoser-1

With the bottle containing the required suspension or lubricant:

- Swing the dosing nozzle on TegraDoser aside to gain access to the connector.

- Insert TegraDoser into the recess on the side of TegraForce.
- Press downwards until TegraDoser clicks into position.
- Swing the dosing nozzle back into position under the body of TegraDoser.

Removing TegraDoser-1


- Swing the dosing nozzle outwards away from TegraForce until it is clear of the connector.

- Grip the sides of the bottle and pull upwards until the TegraDoser comes clear of TegraForce.

Configuring TegraDoser-1

Each TegraDoser can have a unique ID. This is achieved using the small rotary switch on the side of TegraDoser-1.

- Use the small flat blade screwdriver included with the TegraDoser-1 to rotate the arrow to the desired ID number.
- Mark the bottle using the labels supplied with the TegraDoser-1. Show clearly the ID address of the TegraDoser and the type of suspension or lubricant present in the bottle.

Important!

The ID number of any TegraDoser should be unique. None of the other TegraDosers used on this machine should have the same ID selected as the complete set of preparation parameters of one preparation step can be saved together with the bottle ID.

If the new DiaPro suspensions combining both diamond suspension and lubricant are used, only one liquid has to be applied during each preparation step. This is the ideal application for TegraDoser-1.

By filling several TegraDosers with unique ID's with a different suspension each, a complete preparation method can be carried out without reprogramming of the equipment.

As soon as one TegraDoser-1 is replaced with another, the preparation parameters automatically change to the parameters used last time the TegraDoser-1 with this ID was used. Thus an entire preparation method can be carried out just by exchanging TegraDosers as the correct parameters are set automatically.

Setting Preparation Parameters When a TegraDoser-1 has been fitted to the TegraPol / TegraForce system, the dosing level must be set. This is done using the control panel software:

Note:
 The dosing parameters are only displayed when the water parameter is set to OFF.


↓ With the SINGLE SAMPLES menu already displayed:

SINGLE SAMPLES

RPM Disc/Holder: 150/100
 Time: 1m00s
 Force/Dir.: 30N / 
 Water: OFF
 Doser No./Level: 1 / 

↓ Within this menu, there are 2 new preparation parameters, Doser No. and Level

Options	Settings	Comment/Change Increment
Water	ON or OFF	Must be set to OFF to use dosing feature (is automatically set to OFF when a TegraDoser-1 is inserted)
Doser No.	0 – 9	<i>Information only.</i> indicates ID setting of rotary switch on TegraDoser
Level	0 – 20	1

 Use the Arrow buttons ▲▼ to select the desired dosing level and press ENTER ↵ to confirm.

The TegraPol software is now ready to begin the preparation process.

Note

All preparation parameters set in this menu are remembered by the firmware in the TegraPol together with the bottle ID of the TegraDoser. When the TegraDoser is used on another machine, a different set of parameters can be set as it is the TegraPol that records the information, not the TegraDoser.

Running an Entire Preparation Method by Changing TegraDoser-1

As explained above the setting of the ID on the TegraDoser can be used to run an entire preparation method. This is shown in the following example:

Grinding

- ↓ The first preparation step usually is a grinding step where no abrasive is used but only water is applied as a lubricant. Here either the TegraDoser-1 is removed from the TegraForce or the parameter Water is set to: ON. This also changes the parameters Speed, Time and Force to what was previously used in a grinding step.

SINGLE SAMPLES	
RPM Disc/Holder:	300/150
Time:	1 m 00 s
Force/Dir.:	30N / 
Water:	<input checked="" type="checkbox"/> ON

Fine grinding

- ↓ The second step is the fine grinding step, here a TegraDoser-1 with the ID 1 is used, containing the DiaPro suspension for fine grinding on MD-Allegro/Largo.

SINGLE SAMPLES	
RPM Disc/Holder:	150/100
Time:	1 m 00 s
Force/Dir.:	30N / 
Water:	OFF
Doser No./Level:	1 / 

Polishing

- ↓ The third step is the first polishing step, here a TegraDoser-1 with the ID 2 is used, containing a DiaPro suspension for polishing.

SINGLE SAMPLES	
RPM Disc/Holder:	150/100
Time:	4 m 00 s
Force/Dir.:	20N / 
Water:	OFF
Doser No./Level:	2 / 

Final polishing

- ↓ The last step is the final polishing step, here a TegraDoser-1 with the ID 3 is used, containing a DiaPro suspension for final polishing.

SINGLE SAMPLES	
RPM Disc/Holder:	150/100
Time:	2m 30s
Force/Dir.:	20N / 
Water:	OFF
Doser No./Level:	3/10

Simply by changing from one TegraDoser-1 to the next with a unique ID setting all the preparation parameters have been changed automatically.

Cleaning TegraDoser-1

Under normal operating conditions, it is recommended that the contents of the TegraDoser bottles are not changed. It is better to provide one doser per lubricant / suspension. However, if for any reason the contents of a bottle must be changed or will not be in use for a longer period of time, the following procedure can be used to clean the TegraDoser:

- Ensure that no preparation process is active.
- Remove the TegraDoser from its mounting on the TegraForce.
- Hold the TegraDoser vertically and open the snap lid on the top of the bottle.
- Empty the bottle of any remaining suspension or lubricant through the hole in the top.
- Fill the bottle with lukewarm water.
- Mount the TegraDoser on the TegraForce.
- Position a bottle under the outlet point of the TegraDoser. Press and hold  or  until the water has run through the nozzle for approximately 1 min.
- Remove the TegraDoser from the TegraForce and empty out the water remaining in the bottle.
- Now repeat the rinsing process with the bottle fitted on the TegraForce using clean tap water.
- Remove the TegraDoser from the TegraForce and empty out the clean water remaining in the bottle.
- Fit the empty TegraDoser to the TegraForce. Press and hold  or  to empty the remaining water present in the nozzle.
- The TegraDoser is now ready to be filled with a new suspension.

8. Accessories

TegraPol-21/ -25

Specification	Cat. No:
<i>Polishing Discs</i> 250 mm dia., aluminium Disc for magnetic fixation, 250 mm dia.	02426907 02426919
Lid for TegraPol-21, -25	05536901
<i>Struers Cooling Unit, System 3</i> With 50 l tank, small pump, Cooli-1, 50 l static filter with filter paper. 1 x 100 V / 50 Hz 1 x 220-240 V / 50 Hz 1 x 100-120 V / 60 Hz 1 x 220-240 V / 60 Hz 1 x 100-120 V / 50-60 Hz CSA	05766516 05766522 05766523 05766524 05766616

TegraPol-31/ -35

Specification	Cat. No:
<i>Polishing Discs</i> 300 mm dia., aluminium Disc for magnetic fixation, 300 mm dia.	02426906 02426918
Lid for TegraPol-31, -35	05556901
<i>Struers Cooling Unit, System 3</i> With 50 l tank, small pump, Cooli-1, 50 l static filter with filter paper. 1 x 100 V / 50 Hz 1 x 220-240 V / 50 Hz 1 x 100-120 V / 60 Hz 1 x 220-240 V / 60 Hz 1 x 100-120 V / 50-60 Hz CSA	05766516 05766522 05766523 05766524 05766616

TegraPol-21/ -25/ -31/ -35
TegraForce-3/ -5
Instruction Manual

TegraForce-3/ -5

Specification	Cat. No:
<i>Levelling device</i> for levelling of specimens in specimen holders (\varnothing 140 mm)	05256903
<i>Levelling device, Uniforce</i> for levelling of specimens in specimen holders \varnothing 140, 160 and 200 mm. With pressure foot for positioning of specimens	04886102
for levelling of specimens in specimen holders \varnothing 140, 160 and 200 mm.	04886101
<i>Dosing unit TegraDoser-1</i> For water based liquids For alcoholic liquids	05616104 05616204
<i>Specimen mover plates</i> For single specimens. For TegraForce-5. 4 mm thick. Coupling (04796902) is ordered separately. Without holes. For 6 specimens 25 mm dia. For 6 specimens 1" dia. For 6 specimens 30 mm dia. For 6 specimens 1.25" dia. For 6 specimens 1.5" dia. For 6 specimens 40 mm dia. For 3 specimens 50 mm dia. For Accustop 30	04796924 04796925 04796926 04796927 04796928 04796929 04796930 04796931 04796932

Specimen Holders

Please refer to the [Struers Specimen Holders brochure](#) for details of the range available.

9. Trouble-shooting

Error message	Explanation	Action required
TegraPol		
NVRAM initialised	Default settings have been set for all non-volatile parameters.	None, message will be removed after a few seconds.
Frequency inverter communication error		Turn device off and on, and try again. If problem persists call a Struers service technician.
Frequency inverter heat sink temp. too high		Turn device off and let it cool down for a few minutes. Turn device on and try again. If problem persists call a Struers service technician.
Frequency inverter internal temp. too high		Turn device off and let it cool down for a few minutes. Turn device on and try again. If problem persists call a Struers service technician.
Frequency inverter unknown error: 1	Internal frequency inverter error.	Turn device off and on, and try again. If problem persists call a Struers service technician.
Frequency inverter type mismatch error	Fabrication error, incorrect motor controller mounted in device.	Call a Struers service technician.
Motor temperature too high		Stop process and let device cool down for a few minutes. Try again. If problem persists call a Struers service technician.
Motor temp. sensor disconnected		Call a Struers service technician.
Cannot read keyboard		Call a Struers service technician.

TegraPol-21/ -25/ -31/ -35
TegraForce-3/ -5
Instruction Manual

Error message	Explanation	Action required
TegraPol with TegraForce		
TegraForce must be down for preparation		Press up/down button before starting process.
TegraForce not locked		Lock TegraForce with handle before starting process.
TegraForce must be up for manual preparation		Press up/down button before starting manual process.
TegraForce must be up when cleaning		Press up/down button before starting cleaning.
Check water expeller on TegraForce		Press the release outlet valve to drain the water/oil filter.
TegraPol with TegraForce and TegraDoser-5		
TegraDoser-5 nozzle block not mounted	Nozzle block missing, or not correctly seated.	Make sure nozzle block is mounted securely.
Method not valid anymore	Method has been edited on TegraDoser-5.	Re-read method by going back to method list and select method again.
Tube cleaning in progress		Wait for cleaning process to finish.
Suspension not applicable	The suspension specified in the method has not been configured on the TegraDoser-5.	Configure the suspension in one of the pump positions, or select another suspension in the method.
Lubricant not applicable	The lubricant specified in the method has not been configured on the TegraDoser-5.	Configure the lubricant in one of the pump positions, or select another lubricant in the method.

Machine problems

	Explanation	Action required
Noise when the machine starts or the machine will not pull.	The belt is not tight enough.	Call a Struers service technician. The belt must be tightened.
Functions are not carried out.	Fuse at the rear of TegraPol blown.	Replace the fuse.
Machine not operating		
Water not draining away.	Drain hose squeezed.	Straighten the hose.
	Drain hose clogged.	Clean the hose.
	Drain hose does not slope downward.	Adjust the hose to an even slope.
Cooling water stops.	Wrong software setting.	Check software setting.
	Water tap on mains closed.	Open for water.
	Built-in water tap closed.	Open for water.
	Built-in water tap blocked.	Clean water tap.
	Filter at the water inlet blocked.	Clean filter.
Cooling water drips after stop.	Defect in the solenoid valve.	Call a Struers service technician. The solenoid valve must be replaced.
The disc vibrates.	Dirt on the underside of the disc or on the turntable.	Clean the contact face between the disc and turntable.
Continuous, irregular wear on a grinding/polishing surface.	Worn coupling on either specimen holder/mover plate or the specimen mover head of the TegraForce.	Contact a Struers Service Technician to check the couplings.

10. Maintenance

Daily Service


- Clean all accessible surfaces with a moist cloth.
- Check and refill TegraDoser-1.

Weekly Service

- Remove the grinding/polishing disc and remove all dirt from the tray and the drain tube.
- Clean the pressure feet and pistons applying the force on the specimens and specimen holder.
- Press the release outlet valve on TegraForce to drain the water/oil filter (please refer to the section on [Emptying the Water / Oil Filter](#)).

From the Main Menu, select Configuration and press ENTER ↵ to activate the CONFIGURATION Menu.

Select Cleaning and press START ↵ to move the pressure feet or the centre piston down.


Clean the pressure feet and the pistons or the centre piston using a dry cloth.

Press STOP ⏏ to move the pressure feet or the centre piston up again.


Service Information

Struers recommends that a regular service check be carried out after every 1500 hours of use. Information on total operation time and servicing of the machine is displayed on the screen at start-up:

A pop-up message will appear after 1,000 hours operation time to remind the user that a service check should be scheduled.

After the 1,500 hours operation time has been exceeded the pop-up message will change to alert the user that the recommended service interval has been exceeded.

- Contact a Struers Service Technician to service the machine.

11. Technical Data

Subject		Specifications	
		TegraPol-21	TegraPol-25
Disc	Diameter	250 mm / 10"	
	Speed	50 - 300 rpm variable	40-600 rpm variable
	Rotational direction	Anti-clockwise	
	No. of motors	1	
	Motor	370 W / 0.50 HP	550 W / 0.75 HP
	Torque at disc		
	Cont. at <300 rpm	11.8 Nm / 8.7 ft-lbf	17.5 Nm / 12.9 ft-lbf
	Cont. at 600 rpm	-	8.7 Nm / 6.4 ft-lbf
	Max.	> 20 Nm / 14.8 ft-lbf	> 30 Nm / 22.2 ft-lbf
Connected equipment	TegraForce-3/-5	1	
	TegraDoser-5	1	
Software and electronics	Controls	Touch pads	
	Memory	FLASH-ROM / RAM / NV-RAM	
	LC display	128x64 dots with back light	
Safety		Please refer to the Declaration of Conformity	
Noise level	At idle running, at a distance 1.0 m / 39.4" from the machine	54 dBA	54 dBA
	During preparation	62 dBA	62 dBA
Surrounding temperature		5-40°C / 41-104°F	
Humidity	Non condensing	0-95% RH	

TegraPol-21/ -25/ -31/ -35
TegraForce-3/ -5
Instruction Manual

Subject		Specifications	
		TegraPol-21	TegraPol-25
Supply	Voltage / frequency	200-240 V / 50-60 Hz	
	Power inlet	1-phase (N+L1+PE) or 2-phase (L1+L2+PE) The electrical installation must comply with "Installation Category II".	
	Power, nominal load	440 W	700 W
	Power, idle	16 W	16 W
	Current, nom.	2.2 A	3.5 A
	Current, max.	5.6 A	6.9 A
	Pressure for tap water	1-10 bar / 14.5-145 psi	
	Water inlet	1/2" or 3/4"	
	Water outlet	ø32 mm / 1 1/4"	
Dimensions and weight	Width	430 mm / 16.9"	
	Depth	715 mm / 28.1"	
	Height	325 mm / 12.8"	
	Weight	33 kg / 72.8 lbs	35 kg / 77.2 lbs

TegraPol-21/ -25/ -31/ -35
TegraForce-3/ -5
Instruction Manual

Subject		Specifications	
		TegraPol-31	TegraPol-35
Disc	Diameter	300 mm / 12"	
	Speed	50 - 300 rpm variable	40-600 rpm variable
	Rotational direction	Anti-clockwise	
	No. of motors	1	
	Motor	550 W / 0.75 HP	750 W / 1.0 HP
	Torque at disc		
Cont. at <300 rpm	17.5 Nm / 12.9 ft-lbf	23.8 Nm / 17.6 ft-lbf	
Cont. at 600 rpm	-	11.9 Nm / 8.8 ft-lbf	
Max.	> 30 Nm / 22.2 ft-lbf	> 40 Nm / 29.6 ft-lbf	
Connected equipment	TegraForce-3/-5	1	
	TegraDoser-5	1	
Software and electronics	Controls	Touch pads	
	Memory	FLASH-ROM / RAM / NV-RAM	
	LC display	128x64 dots with back light	
Safety		Please refer to the Declaration of Conformity	
Noise level	At idle running, at a distance 1.0 m / 39.4" from the machine	56 dBA	56 dBA
	During preparation	60 dBA	60 dBA
Surrounding temperature		5-40°C / 41-104°F	
Humidity	Non condensing	0-95% RH	

TegraPol-21/ -25/ -31/ -35
TegraForce-3/ -5
Instruction Manual

Subject		Specifications	
		TegraPol-31	TegraPol-35
Supply	Voltage / frequency	200-240 V / 50-60 Hz	
	Power inlet	1-phase (N+L1+PE) or 2-phase (L1+L2+PE) The electrical installation must comply with "Installation Category II".	
	Power, nominal load	700 W	920 W
	Power, idle	16 W	16 W
	Current, nom.	3.5 A	4.6 A
	Current, max.	6.9 A	8.0 A
	Pressure for tap water	1-10 bar / 14.5-145 psi	
	Water inlet	1/2" or 3/4"	
	Water outlet	ø32 mm / 1 1/4"	
Dimensions and weight	Width	500 mm / 19.7"	
	Depth	725 mm / 29.5"	
	Height	330 mm / 13.0"	
	Weight	43 kg / 94.8 lbs	45 kg / 99.2 lbs

TegraPol-21/ -25/ -31/ -35
TegraForce-3/ -5
Instruction Manual

Subject		Specifications	
		TegraForce-3	TegraForce-5
Specimen mover	Speed	50 - 150 rpm variable	
	Rotational direction	C.W. + C.C.W	C.W. + C.C.W
	Motor	115 W	
	Torque max.	9.0 Nm	
Connected equipment	TegraDoser 1		
Safety		Please refer to the Declaration of Conformity	
Noise level	At idle running, at a distance 1.0 m / 39.4" from the machine	56 dB (A)	56 dB (A)
	During preparation	62 dB (A)	62 dB (A)
Surrounding temperature		5-40°C / 41-104°F	
Humidity	Non condensing	0-95% RH	
Supply	Voltage	24 Vdc from TegraPol	
	Current, I_{nom} / I_{max}	2.0 A / 3.0 A	
	Air inlet	ø 5 mm + ø 8 mm	
	Air pressure	5-10 bar / 72-145 psi	
	Air quality	The air supplied must be of the quality 5.6.4, specified in the standard ISO 8573-1.	
Dimensions and weight	Width	215 mm / 8.5 "	180 mm / 7.1"
	Depth	388 mm / 15.3"	
	Height	572 mm / 22.5 "	
	Weight	17 kg / 37.5 lbs.	18 kg / 39.7 lbs.

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
Spare Parts and Diagrams

TegraPol-21/-25/-31/-35
TegraForce-3/-5
Spare Parts and Diagrams

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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TegraPol-21/-25

Drawings

TegraPol-21/-25, complete.....	15530001B
Bottom, complete	15530010O
Cabinet, assembled	15550020M
Motor complete with sledge	15550030G
Bearing construction, assembled	15530040F
Disc assembled with QUAD-rings	15530050C
Front plate, assembled.....	15550070D
Water cock, assembled.....	14920045B

Diagrams

Circuit Diagram	15533100G
-----------------------	-----------

TegraPol-31/-35

Drawings

TegraPol-31/-35, complete.....	15550001B
Bottom, complete	15550010P
Cabinet, assembled	15550020M
Motor complete with sledge	15550030G
Bearing construction, assembled	15550040D
Disc assembled with QUAD-rings	15550050B
Front plate, assembled.....	15550070D
Water cock, assembled.....	14920045B

Diagrams

Circuit Diagram	15553100F
-----------------------	-----------

TegraDoser-1

Drawing

TegraDoser-1, complete	15610001D
Cover + bottle, complete without print.....	15610030H

Diagram

Circuit Diagram	15613100A
-----------------------	-----------

TegraForce-3

Drawings

TegraForce-3, complete	15580001H
Chassis, complete	15580005L
Head, complete	15580010C
Coupling, complete	15580015Ò
Pressure control valve, complete	15580030G

Diagrams

Circuit Diagram, TegraForce-1/-3/-5	15573100D
---	-----------

TegraForce-5

Drawings

TegraForce-5, complete	15590001G
Chassis, complete	15590005S
Head, complete	15590070E

Some of the drawings may contain position numbers
not used in connection with this manual.

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











Ver.C: (FTH 04-07-2002) - Color codes corrected, - Brake resistor added Ver.D: (FTH 09-02-2009) - Working Light added - Emergency Stop functionality added - TF-5 now on a separate diagram		STRUERS A/S Pedersstrupvej 84 DK-2750 Ballerup Danmark +45 44 600 800	
TegraForce-1/3, Circuit Diagram			
Size A2	CAGE Code <Cage Code>	DWG NO	Rev 1
Monday, February 09, 2009	Scale	FTH / FTH	Sheet 15573100 1 of 1

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
is in conformity with the provisions of the following directives:

Safety of Machinery 2006/42/EC according to the following standard(s):
EN ISO 12100-1:2005, EN ISO 12100-2:2005, EN ISO 14121-1:2007, EN 60204-1:2006/A1:2009,
EN 61010:2010.

EMC-Directive 2004/108/EC according to the following standard(s):
EN 61000-6-4:2007, EN 61000-6-1:2007, EN 61000-6-2:2005, EN 61000-6-3:2007,
EN 61000-6-4:2002, EN 61326-1:2006.

Low Voltage Directive 2006/95/EC according to the following standard(s):
EN 60204-1:2006/A1:2009, EN 61010:2010.

**Supplementary
Information** The equipment also complies with the following standards:
UL3101-1:1993, NFPA70:2011, CAN-CSA 22.2 No. 1010-010-30,
FCC 47 CFR Part 15 Class A, AS/NZS 2064.1/2.

The above has been declared according to the global method, module A

Date: 31.01.2011

Christian Skjold Heyde,
Vice President, R&D and Production, Struers A/S

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Christian Skjold Heyde,
Vice- President, R& D et Production, Struers A/S

< 15557005 >

The TegraSystem is designed as a flexible modular system, which you can combine to fulfil your specific needs. Use a single element or several. Expand and adapt the modules to meet your developing needs.

With the TegraSystem you will gain the benefits of efficiency, reproducibility and network integration.

And you will enjoy the reliability of proven technology from the market leader.

High efficiency and reproducibility is ensured due to precise control of all preparation parameters

Consumable costs are reduced and specimen throughput is increased due to improved efficiency

Easy to operate with large graphical display and pre-installed preparation methods

LED light illuminates the preparation area

The emergency switch stops all moving parts in case of an emergency

Networking allows you to share your methods and improve efficiency and standardization

Preparation methods are automatically recalculated to suit different system sizes

Only one power, water and air connection is necessary to supply the whole system

Durable and proven technology ensures years of trouble-free operation

The modular and flexible system

The Tegra family

With the TegraSystem you can choose between 2 grinding and polishing machines, 1 specimen mover and 2 dosing systems. The machines can be combined in several ways to suit your specific needs.

Reproducibility is the key

In all materialographic specimen preparation, being able to reproduce preparation results that are exactly the same every time is crucial to success. Many factors have an influence on your ability to achieve the necessary reproducibility, and many of these can be controlled by employing automated preparation equipment like the TegraSystem.

Ensuring reproducibility and efficiency

Your workload and need for flexibility should determine the size of your TegraSystem. However, no matter how many or few specimens you need to prepare, the investment in a TegraForce specimen mover will be beneficial. The operator avoids the often tedious manual preparation tasks, and since the TegraForce specimen mover controls the force applied precisely, specimens are prepared exactly the same way every time, independent of operator skills. For labs with many specimens, specimen movers are the only way to ensure cost-effectiveness since several specimens can be prepared simultaneously.


Adding a TegraDoser-5 dosing system further improves reproducibility, efficiency and cost-effectiveness. With the ability to store and recall preparation methods, this dosing system yields exactly the same preparation results every time the preparation procedure is performed.

The high level of automation also allows for unattended preparation, thus improving the lab efficiency by giving the operator free time for other tasks. Cost-effectiveness is further increased by automating and controlling consumable consumption.

The innovative approach

The modular concept lets your TegraSystem grow and develop with your needs and challenges. For the lab with few specimens or repetitive preparation of the same material every day, the TegraDoser-1 is an innovative approach to automated preparation.

TegraDoser-1 is a single dosing unit, which is attached directly onto the TegraForce specimen mover. The dosing unit not only controls the dosing of diamond suspension or lubricant precisely, but also controls the preparation parameters for that specific preparation step.

With 2-3 TegraDoser-1 units it is possible to run an entire preparation method simply by exchanging the dosing unit. Reproducibility is ensured since the preparation parameters are exactly the same every time. The necessary investment is low in relation to the benefits achieved.

Combined with the all-in-one diamond suspensions from Struers, this solution allows for unattended preparation and ensures high quality and cost-effective preparation results.

As you will see on the following pages, the TegraSystem holds unique advantages for fast, efficient and cost-effective grinding and polishing in a safe and healthy working environment.

The TegraPol grinding and polishing machines

The TegraPol grinding and polishing machines are a series of sturdy, reliable and powerful machines, based on proven technology and designed for preparation of all types of materialographic specimens.

- Can be used for manual preparation or be equipped with specimen mover and dosing system.
- Available in 2 speed versions for standard specimen preparation and special applications
- Easy and straightforward operation
- TegraPol automatically detects connected units, no configuration is necessary
- TegraPol power supplies all connected units
- Sturdy and powerful design ensuring years of trouble-free operation

The 200 mm TegraPol-11 and -15


If you are preparing a limited number of specimens and still have high requirements to preparation quality, the 200 mm single disc machines are the perfect solution. TegraPol-11 is a unit covering all basic requirements with a speed of 50-300 rpm, whereas TegraPol-15 has a higher speed range from 40-600 rpm.

Added equipment is automatically detected

As soon as a TegraForce or a TegraDoser is connected, the TegraPol automatically detects which unit has been connected and no configuration is required. The operation of the whole system is carried out on the TegraPol. A large graphical display ensures easy and straightforward operation.

In the TegraSystem, the TegraPol is the only machine connected to the power mains. All other units are supplied directly from the TegraPol. The advantage here is of course the limited need for power outlets at your facilities. The power supply for all connected equipment is only 24 V.

For TegraPol-11 / -15 a disposable bowl liner is available to facilitate cleaning of the machine considerably.

The oval opening on the TegraPol machines grants easy access when exchanging consumables. At the same time it allows to place the specimen holder or specimen mover disc in the best possible position over the preparation disc to achieve optimum preparation results

Emergency stop

All TegraPol machines are equipped with an emergency stop, for immediate stopping of all moving parts in case of an emergency.

Environmentally-friendly recirculation system

A recirculation cooling unit can be connected to the TegraPol grinding/polishing machines. Controlled from the TegraPol, it is started and stopped automatically. Polluted debris from the preparation can be collected safely and the cooling water is utilised repeatedly. This reduces pollution and waste of water and eliminates the risk of corrosion of specimens by ordinary tap water.

Recirculation System

The TegraForce specimen mover

Mounting a specimen mover on the TegraPol is the first step towards automatic preparation. The precise control of all preparation parameters including the force is necessary to ensure reproducible preparation results. This is only obtained by electronic force control and, on TegraForce; the force on the specimens can be exerted very precisely, and the advantages are many:

- High reproducibility through precise control of preparation parameters
- Preparation of single specimens with individual force applied – flexibility is maintained and there is no need for dummy specimens

Preparation of single specimens

When single specimens are inserted in a specimen mover plate, the force is applied individually to each specimen. With this outstanding feature there is no need for dummy specimens. If your workflow demands short response-times you will benefit greatly from the individually applied force and maintain maximum flexibility. The preparation of individual specimens starts with grinding using fine grain sizes. There is no need for coarse pre-grinding since each specimen is levelled and ground independently. Very low force (10 N) can be applied on the individual specimen allowing for preparation of very fragile materials.

TegraForce-1

TegraForce-1 is designed to be mounted on a TegraPol-11 or TegraPol-15 and is suitable for preparation of both single specimens and specimens clamped in a specimen holder. 1, 2 or 3 specimens of up to 40 mm/1½" diameter can be fine ground and polished individually. For high removal rate or when plane grinding is required, the 120 mm diameter specimen holder where 3 specimens are prepared simultaneously is used.

To avoid damage when the preparation is started, the specimen holder or the pressure feet are lowered at reduced force. After a few seconds, the force applied to the specimens is slowly raised to the pre-set value.

On TegraForce the force can be reduced automatically towards the end of the preparation, resulting in lower material removal and smaller scratches. This saves time on the following preparation steps, since the amount of deformation to be removed is minimised.

TegraForce is equipped with LED illumination. While providing extra light to the preparation area, LED's grant very long life-time and minimum power consumption, making them most environmentally-friendly.

TegraForce-1

Accessories

Struers offers a wide range of accessories to the TegraForce-1 specimen mover. For TegraForce-1 both specimen mover plates and specimen holders are available, enabling the preparation of many different sizes and shapes of specimens. For the manufacture of your own special specimen holder designs, Struers can deliver specimen mover plates without holes.

LED illumination

The TegraDoser-5 dosing system

Adding an automatic dosing system to your TegraSystem gives you complete control of the entire preparation process and thus the quality of your preparation results and the costs involved.

- Built-in database with 10 preprogrammed Metalog Guide methods
- Built-in database with up to 200 methods leaves you in complete control
- Easy programming and large graphic display
- Excellent reproducibility
- Cost saving specimen preparation
- High-quality peristaltic pumps guarantee exactly the same dosing every time
- No vaporization or spray mist
- Automatic recalculation of method parameters when exchanged between different machines
- Automatic dosing of Oxide Polishing Suspensions
- Configuration according to your needs

Control the quality of your results

Manual dosing always implies certain variations in the amounts of abrasives and lubricants being dosed, and in the intervals between each dosing. TegraDoser-5 ensures that exactly the same amount is dosed each time – no matter how many times a preparation sequence is carried out.

But TegraDoser-5 is much more than just a dosing system. The built-in database can contain up to 200 preparation methods that can be stored and recalled at any time. This gives you full control of all preparation parameters and methods, regardless of which technician, machine or laboratory performs the job. Exactly repeatable settings give reliably reproducible results and frees resources from routine checks and repeats.

On a large graphical display preparation methods can easily be programmed and edited, whereas the method is selected and started from the TegraPol.

No vaporization or spray mist

The peristaltic pumps on the TegraDoser transport both suspensions and lubricants to the polishing cloth. In other dosing systems the suspension is vaporized and sprayed onto the cloth, resulting in an unhealthy spray mist. With the TegraDoser the liquid is dosed exactly where you want it on the cloth, improving preparation quality and improving the working environment.

Control your consumables costs

With the help of TegraDoser the dosing of suspensions and lubricants can be adjusted to the exact amount necessary for the preparation method in question. This will result in substantial savings on the consumables account. And TegraDoser also enables the operator to attend to other tasks during the preparation process – the dosing is automatically taken care of.


Automatic recalculation of preparation methods

A preparation method always needs to be adapted to the disc size of the grinding and polishing machine employed. As an example, the amount of suspension and lubricant that has to be applied on a 200 mm disc system is much less than the amount needed for a 300 mm disc system, simply because the surface area is smaller. The TegraSystem takes complete care of these required adjustments. Since the TegraSystem automatically detects the connected machines and units it can also automatically recalculate preparation methods to fit all system sizes. This ensures that preparation quality and reproducibility is maintained when preparation methods are being exchanged between machines, colleagues and labs.


Automatic dosing of Oxide Polishing Suspensions

The use of peristaltic pumps also makes it possible to include the dosing of oxide polishing suspensions in an automatic dosing unit. No other separate dosing units offer this feature. In addition, when working with oxide polishing suspensions, water is automatically applied onto the polishing cloth before and after the preparation. This ensures that the cloth is prepared prior to the polishing and that both specimens and cloth are cleaned afterwards.

TegraDoser-5 can be equipped with up to 7 peristaltic pumps, 6 for diamond suspensions or lubricants and 1 for Oxide Polishing Suspensions.

Choose between 3, 4, 5 or 6 pumps for diamond suspensions/lubricant both with and without the pump for Oxide Polishing Suspensions

Data management and method details on TegraDoser-5

The TegraDoser-1 dosing system

If your specimen preparation only involves one method or a small number of specimens, the TegraDoser-1 concept can increase the automation level drastically at a low investment cost.

- High reproducibility and cost-efficiency are ensured due to precise dosing
- Individual ID settings let you dedicate one dosing unit for each preparation step
- Entire preparation methods can be run just by exchanging dosing units
- Especially well-suited for Struers all-in-one diamond suspension products allowing unattended specimen preparation

The single dosing unit TegraDoser-1 is attached directly onto the TegraForce and its built-in peristaltic pump enables precise control of the applied amount of suspension – your guarantee for constant dosing levels and high reproducibility.

The dosing levels are equivalent to the dosing levels of TegraDoser-5 and can be taken over directly. With the peristaltic pump the dosing levels are absolutely reproducible and all types of lubricants and diamond suspensions can be used.

Run an entire method with no programming at all

TegraDoser-1 is equipped with a small switch, which allows 10 different ID settings. The TegraPol recognizes the ID setting of each TegraDoser-1 and automatically recalls all the parameter settings related to the specific preparation step, which was run the last time the TegraDoser-1 in question was used. In this way, you can dedicate a TegraDoser-1 unit to each preparation step, and an entire preparation method can be run without any programming, simply by exchanging TegraDoser-1 units. This not only facilitates operation it also provides higher security because of increased reproducibility.

The winning combination of DiaPro and TegraDoser-1 lets you control the entire preparation method just by changing the dosing unit

Innovative all-in-one diamond products

Struers has developed innovative all-in-one DiaPro and DiaDuo diamond suspensions which combine diamond suspension and cooling lubricant into one product. This not only facilitates easy-of-use but also ensures that the consumables costs are kept under control. These products are extremely well-suited for the TegraDoser-1 dosing system as they will allow for unattended specimen preparation.

For more information about Struers diamond products, please ask for separate product brochures.

The StruersLAN networking module

The TegraDoser-5 can be equipped with a StruersLAN module. This allows you to connect your TegraSystem to your company network or to a single PC.

- Easier programming of preparation methods from any connected PC in the company
- Preparation methods can be exchanged between machines, colleagues and labs
- Standardisation and efficiency is improved instantly
- Preparation methods and consumables lists can be downloaded from the Internet and used immediately

Networking simplifies programming

The StruersLAN module is a network card, which works according to the Ethernet standard, used worldwide in data communication. By providing the TegraDoser-5 with the StruersLAN module, the TegraSystem can be connected to a single PC or directly to the company's LAN network (Local Area Network), allowing the TegraSystem to be accessed from all the connected PCs.

In both cases programming of preparation methods is made extremely easy, because a complete keyboard is available. Also copy and paste functions, rename, delete, insert, - in short - all the functions usually available on a PC are available here too. The program is based on the same functions as the Windows™ Explorer.

Networking lets you share your knowledge

Do you want to share your preparation method with a colleague, another lab in your company or even another company? With the TegraSystem this is possible.

All TegraSystem modules connected to the same network can be accessed by PCs connected, and preparation methods can easily be copied directly from one machine to another.

Differences in machine size and set-up are not a problem since the TegraSystem automatically recalculates all preparation parameters when a method is transferred from one system to another.

Even when the TegraSystems are not connected to the same network, preparation methods can be exchanged. Simply export the preparation method, send it attached to an e-mail and import it into another TegraSystem without any problems.

Networking improves standardisation and efficiency

The networking possibilities greatly improve visibility and efficiency. Preparation methods can easily be made accessible on all the TegraSystems in the lab. Redundant work performed in different labs within the company can be eliminated. A central lab can standardise company methods and easily transfer them to TegraSystems in all satellite labs, ensuring consistent settings throughout the organisation. As soon as a method is approved, it can be transferred and made accessible instantly. Even specimen preparation performed at sub-suppliers can be standardised. The possibilities are numerous as required to match the different needs in modern companies.

Download of methods from the Internet



With a TegraSystem and StruersLAN networking possibilities, you benefit from the Internet tools provided by Struers. Hundreds of different preparation methods can be downloaded directly from Struers e-Metalog on the Internet in a format that is compatible with the TegraSystem. The methods can be employed instantly or adjusted to fit your specific application.

Also available from our homepage is the consumables table used by the TegraSystem. Whenever we release new grinding discs, polishing cloths, diamond suspension etc. the list will be updated and is immediately available for download. The advantage of this list lies in the built-in interaction between the consumables. If you e.g. select SiC as abrasive, the lubricant is automatically set to water and it will not be possible to select a suspension or a lubricant.

Full network integration

State-of-the-art communications technology gives you the potential to synchronise machine settings and preparation methods throughout your laboratory network. You can implement changes the moment they are approved. Using a central database simplifies the standardisation processes and keeps you in full control.

The flexibility of the TegraSystem allows you to tailor the system to meet your exact demands. Use a single element or several. Grow and adapt the modules to meet your developing needs. Gain the benefits of reproducibility and efficiency – in other words, let the TegraSystem revolutionize your world of sample preparation.

Hundreds of different preparation methods can be downloaded directly from Struers e-Metalog on the Internet.

State-of-the-art communications technology gives you the potential to synchronise machine settings and preparation methods throughout your laboratory network.

Grinding/Polishing machine	TegraPol-15	TegraPol-11
Disc diameter	200 mm / 8"	
Disc speed	40 - 600 rpm	50 - 300 rpm
Rotational direction	CCW	
Power supply	1/3 200-240 V / 50-60 Hz	
Motor Power continuous / max	370 / 555 W	180 / 270 W
Automatic torque control	Yes	
Emergency stop	Yes	
Touch Pad control	Yes	
Graphical display for easy selection of preparation methods	Yes	
Soft start with pre-dosing	Yes	
Controlled acceleration and deceleration	Yes	
High speed drying of preparation disc, Spin Function @ 600 rpm	Yes	
Max. torque at disc	>20 Nm	>10 Nm
Compatible with Struers MD-system	Yes	
Automatic detection of connected equipment	Yes	
Connection of a polishing head	TegraForce-1	
Maintenance management	Diagnostic and time used	
Service features	Machine test and self diagnostic	
Width, Depth, Height and Weight	410 mm / 670 mm / 315 mm / 25 kg	410 mm / 670 mm / 315 mm / 23 kg

Polishing head	TegraForce-1
Single specimens	Yes, 1 - 3 specimens
Specimen holder	Yes, 3 specimens
Rotational speed	50 - 150 rpm
Rotational direction	CCW + CW
Motor Power continuous / max (power supplied from TegraPol)	35 / 53 W
Force per specimen, single specimens / specimen holder	10 - 50 N / 30 - 150
LED lighting	Yes
Lifting capacity for specimen holders	-
Connection of dosing system	TegraDoser-5 or TegraDoser-1
Identical Start and Stop position of polishing head	Yes
Soft start function	Yes
Force reduction	Yes
Automatic lift-off after preparation	No
Locking of polishing head	Manual
Push button for quick release of specimen holders	No
Guide for positioning of polishing head	Yes
Fixed stop position of polishing head	Yes
Service features	Machine test and self diagnostic
Width, Depth, Height and Weight	140 mm / 311 mm / 445 mm / 7 kg

Dosing unit	TegraDoser-5	TegraDoser-1
Precise dosing with peristaltic pumps	Yes	
Number of pumps	Up to 7 pumps	Single pump
Handles lubricant, suspension and all-in-one suspensions	Yes	
Handles OP-suspension	Yes	-
Integrated automatic cleaning function	Yes	-
Graphical display for programming of preparation methods	Yes	-
Integrated database	10 Metalog Guide methods and 200 user methods	-
ID switch for parameter step control	-	Yes
Download of preparation methods from Internet (Struers e-Metalog)	Yes	-
Possibility to program methods via PC	Yes	-
Access to consumables list	Yes	-
Width, Depth, Height and Weight	200 mm / 210 mm / 380 mm / 8.5 kg (excl. bottle tray)	120 mm / 50 mm / 240 mm / 550 g

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CANADA

Struers Ltd.
7275 West Credit Avenue
Mississauga, Ontario L5N 5M9
Phone +1 905-814-8855
Fax +1 905-814-1440
info@struers.com

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Struers' equipment is in conformity with the provisions of the applicable International Directives and their appurtenant Standards. (Please contact your local supplier for details)

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice.

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Struers Inc.
Struers GmbH