Transpol-2

Instruction Manual
Gebrauchsanweisung
Mode d'emploi

Manual no.: 15957001

Date of release 07.08.2007

Instruction Manual

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the *Serial No.* and *Voltage* on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2007.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Transpol-2 Safety Precaution Sheet

To be read carefully before use

1. The operator should be fully aware of the use of the apparatus according to the Instruction Manual.
2. If grinding accessories not recommended by Struers are used, they must be approved of for speeds above 8500 rpm. Do not use grinding stones with diameter larger than 10 mm.
3. Make sure that the current-carrying parts are kept in order. If any damage to a cable has been observed it must be corrected immediately.
4. Use an isolating transformer in wet and moist rooms.
5. Always use safety goggles when working with Transpol-2.
6. Always protect the motor and control unit from rain and moisture.
7. The grinding unit must not be used without reduction gear.
8. Always clamp the grinding tools carefully.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

Disposal

Equipment marked with a WEEE symbol  contain electrical and electronic components and must not be disposed of as general waste.

Please contact your local authorities for information on the correct method of disposal in accordance with national legislation.

| Table of Contents | Page |
|------------------------------------|------|
| 1. Getting Started | |
| Checking Contents of Packing | 2 |
| Mounting of Motor | 2 |
| 2. Operation | |
| Getting Acquainted..... | 3 |
| Coarse Grinding | 5 |
| Grinding and Fine Grinding | 5 |
| Grinding Time | 6 |
| Grinding Speed | 6 |
| Grinding Pressure | 6 |
| Practical Hints | 6 |
| Polishing..... | 7 |
| Overload of the Motor | 7 |
| 3. Maintenance | 8 |
| 4. Technical Data | 8 |

1. Getting Started

Checking Contents of Packing

In the packing case you should find the following parts:

- 1 Transpol-2
- 1 Carrying case with shoulder strap
- 1 Mains cable
- 1 DC motor
- 1 Reduction gear
- 1 Straight handle with holder for flapper wheel
- 1 Flapper wheel, grit 80
- 1 Right angle handle with quick coupling for rubber discs
- 4 Rubber discs
- 1 Felt cone
- 1 Square box for storing polishing cloths and felt cone in use
- 3 Spanners for tightening grinding tool
- 1 Spanner
- 1 Lock pin for use when clamping flapper wheel
- 1 Pair of safety goggles
- 2 Fuses
- 1 Set of Instruction Manuals

Mounting of Motor

- Place the plug of the motor cable in the socket on the front plate and tighten the nut.
- Screw motor and reduction gear carefully together.

If the parts cannot be screwed correctly together it is because the driver of the motor is not correctly positioned. By turning the motor shaft the parts can be made to line up and then be tightened up with the spanner.

Never use force!

IMPORTANT


The reduction gear must remain mounted on the motor, whether the straight handle or the right angle handle is to be used.

2. Operation

Getting Acquainted


Take a moment to familiarise yourself with the location and names of all the Transpol-2 components:

Front Panel

- ① Main switch
- ② Overload indicator
- ③ Start/ Stop button
- ④ Speed display
- ⑤ Change direction button and direction indicators
- ⑥ Speed control
- ⑦ Socket for motor cable

Transpol-2
Instruction Manual

- ① Front Panel
- ② Carrying case with shoulder strap
- ③ Box for storing polishing cloths and felt cone in use
- ④ Boxes for storing new grinding discs and polishing cloths (consumables are ordered separately)
- ⑤ Felt cone
- ⑥ Right angle handle with quick coupling for rubber discs
- ⑦ DC motor
- ⑧ Reduction gear
- ⑨ Straight handle with 6 mm chuck for flapper wheel
- ⑩ Flapper wheel
- ⑪ Rubber disc
- ⑫ Spanner and lock pin for tightening grinding tool

Coarse Grinding

- Screw the straight handle on the reduction gear in the same way as the reduction gear was mounted on the motor.
- Place the flapper wheel in the handle, put the lock pin in the hole in the middle of the handle and fasten the flapper wheel in the handle with the spanner.
- Switch the mains switch to ON. The display will light up.
- Press the Start/ Stop button and the motor starts rotating.
- Adjust the speed as required. Coarse grinding with flapper wheel requires the highest speed.

Grinding and Fine Grinding

- Dismount the straight handle (the reduction gear must remain mounted on the motor).
- Mount the right angle handle and tighten it with the spanner in the same way as the reduction gear was mounted on the motor.
- Place the rubber disc in the quick coupling device of the right angle handle.
- Hold back on the angle head with your thumb and at the same time turn the rubber wheel and press it into place.
- Take a piece of grinding paper with the appropriate grit size. Remove the cover paper at the rear and press the grinding paper against the rubber disc with your thumb.

Normally, grinding is carried out in 4 steps (grit 60, 120, 240, 500).

Transpol-2
Instruction Manual

| | |
|--------------------------|--|
| <i>Grinding Time</i> | As a rule of thumb, one grinding step should be twice as long as it takes to remove the scratches from the previous step. |
| <i>Grinding Speed</i> | The speed should normally be highest during coarse grinding and then be reduced with the smaller grit sizes. |
| <i>Grinding Pressure</i> | Depends on the roughness of the surface and the condition of the grinding paper. An unprepared surface and a new piece of grinding paper requires only little pressure. |
| <i>Practical Hints</i> | <p>Do not use the whole surface of the grinding paper, just grind with the outer edge of the paper. This gives better grinding and a more uniform scratch pattern.</p> <p>The right angle handle, gear and motor may become warm to hold, but this does not mean that anything is wrong. As the reduction gear throttles down the rotational speed some energy will be lost as heat.</p> <p>The grinding papers should be kept separate from the polishing cloths in order to prevent the grinding particles contaminating the polishing cloths.</p> |

Polishing

Polishing is carried out with DP-paste, DP-spray or DP-stick on a polishing cloth in 2, 3 or 4 steps, dependant on the material, e.g. 15 µm, 6 µm, 3 µm and finally 1 µm. If the surface has to be especially fine ¼ µm diamonds may be used for final polishing, but normally the preparation can be finished at 3 µm.

| Cloth | Diamond Type |
|-------------------|-------------------------|
| DP-Dur/ DP-Dac | 15 µm and 6 µm diamonds |
| DP-Mol | 3 µm and 1 µm diamonds |
| DP-Nap | ¼ µm diamonds |

- Mount the cloth on the rubber disc by removing the cover paper and pressing the cloth against the rubber disc with your thumb.
- Apply Diamond Paste or Spray on the cloth.

Before polishing and between the individual polishing steps the area to be prepared must be cleaned with water or alcohol.

During polishing the motor speed must be relatively low, in order not to spread the diamonds outside the preparation area. Add a few drops of lubricant, blue lubricant for hard materials like steel, red lubricant for soft materials like aluminium and brass.

Overload of the Motor

If the motor becomes overloaded the inner fuse will switch off the power for the motor which will then stop. At the same time the overload indicator will flash. The apparatus can be re-started by first setting the Main switch to OFF and then again to ON.

3. Maintenance

- Clean Transpol-2 with a soft, slightly damp cloth.

4. Technical Data

| Subject | Specifications | | |
|--|---|--------|--------|
| | | Metric | US |
| Max. power consumption | 64 W | | |
| Rotating speed | <i>Motor including gear: about 100-7000 rpm</i> | | |
| Power Supply | 50/60 Hz / 230 V 50/60 Hz / 120 | | |
| Dimensions and Weight (complete equipment in carrying case) | Width | 230 mm | 9.1" |
| | Depth | 150 mm | 5.7" |
| | Height | 250 mm | 9.8" |
| | Length of flexible cable: up to | 2 m | 6.6' |
| | Weight (incl. accessories) | 5 kg | 11 lbs |

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
Spare parts list

| Pos. | Description | Cat no. |
|------|-----------------------------------|------------|
| | Carrying case with shoulder strap | 437 MP 039 |
| | Carbon brushes for motor | 437 MP 040 |
| | Motor cable NCL-26CD | 437 MP 030 |

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
Herewith declares that

| | |
|----------------------|--------------------------------|
| <i>Product Name:</i> | Transpol-2 |
| <i>Type No:</i> | 595 |
| <i>Machine Type:</i> | Grinding and polishing machine |

is in conformity with the provisions of the following directives:

Safety of Machinery 98/37/EEC according to the following standard(s):
EN ISO 12100-1:2003, EN ISO 12100-2:2003, EN60204-1:1997.

EMC-Directive 89/336/EEC and 92/31/EEC according to the following standard(s):
EN50081-1:1992, EN50082-1:1997.

Low Voltage Directive 73/23/EEC according to the following standard(s):
EN60204-1:1997.

Supplementary Information The equipment complies with the American standards:
UL508.

The above has been declared according to the global method, module A

Date: 07.08.2007

Christian Skjold Heyde,
Vice President, R&D and Production, Struers A/S

Hersteller Struers A/S
Pederstrupvej 84
DK-2750 Ballerup, Danmark
Telefon +45 44 600 800

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Date: 07.08.2007

Christian Skjold Heyde,
Vice- President, R & D et Production, Struers A/S

On-site materialographic preparation and examination

Non-destructive preparation

Conventional materialographic preparation entails cutting a sample from a material. This method is not always appropriate as it destroys the component which is being examined.

Many applications require the ability to perform a materialographic preparation and examination without ensuing damage to the subject. The non-destructive method is used in these cases.

The usual preparation steps such as grinding, polishing and etching are performed on the part to be examined in a limited area where no damage can be done, whereupon the surface can be analysed.

On-site preparation is particularly suitable for quality checks in power stations, oil platforms, bridges, aircraft etc. The non-destructive preparation method enables the necessary inspection to be carried out on site.

Struers manufactures a complete range of portable equipment for non-destructive preparation, from basic grinding through mechanical and electrolytic polishing to electrolytic etching. The range also includes a portable microscope for on-site use and methods for producing replicas of the prepared surface.

Using replica methods makes it possible to perform the microscopic examination in the laboratory - under perfect working conditions.

The entire Struers non-destructive preparation range has been designed for field work. The equipment is thus as compact and as light as possible. It is easy to transport and can be used literally under any conditions.

Transpol-2

All the necessary accessories are contained in Transpol-2's carrying case.

Transpol-2 - Mechanical grinding and polishing under all working conditions

Transpol-2 is a portable grinding and polishing apparatus. It is extremely compact and is designed for field work under even the most difficult conditions. Transpol-2 fits into a practical shoulder bag with all the necessary accessories and weighs an approximate total of only 5 kg. Transpol-2 consists primarily of 2 parts: a control unit and a grinding and polishing unit. The control unit comprises the power supply and speed regulator unit (100-7000 rpm). The grinding and polishing unit is connected to the control unit by a soft, flexible cable. It consists of a direct current motor with a holder with an easily exchangeable rubber disc, upon which the grinding paper and polishing cloth are mounted.

Transpol-2 is available with two types of holder: one for straight mounting and one for right-angle mounting, thus making it possible to use Transpol-2 in even the most inaccessible places.

Grinding

Transpol-2 uses 32 mm diameter SiC grinding paper discs in grit sizes 60, 120, 240 and 500#. Transpol-2 also includes a special flapper wheel for grinding of extremely uneven surfaces.

Polishing

Transpol-2 uses DP-Dac, DP-Dur, DP-Mol, DP-Nap and OP-Felt polishing cloths, all 32 mm diameter, and DP-Paste, DP-Stick or DP-Spray. Transpol-2 is furthermore delivered with a special felt cone for polishing of uneven and curved surfaces.

Movipol-3

Due to the special design of the polishing pencil, electrolytic polishing with Movipol-3 is quick and uncomplicated

Movipol-3 - Electrolytic polishing and etching in the field

Movipol-3 is a portable electrolytic metal polishing and etching apparatus. It is extremely compact and robust and can be used anywhere. Electrolytic preparation is a particularly fast and efficient method of non-destructive metallographic preparation. It can be applied directly to critical surfaces and is widely used for metal safety inspection, especially for critical parts of larger units, e.g. for welds and other joints. The method effectively reveals changes in the microstructure of the metal and allows to prevent possible damage because of cracks and leaks.

Automatic polishing and etching

Movipol-3 consists of a compact unit comprising the power supply, control unit, pump motor, electrolyte container and polishing unit. When the polishing unit is pressed against the metal surface, an electrolytic reaction is started. After 5-10 seconds of use, the prepared area is as reflective as a mirror and etching can commence.

PSM-2 Transcopy

The entire process is automatic: after completion of the polishing, Movipol-3 proceeds to the etching process. When this is finished, the unit emits an audio signal. In less than one minute a surface is obtained which is ready for analysis, either on-site or by means of replica for laboratory examination.

Complete Freedom

Movipol-3 is equipped with rechargeable batteries and supply transformer. In other words, the apparatus is as suitable for stationary laboratory work as it is for field work, depending on the requirements.

PSM-2 and Transcopy constitute an indispensable tool for non-destructive materialography.

PSM-2 - Fast and easy microscopy outside the laboratory

PSM-2 is a small, portable microscope for use in the field. It is battery-powered and provides magnifications of 100 x, 200 x and 400 x. PSM-2 comes in a shoulder bag which also contains the batteries.

PSM-2 is particularly suitable for preliminary examination of the prepared surface. A replica is recommended if a more detailed analysis is required. PSM-2 is highly suitable for preliminary replica examinations when used together with Transpol-2 and Movipol-3.

Transcopy Replica Foil


Transcopy Replica Foil is used to make a replica of any polished and etched surface. It is especially used if on-site analysis using a portable microscope is not sufficient or if micrographs of the prepared structure are required. Transcopy Replica Foil creates a permanent replication of microstructures, cracks and defects for future analysis and comparison in the laboratory.

Function

Transcopy Replica Foil consists of reflecting plastic film with a self-adhesive back. A replica is made by applying Transcopy fluid to the film and pressing it firmly onto the prepared and etched surface. A slightly over etched structure results in a higher contrast of the replica. After 4-5 minutes the foil is removed from the surface. By removing the cover paper, the replica can be adhered to a plain plate and then examined under an optical microscope.

Transcopy

With Transcopy, Replica foil, replicas can be taken on the spot for subsequent examination in the laboratory.

RepliSet

The **RepliSet** system is designed to produce an exact 3D copy of a surface. It is used for non-destructive testing and field applications allowing the structure or irregularities on critical components to be examined and measured under laboratory conditions.

RepliSet is a specially formulated fast curing two-part silicone rubber with a good releasing ability for flexible high-resolution 3D replicas, which behave like a metallic surface when examined in an optical microscope. RepliSet compounds are supplied in cartridges and are dispensed using a hand-operated dispensing gun. The cartridges contain both polymer and curing agent, which are automatically mixed in a disposable static-mixing nozzle during application to the surface. The application system offers superior and fast results, regardless of the conditions.

RepliFix is a less advanced parallel to RepliSet. RepliFix and RepliSet are designed to bond together. The two components are mixed and applied by hand. RepliFix is used as support for RepliSet or as a stand alone product for moulding of surface shape for low tech applications.

A backing slide bonds to the RepliSet or RepliFix replica. The backing slide serves to maintain the original profile and ensures a flat back to the replica.


A replica can be taken from all metallic materials and most other solid materials like ceramics, plastics and glass. There are no size, shape or thickness limitations on the replicas that can be made and it is even possible to take replicas from otherwise inaccessible surfaces. The RepliSet system can produce replicas, which at the same time are dimensionally correct, with ultra fine detail reproduction and with a flat back.

An operator dispenses RepliSet onto a sample, and then peels off the cured replica.

*Pure copper.
Sample etched with cupric chloride and ammonia.
Magnification 100 x*

Sample Replica

Further to examination by optical microscopes, 3D examination can be carried out using non-contact measuring instruments such as laser measuring equipment or measuring projectors with 2D or 3D facilities. A replica of a cavity, for instance inner holes, can be examined using measurement and profile projectors. The replicas are suitable for 3D examination by SEM.

Metallographic applications

Typical applications are on-site non-destructive testing in connection with quality control, inspection and maintenance of power plants, oil platforms, bridges, aircraft, etc.

Water or high temperatures do not affect RepliSet so the application range is very broad.

Typical tasks are:

- Quality control of edges, corners, heights, angles, surface finish and other dimensions.
- Inspection of internal surfaces such as bolt hole threads.
- Detection and monitoring of pitting, corrosion, cracking, creep and wear.
- Metrology.
- Forensic investigation.

RepliSet 50 ml system.

Case for RepliSet

The RepliSet Case is designed for transportation and use of the 50 ml RepliSet system. It is made of aluminium and is at the same time elegant and sturdy. It can be carried as hand luggage on flights and is compact and sturdy to such an extent, that it can be taken to locations with narrow or difficult access.

The contents is either fixed by straps in the lid or placed in compartments in the two detachable foam rubber inserts. Each item has its fixed position. The user has access to all that is needed to perform a regular replication by just opening the lid of the RepliSet Case. The lower insert carries a small stock of consumables.

Transpol-2

Technical Data

| | |
|-----------------------------------|-----------------------------------|
| Voltages | 230 V, 50/60 Hz; 120 V, 50/60 Hz; |
| Speed | 100-7000 rpm |
| Dimensions (without shoulder bag) | Width 180 mm |
| | Height 75 mm |
| | Depth 150 mm |
| Total weight incl. accessories | 5 kg |

Specifications

| | |
|--|----------|
| Transpol-2, complete with soft cable, DC motor 230 or 120 V, straight handle, right angle handle, flapper wheel, felt cone, 4 rubber discs | 059561xx |
| Rubber disc, 30 mm dia. | 40800055 |
| SiC grinding papers, adhesive, bundle of 100 | |
| Grit 60 | 40400049 |
| Grit 120 | 40400129 |
| Grit 240 | 40400130 |
| Grit 500 | 40400131 |
| DP cloths, adhesive, dia. 32 mm, bundle of 25 | |
| DP-Dur | 40500040 |
| DP-Dac | 40500145 |
| DP-Mol | 40500041 |
| DP-Nap | 40500042 |
| OP-Felt | 40500043 |
| Flapper wheel, 80 grit, 40 mm dia. x 10 mm | 40800053 |

Movipol-3

Technical Data

| | |
|--|-------------------------------|
| Voltage supply to transformer | 220 V 50/60 Hz, 115V 50/60 Hz |
| Voltage supply from transformer | 42 V AC |
| Input fuse | 5 A |
| Battery fuse | 10 A |
| Thermal overload protection | |
| Polishing power | Max. 140 VA |
| Polishing area | About 9 mm dia |
| Dimensions | Width 410 mm |
| | Height 280 mm |
| | Depth 200 mm |
| Weight with batteries, but without electrolyte | Net 8.80 kg |

Specifications

| | |
|---|----------|
| Movipol-3, complete with batteries and supply transformer, 220 V or 115 V | 039265xx |
| 42 V charging/supply transformer for continuous use for Movipol 3, 230 V or 115 V | 044462xx |
| Set of extra batteries (9 pcs.) | 03926901 |
| Polishing Chambers for Movipol-3, Flexible type, 10 pcs. | 03926904 |

PSM-2

Technical Data

| | | | |
|----------------------|----------------------|-----------|----------------------|
| Batteries | 4 x 1.5 V, type LR20 | Length | 210 mm (8.3") |
| Lamp (threaded base) | 5 V, 6 W | Diameters | 15-60 mm (0.6"-2.4") |
| | | Weight | 600 g |

Specifications

| | |
|---|----------|
| Portable microscope PSM-2, complete with 10 x eyepiece, 10 x objective, lamp housing and carrying case, without batteries | 04286101 |
| Objective, magnification 20 x | 04286901 |
| Objective, magnification 40 x | 04286902 |
| Measuring set consisting of measuring eye piece and plate micrometer | 04286903 |

Transcopy

Transcopy Kit

| | |
|--|----------|
| Set consisting of 40 ml Transcopy Liquid, 50 replica foils 20 x 30 mm (0.8" x 1.2"), 1 pipette, spray nozzle and 50 microscope slides 25 x 75 mm (1" x 3") | 40900090 |
| Transcopy Replica Foils, 50 pcs. | 40900091 |
| Transcopy Liquid, 40 ml | 40900092 |

RepliSet

Technical Data

| | |
|--|---|
| Resolution of cu red replica | Down to 0.1 micron |
| Shrinkage | Negligible |
| Tear Strength | 15-20 kN/m ² |
| Temperature range for the surface to be examined | -10°C to +180°C |
| Life span of the finished replicas is practically indefinite provided they are stored according to the instructions. | |
| Content in static-mixing nozzle | 1.1 ml in nozzle for 50 ml cartridge 9.3 ml in nozzle for 265 ml cartridge |

Specifications

Replication system for non-destructive testing of a microstructure or a 3D structure. Fast curing two-part silicon rubber compound for flexible high-resolution 3D replicas. For the 50 ml system, the hand-operated dispensing gun (40900066) and the static mixing nozzles (40900088) are used in combination with the 50 ml cartridges. For the 265 ml system, the hand-operated dispensing gun (40900065) and the static-mixing nozzles (40900056) are used in combination with the 265 ml cartridges.

Struers A/S
 Pederstrupvej 84
 DK-2750 Ballerup, Denmark
 Phone +45 44 600 800
 Fax +45 44 600 801
 struers@struers.dk

Specifications

Cat. no.

RepliSet-F1

Particularly useful for replicating horizontal or sloping surfaces in low temperature conditions or where rapid results are required. Fluid rapid curing compound with working life of 0.5-1 min. and curing time of 4 min. at 25°C
 1 cartridge of 50 ml
 5 cartridges of 50 ml
 2 cartridges of 265 ml

40900069
 40900047
 40900051

RepliSet-F5

General purpose material. Particularly useful for replicating horizontal or sloping surfaces in normal or high temperature conditions. Fluid fast curing compound with working life of 5 min. and curing time of 18 min. at 25°C
 1 cartridge of 50 ml
 5 cartridges of 50 ml
 2 cartridges of 265 ml

40900068
 40900046
 40900050

RepliSet-T1

Particularly useful for replicating vertical or overhead surfaces in low temperature conditions or where rapid results are required. Thixotropic rapid curing compound with working life of 0.5-1 min. and curing time of 4 min. at 25°C
 1 cartridge of 50 ml
 5 cartridges of 50 ml
 2 cartridges of 265 ml

40900071
 40900049
 40900053

RepliSet-T3

General purpose material. Particularly useful for replicating vertical or overhead surfaces in normal or high temperature conditions. Thixotropic fast curing compound with working life of 3 min. and curing time of 10 min. at 25°C
 1 cartridge of 50 ml
 5 cartridges of 50 ml
 2 cartridges of 265 ml

40900070
 40900048
 40900052

RepliSet-GF1

Replication system especially for comparatory macroscopy and metrology. Particularly useful for replicating horizontal or sloping surfaces and filling holes. Fluid rapid curing compound with working life of 0.5 - 1 min. and curing time of 4 min. at 25°C.
 1 cartridge of 50 ml
 5 cartridges of 50 ml

40900078
 40900076

RepliSet-GT1

Replication system especially for comparator macroscopy and metrology. Particularly useful for replicating vertical or overhead surfaces. Thixotropic rapid curing compound with working life of 0.5 - 1 min. and curing time of 4 min. at 25°C.
 1 cartridge of 50 ml
 5 cartridges of 50 ml

40900079
 40900077

RepliFix

Specially formulated hand mixed fast curing two-part silicone rubber. Bonds to RepliSet. Particularly useful in combination with RepliSet for producing a rigid backing. It can be used directly for moulding of surface shape for profile measurement

RepliFix-2

For low temperature conditions or where rapid results are required. Working life of 2-3 min. and curing time of 10 min. at 25°C. Net 500 g

40900084

RepliFix-20

For high temperature conditions or for taking replicas of complicated geometry or large areas. Working life of 20 min. and curing time of 60 min. at 25°C. Net 500 g

40900086

Accessories

Hand-operated dispensing gun

For 50 ml cartridges of RepliSet
 For 265 ml cartridges of RepliSet

40900066
 40900065

Static-mixing nozzles

For 50 ml cartridges. 35 pcs.
 For 265 ml cartridges. 10 pcs.

40900088
 40900056

Nozzle Tips

For replicating flat surfaces. Fishtail spreaders, 10 mm width.
 To be mounted on 50 ml static mixing nozzle (40900088). 30 pcs.

40900089

For replicating small holes. Luer needle, 1 mm dia., 30 mm long. To be mounted on 50 ml static-mixing nozzle (40900088). 10 pcs.

40900060

For replicating larger holes. Flexible hose, 6 mm dia., 100 mm long. To be mounted on 50 ml static-mixing nozzle (40900088) 10 pcs.

40900061

Backing Slides

A flexible plastic slide, which bonds to the replica and ensures a flat back to the replica. For levelling of replicas to assist microscopic examination, as dimensional support for metrology and for well-ordered labelling, transport and storage of RepliSet replicas 26 x 76 x 1 mm, 50 pcs.

40900087

Backing paper

For RepliSet replication system. Bonds to the replica and facilitates labelling, handling and the levelling of replicas to assist microscopic examination. 60 x 70 mm. 100 pcs.
 A4 (210 x 297 mm), for cutting up to the required size. 10 pcs.

40900062
 40900063

Case for RepliSet 50 ml system

Aluminium case with room for all necessities for field applications.
 The contents of the RepliSet Case is ordered separately

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CANADA

Struers Ltd.
 7275 West Credit Avenue
 Mississauga, Ontario L5N 5M9
 Phone +1 905-814-8855
 Fax +1 905-814-1440
 info@struers.com

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Struers Limited
Struers Inc.
Struers GmbH