Unitom-2

Instruction Manual

Manual No.: 15377001

Date of Release FGF€20FG

Table of Contents	Page
User's Guide	1
Reference Guide.....	12
Quick Reference Guide	25

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the *Serial No.* and *Voltage* on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

Original instructions. The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2012.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Unitom-2 Safety Precaution Sheet

To be read carefully before use

1. The operator(s) should be fully instructed in the use of the machine and its cut-off wheels according to the Instruction Manual and the instructions for the cut-off wheels.
2. The machine must be placed on a safe and stable support table. All safety functions and guards of the machine must be in working order.
3. Use only intact cut-off wheels. The cut-off wheels must be approved for min. 50 m/s.
4. Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
5. The workpiece must be securely fixed in the quick-clamping device or the like. Large or sharp workpieces must be handled in a safe way.
6. The cutting handle should be lowered slowly and carefully, in order to avoid breaking the cut-off wheel. The cutting handle automatically returns to its top position when you let go the handle.
7. Struers recommends the use of exhaustion as the materials to be cut may emit harmful gasses or dust.
8. Use of working gloves is recommended as workpieces may be both very hot and produce sharp edges.
9. Use of safety shoes is recommended when handling large or heavy workpieces.
10. Use of safety goggles is recommended when using the flushing hose.

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

User's Guide

Table of Contents	Page
1. Getting Started	
Checking the Contents of Packing	2
Unpacking Unitom-2.....	2
Placing Unitom-2.....	2
Getting Acquainted with Unitom-2.....	3
Noise Level	4
Supplying Power	4
Direction of the Cut-off Wheel.....	4
Connecting a Cooli Cooling Unit	5
2. Basic Operations	
Using the Controls.....	6
Front Panel Controls of Unitom-2	6
Front Panel Controls	7
Cooling Valve Positions	8
Changing the Cut-off Wheel.....	9
Clamping the Workpiece	9
Starting/Stopping the Cutting Process	10
Flushing the Cutting Chamber	11

1. Getting Started

Checking the Contents of Packing

In the packing box you should find the following parts:

- 1 Fork spanner, 30 mm, for changing the cut-off wheel
- 1 Grease gun, for lubrication of tilting frame bearings
- 2 Cover plugs for threads for cutting tables
- 1 Drain angle, 45⁰
- 1 Drain angle, 90⁰
- 1 Outlet hose 2 m
for connection to external cooling unit
- 2 Hose clamps, 35-60 mm
- 1 Set of Instruction Manuals

Unpacking Unitom-2

- Remove the wooden block holding the spindle using a 17 mm spanner.
- Insert the two cover plugs delivered with the machine in the threaded holes.

Placing Unitom-2


Unitom-2 should be placed on the optional Struers Table Unit, or a table which is able to carry a weight of min. 200 kg.

Adjust the table height using the screws in the table legs. It is important that the table and Unitom-2 are level. Place Unitom-2 on the table. If the machine is not horizontal, you must adjust the height for each corner, as follows:

- Loosen the upper 13 mm nut under the corner screw.
- Adjust the lower 8 mm nut with a spanner.
- Fasten the upper 13 mm nut again.
- Repeat the adjustment for each of the other three corners until Unitom-2 is levelled.

**Getting Acquainted
with Unitom-2**

Take a moment to familiarise yourself with the location and names of the Unitom-2 components.

- ① Front panel controls
- ② Cutting Table
- ③ Quick-clamping device
- ④ Table Unit (optional)
- ⑤ Recirculation cooling unit
- ⑥ Cutting handle
- ⑦ Side plate
- ⑧ Main switch
- ⑨ Protection guard with window
- ⑩ Cooling valve

Noise Level

Approx. 76 dB (A) measured at idle running, at a distance of 1.0 m/39.4" from the machine.

Supplying Power

IMPORTANT
 Check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.

- Unitom-2 is factory mounted with an electric cable. Mount a plug on the cable according to the following:

50 hz	60 hz
Phase Black and brown Earth Yellow/green	L ₁ Black L ₂ Red L ₃ Orange Earth Green

Direction of the Cut-off Wheel

Check that the cut-off wheel rotates in the direction indicated on the label attached to the guard of Unitom-2. If the direction of rotation is incorrect:

50 hz	60 hz
Switch two of the phases	Switch phases L ₁ and L ₂

Connecting a Cooli Cooling Unit



To ensure optimal cooling, Unitom can be fitted with a Cooli unit. The Cooling System 5 is a Cooli configuration designed for use with Unitom.

Note:

Cooling System 5 includes a static filter.
For intensive use, and for materials generating a lot of swarf, a bandfilter is recommended.

Note

Before connecting the cooling unit to the Unitom, follow the instructions in the Cooli Instruction Manual to prepare it for use.

To connect the Unitom to a cooling unit:

- Plug the Cooli control unit's communication cable into the Unitom's control socket.
- Connect the quick coupling on the Unitom's cooling fluid inlet tube to the cooling unit's pump outlet (A).

- Insert the drain pipe ① in the water outlet on the back of Unitom-2 and mount the 90° elbow pipe ②. Lubricate the sealing ring with grease or soap to facilitate insertion. (Use the other elbow pipe if more suitable).
- Strip the steel spring from approx. 3 cm of the outlet hose and cut. Bend the cut end towards the centre of the hose. Mount the outlet hose ③ onto the elbow pipe and clamp the stripped section using a hose clamp.
- Check that the outlet hose slopes downwards when connected. If necessary adjust the length of the hose.
- Insert the open end of the hose into mounting hole in the bracket on top of the Cooli filter unit (B).
- Connect the cooling unit to the mains power supply.


IMPORTANT

Before connecting, check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.


2. Basic Operations

Using the Controls

Front Panel Controls of Unitom-2





Location of the Main Switch




*Unitom-2
Instruction Manual*

Front Panel Controls


Name	Key	Function	Name	Key	Function
① MAIN SWITCH		The main switch is located to the right hand side of Unitom-2. Turn clockwise to switch on the power.	⑤ FLUSH		Lights when cooling valve is in flushing position.
② POWER		Lights when main power is switched on.	⑥ OVERLOAD		Lights when emergency stop has been pressed or the cutting motor is overloaded.
③ START		Starts the machine. The cut-off wheel starts rotating and the cooling water is turned on.	⑦ EMERGENCY STOP		- Push the red button to stop. - Pull the red button to release.
④ STOP		Stops the machine. The cut-off wheel stops rotating and the cooling water is turned off.			

© COOLING VALVE: See next page!


Cooling Valve Positions

<p>⑧ COOLING VALVE</p>	<p><i>The cooling valve for cooling/cleaning is located in the back of the cutting chamber to the left.</i></p>
	<p>Turn the knob counter-clockwise to vertical position for cooling of cut-off wheel.</p>
	<p>Turn the knob to regulate the flow while flushing the cutting chamber.</p>
	<p>Turn the knob clockwise to vertical position for flushing the cutting chamber with maximum flow.</p>

Location of Cooling Valve for cooling/cleaning

Changing the Cut-off Wheel

- Open the protection guard and swing the right cooling tube to the right.
- Press the black locking knob on the left-hand side of the cut-off wheel, turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner. Remove the flange and the old cut-off wheel.

IMPORTANT

The spindle at Unitom-2 is right-hand threaded.

- Mount the new cut-off wheel.

IMPORTANT

Conventional cut-off wheels based on Al_2O_3/SiC abrasives should be placed between two cardboard washers, to protect the cut-off wheel. For maximum precision with diamond or CBN cut-off wheels, do not use cardboard washers.

- Mount the flange and the nut. Tighten carefully.
- Swing the right cooling tube back into position.

Clamping the Workpiece

- Place the workpiece between the clamp and the back stop at the right hand cutting table.
- Push the clamp towards the workpiece and lock the quick-clamping device with the locking handle.

IMPORTANT

It is very important that the workpiece is firmly and immovably fixed in the quick-clamping device.


Starting/Stopping the Cutting Process

- Carefully close the protection guard.
- Press the button START  to start the machine. The cut-off wheel starts rotating and the cooling water is turned on.
- Carefully lower the cut-off wheel by drawing the cutting handle downwards until the cut-off wheel touches the workpiece.
- Let the cut-off wheel work itself a little into the workpiece. Then increase the force and begin the cutting.
- When the cut-off wheel has almost cut through the workpiece the cutting force should be reduced.
- When the cut-off wheel has cut through the workpiece the cutting handle is returned to its top position.
- Press the button STOP  to stop the cut-off wheel and cooling water.
- Open the protection guard and release the workpiece from the quick-clamping device.

IMPORTANT

Unitom-2 is equipped with a protection guard and a guard lock. The lock prevents access to the cutting chamber when the cut-off wheel is running. After pressing STOP  the protection guard cannot be opened for about 5 seconds. This safety function ensures that the cut-off wheel has stopped before you access the cutting chamber.

Flushing the Cutting Chamber

To ensure a longer lifetime of your Unitom-2, Struers strongly recommends that you clean the cutting chamber with the cleaning hose every day.


- Point the cleaning hose towards the bottom of the cutting chamber.
- Turn the cooling valve to cleaning position (see sketch in table) and adjust the flow as needed.
- Clean the cutting chamber thoroughly.

NB!


Do not clean the protection guard with the cleaning hose; use a damp cloth instead.

Cooling Valve


The cooling valve for cooling/cleaning is located in the back of the cutting chamber to the left.

Turn the knob counter-clockwise to vertical position for cooling of cut-off wheel.

Turn the knob to regulate the flow while flushing the cutting chamber.


Turn the knob clockwise to vertical position for flushing the cutting chamber with maximum flow.

- Turn the cooling valve to cooling position.
- Leave the protection guard open to let the cutting chamber dry completely.

AVOID RUST!

Leave the protection guard open to let the cutting chamber dry completely.

Reference Guide

Table of Contents	Page
1. Advanced Operations	
Adjusting the Back Stop	13
Clamping Irregular Workpieces	13
Optimising the Cutting Results	14
Connection to External Exhaustion System	15
Mounting an Extra Cutting Table/Quick-Clamping Device	15
Mounting Cutting Table 2.....	15
Mounting the Quick-Clamping Device.....	15
Explanation of Safety Factors	15
2. Accessories	16
3. Consumables	16
4. Trouble-Shooting	17
5. Maintenance	
Daily Service	20
Weekly Service	20
Monthly service	21
Replacing the Cooling Water	21
Lubrication of Bearings for Tilting Frame	21
Maintenance of Cutting Tables	22
Maintenance of Cut-off Wheels	22
Storing of Bakelite Bonded Al ₂ O ₃ Cut-off Wheels	22
Maintenance of Diamond and CBN Cut-off Wheels.....	22
Maintenance of Clamping Devices.....	22
6. Technical Data	23

1. Advanced Operations

Adjusting the Back Stop

Adjust the back stop of the quick-clamping device to align the centre of the workpiece beneath the centre of the cut-off wheel:

- Loosen the back stop.
- Insert the workpiece and adjust the back stop.
- Tighten the back stop.
- If you also need to adjust the quick-clamping device, place it parallel to the back stop and tighten it.

Clamping Irregular Workpieces

Irregular workpieces without plane clamping surfaces must be clamped using special clamping tools, as the workpieces must not move during the cutting. This could result in damage to the cut-off wheel or to the workpiece itself. Use the T-slots to mount the special clamping tools. Struers offers a kit of Clamping Tools (See Accessories).

To achieve faster cutting, orientate the workpiece so that the wheel will cut the smallest possible cross-section.

Optimising the Cutting Results The following table shows possible answers to a number of common questions:

Optimising the Cutting Results	
Question	Answer
How can I avoid discoloration or burning of the sample?	Apply a lower cutting force
	Change the cut-off wheel as the hardness of the present cut-off wheel may be inappropriate for the hardness of the workpiece*)
How can I avoid burrs?	Use a softer cut-off wheel*)
	In case extra cutting table (optional) and extra quick-clamping device (optional) is mounted: Clamp the workpiece securely at the right hand clamping device Tighten the left hand clamping device just enough to prevent the workpiece to shifting on being cut
How can I avoid the cut-off wheels wearing too quickly?	Apply a lower cutting force or use a harder cut-off wheel*)
How can I achieve faster cutting?	Place the workpiece in a direction which allows the cut-off wheel to cut the smallest possible cross-section. Apply a higher cutting force

*) Please refer to the Selection Guide in the [Struers Cut-off Wheels brochure](#).

Connection to External Exhaustion System

Struers recommends the use of exhaustion as the materials to be cut may emit harmful gasses or dust.

On the back of Unitom-2 you will find a joint for a 50 mm dia. exhaustion hose.

- Remove the red cap.
- Mount an exhaustion hose from your local exhaustion system onto the joint.

Mounting an Extra Cutting Table/Quick-Clamping Device

Unitom-2 is delivered with one cutting table with a quick-clamping device, but it is often useful to have an extra cutting table for irregular workpieces.

Mounting Cutting Table 2

- Remove the 4 plastic knobs in the bottom of the cutting chamber.
- The second cutting table must be mirrored compared to the first cutting table.
- Adjust the second cutting table to parallel position and a distance of 6 mm compared to the original cutting table.

Mounting the Quick-Clamping Device

- Clean the cutting table with the flush hose.
- Place the back stop parallel to the back stop of the original cutting table. Use a steel ruler to adjust the position.
- Tighten the back stop.
- Place the quick-clamping device parallel to the quick-clamping device of the original cutting table and tighten it.

Please refer to the [Struers Clamping Tools brochure](#) for details of the range available.

Explanation of Safety Factors

The protection guard has a safety switch to prevent the cut-off wheel from starting while the guard is open. Furthermore, a locking mechanism prevents the opening of the protection guard before the cut-off wheel is standing absolutely still.

2. Accessories

Clamping Tools

Please refer to the [Struers Clamping Tools brochure](#) for details of the range available.

Specification	Cat. No:
Table Unit For Unitom machines. With compartment for recirculation cooling unit.	05116916
Cutting Table with 10 mm T-slots Additional table for support of large workpieces. With exchangeable stainless steel bands. 140 x 315 mm (5.5" x 12.4")	05336904
T-slot cleaner T-slot cleaner for 10 and 12 mm T-slots	05486910
Cooli System 5 Recirculation System With 100 l tank, large pump, Cooli-1 filter bag and adapter plate:	
1 x 100 V / 50 Hz	05766816
1 x 220 V / 50 Hz	05766822
1 x 120 V / 60 Hz	05766823
1 x 240 V / 60 Hz	05766824
1 x 100-120 V / 50-60 Hz CSA	05767016

3. Consumables

Cut-off Wheels

Please refer to the Selection Guide in the [Struers Cut-off Wheels brochure](#).

IMPORTANT
The max. rotational speed of Unitom-2 is 2775 rpm when running idle. The Struers EXO range of cut-off wheels is marked max. 2720 rpm. However, the rotational speed will always be within an acceptable $\pm 5\%$ tolerance.

Other Consumables

Specification	Cat. No:
Corrozip Environment friendly. To protect the machine from corrosion and to improve cutting and cooling qualities. For Recirculation Cooling Unit.	
1 l	449900045
5 l	449900046

4. Trouble-Shooting

Recirculation Cooling Unit

For Messages and Errors relating to the Recirculation Cooling Unit, please refer to the Cooli Instruction Manual.

Error	Explanation	Action
Indications on the Front Panel		
OVERLOAD	Emergency stop active.	Solve the problem causing the emergency stop, then release the emergency stop.
	Cutting motor overloaded due to heavy use.	Open the protection guard and let the cutting motor cool down for about 5 - 10 min.
FLUSH	Lights when cooling valve is in flushing position. Cooling water is stopped.	Turn the knob on the cooling valve counter-clockwise to vertical position for cooling of cut-off wheel.
Machine Problems		
The cut-off wheel does not rotate.	The protection guard open.	Press the safety guard completely down. If this does not help: Call a Struers service technician.
Chamber light does not work.	Fuse F2 or the bulb may have blown.	Replace fuse F2 or the bulb (NB 230V). The F2 fuse is placed behind the back plate.
Water leaking.	Leak in a hose of Recirculation Cooling Unit.	Check the hose and tighten the hose clamp.
	Water overflow in the cooling water tank.	Remove the excess water in the tank.
Workpieces or cutting chamber rusty.	Insufficient additive for cooling fluid.	Add Struers Additive for cooling fluid to the cooling water, using the correct concentration. Check with a refractometer. Follow the instructions in the Section Maintenance.
	The machine is left with closed protection guard.	Leave the protection guard open to let the cutting chamber dry.
Protection guard window blurred.	Insufficient cleaning.	Clean with mild soapy water (do not use the cleaning hose).
The protection guard will not open.	The safety lock has stuck.	Open the lock with a thin flat screw driver. Clean and lubricate the lock. Follow the Weekly service instructions in Section Maintenance to avoid the situation.
Quick-clamping device unable to hold the workpiece.	The quick-clamping device is not in balance.	Adjust the screw underneath the clamping column. Use a 3 mm hexagon key.
	Clamping heart worn.	Call Struers service technician.

*Unitom-2
Instruction Manual*

Error	Explanation	Action
Cutting Problems		
Discoloration or burning of the workpiece.	The hardness of the cut-off wheel is inappropriate for the hardness / dimensions of the workpiece.	See Section Consumables, Cut-off Wheels.
	Inadequate cooling.	- Check that there is enough water in the recirculation cooling unit. - Check the cooling water hoses.
	The force on the cut-off-wheel is too high.	Apply a lower force on the cut-off wheel.
Unwanted burrs.	Wheel too hard.	See Section Consumables, Cut-off Wheels.
	Too high force on the cut-off wheel near the end of the operation.	Reduce the cutting force near the end of the operation.
	Lack of support.	Support with two cutting tables (optional).
The cutting quality differs.	Cooling water hose clogged.	Clean the cooling water hose and the cooling tube. Check the water flow by turning the cooling valve to cleaning position.
	Insufficient cooling water.	Refill tank with water. Remember Struers Additive for cooling fluid.
The cut bends to a side.	The force on the cut-off-wheel is too high.	Apply a lower force on the cut-off wheel.
The cut-off wheel breaks.	Incorrect mounting of the cut-off wheel.	- Check that the hole has the correct diameter. - Check cardboard washer on both sides of the cut-off wheel. The nut must be tightened properly.
	Incorrect clamping of the workpiece.	Make sure that only one of the quick-clamping devices is tight. The other device (optional) should only press lightly. Use support tools (optional) if the geometry of the workpiece makes support necessary.
	Not enough support of the workpiece.	Support with two cutting tables (optional).
	Wheel too hard.	See Section Consumables, Cut-off Wheels.
	The force on the cut-off-wheel is too high.	Apply a lower force on the cut-off wheel.
	Inadequate cooling.	- Check that there is enough water in the recirculation cooling unit. - Check the cooling water hoses.

Unitom-2
Instruction Manual

Error	Explanation	Action
The cut-off wheel wears down too quickly.	The force on the cut-off-wheel is too high.	Apply a lower force on the cut-off wheel.
	The cut-off wheel is too soft for the task.	See Section Consumables, Cut-off Wheels.
	Unitom-2 vibrates (worn bearings).	Call Struers service technician.
The cut-off wheel does not cut through the workpiece.	Incorrect choice of cut-off wheel.	See Section Consumables, Cut-off Wheels.
	Cut-off wheel worn.	Replace the cut-off wheel.
	Diameter of cut-off wheel too small.	Change to a cut-off wheel with another diameter, e.g. $\varnothing 350$ mm.
	The cut-off wheel gets caught in the workpiece.	Support the workpiece and clamp it on both sides of the cut-off wheel in such a way to allow the cut to stay open.
	The stop position is set too low.	See Section Basic Operation, Stopping the cutting.
The workpiece breaks when clamped.	The workpiece is brittle.	Place the workpiece between two polystyrene plates. NB! Always cut brittle workpieces very carefully.
The sample is corroded.	The sample is not resistant to water.	Use a neutral liquid as cooling fluid or cut without using cooling fluid at all. DO NOT USE AN INFLAMMABLE LIQUID (You may lift the pump carefully out of the recirculation cooling unit).
	The sample has been left in the cutting chamber for too long.	Leave the protection guard open, when you leave the machine.
	Insufficient additive for cooling fluid.	Add Struers Additive for cooling fluid to the cooling water in the correct concentration. Check with a refractometer. See Section Maintenance.

5. Maintenance

Recirculation Cooling Unit

For Maintenance of the Recirculation Cooling Unit please refer to the Cooli Instruction Manual.

To ensure a longer lifetime for your Unitom Struers strongly recommends daily cleaning of the cutting chamber. Clean the cutting chamber thoroughly if the Unitom is not to be used for a longer period of time.

Service and greasing of the quick-clamping devices is part of a regular Struers Service Schedule. Struers recommends that a regular service check be carried out after every 1500 hours of use. Please contact a Struers Service Technician to service the machine.

Daily Service

- Clean the cutting chamber, especially the cutting table with the T-slots.
- Clean the protection guard window with a damp cloth. Do not use the cleaning hose.

IMPORTANT

Do **not** clean the lamp glass with alcohol.
Use a damp cloth.

Cooling Unit

- Check the level of the cooling water after 8 hours use or at least every week.

Weekly Service

- Clean the cutting chamber thoroughly.


Monthly service

Replacing the Cooling Water

- Replace the cooling water in the Recirculation Cooling Unit at least once a month.

Lubrication of Bearings for Tilting Frame

- Press emergency stop.
- Open the protection guard.
- Remove the red protection caps from each of the two grease nipples, located as indicated with a white circle in the illustration.

- Lubricate in the two grease nipples with the grease gun and move the spindle/cutting motor up and down by hand, until grease comes out at the shaft. Use water-repelling grease, e.g. Shell Albida Grease LX.
- Clean superfluous grease away and mount the protection caps at the grease nipples.
- Release the emergency stop.

Maintenance of Cutting Tables

The stainless steel rails (available as spare parts) should be replaced in case they have become worn or damaged.

Maintenance of Cut-off Wheels

Storing of Bakelite Bonded Al_2O_3 and SiC Cut-off Wheels

These cut-off wheels are sensitive to humidity. Therefore, do not mix new, dry cut-off wheels with used humid ones. Store the cut-off wheels in a dry place, horizontally on a plane support.

Maintenance of Diamond and CBN Cut-off Wheels

The precision of diamond and CBN cut-off wheels (and thus the cut) depends on how carefully the following instructions are observed:

- Never expose the cut-off wheel to overload, such as heavy mechanical load, or heat.
- Store the cut-off wheel in a dry place, horizontally on a plane support, preferably under light pressure.
- A clean and dry cut-off wheel does not corrode. Therefore, clean and dry the cut-off wheel before storing. If possible, use ordinary detergents for the cleaning.
- Regular dressing of the Diamond and CBN cut-off wheel is also part of the general maintenance.

Maintenance of Clamping Devices

Important

It is recommended to thoroughly clean and lubricate the Quick Clamping Device and Vertical Clamping Device at regular intervals.

6. Technical Data

Subject		Specifications	
		Metric/International	US
Cutting	Three-phased cutting motor Rotational speed (running idle)	4.7kW 2775rpm	6.4HP 2775rpm
Cut-off Wheel	Diameter x thickness	ø350x2.5x32mm ø300x2x32mm	13.8x0.1x1.25" 11.8x0.08x1.25"
Cutting Chamber Geometry	<i>Max. workpiece dimensions:</i> Thickness Width Length <i>Max cutting capacity:</i> Max. cutting thickness Max. length of cut (at a workpiece thickness of:) (The max. cutting capacity might be less, depending on the workpiece material and the cut-off wheel).	150mm 415mm 630mm 120mm 200mm (90mm)	5.9" 16.3" 24.8" 4.7" 7.9" (3.5")
Cutting Table Dimensions	Width Length	286mm 316mm	11.26" 12.44"
Additional Cutting Table	An optional table may be mounted	140x315mm	5.5x12.4"
Dimensions and Weight	<i>Unitom-2</i> Width Depth Height Weight <i>Optional Table Unit</i> Width: Depth: Height:	700mm 750mm 660mm 160kg 900mm 750mm 800mm	27.5" 29.5" 26" 353lbs 35.4" 29.5" 31.5"

*Unitom-2
Instruction Manual*



Subject	Specifications	
Electrical Data	Voltage / frequency:	Max. Load
	3 x 200 V / 50 Hz	21.3A
	3 x 200-210 V / 60 Hz	21.3A
	3 x 220-230 V / 50 Hz	19.8A
	3 x 220-240 V / 60 Hz	19.8A
	3 x 380-415 V / 50 Hz	11.8A
	3 x 460-480 V / 60 Hz	11.8A
Residual Current Circuit Breaker	type A, 30 mA (or higher) is recommended.	
Safety standards	Please refer to the Declaration of Conformity	
Environment	Noise level	Approx. 76dB(A) during cutting, at a distance of 1.0 m/39.4" from the machine.
	Dynamic Balancing	Max. admissible unbalance U_{per} according to ISO 1940/1, Balance Quality Grade G6.03: $U_{per} = 2 \text{ gmm/kg} = 110 \text{ gmm}$.

Quick Reference Guide


Clamping the Workpiece

- Place the workpiece between the clamp and the back stop at the right hand cutting table.
- Push the clamp towards the workpiece and lock the quick-clamping device with the locking handle.

Starting/Stopping the Cutting

- Carefully close the protecting guard.
- Press **START** . The cut-off wheel starts rotating and the cooling water starts running.
- Lower the cut-off wheel onto the workpiece by pulling the cutting handle slowly until contact with the workpiece. Do not apply too heavy cutting force.
- Let the cut-off wheel work itself a little into the workpiece. Increase the pressure and keep it constant while cutting.
- When the cut-off wheel is almost through the workpiece the cutting pressure should be reduced.
- After cutting, return the cut-off wheel to top position. Stop the cut-off wheel and the cooling water by pressing **STOP** . After pressing **STOP** , the protection guard will be locked for about 5 seconds for safety reasons.
- Open the protection guard and release the clamping device.
- Remove the workpiece.

Changing the Cut-off Wheel

- Open the protection guard and swing the right cooling tube to the right.
- Press the black locking knob and turn the cut-off wheel until the spindle lock clicks.
- Remove the nut with a spanner. Remove the flange and the old cut-off wheel.
- Mount the new cut-off wheel. Mount the flange and the nut. Tighten carefully.
- Swing the right cooling tube back into position.

Flushing the Cutting Chamber

- Point the cleaning hose towards the bottom of the cutting chamber.
- Turn the cooling valve clockwise to flushing position.
- Clean the cutting chamber.
- Turn the cooling valve counter-clockwise to cooling position.
- Always leave the protection guard open to let the cutting chamber dry.

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Date of Release FGF0.20FG

Unitom-2
Spare Parts and Diagrams

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2012.

Struers A/S

Pederstrupvej 84

DK-2750 Ballerup

Denmark

Telephone +45 44 600 800

Telefax: +45 44 600 801

Spare Parts and Diagrams

Table of contents

Drawing

Unitom-2

Drawings

Unitom-2, complete	15370001W
Cutting Table, complete	15330003B
Base, assembled.....	15370011B
Cover, assembled	15370050B
Assembled spare parts for cover	15370042A
Lamp, assembled.....	15370044B
Water Valve, complete	15330025F
Spindle unit, complete.....	15370015E
High power plate, rear side	15370030F

Diagrams

Main Circuit Diagram.....	15373100D
Power Circuit Diagram	15373110C
Control In/Out Circuit Diagram	15373120C
Circuit Diagram, Logic control circuit.....	15373200D
Wiring of Power connection diagram	15373420E
Wiring diagram, Back-Panel, Control	15373410E

Some of the drawings may contain position numbers
not used in connection with this manual.




Unitom-2
Spare Parts and Diagrams

Spare Parts list, Unitom-2


Drawing	Pos.	Spare Part	Cat no:
15370001		Unitom-2, complete	
	170	50Hz timing belt 600-8MP-30mm	2JT20600
	170	60Hz timing belt 560-8MP-30mm	2JT20560
	190	Gas spring 700N Lift_O-Mat S100 with bellows	RYS00190
	310	1.5 AT Fuse Mini TR5	2FU32950
	310	Unitom-2 PCB, complete	15373001
	382	Gas Spring 400 N CAMLOC SL=150 (2 pcs.)	RYS00370
	660	Fixed flange with key	R4981210
	661	Cardboard blotter ø32/105 f. xxUNI (25 pcs)	15332902
	665	Key A 4x4x10	2ZF10410
	670	Loose flange	R4981209
	680	Nut M20	2TA10200
	690	Disc ø21/37/3 A2 DIN 125	2ZA10020
	770	Transparent silicon sealant for window glass	2LL00330
	800	Cooling pump (3x220-240V/50Hz)	381MP408A
	800	Cooling pump (3x380-440V/50Hz)	381MP408B
	800	Cooling pump (3x200-220V/60Hz)	381MP409A
	800	Cooling pump (3x380-480V/60Hz)	381MP409B
	801	Sieve for 2AB4 (ø99)	14980074
	15330003		Cutting table, complete
10		Cutting table, complete	05336904
	20	Replacement steel bands, 3 pcs	05336905
15370011		Base, assembled	
	20	Machine supports	14980015
	40	Gasket for flange	15330150
	45	Flange for exhaust (not shown)	13811000
	200	Water Valve, complete without knob	15330025
	210	Knob for Cooling Water valve with grease	15332904
	220	Safety Switch, T3C 236-O2Z	15330120
	260	Sealing ring JS16	2IF03713
290	Flushing hose, complete	15330075	
15370050		Cover, assembled	
	20	Window with serigraphy, right	15371401
	30	Window with serigraphy, left	15371402
	90	Lamp, assembled	15370044
	120	Assembled parts for cover	15370042
130	Universal Silicone 512. Transparent	2LL00330	

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























Pos. 30 High vacuum grease
 * 50H seal
 > wire length 600 mm.

D	22.05.2007	Solid works drawing, translated to english.	SPE	22.05.2007	CJE
C	03.12.1999	100M seal ->50H	BMJ	03.12.1999	JTV
B	17.09.1997	Pos.80, 90, 100 tilf. seal + led.lgd., fedt tilf.	BMJ	17.09.1997	JTV
A	16.01.1997		BMJ		
Rev	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init

Pederstrupvej 84
 DK-2750 Ballerup
 Copenhagen
 Denmark
 Phone: +45 44600 800
 Fax: +45 44600 804

Material:

Scale:

1:1

Format:

A4

Tolerance: DS/ISO 2768 - mK

Surface treat.: None


ID:

Description:


15330025 Valve body, assembled

Rev:


D



SECTION C-C
SCALE 1 : 1


DETAIL A
SCALE 2 : 2.5

DETAIL B
SCALE 2 : 2.5


E	24.05.2007	SW drawing, translated to english	SPE	24.05.2007	C.JE
D	28.11.2005	Tekst tilføjet for dæksel ved detalie af leje	BRY	28.11.2005	JTV
C					
B					
A	04.06.1997		SOR	02.07.1997	JTV
Revision	Crea. date	Revision description	Draw. Init	Appr. date	Appr. Init
	dd-mm-yy			dd-mm-yy	
		Material:	Scale: 1:2.5	Format: A2	Tolerance: DS/ISO 2768- mK
		ID:	Description:	Surface treat.: None	Rev:
					15370015 Tilting arrangement for Uinitom-2
					E

The back side of the sealing must be cut


F	2012.05.25	Pos.140 2XK46199-> Pos.130 2XK46198	SPE	2012.05.25	
A	10.06.1997		BMJ	02.07.1997	JTV
Revision	Crea. date yyy-mm-dd	Revision description	Draw. Init	Appr. date yyy-mm-dd	Appr. Init
		Material:	Scale: 1:2	Format: A3	Tolerance: DS/ISO 2768- mK Weight : g
ID: 15370030 Description: Electrical back plate, assembled		Rev: F			

Pederstrupvej 84
 DK-2750 Ballerup/Copenhagen
 Denmark
 Phone: +45 44 600 800
 Fax: +45 44 600 804

M1
CUTTING MOTOR

M2
RECIRCULATION MOTOR

CUTTING MOTOR M1,
NOT IN OVERLOAD.

RECIRCULATING PUMP
NOT IN OVERLOAD

RECIRCULATING PUMP ON.

CUTTING MOTOR RUN

SAFETY GUARD LOCK

START
PROCES ON

STOP
PROCES OFF

WATER VALVE IN
COOLING POSITION.

WATER VALVE IN
SLICING POSITION.

EMERGENCY STOP
AND
OVERLOAD AT

RECIRCULATION -
OR
CUTTING MOTOR

SAFETY GUARD
CLOSED.


SAFETY GUARD
RELEASE

B01
NAMUR

B02
NAMUR

S4
SAFETY SWITCH

YM1
AKTUATOR SWITCH

15373110 POWER CIRCUIT DIAGRAM


15373200 PCB CIRCUIT DIAGRAM

15373410 WIRING DIAGRAM, CONTROL, BACK-PANEL


15373420 WIRING DIAGRAM, POWER, BACK-PANEL

15373430 WIRING DIAGRAM, CONTROL, FRONT-PANEL

Rev. C: Test change 1 YW1: 5V + 5V F1: F2, G, R+D.3 F2: F2, G, R+D.2 B1: B3, G, R+D.2 S2: S4, G, R+D.2 G: G, R+D.2 Rev. D: (F1) 13, G, R+D Ref. G1 changed to QS1 Colors on BQ14802 changed G.R. = Gnd Referenced Friday, December 17, 2004	Rev. A: COLOR CHANGE. YW1: 5V + 5V F1: F2, G, R+D.3 F2: F2, G, R+D.2 B1: B3, G, R+D.2 S2: S4, G, R+D.2 G: G, R+D.2 Rev. B: X2.4: BU + X2.5 X2.5: RD + X2.4	STROBERE AG MAIN CIRCUIT ASSEMBLY: 178 GRC-21113 PRODUKTION LEHRMANN PRODUKTION: 2004: 28000 MAIN CIRCUIT DIAGRAM: UNITOM - 2: Size: A1 SCALE: 12 mg 97 / Dan Molding Sheet: 1 of 3
---	---	---

Rev.: A. Text change	STRUERS A/S VALHOESALLE 176 DK-2610 ROEDOVRE DENMARK PHONE: +45 3670 3500		
Rev.: B.: Text change ! F1 -> F2, G.R.->J,4 F2 -> F3, G.R.->J,3	UNITOM -2: POWER CIRCUIT DIAGRAM		
Rev.: C. Ref. Q1 changed to QS1	Size A2	CAGE Code	DWG NO 15373110
Monday, December 13, 2004	Scale	D 16.06.97 / DEM	Sheet 1 of 1

SEE
DRAWING
NO.
15373110

Rev.: B.: Text change 1
S1 -> S3, G.R.=E.2
S2 -> S4, G.R.=F.2,G.2
HL** -> Wiped off.

Rev.: C.: Colors on
BQ1 & BQ2 changed

Rev.: A.: COLOR
CHANGE:
YM1: 12V -> 24V,
YM1: RD -> BK, WH -> BN,
BN -> VT,
YM1: GN -> BU, BU -> RD,
X5:1 -> -5, -5 -> -1.

Friday, December 17, 2004


STRUJERS A/S
VALHOESALLE 176
DK-2810 ROEDDovre
DENMARK
PHONE: +45 3670 3500

Control In / Out Circuit Diagram UNITOM - 2:


Size	A2	CAGE Code	DWG NO	Rev
			15373120	I C
Scale		D.16.06.97 / DEM	Sheet	1 of 1

X1:

Power Supply
24V AC
Chassis


POWER ON
GREEN

X4:

Signal In:
SAFETY
GUARD CLOSED
ACTUATOR:
"SAFETY GUARD"
IN POS. LOCK
OVERLOAD, EMERG.
CUT. MOTOR
REC. MOTOR
EMERGENCY STOP

X2:

24V AC,
TO CONTROL CIRCUIT I:
IN I 24V AC,
FROM CONTROL CIRCUIT.
CUTTING MOTOR ON.
IN II 24V AC,
FROM CONTROL CIRCUIT.
RECIRCULATING PUMP ON.
N24V

ALARM FOR
OVERLOAD OR EMERGENCY STOP


SLICING ON

X5:

Signal In:
COOLING
VALVE READY
SLICING
VALVE READY

X3:

+24V DC OUT.
"SAFETY GUARD"
RELEASE
AKTUATOR:
"SAFETY GUARD"
LOCK
GND-P



Note A:
 Cooling "on board"!
 R15 24V - 9K1
 C21 100uF/450V
 R22 56K - 402K
 Rev: B New SIL 1620V
 Rev: C Label "START" Change
 Rev: D of 17th Dec 2004
 (Rev. changed from 242 to 220R)

Note B:
 Cooling "on board"!
 Cu area 25cm x
 2.5cm to GND.

Rev: A RE 9522 - 000	Rev: B 03	Rev: C 01	Rev: D 01
R15 24V - 9K1	Varidex 100 178	DA 212-Resistor	DA 212-Resistor
C21 100uF/450V	DA 212-Resistor	DA 212-Resistor	DA 212-Resistor
R22 56K - 402K	DA 212-Resistor	DA 212-Resistor	DA 212-Resistor
Rev: B New SIL 1620V	DA 212-Resistor	DA 212-Resistor	DA 212-Resistor
Rev: C Label "START" Change	DA 212-Resistor	DA 212-Resistor	DA 212-Resistor
Rev: D of 17th Dec 2004	DA 212-Resistor	DA 212-Resistor	DA 212-Resistor
(Rev. changed from 242 to 220R)	DA 212-Resistor	DA 212-Resistor	DA 212-Resistor
Circuit Diagram	DA 212-Resistor	DA 212-Resistor	DA 212-Resistor

Unimot - 2
 Logik Control

1573200
 Tuesday, August 10, 2004
 6:17 Dec 96 / OEM

SEE CONNECTIONS DIAGRAM 14983260


Must be short circuit protected with external fuses according to local regulations.
 I_{k,min.} = 145 Amp. in 5 sec.
 Fuses max.:
 200 - 240 Volt. 25 Amp. type gl.
 380 - 480 Volt. 20 Amp. type gl.

WIRE-GAUGE
 200-220V - 3,3mm²/12AWG
 380-480V - 2,0mm²/14AWG







COLOR CODES:
 BK = BLACK
 BN = BROWN
 RD = RED
 OG = ORANGE
 YE = YELLOW
 GN = GREEN
 BU = BLUE
 VT = VIOLET
 GY = GREY
 WH = WHITE

Rev.: E: Doc name change (402 13307)	Rev. C: Text Change, 380-480V 25Amp => 20Amp. 2,0mm ² /14AWG => 2,0mm ² /14AWG. Rev.: B: Add Text "k min." M1 Thermo Switch Symbol NO->NC. Rev.: A: Text change	STRØJERS A/S VALHOEJSALLE 176 DK-2810 ROEDOVRE DENMARK PHONE: +45 3670 3500
G.R. = Grid Reference!	Size: A2	CAGE Code: DWG NO: 15373420
Thursday, September 13, 2007	Scale: 20 JUN 87 /	Sheet: 1 of 1


WIRING DIAGRAM:
UNITOM - 2 (without Cooli conn and safety lock):
 Wiring of Power connection.

SEE CONNECTIONS DIAGRAM 14993200 14993200

TO THE FRONT PANEL:

ALL WIRINGS 0,75mm². - EXCEPT OTHERWISE MARKED.

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Rev.: E.: Doc name change! (AKR 13.9.07)

STRØJERS A/S
 VALHOEJSALLE 176
 DK-2810 ROEDOVRE
 DENMARK
 PHONE: +45 3670 3500

WIRING DIAGRAM:
 UNITOM - 2 (without Cooli conn. and safety lock):
 BACK-PANEL, CONTROL.

G.R. = Grid Reference	Size A2	CAGE Code	DWG NO	15373410	Rev. I E
Thursday, September 13, 2007	Scale	19 JUN. 97 / ZJEM	Sheet 1 of 1		

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Herewith declares that

<i>Product Name:</i>	Unitom-2
<i>Type No:</i>	537
<i>Machine Type:</i>	Cut-off machine

is in conformity with the provisions of the following directives:

Safety of Machinery 2006/42/EC according to the following standard(s):
EN ISO 12100:2010, EN ISO 13850:2006, EN60204-1:2006/AC2010.

EMC-Directive 2004/108/EC according to the following standard(s):
EN 61000-6-1:2007, EN 61000-6-3:2007+A1:2011, EN 61000-6-4:2007/A1:2011.

Low Voltage Directive 2006/95/EC according to the following standard(s):
EN60204-1:2006/AC2010.

RoHS 2011/65/EU according to the following standard(s):
EN 50581:2012.

Supplementary Information The equipment complies with the American standards:
UL508.

The above has been declared according to the global method, module A

Date: 21.01.2013

Christian Skjold Heyde,
Vice President, R & D and Production, Struers A/S

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
Unitom-2

**Manual cutting
machine for large
and irregular
workpieces**


The Unitom-2 is a powerful cutting machine, designed for 350 mm cut-off wheels. Unitom-2 is based on proven, superior cutting technology. The concept of high quality and high capacity combined with ease-of-use has been extremely successful.

Features

- Strong motor for high cutting capacity
- Handles large and irregular workpieces
- Easy operation and maximum safety
- The cutting chamber has a built-in, adjustable flush hose making daily cleaning easy and efficient

A sturdy workhorse

The cutting table is provided with 10 mm T-slots. This allows difficult workpieces to be cut using various clamping tools, here 10 mm turntable

The Unitom-2 is a powerful manual cutting machine for a wide range of applications.

Unitom-2 will enable the user to cut a large number of samples, hour after hour. All parts of the machine subject to stress have been put through a careful life-time test. The rotating parts of the machine have been dynamically balanced to ensure stability, and to maximize the working life of your Unitom-2.

All metal parts in the cutting chamber are made of either stainless steel or surface protected materials, practically eliminating corrosion. Along the travel of the cut-off wheel the cutting table is equipped with exchangeable stainless steel bars.

Specifications

- Cutting capacity: Up to 120 mm (4.7") diameter
- Cutting power: 4.7 kW (6.4 HP)
- Cut-off wheel: Max 350 mm diameter (13.7")

User-friendly design

Unitom-2 is very easy to operate making it suitable for laboratory as well as production use.

Easy access

Extra large windows and fully, indirect lightning, give an excellent, non-glare view of the cutting chamber. A balanced mechanism ensures effortless opening and closing of the protection guard.

Easy cleaning

A powerful lamp lights up the cutting chamber, and a flush hose with adjustable

Built-in convenience: Flush hose for cleaning of the cutting chamber

flow takes care of fast and efficient daily cleaning.

Simple installation

Installation and levelling of the machine are made simple, by virtue of four height adjustable vibration dampers.

Cutting of large workpieces

If the user wants to cut large or protruding workpieces, it is possible to replace one or both of the side walls on the machine with optional rubber curtains.

Safety features

Unitom-2 complies with international safety standards, and safety features include emergency stop and cut-off switch for cutting motor and all machine functions. The protection guard cannot be opened before the cut-off wheel has

Built-in safety: Cut-off switch for cutting motor etc.


A third additional cutting table can be added

stopped. The windows of the protection guard are made of impact resistant material. Unitom-2 is prepared for connection to an exhaust system, as some cutting materials may emit harmful gases or irritating fumes.

Accessories

An extensive range of optional accessories is available for Unitom-2. A combined table and cabinet is available for bringing the cutting machine into an ideal working height. Quick-clamping tools and vertical clamping tools ensure full flexibility in securing workpieces. Please see the back of the brochure for more information about accessories.


Vertical clamping tools

Quick clamping tool for T-slots 10 mm

Quick clamping tool with optional jaw blocks


Technical Data

Subject	Specification	Metric/ International	US
CUTTING SPECIFICATIONS			
Workpiece Dimensions (Max)	<i>Workpiece inside cutting chamber:</i>		
	Height	120 mm	5.9"
	Width	415 mm	16.3"
	Depth	630 mm	24.8"
	<i>Workpiece protruding cutting chamber:</i>		
	Height	120 mm	5.9"
Cutting Capacity*	Max. cutting diameter	120 mm	4.7"
	<i>Rectangular sample**:</i>		
	Height	90 mm	3.5"
	Depth	200 mm	7.9"
PHYSICAL SPECIFICATIONS			
Cutting Motor	Cutting power	4.7 kW	6.4 HP
	Max power (depending on voltage)	5.7 - 6.7 kW	7.6 - 9.0 HP
Cut-off Wheel	Diameter x Thickness x Centre-hole	350 x 2.5 x 32 mm	13.7 x 0.09 x 1.25"
	Rotational speed (running idle)	2775 rpm	2775 rpm
Cutting table	Width	286 mm	11.26"
	Depth	316 mm	12.44"
	T-slots	10 mm	0.39"
Dimensions and Weight	Height	660 mm	26"
	Width	700 mm	27.5"
	Depth	750 mm	29.5"
	Weight	162 kg	353 lbs
	<i>Optional Table Unit (LABUN):</i>		
	Height	795 mm	31.3"
Recirculation Cooling Unit	Tank volume	65 l	17 gallons
	Approx. flow	12 l/min	3.2 gallons/min
ENVIRONMENTS			
	Noise level	Approx. 73 dB(A) measured at idle running, at a distance of 1.0 m/39.4" from the machine.	

* Projected cutting capacity with new cut-off wheel. Actual cutting capacity depends on sample material, cut-off wheel and clamping technique

** Example of cutting capacity. Please refer to Cutting Diagram above

Struers cut-off wheels are especially designed for Struers machines, taking into consideration the most recent developments in wet cutting techniques.

Cut-off wheels

Unitom-2 is designed for cut-off wheels of 300 and 350 mm diameter. For normal use, the 300 mm wheels (with a thickness of 2 mm) are recommended. For large workpieces, the 350 mm wheels (with a thickness of 2.5 mm) are used, allowing Unitom-2 to cut samples up to 120 mm diameter. A wide range of Struers cut-off wheels is available for cutting any material.

Struers' equipment is in conformity with the provisions of the applicable International Directives and their appurtenant Standards. (Please contact your local supplier for details)

Specifications

Unitom-2

Manual cutting machine.
4.7 kW (6.4 HP) motor. For 350/300 mm (14/12") cut-off wheels. Including two cutting tables with 10 mm T-slots (05336904). Recirculation cooling unit (05766xxx), table unit (05116916), quick clamping tools (05876909+0576910) and vertical clamping tool (05116905) has to be ordered separately.

3 x 200 V / 50 Hz	05376729
3 x 200-210 V / 60 Hz	05376730
3 x 220-230 V / 50 Hz	05376735
3 x 220-240 V / 60 Hz	05376736
3 x 380-415 V / 50 Hz	05376746
3 x 460-480 V / 60 Hz	05376754

Cat. no.

Accessories

Table Unit

For support of Unitom-2. With compartment for recirculation cooling unit.
Width: 900 mm, depth: 750 mm, height: 800 mm.

05116916

Cooling System 3

With 50 l tank, small pump, Cooli-1, 50 l static filter with filter paper.

05766xxx

Cutting table for 10 mm T-slots

Additional table for support of large workpieces. With exchangeable stainless steel bands. 140 x 315 mm (5.5" x 12.4")

05336904

Replacement Stainless Steel Bands

For Cutting Table with 10 mm T-slots (05336904). Set of three.

05336905

Flange Set, 119 mm dia.

Alternative to 100 mm standard flange. For use with 350 mm cut-off wheels on Unitom-2.

04986901

Clamping tools

Vertical clamping tool for 10 mm T-slots

For clamping irregular shaped workpieces. Complete with operating key and one exchangeable flat swivel shoe. Max clamping height: 58 mm.

05116905

Quick-clamping tool for 10 mm T-slots, left

For securing the workpiece. Complete with back stop. To be mounted on a cutting table placed at the left hand side of the cut-off wheel.

05876909

Quick-clamping tool for 10 mm T-slots, right

For securing the workpiece. Complete with back stop. To be mounted on a cutting table placed at the right hand side of the cut-off wheel.

05876910

Riser Block for 10 mm T-slots

Increases the clamping height of the 10 mm Vertical clamping tool (05116905) by 60 mm.

05116910

Swivel Shoes for Vertical clamping tool

Set of 4 multi-shaped swivel shoes for use with 10 mm Vertical clamping tool (05116905).

05116911

Rubber-coated jaw blocks for 10 mm T-slots

Set of 2 rubber-coated jaw blocks for optimum clamping of e.g. thin-walled tubes.
For use with Quick Clamping tool for 10 mm T-slots, left (05876909) and right (05876910).

05876915

Grooved jaw blocks for 10 mm T-slots

Set of 2 jaw blocks for optimum clamping of regular workpieces.
For use with Quick clamping tool for 10 mm T-slots, left (05876909) and right (05876910).

05876917

Prism type jaw blocks for 10 mm T-slots

Set of 2 jaw blocks for optimum clamping of round and oval workpieces.
For use with Quick Clamping tool for 10 mm T-slots, left (05876909) and right (05876910).

05876913

T-slot turntable for 10 mm T-slots

0-90 degree turntable for angular cutting. For use with various 10 mm T-slot clamping tools. Table height 42 mm.

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CANADA Struers Ltd.

7275 West Credit Avenue
Mississauga, Ontario L5N 5M9
Phone +1 905-814-8855
Fax +1 905-814-1440
info@struers.com

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Struers Inc.
Struers GmbH