

AbraPlan-20



Instruction Manual

Manual No.: 15897001

Date of Release 2H€ .201H



AbraPlan-20
Instruction Manual

Table of Contents	Page
User's Guide	1
Reference Guide.....	45
Quick Reference Guide	56

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the *Serial No.* and *Voltage* on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

Original instructions. The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 201H

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801



AbraPlan-20

Safety Precaution Sheet

To be read carefully before use

1. The operator should be fully aware of the use of the machine according to the Instruction Manual.
2. The machine must be installed in compliance with local safety regulations.
3. Be aware that the machine's centre of gravity is located in the upper half of the machine.
4. Before lifting the machine, ensure that the supplied lifting boom is securely bolted to the machine.
5. When lifting the machine using a forklift, lift from front - never lift the machine from the side or the rear.
6. When lifting the machine using lifting straps, ensure that the straps are crossed and do not press on the sides of the machine.
7. The machine must be placed on a safe and stable support, which is capable of bearing the weight of this machine. Before using the machine, it must be levelled using the adjustable legs.
8. Be sure that the actual voltage corresponds to the voltage stated on the side of the machine and that the connections comply with local regulations. The machine must be earthed.
9. Be aware that when the machine is connected to a compressed air supply the specimen holder arm moves upwards.
10. Make sure that the specimens in the specimen holder are securely fixed, and ensure that the securing screws are not sticking out.

AbraPlan-20
Instruction Manual

- 11.** If you observe malfunctions or hear unusual noises - stop the machine and call technical service.
- 12.** The machine must be disconnected from the mains supply prior to any service.
- 13.** To achieve maximum safety and lifetime of the machine, use only original Struers consumables.
- 14.** Use of working gloves is recommended when changing the grinding stone/ diamond grinding disc.

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.)

User's Guide

Table of Contents	Page
1. Getting Started	
Checking the Contents of the Packing Box.....	2
Unpacking and Placing AbraPlan-20	2
Getting Acquainted with AbraPlan-20	3
Noise Level	3
Supplying Power	4
Direction of Rotation	4
Supplying Compressed Air.....	4
Connection to an External Exhaust System.....	4
Connecting the Recirculation Cooling Unit (accessory)	5
Mounting the Grinding Stone/ Diamond Grinding Disc.....	6
Flushing Unit	13
Cleaning the Inside of the Safety Guard	13
2. Basic Operations	
Front Panel.....	14
Front Panel Controls	15
Software Settings	18
Setting the Language.....	19
Reading the Display	21
Changing/Editing Values	23
Numeric Values.....	23
Text Values	24
Setting up the Software	26
Inserting/Removing the Specimen Holder.....	28
Inserting the Specimen Holder.....	28
Removing the Specimen Holder	28
Grinding Setup	29
Setting the Process Time.....	30
Setting Removal.....	31
Setting Removal/Time.....	32
Cooling Water	32
Starting the Preparation Process	33
Stopping the Preparation Process	34
Dressing Functions and Changing the Grinding Stone	35
Setting the Dressing Parameters	40
Dressing the Diamond Grinding Disc	42
3. Maintenance	
Daily Service	43
Weekly Service	43
Checking the Recirculation Cooling Unit.....	44
Monthly Service.....	44
Replacing the Cooling Water	44
Yearly Service	44
Inspection of Cover	44

1. Getting Started

Checking the Contents of the Packing Box

AbraPlan-20

In the packing box you should find the following parts:

- 1 AbraPlan-20
(Automatic, high-capacity machine for fast and efficient plane grinding)
- 1 Outlet hose 0.8 m for connection to internal cooling unit (factory mounted)
- 1 Outlet hose 2.5 m for connection to external cooling unit
- 1 Inlet hose (factory mounted)
- 2 Hose clamps, 17 mm
- 1 Drain angle, 90°
- 1 Hose clamp, 35-60 mm
- 1 Hose for compressed air
- 2 Hose clamps, 12 mm
- 1 Hose connection for compressed air
- 1 Rubber disc 350 mm
- 1 Flange
- 1 Bolt M12
- 1 Allen key 8 mm
- 1 Set of Instruction Manuals

Unpacking and Placing AbraPlan-20

AbraPlan-20 should be placed directly on a plane and horizontal floor, capable of bearing the weight of this machine, please refer to the Technical Data section.

- Unscrew the nuts from the transport brackets fixing the machine to the pallet.
- Lift the machine from the pallet by means of a forklift truck from the front, and place in a suitable location.
- Remove the safety-springs from the front crossbar, and remove the bar.


Important

If necessary, turn the adjustable feet so that the machine stands firmly and is level.


**Getting Acquainted with
 AbraPlan-20**

Front view

Take a moment to familiarise yourself with the location and names of the AbraPlan-20 components.



Side views




- ⑧ Mains power switch
- ⑨ Electrical connections
- ⑩ Compressed air inlet
- ⑪ Connection to exhaust
- ⑫ Connector for Recirculation Cooling Unit
- ⑬ Opening for water outlet hose

Noise Level

Approx. 68 dB (A) measured during idle running, at the operator's position in front of the machine.

Supplying Power



Direction of Rotation

IMPORTANT

Check that the mains supply voltage corresponds to the voltage stated on the Type Plate (located under the mains switch on the side of the machine).
If the machine is already connected to a mains supply, disconnect this supply before removing the Allen screws securing the electrical panel.

- Remove the 4 Allen screws, on the right hand side of the machine, securing the electrical panel, and let the electrical panel rest on the tabs.
- Lead the cable through the conduit in the panel and connect the 3 phases and earth according to local regulations.

Check that when the power is turned on the grinding stone/ diamond grinding disc rotates counter-clockwise.
If this is not the case, switch off AbraPlan-20 and unplug the machine.

- Change two of the phases.
- Repeat the rotation check.

Supplying Compressed Air

- Connect the compressed air supply to the inlet, located at the rear of the left-hand side of the machine, using the air hose and the hose connection delivered with the machine.
- Fasten the air hose with a hose clamp.

The pressure supply should be 6-10 bar and can be supplied either from a central compressor, a portable compressor with a compressed air reservoir, or a compressed-air bottle. A capacity of 20 l/min at atmospheric pressure is sufficient.

Connection to an External Exhaust System

- Connect a 50 mm pipe to the outlet at the rear of the machine, on the left and connect to the exhaust system.
Recommended capacity for exhaust system: 180m³/h at 0mm water gauge.


Connecting the Recirculation Cooling Unit (accessory)

Please refer to the manual for the Recirculation Cooling Unit for details.

IMPORTANT

Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the Additive container). Remember to top up with Struers Additive each time you refill with water.


- Lead the water outlet tube down into the Recirculation Cooling Unit. Whenever the recirculation tank is removed from the compartment underneath the AbraPlan-20, place the tube in the tube holder to avoid dripping. Remember to put the tube back into the tank before starting the machine.



Tube holder




IMPORTANT

When AbraPlan-20 is connected to an external cooling unit the short outlet tube has to be replaced with the long tube. The long tube is then guided through the outlet hole on the left-hand side of the machine (see picture) and connected to the cooling unit. Shorten the hose to the required length. The inlet hose is led underneath the machine from the pump to the water inlet.







Outlet tube

**Mounting the Grinding Stone/
 Diamond Grinding Disc**

GRINDING SETUP	
Grinding mode:	Removal
Stock removal:	200 µm
Force:	250 N
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;"> F1  <small>Water ON</small> </div> <div style="text-align: center;"> F2  <small>One dress</small> </div> <div style="text-align: center;"> F3  <small>Dressing setup</small> </div> </div>	


- From the grinding set-up screen, press **F3** to go to the Dressing Set-up.

DRESSING	
Dresser step:	40 µm
Dresser speed [1-low...100-high]:	30
Automatic dressing during process:	Yes
Automatic dressing after process:	Yes
Dressing mode:	Removal
Dresser sensitivity:	60 %
Dressing during grinding:	No
Remaining height of stone:	9.89 mm
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;"> F1  <small>Finding TOP of stone</small> </div> <div style="text-align: center;"> F2  <small>One dress</small> </div> <div style="text-align: center;"> F3  <small>Up</small> </div> <div style="text-align: center;"> F4  <small>Change Grinding Stone</small> </div> </div>	

- Press **F4** Change Grinding Stone to move the dresser to the top position.


AbraPlan-20
Instruction Manual

- ① Flushing unit (head)
- ② Dressing arm
- ③ Stone Guard
- ④ Securing finger screws



- Ensure that the sample motor is fully raised, and lift the safety guard to gain access to the grinding area.
- Lift the flushing unit head and, using its built-in magnet, attach it to the back of the grinding chamber.
- Unscrew the 2 finger screws.
- Carefully lift the stone guard up and away to the front of the machine.

- Assemble the grinding stone/ diamond grinding disc on the motor flange as illustrated:



- ① Motor Flange
- ② Rubber disk
- ③ Grinding stone/ diamond grinding disc
- ④ Securing flange and cardboard washer
- ⑤ Securing bolt

IMPORTANT

Make sure that the grinding stone/diamond grinding disc is intact. The stone/disc must be dry when mounted and the flange should be clean and smooth.

- Mount the bolt and fasten firmly using an 8 mm Allen key.



IMPORTANT

Do not over-tighten the securing bolt as this may damage the grinding stone/diamond grinding disc. The bolt should be tightened with a force of minimum 8 Nm (5.9 lbf-ft), maximum 10 Nm (7.4 lbf-ft).

Warning

Take care of sharp or rough edges whilst securing/ removing the grinding stone/diamond grinding disc.

- Replace the stone cover and secure it with the 2 nuts.
- Re-attach the flushing unit head.
- Lower the safety guard and press F1 to configure the automatic dresser.




- Press F1 if you have installed a grinding stone or press F2 for a diamond grinding disc.


IMPORTANT

Pressing F1 when a diamond grinding disc has been inserted can result in severe damage to both the dressing tool and the diamond grinding disc.

When a grinding stone has been inserted the following screen will appear:



- Press F1 to continue and start the calibration process. The dresser checks the height of the grinding stone in two positions, one close to the centre, and the other at the periphery. Wherever the stone is highest, the following fine adjustment will start.



- Press F1 to continue with the fine adjustment.

FINE ADJUSTMENT OF DRESSER

Vertical dresser position (encoder units): -5000
Horizontal dresser position (encoder units): +350

Close the hood and press F1 to continue. The stone will start rotating.

Press Esc if you do not want to use the guide for the following procedure.



- Press F1 to continue with the fine adjustment.

FINE ADJUSTMENT OF DRESSER

Vertical dresser position (encoder units): -5000
Horizontal dresser position (encoder units): +350

Repeat pressing Enter/(F4) until the dresser has reached the top of the stone.

Each Enter/(F4) activation moves dresser down by 40 μm .


Each F3 activation moves the dresser up by 20 μm .

Turn the knob to make a small left-right movement between each Enter activation for better recognition of the first contact of the dresser tip with the stone.

Press F1 as soon as the dresser tip touches the stone, or press Esc to interrupt this process.




- Press Enter or F4 to move the dresser downwards in steps of 40 μm .
- Turn the knob to move the dresser sideways to ensure that the dresser has touched the grinding stone.
- Press F1 to finish the fine adjustment.



- Press F1 again to store the current dresser position as the top of the grinding stone.

Flushing Unit



When attached to its mounting (as illustrated), the flushing unit supplies water/coolant to the grinding stone/ diamond grinding disc during the grinding and dressing processes. However, the flushing unit can be removed from the mounting and used to hose down the grinding area and samples etc. To do this:

- Squeeze in and hold the clamp buttons (illustrated by arrows) to cut-off the water/coolant flow.
- Press **F1** to start the pump.
- Lift the flushing unit free of its mounting and pull out the tube.
- Direct the flushing unit in the desired direction and release the clamp buttons.
- Press **F1** to stop the pump.
- Push the tube back into the opening at the back of the grinding chamber and reinsert the flushing unit in its mounting.

If the grinding/dressing process is not active, a water/coolant flow can be achieved by pressing **F1** when the GRINDING or DRESSING Menus are displayed.

Cleaning the Inside of the Safety Guard




To clean the inside of the safety guard use the built-in magnet of the flushing unit head to attach it to the back of the grinding chamber.










- Lift the flushing unit free of its mounting and pull out the tube.
- Use the magnet to attach the flushing unit to the back of the grinding chamber.
- Close the safety guard.
- Press **F1** to start the pump.
- When the safety guard has been flushed sufficiently press **F1** to stop the pump.
- Open the safety guard.
- Push the tube back into the opening at the back of the grinding chamber and reinsert the flushing unit in its mounting.

2. Basic Operations


Front Panel




Front Panel Controls

Name	Key	Function	Name	Key	Function
FUNCTION KEY		Controls for various purposes. See the bottom line of the individual screens.	Push/Turn Knob		Used for entering and changing parameters. Combined cursor and enter key.
FUNCTION KEY		Controls for various purposes. See the bottom line of the individual screens.	START		Starts the preparation process
FUNCTION KEY		Controls for various purposes. See the bottom line of the individual screens.	STOP		Stops the preparation process.
FUNCTION KEY		Controls for various purposes. See the bottom line of the individual screens.	EMERGENCY STOP		The EMERGENCY STOP is located on the front of the machine. - Push the red button to stop. - Pull the red button to release.
Esc		Leaves the present menu or aborts functions/changes.	MAIN SWITCH		The main switch is located on the right side of the machine.


Switch on the power at the main switch located at the right hand side of the machine. The following display will appear briefly:



AbraPlan-20 will then search for the reference position of the dresser, showing the following screen:



When the reference position has been found, the next screen will appear for a short while:




Software Settings

The display will then change to the last screen shown before AbraPlan-20 was switched off.

When switching on AbraPlan-20 for the first time, a message will appear to request selection of the language of your choice.

The MAIN MENU display will then appear. If the heading in the display is different, press **Esc**, until the MAIN MENU appears.

The MAIN MENU is the highest level in the menu structure. From this menu, you can enter the configuration menu and grinding process menu.



Setting the Language



To change the language, carry out the following steps:




Turn knob to select CONFIGURATION.




Push knob to activate the CONFIGURATION Menu.




 Turn knob to select Language.




 Push knob to activate the LANGUAGE pop-up menu.



 Turn knob to select the language you prefer.




 Push knob to accept the language.

The CONFIGURATION Menu now appears in the language you have chosen.

 Press **Esc** to return to the MAIN MENU.


Reading the Display

The display is primarily divided into 3 areas. The position of these areas and the information they contain are explained in the illustration below, which uses the CONFIGURATION Menu as an example:



- A** Heading: this tells you where you are in the software.
- B** Information fields: these will either be numerical values or text fields providing information associated with the process shown in the heading. The inverted text shows the cursor position.
- C** Function key options: the functions of these change with the window displayed.

During the grinding process the screen could look as follows:




- A Remaining stock removal
- B Force applied on specimen holder
- C Load on main motor
- D Function(s) selectable during process

Please Note
The sample screens in this Instruction Manual show a number of possible texts. The actual screen displayed may differ from the samples in this manual.




Changing/Editing Values

Numeric Values


Depending on the type of value, there are two different ways of editing.

 Turn knob to select the value to be changed, e.g. *Stock removal*:







GRINDING SETUP		
Grinding mode:	Removal/Time	
Stock removal:	200 μm	
Force:	250 N	
Final time:	0:10 min.	
Final force:	100 N	
F1  Water ON	F2  One dress	F3  Dressing setup




 Push knob to edit value.

A scroll box appears around the value.




GRINDING SETUP		
Grinding mode:	Removal/Time	
Stock removal:	200 	
Force:	250 N	
Final time:	0:10 min.	
Final force:	100 N	
F1  Water ON	F2  One dress	F3  Dressing setup




 Turn knob clockwise to increase, or counter-clockwise to decrease the numeric value.



 Push knob to accept the new value. (Pressing **Esc**, aborts the changes, preserving the original value.)

Text Values


 Turn knob to select the text value to be changed, e.g. *Language*



CONFIGURATION	
Display contrast:	40
Units:	µm
Language:	English
Acoustic signal:	Yes
Water when grinding:	Yes
Type of grinding disc:	Stone

F1 / Default value



 Push knob to edit the value.



If there are only two options, then a toggle function is active, e.g. Yes / No.

If there are more than two options, a pop-up menu appears:

SELECT LANGUAGE	
English	▲
Deutsch	
Français	
Español	
ニホンコ*	▼



AbraPlan-20
Instruction Manual



Turn knob to select/toggle the correct option.




Push knob to accept the new value. (Pressing **Esc**, aborts the changes, preserving the original setting.)

Setting up the Software


Before operating the AbraPlan-20, check/edit the software configuration values:



Turn knob to select CONFIGURATION.



Push knob to activate the CONFIGURATION Menu.



Turn knob to select desired parameter.



Push knob to edit the parameter setting.



The settings possible and the default settings are shown below:

	Possibilities	Increment	Default
Display contrast	10 - 100	1	40
Units	µm / mils		µm
Language	English German French Spanish Japanese Chinese		English
Water when grinding	Yes / No		Yes
Keyboard sound	Yes / No		No
Type of grinding disc	Stone / diamond grinding disc		Stone




Repeat the process until all values are correct.



Press **Esc** to return to the MAIN MENU.

Inserting/Removing the Specimen Holder



Inserting the Specimen Holder


- Position the specimen holder under the quick coupling, and support it with your fingertips.
- Press and hold the flange of the column down with the heel of your hand while guiding the pressure tap of the specimen holder into the coupling.
- Turn the specimen holder until the three pins engage with the corresponding holes. Check that the specimen holder is in position by trying to turn it to the side.
- Release the flange.


Removing the Specimen Holder


- Place your hands as shown in the illustration above.
- Press upwards with your fingers to lift the specimen holder slightly. At the same time, press and hold down the flange with the heel of your hand.
- Using your fingers to support the specimen holder; lower it free of the coupling.
- Release the flange and completely remove the specimen holder.

Grinding Setup

Before starting the grinding process, the grinding time and force setting should be checked/ changed. This is done using the control panel software:

- ↓
If the MAIN MENU is not currently displayed, press **Esc** twice
-  Turn knob to select *Grinding*:




- ↓
-  Push knob to enter the GRINDING MENU.



Setting the Process Time



Turn knob to select *Grinding mode*.



Push knob to display square brackets [] and turn knob to change to *Time* setting.

Note that stock removal is not available when time is selected.





Push knob to accept the new setting.


Then adjust the process time and the force value to the correct settings.




Setting Removal

 Turn knob to select *Grinding mode*.



 Push knob to display square brackets [] and turn knob to change to *Removal* setting.



 Push knob to accept the new setting.






The Removal process has a 15 minutes time-out. If the prescribed amount of material has not been removed within this time, the process stops automatically. Then adjust the stock removal and the force value to the correct settings.

Setting Removal/Time



Turn knob to select *Grinding mode*.



GRINDING SETUP	
Grinding mode:	Removal/Time
Stock removal:	200 μm
Force:	250 N
Final time:	0:05 min.
Final force:	100 N
F1  F2  F3 	



Push knob to display square brackets [] and turn knob to change to *Removal/Time*.



Push knob to accept the new setting.



Then adjust the process time, the stock removal and the force values to the correct settings.

Note

Removal/Time is used when absolute planeness is required. First the required amount of material is removed, then the stone is dressed and the samples are ground again for a very short time. This ensures a maximum planeness after the required amount of material has been removed.

Cooling Water

The cooling water pump will automatically start when the preparation process is started. To stop the pump, press **F1**.

Please Note

All the value and function settings defined in the GRINDING SETUP menu are saved in the AbraPlan-20's memory. Therefore, if the power to the machine is interrupted, these values are remembered.

Starting the Preparation Process

- Insert the specimen holder.
- Lower the safety guard.
- Using the software display, enter the GRINDING menu and select *Removal*, *Time* or *Removal/Time* mode.

AbraPlan-20 has 3 different grinding modes:


Removal: Select *Removal* to remove a specific amount of material from the specimens.

Time: Select *Time* to always use the same grinding time.

Removal/Time: For maximum planeness a combination of removal and time can be used. Firstly a specified amount of material is removed. The stone is then dressed. A very short grinding step on the newly dressed, plane stone is then carried out.

- Check the correct settings for *Removal*, *Time* and *Force*.
- Start the grinding process by pressing the Start button.


The software display will appear as:



Please Note
Although the force setting can be changed during the grinding process, this is not recommended when using Stock removal.




**Stopping the
Preparation Process**

When the time has elapsed or the specified material removal has been obtained, the grinding stone will automatically stop rotating and the specimen holder will return to its initial position.







Dressing Functions and Changing the Grinding Stone

AbraPlan-20 is fitted with a diamond tool for automatic dressing of the grinding stone. It is important to dress the grinding stone at regular intervals to keep the stone plane and sharp. It is recommended that the automatic dressing function is enabled.

GRINDING SETUP	
Grinding mode:	Time
Process time:	0:30 min.
Force:	250 N
<p> F1  Water ON F2  One dress F3  Dressing setup </p>	


To exchange or insert a grinding stone:

- Press **F3: Dressing setup** from the GRINDING SETUP menu:


DRESSING	
Dresser step:	20 µm
Dresser speed [1-low...100-high]:	30
Automatic dressing during process:	Yes
Automatic dressing after process:	Yes
Dressing mode:	Removal
Dresser sensitivity:	60 %
Dressing during grinding:	No
Remaining height of stone:	9.89 mm
<p> F1  Finding TOP of stone F2  One dress F3  Up F4  Change Grinding Stone </p>	

- In the DRESSING menu, press **F4: Change Grinding Stone** to be guided through the changing sequence.

- The dresser is moved into the topmost position and the following screen is displayed:




- Press F1 and the next screen appears:




- Select the inserted grinding media by pressing F1 or F2.

Important

The dressing arm with the diamond tool should never be used on a diamond grinding disc since this would destroy the disc or the diamond tool. Therefore always select the correct media when inserting or changing the grinding media.



- Press F1 to continue and start the calibration process. The dresser checks the height of the grinding stone in two positions, one close to the centre, and the other at the periphery. Wherever the stone is highest, the following fine adjustment will start.



- Press F1 to start fine adjustment. The dresser will retract slightly to avoid any possible damage to the stone. The stone will begin to rotate.

FINE ADJUSTMENT OF DRESSER

Vertical dresser position [encoder units]: -5000
Horizontal dresser position [encoder units]: +350

Close the hood and press F1 to continue. The stone will start rotating.

Press Esc if you do not want to use the guide for the following procedure.



- Press F1 to continue with the fine adjustment.

FINE ADJUSTMENT OF DRESSER

Vertical dresser position [encoder units]: -5000
Horizontal dresser position [encoder units]: +350

Repeat pressing Enter/(F4) until the dresser has reached the top of the stone.

Each Enter/(F4) activation moves dresser down by 40 μm .

Each F3 activation moves the dresser up by 20 μm .


Turn the knob to make a small left-right movement between each Enter activation for better recognition of the first contact of the dresser tip with the stone.

Press F1 as soon as the dresser tip touches the stone, or press Esc to interrupt this process.



- Press Enter or F4 to move the dresser downwards in steps of 40 μm . Turn the knob to move the dresser sideways to ensure that the dresser has touched the grinding stone.

- Press F1 to finish the fine adjustment.



- Press F1 again to store the current dresser position as the top of the grinding stone.

Setting the Dressing Parameters

Several parameters can be set for the best possible utilization of the dresser.

DRESSING	
Dresser step:	20 μm
Dresser speed [1–low...100–high]:	30
Automatic dressing during process:	Yes
Automatic dressing after process:	Yes
Dressing mode:	Removal
Dresser sensitivity:	60 %
Dressing during grinding:	No
Remaining height of stone:	9.89 mm
F1 Finding TOP of stone F2 One dress F3 Up F4 Change Grinding Stone	

Dresser step:

The distance the dresser is moved down for every step.
Can be set between 20 μm and 100 μm in steps of 20 μm . Always try to use the smallest possible step size, this will keep the consumables cost down. However, enough material must be removed from the stone to ensure the best possible grinding result.

Dresser speed:

The speed of the dresser when moving across the stone.
Can be set between 1 and 100.
A speed of 30 is recommended for most types of grinding stones. But depending also on the dresser step size, adjustments can be made to both higher and lower values.

Dresser speed can be adjusted depending on the dresser step size:
Reduce speed for hard grinding stones or large dresser steps
Increase speed for soft grinding stones and small dresser steps.

Automatic dressing during process:	Should be set to <i>Yes</i> to enable automatic dressing while grinding, especially when removal rate mode is used.
Automatic dressing after process:	Can be selected to automatically dress the stone after every process. Should be set to <i>Yes</i> when time mode is used to ensure a sharp grinding stone when starting the process. Can be set to <i>No</i> when removal mode is used and <i>Automatic dressing during process</i> is set to <i>Yes</i> .
Dressing mode:	Dressing mode can be set to either <i>Removal</i> or <i>Time</i> , depending on the grinding mode. When the Grinding mode is set to <i>Time</i> , the Dressing mode is also set to <i>Time</i> . When the Grinding mode is set to <i>Removal</i> , both <i>Time</i> and <i>Removal</i> can be selected. Removal:- dressing will automatically start when the removal rate decreases to a certain level. Time:- the stone will be dressed in regular intervals.
Dresser sensitivity:	When the dressing mode is set to <i>Removal</i> , <i>Dresser Sensitivity</i> can be adjusted. The sensitivity of the dresser during the grinding process can be set to values between 20 – 100 % in steps of 20 %. High sensitivity means that the stone is dressed as soon as the removal rate decreases, low sensitivity allows for a higher decrease in removal before the stone is dressed. High sensitivity provides the shortest possible grinding times whereas low sensitivity gives a longer lifetime of the grinding stone.
Dressing interval:	When the dressing mode is set to <i>Time</i> , <i>Dressing interval</i> can be adjusted. The interval can be set to between 0:10 and 5:00 minutes in steps of 10 seconds.
Dressing during grinding:	This function allows dressing while the specimen holder is still on the grinding stone. It should be used when grinding very hard materials where frequent dressing is necessary to keep the grinding time at a minimum.

Remaining height of stone:	Shows how much grinding stone is left.
Function buttons	F1: Function to automatically find the top of the grinding stone F2: Carries out a single dressing of the grinding stone F3: Moves the dresser upwards F4: Starts sequence to change the grinding stone
Dressing the Diamond Grinding Disc	In order to dress a diamond grinding disc, mount 3 aluminium oxide dressing sticks in a sample holder and grind for a few seconds.

3. Maintenance

Daily Service

- Remove the filter paper in the static filter of the Recirculation Cooling Unit.
- Refill the Recirculation Cooling Unit.

IMPORTANT


Always use the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to top up with Struers Additive each time you refill with water.

Weekly Service

AbraPlan-20 and the Recirculation Cooling Unit should be cleaned regularly, in order to avoid damaging effects to the machine and the specimens from abrasive grains or metal particles.

Cleaning the Drain Pan

- Remove Stone guard as described earlier.
- Remove the drain pan cover by unscrewing the 2 finger screws (located on the right-hand side).



- Remove any material that may have accumulated on the bottom of the drain pan.
- Replace the drain pan cover and the stone guard.

*Checking the Recirculation
Cooling Unit*

The cooling unit should be checked for cooling water after 8 hours use or at least every week. The unit must be refilled if the flushing pump cannot reach the cooling water or if the cooling water is too contaminated.

Please refer to the manual supplied with the Struers Cooling Units for instructions.

Monthly Service

When cleaning the Recirculation Cooling Unit empty the tank of cooling water and clean the tank and the static filter mounted under the discharge branch.

Replacing the Cooling Water

- Replace the cooling water in the Recirculation Cooling Unit at least once a month.

Please refer to the manual supplied with the Struers Cooling Units for filling instructions.

Yearly Service

Inspection of Cover

- Visually inspect the cover and the glass for signs of wear or damage.


Important

Carry out inspection at more regular intervals if AbraPlan-20 is used for more than one 7 hour shift a day.

Struers recommends that the PETG glass in the cover is replaced after 5 years of routine operation.

The cover should be replaced immediately if it has been weakened by collision with projectile objects or if there are visible signs of deterioration as a result of using a cooling fluid other than those produced by Struers.

A label on the cover indicates when the cover glass is due to be replaced.



Reference Guide

Table of Contents	Page
1. Accessories and Consumables	46
Service Information	47
2. Struers Metalog Guide™	48
3. Trouble-Shooting	49
4. Technical Data	54

1. Accessories and Consumables

Please refer to the [AbraPol-20 Brochure](#) and the [Consumables Catalogue](#) for details of the range available.

Remember...

Struers offers a comprehensive range of consumables for grinding and polishing.

Service Information

Struers recommends that a regular service check be carried out on a yearly basis or after every 1500 hours of use. Information on total operation time and servicing of the machine is displayed on the screen at start-up:



A pop-up message will appear after 1,000 hours operation time to remind the user that a service check should be scheduled.

After the 1,500 hours operation time has been exceeded the pop-up message will change to alert the user that the recommended service interval has been exceeded.



- Contact a Struers Service Technician to service the machine.

2. Struers Metalog Guide™

In Struers Metalog Guide™ you will find a detailed description of grinding/polishing methods for automated mechanical specimen preparation.

Struers Metalog Guide™ offers preparation methods for the most common materials, based on a simple analysis of two key properties: hardness and ductility. Finding the right method is easy, including choice of consumables. Always consult Struers Metalog Guide™ on the Struers website for the correct preparation method for the actual specimens.

Metalog Guide™

A complete guide to materialographic specimen preparation.

www.struers.com/KNOWLEDGE/Metalog_Guide.

3. Trouble-Shooting

Troubleshooting the AbraPlan-20 is achieved using popup messages displayed in the LCD display on the console unit. These messages are divided into several categories and are listed below with decreasing severity:

Errors

Process cannot continue before an authorised technician has rectified the error. Turn off the unit at the main switch immediately. Do not attempt to operate the unit before a technician has rectified problem.

Warnings

Warnings must be rectified before process can continue.


Messages

Messages are intended to inform the operator of the machine's progress and advise about minor operational errors.

When a popup message is displayed, it will have the following format:

- a heading showing one of the categories listed above.
- an information field providing a message or instruction.

An example of a popup is shown below:



*AbraPlan-20
Instruction Manual*

The information in the following table is divided by category.

Error message	Explanation	Action required
Error		
Main supply voltage is too low! Please restart the machine. (#18)	No power to circuitry in AbraPlan-20	Contact a Struers Service Technician.
15V DC in PCB missing! Please restart the machine. (#19)	No power to circuitry in AbraPlan-20	Contact a Struers Service Technician.
24V DC in PCB missing! Please restart the machine. (#24)	No power to circuitry in AbraPlan-20	Contact a Struers Service Technician.
No RS 485 communication. Please restart the machine. Call service if error persists. (#45)	Communication problem between elements inside AbraPlan-20.	Restart the machine. Contact Technical Support if error persists.
Machine and console software are not compatible, please upgrade. (#43)	Incompatible software in AbraPlan-20.	Contact a Struers Service Technician.
There is problem with contactors K7 or K8. Please call technical service. (#46)	Contacto(r)s are wrong.	Contact a Struers Service Technician.
LIN bus: X SMU is off line. Please restart the machine. (#1)	Dresser X motor is not responding.	Restart the machine. Contact a Struers Service Technician if error persists.
LIN bus: Y SMU is off line. Please restart the machine. (#2)	Dresser Y motor is not responding.	Restart the machine. Contact a Struers Service Technician if error persists.
LIN bus: X and Y SMU are off line. Please restart the machine. (#3)	Both X and Y motors are not responding.	Restart the machine. Contact a Struers Service Technician if error persists.

*AbraPlan-20
Instruction Manual*

Error message	Explanation	Action required
Warning		
The air pressure is too low (#27)	There may be a leak in the hose or the compressor maybe defective.	Check the compressed air system for possible cause.
Disc motor is overloaded (#16)	Load pressure is too high.	Reduce force value.
Emergency Stop is active	The Emergency Stop is pushed in. This is displayed until action is taken.	Release the Emergency Stop.
The removal rate is too low. The time limit was exceeded! (#22)	The stone is not dressed regularly.	Allow for automatic dressing during the process.
	The stone is not suited for the material to be ground.	Replace the stone with one suited for that application.
The sample mover motor is unable to move upwards after the process! (#23)	There may be a problem with the compressed air system.	Check the compressed air system.
	There may be an internal electrical problem.	Contact a Struers Service Technician.
Sample motor is overloaded!(#17)	Loading on the motor is too high.	Reduce the force value.

*AbraPlan-20
Instruction Manual*

Message	Explanation	Action required
Messages		
The grinding stone must be replaced ! (#26)	Message received during the dressing process.	Stone is worn and must be replaced.
The process is already stopping. (#15)	Message if the Stop button is pushed when the process is already finished.	
Process in progress (#12)	Message if a button is pressed while the grinding process is running.	
The process is finished. (#14)	Message at the end of the process.	
The process was stopped through the emergency stop. Press F1 to raise the sample mover motor.	The emergency stop has been pressed and the sample motor must be raised manually.	Press F1 to raise the sample mover motor.
Safety guard not closed! (#11)	The grinding process can not start because the safety guard is not closed.	Lower the safety guard and start the process.

AbraPlan-20
Instruction Manual

	Explanation	Action
Physical Observations/Problems		
No material removed	The grinding stone/diamond grinding disc is covered with material.	Dress the stone/disc. Replenish the cooling water volume.
	Insufficient grinding force.	Regulate the grinding force.
The specimen heats up	No cooling water.	Replenish the volume of cooling water in the recirculation unit.
		Cooling water pump is blocked.
Un-plane specimens	Stone not dressed.	Dress the stone. Select Removal/Time as the correct grinding mode
	Too few specimens in the specimen holder. Badly centred large specimen or specimen with too small a distribution in one direction.	Put blank specimen/s in the specimen holder.
Squeaking noise	V-belt slides.	Please contact a Struers Service Technician.
Hissing noise when the machine is in operation and/or is switched off	Leak in the air system.	Tighten the fittings and/or replace the defective air tubing.
Violent vibrations when the machine is running idle	The grinding stone is defective and out-of-balance.	Exchange the stone. Before doing so, try to turn the stone in relation to the turntable.
Machine very noisy when running idle	Axial bearing defective.	Please contact a Struers Service Technician.
	Spindle or motor bearings defective.	Please contact a Struers Service Technician.
Continuous, irregular wear on a grinding/polishing surface.	Coupling on either the specimen holder/mover plate or the specimen mover head of the polishing machine is worn.	Please contact a Struers Service Technician to replace the coupling.

4. Technical Data

Subject		Specifications	
		Metric/International	US
Grinding Stone/ Disc	Rotational speed	1450 rpm	1450 rpm
	Size	356 mm	14.0"
	Power consumption	4 kW	5.4 Hp
Specimens	Speed	150 rpm	150 rpm
	Direction	CCW	
	Force	50-700 N	10-150 lbf
	Motor power consumption	0.37 kW	0.5 Hp
Software and Electronics	LC Display	320x240 pixels	
	Controls	Touch pads/Push-turn knob	
	Memory	EPROM/RAM/NV-RAM	
Compressed Air	Compressed air supply	6-10 bar	
Dimensions and Weight	Width	840 mm	33.1"
	Depth	980 mm	38.6"
	Height	1560 mm	61.4"
	Weight	400 kg	880 lbs

AbraPlan-20
Instruction Manual

Subject		Specifications			
Electrical Data					
Supply Voltage	<i>Power consumption</i>	4.4 kW			
	<i>No. of phases</i>	3 (3L+PE)			
	<i>Output, main motor</i>	4.0 kW			
	<i>Voltage/frequency:</i>	<i>Max. Load:</i>			
	3 x 200 V / 50Hz	16.9 A			
	3 x 200 - 210 V / 60Hz CSA	15.7 A			
	3 x 220 - 230 V / 50 Hz	16.9 A			
	3 x 220 - 240 V / 60 Hz	15.7 A			
3 x 380 - 415 V / 50 Hz	8.9 A				
3 x 380 - 415 V / 60 Hz	10.3 A				
3 x 460 - 480 V / 60 Hz CSA	8.5 A				
Mains Cable Recommendation	<i>Voltage/ frequency</i>	Min. Fuse	Minimum cable size @ Min. fuse	Max. Fuse	Minimum cable size @ Max. fuse
	3 x 200 V / 50Hz	25	3x2,5mm ² + PE	40	3x2,5mm ² + PE
	3 x 200 - 210 V / 60Hz CSA	25	3xAWG12 + PE	40	3xAWG12 + PE
	3 x 220 - 230 V / 50 Hz	25	3x2,5mm ² + PE	40	3x2,5mm ² + PE
	3 x 220 - 240 V / 60 Hz	25	3xAWG12 + PE	40	3xAWG12 + PE
	3 x 380 - 415 V / 50 Hz	20	3x2,5mm ² + PE	40	3x2,5mm ² + PE
	3 x 380 - 415 V / 60 Hz	20	3xAWG12 + PE	40	3xAWG12 + PE
	3 x 460 - 480 V / 60 Hz CSA	20	3xAWG12 + PE	40	3xAWG12 + PE
Important: Local standards may overrule the recommendations for the main supply cable. If necessary, please contact a qualified electrician to verify which option is suitable for the local installation setup.					
Residual Current Circuit Breaker	type A, 30 mA (or better) is recommended.				
Environment	Safety standards	Please refer to the Declaration of Conformity			
	Noise level(idle)	77dbA			
	Surrounding temperature	5-40°C			
	Humidity	Max. 95%RH			

Quick Reference Guide

- | | |
|---|--|
| Inserting the Specimen Holder | <ul style="list-style-type: none">■ Position the specimen holder under the quick coupling.■ Press and hold the flange of the column down with the heel of your hand while guiding the pressure tap of the specimen holder into the coupling.■ Turn the specimen holder until the three pins engage with the corresponding holes.■ Release the flange. |
| Removing the Specimen Holder | <ul style="list-style-type: none">■ Press upwards with your fingers to lift the specimen holder slightly. At the same time, press and hold down the flange with the heel of your hand.■ Using your fingers to support the specimen holder; lower it free of the coupling.■ Release the flange and completely remove the specimen holder. |
| Starting the Preparation Process | <ul style="list-style-type: none">■ Insert the specimen holder.■ Lower the safety guard.■ Enter the GRINDING menu and set/check the correct removal and/or time and force.■ Start the grinding process. |
| Stopping the Preparation Process | <ul style="list-style-type: none">■ When the time has elapsed, the grinding stone will automatically stop rotating and the specimen holder will return to its initial position. |
| Dressing the Grinding Stone | <ul style="list-style-type: none">■ Press F2 for a single dressing of the grinding stone. |
| Dressing the Diamond Grinding Disc | <ul style="list-style-type: none">■ Mount 3 aluminium oxide dressing sticks in a sample holder and grind it for a few seconds. |

English

Declaration of Conformity

 Struers

**Manufacturer,
responsible for
Technical File**

Struers ApS
Pederstrupvej 84
DK-2750 Ballerup, Denmark
Telephone +45 44 600 800

Herewith declares that

<i>Product Name:</i>	AbraPlan-20
<i>Type No.:</i>	589
<i>Machine Type:</i>	Grinding machine

is in conformity with the provisions of the following directives:

Safety of Machinery 2006/42/ EC according to the following standard(s):
EN ISO 12100:2010, EN ISO 13849-1:2008/AC:2009, EN ISO 13849-2:2012, EN ISO 13850:2008,
EN 60204-1:2006/AC:2010, EN 574:1996+A1:2008; EN 953:1997+A1:2009,
EN 349:1993+A2:2008, EN 1037:1995+A1:2008.

EMC-Directive 2014/30/EU according to the following standard(s):
EN 61000-6-1:2007, EN61000-6-3:2007/A1:2011.

RoHS 2011/65/EU according to the following standard(s):
EN 50581:2012.

**Supplementary
Information** The equipment complies with the American standards:
UL508, NFPA70:2014; NFPA79:2012, FCC 47 CFR part 15.

The above has been declared according to the global method, module A

Date: 23.02.2016


Christian Skjold Heyde,
Vice President, R & D and Production, Struers ApS

Dansk

Overensstemmelseserklæring

 Struers

**Fabrikant,
ansvarlig for Teknisk
Dossier**

Struers ApS
Pederstrupvej 84
DK-2750 Ballerup, Danmark
Telefon 44 600 800

erklærer herved, at

<i>Produktnavn:</i>	AbraPlan-20
<i>Type nr.:</i>	589
<i>Maskintype:</i>	Slibe maskine

er i overensstemmelse med følgende EU-direktiver:

Maskindirektivet 2006/42/EF efter følgende norm(er):
EN ISO 12100:2010, EN ISO 13849-1:2008/AC:2009, EN ISO 13849-2:2012, EN ISO 13850:2008,
EN 60204-1:2006/AC:2010, EN 574:1996+A1:2008; EN 953:1997+A1:2009,
EN 349:1993+A2:2008, EN 1037:1995+A1:2008.

EMC-direktivet 2014/30/EU efter følgende norm(er):
EN 61000-6-1:2007, EN61000-6-3:2007/A1:2011.

RoHS 2011/65/EU efter følgende norm(er):
EN 50581:2012.

**Supplerende
oplysninger** Endvidere overholdes de amerikanske normer:
UL508, NFPA70:2014; NFPA79:2012, FCC 47 CFR part 15.

Ovenstående overensstemmelse(r) er erklæret iflg. den globale metode, modul A

Dato: 23.02.2016


Christian Skjold Heyde,
Vice President, Udvikling og Produktion, Struers ApS



Pederstrupvej 84
DK-2750 Ballerup
Denmark

AbraPlan-20



Spare Parts and Diagrams

Manual No.: 15897001

Date of Release 2H01 .201H



AbraPlan-20
Spare Parts and Diagrams

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manuals may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manuals may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 201H.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Spare Parts and Diagrams

Table of contents

Drawing

AbraPlan-20

Drawings

Sample motor, assembly.....	15890045C
Motor for grindstone	15890021C
AbraPlan-20, complete.....	15890001P
Casing with electrical, assembly	15890006E
Casing with motor, assembly	15890007G
Plate with PCB & pneumatic distribution	15890032D
Stock removal unit, complete	15890083G
Control box, assembly.....	15890082H
Safety guard assembly.....	15890081C
Arm, assembly	15890073F
Flushing gun, assembled	15490009C
Quick-release coupling, complete	15490007C
Box for tub, assembly.....	15890011G
Air connection, assembled	15090032E
Bearing housing, assembled.....	15090040A
Step motor, assembled	15480018J
Cover for grindstone, assembly	15890008C
Casing, assembly	15890010K
Contact box, assembled	15890013J
Main mechanism, assembly.....	15890020G
Dresser, assembly	15890070M
Moving part of dresser, assembly	15890071C

Diagrams

Transformer Connections (2 pages)	15093452A
Air diagram.....	15892000C
Block diagram	15893050C
Circuit diagram main voltage.....	15893100E
CPU Board A1 (5 pages)	15893120D
Wiring diagram (6 pages).....	15893450B
Wiring of Brake Transformer	15893451A
Variants parts	15897600G

Some of the drawings may contain position numbers
not used in connection with this manual.

AbraPlan-20
Spare Parts and Diagrams

The following is a list of the spare parts that may need replacement during the lifetime of the equipment.

To check the availability of other replacement parts, please contact your local Struers Service Technician. It may help identify the part by referral to its position number on the assembly drawings included in this manual.

Spare Part list for AbraPlan-20

Drawing
 15890045

Pos.		Cat no.
	Sample motor, assembly	
20	Coupling	15490410
10	GEAR MOTOR 3X200V 50HZ. painted	15499016
10	GEAR MOTOR 3X220-230V 50HZ. painted	15499017
10	GEARMOTOR. 3X380-415V 50HZ. painted	15499018
10	GEAR MOTOR. 3X200-208V 60HZ.CSA.painted	15499019
10	GEAR MOTOR 3X220-240V 60HZ.CSA.painted	15499020
10	GEAR MOTOR. 3X380-415V 60HZ. painted	15499021
10	GEAR MOTOR 3X460-480V 60HZ.CSA.painted	15499022
60	Hose Adaptor 45° SVAD-P167GT	2NM10437
40	Hose Adaptor. SVNV-M257.50	2NM10472
40	Hose Adaptor. SVNV-M257.50	2NM10472
40	Hose Adaptor. SVNV-M257.50	2NM10472
40	Hose Adaptor. BVND-N027GT.100	2NM11027
40	Hose Adaptor. BVND-N027GT.100	2NM11027
40	Hose Adaptor. BVND-N027GT.100	2NM11027
40	Hose Adaptor. BVND-N027GT.100	2NM11027
50	Flexible Hose PMA PIST-17S.30, 0.4 m	2NU31200

15890021

	Motor for grindstone	
60	V-Belt SPZ/3 ø106	2JE10106
10	Motor 3x220-240VD / 50Hz 4kW	2ME06205
10	Motor 3x220-240VD / 50Hz 4kW	2ME06205
10	Motor 3x220-240VD / 50Hz 4kW	2ME06205
10	Motor 3x380VD/60Hz 4,0kW	2ME06386
10	Motor 3x400VD/50Hz-480V/60 4kW	2ME06405
10	Motor 3x208VD/60Hz 4,0kW CSA	2ME56206
10	Motor 3x480VD/60Hz 4,0kW CSA	2ME56486
30	Hose Adaptor. SVNV-M257.50	2NM10472
30	Hose Adaptor. SVNV-M257.50	2NM10472

Spare Part list for AbraPlan-20

Drawing	Pos.		Cat no.	
15890001		AbraPlan-20, complete		
	140	Down arm, 2pcs	15890930	
	150	Top right arm, welded	15890900	
	160	Top left arm, welded	15890905	
	170	Flange bearing GFM-2528-21, 2pcs	2BG00089	
	180	Flange bearing GFM-2023-07, 2pcs	2BG00088	
	210	Safety guard AbraPlan-20, assy	15890081	
	300	Arm, assembly	15890073	
	340	Flushing gun, complete	15490009	
	360	Rubber	15890508	
	470	AbraPlan-10 Quick-release coupling, complete	15490007	
	15890006		Casing with electrical, assembly	
		110	35A 800V KBPC3508 BRIGDE RECT.	2VB30750
140		Terminal block with spring	2XL00301	
	145	Terminal block double with spring, grey	2XL00331	
15890007		Casing with motor, assembly		
	100	Rubber disc Ø12/Ø26.4, 4pcs	11440069	
	110	Rubber bushing	15090690	
15890032		Plate with PCB & pneumatic distribution		
	30	PCB AbraPlan-20 A2, testet	15893002	
	60	Pressure Regulator, 5-8.5 bar 1/4in	2YR00001	
	70	Gasket, PVC O-1/8	2IF00011	
	80	Throttle-sound absorber. RSS-111-M35-1/8	2YL00035	
	100	Gasket, PVC 1/4"	2IF00012	
	110	Banjo til quick-coupling ø5-1/8	2NF10034	
	120	Banjo screw 1631-03-1/8"	2NF20080	
	130	Quick release angle swivel connector ø5-1/8"	2NF10082	
	150	Magnet vent. 3/2 24V DC 1/8	2YM10030	
	210	3/2 solenoid valve 24VDC	2YM10124	
	220	Sound absorber, SINTER 2931-M5	2YL00015	
	240	Gasket, PVC M5	2IF00010	
	270	Pressure nipple RTU PK3/3	2NF40242	
	290	Neopren nipple ø36/ø47/ø54-2.5	2GK90457	
	370	Air tube ø5/ø3.2 Superflex	2NU12445	

Spare Part list for AbraPlan-20

Drawing	Pos.		Cat no.
15890083		Stock removal unit, complete	
	20	Self-lubricating bearing ø20/28x32	2BG32032
	40	Charnier for potentiometer	15490830
	120	Pin for rate measur. unit	15890800
15890082		Control box, assembly	
	10	Display, 320X240 w. white LED	2HD32024
	40	Main PCB f.AbraPlan-20, tested	15893000
	130	Pushbutton Head RVAT DG stainl.	2SA00400
15890081		Safety guard assembly	
	10	Hood for AbraPlan-20	15890441
	100	Brace of safety guard	15890410
	140	Straight Actuator AZ 17/170-B1	2SS10017
15890073		Arm, assembly	
40	Diamond dresser/CDP8181-18/22	12660212	
15490009		Flushing gun, assembled	
	10	Flushing head	15490535
	20	Tube, flushing guns	15490537
	25	Tube, internal, flushing guns	15490538
	30	Silicone hose ø8/ø12	2NU19208
	70	Push button	15490545
	80	Slide bearing.M.KR. 12x15x8/18x1.5	2BG00120
120	Magnet 10x10x3 VACODYM 351 WZ	2LM00034	
15490007		Quick-release coupling, complete	
	10	Quick-release coupling	15090009

Spare Part list for AbraPlan-20

Drawing
15890011

Pos.		Cat no.
Box for tub, assembly		
20	Tub, assembly	15890057
30	Sealing disc	15490512
60 + 65	Cover for grindstone, assembly	15890008
65	Cover for grindstone, welded	15890511
60	Top of Grindstone cover	15890514
70	Distance Bushing	15490511
110	Disc for stone. Replaced by R5490006	15490006
170	Elbow 87 for hose ø51(2") pipe socket ø50	2NG20587
175	Drain tube, straight ø50x250	2NG25026

15090032

Air connection, assembled		
40	Air filter, air regulation EAW3000-F02D-6	2YF00005
90	Gasket, PVC 1/4"	2IF00012
100	Nipple 2531-1/4-1/8	2NF40041
110	Gasket, PVC O-1/8	2IF00011
120	Banjo screw 1631-03-1/8"	2NF20080
130	Banjo to quick-coupling ø5-1/8	2NF10034
140	PVC-Hose, clear 13/32"-Ø10	2NU19313
145	PVC-hose 10 mm	2NP00010
150	Air tube ø5/ø3.2 Superflex	2NU12445
160	Quick coupling	2NF10024
170	End piece	2NF40071
190	Angle Quick coupling, Ø8-1/4"	2NF10087
200	Distance nipple.2525-1/4-1/4-27	2NF40181

15090040

Bearing housing, assembled		
3	Ball bearing 6208-2RS1 ø40/80	2BK00120
4	Angle contact bearing ø50/ø90x20	2BK30050
5	Nilos-ring 7210AVH	2BK97210
6	Disc spring for ball bearing 79.5x55.5x0.8.	2GF51026
7	Locking ring J80 DIN 472	2ZL20800

Spare Part list for AbraPlan-20

Drawing	Pos.		Cat no.
15480018		Stepmotor, assembled	
	10	Stepper Motor assembl. with plug	15483532
	20	Bushing for magnet, 2LS00050	15480624
	30	Magnet ø6x2.5 NdFeB	2LS00050
	60	Distance piece F-F, M3x25mm	2GZ10325
	70	PCB for magnet SMU, tested	15483005
15890008		Cover for grindstone, assembly	
	20	INA-Sealing ring G 10X17X3	2II01017
	30	O-RING 12.42-1.78 72 NBR 872	2IO17817
	60	Screw with ball and spring. GN615-M10-KN	2TX91019
	90	Nozzle for dresser	15890522
15890010		Casing, assembly	
	200	Key Lock Switch AZM 170-02ZRKA 24V	2SS00007
	305	Neopren bushing ø53/ø64/ø75-2.5	2GK90459
	310	Hose nipple 2601-12-1/4	2NF40087
	320	Gasket, PVC 1/4"	2IF00012
	330	Ball valve MINIBALL ¼ in-¼ in internal	2YH03622
	340	Armed PVC HOSE 1/2" -ø12.5 for water.	2NU29316
	360	GEKA hose connection 1-2	2NF60000
15890013		Contacto box, assembled	
	30	Contacto CA4-5-0, 24V-50/60HZ	2KM04501
	35	Contacto CI4-5-01, 24VDC	2KM04502
	80	Contacto Danf. CI-25A/24VAC	2KM10641
	100	Auxiliary switch block for K1 CB-NO 037H0111	2KH00111
	110	Contacto CI 12 37H0032/13	2KM10232
	120	TRAFO 200-460V/24V+24V/200VA	2MT72034
	125	4.00A T FUSE GLASS 6,3x32 250V	2FU14200

Spare Part list for AbraPlan-20

Drawing

15890020

Pos.		Cat no.
	Main mechanism, assembly	
90	REED-KONTAKT D-A73L	2KR30177
140	Ball bushing KH4060	2BF20040
150	INA-Sealing ring G 40x52x5	2II04052
230	Distance ring-ball bearing	14590017
240	Sensor read disc	15490568
250	V-belt pulley SPZ/3 ø125	2JE10125
270	Adapter 2012/ø35	2JE92035
280	Pressure disc-V-belt	14590018
330	Proximity sensor A01G142	2HQ00023
340	Blocking valve R 1/4"	2YH60004
350	Nipple 2531-1/4-1/8	2NF40041
360	Banjo til quick-coupling ø5-1/8	2NF10034
380	Gasket, PVC O-1/8	2IF00011
390	Gasket, PVC 1/4"	2IF00012
400	Quick-coupling, straight ø5-M5	2NF10011
420	Terminal block with spring	2XL00301
440	V-belt A XPZ/3V 1060mm	2JD01060
450	Air tube ø5/ø3.2 Superflex	2NU12445

15890070

	Dresser, assembly	
20	Stepmotor, assembled	15480018
30	Coupling ROTEX GS12-22ø6,35ø12	2JH00003
60	Spherical ball Bearing 2201	2BK20012
70	Locking ring J32 DIN 472	2ZL20320
110	Locking ring A12 DIN 471	2ZL10120
120	Pressure spring ø25.0 x ø2.0 Lo=195 22830	2GF10250


15890071

	Moving part of dresser, assembly	
30	Radial bearing SMS 777 20-26-15	2BG30088
60	Locking ring A35 DIN 471	2ZL10350
70	Locking ring A26 DIN 471	2ZL10260
80	Cylinder pin, stainless 8m6x30 DIN 7	2ZS01530
100	Cylinder pin, steel 6m6x25 DIN 7	2ZS02455
110	Ball bearing,.61908-2RS1 ø40/62	2BK00118
120	Locking ring J62 DIN 472	2ZL20620
130	Wave spring Ø51X61X0.5 (5 pcs.)	2GF60038
140	Tooth wheel Ulmer T43519	15890188



Spare Part list for AbraPlan-20

Drawing
15890072

Pos.		Cat no.
	Motor with gear	
10	Stepmotor, assembled	15480018
30	Tooth wheel Ulmer T 16868	15890189
	Wireset+Cont.box f.AbraPlan-20	15893590
	Main switch KG32 K300E	2SE20317
	Terminal block with spring	2XL00301
	Terminal block double with spring, grey	2XL00331
	Contact block 1 NC 1/2 typeMTO	2SB10071
	Contact block 1 NO 3/4 typeMTI	2SB10072



Variant
M1: See AbraPlan-20, complete





20 30 Antisize grease
(Parting Lubricant 785FG)


20 Varmes op til 150°C og krympes på



30 Locked by Loctite 2701 (Låses med Loctite 2701)

C	2013.07.31	PMA fittings changed, Pos40,50,60,70	SPE	2013.07.31	
A	28.11.2006		YKJ		
Revision	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
		Material:	Scale: 1:3	Format: A3	Tolerance: DS/ISO 2768- mK Surface treat.: None
		ID:	Description:		Rev:
<small>Pedersstrøget 84 DK-2750 Ballerup/Copenhagen Denmark Phone:+45 44 600 800 Fax: +45 44 600 804</small>		15890045 Sample motor, assembly			C






ISOMETRIC VIEW
SCALE 1:4




C	17.12.2008	Pos.80 15090500->15890540	SPE	17.12.2008	
A	22.12.2006		JFR		
Revision	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
		Material:	Scale: 1:5	Format: A3	Tolerance: DS/ISO 2768- mK Surface treat.: None
ID:	Description:				Rev:
	15890021 Motor for grindstone , assembly				C


Pedersbølvej 84
DK-2750 Ballerup/Copenhagen
Denmark
Phone: +45 44 600 800
Fax: +45 44 600 804



Seen without Top!




SECTION N-N
SCALE 1 : 1



For adjustment of the Small front plate:
Standard view: 1 washer is necessary.
Option 1: The washer can be removed.
Option 2: One more washer can be added (+7mm).


Sheet 1/2

Project	15890001	AbraPlan-20	Scale	1:5	Sheet	1/2
Revision	001	AbraPlan-20	Scale	1:5	Sheet	1/2
Author	Scale	1:5	Sheet	1/2
Checked	Scale	1:5	Sheet	1/2
Drawn	Scale	1:5	Sheet	1/2
Appr.	Scale	1:5	Sheet	1/2
Project	15890001 AbraPlan-20, complete					





POS. NO.
10
20
30
40
50
60
70
80
90
100
110
120
130
140
150
160
200
210
220
230
250
260
270
280
290
300
310
320
330
350
360
370
380
390
400
410
420
430
440
450
460
470
480
490
500
510
520
530
540
550
560
570
580
590
600
610
620
630
640
650
660
670
680
690
700
710
720
730
740
750
760
770
780
790
800
810
820
830
840
850
860
870
880
890
900
910
920
930
940
950
960
970
980
990
1000


E	2013-06-12	DE: Pos. 190 -> pos. 250.	JMO	
A	12-01-06		JFR	12-01-06
Revision	Crea. date	Revision description	Draw. Init	Appr. date
				Appr. Init
			Scale	Forma
			1:5	A1
			Material	Tolerance: DIN/ISO 2768- msk
				Surface treat: None
ID	Description:			Rev.
	15890006 Casing with electrical, assembly			E



SECTION A-A





SECTION E-E



SECTION C-C
SCALE 1 : 2


SECTION D-D



250 All screws pos. 250 must be tightened by torque wrench to 24Nm



G	2012.05.16	2TJ10825 removed	SPE	2012.05.16	JTV
A	11.01.2006		JFR	11.01.2006	FPG
Revision	Crea. date	Revision description	Draw. Init	Appr. date	Appr. Init
			Scale	Forma	Tolerance: ISO/ISO 2768- msk
			1:5	A1	Surface treat: None
ID	Description:				Rev.
	15890007 Casing with motor, assembly				G

POS. NO.	10
	20
	30
	40
	70
	80
	90
	100
	110
	120
	130
	140
	150
	160
	170
	190
	200
	210
	230
	240
	250
	260
	270
	280
	290
	300
	310
	320
	340
	350
	360



- 60 Indstilles paa 3 bar
- 80 AAbnes 1/4 omgang

D	2009.04.16	L=500 til toppen af cyl. -> L=500 til bunden af cyl.	SPE	2009.04.16	
A	22.8.2006		JF	12-02-2007	FPG
Revision	Crea. date www.mm-td	Revision description	Draw. Init	Appr. date www.mm-td	Appr. Init
		Material:	Scale: 1:5	Format: A2	Tolerance: DS/ISO 2768- mK Surface treat.: None
ID:	Description: 15890032 Plate with PCB and pneumatic distribution, assembly				Rev: D



DETAIL A
SCALE 2 : 1


Stock removal unit controlled by making three resistance measurements:

Benchmarks	Terms	Measurement
Between pin 1 and pin 3	-	1 kohm
Between pin 2 and pin 3	pin completely out	1 kohm ± 100 ohm
Between pin 2 and pin 3	pin pushed totally in	170 ohm ± 50 ohm

Aftagningsmåleren kontrolleres ved at lave tre modstandsmålinger:


Målepunkter	Betingelser	Måling
Mellem ben 1 og ben 3	-	1 kohm
Mellem ben 2 og ben 3	Stift helt ude	1 kohm ± 100 ohm
Mellem ben 2 og ben 3	Stift helt inde	170 ohm ± 50 ohm

NB! Piston DO NOT lubricate
NB! Stempel må IKKE smøres



G	2012.10.08	Translated to english	SPE	2012.10.08	JTV
A	03-03-08		BMJ		
Revision	Crea. date yyyy-mm-dd	Revision description	Draw. Init	Appr. date yyyy-mm-dd	Appr. Init
		Material:	Scale: 1:1	Format: A3	Tolerance: DS/ISO 2768- mK Weight : g
ID:		Description: 15890083 Stock Removal unit, assembled			Rev: G

Pederstrupvej 84
DK-2750 Ballerup/Copenhagen
Denmark
Phone: +45 44 600 800
Fax: +45 44 600 804




40 Clamping torque: 2Nm.
Use Loctite 243

30 Round edges against plastic.

C	2010-12-20	Pos 30 exchanged and drawing updated	JTV	2010-12-20	JTV
A	2007-07-22		JFR		
Revision	Crea. date yyyy-mm-dd	Revision description	Draw. Init	Appr. date yyyy-mm-dd	Appr. Init
		Material:	Scale: 1:7	Format: A3	Tolerance: DS/ISO 2768- mK
		ID:	Description:	Surface treat.:	None
15890081 Safety guard AbraPlan-20, assembly					C


Struers
 Pederstrupvej 84
 DK-2750 Ballerup/Copenhagen
 Denmark
 Phone: +45 44 600 800
 Fax: +45 44 600 804



30 20 Locked by Omnifit 200M (Låses med Omnifit 200M)


ISOMETRIC VIEW

F	2009.10.13	P20 15890112->15890783, P40 15890161->12660212	SPE	2009.10.13	
A	2006.08.08		FPG	2006.08.08	FPG
Revision	Crea. date yyyy-mm-dd	Revision description	Draw. Init	Appr. date yyyy-mm-dd	Appr. Init
		Material:	Scale: 1:2	Format: A3	Tolerance: DS/ISO 2768- mK Surface treat.: None
<small>Pederstrupvej 84 DK-2750 Ballerup/Copenhagen Denmark Phone: +45 44 600 800 Fax: +45 44 600 804</small>		ID: Description: 15890073 Arm, assembly	Rev: F		





C: Pos. 10+120 changed to 15490065.
 2009.03.20 SPE
 B. 2 stk pos. 20 -> pos. 25
 010601: BRY/KLB

Matr.:	Overfi.beh.: Nej	Målforshold:	Ikke ang. tol. efter DS/ISO 2768- mK	
Struers	Projektionsmetode 	1:1	Date	Sign.
			Tegn:	141200 BRY
			Kontr.:	141200 KLB
Spulepistol, monteret		Erst.:		
		549-0009 C		



B. Pos. 20 (demperskive) fjernet. 020101 BRY/KLB		C: Pos.10: 15090007 -> 15090009. 051101 BMJ/JTV		
Matr.:	Overf.beh.:	Måforhold:	Ikke ang. tol. efter DS/ISO 2768-	
Struers	Projektionsmetode 	1:1	Date	Sign.
			Tegn:	160800 BRY
Opspændingsenhed, mont. AbraPlan-10		Erst.:	549-0007C	
		kontr. 160800 KLB		








DETAIL Q
SCALE 1 : 1

Degreased with alcohol.
Glued with silicone 515 all the way around.


Affedtes med sprit.
Limes med silicone 515 hele vejen rundt.

Afkortes til 160mm




G	2012-01-24	B.7: Glue regulatory	JJO		
A	20.12.2005		JFR	20.12.2005	FPG
Revision	Crea. date	Revision description	Draw. Init	Appr. date	Appr. Init
	08.01.2005			08.01.2005	
		Material:	Scale: 1:4	Format: A2	Tolerance: DS/ISO 2768: mK
					Surface treat.: None
ID:		Description:			Rev:
		15890011 Box for tub, assembly			G

POS. NO.	10
	20
	30
	40
	50
	60
	65
	70
	80
	90
	100
	110
	120
	130
	170
	175
	180
	190
	191



Pos. 40: Adjust to 6 bar
 Pos. 40: Indstilles til 6 bar
 Pos. 90: Omnifit 50 H seal

E	08.09.2008	Pos.210 added, pos. 110+150+160+170 moved.	SPE	08.09.2008	
A	29-02-08		SPE	02-04-07	JTV
Rev	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
F	 Pederstrupvej 84 DK-2750 Ballerup Copenhagen Denmark Phone :+45 44800 800 Fax : +45 44600 804	Material:	Scale: 1:2	Format: A4	Tolerance: DS/ISO 2768 - mK Surface treat.: None
		ID:	Description: 15090032 Air connection, assembled	Rev: E	



Smøres m. Tectyl

30

40

50

10

70






60

20


2:1

60 Alvanina Grease EP2


A: Vinkelkontakleje vendt korrekt; 250400/KLB				
Matr.:	Overfl.beh.:	Målforhold:	Ikke ang. tol. efter DS/ISO 2768-	
		1:1	Dato	Sign.
			Tegn:	011199 BRY
			Kontr.:	011199 KLB
Lejeopbygning, monteret		Erst.:		
		15090040A		



SECTION B-B








Magnet
Loctite 409Gel




POS. NO.	AMOUNT	DRAW. NO.	NOTE
10	1	15483532 Stepmotor, assembled with plug	
20	1	15480624 Bushing for magnet	
40	2	2TI10303 Msp skrue M3x3 A2	
50	4	2TI10310 Msp skrue M3x10 DIN916 70A2	
60	4	2GZ10325 Afstandsstag 6-KT M3x25	
70	1	15483005 SMU PCB+magnet, testet	
80	4	2TR50308 MC skrue M3x8 A2	

J	15.08.2008	Text for assembling of pos. 50 and 60 added.	SPE	15.08.2008	
A	13-03-08		CJE	27-11-07	AKR
Revision	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init
		Material:	Scale: 1:1	Format: A3	Tolerance: DS/ISO 2768- mK Surface treat.: None
ID:	Description: 15480018 Print og stepmotor, monteret				Rev: J



SECTION A-A
SCALE 1 : 1




SECTION B-B
SCALE 1 : 1

70 Must be sealed to pos. 50 and pos. 90 by Omnifit Seal 50 H


60 80 Adjusted according doc. 15897507 and sealed with Omnifit 50H

C	13.11.2007	Changed pos.10 (15890510→15890514)	JFR	13.11.2007	ppl
B	5.7.2007	Completely changed	JF	5.7.2007	MD
A	20.12.2005		JF	20.12.2005	FBG
Revision	Crea. date	Revision description	Draw. Init	Appr. date	Appr. Init
	ds-omny			ds-omny	
		Material:	Scale: 1:3	Format: A2	Tolerance: DS/ISO 2768- Surface treat.:
ID:	Description:				Rev:
	15890008 Cover for grindstone, assembly				C



30 Pos 30 to be adjusted 50 0/+5mm from the floor to the edge of the cabinet (see section A-A).





K	2013-06-14	F6: Pos. 200 removed. Text added.	JJO		
A	29-02-08		SPE	15.01.2006	FPB
Revision	Crea. date	Revision description	Draw. Init	Appr. date	Appr. Init
		Material	Scale 1:5	Format A1	Tolerance: DS/ISO 2768- m/k Surface treat: None
ID:	Description:				Rev:
	15890010 Casing, assembly				K




POS. NO.	1	2	3	4	5	6	7	8	9	10	11	12
1												
2												
3												
4												
5												
6												
7												
8												
9												
10												
11												
12												
13												
14												
15												
16												
17												
18												
19												
20												
21												
22												
23												
24												
25												

M	2012.05.16	Drawing updated, 2TJ10825 removed	SPE	2012.05.16	JTV
A	08.08.2006		JFR	08.08.2006	FPG
Revision	Crea. date	Revision description	Draw. Init	Appr. date	Appr. Init
	08.08.2006			08.08.2006	
		Material:	Scale:	Format:	Tolerance: DS/ISO 2768: mK
			1:2	A2	Surface treat.: None
ID:	Description:				Rev:
	15890070 Dresser, assembly				M

Produktion B4
 Dec 27th Ballerup/Copenhagen
 Denmark
 Phone: +45 44 600 800
 Fax: +45 44 600 804



- 180 Locked by Loctite 222 (Låses med Loctite 222)
- 80 100 Locked by Loctite 2701 (Låses med Loctite 2701)




ISOMETRIC VIEW

C	2010-04-06	Pos.180 added.	JTV	2010-04-06	
A	2006.08.08		JFR		
Revision	Crea. date yyyy-mm-dd	Revision description	Draw. Init	Appr. date yyyy-mm-dd	Appr. Init
		Material:	Scale: 1:2	Format: A3	Tolerance: DS/ISO 2768- mK Surface treat.: None
ID:	Description: 15890071 Moving part of dresser				Rev: C


Pødenstrupsvej 84
 DK-2750 Ballerup/Copenhagen
 Denmark
 Phone: +45 44 600 800
 Fax: +45 44 600 804

5 4 3 2 1


CONNECTION FOR 200V / 50Hz




CONNECTION FOR 200-210V / 60Hz




CONNECTION FOR 220-230V / 50Hz




CONNECTION FOR 220-240V / 60Hz




CONNECTION FOR 380-415V / 50Hz




CONNECTION FOR 380-415V / 60Hz



CONNECTION FOR 430-460V / 60Hz



CONNECTION FOR 460-480V / 60Hz




COLOR CODES:


- BK = BLACK
- BN = BROWN
- RD = RED
- OG = ORANGE
- YE = YELLOW
- GN = GREEN
- BU = BLUE
- VT = VIOLET
- GY = GREY
- WH = WHITE

Rev. A: Baan PDM		STRUERS A/S VALHOEJS ALLE 1176 DK-2610 ROEDOVRE DENMARK PHONE: + 45 3670 3500		
Transformer connections.				
FILE NAME.: 5093452.DSN PAGE1.SCH	Size A3	CAGE Code <Cage Code>	DWG NO 15093452	Rev A
Thursday, September 21, 2000	Scale	SLN / SLN		Sheet 1 of 3

Primær:




Sekundær:




| = intern forbindelse

Rev. A: Baan PDM	STRUERS A/S VALHØEJS ALLE 1176 DK-2610 ROEDOVRE DENMARK PHONE: + 45 3670 3500			
	Transformator Construction - electrical			
FILE NAME.: 5093452.DSN PAGE2.SCH	Size A4	CAGE Code <Cage Code>	DWG NO 15093452	Rev A
Thursday, September 21, 2000	Scale	 SLN / SLN	Sheet 2 of 3	



C	2010-05-17	BP1 and BP2 added	JTV	2010-05-17	JTV
B	15.5.2008	Counter pressure corrected from 2,8 to 3,0 bar	JTV	15.5.2008	JTV
A	23.8.2006		JF	23.8.2006	FPG
Rev	Crea. date dd-mm-yy	Revision description	Draw. Init	Appr. date dd-mm-yy	Appr. Init

<p>Pederstrupvej 84 DK-2750 Ballerup Copenhagen Denmark Phone: +45 44600 800 Fax: +45 44600 804</p>		Material:	Scale: 1:1	Format: A4	Tolerance: DS/ISO 2768 -	Rev: C
					Surface treat.:	
		ID: Description: 15892000 Air diagram AbraPlan-20				




Note:
Refer to drawing 15893100-dsn
for correct wiring.

- Inputs:**
- Emergency stop
 - Temperature disc motor
 - Temperature sample mover motor
 - Sample mover top position
 - Removal measurement unit installed
 - Absolute pressure button
 - Absolute pressure top
 - Voltage supply 1 (input)
 - Voltage supply 2 (+15V)
 - Voltage supply 3 (+24V)
 - Removal measurement
 - Motor current
 - Speed disc motor

- Outputs:**
- 7 x contactors:
 - Main
 - Sample mover 150 rpm
 - Disc motor 1450 rpm
 - Disc motor star
 - Disc motor delta
 - DC brake disc motor
 - 5 x 24VDC air valve:
 - Sample mover up
 - Air outlet 1
 - Air outlet 2
 - Sample mover down
 - Dresser up

REV A: OH (2006/10/23)		STRUER'S AS	
REV B: OH (2007/01/24)		Pederstrupvej 84	
- Issue FR added		DK-2750 Ballerup	
REV C: OHO (2011/11/03)		Denmark	
- F7, F8 added, wires reconnected		ABRAPLAN 20, BLOCK DIAGRAM	
Size A2	CAGE Code	DWG NO 15893050	Rev C
Scale	OHD		Sheet 1 of 1
Thursday, November 03, 2011			



SEE REQUIREMENTS FOR MAINS CABLE
- note #1

MUST BE EXTERNALLY FUSED
Fmax: See note #2

MOTOR COUPLING FOR M2

SINGLE SPEED MOTOR
Type: SEW R27DT71D4

400V 50Hz ∇

230V 50Hz Δ

ALL OTHERS
AS DELIVERED

COLOR CODE:


- BK = black
- BN = brown
- RD = red
- OG = orange
- YE = yellow
- GN = green
- BU = blue
- VT = violet
- GY = grey
- WH = white
- PK = pink
- BE = beige
- RO = rose

VOTAGE / FREQ. (from nameplate)	note #1 W1 - mains cable	note #2 max. ext. fuse	F1+F2+F3 (fuse size) See note #3	F4+F5 (fuse size) See note #3	F6 (fuse size)	F7+F8 (fuse size) See note #3	M1 Connection
3 x 200V / 50Hz	2,5mm ²	3 x 40AT	3 x 4AT (aM)	2 x 2AT (aM)	4AT	2 x 6AT (aM)	DELTA
3 x 200-210V / 60Hz	AWG = 12	3 x 40AT	3 x 4AT (CC)	2 x 2AT (CC)	4AT	2 x 6AT (CC)	DELTA
3 x 220-230V / 50Hz	2,5mm ²	3 x 40AT	3 x 4AT (aM)	2 x 2AT (aM)	4AT	2 x 6AT (aM)	DELTA
3 x 220-240V / 60Hz	AWG = 12	3 x 40AT	3 x 4AT (CC)	2 x 2AT (CC)	4AT	2 x 6AT (CC)	DELTA
3 x 380-415V / 50Hz	2,5mm ²	3 x 40AT	3 x 4AT (aM)	2 x 1AT (aM)	4AT	2 x 4AT (aM)	STAR
3 x 380-415V / 60Hz	AWG = 12	3 x 40AT	3 x 4AT (CC)	2 x 1AT (CC)	4AT	2 x 4AT (CC)	STAR
3 x 460-480V / 60Hz	AWG = 12	3 x 40AT	3 x 4AT (CC)	2 x 1AT (CC)	4AT	2 x 4AT (CC)	STAR

Note #3:
F1, F2, F3, F4, F5 are time delay fuses
CC...Class-CC characteristic
aM...aM characteristic

REV D: OHO (2011/11/03)
- F7, F8 added, wires reconnected
- F4, F5 value updated
- T2 Temp. monitoring added, new RD110
Rev E: FTH (2011-11-25)
- M1 connection label corrected


REV A: OH (2008/10/23) REV B: OH (2007/01/24) - Fuse F6 added - max. value for ext. fuse changed - variant overview table updated REV C: AKR (2010/02/22) - Core, changed to Star for 3x 200-210V / 60Hz	STRUER A/S Fædstrupvej 84 DK-2750 Ballerup Denmark		
ABRAPLAN 20, CIRCUIT DIAGRAM - MAIN VOLTAGE			
Size A2	CAGE Code	DWG NO	Rev
		15893100	E
Friday, November 25, 2011	Scale	OHO	Sheet 1 of 1



COLOR CODE:


- BK = black
- BK = brown
- RD = red
- OG = orange
- YE = yellow
- GN = green
- BL = blue
- VT = violet
- GY = grey
- WH = white
- PK = pink
- BE = beige
- RO = rose
- GNYE = yellow-green

REV D: CHD (2011/1/103) - see changes on page 4		STRuers A/S Pneumotekvej 44 DK-2730 Sønderlyng Danmark	
REV B: CHD (2007/01/24) - page 5 added - cable W7 changed - cable W7 changed - SMU's changed to embedded types.		ABRAPLAN 20, CIRCUIT DIAGRAM CPU - BOARD, A1	
Rev	Part No	Rev	Part No
A1	15893120	D	
Thursday, November 03, 2011		1 of 5	




COLOR CODE:
 BK = black
 BR = brown
 RD = red
 OG = orange
 YE = yellow
 GN = green
 BU = blue
 VT = violet
 GY = grey
 WH = white
 PK = pink
 BE = beige
 RD = rose
 GNYE = yellow-green

REV D: CHG 020111103 see changes on page 4		REV A: CH (2006/10/08) Interlocks4 REV B: CH (2007/01/24) - SB2 renamed to SB3-1 - SB3 renamed to SB3-2 - keyboard labeled with A7 - see changes on other pages - BR1 renamed to SA1		STRUINS A7 Interlocks4 020111103 Date	
ABRAPLAN 20, CIRCUIT DIAGRAM CPU - BOARD, A1					
Rev	Doc Code	Doc No	Rev	Doc Code	Doc No
A1					15893120
Date		OH		2	of 5



REV D: OHO (2011/11/03)
- see changes on page 4

REV A: OH (2006/10/08)		STRUERS A/S Pederstrupvej 84 DK-2750 Ballerup Denmark	
REV B: OH (2007/01/24)		ABRAPLAN 20, CIRCUIT DIAGRAM EARTH CONNECTIONS	
- page 5 added - see changes on other pages - GNYE4, GNYE7 changed - name of XPE6 changed			
Size	CAGE Code	DWG NO	Rev
A3		15893120	D
Scale		OH	Sheet 3 of 5



D

D

C

C

B

B

A

A

5

4

3

2

1



5

4

3

2


1



COLOR CODE:

- BK = black
- BN = brown
- RD = red
- OG = orange
- YE = yellow
- GN = green
- BU = blue
- VT = violet
- CY = grey
- PK = pink
- BE = beige
- RO = rose
- CNVE = yellow-green

REV: D: CHD (2011/11/03) - T2 Temp. monitoring added, new RD11	REV: A: CH (2006/10/08) - Discrepancy 64 - Discrepancy 65 - Page 5 added, RD22 added - S25 moved to SB3-3, SB3-3 chan. - 24 removed to B1 - New changes on other pages - All changed to RP1	STILERS A21 Discrepancy 64 Discrepancy 65 Discrepancy 66
ABRAPLAN 20, CIRCUIT DIAGRAM MACHINE CONTROL - BOARD		
Page: A1	Page: NO	Rev: D
15893120		
Monday, November 07, 2011	CH	4 of 5



D

D

C

C

B

B

A

A

5

4

3

2

1

5

4

3

2

1

COLOR CODE:
 BK = black
 BN = brown
 RD = red
 OG = orange
 YE = yellow
 GN = green
 BU = blue
 VT = violet
 GR = grey
 PK = pink
 BE = beige
 RO = rose
 GVE = yellow-green


REV: D: CHD (2011/11/03)
 - T2 Temp. monitoring added, new RD14

REV: A: CH (2006/10/08)
 - Discrepancy #4
 - Discrepancy #5
 - Page 5 added, RD22 added
 - 225 moved to SB3-3, SB3-3 chan.
 - 24 removed to B1
 - New changes on other pages
 - All changed to RP1

STILERS A2
 Drawings: A4
 D: 1/12 (12/08)
 Drawn:

**ABRAPLAN 20, CIRCUIT DIAGRAM
 MACHINE CONTROL - BOARD**

REV	DATE	BY	CHK
A1			
Rev	CH		



Color codes:
 BK = BLACK
 BN = BROWN
 RD = RED
 OR = ORANGE
 YE = YELLOW
 GN = GREEN
 BU = BLUE
 VT = VIOLET
 GR = GRAY
 WH = WHITE
 PK = PINK
 GW = GREEN/YELLOW

Wiring Diagram Contactors AbraPlan-20 Page 1/6	Overhaul: (1:2)	Material: 080507 BRV	Rev. order, total, after OS/CS 2706- 185893450-1A
--	---------------------------	--------------------------------	---

from F1, K5 (page 1)

from A2:X13 (page 3) / W31

T3

BK7

BK32

BK16

T2

See drawing 2MT72045

BK17

BK18

BK19

BK20

V1

GNYE4-XPE3

to XPE4
(page 1)

Color codes:

BK = BLACK
BN = BROWN
RD = RED
OC = ORANGE
YE = YELLOW
GN = GREEN
BU = BLUE
VT = VIOLET
GY = GREY
WH = WHITE
PI = PINK
GNYE = GREEN/YELLOW

Matr.:

Overlæb.:

Målfors.:

ikke ang. tol. efter
DS/ISO 2768-

 **Stuers**



1:2

Projektm. metode

Dato Sign.


Tegn. 110507 BRY

Konstr. 110507

Wiring Diagr. Contactor box
AbraPlan-20 Page 2/6

Erst.:


15893450-2A



Color codes:
 BK = BLACK
 BN = BROWN
 RD = RED
 OG = ORANGE
 YE = YELLOW
 GN = GREEN
 BU = BLUE
 VT = VIOLET
 GY = GREY
 WH = WHITE
 PI = PINK
 GNYE = GREEN/YELLOW

Matr.:	Overfl.beth.:	Wårforhold:	ikke ang. tol. efter DS/ISO 2768-
		1:2	Dato: _____
			Sign.: _____
Wiring Diagr. Control Board AbraPlan-20		Teg.: 11D507 ERFY	
Page 3/6		Erst.: _____	
		15893450-3A	


Color codes:
 BK = BLACK
 BN = BROWN
 RD = RED
 OG = ORANGE
 YE = YELLOW
 GN = GREEN
 BU = BLUE
 VT = VIOLET
 GY = GREY
 WH = WHITE
 PI = PINK
 GNYE = GREEN/YELLOW



to A2 (page 3)
 to A1 (page 4)
 to A2 (page 3)
 to A2 (page 3)

B: BUS junction and SMU added, 2009.11.16 SPE



Matr.:	Overfl.beh.:	Målførhold:	Ikke ang. tol. efter DS/ISO 2768-	
			1:2	Date
				Sign.
Wiring Diagram AbraPlan-20		Eerst:		15893450-5B
Page 5/6				





B. Star/Delta added motor M1. 2010-02-19 BRY/JTV

Matr.:	Overf.beh.:	Målførhold:	Ikke ang. tol. efter DS/ISO 2768-	
			1:2	
			Date: 160507 Sign: BRY Kontr.:	
Wiring Diagram AbraPlan-20			Page 6/6	
Erst.:			15893450-6B	

TRAFNO NO. AA-72045 (200-240V)




TRAFNO NO. AA-72045 (400-460V)



COLOR CODES:

- BK = BLACK
- BN = BROWN
- RD = RED
- OG = ORANGE
- YE = YELLOW
- GN = GREEN
- BU = BLUE
- WH = WHITE
- GY = GREY
- VT = VIOLET

Max. 10% overspænding pr. kobling.



REV A: OH (2006/10/23)	STRUERS A/S Pederstrupvej 84 DK-2750 Ballerup Denmark		
	WIRING DIAGRAM: WIRING AT BRAKE TRANSFORMERS 2MT72045B: POWER SUPPLY 200V - 480V 50/60Hz.		
Size A4	CAGE CODE	DWG NO	15893451
Tuesday, October 24, 2006	SCALE	24.10.2006 / OH	Sheet 1 of 1

Overview, variant parts in Abraplan-20:

BRY/JTV 2010-02-22: Star/Delta added.
 JTV 13-02-2007: Fusetype changed in 60Hz versions to Class CC.
 Ver. E: F7, F8 added. F4, F5 values updated. (OHO 2011-11-03)
 Ver. F: Reduced number of M1 motor variants (FTH 2011-11-25)
 Ver. G: M2: 2ME06405 and 2ME56486 merged into one new ICME version of 2ME56486 (FTH 2013-05-07)

Country Nom. voltage/frequency.	Motor M1	Motor M1 data:	Motor M2	Motor M2 data:	Transformer T2		Fuses F1, F2 & F3	Fuses F4 & F5	Fuses F7 & F8	Variant parts	Ordering number
					Type	Connections					
Japan: 3x200V 50Hz.	Item no.: 15499016 (Delta connection) Gear motor 169 rpm. Cable gland 2NM10472 (M25)	200V/50Hz (S1:) kW = 0,37 Amp = 2,3	Item no. 2ME06205 1440 rpm. Star/Delta start	240V/50Hz (Delta) (S1:) kW = 4,0 Amp = 15,4	Item no. 2MT72045 Ulveco nr. AA-72045 See diagram 15893451	200 V. I _{max} . 0,65A Item no. 2x 2XL30402	4 AT (Type: aM) Item no. 2FC10020	2 AT (Type: aM) Item no. 2FC10060	6 AT (Type: aM) Item no. 2FC10060	Belt wheel Item no. 2JE10125 2JE92028	Mains Fuses 25 Amp.
						230 V I _{max} . 0,65A Item no. 2x 2XL30402					
Norway: 3x220-230V / 50 Hz	Item no. 15499018 (Delta connection) Gear motor 169 rpm. Cable gland 2NM10472 (M25)	230V/50Hz (S1:) kW = 0,37 Amp = 2,15	Item no. 2ME56486 1430 rpm. Star/Delta start	400V/50Hz (Delta) (S1:) kW = 4,0 Amp = 9,4		400 V. I _{max} . 0,33A Item no. 1x 2XL30402	4 AT (Type: aM) Item no. 2FC10040	1 AT (Type: aM) Item no. 2FC10010	4 AT (Type: aM) Item no. 2FC10040	Mains Fuses 25 Amp.	
Europe: 3x380-415V / 50Hz.	Item no. 15499018 (Star connection) Gear motor 169 rpm Cable gland 2NM10472 (M25)	415V/50Hz (S1:) kW = 0,37 Amp = 1,24	Item no. 2ME56486 1430 rpm. Star/Delta start	400V/50Hz (Delta) (S1:) kW = 4,0 Amp = 9,4	400 V. I _{max} . 0,33A Item no. 1x 2XL30402	4 AT (Type: aM) Item no. 2FC10040	1 AT (Type: aM) Item no. 2FC10010	4 AT (Type: aM) Item no. 2FC10040	Mains Fuses 20 Amp.		

Overview, variant parts in Abraplan-20:

BRY/JTV 2010-02-22: Star/Delta added.
 JTV 13-02-2007: Fusetype changed in 60Hz versions to Class CC.
 Ver. E: F7, F8 added. F4, F5 values updated. (OHO 2011-11-03)
 Ver. F: Reduced number of M1 motor variants (FTH 2011-11-25)

Ver. G: M2: 2ME06405 and 2ME56486 merged into one new ICME version of 2ME56486 (FTH 2013-05-07)

Country nom. voltage/frequency	Motor M1	Motor M1 data:	Motor M2	Motor M2 data:	Transformer T1		Fuses F1, F2 & F3	Fuses F4 & F5	Fuses F7 & F8	Variant parts	Ordering number		
					Type	Connections							
UL / CSA: USA, Canada, Japan. 3x200-210V / 60Hz.	Item no.: 15499019 (Delta connection) Gear motor 168 rpm. Cable gland 2NM11027 (NPT ½")	208V/60Hz. (Delta) Hp = 0,37 Amp = 2,05	Item no. 2ME56206 1740 rpm. Star/Delta start	208V/60Hz. (Delta) Hp = 5,4 Amp = 15,9	Item no. 2MT72045 Ulveco nr. AA-72045	200 V. Imax. 0,65A Item no. 2x 2XL30402	4 AT (Class: CC)	2 AT (Class: CC)	6 AT (Class: CC)	Belt wheel Item no. 2JE10106 2JE91628	05896130. 3x200-210V/60Hz Mains Fuses 25 Amp.		
UL / CSA: USA, Mexico, S-Korea. 3x220-240V / 60Hz.	Item no. 15499020 (Delta connection) Gear motor 168 rpm. Cable gland 2NM11027 (NPT ½")	240V/60Hz. (Delta) Hp = 0,37 Amp = 1,76	Item no. 2ME06205 1740 rpm. Star/Delta start	220V/50Hz. (Delta) Hp = 5,4 Amp = 15,9		230 V Imax. 0,65A Item no. 2x 2XL30402	Item no. 2FC11040	Item no. 2FC11020	Item no. 2FC11060		05896136. 3x220-240V/60Hz Mains Fuses 25 Amp.		
UL / CSA: Canada, Brazil, S-Korea. 3x380-416V / 60Hz.	Item no. 15499020 (Star connection) Gear motor 168rpm Cable gland 2NM11027 (NPT ½")	415V/60Hz. (Star) Hp = 0,37 Amp = 1,01	Item no. 2ME06386 1740 rpm. Star/Delta start	380V/60Hz. (Delta) Hp = 5,4 Amp = 9,7		See drawing 15893451	400 V. Imax. 0,33A Item no. 1x 2XL30402	4 AT (Class: CC)	1 AT (Class: CC)		4 AT (Class: CC)	Cable gland Item no. 2GK20026 2GK20045	05896147. 3x380-415V60Hz Mains Fuses 20 Amp.
UL / CSA: USA, Canada. 3x460-480V / 60Hz.	Item no. 15499022 (Star connection) Gear motor 168 rpm Cable gland 2NM11027 (NPT ½")	480V/60Hz. (Star) Hp = 0,37 Amp = 0,88	Item no. 2ME56486 1730 rpm. Star/Delta start	480V/60Hz. (Delta) kW = 4.8 Amp = 9,5			460 V. Imax. 0,33A Item no. 1x 2XL30402	Artikel nr.: 2FC11040	Artikel nr.: 2FC11010		Artikel nr.: 2FC11040	05896154. 3x460-480V60Hz Mains Fuses 20 Amp.	



Pederstrupvej 84
DK-2750 Ballerup
Denmark