

CitoPress-10/-20



The fastest way to mounting and increased productive capacity



CitoPress is offering ultra-short mounting times and maximum user-friendliness. Speed is further improved by using the optional dosing system. The unique automatic dosing system enables dosing of a pre-set amount of mounting resin.

Main features

- Low mass mounting units for fast heating
- Automatic detection of cylinder dimension
- Automatic dosing system for faster filling of the cylinders and less spilling of resins
- Turn/push-knob and large display for intuitive operation
- On-screen Hot Mounting Application Guide for convenient information and minimum of errors

Additional features

- Extended temperature range: 80-180°C
- Remote signalling (flashing light when process is ready)
- Selectable Stand-by heat with sleep function
- Software options: Sensitive, Database
- Low build-height for easy visual control
- Triple-threaded, no-jam top closure
- Individually controlled parameters: temperature, force/pressure, heating time, cooling time
- Simple access for easy cleaning
- Top closure elevator for easy positioning
- Powerful heating unit (1000 W)
- Hydraulic system with high accuracy pressure control
- 3 preset cooling levels (Hi/Med/Lo)





CitoPress increases your capacity

CitoPress is a powerful mounting press, giving you high speed and maximum user-friendliness. The result is faster turnaround and increased productive capacity.

CitoPress is provided with a hydraulic system, where the pressure is controlled with very high accuracy, securing fragile specimens.

Available as CitoPress-20 with dual cylinder, and as CitoPress-10 with a single cylinder.

Automatic dosing system

Speed is further improved by using the optional dosing system. The unique automatic dosing system enables dosing of a pre-set amount of mounting resin. The advantages are faster filling of the cylinders and less spilling of resins.

How it works

A dosing unit is linked to a specific resin by a RFID chip, which is recognized by the mounting press. In turn, when using CitoDoser, the last used method is automatically stored and recalled when the relevant Doser is placed on the press.



CitoDoser: Automatic dosing option


With its flat base, CitoDoser is easily stored on the workdesk when not in use.

To be able to link CitoDoser and the resin inside with more than one method, it is necessary to purchase the Database option.

Large display and turn/push-knob

CitoPress-10/-20 feature a large, LCD display, which graphically shows parameter settings and mounting progress. Changing parameters is very fast and easy by using the turn/push-knob. Intuitive operation secures faster turnaround, a minimum of errors and the shortest possible set-up time.

CitoPress mounting times for different resins, compared to competing presses (competing presses = 100%)



CitoPress can reduce mounting times to nearly one third, compared to competing presses. (Applies for a 30 mm dia. cylinder. When using a 40 mm dia. cylinder CitoPress can reduce mounting times by half)



Intuitive operation: Push a single key to activate the automatic dosing system

Method parameters (temperature, pressure, heating time, cooling time) can be controlled individually, ensuring full flexibility in mounting.

The user-friendly turn/push-knob gives access to the following standard options:

Stand-by mode

Used for serial mounting. Maintains a minimum temperature of 20-60°C on the mounting unit, thus cutting down on the heating time. Unlike other presses, stand-by heating can be switched off.


Operation mode

Three different operation levels are available for selecting levels of user access:

- Production mode, gives access to mounting and selection of method only.
- Development mode, gives access to changing method parameters (plus Production mode access).
- The Configuration mode, gives full access, including access to configuration of set-up parameters.



Fast and easy selection: The large display and the turn/push-knob facilitate a fast selection of selection parameters



User-friendly software

On-screen Hot Mounting Application Guide

The built-in Hot Mounting Application Guide takes user-friendliness to a new level. You have the mounting parameters for all Struers resins right on the press. A detailed description of the Struers resins is also available on-line.

Select the Struers resin to be used, and CitoPress automatically checks the cylinder diameter used, adapts the appropriate method, after which it is displayed – ready for use. Everything you need to know about the resins and their usage – right at your fingertips!

Automatic detection of cylinder dimension

When changing the cylinder dimension, the press automatically detects the new size, and adjusts the method parameters accordingly. In this fashion, errors from an incorrect cylinder setting are eliminated.

Remote signalling

The mounting unit tower is equipped with a LED-bar, which lights up when a process is started, and flashes when the process is completed. Thus, even from

across the lab, the operator can tell when his mount is ready. The LED-bar may be supplemented by an acoustic alarm, which is available on both models.

Software options

In addition to the following software options are available:

- Sensitive option:
Dual phase heating and pressure, for fragile/porous specimens
- Database option:
Method database for up to 15 user-defined methods

User-friendly design

Low build height

Low build height means easy visual inspection of the mounting cylinder. This is relevant e.g. if tall, thin samples are to be mounted, and the operator wants to check if the sample still stands after adding resin. Another typical situation would be if multiple samples are to be included in the same mount.

Easy Cleaning

The door to the mounting unit tower swings open for easy access to the mounting unit.

Compact mounting units

CitoPress is equipped with compact, low-mass mounting units significantly reducing heating time. Available cylinder diameters: 25, 30, 40, 50 mm, 1¼" and 1½".

Simple top closure

All presses come with a top-closure elevator.

Any handling problems of the top closure as seen on other presses are solved by the top closure design.

The CitoPress top closure is with triple-thread. Thus, the entire process of mounting the top closure follows with a minimum of effort.



LED bar for remote signalling



Easy access means quick change of mounting units and easy cleaning



Top closure elevator for easy handling



ClaroFast is a transparent hot mounting resin, offering you transparent mounts with a clear view to all sample details

Saves money and resources

The benefit of shorter heating times is reduced energy consumption. Compared to competing presses, less than half the power is required to heat the mounting unit.

Recirculation cooling unit

CitoPress-10 and -20 are prepared for quick connection to a recirculation unit (Cooling System 3).

Three preset cooling levels

All CitoPress models come with three preset cooling levels. The main advantage is water saving, meaning that the MEDIUM mode can be used in most cases, and water usage reduced significantly.

Thermoplastic resins require far longer cooling times than thermosetting. With CitoPress, the cooling level can be set to match the resin. E.g., instead of reducing the cooling level by manually reducing the flow at the water tap, the level can be reduced by selecting the LOW mode.

Specifications

Cat. no:

CitoPress-20

05746127

Advanced, electro-hydraulic hot mounting press with two cylinders, which can be used individually. Built-in Hot Mounting Application Guide with automatic conversion of method parameters. Mounting units are ordered separately. Optional dosing unit CitoDoser and Base Unit (05796101+05796102), recirculation cooling unit (05766xxx), Database (05736900) and Sensitive mode (05736901).

CitoPress-10

05736127

Advanced, electro-hydraulic hot mounting press with one cylinder. Built-in Hot Mounting Application Guide with automatic conversion of method parameters. Mounting units are ordered separately. Optional dosing unit CitoDoser and Base Unit (05796101+05796102), recirculation cooling unit (05766xxx), Database (05736900) and Sensitive mode (05736901).

Mounting units

Consisting of heating/cooling unit, lower ram and top closure with upper ram.

Cylinder diameter 1½"	05786127
Cylinder diameter 1¼"	05786227
Cylinder diameter 25 mm	05786327
Cylinder diameter 30 mm	05786427
Cylinder diameter 40 mm	05786527
Cylinder diameter 50 mm	05786627

CitoDoser

05796101

Semi-automatic dosing unit with electronic ID-tag for automatic recognition by hot mounting press. Requires CitoDoser Base (05796102).

CitoDoser Base Unit

05796102

Base unit for use with CitoDoser (05796101).

CitoPress Database

05736900

Database option for CitoPress-10 and -20. Store up to 15 mounting methods

CitoPress Sensitive

05736901

Sensitive option for CitoPress-10 and -20. Dual phase heating and pressure for fragile/porous specimens.

Cooling system 3

See voltage

With 50 l tank (05766906), small pump (05766xxx), Cooli-1 (05761116) and 50 l static filter (05766908) with filter paper (05766914).

Voltage:	1 x 100 V / 50 Hz	05766516
	1 x 120 V / 60 Hz	05766523
	1 x 220 V / 50 Hz	05766522
	1 x 240 V / 60 Hz	05766524
	1 x 100-120 V / 50/60 Hz CSA	05766616

CitoPress with automatic dosing system for faster filling of the cylinders and less spilling of resins



Technical data

SUBJECT	SPECIFICATION
Mounting unit	
Mounting dia.	25, 30, 40, 50 mm, 1¼", 1½"
No. of mounting units	CitoPress-10: 1 CitoPress-20: 2
Mounting press	
Mounting parameters	
Pressure	50-350 bar / 725-5076 psi in steps of 25 bar / 363 psi
Heating time	1-15 min in steps of 0.5 min
Heating temperature	80-180°C / 176-356°F in steps of 5°C / 41°F
Cooling time	1-15 min in steps of 0.5 min
Cooling rate	High: Full flow Medium: 20% of full flow Low: 3% of full flow
Dosing rate*	20-150%
*) with optional CitoDoser	
Software and electronics	
Controls	Touch pads, turn/push-knob
LCD display with white LED backlight	320x240 dots
Noise levels	
Idle	0 dB (A)
Max.	56 dB(A)
Working environment	
Temperature, operational	5-40°C / 41-104°F
Humidity, non-condensing	0-95% RH
Supply	
Voltage / frequency	100-120V / 50-60Hz, 200-240V / 50-60Hz Auto detection and auto switch over
Power inlet (IEC320-C13)	1-phase (N+L+PE) or 2-phase (L1+L2+PE) <i>The electrical installation must comply with "Installation Category II"</i>
Power, idle	8W
Power, max.	CitoPress-10: 1300W@100-120V, 1300W@200-240V CitoPress-20: 1300W@100-120V, 2300W@200-240V
Current, max.	CitoPress-10: 13A@100-120V, 5.6A@200-240V CitoPress-20: 13A@100-120V, 10A@200-240V
Water pressure (tap water)	1-10 bar / 14.5-145 psi
Water inlet	¾"
Water outlet	ø13 mm
Dimensions and weight	
Width	CitoPress-10: 480 mm / 19", CitoPress-20: 550 mm / 21½"
Depth	560 mm / 22"
Height	450 mm / 17,7" (excl. CitoDoser) 550 mm / 21½" (incl. CitoDoser)
Weight	CitoPress-10: 34 kg / 75 lbs, CitoPress-20: 48 kg / 106 lbs CitoDoser: 3.1 kg / 7 lbs

Struers A/S

Pederstrupvej 84
DK-2750 Ballerup, Denmark
Phone +45 44 600 800
Fax +45 44 600 801
struers@struers.dk

**BELGIQUE (Wallonie)****Struers S.A.S.**

370, rue du Marché Rollay
F- 94507 Champigny
sur Marne Cedex
Téléphone +33 1 5509 1430
Télécopie +33 1 5509 1449
struers@struers.fr

BULGARIA**Struers GmbH**

Branch Office Bulgaria
21 Mur Str.
BG-1680 Sofia
hristo.ivanov@struers.de

CANADA**Struers Ltd.**

7275 West Credit Avenue
Mississauga, Ontario L5N 5M9
Phone +1 905-814-8855
Fax +1 905-814-1440
info@struers.com

CZECH REPUBLIC**Struers GmbH**

Organizační složka
Havlíčková 361
CZ-252 63 Roztoky u Prahy
Phone +420 233 312 625
Fax: +420 233 312 640
david.cernicky@struers.de

CHINA**Struers Ltd.**

Office 702 Hi-Shanghai
No. 970 Dalian Road
Shanghai 200092, P.R. China
Phone +86 (21) 5228 8811
Fax +86 (21) 5228 8821
struers.cn@struers.dk

DEUTSCHLAND**Struers GmbH**

Karl-Arnold-Strasse 13 B
D-47877 Willich
Telefon +49 (0) 2154 486-0
Fax +49 (0) 2154 486-222
verkauf.struers@struers.de

FRANCE**Struers S.A.S.**

370, rue du Marché Rollay
F- 94507 Champigny
sur Marne Cedex
Téléphone +33 1 5509 1430
Télécopie +33 1 5509 1449
struers@struers.fr

HUNGARY**Struers GmbH**

Magyarországi fióktelep
Puskás Tivadar u. 4
H-2040 Budaörs
Phone +36 (23) 428-742
Fax +36 (23) 428-741
zoltan.kiss@struers.de

IRELAND**Struers Ltd.**

Unit 25a
Monkspath Business Park
Solihull B90 4NZ
Phone +44 (0)121 745 8200
Fax +44 (0)121 733 6450
info@struers.co.uk

JAPAN**Marumoto Struers K.K.**

Takara 3rd Building
18-6, Higashi Ueno 1-chome
Taito-ku, Tokyo 110-0015
Phone +81 3 5688 2914
Fax +81 3 5688 2927
struers@struers.co.jp

NEderland/BELGIE**Struers GmbH Nederland**

Electraweg 5
NL-3144 CB Maassluis
Tel. +31 (0) 10 599 72 09
Fax +31 (0) 10 599 72 01
glen.van.vugt@struers.de

POLAND**Struers Sp. z o.o.**

Oddział w Polsce
ul. Lirowa 27
PL-02-387 Warszawa
Phone +48 22 824 52 80
Fax +48 22 882 06 43
grzegorz.uszynski@struers.de

ROMANIA**Struers GmbH**

Sucursala Sibiu
Str.Scoala de Inot, nr. 18
RO-550005 Sibiu
Phone +40 269 244 558
Fax +40 269 244 559
nicu-ovidiu.aunguroancei@struers.de

SCHWEIZ**Struers GmbH**

Zweigniederlassung Schweiz
Weissenbrunnenstrasse 41
CH-8903 Birmensdorf
Telefon +41 44 777 63 07
Fax +41 44 777 63 09
rudolf.weber@struers.de

SINGAPORE**Struers A/S**

627A Aljunied Road,
#07-08 BizTech Centre
Singapore 389842
Phone +65 6299 2268
Fax +65 6299 2661
struers.sg@struers.dk

SWEDEN**Struers A/S**

Ekbacksvägen 22, 3 tr
SE-168 69 Bromma
Telefon +46 (0)8 447 53 90
Telefax +46 (0)8 447 53 99
struers@struers.dk

UNITED KINGDOM**Struers Ltd.**

Unit 25a
Monkspath Business Park
Solihull B90 4NZ
Phone +44 0121 745 8200
Fax +44 0121 733 6450
info@struers.co.uk

USA**Struers Inc.**

24766 Detroit Road
Westlake, OH 44145-1598
Phone +1 440 871 0071
Fax +1 440 871 8188
info@struers.com

ÖSTERREICH**Struers GmbH**

Zweigniederlassung Österreich
Ginzkeyplatz 10
A-5020 Salzburg
Telefon +43 662 625 711
Fax +43 662 625 711 78
stefan.lintschinger@struers.de

Struers' equipment is in conformity with the provisions of the applicable International Directives and their apurtenant Standards. (Please contact your local supplier for details)

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice

www.struers.com