ViaFix

Instructions

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Mode d'emploi

#50300110

2016.11.01

 Struers

The Struers logo consists of a stylized 'S' shape formed by three slanted bars of increasing height from left to right, followed by the word 'Struers' in a bold, sans-serif font.

English

Description

ViaFix is a two component acrylic mounting system, consisting of one powder and one liquid. ViaFix is ideal for filling vias and microscopic voids, i.e. blind vias and microvias. The system cures at room temperature. To obtain the clearest mount the use of a pressure chamber is recommended. ViaFix is affected by alcohol. We recommend using **DiaPro** for preparation. The use of diamond products or lubricants containing alcohols will affect the surface and the polymer will appear darker. We also recommend cleaning the specimen with water and soap. When using cotton wool with alcohol for cleaning, there will be a risk of smearing the resin over the specimen.

Safety Data Sheet

Please download a copy of the **Safety Data Sheet** from the Library section on www.Struers.com and read carefully before use.

Access the Struers Library from the top menu on the website.

Please check the Library periodically for updates to the Safety Data information.

Storage

Acrylic powders absorb moisture. For maximum shelf life, store ViaFix Powder under dry and cool conditions less than 20°C / 68°F.

Acrylic liquids are sensitive to temperature. For maximum shelf life, store ViaFix Liquid under dry and cool conditions between 0°C / 32°F and 20°C / 68°F.

Mounting of Specimens

Measuring of Components

The ratio for ViaFix Powder:ViaFix Liquid is **11:9** by weight and **2:1** by volume. Weighing is recommended to obtain the best mounts and for reproducibility.

Important!

Do NOT pour excess liquid back into the bottle as this will significantly reduce the shelf-life.

Mixing of Components

- Measure out the amount of ViaFix Powder into a mixing cup (max. 55 g / 5 x 20 ml).
- Measure out the amount of ViaFix Liquid into a mixing cup (max. 45 g / 5 x 10 ml).
- Sprinkle the powder into the liquid while stirring.
- Stir thoroughly until the mixture is homogeneous, without introducing air bubbles into the mixture.
Immediately after the components are mixed polymerisation of the resin starts and the viscosity will increase. To get the best filling use the resin as soon as possible after mixing.

Selecting the Mounting Cup

All Struers mounting cups can be used.

Pre-treatment of specimens

Ensure that the specimens are thoroughly cleaned and dry prior to mounting. Each surface must be free from grease, moisture and other contaminants to promote the best possible adhesion of ViaFix to the specimen.

Mounting

- Place the cleaned specimen in the centre of the mounting cup, mix the necessary amount of resin and carefully pour over the specimen.
- Ensure that no air remains trapped within cracks or cavities in the specimen by moving the specimen or by tapping the mounting cup on the table.

Curing

ViaFix is a room temperature curing system. For a 30 mm dia. mount, the curing time is about 20 minutes. It is recommended to always use a pressure chamber in order to obtain a clear mount with no air trapped within the resin. When using a pressure chamber the temperature in the chamber will normally rise during curing, especially if many mounts are cured at the same time. We recommend keeping the mount in the mounting cup until the temperature is below 30°C/86°F. Above this temperature the mount is soft, the surface sticky and is difficult to take out of the cup.

When specimens with thin layers (d 2 mm) are mounted with ViaFix, the mount must be cured at 40°C/104°F.

After Curing

After removal from the mounting cups, the mounts are immediately ready for preparation. The top surface of the mount will be rubbery for approx. 24 hours. If the mounts are prepared as single specimens the top surface must either be ground for approx. 10 s on #220 SiC-paper or covered by a plastic Dymo® or the like to prevent the pressure feet from sticking to the mounts.

NB! Remember that ViaFix is affected by alcohol.

Technical Data

	Components	Ratio by weight	Ratio by volume
ViaFix	ViaFix Powder	11	2
	ViaFix Liquid	9	1

NB For the best mount use weighing method.

Density ViaFix Powder: 0.71 g/ml
 ViaFix Liquid: 0.89 g/ml

<i>Curing</i>	Can be prepared (ground and polished) after:	20 min.
	Peak temperature:	115°C / 240°F
	Pot life:	2 min.

30 mm dia. mounting cup without specimen at 21°C / 70°F.

Properties after Curing Hardness: 83 Shore D
 Density: 1.1 g/ml

Colouring ViaFix can be coloured with EpoDye (fluorescent dye) and AcryDye.

Tips – Some materials have very poor adhesion to acrylics which can result in a gap between the specimen and resin. To minimise this gap the mount can be refilled with ViaFix after the first preparation cured in pressure chamber and prepared again.

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