Accutom-2

Manual No.: 14277001

Date of Release 11.10.1993

Instruction Manual

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Mode d'emploi

Spare Parts and Diagrams

Table of Contents	Page
Instruction Manual	1
Gebrauchsanweisung	29
Mode d'emploi.....	58
Spare Parts and Illustrations.....	87

*Accutom-2
Instruction Manual*

Instruction Manual

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2003.

Struers A/S
Pederstrupvej 84
DK- 2750 Ballerup
Tel 44 600 800
Fax: 44 600 801Denmark

Accutom-2 Safety Precaution Sheet

To be read carefully before use

1. The operator(s) should be fully instructed in the use of the apparatus according to the Instruction Manual. The operator(s) should be fully instructed in the use of cut-off wheels for Accutom-2.
2. Use only intact cut-off wheels. See the manual, section 7.3. If Struers cut-off wheels are not used, it should be noticed whether they contain materials which give off harmful vapors or dust requiring special ventilation.
3. If Struers cut-off wheels are not used, the wheels must be approved for spindle speed of minimum 3000 rpm.
4. Make sure to have suitable ventilation when cutting materials which give off harmful dust particles or vapors. See Safety Data Sheet for the materials in question.
5. Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the cutting fluid in question.
6. Store the cutting fluid and use it as indicated in the instructions for the cutting fluid in question so that harmful bacterial growth and contact with the skin is avoided.
7. Do not start the apparatus until the wheel guard is mounted. See the manual, section 6.6. Do not remove the wheel guard until the cut-off wheel has stopped completely.
8. If you use the Accutom-2 CATAS stand you must check that the insulation rubber cap on the switch is intact before start of the apparatus.
9. Always use safety goggles in connection with manual cutting. Use only 511CA cut-off wheels for manual cutting. See the manual, section 6.9.
10. The work piece must be clamped securely. See the manual, section 6.5.2.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.)

Table of Contents	Page
1. Description	
Application.....	5
Function	5
Construction	5
Accutom-2.....	5
Recirculating cooling unit.....	6
2. Technical Data	
Accutom-2	7
Recirculating Cooling Unit.....	7
3. Accessories and Consumables	8
4. Contents of Packing Box	10
5. Installation	
Placing.....	11
Set-up.....	11
Accutom-2 with recirculating cooling unit	11
Accutom-2 without recirculating cooling unit	11
Connection to water supply	12
Connection, electrical.....	12
Accutom-2.....	12
Recirculating cooling unit.....	12
Stand for rotating/swivelling specimen holder (CATAS).....	12
6. Operation	
Operating Elements	13
Stands	14
Stand for fixed specimen holder (CATYS).....	14
Stand for rotating/swivelling specimen holder (CATAS).....	14
Cut-off Wheels	15
Choice of Cut-off Wheel.....	15
Mounting of cut-off wheel.....	15
Choice of Cutting Speed	16
Specimen Holders.....	16
Choice of Specimen Holder	16
Mounting of specimen and adjustment of specimen holder	17
Preparation of automatic cutting	18
Automatic cutting.....	19
After cutting	20
Manual cutting	21
Mounting of manual cutting table (CATMA).....	21
Cutting.....	21

7. Maintenance

General cleaning	22
Cleaning of waste-tray	22
Storing and maintenance of diamond cut-off wheels	22
Storing of bakelite bonded Al ₂ O ₃ cut-off wheels	23
Cutting in of fuse	23
Maintenance of stands and specimen holders	24
Cutting stands	24
Specimen holders	24
Change of cooling water in the recirculating cooling unit.....	25

8. Service

Cleaning of the cutting table bearings.....	26
Ventilating the hydraulic system	26
Fault in the hydraulic system	27
Changing of gear ratio.....	27
Belt tightening	28

1. Description

Application

Accutom-2 is a laboratory precision cut-off machine for cutting all types of materials of small dimensions. According to accessory it can be used for automatic cutting, fast cutting, cutting with particularly low deformation, cutting with a minimum of burrs, or for manual cutting.

Function

The specimen to be cut is clamped in a specimen holder mounted in a stand. This can move the specimen rotating or swivelling, or it can retain the specimen stationary.

The specimen is moved towards the rotating cut-off wheel by means of a linear movable cutting table at a speed which can be adjusted by a continuously variable hydraulic braking device.

The speed of the cut-off wheel is continuously variable within the ranges 100-1000 rpm and 300-3000 rpm.

Heat which is generated from the cutting is removed with cooling water applied abundantly.

A recirculating cooling unit makes it possible to cool with other agents than tap-water.

Construction

Accutom-2

The spindle for the cut-off wheel is a precision shaft suspended in ball bearings, trued up accurately in relation to the cutting table.

The spindle is operated via 2 V-belts by a DC motor with continuously variable speed. The V-belt discs have 2 gear ratios whereby speed ranges for the spindle of 100-1000 rpm or 300-3000 rpm can be obtained.

The motor is protected with a fuse and a thermal switch.

The cutting table is suspended in linearly ball bearings. The feeding (movement) of the cutting table is done with 2 variable springs. The movement is regulated by the hydraulic braking device and can be varied continuously. The hydraulic pressure in the braking system can be read on a pressure gauge.

The cutting table activates a motor stop which makes automatic stop of the cutting possible. The movement of the cutting table can be monitored on a dial gauge.

Feeding of cooling water takes place through 2 tubes which cool and clean the cut-off wheel on both sides.

The cooling is best carried out when connecting the recirculating cooling unit, which belongs to the Accutom-2. At the same time, this opens up the possibility of using special cooling agents. However, cooling can also be carried out through connection to a water tap.

Recirculating cooling unit

The recirculating cooling unit is made with a pump, the motor of which is supplied with voltage from Accutom-2. This ensures synchronous start and stop. The 12 l vessel functions as sedimentation vessel, as the pump takes the cooling liquid from the upper layers.

The pump is resistant to water, to the current synthetic additives as well as to oil and paraffin oil based cutting liquids.

2. Technical Data

Accutom-2

<i>Voltage</i>	115/220 V, 50-60 Hz
<i>Effect</i>	95 W
<i>Cutting speed</i>	100-1000 rpm, alternatively 300-3000 rpm
<i>Feeding speed</i>	0-3 mm/s
<i>Feeding force</i>	Max. 40 N
<i>Cutting length</i>	Max. 75 mm, measurable with dial gauge max. travel 5 mm, resolution 0.01 mm
<i>Diameters of cut-off wheels</i>	50 to 175 mm
<i>Shaft diameter for cut-off wheel handle</i>	12.7 mm
<i>Water connection</i>	8 mm tube
<i>Dimensions, with splash screen</i>	Width: 420 mm Depth: 350 mm Height: 270 mm
<i>Weight</i>	21 kg

Recirculating Cooling Unit

<i>Voltage</i>	115/220 V, 50-60 Hz, 66VA
<i>Voltage supply</i>	From Accutom-2
<i>Volume</i>	12 l
<i>Dimension</i>	Width: 415 mm Depth: 515 mm Height: 120 mm

3. Accessories and Consumables

Accessories

Specification	Code word
<i>Stand</i> with dovetail clamp for specimen holder and micrometer screw. With rotating or rocking movement of specimen during cutting.	CATAS
<i>Stand</i> with dovetail clamp for specimen holder and micrometer screw.	CATYS
<i>Stand Extension</i> to raise stand (CATYS or CATAS).	CATIF
<i>Specimen Holders for Accutom-2</i> For general use. Vice type with max. opening 40 mm. For round or square specimens. Teardrop type with max. opening 25 mm /Ø 1" For irregular specimens, with 7 screws. with goniometer. For adhering specimens. With ceramic vacuum check for thin sections. for small specimens. Vice type.	CATCH CATEL CATIN CATNI CATPE CATAP CATAN
<i>Base Plate</i> with dovetail. For mounting of other types of specimen holders.	CATFU
<i>Manual Cutting Table</i> complete with splash guard and guide rail.	CATMA
<i>Open splash guard, U-shaped</i> Used together with stand extension CATIF and at long, protruding specimens	427 MP 2

Consumables

Cut-off wheels Application	Diameter x thickness mm	Bond	Abrasive	Code word
General use	125 x 0.35	Soft metal	Diamond	330CA
Ceramics, sintered carbides, hardened steel	125 x 0.5	Bakelite	Diamond	352CA
General use, smaller specimens	100 x 0.3	Soft metal	Diamond	230CA
High precision, costly materials	75 x 0.15	Soft metal	Diamond	100CA ¹
General use, metals	125 x 0.5	Bakelite	Al ₂ O ₃	357CA
Medium hard and hard steel	125 x 0.5	Bakelite	Al ₂ O ₃	356CA
Extremely hard ferrous metals	125 x 0.5	Bakelite	CBN	355CA
Manual cutting	175 x 0.8	Metal	Diamond	551CA

¹ Should be used with flanges Ø 65 mm

Specification	Code word
<i>Additive for Cooling Fluid</i>	ADDUN
<i>Bacterial Cleaning Liquid for Recirculating Cooling Unit</i>	BACTE

4. Contents of Packing Box

Remove the upper clamps on the transportation box with a screwdriver, so that the lid can be removed. Open the plastic bag and remove all loose parts. Take off the lower clamp and remove the casing. Remove the 4 bolts which keep the wooden bottom in place, and the apparatus is ready to be installed.

Check that besides the Accutom-2 basic unit the following loose parts are present:

- 1 Splash guard, plastic
- 1 Bottom part for same
- 1 Flange set, Ø 65 mm
- 1 Flange set, Ø 42 mm
- 1 Screw for spindle shaft
- 1 Dial gauge
- 1 Two-pin spanner
- 1 Hexagon key, 2 mm
- 1 Hexagon key, 2.5 mm
- 1 Hexagon key, 3 mm
- 1 Hexagon key, 4 mm
- 1 Hose connection
- 2 Hose clips S13
- 1 Reinforced PVC tube, inside Ø 5/16", 2 m
- 1 Latex tube, inside Ø 8 mm, 2 m
- 1 O-belt

5. Installation

Placing

Accutom-2 should be placed on a stable table at a suitable working height.

If the precision should be utilized fully, Accutom-2 should not be placed near strongly vibrating or noisy machines.

If Accutom-2 is used without recirculating cooling unit it should be placed near to or less than 1 m from water supply and outlet. The outlet must be placed at least 20 cm below the level of Accutom-2 in order to ensure good drainage of the cooling water.

Set-up

Accutom-2 with recirculating cooling unit

Place the recirculating cooling unit and then Accutom-2 on top of it, with the 4 rubber legs placed in the appropriate slots. Mount the L-outlet tube so that it points vertically downwards.

Accutom-2 is now adjusted sideways so that the outlet tube can move freely without touching the wall of the recirculating cooling unit along the whole travel of the cutting table (fig. 2.6).

The easiest way to obtain this is to place the recirculating cooling vessel against the right wall of the cabinet and then to mount the L-outlet tube.

The recirculating cooling vessel can now be filled with cooling liquid. See chapter 7.7.

Connect the pump to the cooling water supply on the rear side of Accutom-2 as follows:

- Put the nut on the hose and connect the hose to the L-tube section by means of the nut (fig. 2.7).
- The pump is supplied electrically when connected to the rear side of Accutom-2. Unscrew the blind flange and put the plug in (fig. 2.8).

Accutom-2 without recirculating cooling unit

Place Accutom-2 in the chosen place and mount the Latex tube on the outlet branch and lead it to a drain.

If there is no drain a container may be used. However, attention must be drawn to the fact that it must be emptied regularly.

It must be ensured that the height difference is more than 20 cm and that there are not sharp bends as well as that it slopes towards the drain.

Connection to water supply

- Dismount the L-tube section which is meant for use together with the recirculating cooling unit.
- Place the pressure reducing section in the pipe branch and mount it on the inlet valve of Accutom-2 instead of the L-tube section (also see fig. 3).
- Place a 13 mm tube clip on the reinforced PVC-tube and place the tube on the tube branch. Tighten the tube clip so that the joint becomes tight.
- Place the other end of the hose on the water tap in a similar way. If necessary, shorten the hose.

IMPORTANT

Always turn off the water supply at the water tap when Accutom-2 is not in use.

Connection, electrical

IMPORTANT

Before connecting the apparatus check that the operating voltage of the apparatus corresponds to the mains voltage.

The operating voltage is indicated on the tag at the mains voltage cable.

Accutom-2

Mount a plug on the mains cable:

Yellow/green: earth
Blue: neutral
Brown: phase
(blue and brown lead are freely interchangeable).

or Green: earth
Black: neutral
White: phase
(Black and white lead are freely interchangeable).

Recirculating cooling unit

The recirculating cooling unit is meant for use together with Accutom-2 and will be supplied from there. The connection takes place by connecting the plug of the pump on the rear side of Accutom-2.

Stand for rotating/swivelling specimen holder (CATAS)

Mount a plug on the mains cable:
Yellow/green: earth
Blue: neutral
Brown: phase
(blue and brown lead are freely interchangeable).

or Green: earth
Black: neutral
White: phase
(Black and white lead are freely interchangeable).

6. Operation

Operating Elements

Knob marked SPEED (fig. 1.1)

Turn the knob slowly clockwise to start the spindle. When the spindle has to stop turn the knob towards OFF.

Lamp marked TABLE STOP OR OVERLOAD (fig. 1.2)

When this lamp lights the motor has stopped. This may be because both the spindle and the cutting table have been stopped by the STOP SCREW and the cutting therefore is finished, or because the motor has run warm due to thermal overload.

IMPORTANT

The motor cannot start until the knob SPEED has been on OFF.

Speed indicator (fig. 1.3)

Indicates the spindle speed in rpm. It has two scales: 0-1000 and 0-3000 depending on which gear ratio has been chosen.

Knob marked FEED (fig. 1.4)

Is used to adjust the speed of the cutting table. When the knob is turned counter-clockwise the speed is increased.

Button marked FORWARD (fig. 1.5)

Press the button for fast forward movement of the cutting table.

Pressure gauge (fig. 1.6)

Measures the actual pressure in the hydraulic braking circuit. Is used to monitor the cutting.

Thermal fuse (fig. 1.7) and lamp (fig. 1.8)

The lamp indicates whether there is power on the apparatus. If the fuse has been cut-out it can be cut-in again by pressing the button.

Stop screw (fig. 1.9)

Stops the cutting table mechanically and switches off the motor electrically. The lamp TABLE STOP OR OVERLOAD (fig. 1.2) lights.

Adjust the screw before the cutting starts so that cutting table and motor are stopped when requested (when the specimen has been cut through, or when the requested cutting depth has been obtained).

Hasp (fig. 1.10)

Retains the cutting table in its starting position.

Rod for dial gauge (fig. 1.11)

Adjusts the working area of the dial gauge. The rod is held by friction. The size of the friction is set with the hexagon-headed screw. The friction is increased by turning the screw clockwise. The friction is adjusted so that the rod can be moved by pushing it manually, without loosening the screw.

Chains (fig. 2.1)

The force on the cutting table is increased by pulling the chains backwards. The drive must be symmetric.

*Lock screw on stand CATAS
(fig. 2.2)*

Locks the micrometer screw so that it is retained at a certain measurement, when the shaft of the stand is in motion. It is locked by turning clockwise.

*Lock screw on the stands CATAS
and CATYS (fig. 2.3)
Handle on cutting table (fig. 1.12)*

Locks all movement of the dovetail holder.

Is used to bring the cutting table back in its starting position.

Flange set (fig. 2.4)

Is used for mounting the cut-off wheel. Comprises inner flange, flange, flange screw. They are mounted in that order.

*Valve for regulating cooling liquid
(fig. 2.5)*

Stands

*Stand for fixed specimen holder
(CATYS)*

When cutting uncomplicated specimens where the requirement to surface roughness and cutting speed do not have first priority, this stand can be recommended. The position of the specimen holder can be adjusted horizontally with the micrometer screw. It can be turned in the same level as the cut-off wheel but it must be clamped before cutting.

*Stand for rotating/swivelling
specimen holder (CATAS)*

This stand will always provide a better cutting, less surface roughness as well as higher precision and speed than CATYS. In many cases, especially hard and ductile materials would be impossible to cut without this possibility of moving the specimen.

When cutting with the Al_2O_3 cut-off wheels 356CA and 357CA the rotating movement will reduce the friction of the cut-off wheel, thus providing the possibility of cutting bigger specimens than when using CATYS.

Circular, square or irregularly shaped specimens are best cut with the rotating movement.

Oblong, rectangular and ellipse-shaped specimens are best cut with the swivelling movement.

CATAS can be used as CATYS by just turning off the motor and clamping the holder like on CATYS.

Cut-off Wheels

Choice of Cut-off Wheel

Accutom-2 accepts cut-off wheels with a diameter of 50 mm (smallest flange is 42 mm) up to 175 mm. Hole 12,7 mm.

For cutting metallic materials the fastest and most economical cut is obtained with the Al_2O_3 cut-off wheels 356CA or 357CA. However, the specimens must not contain phases harder than 1700 Vickers, if these phases have to be cut.

Diamond cut-off wheels are used for all harder materials or materials containing very hard phases, carbides, nitrides etc.

For hard ductile materials like sintered carbides or nitrides diamond cut-off wheels with bakelite bonding like 351CA are used.

For hard brittle materials such as geological materials, electronic components or crystals diamond cut-off wheels with metal bonding 100CA, 230CA, 330CA and 551CA are used.

In order to obtain the best possible surface quality as thin cut-off wheels as possible, with a high concentration of diamonds and small grain sizes are used. 100CA, 230CA.

For general use and routine inspection a cut-off wheel with low concentration and coarse grains like 330CA or 351CA are used. These cut-off wheels, and especially 551CA, are well suited for manual cutting.

In order to obtain the highest possible accuracy as big flanges as possible should always be used (65 mm).

Mounting of cut-off wheel

Having chosen a cut-off wheel and a flange these are mounted as follows:

- Turn the right cooling tube to the right in order to facilitate access to the spindle.
- Put on the inner part of the chosen flange (fig. 1.13).

IMPORTANT

The fitting between spindle and inner flange being tight, these 2 surfaces must be absolutely clean.
If the inner flange does not fit directly, under no circumstances try to press it on, as this may damage the spindle. See if there should be small burrs and if so, take them off with grinding paper, grain size 1200.

- Mount the cut-off wheel.
- Put outer flange on (fig. 1.14).
- Mount the flange screw (fig. 1.15) and tighten with a two-pin spanner. It is sufficient just to hold on to the flanges when tightening the flange screw.
- After cutting the cut-off wheel should be dismantled in order to prevent corrosion on the sides of the cut-off wheel. The screw may be a little tight because of self-tightening, therefore a backstop with a hexagonal key at the left side of the spindle may be a possibility.

IMPORTANT

In any case it must be ensured that keys in question are moved right to the bottom of the holes before loosening or tightening.

Choice of Cutting Speed

On Accutom-2 it is possible to use 2 gear ratios between motor and spindle. This results in two ranges of speed of the cut-off wheel, viz. 100-1000 rpm and 300-3000 rpm. The switch is carried out manually by changing belts. (See section 8.3).

As a general rule, the range 300-3000 rpm should always be used, as this gives the fastest cutting, the smallest surface roughness and the best planeness .

The torque of the shaft being three times larger at the range 100-1000 rpm, this range should be used with 356CA and 357CA, and 551CA at manual cutting.

Specimen Holders

Choice of Specimen Holder

The following specimen holders are available as accessories:

- Specimen holder with parallel vice (CATCH). Opening between the screws allows a specimen width of max. 26 mm. Screws of various lengths are included so that it is possible to clamp specimens with lengths of up to 40 mm.
- Specimen holder with drop shaped cut-out (CATEL). Max. specimen width 25 mm. For clamping of circular and rectangular specimens. Clamping with one screw.
- U-shaped specimen holder with 7 screws (CATIN). Max. specimen width 32 mm. For clamping of very irregular specimens.
- Specimen holder with goniometer device (CATNI). 1 x $\pm 15^\circ$ ($\pm 0,1^\circ$), 1 x $\pm 20^\circ$ ($\pm 0,1^\circ$) and 1 x 360° (without dial). Can be displaced horizontally by means of the micrometer screw of the cutting stand. The specimen is affixed on a circular plate, \varnothing 13 mm.

- Specimen holder for affixation of specimen (CATPE). The block has a surface of 40 x 30 mm, suitable for affixation of the specimen with double adhesive tape, thermoplastic or ordinary glue.
- Specimen holder with vacuum holder (CATAP). The vacuum holder has the dimensions 30 x 50 mm. Well suited for vacuum clamping of plane specimens, e.g. thin sections on a glass plate. Without vacuum pump.
- Specimen holder with parallel vice, small (CATAN). As CATCH. However, max. dimensions between the screws 16 mm.
- Dovetail plate for mounting of user's own specimen holder (CATFU). Length 60 mm.

*Mounting of specimen and
adjustment of specimen holder*

The specimen must be clamped/be supported in the specimen holder so that it moves as little as possible during the cutting. Inconvenient movements may lead to damaging of the cut-off wheel and the specimen, and under all circumstances cause a bad cutting quality.

When cutting with the rotating or swivelling cutting stand the specimen should be placed so that by rotation or swivelling it moves as little as possible, i.e. it moves on its own axis. Thus the fastest cutting is obtained, as most of the time the cut-off wheel will be cutting.

The movement of the cutting table is so precise that an asymmetric clamping of the specimen will not affect the cutting quality, it only gives a longer cutting time.

When using the stand for fixed specimen holder the specimen should be placed so that the cut-off wheel cuts through the smallest cross section.

It is advantageous to place the specimen so that as little a variation of the cross section as possible is obtained.

Notice that the cut-off wheel must be able to cut through the specimen without the latter touching the flanges. Normally, this will not provide any problem when using diamond cut-off wheels, as they do not wear.

Depending on the specimen to be cut, the wear can vary with the cut-off wheels 356CA and 357CA. Therefore, the condition of the cutting should be checked at regular intervals.

The specimen should protrude at least 6 mm from the specimen holder which might otherwise come in touch with the flanges.

It is especially important to pay attention to this fact when the rotating/swivelling stand is used.

Preparation of automatic cutting

Having chosen the cut-off wheel, mounted the specimen, adjusted the specimen holder, and levelled the stand, the cutting can now be started according to the following steps:

- Check that there is sufficient cooling agent in the recirculating cooling unit (see section 7.7 for re-filling). If a recirculating cooling unit is not mounted then check if the water supply is opened.
- Adjust the stop screw (fig. 1.9) so that Accutom-2 stops when the specimen has been cut through or the desired cutting depth has been reached.
 - a) Turn on the Accutom by turning the SPEED knob (fig. 1.1) and let the spindle rotate at low rate of revolution.
 - b) Push button FORWARD (fig. 1.5) and hold it until the cutting table stops and the lamp TABLE STOP OR OVERLOAD (fig. 1.2) lights.
 - c) Hold the FORWARD button in and adjust the spacer screw to the desired cutting depth.
 - d) Turn the SPEED knob back to OFF to re-couple the switch.
- Adjust the placing of the cut by turning the micrometer screw. This can be done most easily by moving in the cutting table so that the specimen almost touches the cut-off wheel. If necessary, a cut can be made so that a reference surface is then obtained. If the cut is not all the way through there is furthermore the possibility of measuring accurately the actual width of the cut.

IMPORTANT

Both lock screws (fig. 2.2 and 2.3) must be loose.

- Lock the travel of the specimen holder in horizontal direction by: either locking the micrometer screw (fig. 2.2) if rotating/swivelling cutting is requested or locking the specimen holder tightly (fig. 2.3) if direct cutting is requested.
- Mount the protective screen of clear PVC.
- Start the Accutom-2 by means of the SPEED knob (fig. 1.1) and set to the requested speed.
- Adjust flow of cooling agent (fig. 2.5) to requested quantity, usually maximum.
- Make sure that the discharge is not clogged.
- Press the FORWARD button (fig. 1.5) to place the specimen closely in front of the cut-off wheel.
- Turn the FEED knob a little counter-clockwise (opens the valve) (fig. 1.4) and read out the pressure P_1 on the pressure gauge (fig. 1.6).

- Let the cut-off wheel cut its way a little into the specimen before reading out the actual pressure P_2 . Then open the FEED valve a little more until $\Delta P = P_1 - P_2$ is suitable.
- The cutting now takes place automatically so that the Accutom can be left.

Only if the cross section of the specimen varies greatly, or if cutting is done over a longer distance (30-40 mm) may it be necessary to change P_2 or to adjust the pull from the springs (fig. 2.1).

Automatic cutting

The size ΔP equals $P_1 - P_2$, where P_1 is the pressure in the brake-circuit, read before the cutting, and P_2 is the pressure read during the cutting.

ΔP indicates whether the values for FEED, speed of the cut-off wheel, have been set correctly, as well as if the cut-off wheel and feeding of cooling agent are correct.

ΔP is the expression for the force which affects the cut-off wheel.

If the value is high it means that the cut-off wheel does not manage to cut itself free. Therefore, the cut-off wheel is exposed to a higher pressure and it may bend, which would result in a poorer cutting quality.

Because of the friction a higher pressure will also mean more generation of heat.

The following ΔP values can be recommended:

ΔP : 0,05-0,15

Very plane parallel cutting ($< 2^\circ$) with very fine surface quality.

ΔP : 0,15-0,25

Good plane parallel cutting ($< 1^\circ$) and acceptable surface quality, i.e. "normal" cutting.

ΔP : 0,25-0,35

It should be considered to change the cutting parameters, as this range will normally give a poor result.

$\Delta P = \text{over } 0,35$

The cutting parameters should be changed: rotating/swivelling cutting, higher cutting speed, lower FEED, softer cut-off wheel, coarser grains in the cut-off wheel.

ΔP is increased by 0.1 if the diamond cut-off wheel is used for very ductile materials (sintered carbides, metals) or if Al_2O_3 cut-off wheels are used for metals with hard phases.

If P_2 , i.e. the pressure of the brake-circuit during the cutting, moves towards 0, the table is not moving. I.e. the pressure from the springs rests on the cut-off wheel, however, without any cutting. This is because the cut-off wheel does not cut itself free during the cutting. Therefore, one has the choice of cleaning the cut-off wheel often (see section 7.3), or to change to a softer binding.

ΔP may change during the cutting, especially if the specimen has a very varying cross section, or if the cutting is done over a length of 30-40 mm. In such cases, therefore, it may be an advantage to check ΔP by reading out P_2 .

After cutting

The cut specimen usually falls down on the cutting table, but with very thin specimens and low cutting speed the specimen may stick on the side of the cut-off wheel.

If the spacer screw has been adjusted correctly, the spindle will have stopped and the lamp TABLE STOP OR OVER LOAD (fig. 1.2) will be lit. Otherwise the motor is switched off by turning the SPEED knob (fig. 1.1) to OFF. Pull the cutting table back and lock it with the hasp (fig. 1.10).

If no more cutting is to be made, the cut-off wheel should be taken off and kept as prescribed in sect. 6.1.

IMPORTANT

Remove the protective screen after use in order to avoid that enclosed damp air causes corrosion.

Notice that the motor cannot rotate until the SPEED knob has been OFF, when the motor has been stopped by the stop screw or by the thermal fuse.

Manual cutting

For the coarse cutting of specimens (e.g. printed circuit boards, thin metal sheets, mineral samples) it would be advantageous to use the manual cutting table.

Mounting of manual cutting table (CATMA)

- Choose a diamond cut-off wheel with a diameter of at least 100 mm, e.g. 330CA, 352CA, 551CA. Use only the big flange set, Ø 65 mm. (Mount as described in sect. 6.3.2).
- Turn the cooling tubes against the cut-off wheel.
- Mount the 2 screws so that the screw heads are approx. 5 mm above the bottom position.
- Place the cutting table so that the cut-off wheel goes through the cut-out of the cutting table and the cut-outs fit around the screws.
- Push the table to the left so that it makes full contact with the side of the Accutom and tighten the screws.
- Mount the guide rail.

Loosen the screws for the guide rail approx. 5 mm, place them in the cut-outs of the cutting table and tighten.

IMPORTANT

The gear ratio between motor and spindle must be in the low range, see sect. 7.3. if necessary.
Use the protective goggles for protection against splashing of cooling agent and possible burrs from the cutting.

Cutting

Adjust the guide rail to obtain the requested cutting width.
Set the cut-off wheel to max. speed and move the specimen into the cut-off wheel under light pressure.

IMPORTANT

Do not turn the specimen during the cutting as this might cause damage to the cut-off wheel.

Do not press hard, only a light pressure is necessary. If the cut-off wheel does not cut satisfactorily it should be dressed (sect. 7.3), or it should be replaced by a coarser type, possibly 551CA.

7. Maintenance

General cleaning

Clean the painted surfaces with ordinary detergents and a soft cloth.

IMPORTANT

Never use alcohol or other solvents, as the long-range effects will damage the lacquer.

Cleaning of waste-tray

Ordinary water, possibly with detergent added, and a soft brush, is used for cleaning.

It is an advantage to clean the waste tray immediately after the cutting, as dried-up scrap may be difficult to remove.

IMPORTANT

If the Accutom is connected to the recirculating cooling unit the Latex tube, which is included in the package, is placed on the outlet branch and led to a suitable discharge. In this way pollution of the cooling agent and running over of the recirculating cooling unit are avoided.

Storing and maintenance of diamond cut-off wheels

The precision of diamond cut-off wheels and thus the cut depends on how carefully the following is observed:

- Never expose the cut-off wheel to overload, such as heavy mechanical load or heat.
- Store the cut-off wheel in a dry place, horizontally on a plane support, preferably under light pressure.
- A clean and dry cut-off wheel does not corrode. Therefore, clean and dry the cut-off wheel before storing. If possible, use ordinary detergents for the cleaning.

- Regular dressing of the cut-off wheel is also part of the general maintenance.
A newly dressed cut-off wheel will have optimum cutting qualities. A badly maintained and poorly dressed cut-off wheel requires a higher cutting pressure and will, therefore, generate more frictional heat. It may also bend which will result in an oblique cut.
A combination of both might result in damaging the cut-off wheel. Actually, a poorly dressed cut-off wheel is the most frequent reason for damage of the cut-off wheel.
The cut-off wheel is dressed as follows:
Let the cut-off wheel rotate on min. revolutions. Now dress the cut-off wheel by moving the moistened dressing stone 1 or 2 times all over the diamond layer. This can be done manually, but it is safer (especially at the thin cut-off wheels) to clamp the dressing stone in the specimen holder and cut it in the normal way.
After the dressing the cut-off wheel is cleaned by applying plenty of cooling water.

Storing of bakelite bonded Al_2O_3 cut-off wheels

These cut-off wheels are sensitive to humidity. Therefore, do not mix new and old (dry and humid) cut-off wheels. It is best to store them in a dry place, horizontally on a plane support.

Cutting in of fuse

The lamp (fig. 1.8) indicates whether there is power on the apparatus. A constant heavy load of the motor will cut-out the fuse on the front plate (fig. 1.7). It can be cut-in again by pressing the button.

**Maintenance of stands and
specimen holders**

Cutting stands

They should be kept clean of cutting scrap and be regularly lubricated in movable parts with 1 or 2 drops acid free oil. This is especially important for the O-rings in the connecting lever for the swivelling movement (fig. 7).

IMPORTANT

By use of stands with micrometer screw the latter should be protected against corrosion by lubricating with a thin coat of acid-free oil. Turn the micrometer screw quite back before lubrication in order to make the protection as effective as possible.

Every second month, the internal parts of the micrometer screw will have to be lubricated with 1 or 2 drops of acid-free oil, poured in one of the two holes at the end of the micrometer screw.

Specimen holders

Here all movable parts, dove tail feeds and screws should be kept clean and be lubricated at regular intervals. Use acid-free oil.

Change of cooling water in the recirculating cooling unit

The cooling liquid should be changed regularly in order to avoid growth of bacteria, which would disintegrate the liquid and thus make it change its qualities. Change of liquid would depend on such factors as the quality of the water, the ambient temperature and the frequency of use. However, we suggest that change of liquid at intervals of 3-4 weeks would be a good starting point.

The liquid is changed as follows:

- Disinfect the cooling liquid the day before by adding 0,2% corresponding to 24 ml bacterial disinfectant BACTE to the recirculating cooling unit.
Hereby the recirculating cooling unit is disinfected and new attacks are prevented.
- Throw out the cooling liquid.
The easiest way to do this is to let the pump almost empty the unit after removing the pump tube from the inlet valve of Accutom-2 and lead it to the outlet. Thus the unit is easier to carry.

IMPORTANT

Cooling liquids may contain agents which makes it illegal to pour into ordinary drains. Inquiries should be made with the local authorities.

- Turn the L-outlet (fig. 2.6) piece of Accutom-2, so that it points upwards.
- Remove the pump - use handle (fig. 2.11).
- Remove the front plate of the recirculating cooling unit and take the vessel out and clean it of sedimented cutting scrap.
- Re-mount the vessel and fill it with water, 11.5 l, then add additive and stir to spread the additive.
When using ADDUN apply from 2-6% = 230-690 ml. 2% is used if the cut-off wheel has to cut freely, typically when cutting metals.
6% is used if a lubricating effect is desired, which gives the cutting a low surface roughness.
- Replace the pump.

The pump unit must be full of liquid before starting. This will occur naturally if the pump is lowered slowly into place in the liquid in the recirculation unit.

It could still be necessary to remove any air which may be trapped in the pump unit using the enclosed 20 ml syringe.

- Fill the syringe slowly by drawing up the liquid through the pump's tubing.
- Then empty the syringe quickly, again through the pump's tubing.
- If necessary this operation should be repeated.

8. Service

Cleaning of the cutting table bearings

If the cutting table does not move regularly at low speeds it is necessary to clean the bearings.

- Remove the 4 internal hexagon headed screws, and dismount the aluminium plate, the water tray, the spacer tubes, the O-rings and the lower cover plate. Notice the sequence of mounting.
- Clean the bars with a lint-free cloth without solvents. Also remove any surrounding dirt so that the bearings do not quickly become dirty again.
Rub the bars with a thin layer of acid free oil.
- Check that the cutting table can now move at low speeds.
- Mount in reverse sequence.
- Align the cutting table in relation to the division plate.
- Fill up the holes in the hexagon headed screws with silicone matter.

Ventilating the hydraulic system

Irregular, or no movement at all of the cutting table may be due to air in the hydraulic damping circuit. Usually, this was caused during the transportation of the apparatus.

Procedure

- Take hold of the left side of the apparatus so that it forms an angle of 60° with the support.
- Turn off the FEED valve.
- Press the FORWARD button and move the cutting table back and forward 5 to 10 times in all its travel.
- Repeat the process, this time with the FEED valve opened completely, until the table moves satisfactorily.

IMPORTANT

Watch the pressure gauge. The pressure must not exceed max. indication.

Remember to pull the table slowly back in its starting position.

Fault in the hydraulic system

If the air in the hydraulic system is not due to the transportation, it may have the following causes:

- The cutting table has been moved forward. The FORWARD knob has not been used.
This may have changed the deflection of the pressure gauge, or resulted in leaking hose connections.
- The oil reservoir is empty.
This is due to a leakage which must be found before the reservoir is filled.
- The washers in the damping cylinder are worn. When changing them it should be noticed that the plunger shaft has floating suspension which should be maintained when remounting.

Changing of gear ratio

IMPORTANT

Disconnect Accutom-2 from the mains
by taking out the plug from the socket.

The belt pulleys of Accutom-2 have been designed so that it is possible to change the gear ratio between 1:1 and 3:1. At the gear ratio 3:1 both V-belts are used, but at the gear ratio 1:1 the O-belt is used (See also fig. 8).

Change of the gear ratio is done as follows:

- Switch off the Accutom from the mains and remove the motor screen.
- Force the belt over the edge of the groove in the direction in which the belt is to be moved, and at the same time turn the spindle round so that the belt moves.
In this way the belts are moved one by one.
Check the belt tightening (sect. 8.4).
- Mount the motor screen.

Belt tightening

If the belts can move more than approx. 5 mm inward when pressing in the middle between the belt pulleys, the belts must be tightened.

Loose belts will result in vibrations which transmit through the spindle to the cut-off wheel, which may result in a poorer cutting quality.

The belts are tightened as follows:

- Loosen the 4 screws on the black motor plate.
- Pull the motor back until the belts are correctly tightened, i.e. they can move approx. 5 mm in the middle between the belt pulleys.

IMPORTANT

Never tighten the belts harder, as this would result in a deflection of the spindle arrangement, causing the cut-off wheel to become oblique in relation to the cutting table.

- Check that the belt pulleys are parallel to each other and tighten the screws.
- Adjust the connection between motor and spindle so that the latter is tightened lightly.

IMPORTANT

If it is tightened too much it will, like when tightening the belts, result in a deflection of the spindle.

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
Some of the drawings may contain position numbers
not used in connection with this manual.

Drawing

Fig. 1

Pos. Description

9.1 Accutom-2, front


Drawing

Pos. Description

Fig. 2

9.2 Accutom-2, back


Drawing

Fig. 3
Drawing 14270164

Pos. Description

9.3 Inlet valve

Drawing

Pos. Description

Fig. 4-6

Fig. 4-6

Mounting of CATCH


Fig. 4


Fig. 6



Fig. 5

Drawing

Fig. 7
Drawing 314 M 385R

Pos. Description

9.7 Rocking lever


Drawing

Fig. 8
 Drawing 14270159


Pos.	Description	Spare Part No.
9.8	Belt pulley	

Belts, set of 2 (100-1000 rpm)	314 MP 064
Belt (300-3000 rpm)	427 MP 078

100 - 1000 R.P.M.

300 - 3000 R.P.M.


Accutom-2
Spare Parts and Illustrations

Drawing

Fig. 9
Drawing 14270201A

Pos.	Description	Spare Part No.
9.9	Cabinet	
1	Cabinet	14279023
2	Bracket for dial gauge.....	427 MP 005
3	Speed control PCB	314 MP 053
5	Bottom plate.....	427 MP 080
6	Intermediate wall.....	427 MP 007
7	Speed-knob	314 MP 007
8	Feed-knob.....	314 MP 007
9	Lamp, complete	314 MP 008
10	Stop-lamp, complete.....	314 MP 011
11	Fuse (please state voltage and serial number)	427 MP 086
12	Speed indicator.....	14279026
13	Pressure gauge 0-1.6 bar, with holder	14279027
14	Shield	427 MP 033
15	Transparent PVC shield.....	427 MP 034
16	Jet-pipe block housing	427 MP 035
17	Jet-pipe, left	427 MP 037
18	Jet-pipe, right	427 MP 038
19	Dial gauge.....	398 MP 007
	Outlet Latex tubing, 2 m.....	427 MP 041
	Set of Hexagon keys, 2, 2.5, 3 and 4 mm	427 MP 042
	Two-pin spanner for flange screw	427 MP 085
	Gear set 3000/1000 rpm.....	427 MP 001
	Shield, U-shaped	427 MP 002
	Block with cooling water jet tubes.....	427 MP 003
	Hydraulic oil, 0.5 l	314 MP 137
20	Nozzle for jet pipe	427 MP 089
	Set of sealings for jet-pipe	427 MP 039
	Straight fitting for cooling water	427 MP 040

(Fig. 9, Cabinet, Drawing 14270201A)


Drawing

Fig. 10
Drawing 14270210

Pos.	Description	Spare Part No.
9.10 Bracket		
1	Bracket for front and back plate.....	14279029
	Bracket for intermediate wall	427 MP 009
2	Legs, 4 pcs, with 4 screws M4 x 25 mm.....	314 MP 002
3	Spindle housing	314 MP 074
4	Spindle	427 MP 010
5	Ball bearings SKF 6004-2RS, 2 pcs.	314 MP 071
6	Spring ring K6004, 4 pcs.	314 MP 072
7	Seeger ring 142 x 1.75 DIN 472	427 MP 011
	Spindle with housing, complete	427 MP 062
	Inner flange, Ø 42 mm	314 MP 075
	Outer flange, Ø 42 mm	314 MP 076
	Inner flange, Ø 65 mm	314 MP 138
	Outer flange, Ø 65 mm	314 MP 139
	Flange screw	427 MP 012
8	Spindle pulley, Ø 60 mm.....	427 MP 013
10	Terminal strip 3/12 6ESDF	427 MP 083
11	Aligning rod, threaded.....	314 MP 066
12	Smoothing filter.....	427 MP 082


(Fig. 10, Bracket, Drawing 14270210)

Drawing

Fig. 11
Drawing 14270203

Pos.	Description	Spare Part No.
9.11	Valve plate	
1	Front plate.....	427 MP 079
4	Feed valve, complete.....	314 MP 057
3	Oil reservoir, 30 ml.....	314 MP 056
2	By-pass valve, complete.....	314 MP 012
5	PCB for instruments, 110 V or 220 V 50-60 Hz (Please state voltage).....	314 MP 054

Accutom-2
Spare Parts and Illustrations


Drawing

Fig. 12
Drawing 14270211A

Pos. Description Spare Part No.

9.12 Accutom-2 back, details		
Pos.	Description	Spare Part No.
1	Back plate	14279028
2	Copper tube for cooling water	427 MP 018
3	Plastic cap	427 MP 019
4	Ball valve 1/4"	383 MP 020
5	Angle fitting for valve	427 MP 020
6	Fuse holder	381 MP 220
7	Cable entry PG9	427 MP 021
8	Female plug for pump	427 MP 022
9	Terminal strip 20/2	427 MP 023
	Fuse 1 AT	427 MP 024
10	Noise filter RD 52-2	314 MP 058
11	2 screws and 4 gaskets for copper tube	427 MP 087



(Fig. 12, Accutom-2 back, details, Drawing 14270211A)

Drawing

Fig. 13
Drawing 14270205

Pos.	Description	Spare Part No.
9.13 Motor		
1	Motor 220 V, 110 V (Please state voltage and frequency)	314 MP 059
2	Pulley for motor	427 MP 025
3	Plate for motor	427 MP 026
4	Angle piece for motor	314 MP 065
5	Noise filter, complete	314 MP 060


Drawing

Fig. 14
Drawing 14270206

Pos.	Description	Spare Part No.
9.14 Table		
1	Bearing housing.....	427 MP 027
2	Table bearing support.....	314 MP 036
3	Table bearing A	314 MP 037
4	Table bearing B	427 MP 028
5	Shaft for table bearing	314 MP 040
6	Shaft support frame	427 MP 029
7	Back support for table.....	314 MP 032
8	Front bridge for table	314 MP 031
9	Handle for table	314 MP 030
10	Stop adjustment screw	314 MP 145
11	Stop pin for dial gauge.....	314 MP 034
12	Cylinder fixture.....	314 MP 043
13	Spring-set	314 MP 041
14	Spring fixture.....	314 MP 042
15	Table plate	427 MP 030
16	Distance ring for table, 4 pcs.....	314 MP 038
17	O-ring 8 x 3 mm, 4 pcs.....	314 MP 039
18	Lower PVC-cover.....	314 MP 021
19	PVC-tray	427 MP 031
20	Pin for splash-guard.....	314 MP 019
	Angular outlet tube.....	427 MP 032
21	Non-return valve, complete	314 MP 143
22	Table stop-switch.....	314 MP 144
	Cylinder, complete	314 MP 142
23	Cylinder housing, with 2 fittings	314 MP 047
24	Piston rod, complete	314 MP 141
25	Cylinder bottom piece, complete	314 MP 048
26	Cylinder top, complete	314 MP 049
<i>Recirculating cooling unit</i>		
	Chassis	14279032
	Front plate.....	14279035
	Reservoir	427 MP 045
	Tray for pump	427 MP 046
	Lid with knob.....	14279034
	Pump 110 V or 220 V, 50-60 Hz (Please state voltage).....	427 MP 048
	Male plug for pump	427 MP 049

(Fig. 14, Table, Drawing 14270206)

Drawing

Fig. 15
Drawing 14270207


Pos.

Description

Spare Part No.

9.15 Manual cuttings table


1	Base plate	427 MP 050
2	Splash guard	427 MP 051
3	Tightening piece, narrow	427 MP 052
4	Tightening piece, wide	427 MP 053
5	Distance piece for screw	427 MP 054
6	Handle for splash guard	427 MP 055
7	Block for guiding rail	427 MP 056
8	Guiding rail	427 MP 057
9	Tightening screw for guiding rail	427 MP 058 </td
10	Guiding shaft	427 MP 059
11	Screws for guiding shaft, 2 pcs	427 MP 060

Accutom-2
Spare Parts and Illustrations

Drawing	Pos.	Description	Spare Part No.
Fig. 16 Drawing 14270208A	9.16	Stand, with motor	
CATAS stand	1	Housing.....	314 MP 099
	2	Lid	427 MP 061
	3	Bottom plate.....	314 MP 082
	4	Dove tail clamp, complete	314 MP 092
	5	Fixture for micrometer screw	427 MP 063
	6	Lock for micrometer screw	427 MP 064
	7	Brass ring for stand.....	427 MP 065
	8	Micrometer screw, complete.....	427 MP 066
	9	Tightening piece for tube	427 MP 067
	10	Knob for 427 MP 067.....	427 MP 068
	11	Motor pulley, small.....	314 MP 095
	12	Micrometer pulley, large	314 MP 096
	13	Motor 220 V or 110 V (Please state voltage).....	314 MP 094
	14	Switch for motor, with cap.....	314 MP 150
	15	Cable entry	427 MP 069
	16	Crank with screws.....	314 MP 097
	17	Bearings with springs, set of 2.....	314 MP 090
	18	Bearing housing, set of 2.....	314 MP 083
		Eccentric screw.....	427 MP 070
	19	Threaded rod for housing	314 MP 151
		Hexagon key 4 x 200 mm	314 MP 088
	20	Sealing for lid	427 MP 075
	21	Sealing for knob.....	427 MP 076
	22	Sealing for micrometer screw	427 MP 077
	23	Ball belt	314 MP 098
CATYS stand	1	Housing.....	314 MP 081
	2	Lid	314 MP 080
	3	Bottom	314 MP 082
	4	Dove tail clamp, complete	314 MP 140
	5	Locking device	314 MP 079
	6	Micrometer screw	314 MP 086
	7	Coupling, complete	314 MP 087
	8	Supporting tube	314 MP 084
		Eccentric screw.....	427 MP 070
		Hexagon key 4 x 200 mm	314 MP 088
	9	Sealing for lid	427 MP 075
	10	Sealing for knob.....	427 MP 076
	12	Knob.....	427 MP 068
	11	Sealing for micrometer screw	427 MP 077
CATIN sample holder		Set of screws	427 MP 072
CATCH Sample holder		Sample holder CATCH	04276908
		Set of screws	427 MP 073
CATAP sample holder		Frame for vacuum chuck	314 MP 126
		Base part with valve and dovetailed plate	314 MP 127

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







Herewith declares that

<i>Product Name:</i>	Accutom-2
<i>Type No:</i>	427
<i>Machine Type:</i>	Precision cut-off machine

is in conformity with the provisions of the following directives:


Safety of Machinery 89/392/EEC and 91/368/EEC and 93/44/EEC according to the following standard(s):
EN292-1:1991, EN292-2:1991, EN60204-1:1992.

EMC-Directive 89/336/EEC and 92/31/EEC according to the following standard(s):
EN50081-1:1991, EN50082-1:1992.

Low Voltage Directive 73/23/EEC and 93/68/EEC according to the following standard(s):
EN60204-1:1992.

Supplementary Information The equipment complies with the American standards:
FCC part 15, subpart B and UL508.

The above has been declared according to the global method, module A

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Klaus Kisbøll, Entwicklungsleiter, Struers A/S
Klaus Kisbøll, Chef du service développement, Struers A/S