Accutom-50

Instruction Manual

Manual No.: 15017001

Date of Release 15.11.2012

Accutom-50
Instruction Manual

Table of Contents	Page
User's Guide	1
Reference Guide.....	35
Quick Reference Guide	60

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

Original instructions. The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2012.

Struers A/S
Pederstrupvej 84
DK 2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Accutom-50 Safety Precaution Sheet

To be read carefully before use

1. The operator should be fully instructed in the use of the machine, its cut-off wheels and cup wheels according to the Instruction Manual and the instructions for the cut-off wheels and cup wheels.
2. The machine must be placed on a safe and stable support table.
3. Be sure that the actual voltage corresponds to the voltage stated on the back of the machine. The machine must be earthed.
4. Use only intact cut-off wheels or cup wheels. The cut-off wheels and cup wheels must be approved for min. 5000 rpm. If other cut-off wheels, saw blades or cup wheels are used, make sure that the speed setting of Accutom-50 complies with the max. speed for the cut-off wheels, saw blades or cup wheels.
5. Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
6. The sample must be securely fixed in the specimen holder.
7. Do not touch the sample, the specimen holder head, the cut-off wheel or cup wheel while positioning the sample with the POSITION controls.
8. Never try to open the cover before the cut-off wheel or cup wheel has stopped completely.

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

Disposal

Equipment marked with a WEEE symbol  contain electrical and electronic components and must not be disposed of as general waste.

Please contact your local authorities for information on the correct method of disposal in accordance with national legislation.

User's Guide

Table of Contents	Page
1. Getting Started 3	
Checking the Contents of Packing	3
Placing Accutom-50	3
Getting Acquainted with Accutom-50	3
Supplying Power	4
Changing the voltage setting	4
Recirculation Unit	6
Software Settings	7
Configuration Menu.....	7
Setting the Language.....	8
2. Basic Operations	
Using the Controls.....	9
Front Panel Controls of Accutom-50.....	9
Groups of Keys	9
Acoustic Signals.....	9
Location of Main Switch	9
Front Panel Controls	10
Display	11
Reading the Display	12
Changing/Editing Values.....	13
Numeric Values.....	13
Alphanumeric Values	14
Positioning the Sample	15
Reference Position.....	15
Absolute Position	15
Relative Position	15
Relative Zero.....	15
Stop Position.....	15
Selecting between Cutting and Grinding.....	16
Cutting Display	17
Grinding Display.....	17
Cutting.....	18
Changing the Cut-off Wheel.....	18
Clamping the Sample and Specimen Holder	19
Positioning the Sample	19
Setting the Cutting Parameters	20
Wheel.....	20
Speed.....	20
Feed.....	20
Force.....	20

Cutting Mode	21
Single Cut	21
Multiple Cut	21
Rotation	22
Water	22
Starting the Cutting	22
During Cutting	23
Changing the Feed Speed	23
Retracting the sample	23
Force Limit Reached	23
Stopping the Cutting	24
Automatic Stop	24
Manual Stop	24
Grinding	25
Changing the Cup Wheel	25
Connecting the Vacuum Chuck to the Vacuum System	25
Clamping the Sample and Specimen Holder	26
Positioning the Sample	27
Positioning the Sample	27
Setting the Grinding Parameters	28
Cup Wheel	28
Speed	28
Y-Feed  	28
Sweeps  	28
X-Inc. 	28
Water	28
Rotation	29
Removal Mode	29
Starting the Grinding	31
During Grinding	31
Stopping the Grinding	31
Automatic Stop	31
Manual Stop	31
Continuing the Grinding	31

3. Maintenance

Daily Service	32
Checking the Recirculation Unit	32
Weekly Service	33
Refilling the Cooling Water Tank	33
Emptying and Cleaning the Tank	33
Refilling the Tank	33
Maintenance of Cut-off Wheels	34
Maintenance of Diamond and CBN Cut-off Wheels	34
Storing of Abrasive Cut-off Wheels	34
Maintenance of Cup Wheels	34

1. Getting Started

Checking the Contents of Packing

In the packing box you should find the following parts:


- 1 Accutom-50
- 2 Mains cables
- 1 Specimen Holder with parallel vice
- 1 Flange for cut-off wheels (42 mm dia.)
- 1 Stop pin
- 1 Spanner, 17mm
- 1 Small grate
- 1 Large grate
- 1 Debris collecting tray
- 1 Tube for vacuum connection
- 1 Cooling tube for use with cup-wheels
- 1 Allen key, 2 mm
- 1 Allen key, 2.5 mm
- 1 Allen key, 3 mm
- 1 Allen key, 4 mm
- 1 Allen key 5 mm
- 2 Screws M4x20
- 2 Screws M4x35
- 1 Set of Instruction Manuals

Placing Accutom-50

Accutom-50 should be placed on a stable and plane (tolerance ± 1 mm) table. The table must be able to carry a weight of min. 50 kg.

Getting Acquainted with Accutom-50

Take a moment to familiarise yourself with the location and names of the Accutom-50 components.

- A Front panel/Front panel control(s)
- B Main switch
- C Recirculation Cooling Unit
- D Cut-off wheel/Cup wheel
- E Specimen holder head
- F Screw to vacuum connection

Supplying Power

Always remember to switch the power off when installing electrical equipment!

IMPORTANT
Check that the mains voltage corresponds to the voltage stated on the type plate on the back of the machine.

Changing the voltage setting

The factory setting for Accutom-50 is 240V.

If the factory setting is not the correct setting for your mains supply the setting can be changed from 240V to 220V:

- Pull out the fuse holder from the cable terminal at the back of the machine.
- Pull out the voltage switch and turn it to the correct setting.


Voltage Required	Setting
230 or 240V	240V
200 to 220V	200V

Note: The two additional settings, 110V and 120V are not to be used.

- Reinsert the voltage switch.
- Put the fuse holder back into the cable terminal.

Accutom-50 is shipped with 2 types of Mains cables:

Single-phase Supply

The 2-pin (European Schuko) plug is for use on single-phase connections.

If the plug supplied on this cable is not approved in your country, then the plug must be replaced with an approved plug. The leads must be connected as follows:

Yellow/green: earth
Brown: line (live)
Blue: neutral

3-phase Supply

The 3-pin (North American NEMA) plug is for use on 3-phase power connections.

If the plug supplied on this cable is not approved in your country, then the plug must be replaced with an approved plug. The leads must be connected as follows:

Green: earth
Black: line (live)
White: line (live)

Connection to the Machine

Both cables are on the other end equipped with an IEC 320 cable connector that has to be connected to the Accutom.

WARNING!

The output voltage from this cable is 200 – 240V and not 110V.
DO NOT use this cable to connect equipment that use a 110V power supply. Failure to adhere to this may result in material damage.

Recirculation Unit

- Pull out the recirculation tank.
- Fill the tank with 3.88 l of water and 120 ml Struers Additive. The water level should be 5 mm below the edge of the front hole in the tank cover.

IMPORTANT

Always ensure that there is sufficient water in the tank as the recirculation pump will be damaged if it is running dry.

- Check that the cover is fitted securely into place in the recirculation tank and push the drawer with the recirculation tank into place again.
- Check that the end of the inlet tube has dropped into position.

Note:

Changing of cooling water should be done at least once a month.

IMPORTANT

Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to add Struers Additive each time you refill with water.
Do not use any oil, petrol, or turpentine-based additives, only Struers additive.




Software Settings
Configuration Menu

Switch on the power at the main switch located at the back of the machine. The following display will appear briefly:

Afterwards the display will change to the same screen which was shown before Accutom-50 was switched off, usually a cutting or grinding method. When switching Accutom-50 on for the first time, the display to appear should be the MAIN MENU. If the heading in the display is different, press Esc, until the MAIN MENU appears. (A long beep can be heard)

The MAIN MENU is the highest level in the menu structure. From here you can go to configuration, manual functions and cutting or grinding methods.


Setting the Language

 Press MENU UP/DOWN  to select OPTIONS.

↓

 Press ENTER  to activate the OPTIONS menu.

↓


 Press ENTER  to activate the LANGUAGE Menu.










↓

 Press MENU UP/DOWN  to select the Language you prefer.


↓

 Press ENTER  to accept the language. The OPTIONS Menu now appears in the language you have chosen.

Front Panel Controls

Name	Key	Function	Name	Key	Function
① MAIN SWITCH		The main switch is located on the back of the machine.	⑦ ESC		Leaves the present menu or specimen holder position screen.
② START		Starts the cutting or grinding process according to the pre-set method. The recirculation water, if selected, is turned on.	⑧ FAST POSITION LEFT		Changes to POSITION menu or moves the specimen holder to the left in the X-direction in steps of 100 µm. Keep the key pressed to increase the speed.
③ STOP		Stops the cutting or grinding process. The recirculation water, if selected, is turned off.	⑨ FAST POSITION RIGHT		Changes to POSITION menu or moves the specimen holder to the right in the X-direction in steps of 100 µm. Keep the key pressed to increase the speed.
④ FUNCTION KEYS		Controls for various purposes. See the bottom of the individual screens.	⑩ POSITION LEFT		Changes to POSITION menu or moves the specimen holder slowly to the left in the X-direction in steps of 5 µm. Keep the key pressed to increase the speed.
⑤ MENU		Scrolls up (▲) or down (▼) in the menu tree structure of Accutom-50. When setting a parameter the value is increased (▲) or decreased (▼).	⑪ POSITION RIGHT		Changes to POSITION menu or moves the specimen holder slowly to the right in the X-direction in steps of 5 µm. Keep the key pressed to increase the speed.
⑥ ENTER		Selects a marked parameter value or chooses a menu.	⑫ POSITION UP/DOWN		Changes to POSITION menu or moves the specimen holder up- or downwards in the Y-direction in steps of 100 µm. Keep the key pressed to increase the speed.

Display

- A Heading
- B Inverted text. Cursor
- C Function key options
- D Arrow indicates, that there are more lines in the picture

Please Note

The examples of display screens in this Instruction Manual show a number of possible texts. The actual display screen may differ from the examples in the Instruction Manual.

Reading the Display

The display can show various kinds of information, for example about the cutting or grinding method or about the sample position. A screen for a cutting or grinding method could look as the following example:

```

Cutting Method:  1. EMPTY METHOD

Wheel   :      B0D15   Speed  : 2700 rpm
Feed    : 0.100 mm/s   Force   :  MEDIUM
Rotation:      OFF
Water   :      ON

F1:MULTI. CUT
    
```

Pressing one of the X-axis positioning keys will change the screen to the following:

```

X-POSITIONS  ⬆⬇⬆

Absolute Position:      15.255 mm
Relative Position:      5.000 mm

F1:RESET  F2:ADD  F3:MOVE TO REL. ZERO
    
```

F2:ADD is only available in cutting methods.

Pressing one of the Y-axis positioning keys will change the screen to the following:

```

Y-POSITIONS  ⬆⬇⬆

Absolute Position:      55.7 mm
Relative Position:      0.0 mm
Stop Position:          20.0 mm

F1:RESET  F2:SET STOP  F3:MOVE TO REL.ZERO
    
```

Press Esc, to return to the cutting or grinding method screen.

The “heavy” arrows on the display indicate the direction in which the sample holder is moved, for example:

X-POSITIONS ⬆⬇⬆
Y-POSITIONS ⬆⬇⬆


Changing/Editing Values

Depending on the type of value, there are two different ways of editing.

Numeric Values

```
Cutting Method:  1. EMPTY METHOD
Wheel   :      B0D15   Speed   : 2700 rpm
Feed    : 0.100 mm/s   Force    : MEDIUM
Rotation:      OFF
Water   :      ON
F1:MULTI. CUT
```

▲ Press MENU UP/DOWN ▲▼ to select the numeric value you want to change, e.g. Feed:

```
Cutting Method:  1. EMPTY METHOD
Wheel   :      B0D15   Speed   : 2700 rpm
Feed    : 0.100 mm/s   Force    : MEDIUM
Rotation:      OFF
Water   :      ON
F1:MULTI. CUT
```

 Press ENTER ↵, to edit the value.


Two square brackets [] appear around the value.

```
Cutting Method:  1. EMPTY METHOD
Wheel   :      B0D15   Speed   : 2700 rpm
Feed    [0.100]mm/s   Force    : MEDIUM
Rotation:      OFF
Water   :      ON
F1:MULTI. CUT
```

▲ Press MENU UP/DOWN ▲▼ to increase or decrease the numeric value.

 Press ENTER ↵, to accept the new value.

or

 Press Esc, to keep the original value.

Alphanumeric Values

```
Cutting Method: 1. EMPTY METHOD
Wheel : B0D15 Speed : 2700 rpm
Feed : 0.100 mm/s Force : MEDIUM
Rotation: OFF
Water : ON
F1:MULTI. CUT
```

Press MENU UP/DOWN to select the alphanumeric value you want to change, e.g. Cut-off wheel:

Press ENTER ↵, to edit the value.

A submenu appears.

```
CUT-OFF WHEELS +
M0D08 M1D08 M0D10 M1D10 M0D13
M1D13 B0D13 B0C13 50A13 30A13
SAW13 M0D15 M1D15 E0D15 B0D15
B0C15 50A15 40A15 10S15 30A15
USER 1 USER 2 USER 3 USER 4 USER 5
USER 6 USER 7 USER 8 USER 9 USER10
```

Press MENU UP/DOWN to select the correct cut-off wheel.

Press ENTER ↵, to accept the new value and to return to the previous screen.

or

Press Esc, to keep the original value and to return to the previous screen.

Positioning the Sample

After clamping the sample in the specimen holder, the holder is placed in the specimen holder head and fixed there. To move the specimen holder and the sample, the POSITION keys ▼▲▶▶◀◀ are used. The display shows the position in either X- or Y- direction.

Reference Position

Every time the power is switched on while the cover is closed, Accutom-50 checks its own reference position. The sample holder head will be moved back as far as possible, to the reference position (X=0.000 Y=0.0), and thereafter it will return to the position where it was before the power was switched on.

Absolute Position

The absolute position shows the total distance the specimen holder has travelled from the reference position.

Relative Position

The relative position equals that of the absolute position until it is set to zero at a desired point. By setting it to zero, calculation of the sample movement close to the cut-off wheel or cup wheel is made easier. The screen value always relates to the distance the specimen holder has travelled since being set to zero.

Relative Zero

The relative zero position is the point where the relative position in either the X- or Y- position was set to zero. Having completed the cutting process, the sample holder automatically returns to this point. Pressing F3 when in the X-or Y- position screen also returns the sample to the relative zero position.

Stop Position

A stop position can be set to halt the cutting process at a precise point. After reaching this point the sample will be retracted and returned to the relative zero position.

Compensating for Wheel Wear


Please make sure when setting the stop position to compensate for possible wear of the cut-off wheel. This is especially important when using Al₂O₃ or SiC wheels.

Selecting between Cutting and Grinding

To select between cutting and grinding you must be in the Main Menu. As Accutom-50 always starts up in the same state as when the voltage was switched of last time, it might be necessary to go back to the Main Menu.

To go back to the Main Menu:

Esc Press Esc, until the MAIN MENU display appears.
(A long beep can be heard).

 Press MENU UP/DOWN   to select between CUTTING METHODS and GRINDING METHODS.

 Press ENTER , to accept the selection.

If you selected CUTTING METHODS, the screen will appear as follows.

Cutting Display

SELECT CUTTING METHOD ↓			
1. EMPTY METHOD	2. EMPTY METHOD		
3. EMPTY METHOD	4. EMPTY METHOD		
5. EMPTY METHOD	6. EMPTY METHOD		
7. EMPTY METHOD	8. EMPTY METHOD		
9. EMPTY METHOD	10. EMPTY METHOD		
11. EMPTY METHOD	12. EMPTY METHOD		
F1: COPY	F2: INSERT	F3: RESET	F4: RENAME


↓ If you selected GRINDING METHODS, the screen will look as follows.

Grinding Display

SELECT GRINDING METHOD			
1. EMPTY METHOD	2. EMPTY METHOD		
3. EMPTY METHOD	4. EMPTY METHOD		
5. EMPTY METHOD	6. EMPTY METHOD		
7. EMPTY METHOD	8. EMPTY METHOD		
9. EMPTY METHOD	10. EMPTY METHOD		
F1: COPY	F2: INSERT	F3: RESET	F4: RENAME

Cutting

Changing the Cut-off Wheel

- Open the cover and swing the right cooling tube in the air.
- Insert the stop pin in the hole of the inner flange.
- Use the spanner (17 mm) to loosen the flange screw.
- Remove the outer flange and the wheel.

IMPORTANT

The tolerance between the spindle and inner flange is very small which means that the two surfaces must be absolutely clean. Never try to squeeze the cut-off wheel on as this may damage the spindle. If there are any small burrs remove them with grinding paper (grit size 1200).

- If a cup wheel has been used before remember to exchange the cooling tube for use with the cup-wheel with the cooling tubes for cutting. This is done by removing the end screw and sliding the appropriate parts off.
- Mount the new cut-off wheel and remount the outer flange, with the machined face towards the inner flange.
- Insert the locking pin in the hole in the inner flange.
- Gently fasten the flange screw using the spanner.
- Swing the right cooling tube back in its place.

Clamping the Sample and Specimen Holder


- Clamp the sample securely in the specimen holder using the appropriate Allen key.
- Fasten the specimen holder in the specimen holder head by pushing the specimen holder into the dovetail fixture and tightening the screw with an Allen key.
- When cutting with rotation or oscillation, the sample and the specimen holder should be clamped so that they rotate evenly around the centre of the sample. This way the fastest cutting is obtained as the cut-off wheel will be cutting most of the time and the possibility of damaging the cut-off wheel is limited

Positioning the Sample

- Move the sample into the correct start position, close to the cut-off wheel by using the POSITION keys.
- Reset the relative position in both X- and Y-direction by pressing F1 when in the respective screens.
- Set the stop position to define the length of the cut:

There are two ways of setting the stop position:

1. If you know the size of your sample:
 - Press ENTER ↵ and use MENU UP/DOWN ⏮⏭⏮⏭ to set the stop value.
 - Press ENTER ↵ again to accept the value or Esc to cancel.
2. If you do not know the exact size of your sample:
 - Position the sample using the POSITION key ► so you can move it along the cut-off wheel.
 - Move the sample along the cut-off wheel to the required stop position using the POSITION key ▲.
 - Press F2:SET STOP to record the stop position.
 - Press F3:MOVE TO REL.ZERO to move the sample back into the initial position.
 - Position the sample correctly in front of the cut-off wheel using the POSITION key ◀. If possible use F3:MOVE TO REL.ZERO.

Setting the Cutting Parameters All cutting parameters can be edited independently of each other.

See section on Changing/Editing values for details on how to change values.

Method Screen

Cutting Method: 1. EMPTY METHOD					
Wheel	: B0D15	Speed	: 2700 rpm		
Feed	: 0.100 mm/s	Force	: MEDIUM		
Slices	: 2	Thickness	: 1.000 mm		
Drift Correction Value (δ)		:	0.000 mm		
Rotation	: CONT.	Speed	: 1	Dir.:	> <
Water	:	ON			
F1: SINGLE CUT		F3: TEST ROTATION		F4: SAVE	

Wheel All Struers cut-off wheels are listed with their default values for force limit, wheel thickness and speed. Also eight undefined wheels are included for your own choice of cut-off wheels. (See also configuration of cut-off wheels).

Speed The speed of the cut-off wheel can be set between 300 and 5000 rpm in steps of 100 rpm. The default values for all Struers cut-off wheels are already saved together with the wheel definitions. Whenever a different cut-off wheel is selected the recommended speed will automatically be inserted.

Feed The feed speed can be set between 0.005 mm/s and 3.000 mm/s. (See the section on Cutting Parameters in chapter 3. Consumables for recommended settings).

Force The force limit is a protection for the cut-off wheel and/or the samples to be cut. Depending on the thickness or strength of the cut-off wheels the force limit can be set to three different levels, LOW, MEDIUM and HIGH. The default values for all Struers cut-off wheels are already saved together with the wheel definitions. Whenever a different cut-off wheel is selected, the recommended force limit will automatically be inserted.

Cutting Mode	Accutom-50 can be programmed to automatically cut several slices from the same sample. Pressing F1 toggles between SINGLE CUT and MULTIPLE CUT.
<i>Single Cut</i>	The normal cutting operation, where only one single cut is carried out.
<i>Multiple Cut</i>	A series of cuts is carried out.
Slices	If Multiple Cut has been selected, the desired number of slices to be cut from the sample can be set.
Thickness	<p>The thickness of the slices to be cut is defined.</p> <p>After start is pressed, Accutom-50 will perform one straight cut to have a reference surface. Then the sample is retracted and moved to the left in the X-direction, according to the set thickness of the slice. The thickness of the cut-off wheel will automatically be compensated for. Then the next cut will be performed and the above will be repeated until the defined number of slices has been obtained. <i>Note!</i> The defined number of slices includes the reference slice.</p> <p>If the programmed number and thickness of slices exceeds the possible number of slices, Accutom-50 will automatically tell you how many slices can be obtained.</p>
Drift Correction Value (δ):	<p>The nominal thickness values for all Struers cut-off wheels are already saved within the wheel definitions. Whenever a cut-off wheel is selected, that particular wheel-thickness will automatically be used. For user defined wheels though, the thickness has to be entered manually when configuring the wheels in the configuration menu.</p> <p>Accutom-50 automatically compensates for the thickness of the cut-off wheel when multiple cut is used. However, due to differences in feed speed and wheel speed between different methods where the same wheel is used, additional compensation might be necessary:</p> <ul style="list-style-type: none">■ After cutting a few test slices, the obtained thickness is measured, and the deviation from the pre-set thickness has to be entered as drift correction value, delta(δ).


Rotation

Sample rotation can be an advantage when cutting large, very hard, coated or very long samples. On Accutom-50, three different settings are possible.

Rotation: Off


The sample does not rotate.

Rotation: Continuous

The sample is rotating around its centre.
The speed can be set to three different levels, 1, 2 or 3.
The direction of rotation can be set to either counter- or co-rotation.
Counter-rotation is recommended.
Press F3 to test sample rotation with the set parameters.

Rotation: Oscillating

The sample is rocking/oscillating around its centre.
The speed can be set to three different levels, 1, 2 or 3.
The angle can be adjusted from 10° to 400°.
Press F3 to test sample rotation with the set parameters.

Note

When the sample is moved in the x- or y- direction, with Rotation set to Cont. or Osc. a warning message, "Rotation Mode!", will briefly flash on the screen as a reminder that rotation mode is selected.

Water

The cooling water can be set to either on or off. For all normal cutting operations the setting should be on.

Starting the Cutting

Before you start cutting make sure that the cutting chamber is clean to ensure a free flow of cooling water. The debris collecting tray must be free of removed material and other fragments.

IMPORTANT

Especially when using Multiple Cut, the debris collecting tray must be cleaned before Accutom-50 is started to avoid a build-up of removed material. This might result in a reduction of the cooling water flow, which might damage the sample or the cut-off wheel.

- Position the sample correctly.
- Set the correct cutting parameters.
- Make sure that the correct cut-off wheel is mounted.
- Close the cover of the machine.
- Press START .

During Cutting

After Start has been pressed, the display changes to the following:

Feed Speed		CUTTING STARTED	
			
0%	100%		
Pre-set Feed Speed	:	0.100	mm/s
Actual Feed Speed	:	0.100	mm/s
Remaining Process Time:		3	min

Both the pre-set and the actual feed speed are displayed. The horizontal bar is used to display the feed speed graphically.

Changing the Feed Speed

During the cutting process the feed speed can be changed. Simply press Enter ↵ and change the feed speed to the correct value. Press Enter ↵ again to confirm the change or Esc to cancel.

Retracting the sample

The only possibility of moving the sample after the cutting process has been started is to retract the sample from the cut-off wheel. This can be done by pressing the POSITION DOWN ▼ key. The forward movement of the sample is stopped and the sample moved backwards until the key is released again. Then the sample will start again to move forward with the pre-set feed speed. If the pre-set feed speed can not be achieved because the force limit is reached, Accutom-50 automatically reduces the feed speed to the maximum possible speed. This value is displayed, and after the cut is finished, the following message is shown on the display:

Force Limit Reached

CUTTING FINISHED	
Force limit was reached during cutting feed speed reduced to:	
0.070 mm/sec	

For similar samples to be cut afterwards, the feed speed should be reduced to the new value or below.

Stopping the Cutting



Automatic Stop

- Accutom-50 automatically stops the cutting process at the pre-set stop position. (Make sure to compensate for possible wheel wear when setting the stop position).
- The sample is then retracted and the cut-off wheel is stopped.

Manual Stop

- The cutting process can be stopped at any time during operation by pressing the STOP Ⓢ key. The sample remains in its actual position and the process stops there.
- To avoid stopping the process while the cut-off wheel still is in the sample, press the POSITION key ▼ for the Y-direction and then press F2:SET STOP. The sample is retracted immediately and the process will be stopped when the relative zero position is reached.
- The display changes back to what it was before cutting was started.

Grinding Changing the Cup Wheel

- Open the cover and swing the right cooling tube in the air.
- Insert the stop pin in the hole of the inner flange.
- Use the spanner (17 mm) to loosen the flange screw.
- Remove the outer flange and the cup wheel.

IMPORTANT

The tolerance between the spindle and inner flange is very small which means that the two surfaces must be absolutely clean. Never try to squeeze the cup wheel on as this may damage the spindle. If there are any small burrs remove them with grinding paper (grain size 1200).

- If a cut-off wheel has been used before, remove also the inner flange. Please remember also to remove the left cooling tube and to put on the cooling tube for use with the cup wheel.
- Mount the inner flange if necessary, (only applicable if a cutting operation has been performed previously), the new cup wheel and the outer flange.
Note: When mounting the inner flange, the notch on the flange must meet with the guiding pin on the spindle.
- Insert the locking pin in the hole in the inner flange.
- Fasten the flange screw gently using the spanner.
- Swing the cooling tube back into its proper position.

Connecting the Vacuum Chuck to the Vacuum System

Before the vacuum chuck is used for the first time, it has to be connected to a vacuum pump. Accutom-50 is prepared for the connection through a lead-in bushing into the cutting chamber.

- Remove the long vacuum tube from the vacuum chuck and replace it with the short tube delivered together with Accutom-50.
- Remove the screw on the left hand side inside the cutting chamber (see Getting Acquainted with Accutom-50) and connect the other end of the short tube to this outlet.
- Remove the screw on the left hand side of the back of Accutom-50 and connect the long tube removed from the vacuum chuck to this inlet.
- Connect the other end of the long tube to a vacuum pump.


*Clamping the Sample and
Specimen Holder*

- Clamp the sample securely in the specimen holder using the appropriate Allen key or place the sample mounted on a glass plate on the vacuum holder.
- Fasten the specimen holder in the specimen holder head by pushing the specimen holder into the dovetail fixture and tighten the screw with an Allen key.
- When grinding with rotation or oscillation, the sample and the specimen holder should be clamped so that they rotate evenly around the centre of the sample. This way the fastest grinding is obtained as the cup wheel will be grinding most of the time and the possibility of damaging the cup wheel is limited.

IMPORTANT

Do not use rotation when working with the vacuum holder. The vacuum tube will be wrapped around the holder. Use oscillation instead.

Positioning the Sample

- Move the sample into the correct start position, close to the cup wheel by using the position keys, as follows:
 - Move the sample in the X-direction until it slightly touches the cup wheel grinding surface. *NB:* The wheel should still be able to be freely turned manually.
 - Move the sample back in the Y-direction until it is just clear of the cup wheel grinding surface.
- Reset the relative position in both X- and Y-direction by pressing F1 when in the respective screens.
- Set the stop position to define when the sample has passed all the way over the face of the cup wheel and the back movement should start.

There are two ways of setting the stop position.

1. If you do know the size of your sample:
 - Press ENTER ↵ and use MENU UP/DOWN ⏮▲⏭▼ to set the stop value.
 - Press ENTER ↵ again to accept the value or Esc to cancel.
2. If you do not know the exact size of your sample:
 - Position the sample using the POSITION key ► so you can move it along the cup wheel.
 - Move the sample along the cup wheel to the required stop position using the POSITION key ▲.
 - Press F2:SET STOP to record the stop position.
 - Press F3:MOVE TO REL.ZERO to move the sample back into the initial position.
 - Position the sample correctly in front of the cup wheel using the POSITION key ◀. If possible use F3:MOVE TO REL.ZERO.

Setting the Grinding Parameters All grinding parameters can be edited independently of each other.

See section Changing/Editing values for details on how to change values.

Method Screen A

Grinding Method: 1. EMPTY METHOD			
Cup Wheel :	MOP10	Speed :	4000 rpm
Y-Feed \updownarrow :	0.5 mm/s	Sweeps \updownarrow :	1
X-Inc. \updownarrow :	0.100 mm	Water :	ON
Rotation :	OFF		
Material to be removed		:	0.100 mm
F1:	RECALL REL. X POS.		
F2:	CHANGE TO STOP POS. MODE		

Cup Wheel

The Struers cup wheels are listed with their default values for speed. Also eight undefined wheels are included for your own choice of cup wheels. (See also configuration of cup wheels).

Speed

The speed of the cup wheel can be set between 300 and 5000 rpm in steps of 100 rpm. The default values for all Struers cup wheels are already saved together with the wheel definitions. Whenever a different cup wheel is selected the recommended speed will automatically be inserted.

Y-Feed \updownarrow

The movement of the sample along the Y-axis (parallel to the cup wheel).

Y-feed speed can be set between 0.1 mm/s and 10.0 mm/s (see Section Grinding Parameters in Chapter 3, Consumables for recommended setting).

Sweeps \updownarrow

The number of sweeps performed between any feed of the sample in X-direction. Every sweep consists of both a forth and a back movement of the sample. The number of sweeps can be set between 1 and 10.

X-Inc. \updownarrow

The feed of the sample towards the cup-wheel. The X-increment can be set in 0.005 mm steps between 0.005 and 1.000 mm.

Water

The cooling water can be set to either on or off. For all normal grinding operations the setting should be on.

Rotation

Sample rotation can be an advantage to get a more uniform scratch pattern on the surface and a better planeness. On Accutom-50, three different settings are possible.

Rotation: Off

The sample does not rotate.

Rotation: Continuous

The sample is rotating around its centre.

Press F3 to test sample rotation with the set parameters.

Rotation: Oscillating

The sample is rocking/oscillating around its centre.
The angle can be adjusted from 10° to 400°.

Press F3 to test sample rotation with the set parameters.

Removal Mode

On Accutom-50 two different ways of removing a certain amount of material can be used.

Material to be removed

If you know the exact position of a component underneath the sample surface, for example 0.125 mm, this distance can be pre-set.

- Move the sample as close to the cup wheel as possible, without getting in contact with each other using the POSITION ◀ key.
- Press the POSITION◀ key to move the sample slowly against the cup wheel just until there is just contact between sample and wheel.
- Press F1:RESET to set the relative X-position to zero.
- After having defined the amount of material to be removed, move the sample slightly back from the cup wheel, in the Y-direction, press START ◇, and after completing the grinding process Accutom-50 will stop precisely at the pre-defined depth.

Method Screen B

```
Grinding Method: 1. EMPTY METHOD
Cup Wheel : MOP10 Speed :4000 rpm
Y-Feed  0.5 mm/s Sweeps: 1
X-Inc.  0.100 mm Water : ON
Rotation : OFF
Stop grinding at rel. X pos.: - 0.200 mm
F1:RECALL REL. X POS.
F2:CHANGE TO REMOVAL MODE
```

Stop grinding at rel. X pos.

A function mainly used for mineralogy.

- Press F2:CHANGE TO REMOVAL MODE. Grind the ceramic plate of the vacuum holder so that it is plane parallel to the cup wheel.
- After the grinding press F1:RESET to set the relative X-position to zero and press F2: CHANGE TO STOP POS. MODE.
- Define the desired final thickness of the glass plate you want to grind by setting the “Stop grinding at position:” to the required value.
- Move the holder away from the cup wheel to allow the glass plate to be inserted.
- Place the glass plate on the vacuum holder.
- Reposition the glass plate close to the cup wheel.
- Press START \diamond to grind the glass down to the pre-set thickness.
- Press F1:RESET to set the relative X-position to zero.
- Move the holder with the glass plate away from the cup wheel to allow the sample to be inserted.
- Mount the sample on the glass plate and define the thickness of the sample by resetting the Stop grinding at position: to the new value.
- Reposition the sample close to the cup wheel.
- Press START \diamond to grind the sample down to the correct thickness.

Hint

To facilitate the production of thin sections, using permanent start positions, Accutom-50 can recall the original position settings.

When changing back to a grinding method which has been used before you will see F1:RECALL REL. X. POS Pressing F1 will recall the original relative positions. To move the sample to these positions, go to the X and Y position screen and press F3:MOVE TO REL. ZERO

Starting the Grinding

Before you start grinding make sure that the grinding chamber is clean to ensure a free flow of cooling water. The grate must not be clogged by debris.

- Position the sample correctly
- Set the appropriate grinding parameters.
- Make sure that the correct cup wheel is mounted
- Close the cover of the machine.
- Press START \diamond .

During Grinding

After START \diamond has been pressed, the display changes to the following:

GRINDING STARTED		
Pre-set Removal	:	0.100 mm
Actual Removal	:	0.100 mm
Remaining Process Time:		3 min
Y-Feed Speed \updownarrow	:	10.0 mm/s
Wheel speed	:	4000 rpm

Both the pre-set and the actual removal or thickness are displayed.

Stopping the Grinding

Automatic Stop

- Accutom-50 automatically stops the grinding process at the pre-set stop position or after the pre-defined amount of material has been removed.
- The sample is then retracted and the cup wheel is stopped.

Manual Stop

- The grinding process can be stopped at any time during the operation by pressing the STOP ∇ key. The sample remains in its actual position and the process stops there.
- The display changes back to what it was before grinding was started.

Continuing the Grinding

If the process has been stopped manually before the specified amount of material was removed the process can be continued by pressing START \diamond again. The process will start again and only the remaining amount of material will be removed.

To start all over again press ENTER \leftarrow twice to reset the removal value to the original setting.

3. Maintenance

Daily Service

- Clean the cutting chamber with a damp cloth. Do not use tap water as you risk overflow in the recirculation cooling water tank. Remember to remove all dirt particles from the grate.
- Clean the specimen holder head and the clamps for the dovetail feed.
- Clean the flanges.
- Clean the transparent cover with a damp cloth.

WARNING!

Do not use alcohol, acetone or similar solvents.

Checking the Recirculation Unit

The cooling unit should be checked for cooling water after 8 hours use or at least every week. The unit must be refilled if the flushing pump cannot reach the cooling water. Remember to add Struers Additive for Cooling Fluid: One part of Additive for 33 parts of water.

To check the concentration of additive, use a refractometer. Concentration = Brix value. The concentration of additive should always lie between 2.7 and 3.3 %. Add Struers Additive for Cooling Fluid if the concentration is too low.

IMPORTANT



Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to add Struers Additive each time you refill with water.

Weekly Service

- Clean the specimen holder(s): movable parts, dovetail feeds and screws. Lubricate with acid free oil.
- Clean the cutting chamber and cover thoroughly. Remember to remove all dirt particles from the grate.

Refilling the Cooling Water Tank

Emptying and Cleaning the Tank

Replace the cooling water in the Recirculation Cooling Unit at least once a month.

- Carefully pull out the drawer and lower the hinged front plate.
- Carefully pull out the flat tank.
- Remove the cover plate and empty the tank into a drain approved for waste chemicals.
- Clean the tank, the dividers and the cover with tap water. Place the dividers correctly in the tank and replace the cover plate.
- Press the cover plate into the tank and make sure it is seated firmly.

IMPORTANT

The container and the cover plate have to be placed correctly, otherwise the pump cannot suck up the water.

IMPORTANT

Flush the recirculation system with clean water if Accutom is not to be used over longer periods of time. This will prevent any dried residue of cutting material from damaging the inside of the pump.

Refilling the Tank

- Refill the tank with 3.88 l of water and 120 ml Struers Additive . The water level should be 5 mm below the edge of the front hole in the tank cover.
- Push the drawer with the recirculation tank back into place.

IMPORTANT

Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to add Struers Additive each time you refill with water.

Maintenance of Cut-off Wheels

Maintenance of Diamond and CBN Cut-off Wheels

The precision of diamond and CBN cut-off wheels and thus the cut depends on how carefully the following instructions are observed:

- Never expose the cut-off wheel to overload, such as a heavy mechanical load, or heat.
- Store the cut-off wheel in a dry place, horizontally on a plane support, preferably under light pressure.
- A clean and dry cut-off wheel does not corrode. Therefore, clean and dry the cut-off wheel before storing. If possible, use ordinary detergents for cleaning.
- Regular dressing of the cut-off wheel is also part of the general maintenance (see **ADVANCED OPERATIONS**).

Storing of Abrasive Cut-off Wheels

These cut-off wheels are sensitive to humidity. Therefore, do not mix new, dry cut-off wheels and used, humid ones. Store the cut-off wheels in a dry place, horizontally on a plane support.

Maintenance of Cup Wheels

The precision of diamond cup wheels and thus the grinding depends on how carefully the following instructions are observed:

- Never expose the cup wheel to overload, such as heavy mechanical load, or heat.
- Store the cup wheel in a dry place, horizontally on a plane support.
- A clean and dry cup wheel does not corrode. Therefore, clean and dry the cup wheel before storing. If possible, use ordinary detergents for the cleaning.
- Regular dressing of the cup wheel is also part of the general maintenance (see **ADVANCED OPERATIONS**).


Reference Guide



Table of Contents	Page
1. Advanced Operations	
Configuration Menu	36
Cut-off Wheels	37
Cup Wheels	38
Options.....	39
Manual Functions.....	40
Sample Rotation	41
Recirculation Pump.....	42
Dressing.....	43
Method Options.....	46
Saving a Method.....	46
Copying a Method.....	47
Inserting a Method	48
Resetting a Method.....	49
Editing Names.....	50
Name Editing Principles.....	51
2. Accessories	52
3. Consumables	
Cut-off Wheels	53
Optimising the Cutting Results.....	53
Cutting Parameters	53
Common Cutting Questions.....	54
Cup Wheels.....	54
Optimising the Grinding Results.....	55
Grinding Parameters.....	55
Consumables	55
4. Trouble-Shooting	56
5. Technical Data	58
6. Menu Overview	59

Cut-off Wheels

If you want to use other than Struers cut-off wheels, these must be configured to contain the correct values for force limit, wheel thickness and speed.

To configure your own cut-off wheels,

 Press MENU UP/DOWN ▲▼ to select CUT-OFF WHEELS.

↓
 Press ENTER  to activate the CUT-OFF WHEELS MENU.



↓

CONFIGURATION OF CUT-OFF WHEELS	
Cut-off Wheel :	USER2
Force Limit :	MEDIUM
Wheel Thickness:	0.50 mm
Speed :	3000 rpm
F1:PREV. WHEEL F2:NEXT WHEEL F4:RENAME	

↓
 Press MENU UP/DOWN ▲▼ to select the value to be changed.

↓
 Press ENTER  to edit the value.

After changing the value,

↓
 Press ENTER  to accept the new value and to return to the previous screen.


↓
F4 Press F4 to rename the cut-off wheel.
(see Editing Names for details).

Cup Wheels

If you want to use non-Struers cup wheels, be sure to configure the correct values for the speed.


Press MENU UP/DOWN ▲▼ to select CUP WHEELS.


Press ENTER  to activate the cup wheels menu.

```
CONFIGURATION OF CUP WHEELS
Cup Wheel      :                USER 1
Cup Wheel Speed:                4000 rpm
F1:PREV. WHEEL F2:NEXT WHEEL F4:RENAME
```

Press ENTER  to edit the value.

After changing the value,

Press ENTER  to accept the new value and to return to the previous screen.

Press F4 to rename the cup wheel. (see Editing Names for details).

Options

In the options menu various items can be configured.

Language

See how to change the language in the Getting Started section of this Instruction Manual.

X-retract 

After the stop position is reached, the sample is normally moved straight back to the relative zero position. To avoid any contact with the cut-off wheel during that movement, the sample can be retracted from the wheel before it is repositioned by setting X-retract to ON. After reaching the relative zero position, the sample is then moved back into the original position.

IMPORTANT

Make sure that the Y-stop position is set correctly. If the sample is not cut through before the sample is retracted the cut-off wheel might be damaged.

Contrast

The contrast setting of the display can be adjusted.

No. of Final Grinding Sweeps

When grinding, the number of sweeps carried out before the sample is moved towards the cup wheel can be adjusted. To save time, very often only a small number of sweeps is set during normal operation. However, to achieve a better surface finish, the number of final sweeps can be set to a maximum of 15.




F1: DEFAULT VALUE

Pressing F1 will reset the inverted option to the factory adjusted default value.


Manual Functions


Several functions on Accutom-50 can also be carried out manually. All of these functions are available from the Manual Functions Menu.

Press Esc until you reach the Main Menu.

Sample Rotation

 Press MENU UP/DOWN  to select SAMPLE ROTATION.

↓
 Press ENTER  to activate the SAMPLE ROTATION menu.

↓

SAMPLE ROTATION	
Sample Rotation:	CONTINUOUS
Speed Level :	1
Direction :	> <
F1:START ROTATION	

↓
 Press MENU UP/DOWN  to select the value to be changed.

↓
 Press ENTER  to edit the value.


After changing the value,

↓
 Press ENTER  to accept the new value and to return to the previous screen.

↓
F1 Press F1 to start rotation and check that everything has been aligned correctly. Press F1 again to stop rotation.


Recirculation Pump

 From the MANUAL FUNCTIONS MENU press MENU UP/DOWN  to select RECIRCULATION PUMP.

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
Dressing

Diamond and CBN cut-off wheels and cup wheels may have to be dressed to open up the wheel again and achieve maximum effect after cutting or grinding of ductile materials. This can be done using the following routine:

Dressing of Cut-off Wheels

- Mount the cut-off wheel or cup wheel to be dressed.
- Clamp the dressing stick in the sample holder and position it correctly.
 - The dressing stick has to be in front of the wheel so that the wheel is cutting slices of the dressing stick.

 ▲ From the MANUAL FUNCTIONS MENU Press MENU
 ▼ UP/DOWN  ▲▼ to select
DRESSING OF CUT-OFF WHEEL.


 Press ENTER ↵ to activate the dressing option.

DRESSING OF CUT-OFF WHEEL	
Feed Speed :	0.5 mm/s
Wheel Speed :	5000 rpm
Cooling Water:	ON
Y-Distance \blacktriangle \blacktriangledown :	10.0 mm
F1: START DRESSING	

↓
 Press MENU UP/DOWN  \blacktriangle \blacktriangledown to select the value to be changed.

↓
 Press ENTER \blacktriangleleft to edit the value.


↓
 Press ENTER \blacktriangleleft to accept the new value and to return to the previous screen.

↓
 Press F1 to start the dressing operation.

Dressing of Cup Wheels


- Mount the cup wheel to be dressed.
- Clamp the dressing stick in the sample holder and position it correctly.
 - The dressing stick has to be positioned so that the cup wheel is grinding a piece of the end of the stick.

 Press MENU UP/DOWN ▲▼ to select DRESSING OF CUP WHEEL.

↓

 Press ENTER  to activate the dressing option.

↓

 Press MENU UP/DOWN ▲▼ to select the value to be changed.

↓

 Press ENTER  to edit the value.

↓

 Press ENTER  to accept the new value and to return to the previous screen.

↓

 Press F1 to start the dressing operation.

Method Options

Saving a Method

While working with a specific method, you can save the changes you have made in the database.

- As soon as a parameter is changed in a cutting or grinding method, F4:SAVE will be shown on the bottom line of the display.

```
Cutting Method:  1. EMPTY METHOD
Wheel   :      B0D15   Speed  : 2700 rpm
Feed    : 0.200 mm/s   Force   : MEDIUM
Rotation:      OFF
Water   :      ON
F1:MULTI. CUT                               F4:SAVE
```

F4 Press F4:SAVE to save the changes after all necessary parameters have been altered.

```
SAVE METHOD
Save changes in method
1. EMPTY METHOD ?
ENTER: YES                               ESC: NO
```

 Press ENTER ↵ to save the modified method.

Unless you have already named the method, the name will automatically change from EMPTY METHOD to UNNAMED METHOD. This will always show that at least one parameter has been changed compared to a default method.

Important


When saving changes, the original method will be overwritten. If you want to keep the original method, you should make a copy of the method with a new name, thus making the changes in a copy rather than changing the original method. See Copying a Method.

Copying a Method

Copying a method is a shortcut to creating a new cutting or grinding method on the basis of an existing one.

- Select the cutting or grinding method you want to copy from.

SELECT CUTTING METHOD ↓		
1. UNNAMED METHOD	2. EMPTY METHOD	
3. EMPTY METHOD	4. EMPTY METHOD	
5. EMPTY METHOD	6. EMPTY METHOD	
7. EMPTY METHOD	8. EMPTY METHOD	
9. EMPTY METHOD	10. EMPTY METHOD	
11. EMPTY METHOD	12. EMPTY METHOD	
F1: COPY	F3: RESET	F4: RENAME


↓
 Press ENTER ↵ if you would like to see the method parameters.

↓
 Press Esc to return to the above display.

↓
 Press F1: COPY to copy the method to the buffer.



↓

COPY METHOD	
Do you want to copy method	
1. UNNAMED METHOD to buffer?	
ENTER: YES	ESC: NO

↓
 Press ENTER ↵ to accept copying the method to the buffer.



Inserting a Method

If you want to insert the method in a different method:


 Press MENU UP/DOWN ▲▼ to select the method, in which you want to insert the method.


If it is not an EMPTY METHOD:

 Press ENTER  to see the method.

 Press Esc to return to the above display.

 Press F2:INSERT to insert the method from the buffer.

```
INSERT METHOD
-----
Do you want to copy from buffer to
method 4. EMPTY METHOD      ?
ENTER: YES                      ESC: NO
```

 Press ENTER  to accept insertion of the method.

Resetting a Method

If a method is not in use any more, it should be reset.

The parameters will change to default values, which can easily be changed to a new method later on.

The method name will change to: EMPTY METHOD, telling you, that you can copy to or modify in this method without replacing any valuable data.

- Select the method you want to reset.

SELECT CUTTING METHOD A	
1.UNNAMED METHOD	2.EMPTY METHOD
3.EMPTY METHOD	4. UNNAMED METHOD
5.EMPTY METHOD	6.EMPTY METHOD
7.EMPTY METHOD	8.EMPTY METHOD
9.EMPTY METHOD	10.EMPTY METHOD
11.EMPTY METHOD	12.EMPTY METHOD
F1: COPY	F2: INSERT
F3: RESET	F4: RENAME

↓
 Press ENTER ↵ to see the method and to make sure that you do not delete valuable data.

↓
 Press Esc to return to the above display.

↓
 Press F3:RESET.

↓






RESET METHOD	
Do you want to reset this method ?	
4. UNNAMED METHOD	
ENTER: YES	ESC: NO

↓
 Press ENTER ↵ to reset the method.

Editing Names

The names for the methods, cut-off wheels and cup wheels can be edited and changed to suit your preference.

- Select the method, cut-off wheel or cup wheel you want to rename.


Name Editing Principles


- Place the main cursor on the character you want to change, using F1:LEFT or F4:RIGHT. Use MENU DOWN  to move to the character set in the next line. An auxiliary cursor in the text line shows the position in the method name.
 - Write the new name using the following keys:
 - F1** Moves the main cursor to the left
 - F2** Inserts a space in the text
 - F3** Deletes one character to the left in the text
 - F4** Moves the main cursor to the right
- ↵ ENTER places the new character in the method name and moves the auxiliary cursor in the name to the right. Repeat the procedure for each character.

↓ Write the new name using the above name editing procedures.

Esc Press Esc to leave the editor again.

↓

 Press ENTER ↵ to accept the new name.

2. Accessories

Flange Sets

Specification	Cat. No:
For cut-off wheels. 65 mm dia. For use on Accutom-50 to achieve highest precision.	04946902
For cup wheels.	05016901

Specimen Holders

Specification	Cat. No:
For general use. Vice type with max. opening 60 mm	04946903
For round or square specimens. Teardrop type with max. opening \varnothing 40 mm / \varnothing 1 1/2"	04946904
For irregular specimens, with 7 screws. Max. width 40 mm / 1 1/2"	04946905
With goniometer	04276911
For adhering specimens	04276912
With ceramic vacuum chuck for thin sections	04276913
For small specimens. Vice type	04276915
Double parallel vice	04946909
<i>Joints to be mounted between the dovetail and the specimen holder</i>	
Tilting joint with max. angle $\pm 10^\circ$	04946906
Angling joint with max. angle $+30/-90^\circ$	04946908
<i>Base Plate</i>	
With dovetail. For mounting of other types of specimen holders	04276914

3. Consumables

Cut-off Wheels

Please refer to the Selection Guide in the Struers [Cut-off Wheels brochure](#).

Accutom-50 accepts cut-off wheels with a diameter of 75 mm (smallest flange is \varnothing 42 mm) up to 152 mm. Hole 12.7 mm.

IMPORTANT!
Always use large flanges (65 mm) for the highest possible accuracy.

Optimising the Cutting Results Cutting Parameters

Use the following table to select proper wheel and cutting parameters according to the sample material.

Recommended Cutting Parameters				
Material	Hardness	Force limit	Feed speed [mm/s]	Wheel speed [rpm]*)
Ceramics, minerals and crystals	> HV 800	LOW	0.005-0.15	5000
		LOW	0.005-0.20	4000
		HIGH	0.005-0.30	3200
		HIGH	0.005-0.30	2700
Sintered carbides and hard ceramics	> HV 800	MEDIUM	0.005-0.25	3200
		MEDIUM	0.005-0.25	2700
Extremely hard ferrous metals	> HV 500	MEDIUM	0.005-0.25	5000
Hard and very hard ferrous metals	HV 350-800	MEDIUM	0.05-0.30	1000-5000
		MEDIUM	0.05-0.30	1000-5000
Hard and very hard ferrous metals with larger dimensions	HV 350-800	MEDIUM	0.05-0.30	1000-5000
Soft and medium soft metals	HV 30-350	MEDIUM	0.05-0.30	1000-5000
		MEDIUM	0.05-0.30	1000-5000
Soft and ductile non ferrous metals	HV 70-400	MEDIUM	0.05-0.30	1000-5000
Plastics and very soft metals	< HV 100	MEDIUM	0.05-0.30	max. 1200

*) To obtain low wheel wear and better surface quality, always use highest recommended wheel speed.

Common Cutting Questions

The following table shows the possible answers to a number of common cutting questions:

Optimising the Cutting Results	
Objective	How to achieve it
Better surface quality	Use lowest recommended feed speed and no specimen holder rotation.
Lower wheel wear	Use the lowest recommended feed speed and no specimen holder rotation. This is especially important when using resin bonded wheels and all abrasive cut-off wheels.
Problems with abrasive cut-off wheels?	Abrasive cut-off wheels should not be used outside their recommended feed speed range. At lower than recommended feed speeds they will produce irregularly cut surfaces. At higher feed speeds excessive wheel wear will occur, along with increased risk of wheel breakage.
Flatter samples	Use primarily low feed speeds, highest recommended wheel speed, largest possible flanges and no specimen holder rotation. The initial kerf is especially critical. If the initial feed speed is too high the wheel will bend and start cutting at an angle. Such a cut will never end up flat.
Better parallelism	Use the lowest recommended feed speed.
Faster cutting	Orientate the sample so that the wheel will cut the smallest possible cross-section and then use maximum recommended feed speed.

Cup Wheels

Material	Hardness	Size	Bond	Abrasive	Grit	Cat. No:
Ceramics, minerals and crystals	> HV 800	ø100	Soft metal	Diamond	#150	M0P10 40800082
Sintered carbides, hard ceramics and hard composites	> HV 600	ø100	Resin	Diamond	#300	B0P10 40800083

Optimising the Grinding Results

Grinding Parameters

Use the following table to select proper wheel and grinding parameters according to the sample material.

Recommended Grinding Parameters							
Material	Hardness	Wheel	Precision	Feed speed [mm/s]	X-increment	No. of final sweeps	Wheel speed [rpm]
Ceramics, minerals and crystals	> HV 800	M0P10	High	0.1-0.2	5-10 µm	10	4000
			Medium	0.2-4.0	10-20 µm	5	
			Low	4.0-8.0	20-30 µm	2	
Sintered carbides, hard ceramics and hard composites	> HV 600	B0P10	High	0.1-0.3	5-10 µm	10	4000
			Medium	0.3-0.5	10-20 µm	5	
			Low	0.5-1.0	20-30 µm	2	

Consumables

Specification	Cat. No:
<p><i>Corrozip</i> <i>Additive for Cooling Fluid</i> Environment friendly. To protect the machine from corrosion and to improve cutting and cooling qualities. 1 l 5 l</p>	<p>49900045 49900046</p>
<p><i>Cutting Fluid</i> Water free Cutting Fluid for cutting of water-sensitive materials, 5 l</p>	<p>49900030</p>

4. Trouble-Shooting

Error Messages		
Display Message	Explanation	Action
Searching for XY-reference position	The machine searches for the initial reference positions	Wait for the search process to finish
Reference pos. not found, close cover	Accutom-50 was started with the cover open and did not search for the reference position	Close the cover, wait for the reference search to finish and press START ◊ again
Cover open!	You are trying to start the cutting/ grinding process or a dressing while the cover is open	Close the cover and press START ◊ again
Process stopped by open cover	You have opened the cover during a cutting/grinding process or a dressing process	Close the cover and restart the interrupted process
Process running	You are trying to start two different processes at the same time	Wait for the previous process to finish and try again
Y-stop position outside cutting range Y-stop position outside grinding range	The set stop position is outside the maximum range in the Y-direction. (max. 105.0 mm)	Adjust Y-stop position.
Y-stop position less than Y-relative	You try to start a cutting process when the relative Y-position is larger than or equal to the Y-stop	Position the specimen holder correctly and/or change the Y-stop position.
Feed less than 90% of pre-set value	The force has exceeded the force limit (LOW, MEDIUM, HIGH) and the feed speed is automatically reduced	Decrease the feed speed or increase the force limit next time you cut the same material
Cutting stopped, feed too low	The feed speed has been reduced to less than 10% due to the force exceeding the force limit	Decrease the feed speed, increase the force limit or change the cut-off wheel and re-start the cutting process
X-relative must be negative before start	In connection with the grinding to a STOP POSITION the relative X-position is positive	Position your sample correctly by moving the sample holder to the right
X-relative must be less than stop pos.	In connection with the grinding to a STOP POSITION the relative X-position is larger than the pre-set stop position	Position your sample correctly by moving the sample holder to the right
Too many cuts, only xx cuts possible	The no. of cuts and/or the slice thickness exceed the no. of possible cuts	Re-position your sample or press START ◊ again to allow for automatic adjustment of the no. of slices
Removal too large, move X-table	In connection with the grinding mode REMOVAL, the removal is set higher than the possible movement of the table	Reduce the amount of material to be removed or clamp the sample differently

*Accutom-50
Instruction Manual*


Error Messages		
Display Message	Explanation	Action
Cutting or grinding method not selected	START ◊ was pressed before a cutting or grinding method was selected	Select a method and press START ◊ again
Dressing range exceeded	The parameter Y-distance is set to a value larger than the possible movement of the table	Adjust the Y-distance and press START ◊ again
X-MIN sensor not activated	The inductive sensor cannot detect the end position of the X-table	Check for mechanical blockage of the table. If that is not the case, call a Struers service technician
X-MIN sensor not deactivated	The inductive sensor for the X-position cannot be deactivated	Check for mechanical blockage of the table. If that is not the case, call a Struers service technician
Encoder error	The Y-table does not move, or the encoder is defective	Check for mechanical blockage of the table. If that is not the case, call a Struers service technician
Cutting/grinding motor overloaded	The main motor has been exposed to a high load for a long time	Wait until the motor has cooled down. "Cutting motor ready after overload" will be shown on the display. Adjust process parameters to avoid repetitions
Machine Problems	Explanation	Action
The recirculation pump is not working correctly	The water level in the recirculation tank is too low	Check that there is sufficient water in the recirculation tank
	The water outlet is clogged	Remove the cooling tubes and run the recirculation pump. Flush the cooling tubes with clean water before replacing
	The recirculation pump is contaminated with cutting residue	Flush the pump with clean water by gently pressing fresh water into the inlet tube

5. Technical Data

Subject	Specifications	
Cutting	<p><i>Cutting Speed:</i> 300-5000 rpm (adjustable in steps of 100 rpm) <i>Feed Speed:</i> 0.005-3.000 mm/s (adjustable in steps of 0.005 mm/s) <i>Max. Positioning Speed:</i> Y = 13 mm/s, X = 13 mm/s <i>Force Limits:</i> LOW: approx. 20N; MEDIUM: approx. 40N; HIGH: approx. 60N</p>	
Positioning Range	<p>Y direction: 105 mm (precision 0.1 mm) X direction: 60 mm (precision 0.005 mm)</p>	
Sample Size	<p><i>Max. length of Cut-off Sample:</i> 30 mm, 140 mm at $\varnothing 20$ mm <i>Max. length of Sample to be clamped:</i> 225 mm <i>Max. Sample Cross Section:</i> $\varnothing 127$mm cut-off wheel and $\varnothing 42$mm flange: $\varnothing 40$mm (without rotation) $\varnothing 80$mm (with rotation) $\varnothing 152$mm cut-off wheel and $\varnothing 42$mm flange: $\varnothing 50$mm (without rotation) $\varnothing 100$mm (with rotation)</p>	
Cut-off Wheels	<p><i>Diameter:</i> $\varnothing 75$ to $\varnothing 152$ mm / 3 to 6" <i>Max. Thickness:</i> 2 mm / 0.08" <i>Hole:</i> $\varnothing 12.7$ mm / 0.5"</p>	
Cup Wheels	<p><i>Diameter:</i> $\varnothing 100$ mm / 4" <i>Hole:</i> $\varnothing 12.7$ mm / 0.5"</p>	
Recirculation Cooling Unit	<p><i>Contents:</i> 4 l <i>Flow:</i> 800 ml/min</p>	
Motor	550 W at 5000 rpm Continuous Torque 0.8 Nm	
Safety Standard	Please refer to the Declaration of Conformity	
Noise Level	Approx. 65 dB (A) measured at idle running, at a distance of 1.0 m / 39.4" from the machine.	
Altitude	Min. 50 m. below water level. Max 3000 m above water level	
Surrounding temperature	5-40°C/41-104°F	
Humidity	0-95% RH non condensing	
Software and Electronics	<p><i>Display:</i> 8 x 40 characters <i>Controls:</i> touch pad <i>Database:</i> 20 cutting methods + 10 grinding methods <i>PCB fuses:</i> 3.15A , 1,6A SUB mini fuses</p>	
Power Supply	50-60Hz (max. load 2.7A) 1/3 x 220 - 240V	
Dimensions and Weight	Width	510 mm / 20"
	Depth	700 mm / 27.6"
	Height	270 mm / 10.6"
	Weight	45 kg / 99 lbs

6. Menu Overview

Accutom-50
Menu Structure

Quick Reference Guide

Cutting with Accutom-50

Loading a Cutting Method

- Press Esc until the Main Menu appears.
- Select Cutting Methods and press Enter ↵.
- Select the method you want to use and press Enter ↵.

Positioning the Sample

- Press one of the positioning keys ▼ ▲ to change to the Y-Positions screen.
- Place the specimen holder with the sample in the specimen holder head and clamp it.

The “heavy” arrows on the display indicate the direction in which the sample holder is moved, for example:

X-POSITIONS	➡➡
Y-POSITIONS	⬆️⬇️

If you already know the cutting length (i.e. the diameter of the sample):

- Position the sample correctly in front of the cut-off wheel using the positioning keys.
- Press one of the positioning keys ▼ ▲ to change to the Y-Positions screen.
- Press F1:RESET to reset the relative Y-position.
- Press Enter ↵ and use the MENU arrows ⬅️➡️ to adjust the Y-stop value. Press Enter ↵ again to accept the new value, or Esc. to cancel.

If you do not know the cutting length in advance:

- Position the sample correctly in front of the cut-off wheel using the positioning keys.
- Press F1:RESET to reset the relative X-position.
- Press one of the positioning keys ▼ ▲ to change to the Y-Positions screen.
- Press F1:RESET to reset the relative Y-position.
- Move the sample to the right using the positioning key ▶️ so that it can pass along the cut-off wheel.
- Move the sample along the cut-off wheel to the desired stop position using the positioning key ▲. Press F2:SET STOP to record the stop position.
- Press F3:MOVE TO REL: ZERO to move the sample back to the initial Y-position.
- Press one of the positioning keys ⬅️➡️ or ⬆️⬇️ to change to the X-Positions screen.
- Press F3:MOVE TO REL: ZERO to move the sample back to the initial X-position.

Starting the Process

- Press START ⬠ to start the cutting process.

Grinding with Accutom-50

Loading a Grinding Method

- Press Esc until the Main Menu appears.
- Select Grinding Methods and press Enter ↵.
- Select the method you want to use and press Enter ↵.

Positioning the Sample

- Place the specimen holder with the sample in the specimen holder head and clamp it.

Pre-set Amount of Material to be Removed

If you know the exact position of a component underneath the sample surface, for example 0.125 mm, this distance can be pre-set:

- Move the sample slowly against the cup wheel just until there is contact between sample and wheel.
- Press F1:RESET to set the relative X-position to zero.
- Define the amount of material to be removed, then move the sample away from the cup wheel in Y-direction.

Stop Grinding at Position

- Grind the ceramic plate of the vacuum holder so that it is plane parallel to the cup wheel.
- After the grinding press F1:RESET to set the relative X-position to zero.
- Define the thickness of the glass plate you want to grind by setting the Stop grinding at position: to the wanted value.
- Move the holder away from the cup wheel to allow the glass plate to be inserted.
- Place the glass plate on the vacuum holder.
- Reposition the glass plate close to the cup wheel.
- Press START ◊ to grind the glass down to the pre-set thickness.
- Press F1:RESET to set the relative X-position to zero.
- Move the holder with the glass plate away from the cup wheel to allow the sample to be inserted.
- Mount the sample on the glass plate and define the thickness of the sample by resetting the Stop grinding at position: to the new value.
- Reposition the sample close to the cup wheel.

Starting the Process

- Press START ◊ to start the grinding process.

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Date of Release 15.11.2012

Accutom-50
Spare Parts and Diagrams

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2012.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Telefax: +45 44 600 801

Spare Parts and Diagrams

	Table of contents	Drawing
Drawings	Accutom-50, complete.....	15010001A
	Cabinet	15010021D
	Water Cock.....	14940044B
	Recirculation Tank	14940041G
	Lid for Recirculation Tank, complete	15010427
	Bottom	15010051N
	Back plate.....	15010054E
	Various Parts	
Diagrams	Block and Wiring Diagram.....	15013052D
	Wiring Diagram	15013452D





Some of the drawings may contain position numbers
not used in connection with this manual.

Spare Part List for Accutom-50





Drawing	Pos.	Spare Part	Cat no.
15010001		Accutom-50, Complete	
	130	Flange ø42, back	14940527
	140	Flange ø42, front	14940529
	150	Flange screw M5	427MP012
	160	Backstop with 1 hole	14278017
	170	Backstop with 2 holes	14278016
	180	Pressure spring	R2600718
15010021		Cabinet	
	15	Seal for cover	R4940218
	25	Safety switch	2SS00012
14940044	70	Cover	R4940206
		Water cock	
	10	Pipe branch	14940420
	15	O-ring 7.65-1.78	RIO17812
	20	O-ring 10.82-1.78	130MP057
	30	Socket for water nozzle	R4940422
	40	End piece for water nozzle	14940428
60	Water nozzle	R4940421	
14940041	70	Stopper for water nozzle	R4940429
		Recirculation Tank	
	10	Bottom plate	14949005
	20	Fitting for Drawer	2JX04230
	50	Front plate for Recirculation tank	14949006
	60	Foil for recirculation drawer	15010315
15010427	90	Tension spring	R4940416
		Lid for Recirculation Tank, complete	
	10	Lid for recirculation tank, complete	15010427
15010051	20	Sealing list 10x4	RIP00401
	30	Sealing profile	RIP00203
15010051		Bottom	
	20	Rubber foot, grey	R4910140
	910	Suction hose ø8/ø12 (3 Pcs)	14942908

Spare Part List for Accutom-50

Drawing	Pos.	Spare Part	Cat no.
15010054		Back Plate	
	50	Fuse glass tube 10A T (10 pcs.)	392MP075
		Various Parts	
		Grate	14940220
		Grate, large	R4940222
		Backstop spanner for cut-off wheel mounting	14940701
		2 pin spanner for cut-off wheel mounting	427MP085
		Allen key 2 (10 pcs.)	RGR00020
		Allen key 2.5 (10 pcs.)	RGR00025
		Allen key 3 (10 pcs.)	RGR00030
		Allen key 5 (5 pcs.)	RGR00050
		Cooling tube for Cup Wheel	15010406
		Connection hose for vacuum	15010402

A: Pos 90 fjermet, Pos 200 + 210 tilf. PMP 040789			
Matrik:	Overf. beh.:	Målestok:	Ikke ang. tol. efter DS/ISO 2768-
Struers	Projeksjonsmetode:	1:2	Dato: 201295 Sign.: BMJ
Accutom-50, komplet			Erst.: 15010001A

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





















*Note:



Old Version	New Version	Signal	Connect to
GN	GN	Channel A	YE
YE	WR	Channel B	WR
BN	RD	VCC	RD
WH	BK	GND	BK

Part No.	Part Name	Quantity	Unit
ACCUTOM 5/50	ACCUTOM 5/50	1	PCB
15013452	15013452	1	PCB



Q1 POWER ENTRY MODULE

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T1 TRANSFORMATOR

KEYBOARD


X-MIN. LIMIT.:


EXTRA INPUT.

EXTRA INPUTS.

SAFETY GUARD CLOSED



CUTTING MOTOR READY FOR RUN


ROTARY OPTICAL ENCODERS AT X-POS.

*Note

Z3


CUTTING MOTOR M1

MB1 COOLING-FAN


EXTRA +5V DC OUTPUT FOR COOLING FAN


YO1 OSCILLATING PUMP

COOLING WATER.


CONTACTOR

CUTTING MOTOR.


SERVICE COM.

M4 STEPMOTOR


DRIVE AT SAMPLE.

OPTION AT ACCUTOM 5.


M3 STEPMOTOR

DRIVE AT X-POS.

M2 STEPMOTOR

DRIVE AT Y-POS.

LCD 240 X 56 DOT OPTREX 50316 CCFT

*Note:

Old Version	New Version	Signal	Connect to
GN	GN	Channel A	YE
YE	WH	Channel B	WH
BN	RD	VCC	RD
WH	BK	GND	BK

Rev B : Motor Control, S1, S2, S3 added.
Rev C : Motor Control, Connections changed
Rev D : S1 changed to S2

English

Declaration of Conformity

 Struers

Manufacturer Struers A/S
Pederstrupvej 84
DK-2750 Ballerup, Denmark
Telephone +45 44 600 800

Herewith declares that

<i>Product Name:</i>	Accutom-50
<i>Type No:</i>	501
<i>Machine Type:</i>	Precision cut-off machine

is in conformity with the provisions of the following directives:

Safety of Machinery 2006/42/EC according to the following standard(s):
EN ISO 12100:2011, EN60204-1:2006/AC:2010.

Low Voltage Directive 2006/95/EC according to the following standard(s):
EN60204-1:2006/AC:2010.

EMC-Directive 2004/108/EC according to the following standard(s):
EN61000-6-1:2007, EN61000-6-3:2007/A1:2011.

RoHS 2011/64/EC according to the following standard(s):
EN 50581:2012.

Supplementary Information The equipment complies with the American standards:
FCC part 15, subpart B and UL508.

The above has been declared according to the global method, module A

Date: 14.01.2014

Christian Skjold Heyde,
Vice President, R & D and Production, Struers A/S

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... and grinding machine

Common features of Accutom-5/-50

- Positioning accuracy of 5 μm
- Pre-set, constant feed speed (0.005 to 3 mm/sec)
- Adjustable force limit
- Sample rotation or oscillation
- Variable cut-off wheel speed up to 3,000 rpm
- Menu control and digital read-out
- Built-in cut-off wheel list
- Database with user-definable cutting methods
- Built-in recirculation cooling unit
- Enclosed cutting chamber with safety switch

Common application areas:

Precise and deformation-free cutting of:

- Metals
- Electronic components
- Ceramics
- Composites
- Crystals
- Sintered carbides
- Minerals
- Fibre reinforced materials
- Biomaterials (teeth and bones)

Additional features of Accutom-50

- Serial cutting of slices (Multi-cut)
- Grinding of thin sections
- Target grinding in components
- Variable cut-off wheel speed up to 5,000 rpm
- Database for both cutting and grinding methods
- Improved menu control

Additional application areas:

- Precise serial cutting of all materials
- Thin sections of minerals
- Thin sections of ceramics
- Target grinding in components

Oscillating movement is used for cutting oblong samples, thus reducing the cutting time considerably

Accutom-50 is based on the known and well tested technology of Accutom-5. Sample positioning within 5 µm combined with constant feed speed and sample rotation or oscillation represents the utmost in precision cutting. The addition of fully automatic serial cutting and precision grinding facilities, together with even higher wheel speeds, makes Accutom-50 the most universal tool in precision manufacturing of samples.

Working principle

Both Accutom-5 and Accutom-50 are built on the same, precisely manufactured, cast frame. The use of step motors to control the movement of the sample in both X- and Y-directions enables very precise positioning.

Positioning accuracy of 5 µm

The high resolution of the step motor makes it possible to position the sample with an extraordinary precision of 5 µm. The total movement of the sample in the X-direction is 60 mm.

Pre-set, constant feed speed

The feed speed can be pre-set and controlled in the range of 0.005 to 3 mm/sec. This wide range allows even the most delicate materials to be cut. Both the pre-set and the actual speeds are shown on the display. If the preset force limit is reached, the feed speed will be reduced and the process stopped if the speed has fallen below a specified limit.

Adjustable force limit

To protect thin cut-off wheels or fragile samples, the force limit can be set to three different levels: Low, Medium and High.

Sample rotation or oscillation

Sample rotation or oscillation reduces the contact area between cut-off wheel and sample to a minimum. This results in a decrease in frictional heat and allows fast cutting of extremely hard materials.

Variable cut-off wheel speed up to 3000 rpm

The motor of Accutom-5 has variable cut-off wheel speeds from 300 up to 3000 rpm in steps of 100 rpm. This allows both high and low speed cutting and enables a more accurate choice of wheel speed for different cut-off wheel diameters.

Menu control and digital readout

Accutom-5 is operated from a touch pad on the ergonomically designed, slanting front. All parameters, both for positioning and cutting, are presented on the LC-Display.

Built-in cut-off wheel list

A list of all Struers cut-off wheels is stored in an internal database together with their name, respective thickness and force limit.

Accutom-50 can be programmed to automatically cut series of slices from the same samples

Database with user-definable cutting methods

After establishing a method for a specific material, all settings can be saved in a database and recalled at any time. Accutom-5 has a storage capacity of ten methods.

Built-in recirculation cooling unit

With a capacity of four litres, and a very high through-put, the recirculation cooling unit ensures a constant low temperature during the cut, and an efficient removal of cutting debris.


Enclosed cutting chamber with safety switch

The fully transparent cover is equipped with a safety switch to prevent the cut-off wheel from starting whilst the cover is open. The sample can be inserted and positioned but the cut-off wheel itself can only be activated when the cover has been closed properly, thus providing maximum safety for the operator.


Accutom-50 provides the following additional features:

Serial cutting of slices (Multi-cut)

Accutom-50 can be programmed to automatically cut series of slices from the same sample. The advanced software will also compensate for the wheel thickness and for slight drift of the cut-off wheel. This results in very precise sections with an identical thickness, as required in failure analysis or the examination of biomaterials.

After the process has been started, the display changes and the relevant parameters are shown

Grinding on Accutom-50 can be carried out in two different ways: The grinding operation can be stopped at a specified position, or a well-defined amount of material can be removed

Grinding of thin sections

The thickness ground during each different operation associated with a thin section can be easily adjusted or even saved in different grinding methods. After grinding the glass plate to the precise thickness, the position of the sample holder is set to zero. Thereafter the thickness of the thin section is entered directly.

Target grinding in components

When the precise position of a component within a sample is known, Accutom-50 can be programmed to grind down to precisely that position. The position can be adjusted in steps of 5 µm.

Variable cut-off wheel speed up to 5000 rpm

The variable cut-off wheel speed can be adjusted from 300 up to 5000 rpm in steps of 100 rpm. This grants even higher precision as the drift of the cut-off wheel is reduced because of the higher speed. Also, the feed speed can be increased without experiencing any deterioration in the high quality of the cut surface.

Database for both cutting and grinding methods

Accutom-50 is equipped with a database capable of storing 20 cutting and 10 grinding methods.

Improved menu control

Subsequent menus for the selection of cutting or grinding parameters greatly facilitate the operation. Together with the larger display the new menu control improves the user-friendliness of Accutom-50.

Cut-off wheels

As Accutom-50 works with programmable constant feed speeds, all types of cut-off wheels, including abrasive wheels, can be used.


Struers offer a large variety of Diamond, CBN, Al₂O₃ and SiC cut-off wheels covering all possible applications. Also, saw blades for the cutting of very soft metals or plastics can be used on Accutom-5/50 due to the cover with the safety switch.

Cup wheels

Two different types of cup wheels, resin and metal bonded, are available for grinding on Accutom-50.

Specimen holders

All specimen holders for Accutom-5/50 are provided with a dovetail plate. This means that the samples can be clamped into the sample holder outside the cutting chamber, and then easily inserted into the dovetail fixture of Accutom-5/50.

A wide range of specimen holders can be used on Accutom machines

Technical data

Cutting	Accutom-5	Accutom-50
	Cutting speed: 300-3000 rpm (adjustable in steps of 100 rpm) Feed speed: 0.005-3.000 mm/s (adjustable in steps of 0.005 mm/s) Max.pos.speed: Y=13 mm/s X=10 mm/s Force limits: low: approx. 20 N medium: approx. 40 N; high: approx. 60 N	Cutting speed: 300-5000 rpm (adjustable in steps of 100 rpm) Feed speed: 0.005-3.000 mm/s (adjustable in steps of 0.005 mm/s) Max.pos.speed: Y=13 mm/s X=10 mm/s Force limits: low: approx. 20 N medium: approx. 40 N; high: approx. 60 N
Pos. range	Y direction: X direction:	105 mm (precision 0.1 mm), 60 mm (precision 0.005 mm)
Sample size	Max length of cut-off sample: Max length of sample to be clamped: Max sample cross section: ø127 mm cut-off wheel and ø42 mm flange: ø152 mm cut-off wheel and ø42 mm flange:	30 mm, 140 mm at ø20 mm 225 mm ø40 mm (without rotation) ø80 mm (with rotation) ø50 mm (without rotation) ø100 mm (with rotation)
Recirculation cooling unit	Contents: 4 l Flow: 800 ml/min	Contents: 4 l Flow: 800 ml/min
Motor	370 W at 3000 rpm	550 W at 5000 rpm
Software and electronics	Display: 4 x 40 characters Controls: touch pad Database: 10 cutting methods	Display: 4 x 40 characters Controls: touch pad Database: 20 cutting and 10 grinding methods
Power supply	1-3 x 200-240 V / 50/60 Hz	
Dimensions and weight	Width: 510 mm / 20", Depth: 700 mm / 27.6", Height: 270 mm / 10.6", Weight: 45 kg / 99 lbs	

Specifications

Specifications	Cat. no.
Accutom-5 Precision cut-off machine with variable speed. Automatic feeding with electronic control of feed speed. Motorised positioning system and digital read-out, with built-in sample rotation and oscillation. Complete with recirculation cooling unit, flange sets 42 mm dia. and specimen holder (04946903). Cut-off wheels, additive for cooling fluid and additional specimen holders are ordered separately.	04946227
Accutom-5 Precision cut-off machine with variable speed. Automatic feeding with electronic control of feed speed. Motorised positioning system and digital read-out. Complete with recirculation cooling unit, flange sets 42 mm dia. and specimen holder (04946903). Without built-in sample rotation and oscillation. Cut-off wheels, additive for cooling fluid and additional specimen holders are ordered separately.	04946127
Accutom-50 Precision cut-off and grinding machine with variable speed. Equipped with Multi-cut feature. Automatic feeding with electronic control of feed speed. Motorised positioning system and digital read-out, with built-in sample rotation and oscillation. Complete with recirculation cooling unit, flange sets 42 mm dia. and specimen holder (04946903). Cut-off wheels, additive for cooling fluid, cup wheels, additional specimen holders, and the flange set for cup-wheels are ordered separately.	05016227
Accessories	
Specimen holders	
Three jaw-chuck	05036919
For general use. Vice type with max. opening 60 mm	04946903
For round or square specimens. Teardrop type with max. opening ø 40 mm / ø 1½"	04946904
For irregular specimens, with 7 screws. Max. width 40 mm / 1½"	04946905
With goniometer	04276911
For adhering specimens	04276912
With ceramic vacuum chuck for thin sections	04276913
For small specimens. Vice type	04276915
For long specimens. With double parallel vice	04946909
For cylindrical specimens, max. dia. 40 mm	04946907
Joints to be mounted between the dovetail and the specimen holder	
Tilting joint with max. angle ±10°	04946906
Angling joint with max. angle +30/-90°	04946908
Base Plate	
With dovetail. For mounting of other types of specimen holders	04276914
Flange Sets	
For cut-off wheels, 65 mm dia., to achieve highest precision	04946902
For cup wheels, only to be used on Accutom-50	05016901
Additive for Cooling Fluid	
For Recirculation Cooling Unit. To protect the machine from corrosion, and to improve cutting and cooling qualities	1 l 49900045 5 l 49900046
See special brochure for more information about cut-off wheels.	

Struers' equipment is in conformity with the provisions of the applicable International Directives and their appurtenant Standards. (Please contact your local supplier for details)

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice.

Struers A/S
 Pederstrupvej 84
 DK-2750 Ballerup, Denmark
 Phone +45 44 600 800
 Fax +45 44 600 801
 struers@struers.dk
 www.struers.com

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CANADA

Struers Ltd.
 7275 West Credit Avenue
 Mississauga, Ontario L5N 5M9
 Phone +1 905-814-8855
 Fax +1 905-814-1440
 info@struers.com

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Struers Inc.