Discotom-6

Instruction Manual

Manual No.: 15677001

Date of Release 09.09.2013

Discotom-6
Instruction Manual

Table of Contents	Page
User's Guide	1
Reference Guide.....	18
Quick Reference	31

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

Original Instructions. The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2013.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Discotom-6 Safety Precaution Sheet

To be read carefully before use

- 1.** The operator(s) should be fully instructed in the use of the machine and its cut-off wheels according to the Instruction Manual and the instructions for the cut-off wheels.
- 2.** The machine must be placed on a safe and stable support table. All safety functions and guards of the machine must be in working order.
- 3.** Use only intact cut-off wheels. The cut-off wheels must be approved for min. 50 m/s.
- 4.** Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
- 5.** The workpiece must be securely fixed in the quick-clamping device or similar. Large or sharp workpieces must be handled in a safe way.
- 6.** The cutting handle should be lowered slowly and carefully, in order to avoid breaking the cut-off wheel.
- 7.** Struers recommends the use of an exhaust system as the materials being cut may emit harmful gasses or dust.
- 8.** The rubber curtains on the left hand side make it possible to cut long workpieces with a height of up to 80mm.
- 9.** Do not put your hand through the rubber curtain during operation of the machine.
- 10.** Use of safety goggles is recommended when using the flushing gun.

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.
The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.)

User's Guide

Table of Contents	Page
1. Getting Started	
Checking the Contents	2
Getting Acquainted with Discotom-6	3
Discotom-6 Rear View	4
Noise Level	4
Power Supply	4
Direction of the Cut-off Wheel.....	4
Connecting a Recirculation Cooling Unit.....	5
Connection to an External Exhaust System.....	6
2. Basic Operations	
Using the Controls.....	7
Front Panel Controls of Discotom-6.....	7
Front Panel Controls	8
Cooling Valve Positions	9
Cutting Position.....	9
Cleaning Position.....	9
Light Diode and Beep Signals.....	10
Fitting or Changing the Cut-off Wheel.....	11
Clamping the Workpiece	11
Starting/Stopping the Cutting Process	12
Automatic Cutting (AUTO mode)	12
Manual Cutting.....	14
Combining Manual and Auto Operation	15
Cleaning	15
Flushing the Cutting Chamber	15
3. Routine Maintenance	
Daily Service	16
Weekly Service	16
Monthly service	16
Cleaning the Flushing Gun	16
Replacing the Cooling Water	16
Yearly Service	17
Inspection of Guard.....	17

1. Getting Started

Checking the Contents

In the packing box you should find the following parts:

- 1 Fork spanner (24mm), for changing the cut-off wheel
- 1 Connector pipe for water outlet
- 3 Elbow pipes for water outlet
- 1 Outlet hose, 2 m
- 3 Hose clamps
- 1 Grease for Y-table spindle
- 1 Oil for maintenance of cutting table
- 1 Set of Instruction Manuals

Unpacking Discotom-6

Remove the screws underneath the pallet to free the Discotom-6.


Placing Discotom-6

Discotom-6 should be placed on a table that is able to carry a weight of min. 140 kg.

Discotom-6
Instruction Manual

**Getting Acquainted
with Discotom-6**


Take a moment to familiarise yourself with the location and names of the Discotom-6 components.

(Right side view)


- ① Cutting handle
- ② Protection guard, with window
- ③ Instrument panel
- ④ Emergency Stop
- ⑤ Flushing Gun
- ⑥ Cooling valve
- ⑦ Cutting table

- ⑧ Main power switch
- ⑨ Hours counter
- ⑩ Cut-off wheel locking handle

*Discotom-6
Instruction Manual*

Discotom-6 Rear View

- ① Upper back plate
- ② Cooli Unit connection
- ③ Electrical cable for power supply
- ④ Cut-off wheel locking handle
- ⑤ Exhaust hose
- ⑥ Water outlet to Flushing Gun
- ⑦ Flushing Gun
- ⑧ Drainage hose
- ⑨ Water outlet
- ⑩ Water inlet

Noise Level

Approx. 67 dB (A) measured at idle running, at a distance of 1.0 m/39.4" from the machine.

Power Supply

- First check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.
- Discotom-6 is factory mounted with an electric cable. Mount a plug on the cable:
 - Black and brown = phase
 - Yellow/green = earth

Direction of the Cut-off Wheel

Check that when the power is turned on, the cut-off wheel rotates in the direction indicated by the moulded arrow in the guard of Discotom-6. If the direction of rotation is incorrect, switch the two phases.

Connecting a Recirculation Cooling Unit

To ensure optimal cooling, Discotom-6 can be fitted with a Cooli unit. Cooling System 5 is a Cooli configuration designed for use with Discotom.

Note:



Cooling System 5 includes a static filter.
For intensive use, and for materials generating a lot of swarf, a bandfilter such as Coolimat-200 is recommended.

Note

Before connecting the cooling unit to the Discotom, follow the instructions in the Cooling Units Instruction Manual to prepare it for use.

A

B

To connect the Discotom-6 to a Recirculation Cooling Unit:

- Plug the Cooli control unit's communication cable into the Discotom's control socket.
- Connect the water inlet hose to the Cooli pump using the quick coupling (A).
- Connect the other end of the hose to the quick coupling of the Discotom water inlet.

- Insert the drain pipe ① in the water outlet on the back of Discotom and mount the 90° elbow pipe ②. Lubricate the sealing ring with grease or soap to facilitate insertion. (Use the other elbow pipe if more suitable).
- Strip the steel spring from approx. 3 cm of the outlet hose and cut. Bend the cut end towards the centre of the hose. Mount the outlet hose ③ onto the elbow pipe and clamp the stripped section using a hose clamp.
- Check that the outlet hose slopes downwards when connected. If necessary adjust the length of the hose.
- Insert the open end of the hose into mounting hole in the bracket on top of the Cooli filter unit (B).
- Connect the cooling unit to the mains power supply.

IMPORTANT

Before connecting, check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.

Drainage Hose

- Place the drainage hose in the tank of the Cooling Unit or lead to drain.

**Connection to an External
Exhaust System**


Struers recommends the use of an exhaust as workpieces may emit harmful gasses or dust when cut. The exhaust system will also reduce the level of water condensation on the sides of the protection guard.

- On the back of Discotom-6 you will find a joint for a 50 mm (approx. 2") dia. exhaust hose.
- Mount an exhaust hose from your local exhaust system onto the joint.

2. Basic Operations

Using the Controls








*Front Panel Controls
of Discotom-6*

Control name	Symbol
① Forwards button.....	▲
② Feed Speed Reduced indicator light.....	△
③ Backwards button	▼
④ Feed Speed potentiometer	
⑤ Load bar	
⑥ Emergency stop/overload indicator light	▽
⑦ Auto mode indicator light	Auto
⑧ Power indicator light	Power
⑨ Auto key	Auto
⑩ Start key	◊
⑪ Stop key	▽
⑫ Emergency stop button.....	

Discotom-6
Instruction Manual

Front Panel Controls


Name	Key	Function	Name	Key	Function
MAIN SWITCH		The main power switch is located on the right hand side of Discotom-6. Turn clockwise to switch on the power.	① FORWARDS		Moves cutting table forward.
⑧ POWER INDICATOR		Lights when main power is switched on.	③ BACKWARDS		Moves cutting table backwards.
⑩ START		Starts the machine. The cut-off wheel starts rotating and the cooling water is turned on.	④ FEED SPEED		Regulates the feed speed of the cutting table. The speed can be set from 0.1-2.5 mm/sec.
⑪ STOP		Stops the machine. The cut-off wheel stops rotating and the cooling water is turned off.	② FEED SPEED REDUCED		Lights when the feed speed has been set too high and is reduced by the machine.
⑨ AUTO		Select or deselect automatic cutting mode	⑤ LOAD BAR		Reflects the present motor load status.
⑦ AUTO INDICATOR		Lights when automatic cutting mode selected.	⑫ EMERGENCY STOP		Push the red button to stop. Pull the red button to release.
			⑥ EMERG. STOP OVERLOAD		Lights when emergency stop has been pressed or the cutting motor is overloaded.

*Discotom-6
Instruction Manual*

Cooling Valve Positions


The flushing gun, and a valve for cooling and cleaning, are located on the left hand side of the cutting machine. During cutting the valve is turned clockwise (mark at the top) and the water is directed to cool the cut-off wheel. For cleaning the valve is turned counter-clockwise (mark at the bottom) to flush the cutting chamber.

Cutting Position

Turn the knob clockwise to vertical position to cool the cut-off wheel.

Cleaning Position

Turn the knob counter-clockwise to vertical position for flushing the cutting chamber.

Light Diode and Beep Signals

Discotom-6 is provided with various light diodes indicating the status of the machine. Beep signals will sound if a key is pressed when the function cannot be activated.

Status	Indication	Comments
Green POWER indicator on	The main power switch has been turned on.	
White AUTO indicator on	AUTO mode is active.	
🔔 Beep signal	You pressed AUTO, START ⌂, FORWARDS ▲ or BACKWARDS ▼ while the motor was running.	
2 long 🔔 Beep signals	You pressed START ⌂ while the protection guard was open.	
3 long 🔔 Beep signals	You pressed START ⌂ for cutting while cooling valve isn't in cutting position.	Turn cooling valve to cutting position.
Red ⚡ indicator on	The emergency stop has been pressed.	Release the emergency stop to shut off the ⚡ indicator.
Red ⚡ indicator on and 🔔 beep signal	The emergency stop has been pressed and you tried to press another key.	
Yellow △ FEED SPEED REDUCED indicator blinking	In AUTO mode: the feed speed is too high.	The machine automatically turns the feed speed down.
	An AUTO cut has finished. Feed speed reduced.	Turn the feed speed down or press any key.
Red ⚡ and yellow FEED SPEED REDUCED indicators blinking	Cutting in AUTO mode not possible as the cutting table does not move.	Press any key to stop the indicators. Material may be too hard for the cut-off wheel chosen or cutting table may be jammed. See Trouble Shooting section at back of manual.
Red ⚡ indicator and 1 red cell on LOAD BAR	Pump motor overloaded.	Stop the machine. Push the red button on the back of the machine to reset the safety cut-out.
Red ⚡ indicator and 2 red cells on LOAD BAR	Main motor overloaded. Process halted.	The indicators will disappear when the motor has cooled down.
Red ⚡ indicator, 2 red cells on LOAD BAR and 🔔 beep signal		You tried to press a key before the motor had cooled down.

Fitting or Changing the Cut-off Wheel

Note: The spindle for Discotom-6 is left-hand threaded. If Discotom-6 has been used for manual cutting, fasten the cutting handle with the black handle at the right hand side of the machine before changing the cut-off wheel.

- Press and hold down the black locking knob on the right-hand side of the cut-off wheel, turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner. Remove the flange and the old cut-off wheel.
- Mount the new cut-off wheel.

IMPORTANT

Conventional cut-off wheels based on Al₂O₃/SiC abrasives should be placed between two cardboard washers, to protect the cut-off wheel.
For maximum precision with diamond or CBN cut-off wheels, do not use cardboard washers.

- Mount the flange and the nut. Tighten carefully and release the locking knob.

Clamping the Workpiece

- Place the workpiece between the quick clamping devices and the back stops.
- Push the clamps towards the workpiece and lock the quick-clamping devices with the locking handle.

IMPORTANT

The workpiece must be firmly held in the quick-clamping devices.

- Carefully close the protection guard.

**Starting/Stopping
the Cutting Process**
*Automatic Cutting
(AUTO mode)*

When the machine is used in AUTO mode the cut-off wheel is stationary and the cutting table moves. In MANUAL mode the cutting table is stationary and the operator moves the cut-off wheel.

- Clamp the workpiece.
- Loosen the black handle situated on the right hand side of the cutting chamber.
- Lower the cut-off wheel by drawing the cutting handle downwards until the cut-off wheel is positioned ready to cut the workpiece. Fasten the black handle again.
- Press Δ FORWARDS to position the cutting table and the workpiece just in front of the cut-off wheel.
- Close the protection guard.
- Set the desired feed speed
- Press the AUTO mode key to select AUTO mode.
- Press START \diamond to start the machine. The cut-off wheel starts rotating and the cooling water is turned on. The cutting table moves forward toward the cut-off wheel at the pre-set feed speed.
- If the machine has been stopped with the STOP \ominus key, the cutting table will stay in position. To return the cutting table to its start position press BACKWARDS ∇ once. Struers recommends manually lifting the cut-off wheel away before the cutting table is moved back. Please note that the protection guard has to be closed in order to carry out this operation. If the protection guard has been opened press and hold the BACKWARDS key ∇ to move the cutting table to its start position.
- If the cutting table reaches the back of the cutting chamber during the cutting process, the cut-off wheel automatically stops and the cutting table returns to its start position.

- To set another stop position:
 - Press FORWARDS \blacktriangle to move the cutting table to the desired position.
 - Press and hold the STOP $\bigcirc\!\!\!/\!\!\!\bigcirc$ key and then press once on FORWARDS \blacktriangle to set this position as stop position Then release the STOP $\bigcirc\!\!\!/\!\!\!\bigcirc$ key.
When the cutting table now reaches this position during cutting, the cut-off wheel will stop and the cutting table will return to its start position.
The set stop position will be cancelled when the power is switched off.

Note: the cutting table can be positioned with FORWARDS \blacktriangle and BACKWARDS \blacktriangledown . When the keys are held down, the cutting table moves with increased speed of 20 mm / sec.

Reduced Feed Speed

If the feed speed is set too high, the yellow FEED SPEED REDUCED indicator (\triangle) above the feed speed control flashes and the feed speed is automatically reduced by 10%. Discotom-6 can repeat this reduction up to 5 times.
Once the workpiece has been cut, the operator can reduce the feed speed using the feed speed control until the overload indicator stops flashing. Discotom-6 will then be able to cut another similar workpiece without feed speed reduction.

Manual Cutting

- Clamp the workpiece.
- Loosen the small black handle on the right side of the cutting chamber. The cutting handle is now easily moved up and down.
- Position the cutting table and the workpiece under the cut-off wheel with the **FORWARDS ▲** and **BACKWARDS ▼** keys
- Close the protection guard.
- Make sure that AUTO MODE is de-selected.
- Press START ◇ to start the machine. The cut-off wheel starts rotating and the cooling water is turned on.
- Pull down the cutting handle and let the cut-off wheel work itself into the workpiece. Increase the pressure and begin cutting.
- When the cut-off wheel has almost cut through the workpiece reduce the pressure.
- When the cut-off wheel has cut through the workpiece push the cutting handle back to its top position.
- Press the ☓ STOP key to stop the machine

Remember

When cutting manually, the motor load bar should be used to monitor the force on the workpiece. Try to keep the motor load in the green or the beginning of the yellow area.

Combining Manual and Auto Operation

Manual and Automatic cutting mode may be used in combination. You may start in manual mode and cut into the workpiece. If you lock this position you can continue in automatic mode, and the workpiece will move towards the cut-off wheel.

Cleaning

Flushing the Cutting Chamber

To ensure a longer lifetime for your Discotom-6, Struers strongly recommends that you clean the cutting chamber with the flushing gun every day.

- Press STOP  to stop the cut-off wheel and cooling water.
- Open the protection guard and release the workpiece from the quick-clamping devices.
- Take the flushing gun from the left hand side of the cutting machine and point it towards the bottom of the cutting chamber.
- Turn the cooling valve to cleaning position (see sketch at left)
- Clean the cutting chamber thoroughly.


Note:

Do not clean the protection guard with the flushing gun;
use a damp cloth instead.

- Turn the cooling valve to cooling position.
- Leave the protection guard open to let the cutting chamber dry completely.

AVOID RUST!

Leave the protection guard open to let the cutting chamber dry completely.

3. Routine Maintenance

Accumulated dirt and swarf can restrict or cause damage to the movement of the cutting table. To ensure a longer lifetime for your Discotom Struers strongly recommends daily cleaning of the cutting chamber.

Clean the cutting chamber thoroughly if the Discotom is not to be used for a longer period of time.

Recirculation Cooling Unit

For Maintenance of the Recirculation Cooling Unit please refer to the Cooli Instruction Manual.

Daily Service

- Clean the cutting chamber, especially the cutting table with the T-slots.
- Clean the protection guard window with a damp cloth. Do not use the flushing gun.
- Do **not** clean the lamp glass with alcohol. Use a damp cloth.

Weekly Service

- Clean the cutting chamber thoroughly:
Move the cutting table forwards and backwards to access the whole of the cutting chamber.
- Clean along the length of the guide shafts with the flushing gun and a brush to remove accumulated swarf.
- Clean under the cutting table with the flushing gun and a brush to remove accumulated swarf.
- Check the level of the cooling water after 8 hours use or at least every week.

Monthly service

Cleaning the Flushing Gun

- Remove flushing gun from the hose by squeezing the clamp.
- Soak in warm soapy water to loosen dirt and residue.
- Press the handle and clean the barrel with compressed air or a pipe cleaner. (The nozzle can be removed for better access).

Replacing the Cooling Water

- Replace the cooling water in the Recirculation Cooling Unit at least once a month.

Discotom-6
Instruction Manual

Yearly Service
Inspection of Guard

- Visually inspect the guard and the glass for signs of wear or damage.


Important

Carry out inspection at more regular intervals if Discotom-6 is used for more than one 7 hour shift a day.

Struers recommends that the PETG glass is replaced after 5 years of routine operation.

The guard should be replaced immediately if it has been weakened by collision with projectile objects or if there are visible signs of deterioration as a result of using a cooling fluid other than those produced by Struers.

A label on the guard indicates when the glass is due to be replaced.

Reference Guide

Table of Contents	Page
1. Advanced Operations	
Adjusting the Back Stop	19
Clamping Irregular Workpieces	19
Mounting an Adjustable Cutting Table	19
Optimising the Cutting Results	20
Explanation of Safety Factors	20
2. Accessories	21
3. Consumables	21
4. Trouble-Shooting	22
5. Maintenance	
Hours Counter	25
Maintenance of Cutting Tables	26
Lubricating the Cutting Table	26
Maintenance of Cut-off Wheels	28
Storing Bakelite Bonded Al ₂ O ₃ Cut-off Wheels	28
Maintenance of Diamond and CBN Cut-off Wheels	28
Maintenance of Clamping Devices	28
6. Technical Data	29

1. Advanced Operations

Adjusting the Back Stop

- Use a 13 mm spanner to loosen the back stop.
- Place the workpiece. Adjust the back stop so that the centre of the workpiece is under the centre of the cut-off wheel:
- Tighten the back stop.
- If you also need to adjust the clamping devices, loosen them and position them parallel to the back stop. Tighten the clamping devices.

Clamping Irregular Workpieces

Irregular workpieces without plane clamping surfaces must be clamped using special clamping tools, as the workpieces must not move during cutting. This could result in damage to the cut-off wheel or to the workpiece itself. For this purpose a vertical clamping device with one flat clamping shoe is available. (See Accessories).

Mounting an Adjustable Cutting Table

If you need to cut plane parallel sections, an adjustable table may be mounted on top of the existing cutting table. The adjustable table positions the workpiece at a right angle to the cut-off wheel and is movable up to 60mm. (See Accessories).

- Dismount the clamping tools, if mounted.
- Move the cutting handle to the top position and the cutting table to the back in order to create working space.
- Move the centre part of the adjustable cutting table to the left by turning the handle until the holes for the Allen screws are visible.
- Loosely mount two T-nuts underneath the adjustable cutting table using Allen screws.
- Slide the adjustable cutting table onto the cutting table:
 - Position the adjustable cutting table on the left-hand side of the cutting table with the spindle on the left.
 - Open the hinged protection guard to let the spindle pass.
 - The guide pins and T-nuts on the bottom of the adjustable cutting table are slid into the third T-slot from the left; this automatically aligns the adjustable cutting table.
- Tighten the Allen screws for the T-nuts.
- Mount clamping tools using the T-slots in the adjustable cutting table.

Optimising the Cutting Results

The following table shows answers to a number of common questions about how to achieve better quality cuts:

Question	Answer
How can I avoid discolouration or burning of the sample?	Reduce the feed speed Change the cut-off wheel, as the hardness of the present cut-off wheel may be inappropriate for the hardness of the workpiece*
How can I avoid burrs?	Use a softer cut-off wheel* Clamp the workpiece securely with the quick-clamping devices.
How can I prevent the cut-off wheel wearing out too quickly?	Reduce the feed speed or use a harder cut-off wheel*
How can I cut more quickly?	Place the workpiece in a position which allows the cut-off wheel to cut the smallest possible cross-section. Increase the feed speed.

- *) Please refer to the Selection Guide in the Struers [*Cut-off Wheels brochure*](#).

Explanation of Safety Factors

The protection guard has a safety switch to prevent the cut-off wheel from starting while the guard is open. Furthermore, a locking mechanism prevents the operator from opening the protection guard until the cut-off wheel stops spinning.

2. Accessories

Please refer to the *Discotom-6 brochure* for details of the range available.

Clamping Tools

Please refer to the *Struers Clamping Tools brochure* for details of the range available.

Other Accessories

Specification	Cat. No:
<i>T-slot cleaner</i> T-slot cleaner for 10 and 12 mm T-slots	05486910

3. Consumables

Cut-off Wheels

Please refer to the Selection Guide in the *Struers Cut-off Wheels brochure*.

Other Consumables

Specification	Cat. No:
<i>Corrozip</i> Environment friendly. To protect the machine from corrosion and to improve cutting and cooling qualities. For Recirculation Cooling Unit. 1 l 5 l	449900045 449900046
Grease for maintenance/lubrication of the spindle (GadusS2) 100 ml	16080802
Oil for maintenance of cutting table (Shell Tellus S100) 50 ml	16080845

4. Trouble-Shooting

Recirculation Cooling Unit

For Messages and Errors relating to the Recirculation Cooling Unit, please refer to the Instruction Manual.

Error	Explanation	Action
Indications on the Front Panel		
FEED SPEED REDUCED  and  indicator lights come on.	Emergency stop active.	Solve the problem causing the emergency stop, for example, by reducing the feed speed, or using a different cut-off wheel, then release the emergency stop.
	Cutting motor overloaded due to heavy use.	Open the protection guard and let the cutting motor cool down for about 5 - 10 min.
Machine Problems		
The cut-off wheel does not rotate.	The protection guard may be open.	Press the safety guard completely down. If this does not help, call a Struers service technician.
Chamber light does not work.	Lamp needs to be replaced.	Remove the plastic cap to gain access to the fluorescent lamp. Pull out the lamp and replace it.
Water leaking.	Leak in a Recirculation Cooling Unit hose.	Check the hose and tighten the hose clamp.
	Water overflow in the cooling water tank.	Remove the excess water in the tank.
Workpieces or cutting chamber rusty.	Insufficient additive in cooling fluid.	Use Struers Additive for cooling fluid in the cooling water, at the correct concentration. Check with a refractometer. Follow the instructions in the 'Maintenance' Section.
	The protection guard is left closed after use.	Leave the protection guard open to let the cutting chamber dry.
Protection guard window blurred.	Insufficient cleaning.	Clean with mild soapy water (do not use the flushing gun).
Quick-clamping device unable to hold the workpiece.	The quick-clamping device is not balanced.	Adjust the screw underneath the clamping column. Use a 3 mm Allen key.
	Clamp worn.	Call a Struers service technician.
Flushing gun blocked	Build up of dirt in the gun.	Clean flushing gun. Follow the instructions in the 'Maintenance' Section.

Discotom-6
Instruction Manual

Error	Explanation	Action
Cutting Problems		
Discoloration or burning of the workpiece.	The hardness of the cut-off wheel is inappropriate for the hardness / dimensions of the workpiece.	See Consumables section, under heading 'Cut-off Wheels'.
	Inadequate cooling.	<ul style="list-style-type: none"> - Check that there is sufficient water in the Recirculation cooling unit. - Check the cooling water hoses.
Unwanted burrs	Wheel is too hard.	See Consumables section, 'Cut-off Wheels'.
	Feed speed too high at the end of the operation.	Reduce the feed speed near the end of the operation.
	Lack of support for the workpiece.	If possible support the workpiece on both sides.
Cutting quality differs.	Cooling water hose clogged.	Clean the cooling water hose and the cooling tube. Check the water flow by turning the cooling valve to cleaning position.
Cut-off wheel breaks.	Incorrect mounting of the cut-off wheel.	<ul style="list-style-type: none"> - Check that the hole has the correct diameter. - Check cardboard washer on both sides of the cut-off wheels and replace if worn. The nut must be tightened properly.
	Incorrect clamping of the workpiece.	Make sure that only the left quick-clamping device is tight. The other device should only press lightly.
	Wheel is too hard.	See Consumables section, 'Cut-off Wheels'.
	Feed speed is set too high.	Reduce the feed speed.
	Inadequate cooling.	Check that there is enough water in the cooling unit. Check the cooling water hoses.

Discotom-6
Instruction Manual

Error	Explanation	Action
The cut-off wheel wears down too quickly.	The feed speed is too high.	Reduce the feed speed.
	Insufficient cooling.	- Check that there is enough water in the Recirculation cooling unit. - Check the cooling water hoses.
	The cut-off wheel is too soft for the task.	See Consumables section, 'Cut-off Wheels'.
The cut-off wheel does not cut through the workpiece.	Incorrect choice of cut-off wheel.	See Consumables section, 'Cut-off Wheels'.
	Cut-off wheel worn.	Replace the cut-off wheel.
	The cut-off wheel gets caught in the workpiece.	Support the workpiece and clamp it on both sides of the cut-off wheel so that the cut stays open.
The workpiece breaks when clamped.	The workpiece is brittle.	Place the workpiece between two polystyrene plates. Note! Always cut brittle workpieces very carefully.
The sample is corroded	The sample has been left in the cutting chamber for too long.	Remove the sample directly after cutting. Leave the protection guard open when you leave the machine.
	Insufficient additive for cooling fluid.	Use Struers Additive for cooling fluid in the cooling water at the correct concentration. Check with a refractometer. See Maintenance section.

5. Maintenance

Service and greasing of the quick-clamping devices is part of a regular Struers Service Schedule. Struers recommends that a regular service check be carried out after every 1500 hours of use. Please contact a Struers Service Technician to service the machine.

Hours Counter

The hours counter at the right hand side of Discotom-6 records the number hours that the cutting motor has been in operation. To ensure trouble-free operation of Discotom-6, service provided by a Struers Service Technician is recommended for every 1500 hours of operation. The hours counter can only be reset by a Struers Service Technician.

Maintenance of Cutting Tables

The stainless steel bands should be replaced if they become worn or damaged. The bands are available as spare parts.

To allow humidity to escape from the cutting table and chamber, it is recommended to leave the hood open when the machine is not in use.

Lubricating the Cutting Table


To maintain the optimum performance of Discotom-6, lubricate the cutting table at regular intervals (approx every 50 hours).

Check the service counter on the right hand side of the machine regularly to monitor the actual number of usage hours.

After lubrication, note the date and number of service hours on the Maintenance Log table.

Spindle

- Move the cutting table to its rear position using the ▼ BACKWARDS key.


- Place the grease gun on the nipple in front of the spindle for the cutting table and push twice to grease the spindle.

A grease gun with grease for lubricating the Y-table spindle is supplied with the Discotom-6.

When all the grease has been used, refill the grease gun.

Slide Bearings

- Remove the stainless steel bands located over the guide shafts (see illustration).

- Fill the holes with oil eg. Shell TELLUS oil S100.
Use approx. 12.5ml oil for each side or until the reservoirs are full.
A bottle with oil for lubricating the slide bearings is supplied with the Discotom-6. When all the oil has been used, refill with eg. Shell TELLUS oil S100.

Important

If the oil in the reservoir appears milky;
this indicates that water has entered the oil reservoir.
Call a Struers Service Technician to clean the reservoir.

- Check that the O-rings are in position and remount the stainless steel bands.

Maintenance of Cut-off Wheels

Storing Bakelite Bonded Al₂O₃

Cut-off Wheels

*Maintenance of Diamond and
CBN Cut-off Wheels*

These cut-off wheels are sensitive to humidity. Therefore, do not mix new, dry cut-off wheels with used humid ones. Store the cut-off wheels in a dry place, horizontally on a plane support.

The precision of diamond and CBN cut-off wheels (and thus the cut) depends on how carefully the following instructions are observed:

- Never expose the cut-off wheel to a heavy mechanical load, or heat.
- Store the cut-off wheel in a dry place, horizontally on a plane support, preferably under light pressure.
- A clean and dry cut-off wheel does not corrode. Therefore, clean and dry the cut-off wheel before storing. If possible, use ordinary detergents for the cleaning.
- Regular dressing of the cut-off wheel is also part of the general maintenance.

Maintenance of Clamping Devices

Important

It is recommended to thoroughly clean and lubricate the Quick Clamping Devices at regular intervals.

6. Technical Data

Subject	Specifications		
	Metric/International	US	
Cutting motor	Cutting power constant [S1] Cutting power intermittent [S3] Maximum power	2.2 - 3.0 kW 3.2 - 4.0 kW 4.4 - 6.0 kW	3.0 - 4.0 HP 4.3 - 5.4 HP 6.0 - 8.0 HP
	Rotational speed (running idle)	2850 rpm at 50Hz 3420 rpm at 60 Hz	2850 rpm at 50Hz 3420 rpm at 60 Hz
Cut-off Wheel	Diameter x thickness	Ø250x1.5x32mm	9.8 x 0.06 x1.25"
	Max. cutting diameter: Max. length of cut: (at a workpiece height of 65mm and width of 205mm) The max. cutting capacity might be less, depending on the workpiece material and the cut-off wheel	80mm 140mm	3.1" 5.6"
Cutting Table Dimensions	Width Length	386mm 225mm	15.19" 8.85"
Fume Extraction	Dia. for connecting tube Recommended capacity: at 0mm /0" water gauge	50 mm 240 m ³ /h	2.0" 8475 ft ³ /h
Dimensions and Weight	Discotom-6 Width Depth Height Weight <i>Optional Table Unit</i> Width: Depth: Height:	794mm 640mm 458mm (123kg) 900mm 750mm 800mm	31.26" 25.20" 18.03" 272lbs 35.4" 29.5" 31.5"

Discotom-6
Instruction Manual

Subject	Specifications				
Electrical Data	Voltage/frequency: 3 x 20050Hz 3x200-210/60Hz 3 x 220-230V/50Hz 3 x 220-240V/60Hz 3 x 380-415V/50Hz 3x380-415/60Hz 3 x 460-480V/60Hz	Max. Load: 21.3A 21.3A 19.8A 19.8A 11.8A 11.8A 11.8A	Fuse		
			min.	max.	
<i>Residual Current Circuit Breaker</i>	type A, 30 mA (or higher) is recommended.				
Safety standards	Please refer to the Declaration of Conformity				
Environment	Noise level	Approx. 67dB(A) during cutting, at a distance of 1.0 m/39.4"			
	Dynamic Balancing	from the machine. Max. admissible unbalance U_{per} according to ISO 1940/1, Balance Quality Grade G6.03: $U_{per} = 2 \text{ gmm/kg} = 110 \text{ gmm}$.			

Quick Reference

Clamping the Workpiece

- Place the workpiece between the clamps and the back stop.
- Push the clamps towards the workpiece and lock the quick-clamping device with the locking handle.

Auto Cutting

- Clamp the workpiece.
- Position the cut-off wheel so that it can cut the sample and lock it in place with the black cut-off wheel locking handle on the right of the machine.
- Position the cutting table with \blacktriangle FORWARDS and \blacktriangledown BACKWARDS.
- Carefully close the protecting guard.
- Set feed speed and Press AUTO.
- Press START \diamond to begin cutting.
- After cutting, the cutting table automatically returns to its start position.
- Open the protection guard and release the clamping device.
- Remove the workpiece.

Manual Cutting

- Loosen the black handle on the right hand side of Discotom-6 so that the cut-off wheel moves freely
- Clamp the workpiece.
- Position the cutting table under the cut-off wheel.
- If AUTO mode is selected, de-select it.
- Press START \diamond . Lower the cut-off wheel until it contacts the workpiece and apply pressure until the sample is cut.
- Keep the Load Bar in the green, or first part of the yellow area. Relax pressure near the end of the cut.
- Open the protection guard and release the clamping device.
- Remove the workpiece.

Changing the Cut-off Wheel

- Fasten the black handle to stop the cutting handle from moving.
- Press the black locking knob on the right-hand side of the cut-off wheel, turning the cut-off wheel until the spindle lock clicks.
- Remove the nut with a fork spanner. Remove the flange and the old cut-off wheel.
- Mount the new cut-off wheel.
- Mount the flange and the nut. Tighten carefully.

Flushing the Cutting Chamber

- Pull out the flushing gun and point towards the bottom of the cutting chamber.
- Turn the cooling valve to flushing position.
- Clean the cutting chamber.
- Turn the cooling valve to cooling position.
- Always leave the protection guard open to let the cutting chamber dry.

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Date of Release €J.€J.20FH

Discotom-6
Spare Parts and Diagrams

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 201H

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Telefax: +45 44 600 801

Spare Parts and Diagrams

Table of contents	Drawing
Discotom-6, complete	15670001S
Base and back plates, complete	15670015T
Base plate, complete.....	15670010Q
Stepper motor and gear potentiometer	15670042G
Cutting spindle, complete	15670050K
Cutting wheel guard, complete.....	15110055B
Water valve, house	15110065H
Cover, complete	15670030G
Lamp, complete.....	15670012C
Cabinet, complete	15670020B
Front plate, complete	15670025D
Quick Clamping Tool for 10 mm T-slots, Left	15870040K
Quick Clamping Tool for 10 mm T-slots, Right	15870041L
Stopper for 10mm T-slots.....	15870045C

Diagrams	
Wiring diagram (3 pages)	15673450K
Circuit diagram (2 pages)	15673120K

Some of the drawings may contain position numbers
not used in connection with this manual.

Discotom-6
Spare Parts and Diagrams

The following is a list of the spare parts that may need replacement during the lifetime of the equipment.

To check the availability of other replacement parts, please contact your local Struers Service Technician. It may help identify the part by referral to its position number on the assembly drawings included in this manual.

Spare part list for Discotom-6

Drawing	Pos.	Spare Part	Cat no:
15670001		Discotom-6, complete	
	210	O-ring 6.30-2.40 72 NBR 872	2IO24008
	335	50Hz coil	2MO01614
		60Hz coil	2MO01623
	465	W1 net cable 50Hz.	15113501
		W1 net cable 60Hz.	15113500
	470	Red Push Button For Thermal Relay.	2SA00301
15670015		Base and back plates, complete	
	0020	Motor 3x200V / 50Hz. 3x200-240V / 60 Hz.	2ME51140
		Motor 3x220-415V / 50 Hz	2ME51141
		Motor 3x460-480V / 60 Hz.	2ME51142
	0040	Key A 8X7X28	2ZF10828
	0190	Counter nut M25	2TD20026
	0200	Flex hose PMA PIST-17S.30	2NU31200
	0280	Cardboard flange Ø33/85	50000028
	0290	Flange	14150030
	0300	Washer for spindle	14150012
	0420	Edge list, black	2IP28459
	0430	Sealing band /31018 black	2IP30001
	0440	Sealing profile	2IP00203
	0560	Flushing gun	15670007
	0800	Thermo relay 0.18-0.27A	2FS10021
		Thermo relay 0.4-0.6A	2FS10023
		Thermo relay 0.6-0.9A	2FS10024
	0810	0.25 AT fuse, glass tube	2FU13500

Discotom-6
Spare Parts and Diagrams

Spare part list for Discotom-6

Drawing	Pos.	Spare Part	Cat no:
15670010		Base plate, complete	
	20	Buffer AA-315 ø30x15 m. M8x10	2GS10315
	60	O-RING 50.00-4.00 NBR70	2IO40080
	390	Gas spring	15670046
	630	Transformer	2MT70711
	690	5.00A T Fuse glass tube	2FU14250
		Replacement Stainless Steel Bands	
		1 pc. 40 mm steel band	05676901
		Set of 4 pcs. 40 mm steel bands	05116907
15670042		Stepper motor and gear potentiometer	
	0180	Stepper motor no.: 6600R174	2MI06674
	0260	Buffer ø30x40	2GS03040
15670050		Cutting spindle, complete	
	10	Spindle shaft	15670110
	20	Ball bearing, 6007 2RS1 ø35/ø62	2BK00110
	30	Locking ring J62 DIN 472	2ZL20620
	40	O-ring 31.42-1.78 NBR70	2IO17829
	50	Seal ring	15110158
	60+70	Seal ring	15112927
	80	V-Ring V-40A	2IV10040
	90	Lock nut	15110152
	100	Cylinder pin, hardened. ST.6m6x14 D6325	2ZS03454
	110	Spindle house painted	15679003
	120	Spring	2GF51020
	140	Lock ring	2ZL10350
	150	Coupling	15670070
	160	Disc M10	2ZA20010
	170	ULS screw	2TR80812
	180	Key 8x7x20	2ZF10820
	185	Ring for sensor	15670411
	186	Screw MSP M5x12	2TI10512
	200	Namur MA516NF M5x16 sensor	2HQ00517
15110055		Cutting wheel guard, complete	
	20	Hose nipple 2601-12-3/8	2NF40088
	30 →	Locking pin, complete	15112909
	90	Grip M5	375MP085

Discotom-6
Spare Parts and Diagrams


Spare part list for Discotom-6

Drawing	Pos.	Spare Part	Cat no:
15110065		Water valve, house	
	20	Inductive sensor	2HQ00023
	30	O-ring 25.00-2.00	2IO20055
15670030		Cover, complete	
	20	Rubber curtain, sliced	15110665
	30	Rubber curtain	15110666
	60	Sealing profile EPMD (1M)	2IP21038
15670012		Lamp, complete	
	0010	Lamp cover	15119001
	0050	Bulb	2HG10827
	0060	2uF 250V 50HZ condensator	2CR12301
15670020		Cabinet, complete	
	90	Seal profile 0,76 meter	2IP00205
15670025		Front plate, complete	
	20	Front plate foil	15670400
	60	PCB Discotom-6 Contr. board	15673380
		Accessories	
		Pipe for drain	15110997
		Spanner for cutting wheel 24 mm	2GR00224
		Socket spanner 13 mm	2GR00313
		Grease gun (GadusS2 V22000)	16080802
		Lubrication oil (Shell Tellus S100)	16080845
		Transportation crate	1511EM00
		Wiring	
		A2 Brake circuit PCB	15670047
		Contactor	2KM10640
		Contactor K2-09A01 24V	2KM09024
		Pot. Meter 5K 0,5W	2RP50120
		W5 cable	15673505
		W9 cable	15673509
		w10 cable	15673510
		W18 lamp cable	15673518

Spare part list for Discotom-6

Drawing	Pos.	Spare Part	Cat no:
Vertical clamping system for 10 mm T-slots (TREVS)			
	1	Operating key for TREVS	2GR27021
	2	Spring for TREVS	2GR27023
	3	Flat clamping shoe	2GR26150


Vertical clamping system for 10 mm T-slots (TREVS)





Discotom-6
Spare Parts and Diagrams













Spare part list for Discotom-6

Drawing	Pos.	Spare Part	Cat no:
<i>15870040 TENLE</i>		Quick Clamping Tool for 10 mm T-slots, assembled	
<i>15870041 TENRI</i>	10	Handle for Quick Clamping Tool	381MP128
	20	Base Plate 10mm, Right (TENRI)	15872907
		Base Plate 10mm, Left (TENLE)	15872906
	40	House, machined	15872901
	30	Jaw 10mm	15870122
	50	Rod, hexagon	15870125
	60	Clamping ring	15870126
	90	Handle	15870129
	110	Excenter, Right (TENRI)	15870190
	110	Excenter, Left (TENLE)	15870191
	130	Spring, 22x3-20	15870161
	170	Lube Nipple, M6x1 DIN 71412A Zn	2GN90050
	180	Rod Seal, 25x18x4.35	2IT20060
	190	Cap nut M8 A2 DIN15870	2TA30080
	200	T-Nut M8_10 DIN 508 A2	2TF41810
	210	Pin Bolt M8x30 DIN 939, 8.8+Zn	2TH30840
		Quick Clamping tool Wear Kit	15872905
	60	Clamping ring	
	130	Spring, 22x3-20	
	180	Rod Seal, 25x18x4.35	
		Angle pipe wrench 289N 17mm	2GR00317
<i>15870045 TENLE +TENRI</i>		Stopper for 10mm T-slots, assembled	
	10	Stopper 10mm	15870130
	20	T-nut M8 10 DIN 508 A2	2TF41810
	30	Allen screw M8x70 A2 DIN 933	2TZ10870

- Pos 130 Omnitit 50M
- Pos 320 Omnitit 230L
- Pos 610, 620 Omnitit 50M
- Pos 260, 475 Omnitit 100M
- Pos 20 Shaft seal in shaft-end, and V-seal in the free-end to be removed.

22

10 of 10

11-30 Bushings to be lubricated by Chain Flexia 10W 40

125 125 12,5 ml of the oil to be added in both s

130 130 Retinax CS00 (thin layer on the spindle)

(270) 270 To be secured using Omnitil 50H, and 1





530 280 280 Tightening torque = 35NM. Omnitilt 50M

400 400 To be lubricated with Shell Rotella Grease E5

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


1 2 3 4 5 6 7 8

A
B
C
D
E
F

BQ3 + nut secured with Loctite 243.
Note the adjustment!

Montageinstruktion: 15677300

DETAIL B
SCALE 2 : 1

K	2013-08-16	Spindelmontage noter slættet - se altid 15677300.			JTV	2013-08-16	JTV
A	04-11-08				JTV		
Revision	Crea. date yyyy-mm-dd	Revision description			Draw. Init	Appr. date yyyy-mm-dd	Appr. Init
			Material:	Scale:	Format:	Tolerance: DS/ISO 2768- mK	
				1:2	A3		
		ID:	Description:			Weight :	3 g
15670050 Spindelhus samlet							Rev: K

Pedersensgade 84
DK-2750 Ballerup/Copenhagen
Danmark
Phone: +45 44 600 800
Fax: +45 44 600 804

A

B

C

D

E

F

A





B

C

D

E

F

DETAIL A
SCALE 1 : 1


- (20) Loctite 577
- (40) Albida grease
- (90) Loctite 243

B	2011.06.16	Redrawn, washer andnut added	SPE	2011.06.16	
A	1996.01.13		KLB	1997.08.26	FPG
Revision	Crea. date yyyy-mm-dd	Revision description	Draw. Init	Appr. date yyyy-mm-dd	Appr. Init
	Material:	Scale:	Format:	Tolerance: DS/ISO 2768- mK	
		1:2	A3	Weight : g	
	ID: Description:			Rev:	
				B	

A


SECTION D-D

NAMUR BQ1
Afkortes ved montage

B

D

C

10

20

NAMUR BQ2
Afkortes ved montage

D

D

E

60

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40

50


60

70

Loctite 577

30

High Vacuum

F

H 2012-02-21

A 16.01.1997

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BRY

BMJ

2012-02-21

28.09.1997

CJE

FPG

Rev

Draw. Init

Appr. date

yyyy-mm-dd

Appr. Init

Struers

 Pedestrupvej 84
 DK-2750 Ballerup
 Copenhagen
 Denmark
 Phone: +45 44600 800
 Fax: +45 44600 804

Material:

Scale:

1:1

Format:

A4

Tolerance: DS/ISO 2768 - mK

Weight :

g






ID:

Description:

15110065 Ventilhus, monteret

Rev:

H

SECTION A-A

- 60, 100 Loclite 409
- 120, 160 Omnitil 50M after adjustment
- 140 To be glued with grey Vetrosilicone 512.
The upper and lower joint to be filled with black
Elasticon 33N (2LL00332). Curing time 24h
- A) Joints to be sealed with grey Vetrosilicone 512.

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
















1 2 3 4 5 6 7 8

SECTION A-A

Loctite 648 between house and insert,
both parts must be degreased first!

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




Rev. E: FTH 16-05-2006
Ver. 1.0
- 350-415V / 60Hz variant added
Rev. F: FTH 26-03-2007
- Protection against short circuit added (X13-X14)
Rev. G: FTH 30-10-2007
- Thermo relay removed (K2-F11)
Rev. H: FTH 26-03-2008
- Short circuit protection added
Rev. I: FTH 15-05-2009
- Short circuit protection added
Rev. J: FTH 20-05-2009
- Values in table A1 & F11 corrected

STRUTER'S A/S	
Production line 54	
DK-2720 Ballerup	
PHN001-44 44 000 800	
DISCOTOM-6, Circuit diagram	- Safety related control circuits
A1	15673120
J	
FTH/FTH	
2	2

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Herewith declares that

Product Name:	Discotom-6
Type No.:	567
Machine Type:	Cut-off machine

is in conformity with the provisions of the following directives:

Safety of Machinery 2006/42/EC according to the following standard(s):
EN ISO 12100:2010, EN ISO 13850:2006, EN ISO 13849-1:2008/AC:2009, EN ISO 13849-2:2008,
EN60204-1:2006/AC2010

EMC-Directive 2004/108/EC according to the following standard(s):
EN 61000-6-4:2007/A1:2011, EN 61000-6-1:2007, EN 61000-6-3:2007+A1:2011.

Low Voltage Directive 2006/95/EC according to the following standard(s):
EN60204-1:2006/AC2010

RoHS 2011/65/EU according to the following standard(s):
EN 50581:2012.

Supplementary Information The equipment complies with the following standards:
NFPA70:2011, NFPA79:2007, FCC 47 CFR Part 15 Class A

The above has been declared according to the global method, module A

Date: 21.01.2013

Christian Skjold Heyde,
Vice President, R & D and Production, Struers A/S

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Phone +45 44 600 800
Fax +45 44 600 801
struers@struers.dk
www.struers.com

Bankers:
Nordea Bank Denmark A/S
Account no.: 2191-0126305060
VAT no.: 15683309

Letter of conformity: Safety hood/window for the following equipment

AbraPlan-20	PETG thickness: 8 mm
AbraPol-20	PETG thickness: 8 mm
Discotom-5/-6/-10/-60/-65/-100	PETG thickness: 4 mm
Labotom-3	PETG thickness: 4 mm
Labotom-5	PETG thickness: 3 mm
Labotom-15	PETG thickness: 5 mm
Axitom/-5	PETG thickness: 4 mm
Exotom-100/-150	PETG thickness: 4 mm
Magnutom-400/-500	PETG thickness: Front: 10mm, Side: 8mm

The safety hoods/window primary function is to protect against samples not clamped properly.

The safety hood is made from PETG, which is generally used in all Struers machines as "safety glass".

PETG is resistant to scratches and chemicals and provides the highest degree of safety for the operator.

Commodity Supplier: **RIAS A/S**, Industrivej 11, 7000 Roskilde, Denmark

Producer: **NEOPLAST ApS**, Tingbjergvej 4, 4632 Bjæverskov, Denmark

Supplier: **Struers A/S**, Pederstrupvej 84, 2750 Ballerup, Denmark

This is to certify that the product is produced according to regulations from the Commodity Supplier.

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Yours faithfully

Flemming Perret-Gentil
R&D Manager
R&D Department, Equipment

+45 44 600 913 (direct)
+45 27 130 890 (mobile)
flemming.perret.gentil@struers.dk
www.struers.com

Discotom-6

Versatile automatic
and manual
cut-off machine in one

- Automatic and manual cut-off machine in one
- Very high cutting capacity
- Constant feed speed principle ensures minimum deformation
- Automatic feed speed reduction for fast cutting of hard materials
- Quick positioning and freely selectable stop position functions
- Effective cooling system reduces the risk of burns
- Corrosion resistant cutting chamber
- Electric motor brake and safety lock ensure operator safety

Discotom-6 is a high performance table-top cut-off machine with a strong 3.7 kW (5.0 HP) motor. It is designed for cut-off wheels with a diameter of 250 mm (10").

Discotom-6 can be operated both automatically and manually, making it suitable for a wide range of applications.

Automatic Cutting made Easy

Automatic cutting mode is engaged by a simple press key. The cutting table with the clamped workpiece will move forward towards the cut-off wheel at the set feed speed. Before cutting, you can of course position the cut-off wheel to compensate for the wear of the cut-off wheel.

High Cutting Capacity

Manual mode for Fast and Uncomplicated Cutting

When operating the cut-off machine manually, the cutting table with the clamped workpiece is easily positioned under the cut-off wheel, and you pull the cutting handle to begin cutting.

Very High Cutting Capacity

You can cut workpieces with a diameter of up to 80 mm (3.1") on the Discotom-6. In automatic mode though, the cutting table is movable up to 200 mm (7.9") enabling cutting of very long samples. You can for instance cut a length of 140 mm (5.5") at a workpiece thickness of 65 mm (2.6").

The graph shows the cutting capacity in automatic cutting mode with a new cut-off wheel and flexible clamping tools. The real cutting capacity depends on the sample material

Quick positioning and Auto Stop

The cutting table can be moved with increased speed allowing quick positioning and set-up.


Additionally, a freely selectable stop position can be set where the cutting process will be stopped and the cutting table returns to its start position.

Automatic Feed Speed Reduction

When cutting automatically, you choose a feed speed from 0.1 to 2.5 mm/sec. If the feed speed is set too high and the machine cannot cut through the workpiece, Discotom-6 will automatically reduce the feed speed and continue cutting. This may repeat itself up to 5 times and prevents damage or burns on the sample surface.

Effective Cooling of the Cutting Area

A very important design detail is that the water cooling tubes are integrated into the cut-off wheel guard. The cooling water is applied to the cut-off wheel under high pressure and at a right angle. The cooling water thereby "breaks" the thin layer of air surrounding the spinning cut-off wheel, and the result is excellent cooling of the cutting area. Discotom-6 is available with a large 65 l recirculating cooling unit, ensuring that the cooling water does not overheat.

The adjustable cutting table enables parallel sectioning

Monitoring the Load on the Cutting Motor

The front panel includes a three-coloured LED-bar, indicating the load on the cutting motor.

In automatic cutting mode the LED-bar will show the load on the cutting motor and if the red area is reached the feed speed will be reduced. When cutting manually the LED-bar gives a good indication of how much force to put on the workpiece.

Corrosion resistant Cutting Chamber

The cutting chamber, made of corrosion resistant materials, is provided with a 225 mm x 386 mm (8.6" x 15.2") cutting table, which is easily moved up to 200 mm (7.9") by control keys. The cutting table is made of corrosion resistant cast aluminium and is fitted with a top surface

of stainless steel bands. These steel bands are easily exchanged if damaged or worn.

A flush hose makes daily cleaning quick and efficient, and a powerful lamp lights up the cutting chamber allowing you to monitor the cutting process through the large window in the protection guard. An exhaust system may be connected.

Quick and Handy Clamping of the Workpiece

The integrated cutting table has six 10 mm T-slots, facilitating the use of universal clamping tools.

Adding an Adjustable Cutting Table for Parallel Sectioning

An adjustable cutting table is available for mounting on the Discotom-6. It positions the workpiece at a right angle to the cut-off wheel, allowing parallel sections to be cut. The table can move up to 60 mm (2.4").

High Safety Level

Discotom-6 cannot be started while the protection guard is open. During cutting the protection guard is locked, eliminating any risk of injury. A built-in electric, dynamic motor brake stops the cutting motor within a few seconds, ensuring quick access to the cutting chamber after cutting.

The hour counter registers the hours the cutting motor has been in operation

For irregular shaped workpieces the vertical clamping tools are employed

Quick-clamping devices enable easy and fast securing of workpieces

The flush hose facilitates cleaning

Technical Data

Cutting Motor	Cutting power constant [S1] Cutting power intermittent [S3, 15%] Maximum power	2.2 - 3.0 kW 3.2 - 4.0 kW 4.4 - 6.0 kW	3.0 - 4.0 HP 4.3 - 5.4 HP 6.0 - 8.0 HP
Cut-off Wheels	Cut-off wheel size: 250 mm x 1.5 mm x 32 mm (10" x 0.06" x 1.25") Rotational speed (running idle): 2850 rpm at 50 Hz, 3420 rpm at 60 Hz		
Cutting Table	Cutting table dimensions Width: 386 mm Length: 225 mm		
Dimensions and Weight	Discotom-6 without recirculating cooling unit Width: 794 mm (31.3") Depth: 640 mm (25.2") Height: 458 mm (18.1") Weight: 123 kg (272 lbs)		

Struers' equipment is in conformity with the provisions of the applicable International Directives and their appurtenant Standards. (Please contact your local supplier for details)

Specifications

Discotom-6

Automatic cut-off machine. Automatic feed with electronic control of feed speed. Ability for manual operation. 3.7 kW (5.0 HP) motor. For 250 mm (10") cut-off wheels. Cutting table with 10 mm T-slots. Recirculation cooling unit (05766xxx), quick-clamping devices (05876909+05876910) or vertical clamping system (05116905) are ordered separately. Optional table unit (05116916) can be used.

Cat.no	Voltage
05676129	3 x 200 V / 50 Hz
05676130	3 x 200-210 V / 60 Hz CSA
05676135	3 x 220-230 V / 50 Hz
05676136	3 x 220-240 V / 60 Hz CSA
05676146	3 x 380-415 V / 50 Hz
05676147	3 x 380-415 V / 60 Hz
05676154	3 x 460-480 V / 60 Hz CSA

Accessories

Cooling System 3

With 50 l tank, small pump, Cooli-1, 50 l static filter with filter paper 05766xxx

Adjustable Cutting Table

For cutting plane parallel sections. It positions the workpiece at a right angle to the cut-off wheel. Movable up to 60 mm. With 10 mm T-slots and exchangeable stainless steel bands 05116902

Spring Clamp for 10 mm T-slots, Right

For retaining the cut-off specimen. Complete with back stop. To be mounted at the right hand side of the cut-off wheel 05116904

Vertical Clamping System for 10 mm T-slots

For clamping irregularly shaped workpieces on the cutting table. Clamping height adjustable up to 58 mm. Complete with operating key and one flat clamping shoe 05116905

Quick Clamping Tool for 10 mm T-slots, right

For securing the workpiece. To be mounted on the right side of the cut-off wheel. Complete with backstop. Max clamping width: 50 mm 05876910

Quick Clamping Tool for 10mm T-slots, left

For securing the workpiece. To be mounted on the left side of the cut-off wheel. Complete with backstop. Max clamping width: 50 mm 05876909

T-slot Turntable for 10mm T-slots

0-90 degree turntable for angular cutting. For use with various 10 mm T-slot clamping tools. Table height 42 mm. 05876905

Riser Block for 10 mm T-slots

For elevating Vertical Clamping System (05116905) when clamping high workpieces. To increase the clamping height by 60 mm 05116910

Swivel Shoes for Vertical Clamping System

Set of 4 multi-shaped swivel shoes for (05116905) 05116911

Replacement Stainless Steel Band

For cutting table on Discotom-6. 1 pc. of 40 mm x 225 mm 05676901

Struers' products are subject to constant product development.

Therefore, we reserve the right to introduce changes in our products without notice.

A wide range of cut-off wheels is available for Discotom-6. Please refer to the separate brochure.

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Struers S.A.S.
Struers Ltd.
Struers Inc.