Epovac

Instruction Manual
Gebrauchsanweisung
Mode d'emploi

Manual No.: 14357001

Date of Release 19.06.2007
Updated 2016.07.01

Instruction Manual

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the *Serial No.* and *Voltage* on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2007/2016.

Struers
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

Epovac

Safety Precaution Sheet

To be read carefully before use

1. The operator should be fully aware of the use of the apparatus according to the Instruction Manual. The apparatus must be placed in an adequate working position.
2. Never use vacuum chambers other than the ones delivered from Struers.
3. If any cracks or a visible colour change occur in the vacuum chamber it must be replaced immediately.
Replace the vacuum chamber every second year.
4. Do not leave the chamber under vacuum and filled with impregnation and / or gluing material fumes for longer than 1 hour.
5. Clean the vacuum chamber with a soft damp cloth and common household detergent or ethyl alcohol. Do not use Acetone or similar solvents.
6. Use only consumables which are suitable for vacuum impregnation in a PVC chamber.
Impregnation and / or gluing materials must be handled with care. See the Instruction Manual and Material Safety Data Sheet for the resin.

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

User's Guide

Table of Contents	Page
1. Getting Started	
Description	2
Checking the Contents of Packing	2
Consumables kit (EPOSU)	2
Getting Acquainted with Epovac	3
Installation	4
Connection	4
Epovac with Built-in Air Ejector (EPOMA).....	4
Epovac without Built-in Air Ejector (EPOCA)	4
2. Operation	
Preparing the Impregnation.....	5
Check vacuum chamber	5
Impregnation	5
3. Reference Guide	
Hints	6
Hardening.....	6
Gluing.....	6
Maintenance.....	7
Technical Data	7
Air supply (only EPOMA)	7
Vacuum connection (EPOCA)	7
Dimensions	7

1. Getting Started

Description

Epovac is a vacuum impregnation apparatus, especially designed for embedding and impregnation of porous materials. Epovac is suitable for gluing specimens to slides for the preparation of thin sections. The apparatus consists of a sturdy, compact metal cabinet with a cylindrical embedding chamber of PVC. The impregnation liquid, e.g. Struers Epofix, is supplied through a disposable tube. The tube is led through a ball joint and it is therefore very easy to distribute the liquid to the individual mounting cups.

Using $\varnothing 25$ mm mounting cups, 9 specimens can be embedded at the same time. As the embedding chamber is large (diameter 134mm, height 100 mm), also large specimens can be impregnated and embedded.

The bottom of the embedding chamber is made of polypropylene so that spilled impregnation liquid does not stick.

Epovac comes in 2 models, with built-in air ejector (to be supplied with compressed air 4.5-6 bar) or without air ejector (to be connected to an external vacuum pump).

The Epovac model without air ejector can be used with most external vacuum pumps.

Checking the Contents of Packing

In the packing box you should find the following parts:


- 1 Epovac basic unit
- 1 Vacuum chamber, $\varnothing 134$ mm, height 100 mm
- 1 Cover
- 1 Pressure rod
- 2 Snap holders for pressure rod
- 1 Hose clip in holder
- 1 Allen key 2.5 mm
- 1 Allen key 5 mm

Consumables kit (EPOSU)

- 2 Rubber plugs with $\varnothing 4$ mm center hole
- 80 Beakers
- 80 Polyethylene hoses
- 80 Wooden stirrers

**Getting Acquainted
with Epovac**

Take a moment to familiarise yourself with the location and names of the Epovac components.

- ① Vacuum chamber
- ② Seals
- ③ Cover before operation
- ④ Impregnation hose
- ⑤ Hose clip
- ⑥ Pivot joint of the cover
- ⑦ Rubber plug
- ⑧ Vacuum valve
- ⑨ Vacuum meter
- ⑩ Stand

Installation

Epovac does not require any special place, however, the apparatus should not be placed in direct sunlight, as the PVC vacuum chamber is UV sensitive.

NB As impregnation liquids like epoxy are common, the apparatus should only be operated where ample ventilation is possible, e.g. in a fume cupboard.

Connection

Epovac with Built-in Air Ejector (EPOMA)

The built-in air ejector enables a vacuum of 90%. This model must be supplied with compressed air (4.5-6 bar) by means of the hose supplied with the apparatus:

- Insert the connecting tube with fitting in the quick coupling on the rear of the apparatus.
- Shorten the hose to an adequate length.
- Mount the hose on the compressor.

Epovac without Built-in Air Ejector (EPOCA)

Connect the apparatus to a vacuum pump by means of the hose, supplied with the apparatus:

- Mount the connecting tube on the fitting on the rear of the apparatus.
- Mount the hose clip and tighten it.
- Shorten the hose to an adequate length.
- Mount the hose on the vacuum pump.

Take care that the hose does not have any sharp bends as these will delay the application of resin. If a longer hose is required it is recommended to choose a hose with a larger inside diameter, in order to secure a sufficient flow.

2. Operation

Preparing the Impregnation

- Place the clean and dry specimen in a suitable mounting cup. If there is more than one, place them symmetrically around the center of the embedding chamber.

Check vacuum chamber

- Check that the vacuum chamber and the cover before operation are not cracked or have fissures, or they might implode when exposed to vacuum.
- Check that the seals are clean and undamaged.
- Mount the vacuum chamber.
- Mount the cover.
- Press the impregnation hose through the rubber plug. Put the impregnation hose through the hose clip and place the rubber plug in the pivot joint of the cover.
- Tighten the hose clip and open for the compressed air (model with pump) or start the vacuum pump (model without pump). Set the vacuum valve to max. The actual vacuum can be read on the vacuum meter.
NB Before impregnation the specimens should be cleaned and degreased.

For embedding, Struers manufactures a number of mounting cups that can be used without releasing agent.

Impregnation

- When the specimens have been under vacuum for an appropriate time - from a couple of minutes for not very porous specimens up to half an hour for very porous specimens - the impregnation can begin:
- Place the beaker with the ready-mixed impregnation liquid in the stand. Check that the hose is near the bottom of the beaker, or air may be sucked into the system and create splashing in the vacuum chamber.
- Turn the hose so that it points down into one of the mounting cups and open slowly for the hose clip, until the impregnation liquid runs in an appropriate flow.
- When the specimen is covered with impregnation liquid, turn off the flow, move the hose to the next mounting cup and open the hose clip. Go on like this until all mounting cups have been filled up.
- Close the hose clip and turn the vacuum valve to min.
- Remove the vacuum chamber with cover and leave the specimens for hardening.

3. Reference Guide

Hints

If 2 or 3 component resins like Epofix are used, it is important that the resin is not heated in order to obtain quicker hardening. The resin will quickly become highly viscous that results in poor impregnation.

Avoid using larger quantities of embedding material than about 20 ml as the heat generation then may become too large with a bad impregnation as a result.

At too high vacuum some of the components in the epoxy evaporate and it looks as if it boils. This can be avoided by reducing the vacuum with the valve.

Hardening

If Epovac is not going to be used anymore the same day, the specimens can remain in the apparatus for hardening. When the specimens have hardened, the impregnation hose can also be removed without any risk of leaving resin on the hose clamp or in the pivot joint.

If the apparatus is going to be used again, the hose with rubber plug is removed from the cover. Keep the rubber plug, as it can be used again. Take care not to spill impregnation liquid on the pivot joint or the hose clip, especially if it is a 2- or 3-component liquid which hardens.

Any spilled impregnation liquid must be removed at once with a suitable detergent.

Gluing

- Apply a suitable amount of resin to the specimen. Usually one drop is enough.
- Place the specimen with the glued surface up, and place a glass slide upon the specimen.
- Mount the pressure rod with rubber plug in the pivot joint of the lid.
- Start the vacuum pump, open the compressed air, and set the vacuum valve to max.
- Place the pressure rod on the glass slide and rub. When the specimen is firmly stuck to the slide, set the vacuum valve to min.
- Remove the slide with glued sample.

Maintenance

- Clean the vacuum chamber with a soft damp cloth and common household detergent or ethyl alcohol.
Do not use Acetone or similar solvents

Model with built-in air ejector
(EPOMA)

- Drain the air filter periodically by loosening the valve screw underneath the apparatus.
- When no more liquid emerges, retighten the valve screw.

Technical Data

Air supply (only EPOMA)

Compressed air: 4.5-6 bar
Air quality: 0.3 micron filtered, oil free
Consumption of compressed air: 12.5 l/min
Hose connection: $\varnothing 1/4$ "
Vacuum (at compressed air 6 bar): -860 mbar (-655 mm Hg)

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Some of the drawings may contain position numbers
not used in connection with this manual.

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment. The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2007

Struers
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800


Epovac
Spare Parts and Diagrams

Drawing

Pos. Spare Part

Part No.

**Fig. 2 Epovac with
Pressure Rod**



Epovac
Spare Parts and Diagrams

Drawing

Fig. 3 Epovac, outline drawing

Pos.	Spare Part	Part No.
1	Top cabinet	14359011
2	Front plate foil	435MP002
3	Plate, complete	435MP003
4	Pressure gauge with fitting	435MP004
5	Ball valve 1/4"	383MP020
6	Base	14359012
7	Holder for beaker	435MP006
8	Stand	435MP007
9	Hose clip	435MP008
10	Vacuum chamber	435MP009
11	V-sealing ring VA130 NBR60, 2 pcs	435MP010
12	Vacuum sucking tube	
	O-ring 12.42 x 1.78 mm	435MP011


Epovac
Spare Parts and Diagrams

Drawing

Fig. 4 Lid

Pos.	Spare Part	Part No.
	Cover, complete	435 MP 012
1	Cover	435 MP 013
2	Flange for cover	435 MP 014
3	Ball joint GL8 O-ring 16 x 2 mm, NBR72 Nylon disc M4 x 0.5 mm, 3 pcs. Screw wls 4 x 12 mm, 3 pcs.	435 MP 015

Epovac
Spare Parts and Diagrams

Drawing


Fig. 5 EPOMA, Cover

Pos.

Spare Part

Part No.

1	Mounting plate for ejector	435 MP 016
2	Air pressure filter, EAF 1000 M5	435 MP 017
3	Vacuum ejector ZH05BS 06-06	435 MP 018
4	Ball valve 1/4"	385 MP 020
5	Pressure air tubing 1/4" x 3 m	435 MP 019
6	Holder for pressure gauge	435 MP 020

Epovac
Spare Parts and Diagrams

Drawing

Fig. 6 EPOCA, Cover without Ejector


Pos.

Spare Part

Part No.

- 1 Holder for pressure gauge
- 2 Connecting piece
- 3 Vacuum tubing 5/16" x 3 m

- 435 MP 020
- 435 MP 021
- 435 MP 022

Epovac
Spare Parts and Diagrams


Drawing

Fig. 7 Pressure rod

Pos. Spare Part

Part No.

	Pressure rod, complete	435 MP 023
1	Rubber plug with Ø4 mm center hole	435 MP 024
	Clips for pressure rod	435 MP 025


Drawing

Pos. Spare Part

Part No.

Fig. 8 14350011

Pneumatic diagram (EPOMA)

Drawing

Pos. Spare Part

Part No.

Fig. 9 14350016

Pneumatic diagram (EPOCA)