RotoPol-31

RotoPol-35

Manual No.: 15287001

Date of release: 22.08.2008

Instruction Manual
Gebrauchsanweisung
Mode d'emploi

RotoPol-31/-35
Instruction Manual

Table of Contents	Page
User's Guide	1
Reference Guide.....	13
Quick Reference Guide.....	31

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2008.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

RotoPol-31/-35

Safety Precaution Sheet

To be read carefully before use

1. The operator(s) should be fully instructed in the use of the machine, any connected equipment and accessories and the applied consumables according to the relevant Instruction Manuals.
2. The machine must be placed on a safe and stable table with an adequate working height. All functions on the machine and any connected equipment must be in working order.
3. Operators should ensure that the actual voltage corresponds to the voltage on the back of the machine. The machine must be earthed. Follow the local regulations. Always turn the power off and remove the plug or the cable before opening the machine.
4. Connect only to cold water. Make sure that the water connections are leak-proof and that the water outlet is working.
5. Struers recommend that the mains water supply is shut off or disconnected if the machine is to be left unattended.
6. Establish efficient exhaust from the working area if consumables containing alcohol, e.g. blue DP-Lubricant, are being used.
7. Make sure that the grinding or polishing discs are correctly mounted, the driving pins engaged with the holes in the turntable.
8. Keep clear of the rotational disc and the specimen mover during operation. Make sure that any screws for clamping of specimens in specimen holders protrude as little as possible. While grinding manually, be careful not to touch the grinding paper. Do not drop the specimen while grinding.
9. Do not lift the specimen mover while preparing individual specimens. RotoForce-4 or Pedemat must be set to Auto.
10. If you observe malfunctions or hear unusual noises - stop the machine and call technical service.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

Disposal

Equipment marked with a WEEE symbol  contain electrical and electronic components and must not be disposed of as general waste.

Please contact your local authorities for information on the correct method of disposal in accordance with national legislation.

User's Guide

Table of Contents	Page
1. Getting Started	
Checking the Contents of Packing.....	2
RotoPol-31 (ROT31) RotoPol-35 (ROT35)	2
Lid, Optional (ROTLI)	2
Recirculation Cooling Unit, Optional (TREC1)	2
Recirculation Cooling Unit Connector, Optional (ROPRE).....	2
Getting Acquainted with RotoPol-31/-35.....	3
Unpacking RotoPol-31/-35.....	4
Placing RotoPol-31/-35.....	4
Supplying Power	4
RotoPol-31	4
Direction of the Turntable on RotoPol-35	4
RotoPol-35 Single-phase Versions	4
RotoPol-35 Three-phase versions	4
Supplying Water	5
Connection to Water Mains	5
Connection to Water Outlet	5
Software Settings.....	6
Setting the Language	7
2. Basic Operation	
Using the Controls.....	8
The Front Panel of RotoPol-31.....	8
The Front Panel of RotoPol-35.....	8
Front Panel Controls	9
General Use of Keys	10
Acoustic Signals.....	10
Reading the Display	10
Routine Operations.....	11
Mounting a Disc.....	11
Starting the Disc.....	11
Stopping the Disc.....	11
Selecting Water and Speed	11
Settings with a RotoCom Memory Unit	12

1. Getting Started

Checking the Contents of Packing

RotoPol-31 (ROT31)

RotoPol-35 (ROT35)

In the packing box you should find the following parts:

- 1 RotoPol-31/-35
- 1 Pressure hose 3/4" x 2 m
- 1 Filter gasket
- 1 Gasket $\varnothing 11/\varnothing 24$ x 1.5 mm
- 1 Reduction ring with gasket
- 1 Drain hose, $\varnothing 38$ mm x 3 m
- 1 Drain elbow pipe
- 2 Hose clamps
- 1 Protection ring
- 1 Set of Instruction Manuals

Lid, Optional (ROTLI)

- 1 Lid

Recirculation Cooling Unit, Optional (TRECI)


- 1 TRECI complete
- 1 Outlet angle
- 1 Funnel
- 2 Disposable plastic containers (TREPO)

Recirculation Cooling Unit Connector, Optional (ROPRE)

- 1 Connection box
- 1 Hexagon nipple
- 1 Non-return valve with elbow and hexagon nipple
- 1 Reduction piece, large
- 1 Reduction piece, small
- 1 Hose clamp
- 1 Diagram


Getting Acquainted with RotoPol-31/-35

Take a moment to familiarise yourself with the location and names of all the RotoPol components:

RotoPol (RotoPol-31)

- A** Front panel controls
- B** Main switch
- C** Water tap with sprinkler tube
- D** Turntable

Back of RotoPol (RotoPol-31)

- ① Water inlet tube
- ② Water outlet tube
- ③ Valve for disc cooling
- ④ Socket for recirculation cooling unit connection
- ⑤ Socket for specimen mover
- ⑥ Fuse

Unpacking RotoPol-31/-35

Detach RotoPol-31/-35 from the bottom of the packing case by removing the four screws from below.

Placing RotoPol-31/-35

Place RotoPol-31/-35 on a steady table with an adequate working height. The machine must be close to the power supply, water mains and water outlet facilities. If water recirculation is used, there must be room under the table for the Recirculation Cooling Unit (TREC).

Supplying Power

RotoPol-31/-35 is factory mounted with electric cable.

<p>IMPORTANT Check that the mains voltage corresponds to the voltage stated on the type plate on the back of the machine.</p>
--

RotoPol-31

Mount a plug on the cable or connect to mains:

Yellow-green or just green: earth
Other colours: phase

Direction of the Turntable on RotoPol-31

Check that the turntable rotates in the direction indicated by a label on the turntable (anti-clockwise). If the direction of rotation is incorrect, switch two of the phases.

RotoPol-35
Single-phase Versions

Mount a plug on the cable or connect to mains:

Yellow-green: earth
Blue or white: neutral
Other colours: phase

RotoPol-35
Three-phase versions


Mount a plug on the cable or connect to mains:

Yellow-green or just green: earth
Other colours: phase

 <p>DANGER! The machine must be earthed</p>
--

Supplying Water

Water for wet grinding and disc cooling may be supplied from the water mains or a Recirculation Cooling Unit.

RotoPol (RotoPol-31)

- ① Water inlet tube
- ② Water outlet tube

Connection to Water Mains

- Mount the straight end of the pressure hose onto the water inlet tube (upper) in the hose box on the back of RotoPol-31/ -35:
 - Insert the filter gasket in the coupling nut with the flat side against the pressure hose.
 - Tighten the coupling nut completely.
- Mount the other end of the pressure hose on the water mains tap for cold water:
 - Mount the reduction ring with gasket on the water mains tap, if necessary.
 - Introduce the gasket and tighten the coupling nut completely.

IMPORTANT
Only connect to cold water.

Connection to Water Outlet

- Determine the best arrangement for the outlet. If the drain hose is directed parallel to the machine, the drain elbow pipe should be mounted on the water outlet tube (lower) in the hose box on the back of RotoPol-31/-35.

- Locate the end of the drain hose where a short length of the steel reinforcement has been stripped away. Mount the stripped end of the hose onto the water outlet tube, or onto the drain elbow pipe if such is mounted. Use hose clamp for fastening.
- Lead the other end of the drain hose to the water outlet. Arrange the tube so that it slopes throughout its length. Shorten the hose, if necessary. The drain hose should be attached to the water outlet system if possible.

REMEMBER
Avoid sharp bends in the drain hose and make sure that the hose slopes throughout its length.

Software Settings

- Switch on the power at the main switch. The following display will appear briefly:

- ROTOPOL-31 -
VERSION X.Y

- ROTOPOL-35 -
VERSION X.Y

- When the operation menu appears, you are ready to choose the language you prefer in the display.

Setting the Language

WATER OFF 150RPM


Example screen


The MENU  key activates the set-up menu.

DISC COOLING
WHEN WATER OFF


Press MENU  again. The language screen appears.


LANGUAGE
ENGLISH


Press SELECT  to choose between English, German, French or Japanese

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Press SELECT  until you reach the language you prefer.

LANGUE
FRANCAIS

Press STOP  to accept the choice.

The Operation Menu will re-appear in the language you have chosen.

EAU OFF 150 T/M


Remember
to set the same language on all connected RotoSystem modules. Setting the Communication and Disc Cooling: see Advanced Operation Section.

2. Basic Operation

Using the Controls



The Front Panel of RotoPol-31

RotoPol-31

The Front Panel of RotoPol-35

RotoPol-35









① Main Switch: figure A

⑮ Water Tap with Sprinkler Tube: figure B



RotoPol-31/-35
Instruction Manual

Front Panel Controls

Communication is set to STANDARD mode

Name	Key	Function	Name	Key	Function
① MAIN SWITCH		The main switch is located to the right of the cabinet. Turn clockwise to switch on the power.	⑨ SPEED 300 RPM		Select 300 RPM. (RotoPol-31)
② POWER INDICATOR		Lights up when the main power is on.	⑩ SPEED 150 RPM		Select 150 RPM. (RotoPol-31)
③ START		Starts the disc and any connected RotoSystem equipment set in Auto mode.	⑪ MENU		Activates the set-up menu or moves down in the menu structure.
④ STOP		Stops the disc and any connected RotoSystem equipment set in Auto mode. Exits the set-up menu.	⑫ SELECT		Changes values in the set-up menu.
⑤ DISC		Scrolls up in the operation menu.	⑬ SPEED		Increases the speed. (RotoPol-35)
⑥ DISC		Scrolls down in the operation menu.	⑭ SPEED		Decreases the speed. (RotoPol-35)
⑦ WATER ON		Select WATER ON: press the key briefly. <i>Manual Dosing:</i> keep pressing the key. Release the key to stop the water.	⑮ WATER TAP		Adjust the water on the tap.
⑧ WATER OFF		Select WATER OFF.			

General Use of Keys

- WATER ON , WATER OFF  and SPEED  can only be selected freely at User Select.
- MENU , SELECT  and DISC  can only be activated when the machine has been stopped.
- SELECT  can only be activated when MENU  is selected.

Acoustic Signals

- When pressing a key, a short beep indicates that the command has been accepted, whereas a long beep indicates that the key cannot be activated at the moment.

Reading the Display

**WATER ON 300RPM
GRINDING**


- The first line on the screen shows the setting of parameters WATER and SPEED.
- The second line, if any, shows a selected pre-set mode (GRINDING or POLISHING).

Routine Operations

Mounting a Disc

- Place the disc on the turntable and turn it until the three driving pins engage with the holes in the turntable.

Starting the Disc

- Press START . The disc starts rotating. Selected functions are executed and connected RotoSystem modules start operating.


Stopping the Disc

- Press STOP . The disc stops rotating. Selected functions and connected RotoSystem modules stop.

Selecting Water and Speed

You can either select water and speed yourself or press DISC  to select one of the pre-set GRINDING and POLISHING modes.

Setting to User Select

The User Select screen has an empty line at the bottom.


- Press  or  to select WATER ON or WATER OFF.
- Press SPEED  or  to select the speed. (RotoPol-31)
- Press SPEED  or  to adjust the speed. (RotoPol-35)

Press DISC  or  to change between the settings.

Setting to GRINDING


The standard setting GRINDING is pre-set to WATER ON and 300RPM.

Press DISC  or  to change between the settings.

Setting to POLISHING

The standard setting POLISHING is pre-set to WATER OFF and 150RPM.

*Settings with a
RotoCom Memory Unit*

If RotoPol is connected to RotoCom and communication is in ROTOCOM mode, the following screen will appear:

**WATER ON 300RPM
PG-Paper #120**

(Example screen)

- The first line on the screen shows the settings of parameters for WATER and SPEED.

- The second line shows the surface to be used, according to the preparation program selected on RotoCom.

DISC  and  toggles between the various surfaces defined in RotoCom.

STRUERS CONSUMABLES

Struers offers a wide range of:

- Consumables for both plane and fine grinding.
- Polishing cloths, diamond and oxide abrasive materials.

Reference Guide

Table of Contents	Page
1. Advanced Operations	
Selection of Disc	14
Disc Cooling	14
Setting the Disc Cooling	14
Adjustment of Water Flow at Disc Cooling	15
Mounting of Specimen Mover (Optional).....	16
RotoForce-3/4	16
Pedemat/Pedemax-S.....	16
Adjusting Pedemat	17
Mounting of Multidoser (optional)	18
RotoCom Memory Unit (optional).....	18
Setting the Communication	18
Recirculation Cooling Unit (Optional).....	19
Struers Metalog Guide™.....	20
Oxide Polishing	21
2. Overview of Menu Screens in the Set-up Menu	22
3. Accessories and Connected Equipment	23
4. Trouble-shooting	24
5. Maintenance	
Daily Service	26
Weekly Service.....	26
Replacing the Cooling Water.....	26
Checking the Recirculation Cooling Unit.....	26
Replacing the Swivel for Disc Cooling	27
6. Technical Data	28

1. Advanced Operations

Selection of Disc

RotoPol-31/35 is designed for MD-Discs of 300 mm dia., wet grinding discs of 305 mm dia. and polishing discs of 300 mm dia. However, the whole range of Struers discs can be mounted.

- For automatic preparation of specimen holders mounted in specimen movers, only aluminium discs should be used. See the instructions for RotoForce-3 or RotoForce-4.
- When disc cooling is selected, only aluminium discs should be used.

Disc Cooling

Disc cooling should be selected in the setting whenever RotoPol-31/-35 is used intensively for:


- Preparation of specimen holders.
- Polishing time longer than 5 minutes.
- Successive preparation of several specimen holders on the same disc.

Important

The set-up menu setting can only be implemented when the machine has been stopped.

Setting the Disc Cooling

- Press MENU  to move from one screen to another.
- When in the DISC COOLING menu, press SELECT :

- Press STOP  when the desired setting appears, or press MENU  repeatedly until the operation menu appears again.
- WHEN WATER OFF should normally be selected when disc cooling is applicable.
- ALWAYS ON should only be selected for special applications.

- ALWAYS OFF should be selected if disc cooling is not relevant for the application.


Disc cooling will automatically start and stop with RotoPol-31/-35.

IMPORTANT

- Disc cooling will only be active at rotational speeds up to 300 RPM.
- Disc cooling is only efficient when aluminium discs are being used.

Adjustment of Water Flow at Disc Cooling

- Select mode WHEN WATER OFF or ALWAYS ON in the disc cooling menu.
- Select WATER OFF  and SPEED  150, and start RotoPol-31/-35 without a preparation disc.
- Adjust the flow using the valve at the rear of the cabinet so that the cooling water flows from the hole in the centre of the turntable at a height of about 2 mm. This will result in a water flow of about 250 ml/min which will generally be adequate.

RotoPol (RotoPol-31)

③ Valve for Disc Cooling

Disc cooling works with both tap water or by means of recirculation. The cooling capability for Recirculation Cooling Unit (TRECI) will be sufficient for most situations. If, however, the equipment is used very intensively, using recirculation, the temperature of the cooling water may vary too much during working time to be adequate for disc cooling. In this situation, a separate cooling system (available from most laboratory suppliers) should be connected to the Recirculation Cooling Unit.


Mounting of Specimen Mover (Optional)

RotoForce-3/4

RotoPol-31/-35 can be equipped with a *RotoForce-3* or a *RotoForce-4* specimen mover. To mount the specimen mover, see the Instruction Manual of the specimen mover.

Pedemat/Pedemax-S

Older versions of the specimen movers: *Pedemat* (PEDRO) or *Pedemax-S* (PEDSO) can be mounted on RotoPol-31/-35, operating with automatic start and stop. In such a configuration, however, a RotoCom unit cannot be connected.

① Fastening screws.


② Support hole.

- Remove the cap covering the support hole placed on the top of the RotoPol-31/-35 cabinet by lifting it with a knife.
- Remove the two caps covering the fastening screws on the left-hand side of the RotoPol-31/-35 cabinet by lifting them with a knife.
- Remove the black plastic ring from the supporting column of Pedemat/Pedemax-S.
- Lead the supporting column of the specimen mover down into the support hole in RotoPol-31/-35. Turn the column until the cross pin in the bottom of the column gets caught by the guiding slots inside RotoPol-31/-35.
- Tighten the two fastening screws on the left-hand side of the RotoPol-31/-35 cabinet by means of the hexagon key included in the packing case, and mount the two caps again.

Adjusting Pedemat

Adjust the height of the specimen mover in relation to the preparation disc:

- Mount a preparation disc on RotoPol-31/-35.
- Mount a specimen mover plate in the coupling of Pedemat.
- Lower Pedemat into work position.
- Measure the distance between the specimen mover plate and the preparation disc. The correct distance is 1-2 mm. If the distance is not correct, the stop ring on Pedemat can be disengaged by loosening the retaining screw and adjusting the height. Remember to fasten the screw again after the adjustment.

RotoPol (RotoPol-31)

⑤ Socket for Specimen Mover

Establish the electronic connection between RotoPol-31/-35 and the specimen mover: connect the cable from the specimen mover to the socket at the rear of RotoPol-31/-35.

IMPORTANT

The distance between the specimen mover plate and preparation disc should be 1-2 mm. Do not tighten the fastening screws on the left-hand side of the RotoPol-31/-35 cabinet without having a column mounted.

IMPORTANT

If RotoPol-31/-35 has to be relocated, the specimen mover should be removed.

**Mounting of
Multidoser (optional)**

RotoPol-31/-35, combined with a specimen mover, can be equipped with *Multidoser* for automatic dosing of suspensions and lubricants. To connect Multidoser, see the appropriate Instruction Manual.

**RotoCom Memory Unit
(optional)**



RotoPol-31/-35, combined with RotoForce-3 or RotoForce-4 and Multidoser, can be equipped with a *RotoCom* memory unit. To connect RotoCom to RotoPol, see the appropriate Instruction Manual.


Setting the Communication

- Press MENU  to activate the set-up menu.

COMMUNICATION
STANDARD

Move to the COMMUNICATION menu using the MENU  keys.

Press SELECT  to toggle between STANDARD and ROTOCOM:


COMMUNICATION
ROTCOM

- Select STANDARD if RotoPol is used as a single unit or if combined with a RotoForce-3/-4 and/or Multidoser.
- Select ROTOCOM when RotoPol is connected to a RotoCom unit.


- Press STOP  when the desired setting appears, or press MENU  repeatedly until the operation menu appears again.

**Recirculation Cooling Unit
(Optional)**

A Recirculation Cooling Unit (TRECI) can be connected by means of a Recirculation Cooling Unit Connector (ROPRE).

- Connect the cable from the pump in TRECI to ROPRE, according to the diagram of ROPRE.
- Connect the cable with plug from ROPRE to the socket marked Recirculation on the back of RotoPol-31/-35.

RotoPol (RotoPol-31)

④ Socket for Recirculation

- Connect the cable from ROPRE to the mains.
- Check that the pump rotates in the direction indicated on the pump. If the direction is incorrect, switch two of the phases in the connection to the pump in TRECI (3-phase versions only).
- Replace the outlet tube on the pump in TRECI with the non-return valve with elbow and hexagon nipple from ROPRE. Struers recommends that you seal the threads.
- Connect the pressure hose of RotoPol-31/-35 to the hexagon nipple. Remember to mount the gasket.
- Select the reduction piece from ROPRE which fits the drain hose of RotoPol-31/-35. Mount the piece on the outlet angle of the TRECI.
- Shorten the drain hose if necessary. Fix the drain hose on the reduction piece: strip away a short length of the steel reinforcement in the drain hose and mount a hose clamp on the stripped end.
- Mount the outlet angle on TRECI.

REMEMBER

Avoid sharp bends in the drain hose and make sure that the hose slopes throughout its length.

Struers Metalog Guide™

RotoPol-31/-35 is designed for both grinding and polishing. In Struers Metalog Guide™ you will find a detailed description of grinding/polishing for automated mechanical specimen preparation.

Struers Metalog Guide™ offers preparation methods for the most common materials, based on a simple analysis of two key properties: hardness and ductility. Finding the right method, including choice of consumables is easy. Always consult Struers Metalog Guide™ for the correct preparation method for the actual specimens.

Struers Metalog Guide™ contains 6 useful chapters:

Metalogram: a quick and safe guide to the right preparation method.

Metalog Methods: a complete catalogue of preparation methods, based on Struers' extensive materialographic experience, and employing Struers' range of consumables.

Preparation Philosophy: the basics of modern specimen preparation, seen from a professional point of view.

Metalog Process: the materialographic preparation process from start to finish, logically explained.


Metalog Master: a combined trouble-shooting guide and supply of in-depth information on the processes of mechanical preparation, including an expert system for the solving of preparation problems.

Consumables Specification: quick access to the relevant consumables for the chosen preparation methods.

Struers Metalog Guide™

A complete guide to materialographic specimen preparation.
Contact your local dealer for a free copy of Metalog Guide™.

Oxide Polishing


For polishing with Oxide Polishing Suspension: at the end of the polishing process, select WATER ON . Adjust the water tap if necessary. Let the water flow until the oxide polishing suspension has been cleaned away (about 10 seconds).

IMPORTANT

If a Recirculation Cooling Unit is connected, special precautions must be taken:

- The water coming from the built-in water tap cannot be used to wash the oxide polishing suspension away. Use clean water.
- The outflow must not be conducted to the Recirculation Cooling Unit. Direct the outlet to the drain or collect the outflow in a separate container.

2. Overview of Menu Screens in the Set-up Menu



3. Accessories and Connected Equipment

Specification	Code
Accessories	
<i>Wet Grinding Discs with Ring</i>	
250 mm dia., aluminium	ROTIF
250 mm dia., PVC	ROTRO
305 mm dia., aluminium	ROTAL
<i>Polishing Discs</i>	
250 mm dia., aluminium	DEDIF
300 mm dia., aluminium	DEDAL
<i>Discs for Magnetic Fixation</i>	
MD-Disc, 250 mm dia.	DEMIF
MD-Disc, 300 mm dia.	DEMAL
Connected Equipment	
<i>RotoForce-4, Specimen Mover</i>	ROFO4
<i>RotoForce-3, Specimen Mover</i>	ROFO3
<i>Multidoser (6 Bottles)</i> for automatic dosing of 6 liquids (1-5 suspensions, 0-4 lubricants and 1 OP-Suspension)	MULDO
<i>Multidoser Base Unit (3 Bottles)</i> for automatic dosing of 3 liquids (0-3 suspensions, 0-3 lubricants)	MULBA
<i>RotoCom, Memory Unit</i>	ROCOM
<i>Recirculation Cooling Unit</i>	TRECI
<i>Recirculation Cooling Unit Connector</i> for connection of Recirculation Cooling Unit (TRECI) to grinding/polishing machines	ROPRE
Spare Parts	
Swivel	2NF10130

4. Trouble-shooting

Error	Explanation	Action
Error Message		
OVERLOAD	The motor has been overloaded for a period and the thermal switch has stopped the motor.	Wait until the motor has cooled down. When the display returns to the previous message, the machine can be re-started.
Acoustic Signal		
Long beep.	The selected function is protected. It cannot be activated.	Check the setting.
Machine Problems		
Noise when the machine starts or the machine will not pull.	The belt is not tight enough.	Call a Struers service technician. The belt must be tightened.
Disc rotates in the wrong direction. (RotoPol-31)	Wrong connection to the power supply.	Switch two of the phases.
Motor buzzes. (RotoPol-31)	One phase lacking.	Check all fuses of local power supply.
Functions are not carried out.	Fuse at the rear of RotoPol-31/-35 blown.	Replace the fuse.
RotoCom and/or Multidoser do not operate.		
Water not draining away.	Drain hose squeezed.	Straighten the hose.
	Drain hose clogged.	Clean the hose.
	Drain hose does not slope downward.	Adjust the hose to an even slope.
Water dripping underneath the machine.	The water does not drain away.	See above.
	The swivel for disc cooling is worn.	Replace the swivel. See maintenance. In the meantime the software setting for disc cooling should be changed to ALWAYS OFF.
Cooling water stops.	Wrong software setting.	Check software setting.
	Water tap on mains closed.	Open for water.
	Built-in water tap closed.	Open for water.
	Built-in water tap blocked.	Clean water tap.
	Filter at the water inlet blocked.	Clean filter.
Disc cooling stops.	Wrong software setting.	Check software setting.
	Wrong adjustment of water flow at disc cooling.	Adjust the flow at the valve at the rear of the machine.
Cooling water drips after stop.	Defect in the solenoid valve.	Call a Struers service technician. The solenoid valve must be replaced.

*RotoPol-31/-35
Instruction Manual*

Error	Explanation	Action
The disc vibrates.	Dirt on the underside of the disc or on the turntable.	Clean the contact face between the disc and turntable.
JAPANESE TEXTS	RotoPol-31/-35 has the possibility of showing the screen texts in Japanese. When choosing a European language in the LANGUAGE menu, you might accidentally press  MENU at the Japanese option.	To return to English, German or French texts, just follow this procedure: Press  STOP. Press  MENU two times. Press  SELECT until you reach the language you prefer. Press  STOP.

5. Maintenance

Daily Service

Clean all accessible surfaces with a moist cloth.

Weekly Service

Remove the grinding/polishing disc and remove all dirt from the tray and the drain tube.

Replacing the Cooling Water

Only relevant when a Recirculation Cooling Unit is connected:
Replace the cooling water in the Recirculation Cooling Unit at least once a month.


- Place a Disposable Plastic Container (TREPO) in the tank and fold it over the edge.
- Fill the tank with 30 l of water. NB! Too high a level of coolant in the tank might damage the pump.
- Add Struers Additive according to the instructions.
- Put the lid on the tank.
- Place the pump in the lid.

IMPORTANT

Always maintain the correct concentration of Struers Additive in the cooling water (percentage stated on the container of the Additive). Remember to add Struers Additive each time you refill with water.

Checking the Recirculation Cooling Unit

The cooling unit should be checked for cooling water every week and refilled when the flushing pump cannot reach the cooling water. Remember to use Struers Additive (pH should be around 9). Replace if pH is α 8.0. To check the concentration of Additive use a refractometer (0-10 Brix, see diagram.). Concentration = 2 x Brix.

Replacing the Swivel for Disc Cooling


Only relevant when disc cooling is being used:

Struers recommends that the swivel be replaced at least every two years.

- Turn off the power and remove the plug and the cable before opening the machine. Place RotoPol-31/-35 as indicated in suggestion 1 or 2:

Suggestion 1

Place RotoPol-31/-35 on two tables with the rubber feet resting on the tables as shown in the sketch, so as to have access to the underside of RotoPol-31/-35.

- ① Swivel.
- ② Bottom plate.

Alternative Suggestion 2

Remove the specimen mover, if mounted, and carefully tilt RotoPol-31/-35 to one side so as to have access to the base.

- Remove the protection guard for the swivel.
- Remove the tube on the swivel: push the ring and pull off the tube. Pull the tube out of the clips.
- Remove the bottom cover.
- Support the belt wheel and unscrew the swivel with an open-end spanner.
- Fasten the new swivel with a moderate torque and mount the parts again.

6. Technical Data

Subject		Specifications		
		RotoPol-31	RotoPol-35	
Disc	Diameter	200 - 300/305 mm / 8" - 12"		
	Rotational Speed	150/300 rpm *)	40 - 600 rpm	
	Rotational direction	Anti-clockwise		
	Motor	300 W / 0.41 HP at 150 rpm	1000 W / 1.4 HP	
		600 W / 0.82 HP at 300 rpm		
	Torque at disc	Continuous 19 Nm / 14 lbs-ft	Continuous 18 Nm / 13 lbs-ft	
Max. > 30 Nm / 22 lbs-ft		Max. > 20 Nm / 15 lbs-ft		
Water Supply	Tap water or Recirculation Cooling Unit (TREC)			
	Pressure for tap water	1 - 10 bar / 14.5-145 psi		
	Inlet	1/2" or 3/4"		
	Outlet	ø38 mm / 1.5"		
Dimensions and weight	Width	515 mm / 20.3"		
	Depth	720 mm / 28.4"		
	Height	335 mm / 13.2"		
	Weight	42.5 kg / 93.7 lbs	49.0 kg / 108 lbs	

*) Soft start at 300rpm

RotoPol-31/-35
Instruction Manual

Subject		Specifications		
		RotoPol-31		RotoPol-35
	LC display	2 x16 characters		
	Controls	Touch pad		
	Memory	EPROM, RAM and non-volatile RAM		
	Automatic start/stop	With: Recirculation Cooling Unit (TREC1) connected with Recirculation Cooling Unit Connector (ROPRE) Older versions of specimen movers: Pedemat (PEDRO) or Pedemax-S (PEDSO)		
	Communication including automatic start/stop	With: RotoForce-3 RotoForce-4 Multidoser		
	Data exchange	With: RotoCom		
	Power outlet	For: Multidoser RotoCom Recirculation Cooling Connector		
Environment	Safety standards	General 89/392/EEC & 91/386/EEC & 93/44/EEC EN292-1, EN292-2, EN60204-1 (IEC 204-1) EMC 89/336/EEC & 92/231/EEC EN50081-1, EN50082-1 FCC part 15, subpart B		
	Noise Level	Approx. 56 dB(A) at idle running with max. rotational speed, at a distance of 1.0 m / 39.4" from the machine		
	Surrounding temperature	5-40/C / 41-104/F		
	Humidity	0-95% RH non condensing		
Supply Voltage	Power phases	3 (3L + PE)	Single-phase versions: 1 (L + N + PE) Three-phase versions: 3 (3L + PE)	
	Power consumption	Max continuous load 960 W / 1.3 HP	Max. continuous load 1270 W / 1.7 HP	Max. load 2500 W / 3.4 HP

RotoPol-31/-35
Instruction Manual

Subject		Specifications		
		RotoPol-31	RotoPol-35	
	Voltage / frequency:	Current:		
		Max continuous load	Max. continuous load	Max. load
	3x200 V / 50 Hz	5.0 A		
	3x200 - 210 V / 60 Hz	5.0 A		
	3x220 - 230 V / 50 Hz	5.0 A		
	3x220 - 240 V / 60 Hz	5.0 A		
	3x380 - 415 V / 50 Hz	3.0 A		
	3x460 - 480 V / 60 Hz	3.4 A		
	1x220 - 240 V / 50 Hz		5.3 A	12.5A
	1x220 - 240 V / 60 Hz		5.3 A	12.5 A
	3x200 - 210 V / 50 - 60 Hz		6.2 A	12.5 A
	3x220 - 240 V / 60 Hz		5.3 A	12.5 A
	±5% on nominal voltage			

Quick Reference Guide

Mounting a Disc

Place the disc on the turntable and turn it until the three driving pins engage with the holes in the turntable.

Starting the Disc

Press START . The disc starts rotating. Selected functions are executed and connected RotoSystem modules start operating.

Stopping the Disc


Press STOP . The disc stops rotating. Selected functions and connected RotoSystem modules stop.

Grinding/Polishing

- Mount the required surface on a disc.
- Mount the disc on the turntable.
- Select the desired functions, if RotoCom is not connected.
- Press START .

Settings

Setting to User Select

The User Select screen has an empty line at the bottom.


- Press or to select WATER ON or WATER OFF.
- Press SPEED or to select the speed. (RotoPol-31)
- Press SPEED or to adjust the speed. (RotoPol-35)

Press DISC or to change between the settings.

Setting to GRINDING


The standard setting GRINDING is pre-set to WATER ON and 300RPM.

Press DISC or to change between the settings.

Setting to POLISHING

The standard setting POLISHING is pre-set to WATER OFF and 150RPM.

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Spare Parts and Diagrams

RotoPol-31/-35
Spare Parts and Diagrams

**Always state *Serial No* and *Voltage/frequency*
if you have technical questions or when ordering spare parts.**

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
Some of the drawings may contain position numbers
not used in connection with this manual.

RotoPol-31/-35
Spare Parts and Diagrams

Drawing	Pos.	Spare Part	Part No.
14920013		RotoPol-31, complete	
15280001		RotoPol-35, complete	
	120	Cover plate, complete	14920160
	120	O-ring	2IO20060
	130	Relief clamp	14920121
14920038		Bottom Frame, RotoPol-31	
15280020		Bottom Frame, RotoPol-35	
	030	Rubber foot	2GB00001
	260	2.50A T fuse glass tube	2FU14000
	270	Fuse holder/L 2004	2XS12004
	340	Switch, 3-poled	2SE00125
	390	Poly V-belt, L = 1168mm	2JD31168
14920017		Bearing Construction, RotoPol-31 and -35	
	120	O-ring	2IO16005
	130	Turntable mounted with Quad rings	14920050
	150	Swivel	2NF10130
14920029		Cabinet, RotoPol-31	
15280040		Cabinet, RotoPol-35	
	020	Throttle valve	2YH00008
	030	Single banjo	2NF10029
	040	Aluminium ring	2NF20151
	050	Reduction socket	2NF40061
	060	Hose nipple	2NF40087
	070	Solenoid valve	2YM12306
	100	Hose 3/8" (please state length)	2NU29312
	110	Hose clamp S14	2NS11409
	110	Hose clamp S16	2NS11609
	120	Hose piece 510mm with 3" bypass.	14960022
	160	Hose r4 (please state length)	2NU15704
	170	PVC tube	14920028
	180	Worm belt 32-50	2NS23250
	220	Pipe section	14910070
	230	Nut 3", brass	12660049

RotoPol-31/-35
Spare Parts and Diagrams

Drawing	Pos.	Spare Part	Part No.
14920045		Water Tap	
		Water tap, complete	14920045
	090	Spout	2NU17035
		Gaskets for water tap	14920503
		Various parts for RotoPol-31/-35	
	010	Pressure hose for water, white	2NU93020
	020	Connection piece 2" - 3/4"	2NG30013
	030	Gasket r11/r24 x 1.5	13590359
	040	Gasket with filter 3/4"	2IX20410
	050	Elbow 87E with socket joint r40	2NG20440
	060	Hose r38 (please state length)	2NU30238
	070	Worm belt 32-50	2NS23250
	090	Splash ring	14920610

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




* MUST BE SHORT CIRCUIT PROTECTED WITH EXTERNAL FUSES ACCORDING TO LOCAL REGULATIONS
F = 25 Amp MAXIMUM

FOR DETAILED CONNECTION ON PRIMARY SIDE OF T1 SEE OTHER DRAWING.

COLOR CODES:
WH = WHITE
BR = BROWN
GN = GREEN
YL = YELLOW
GY = GREY
RD = RED
BL = BLUE

Rev.A: A1-X8 corrected to A1-X7 Rev.B: Fuses = 16A -> 25Amp.		STRUERS A/S VALHOEJESALLE 176 DK-2610 ROEDOVRE DENMARK PHONE: +45 3670 3500	
RotoPol-31 CIRCUIT DIAGRAM			
Size C	CAGE Code <Cage Code>	DWG NO 14923100	Rev B
Tuesday, June 06, 2000		Scale FTH / FTH	Sheet 0 1-1 of 0

ALL WIRINGS 0,75², EXCEPT OTHERWISE MARKED.

- COLOR CODES:
- BK = BLACK
 - BN = BROWN
 - RD = RED
 - YE = YELLOW
 - GN = GREEN
 - BU = BLUE
 - GY = GREY
 - WH = WHITE

Rev. B: F1= 3,15A=25A Rev. C: W1 CSA Colour coding of conductors changed! BK → CR(BL), YE/GN → GN(GN/YE) (11-02-00DEM) Rev. D: New motor controller A2(4PC) Rev. E: Change in connection A2(4PC) Rev. F: (PTH 24-102099) Wire colors between Z1 and D1 changed		STRUERS A/S Pedersbølvej 84 DK-2750 Ballerup DENMARK PHONE: +45 44 600 800	
Block- & Circuit Diagram:			
RotoPol-35			
Size A2	CAGE CODE	DWG NO 15283100	Rev F
Monday, October 24, 2005		SCALE 03.03.97 / DEM D.E.M.	Sheet 1 of 1

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Herewith declares that

<i>Product Name:</i>	RotoPol-35
<i>Type No:</i>	528
<i>Machine Type:</i>	Grinding and polishing machine

is in conformity with the provisions of the following directives:

Safety of Machinery 2006/42/EC according to the following standard(s):
EN ISO 12100-1:2005, EN ISO 12100-2:2005, EN 60204-1:2006/A1:2009.

EMC-Directive 2004/108/EC according to the following standard(s):
EN61000-6-1:2002, EN61000-6-3:2002.

Low Voltage Directive 2006/95/EC according to the following standard(s):
EN 60204-1:2006/A1:2009.

Supplementary Information The equipment complies with the American standards:
FCC part 15, subpart B and UL508.

The above has been declared according to the global method, module A

Date: 31.01.2011

Christian Skjold Heyde,
Vice President, R&D and Production, Struers A/S

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Christian Skjold Heyde,
Vice- President, R & D et Production, Struers A/S

 **Struers**

**Maximum
Reproducibility
and Flexibility**

RotoSystem
from **Struers**

The RotoSystem

enables you to choose from several grinding/polishing machines, and different specimen movers. You can also add a dosing unit and a memory unit for the control of all preparation parameters. **A maximum of flexibility** is provided, enabling you to find the best solution for your specimen preparation.

At the same time maximum reproducibility is guaranteed, as both the memory unit for all preparation parameters and the automatic doser will secure that the same results are achieved every time a preparation method is carried out.

Even different laboratories within the same group can share the same preparation methods, as the RotoCom unit can be interfaced with a PC for easy up- and downloading of preparation methods.

To demonstrate the potential of RotoSystem we have designed three complete systems, which will cover most applications in a modern laboratory. The RotoSystem 200 is based on the RotoPol-11 with 200 mm discs. RotoSystem 250 is based on a RotoPol-22, where 250

RotoSystem 200, for small laboratories or a limited number of specimens

mm dia. discs are used, and the RotoSystem 300 is built up on a RotoPol-31 with 300 mm dia. disc size.

RotoSystem 200, the small system

The combination of a RotoForce-1 specimen mover on a RotoPol-11 is suitable for the preparation of a limited number of specimens. Either 1-3 single specimens or 3 specimens clamped in a holder can be prepared at a time. If the specimens arrive at the laboratory one at a time or in small quantities, the fastest and most economic preparation will be achieved on the RotoSystem 200. Together with Multidoser for accurate and reproducible dosing of suspensions and lubricant and RotoCom with pre-pro-

grammed preparation methods, perfect results can be obtained every time.

RotoSystem 250, the smart system

Ideal for the lab where the workload is variable. The RotoSystem 250 can be used for the preparation of single specimens as well as for the preparation of specimens clamped in holders. No need to use "dummy" specimens. If only one specimen has to be prepared, simply switch to single specimen mode on the RotoForce-4 and run only one specimen. That means a considerable saving of consumables over the years, as one specimen requires less diamond suspension than six specimens. Next time a whole load of specimens arrives, clamp up to twelve into one specimen holder. Handling is a lot easier than with twelve single specimens. The RotoPol-22 grinder/polisher has two individually controlled discs. While one is in operation, the second disc can be prepared according to the information already sent by RotoCom. RotoCom, with the built-in Metalog Guide methods, automatically sends force values for both single specimens and specimen holders to the specimen mover. The appropriate value for that set-up will then be used.

RotoSystem 300, the powerful system

For the preparation of many or large specimens.

RotoSystem 250, the ideal solution for the lab with a varying workload

RotoSystem 300, the workhorse for the preparation of many or large specimens

The RotoSystem 300 with its 300 mm dia. disc and the strong RotoForce-3 specimen mover is especially suitable for the preparation of a large number of specimens and also for specimens of a large size. The built-in disc cooling of RotoPol-31 is especially valuable in these cases, as the temperature on the preparation disc will stay low and with little variation during the entire preparation process. The amount of lubricant used during preparation will be reduced accordingly.

RotoForce-3 is designed to handle specimen holders of up to 200 mm dia., making it possible to either prepare very large specimens or many small specimens in one of Struers standard specimen holders. Of course, holders for many different sizes and shapes are available or can be manufactured according to your specifications.

System equipment

All three systems are equipped with the Multidoser automatic dosing unit and the RotoCom memory unit. The

use of Multidoser gives you controlled dosing of both diamond and oxide polishing suspensions as well as of lubricants. RotoCom is used to store up to one hundred preparation methods with a maximum of 10 steps each.

All preparation parameters are saved in RotoCom and can be recalled at any time. Now you do not have to wonder whether the preparation has been carried out in exactly the same way as last time, the RotoSystem provides you with this guarantee.

With RotoCom it is even possible to connect the system to a PC. That means you can copy preparation methods and send them to a lab at your other plant, or to your supplier or anybody else who has to use the exact same method.

Flexibility -an Overview of the Different Modules

All modules of the RotoSystem are equipped with advanced microprocessors. They are operated on touch pads and all parameters can be read on the LC-Display and adjusted by menus. The RotoSystem is based on four inter-connectable modules, available in different versions.

Rotopol grinding/polishing machines

200 mm single disc:

RotoPol-11. With two speeds RotoPol-11 is capable of both grinding and polishing in all standard applications.

RotoPol-15. For special requirements the variable speed from 40-600 rpm is the correct choice.

250 mm twin disc:

RotoPol-21. Driven by a single motor allowing for two speeds, is suited for all standard applications.

RotoPol-22. The two discs of RotoPol-22 have separate motors, each

providing two speeds. The discs can be used independently, or alternatively. The next step can be prepared in advance on the disc not operating. When RotoCom is connected the next step can be prepared on the disc not in operation, as the parameters for this step, will be sent from RotoCom as soon as the first step has started.

RotoPol-25. With a speed varying from 40-600 rpm, RotoPol-25 is the perfect choice for special applications.

300 mm single disc:

RotoPol-31 is a sturdy machine,

designed for intensive use. To obtain optimal conditions during preparation, the temperature of the preparation disc can be kept constant by means of a built-in disc cooling system.

RotoPol-35. With a variable speed of 40-600 rpm, RotoPol-35 is ideal for changing requirements when polishing large or many specimens. Disc cooling can be used with speeds up to 300 rpm.

RotoPol-11

RotoPol-22

RotoPol-35

RotoForce specimen movers

RotoForce-1

For a limited number of specimens, either as single specimens or clamped in a holder.

RotoForce-3

The high capacity specimen mover for preparation of specimens clamped in a specimen holder. RotoForce-3 can


handle specimen holders of up to 200 mm dia. where either large specimens or many small specimens can be prepared at a time.

RotoForce-4

The specimen mover for the preparation of both single specimens and specimens clamped in holders. The

touch of a key on the touch pad changes from one function to the other. Single specimens are inserted in a specimen mover plate and the force is applied individually.

If the specimens are clamped in a holder the force is applied centrally on the specimen holder.

RotoForce-1

RotoForce-3

RotoForce-4

The preparation methods from Metalog Guide are included in RotoCom

Multidoser

The dosing unit for 6 or 3 bottles. Automatic dosing, supplying constant amounts of both suspension and lubricant is necessary whenever reproducibility of preparation results is required. With Multidoser the dosing level can be programmed in twenty steps to assure correct dosing for all types of materials and all specimen sizes. Through the use of peristaltic pumps, the correct amount of liquid is applied exactly where you want it, and vaporisation or spray mist is avoided altogether.

RotoCom

- Database for 100 preparation methods with up to 10 steps each
- 10 Metalog Guide Methods included
- Complete list of all consumables included
- Data communication between the Roto Modules
- Interface to PC
- Development of preparation methods on the PC
- Easy transfer of methods to other laboratories

The memory and control unit

If different preparation methods are carried out in a lab, it is important to be able to repeat the same method again and again, using the exact same values every time. With RotoCom, every single parameter in every step for example: preparation time, force, rpm, type of suspension and lubricant, dosing level and even the preparation disc used, can be saved and recalled and run again at any time. There are

10 Metalog Guide methods included in the software. These can be used directly for the preparation of many different materials, or for fast and easy development of new, user specific preparation methods.

RotoCom is the heart of the RotoSystem

100 preparation methods with 10 steps and all parameters like speed, force, abrasive, lubricant, dosing level, the grinding/polishing surface used and the rotational direction can be saved.

The 10 methods from the Metalog Guide are included in RotoCom, which means that the actual specimen preparation can be started immediately after the RotoSystem has been set up.

RotoCom contains a library with all the consumables available from Struers. This library can easily be updated whenever new consumables are introduced or if you want to use consumables of your own.

With its software system

RotoCom takes care of all data transfer to the single Roto Modules. Simply select a preparation method, press the Auto button on RotoCom and all the preparation parameters are sent to the individual units. Force, time and rotation will be sent to the RotoForce. Type and amount of suspension and lubricant are sent to Multidoser and speed and the preparation disc to be used are sent to RotoPol. After placing the correct grinding or polishing disc on

RotoPol and inserting the specimens in RotoForce, just press start and the complete preparation step will be carried out automatically.

RotoCom, independent unit or connected to a PC

While RotoCom has been developed to work totally independently it can also be connected to a PC for even greater ease of operation. With the RotoCom PC Software™ included, some of the functions are made even easier because the PC keyboard can be used for entries. Also the number of methods to be saved is only limited by the amount of disc space on your PC.

Requirements to the PC

- IBM compatible PC with MS-DOS
- Min. 512 Kb free memory
- Optional Microsoft Windows version 3.1 or newer

The export/import feature is one of the major advantages of RotoCom PC Software™

Methods can easily be distributed from one laboratory to another, from Struers to your own laboratory or from you to your supplier, customer or any other affiliated laboratory. This will then guarantee that exactly the same method is carried out in exactly the same way, producing exactly the same results every time.

This is one of the reasons why we can talk about Maximum Reproducibility when working with a RotoSystem.

How to **design your own** System

On the previous pages you have seen three different systems covering most applications. The strength of RotoSystem is however its many different combinations enabling the tailoring of a system catering for your specific needs. On the following diagram the most important topics are mentioned, which are significant for the selection of the various modules.

First select a grinding/ polishing machine, depending on your requirements	Few specimens, small series of specimens	Series of specimens, single specimens	Many or large specimens (200mm holder size)	
	200 mm disc	2 x 250 mm disc	300 mm disc	
2 fixed speeds	RotoPol-11	RotoPol-21	RotoPol-31	150 / 300 rpm
Variable speed	RotoPol-15	RotoPol-25	RotoPol-35	40 – 600 rpm
2 individual discs		RotoPol-22		2 x 150/300 rpm
Then select the appropriate specimen mover				
		*on RotoPol-22, 2 RotoForce-1 can be mounted		
	Few specimens, small series of specimens	RotoForce-1	RotoForce-1*	150 rpm
	Single specimens, specimen holders, counter rotation		RotoForce-4	RotoForce-4
Large or heavy specimens (> 2 kg, incl. holder)		RotoForce-3	RotoForce-3	150 rpm

With the correct grinding/polishing machine and the appropriate specimen mover selected, you only need to add a Multidoser and the RotoCom to have a fully operational RotoSystem.

Expandable Systems



With RotoCom it is possible to connect two independent preparation systems to one RotoCom unit. This can be different combinations of grinding/polishing machines and specimen movers with Multidoser and RotoCom. This offers a whole new range of possibilities. In the following setup, two different ways of working are possible for continuous preparation of many specimens.

1. The preparation of different materials

Both machines are used as totally independent systems. Two different preparation methods are run simultaneously. With this configuration it may be necessary to add another Multidoser to the system. This depends on the number of suspensions and lubricants required on each system. RotoCom can control either one or two Multidoser units.

2. The preparation of many specimens of the same material

If you have two specimen holders containing the same material, for instance one holder is processed on the first step on the first system. Thereafter the second holder is processed while the first holder is taken to the second system for step two. In this case all preparation steps are carried out alternately on the two systems. In this configuration only one Multidoser will be required, dosing on both systems simultaneously.

Accessories

A complete range of accessories for automatic preparation is available

Specimen Holders


Most automatic Struers equipment is designed to work with specimens clamped in specimen holders. Depending on the size and shape of the specimens to be prepared, a large variety of different holders are available.

The program includes holders for mounted or unmounted, round, square or irregular, large or small specimens.

For the manufacturing of your own special specimen holder designs, Struers can deliver holders without holes. All specimen holders are designed using the same coupling, which makes them interchangeable from one Struers machine to another.

Accustop

Accustop is a specimen holder with a wear-resistant ceramic base. It enables controlled material removal to a pre-set depth in the specimen and ensures that the specimen surface is kept totally plane and not angled. Accustop can be used for manual grinding or be inserted in specimen mover plates.

MD-System

Struers MD-System

has revolutionized the way material-

lographic samples are prepared. A single magnetic disc, the MD-Disc, is used to support all the consumables needed for grinding and polishing. Grinding and polishing discs are easily positioned and exchanged due to the magnetic fixation, saving time and operating costs. With the MD-System, sample preparation can be carried out in only 3-5 steps thanks to the outstanding surface planeness achieved, and the extensive use of SiC paper can often be eliminated.

New Grinding Technology

As a complement to the MD-Disc, Struers has engineered a comprehensive line of metal-backed grinding and polishing discs.

The patterned segments of the grinding discs minimize build-up of abraded material, allowing consistently high material removal and optimal surface planeness. Disc flatness is maintained throughout the life of the discs.

MD-Polishing Cloths suitable for all polishing and final polishing steps are available along with adapters and accessories for storage.

For more information please ask for a separate brochure.

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Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice

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Struers Ltd.
Struers Inc.