Tenupol-3

Instruction Manual
Gebrauchsanweisung
Mode d'emploi

Spare Parts and Diagrams

Table of Contents	Page
Instruction Manual.....	1
Gebrauchsanweisung	28
Mode d'emploi.....	59
Spare Parts and Diagrams	89

Instruction Manual

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 1997.

Struers A/S
Valhøjs Allé 176
DK-2610 Rødovre/Copenhagen
Denmark
Telephone +45 36 70 35 00
Fax +45 38 27 27 01

Tenupol-3 Safety Precaution Sheet

To be read carefully before use

1. The operator must be fully instructed in the use of the instrument, according to this manual.
2. The operator must be fully instructed in the use of the electrolytes used with Tenupol-3.
3. Observe the current safety precautions regarding the handling, mixing, filling and disposal of the electrolyte.
Material Safety Data Sheets on request.
4. Check the model plate and the mains voltage selector prior to the power supply connection.
5. Check that the polishing unit is standing firmly on the support in order not to tip over.
6. Do not leave the instrument unattended when filled with electrolyte.
7. The polishing process develops heat. Keep the electrolyte temperature below the flashpoint of the actual electrolyte.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.)

Table of Contents		Page
1. Installation	4
1.1. Packing List	4
1.2. Placing	4
1.3. Mounting and Connection of the Individual Units	5
1.4. Safety Precautions	6
2. Operation	8
2.1. Introduction	8
2.2. Application	8
2.3. Function	8
2.4. Design	9
2.5. Operation of Control Unit	11
2.6. Operating the Polishing Unit	14
2.7. Operation Parameters	16
2.8. Pre-treatment of the Specimen	18
2.9. Thinning, General	19
2.10. Thinning Examples	21
2.11. Electrolytes	21
2.12. Polishing Defects	22
3. Accessories	24
4. Maintenance	25
4.1. Polishing Unit	25
4.2. Control Unit	25
5. Trouble-Shooting	26
6. Technical Data	27

1. Installation

1.1. Packing List

- 1 Control unit with cable for connection to main power supply
- 1 Polishing unit with polishing cell, cover for polishing cell, built-in IR-transmitter/receiver, cooling coil and pump system as well as cable and plug connector for control unit
- 1 Thermometer pocket
- 1 Thermometer +20/-50°C
- 1 Specimen holder for 03 mm specimens
- 1 Jet holders, inside diameter 1 mm, set with 2 pcs.
- 1 PVC-reservoir, insulated
- 1 PVC-reservoir, uninsulated
- 2 Ø8 mm tube for cooling coil, 2x1 m.
- 1 Jet holder with ascending tube
- 1 Return tube

1.2. Placing

Tenupol-3 is a stationary instrument for setting up in a laboratory. This implies that the instrument must be placed on a plane and horizontal table. The table must be resistant to chemicals. It is recommended to place the polishing unit on a piece of impregnated paper or some other easily replaceable protective foil.

Please note that the electrolytes to be used may result in demands for local exhaust or installation in fume cupboard.

The control unit can either be placed next to the polishing unit or on a shelf above the polishing unit. If placed on a table next to the polishing unit the operating convenience can be increased by placing the control unit with the built-in legs in an inclined position.

1.3. Mounting and Connection of the Individual Units

1.3.1. Polishing Unit

Place the base plate with pump, polishing cell etc. on top of the insulated PVC-reservoir. The cooling coil and the pump must be immersed in the reservoir. The pump motor and polishing cell must turn upwards (see fig. 11).

Put the thermometer pocket in the drilled hole (Fig. 2.5) between the polishing cell and the pump motor. Put the thermometer into the pocket.

Mount the two jet holders and the specimen holder as described in sections OPERATION.

Connect the polishing unit to the control unit. Connect the multi plug to the socket on the right hand side of the control unit and fasten the coupling nut of the plug. The fastening is done by hand, without use of spanners or other tools.

Connect the cooling coil to a separate cooling unit, to the water pipe, or to some other cooling arrangement which meets all technical and safety regulations as well as any possible restrictions on consumption of water resources.

Alternatively, the uninsulated vessel can be used as electrolyte reservoir. Thus, a cooling effect can be obtained if this reservoir is placed in an ice bath or another cooling bath.

1.3.2. Control Unit

On the left hand side of the instrument (Fig. 1.3) a mains voltage switch is placed. Check that it is set according to the actual supply voltage.

The power cord supplied with the instrument is mounted in the corresponding socket on the left hand side of the instrument (Fig. 1.4) and then connected to the mains.

On the lower side of the instrument, approx. 5 cm from the front edge, foldable legs (Fig. 1.5) are mounted at both sides. These legs can be used if it is required that the front plate of the control unit should be more inclined.

1.4. Safety Precautions

1.4.1. *Electrolytes*

The electrolytes mentioned in section OPERATION can be used without involving any risk for user or equipment if the safety precautions stated for these electrolytes are followed.

Some of the electrolytes contain perchloric acid. See section OPERATION for further information.

Tenupol-3 is designed for use with electrolytes recommended by Struers. Other electrolytes, e.g. electrolytes containing strong bases or acids, may harm the construction or endanger the safety of the user.

Note that many electrolytes contain alcohol or other inflammable solvents. Make sure that all safety precautions are followed when using such electrolytes.

It is essential that the user is fully trained in the use of Tenupol-3 as well as the electrolytes involved.

When cleaning the instrument after use, make sure that no electrolyte is allowed to dry and/or crystallize inside the instrument or on the polished material.

1.4.2. *Perchloric Acid*

Electrolytes from Struers with the prefix A consist of approx. 11 stock solution, to which from 15 to 90 ml perchloric acid (60%) is added. Thus, the ready-made electrolyte will contain from 1 to 5% perchloric acid.

Electrolytes with such a low content of perchloric acid will not imply considerable elements of risk. Before mixing, however, there will be 60% perchloric acid present. Therefore, the following safety precautions for perchloric acid must be observed:

- A. *Training operators*
- All personnel involved in the mixing, use, storage, transportation, and disposal of the electrolytes or its components, should be thoroughly trained in the perchloric acid precautions.
 - Great importance must be attached to precautions against inhalation of vapors of the solution or its components, against skin contact, mixing and overheating, and concerning storage and disposal methods.
- B. *Mixing the solution*
- Add the perchloric acid to the solvent/water mixture slowly, stirring continuously, in a water bath with temperature control.
 - The mixing must take place in a ventilated chemical fume hood designed for perchloric acid use.
 - The operator(s) must use the listed protective clothing or devices: full face shield or splash goggles, rubber gloves and lab coat or coveralls.
 - Avoid the use of any combustible or carbonaceous containers, reaction vessels, spill pans, storage shelves, or materials of this type when dealing with the acid.
- C. *Storage of the perchloric acid or solution*
- Do not permit any acid to crystallize on bottle necks, caps, or anywhere else.
 - Store in secure, ventilated area with metal, glass, or ceramic spill catch pan.
 - Store away from other chemicals, combustible materials, and organics.
 - Do not permit solutions to dry out under any circumstances.
- D. *Fire and explosion hazards*
- 60% perchloric acid is classified as a Material that will not burn, but caution should be taken against their coming in contact with organic material which could cause an explosion.
 - Fire fighting should be done from a protected location. Extinguish with water spray.
 - Do not produce anhydrous perchloric acid, either from its salts or from aqueous solutions, e.g. by heating with high boiling acids or dehydrating agents, such as sulfuric acid or phosphorous pentoxide. In addition to spontaneous explosion, the anhydrous acid explodes instantaneously in contact with oxidizable organic materials.
 - The use or storage of perchloric acid should be limited to quantities less than 500 g per hood.
- E. *Disposal*
- Follow local regulations when residues are to be disposed of.
 - Dilution and/or neutralization are the normally recommended methods of disposal for the electrolyte.

2. Operation

2.1. Introduction

Tenupol-3 comprises a base plate with a polishing cell, a built-in, infrared detector system, a cooling coil and a pump system, all placed on top of a PVC-reservoir for electrolyte. The electrical parts, pump etc. are connected to a separate control and power supply unit through a single cable.

2.2. Application

Tenupol-3 is an automatic instrument for electrolytic thinning of specimens for Transmission Electron Microscopy and other methods of examination where the thickness of the specimens must correspond to or be smaller than the grain size of the examined material.

2.3. Function

A specimen, normally max. 0.5 mm thick and with a diameter of 3 mm, is placed between two immersed jets. A pump system sends a flow of electrolyte through the jets against the specimen.

An electric circuit is established with a DC-power supply, with a cathode placed in the electrolyte and with the specimen connected as anode.

When the electric circuit is closed, material will be removed electrolytically from the specimen. This results in the forming of a small hole in the middle of the specimen if the correct operating parameters have been used.

The edge of the hole will have a V-shaped cross section. The point of this "V" will have the thickness required for the further structural examinations.

In Tenupol-3 the hole is detected by an infrared detector system, which automatically cuts off the process. The liquid flow rate and the voltage can be controlled from the control unit.

2.4. Design

The instrument comprises the following main parts:

2.4.1. Control Unit

A cable connects the polishing unit to the control unit. The cable contains all connections to the polishing chamber, pump and IR transmitter/receiver. The connection is done with a multi-plug on the right hand side of the control unit. The sturdy steel cabinet contains power supply, electronics and regulations for polishing voltage, photo sensitivity, flow rate and timer. Two indicating instruments indicate voltage and current respectively during the process.

It is possible to compensate for deviations of different electrolytes by setting the instrument in "Adjust Mode". In this mode Flow rate can be optimized to the electrolyte in question.

2.4.2. Polishing Unit

The polishing unit comprises a PVC base plate (Fig. 2.1) on which the pump system (Fig. 2.2), the polishing cell (Fig. 2.3) with jets and specimen holder, a cooling coil (Fig. 2.4) as well as a thermometer (Fig. 2.5) are mounted. Depending on the chosen cooling method the base plate is placed on an insulated or non-insulated PVC-reservoir.

Pump system (Fig. 3)

The electrolyte pump (Fig. 3.1) is a pressure pump made of PVC, which presses the liquid up through a Y-tube connection to the two symmetric polishing chambers in the polishing cell. The pump motor (Fig. 3.2) on top of the base plate is covered by a housing made of anodized aluminium (Fig. 3.3) and PVC (Fig. 3.4). The bearings of the motor are sealed from electrolyte-vapors and splashing by a V-ring seal (Fig. 3.5), which is mounted on the motor shaft.

The polishing cell (Fig. 4)

Electrolyte is pumped up into two symmetrical chambers (Fig. 4.1) in the polishing cell from the above mentioned connection of the Y-tube. The polishing cell is made of PVC. The electrolyte passes from the two chambers through the two jets (Fig. 4.2), hits the specimen in the specimen holder (Fig. 4.3), and returns to the reservoir.

In the middle part of the polishing cell a spring contact (Fig. 4.4) is built-in so that the specimen holder is connected to the polishing circuit.

The IR-transmitter and -receiver (Fig. 4.9) respectively are mounted in the two symmetrical end covers (Fig. 4.5). The IR-detector system is operated from the front plate of the control unit.

Jets (Fig. 5 and 6)

The two jets and jet holders are made of PVC. The standard jets have an inside diameter of 1 mm and are normally used together with the 3 mm or 2.3 mm specimen holders. The jet holders are also available with jets of Ø2.5 mm for use with 10 mm specimen holders or highly viscous electrolytes (see section ACCESSORIES).

The individual jet is wound with a helix of platinum wire which functions as a cathode. The platinum wire is connected to the socket at the side of the jet holder. When the jet holders are mounted in the polishing cell the two banana plugs are plugged into the sockets (Fig. 5.1), whereby the platinum wires are connected to the polishing circuit.

Specimen holders (Fig. 7)

The standard specimen holder which is supplied with the apparatus is made of PVC with carbon reinforced teflon insertions and is designed for thinning of specimens with diameter 3 mm and max. thickness 0.5 mm. The specimen is placed in a platinum fixture (Fig. 7.4), which by means of a stainless steel contact part (Fig. 7.5) is connected to the polishing circuit when lowered into the polishing cell.

One half part of the specimen holder carries the platinum fixture and the contact part. The other half part has a displaceable diaphragm (Fig. 7.3). The displaceable diaphragm implies that the thickness of the specimen (less than 0.5 mm) is not a critical problem when retaining the specimen in the holder.

Specimen holders for Ø2.3 mm specimens and holders with 10 mm inside diameter are available as accessories (See section ACCESSORIES).

*Cooling coil
and thermometer (Fig. 2)*

On the lower side of the base plate the corrosion resistant cooling coil (Fig. 2.4) is fixed so that the hose connectors (Fig. 2.6) are located on the rear edge of the base plate. Depending on the chosen cooling method, either hoses from a cold water tap or a separate recirculating cooling unit is connected.

In the thermometer pocket (Fig. 2.5) a thermometer -50/+20°C is placed. If an operation requires a different range, the thermometer can easily be replaced.

PVC-reservoir

The standard equipment comprises both an insulated and an uninsulated PVC-reservoir. The insulated reservoir is normally used when the electrolyte is cooled by the cooling coil. The uninsulated reservoir is used when cooling is done by an ice bath or the like.

2.5. Operation of Control Unit

When Tenupol-3 has been installed and connected as described in section INSTALLATION the instrument is switched on/off by using the main switch O/I on the top left hand side of the front plate (Fig. 10.1).

When the instrument is switched on the square indicator lamp (Fig. 10.2) is lit.

Note that the main switch together with the Start/Stop button is used to set the instrument in "Adjust Mode". The procedure is described in section OPERATION, Operation Parameters.

2.5.1. Controls

*Switch for voltage range
(Fig. 10.3)*

With the switch, it is possible to work in two voltage ranges. The low range, 5-45 V, is typically used for polishing aluminium, copper and brass. The high range, 15-120 V, is typically used for stainless steels and iron, as well as for pre-thinning and electrolytical blanking in the specimen holder with 10 mm inside diameter.

*Adjusting knob for setting
the polishing voltage (Fig. 10.4)*

When the required voltage range has been selected the polishing voltage can be set by using the double dial on the adjusting knob. Note that the voltmeter only deflects during the polishing process.

*Switch for ammeter range
(Fig. 10.6)*

The ammeter shows the current in the polishing circuit during the process. The dial range suitable for the actual process is selected with the switch, i.e. either 0-0.5 A or 0-5 A.

*Adjusting knob for
flow rate (Fig. 10.8)*

The rate of revolution of the pump motor and consequently also the flow rate of the electrolyte through the jets in the polishing chamber is controlled by the adjusting knob. The arbitrary dial from 1-10 is approximately proportional to the flow, so that 1 is lowest flow and 10 is highest flow.

*Adjusting screw for
flow rate (Fig. 10.9)*

The adjusting screw moves the minimum range for the flow rate adjusting knob. When turning the screw counter-clockwise the whole minimum range is moved downwards. Correspondingly, the whole minimum range is moved upwards when turning clockwise. The procedure is described in detail in section OPERATION, Operation Parameters.

*Adjusting knob for
photo sensitivity (Fig. 10.10)*

With this knob a potentiometric regulation of the IR-photo system sensitivity can be made. Pos. 10 on the scale gives the highest sensitivity causing the smallest hole.

Adjusting screw for photo sensitivity (Fig. 10.11)

With the adjusting screw the dial range for the photo sensitivity is moved. When the screw is turned counter clockwise the intensification of the signal from the photo cell is reduced, and reverse, when it is turned clockwise. The ratio between the luminous intensities in the two outer positions is approx. 1:4. The procedure is described in detail in section OPERATION, Operation Parameters.

Adjusting knob for timer (Fig. 10.12)

When pre-thinning and blanking electrolytically it may be an advantage to perform the polishing with time control. The process is then carried out in the preset span of time. The photo cell system will switch off the process before the set time is out (photo sensitivity higher than 0), if a hole unexpectedly develops in the specimen.

Push button for Start/stop (Fig. 10.13)

The start/stop-button is activated when the process parameters are set, just as the polishing unit must be prepared for a process performance as far as electrolyte, jets, specimen holder etc. are concerned. The process is automatically switched off, either when the IR-detector system registers a hole in the specimen, or if the preset time on the timer is out.

The process can be switched off manually by re-activating the Start/Stop-button. Switching off the process is marked by a sound signal, whether it is switched off manually or automatically.

The push-button is also used together with the main switch (Fig. 10.1) when the instrument is set in "Adjust Mode" for adjusting the flow rate. The procedure is described in section OPERATION, Operation Parameters

2.5.2. Reading of Instruments and Indicators

Voltmeter (Fig. 10.5)

The voltmeter is a double dial instrument with max. reading for 40 V and 120 V respectively. Change-over between the two dial ranges is done automatically when selecting voltage range with the switch (Fig. 10.3). The voltmeter is cut in when the polishing is started on the push button Start/Stop.

Ammeter (Fig. 10.7)

The ammeter is a dial instrument with max. deflection for 5 A. With the selector button (Fig. 10.6) the deflection range can be changed in the scale 1:10 so that full deflection corresponds to 0.5 A. Always start with the setting 5 A if the current is not known in advance in the actual process.

Start/Stop-diode (Fig. 10.14)

This diode indicates whether the IR-detector, depending on the preset photo sensitivity, registers infrared light from the transmitter. The diode can light either green or red, so that:

Red means light transmittance. I.e. there is no specimen holder in the polishing cell or, there is no specimen in the specimen holder. When the diode lights red the polishing cannot be started by activating the start button.

Green means that there is no light transmittance. I.e. the light beam is blocked by the specimen, presuming that the specimen holder is mounted correctly in the polishing cell. The polishing can be started when the diode lights green.

Polishing-diode (Fig. 10.15)

With the color red the diode indicates that the polishing has started, and with the color yellow that the instrument is in "Adjust Mode":

Red indicates that there is a current in the polishing circuit, a current which can be read on the ammeter. Note that the red reading comes a few seconds after activating the start button. Before the polishing process starts the pump works for a short time at full speed in order to fill the liquid system with electrolyte. When this is done the rate of revolution decreases to the preset rate on the "Flow rate" knob, and the polishing circuit is closed.

Note that a polishing current should be read on the ammeter. If the polishing diode lights red, while the ammeter shows 0. there is no current. The reason for this may e.g. be faulty mounting of the specimen in the specimen holder, that the specimen holder is not mounted correctly in the specimen holder, that the electrodes in the jet holders have not been correctly connected, that there is breakage and faults in the wire connections, or other faults in the polishing circuit.

*Lamp for
Thermal Overload (Fig. 10.17)*

A thermal fuse is built into the transformer. This overload fuse cuts out the transformer at 110°C. The lamp for thermal overload indicates the cutting out by lighting yellow. The fuse cuts in the transformer when it has been sufficiently cooled (normally approx. 20 min).

Automatic fuse 2 A (Fig. 10.16)

Tenupol-3 is protected against overload by an automatic fuse in the primary side of the transformer. When overloaded, the signal button of the fuse springs out. After 15-20 sec, when the fuse has become sufficiently cooled off, the signal button can be reset and the instrument is ready for use again. First, examine whether the fuse has been released by faults, which have to be repaired before starting to use the instrument again.

Automatic fuse 5 A (Fig. 13.1)

The polishing circuit is secured against overload and short circuit with a 5 A automatic fuse. The signal button is released in case of overload or short circuit. After 15-20 seconds it can be reset. First, examine whether the fuse has been released by faults, which have to be repaired before starting to use the instrument again.

Sound signal

When the polishing is stopped, either by the IR-detector, or by the timer this is marked by a built-in signal, which gives off a short audible warning signal three times.

2.6. Operating the Polishing Unit

It is assumed in the following that the mounting instructions under point 5.2.1 have been followed.

2.6.1. *Dismantling and Assembling the Specimen Holder*

To dismantle the two main parts of the specimen holder turn them a little in relation to each other. Now, they can easily be dismantled.

One half of the specimen holder (Fig. 7.1) contains a displaceable diaphragm (Fig. 7.3) This diaphragm can now be pressed out from the rear side.

When assembling the specimen holder (with specimen) the procedure is:

- The two halves are clicked together (the click is heard). The displaceable diaphragm must not be inserted.
- Without the displaceable diaphragm inserted, the specimen holder is placed on a table in such a way that the big cut-out where the diaphragm has to be inserted, turns upward. The platinum strip (Fig. 7.4) will be visible through the cut-out.
- The specimen is now placed in the aperture of the platinum strip over the diaphragm opening (Fig. 7.6) in the bottom half (Fig. 7.2).
- Then the displaceable diaphragm is pressed down into the cut-out, until it fits tightly against the specimen.
- Now the specimen holder can be mounted in the polishing cell with the contact piece facing the pump motor. Remember to press the holder completely down into the polishing cell (See Fig. 9).
- When the polishing is finished, the holder is removed from the polishing cell and placed in a beaker with water or water/alcohol mixture. Now the holder can be opened so that the specimen "falls" out into the beaker. The holder may also be taken out of the beaker and rinsed with a washing bottle. In that case, the holder is opened and the specimen taken out with forceps or the like.

2.6.2. *Mounting of Jet Holders*

The jet holders are equipped with a long sleeve which fits into a corresponding groove in the polishing cell. The holders are inserted and taken out in the longitudinal direction of the polishing cell. The procedure is as follows:

- Place one jet holder in the polishing chamber and slide the sleeve into the groove until the holder fits tightly against the back wall of the chamber.
- Mount the other jet holder similarly at the other side of the chamber. The cut-outs at the top of the holders provide the necessary space for the fingers (Fig. 6.1). Do not use force. Slide the jets at a right angle to the polishing cell aperture.
- Connect the two banana plugs in the corresponding sockets in the holders, one plug in each holder (Fig. 5.1).
- Place the protective casing over the polishing cell (see Fig. 12).

2.6.3. Filling and Emptying of Electrolyte

The volume corresponding to 1 liter is marked on the inside of the electrolyte container. There must be at least 0.7 liter in the reservoir, but it is recommended to use 1 liter.

If there is too little fluid, the pump may suck air into the liquid stream, thereby disturbing the polishing. On the other hand, there is no point in using more than one liter of fluid.

Procedure, filling

- Prepare 1 liter electrolyte according to the prescribed regulations for the electrolyte.
- Place the base plate with motor, cooling coil etc. on the electrolyte reservoir which is not being used.
- Pour the prepared electrolyte into the reservoir (insulated/uninsulated), which is to be used during the polishing processes.
- Move the base plate back to the reservoir filled with electrolyte.

Procedure, emptying

- Move the base plate from the reservoir filled with electrolyte to the reservoir which is not being used. This might be filled with water so that the first cleaning with water can start.
- If the electrolyte is to be re-used, it is poured back into its original packaging or into another resistant packaging which is duly marked. Funnel, gloves, ventilation as well as all other prescribed equipment must be used during emptying.
- If the electrolyte is not to be used again, it should be filled into a container suitable for disposal procedures. Please observe the local regulations and instructions for disposal.
- Clean the polishing unit by pumping water through the polishing cell instead of electrolyte. Make three to four successive flushings by changing from one reservoir to another. Between the changes fresh water is poured into the reservoirs.
- After the final flushing, specimen holder and jet holders are taken out and left to dry. The reservoirs and the polishing unit should be dried with a damp cloth.

2.6.4. Cooling of Electrolyte

For two main reasons it is necessary to cool the electrolyte during use, i.e.:

- In most cases the best polishing results are obtained if the electrolyte is cooled down (See section OPERATION, Electrolytes).
- For safety reasons; because of the fire hazard, evaporation of alcohol etc. it is necessary to cool the electrolyte to a temperature below its flash point.

With Tenupol-3 it is possible to cool the electrolyte by using either the built-in cooling coil or by immersing the uninsulated reservoir in an ice bath.

2.7. Operation Parameters

2.7.1. Adjust Mode

For basic adjustment of flow rate, the apparatus should be in "Adjust Mode". In this mode the normal switch-off security systems and the polishing circuit are not in operation.

Set the apparatus in "Adjust Mode" as follows:

- Replace the rear jet holder by the jet holder with the ascending tube.
- Remove the thermometer from the polishing cell.
- Place the return tube in the thermometer pocket and make sure that the ascending tube goes inside the return tube.
- Place a specimen holder, with a specimen inserted, in the polishing cell.
- Adjust the flow rate to minimum (turn the flow rate knob counter-clockwise).
- Switch off the power switch.
- With the Start/Stop-button activated, turn on the power switch.
- Release the Start/Stop-button.
- Now the motor is working with the low flow rate and the light of the polishing diode is yellow. Dependent on other factors the light of the Start/Stop-diode may be either red or green.

Disconnect "Adjust Mode" in the following way:

- Disconnect "Adjust Mode" either by pushing the Start/Stop-button or by switching off the main switch.
- Remove the specimen holder from the polishing cell.
- Replace the jet holder with ascending tube by the holder with cathode.
- Remove the return tube and insert the thermometer.

NB

Clean the jet holder with the ascending tube
and the return tube with alcohol after use.

2.7.2. Basic Setting of Flow Rate

Different electrolytes can have different viscosity and wetting qualities. Therefore the flow rate should be set at zero each time another electrolyte is to be used.

- Set the apparatus in "Adjust Mode" according to section OPERATION, Operation Parameters.
NB If the tube should contain any air bubbles, repeat the start of the apparatus.
- Set the flow rate at zero by turning the adjusting screw (fig. 10.9) with a screw driver. The zero point is found by adjusting the liquid column so that it is just visible at the passage from the jet holder to the ascending tube.
- Now the flow rate has been set at zero on a liquid column of 0 mm.
- Note the settings of the flow rate knob for liquid columns from 0-200 mm with intervals of 20 mm (measure the liquid column with a ruler).
NB If the liquid column has to be measured often, the tube can be marked with a speedmarker for every 20 mm.
- Disconnect "Adjust Mode" according to section OPERATION, Operation Parameters.
NB By going through the above procedure each time another electrolyte is used, the operator will ensure that the same flow rate will always be used during the preparation.

2.7.3. Basic Setting of Photo Sensitivity

On delivery, the gain adjusting screw for the photocell amplifier (Fig. 10.11) is adjusted so that the apparatus when holding water is switched off when a sample has a hole of 0130 μm . The photo sensitivity knob (Fig. 10.10) must be set on 2.5.

Normally no re-adjusting of the photocell amplifier is required. The amplification is reduced when turning counter clockwise, and increased when turning clockwise. A smaller amplification gives a bigger hole and vice versa. There is, however, a limit, how small the hole size can be. This lower limit should be achieved by setting the photo sensitivity on 10.

If a re-adjusting of the photocell amplifier is required the procedure is as follows:

- Choose a material and an electrolyte which notoriously give good results.
- Set the gain adjusting screw for the photocell amplifier in the middle position.
- Set the photo sensitivity on 1.
- Make a sample.
- Measure the hole of the sample (the hole must be near the center). The hole must have an average diameter of about 200 μm . If not, the photo sensitivity must be adjusted until a hole size of about 200 μm diameter can be achieved in a new sample. See section OPERATION.
- Place a sample with a hole of 200 μm diameter in the sample holder.
- Place the sample holder in the apparatus.
- Dismount the wires to the jet holders.
- Turn the gain adjusting screw for the photocell amplifier as much as possible counter clockwise.
- Set the photo sensitivity on 1.
- Start the apparatus with the Start/Stop knob.

- Turn the gain adjusting screw for the photocell amplifier clockwise until the process is switched off.

As an extra check a polishing with the photo sensitivity set on 10 can be made. The preparation of a new sample should then give a hole of 60-100 µm diameter. If not, try to re-adjust the photocell amplifier. Turning the adjusting screw counter clockwise gives a bigger hole, and turning clockwise gives a smaller hole. If a smaller hole cannot be achieved, the limit for what is obtainable has been reached.

2.8. Pre-treatment of the Specimen

The purpose of pre-treatment is to prepare discs of 3 mm (2.3 mm) diameter thickness 0.1 to 0.5 mm. Methods such as spark-machining or mechanical turning and cutting are well known.

The preparation can also be done on the Tenupol-3 (7.3.1 and 7.3.2), the original sample being a cut-off disc, thickness max. 1 mm, diameter max. 21 mm. Such a disc may e.g. be cut on the Discotom cut-off machine in the following way: a plane metal block is clamped as a stop at a distance of 0.4-2 mm from the side of the cut-off wheel. A plane face on the metal, which later is to be thinned, is pressed against the stop block and clamped. When cutting-off in the normal way, the metal between the stop block and the cut-off wheel will be parted off as a disc. The metal stop block takes care of the cooling of the disc and ensures that the cut-off wheel is not bent. The parted-off disc is now ground down on abrasive paper, e.g. on Knuth-Rotor. It is recommended to fix the disc to a plane block of metal by double-adhesive tape. It is ground first on one side, then on the other on the rotating disc, finished with grain 1000 paper and should then be max. 1 mm thick.

When the specimen has been made ready for preparation, it is important to avoid oxidation or any fouling of the specimen as this may cause the specimen not to be satisfactorily polished on Tenupol-3.

If the specimen has been punched out of a foil, it must be fineground on both sides to remove any impurities before preparation on Tenupol-3.

2.9. Thinning, General

Pour about 1 liter of the electrolyte suitable for the metal in question (see section OPERATION, Operation Parameters) into the electrolyte tank. Adjust the temperature. Place the jets in the polishing cell. Usually 1 mm jets are used for 3 mm (2.3 mm) specimens and 2.5 mm jets for the 10 mm specimen holder, but the 2.5 mm jets may also be used for the 3 mm (2.3 mm) specimens if a nearly stagnant electrolyte is used.

The specimen is placed in the holder so that it comes into contact with the platinum. The holder is closed and slid down into the cell. The voltage-source is switched on and set to the value chosen for the polishing. Also the pump speed is set to the chosen value, 1-4: slow, 4-6: middle, 6- 10: fast. The mains switch (start/stop) is closed, whereby the illumination lamp lights up. The switch "start/stop" is closed. The pump now starts. After 3 seconds the polishing chambers are filled and when the pump speed drops the polishing circuit is closed automatically. It is now possible to readjust current, voltage and pump speed. If in doubt how much electrolyte is passing the cell the user may lift the front edge of the support plate a little to watch the flow. The proper moment to stop the process depends on whether pre-thinning, electrolytic blanking or final thinning is going on. The process is terminated either by the photocell at the moment of perforation, by stop after end of timer (NB if you want smallest possible hole do not use the timer), or by switching off the motor and polishing current by the switch "start/stop". The smallest possible hole will always be achieved by turning the adjusting screw maximally clockwise and set the photosensitivity button on 10 on the scale. The specimen holder is at once drawn up and opened in a rinsing bath, whereby the specimen will fall out. If another polishing is not to be done at once, it is recommended to switch off the mains voltage as well.

2.9.1. Pre-thinning

The cut-off and possibly preground plate can be pre-thinned by means of the specimen holder providing bilateral polishing of an area of 10 mm diameter (TETMA). To fit the specimen holder, the plate must not be more than 1 mm thick and not more than 21 mm diameter. The 2.5 mm jets are used (extra equipment TETET). The polishing should not be carried on to perforation. In most cases 5 minutes will be enough. After pre-thinning the thickness should be 0.1-0.5 mm. Time is controlled with the timer.

2.9.2. “Electrolytic Blanking”

When equipped with the 2.5 mm jets the Tenupol-3 may be used for cutting out 3 mm (2,3 mm) specimens electrolytically in the 10 mm specimen holder. The larger specimen, which should be 0.1-0.5 mm thick, is first degreased in alcohol. One side of the specimen is covered completely with a special acid-resistant tape, having an acid-resistant adhesive. On the other side are stuck 1-4 discs of tape of 3 mm (2.3 mm) diameter within a circle of 10 mm diameter. The tape is available as a kit (accessory: TENKI). The tape must be pressed securely against the metal. The specimen is now put into the specimen holder, so that there is a free area, which can be polished away, around the small tape discs. When the holder has been put into the cell the cathode on the side, where the specimen is completely covered, is disconnected by pulling out the banana plug. The specimen is then polished in the usual way until the exposed area has disappeared. The 3 mm (2.3 mm) metal specimens will now be present under the tape discs. They are disengaged, washed and are ready for the final thinning.

The standard supplied jet holders have jets with an inside diameter of 1 mm. For thinning and electrolytic blanking with the 10 mm specimen holder, however, jets with an inside diameter of 2.5 mm have to be used. When using highly viscous electrolytes it may be necessary to use 2.5 mm jets (TETET) also in connection with 2.3 mm and 3 mm specimen holders. It is easy for the user to carry out the change between 1 mm and 2.5 mm jet holders.

2.9.3. Final Thinning

This will normally consist in polishing of 3 mm (2.3 mm) diameter specimens with the 1 mm jets. The polishing is carried on until a small hole appears and will be stopped as soon as the IR-light breaks through the specimen.

By varying the photo sensitivity the hole size can be changed. Position 10 implies maximum sensitivity causing a small hole to be seen at once. The process can be terminated automatically by means of the photocell which is mounted on the polishing cell. The current is broken and a beeper is activated.

As soon as the polishing is finished the specimen holder must be taken out and opened in a small bath of e.g. ethanol or distilled water, kept ready close to the cell. From this the specimen is transferred to a bath of e.g. ethanol by means of a pair of tweezers. It may then be laid on a piece of filter paper, where it dries in a couple of seconds and is ready for examination or storing.

Thinned specimens may be kept under vacuum in a desiccator with silicagel. In most cases they can also be kept in glycerol, which will protect them against the action of the atmosphere.

2.10. Thinning Examples

2.10.1 Thinning in the 10 mm Specimen Holder

Material	Electrolyte	Temperature	Flow rate mm	A	V	Approx. time for removing 0.1 mm in minutes
Brass	D-2	15	140	0.4	10	6
Stainless Steel	A-8	15	20	0.2	40	10

2.10.1. Thinning in the 3 mm Specimen Holder

Material	Electrolyte	Temperature	Flow rate mm	A	V	Approx. time for removing 0.1 mm in seconds
Brass	D-2	15	120	0.06	4	140
Stainless Steel	A-8	13	20	0.07	40	130

NB

A flow rate measured by a liquid column between 20 and 200 mm is recommended for all materials.

2.11. Electrolytes

2.11.1. Electrolytes and their Applications

A-8 Iron, Stainless Steel, Nickel, Cobalt, Zirconium. All-round electrolyte.

D-2 Copper, Copper-Alloys, Brasses.

2.11.2. Deterioration of Electrolytes

The two components are normally supplied separately so that they may be mixed directly before use. In mixed state the electrolytes have a working life of 3 months or perhaps a little more. It should be noted, however, that the use of the electrolyte is the most important cause of deterioration. The electrolytes must be replaced when they contain an excessive amount of metal ions.

2.12. Polishing Defects

2.12.1. Problem Areas

If a sufficient high quality is not obtained in electrolytic polishing it will be necessary to alter the polishing conditions. There are five important parameters:

- the type of electrolyte
- the flow rate of the electrolyte
- the temperature
- the electrical conditions and
- the geometrical conditions

In the Tenupol-3 the geometrical conditions are fixed and can only be altered by changing or altering certain components.

- As to the other parameters the following will apply:
- The chemical composition of the electrolyte is highly important for the quality of the polishing. An unsuitable electrolyte will cause inferior polishing, oxidized or etched surface, pittings or one-sided polishing, i.e. only one side of the specimen is polished, the other is black and oxidized. It is not at all certain that an electrolyte, which in other apparatus (e.g. the Polectrol) will polish a certain material, will also be suitable for thinning it on the Tenupol-3. A-2, for instance, is not suitable for use with Tenupol-3. In the literature numerous examples are cited of electrolytes suitable for thinning. (See e.g. I.S. Brammer and M.A.P. Dewey: Specimen Preparation for Electron Metallography. Blackwell Scientific Publication, Oxford 1966).
- The flow rate determines first of all whether a viscous anodic layer can be maintained during polishing. Jet polishing tends to remove the layer. The proper choice of flow rate is therefore mostly determined by the material to be polished and by the electrolyte. The best flow rate differs much from case to case and it is not possible to give general rules.
- In some cases a lower temperature will give better results. A lower temperature slows down the polishing and provides less etching and oxidation.
- The electrical conditions of course determined whether a polishing is at all obtained. The right conditions will only be present within a certain range of current densities. If a current versus voltage curve is drawn the best polishing conditions will be found where the current is approximately constant within a range of voltages. Moreover, some users have found that a better polish is obtained on some materials, if towards the end of the process the current is applied as a series of shocks.

In Tenupol-3 polishing, the polishing defects may be divided into the following classes:

1. Defective polishing, i.e. wrong polishing conditions, so that a mirror-like finish is absent on one or both sides of the specimen. Pitting may be present.
2. Polishing takes place, but no thin area is present at the edge of the hole.

Re 1:

try to alter the electrical conditions. Perhaps the voltage was too low to reach the polishing range. Pittings may arise if the current was too high. A lowering of the temperature may make the polishing less sensitive to voltage changes. The flow rate should be altered. A too high flow rate may break up the polishing layer, which may be the cause of one side of the specimen not being polished. A slower flow may allow a polishing layer to be built up.

If nothing helps it can be concluded that the material in question cannot be polished by the electrolyte and another electrolyte should be tried.

Re 2:

the hole may have grown too large, so that the thin foil first formed has disappeared. Just like electropolishing will attack the tops of the asperities of the surface, it will attack the edge of the hole as being an asperity. Especially when the jet effect is strong it is very important to stop the process while the hole is small.

Adjustment of the flow rate may be helpful. A vigorous jet may attack a too small area. Another electrical current density should also be tried.

Try increasing the photo sensitivity by turning the knob (fig. 10.10) clockwise.

Sometimes a lower temperature will give a better profile around the hole, possibly on account of a change of viscosity.

If nothing helps it is probable that the electrolyte is not suitable for thinning the material in question, although it will electropolish it. Try another electrolyte.

3. Accessories

Specification	Code
<i>Specimen holders</i> For Ø3 mm specimens	TETRI
For Ø2.3 mm specimens	TETTO
For pre-thinning, diaphragm aperture Ø10 mm	TETMA
<i>Jet holders</i> Ø1 mm (for TETRI and TETTO)	TETOR
Ø2.5mm (for TETMA)	TETET
<i>Tape Kit</i> For electrolytic blanking of 3 mm specimens	TENKI

4. Maintenance

If during the guarantee period, there should be any repairs they must be performed only by Struers or by authorized Struers' dealers. Also after expiration of the period of guarantee it is recommended to let repairs be made by Struers or Struers' dealers or other competent, skilled personnel.

4.1. Polishing Unit

Residuals of the electrolyte used previously may spoil the following polishing. Therefore it pays to rinse the polishing cell and pump thoroughly before filling in a new electrolyte. This may be done as follows: The electrolyte tank is filled with water and this is pumped through the system in the usual way. After this rinsing the pump and polishing cell should be allowed to drip off, in order that the electrolyte is not diluted with water. If copper or copper alloys have been thinned, some copper may have deposited on the cathodes. This may be removed with a few drops of nitric acid before the rinsing.

Tenupol should not be left filled with electrolyte for an unnecessarily long time, as the submerged parts may be attacked.

Especially A-8 will shorten the life of the metal parts. Never turn the polishing unit upside down, especially not when there is electrolyte in the pump, as this may trickle from the pump up to the motor and cause the axle to stick. Take care that the motor housing does not come in contact with the electrolyte. The polishing unit may be cleaned on the outside with a cloth moistened with water. The polishing unit can be flushed with a soft water spray.

4.2. Control Unit

Avoid spilling electrolyte on the cabinet or front plate of control unit. Clean the front plate with a damp cloth after use.

5. Trouble-Shooting

Electrolyte

- Check the age of the mixed electrolyte. The mixture should not be more than 3 months old.
- Check the number of polishings made with the electrolyte. The electrolyte can be worn out by too many polishings.
- Check that the correct combination of material and electrolyte is used.
- Check that the electrolyte is sufficiently cooled during operation.

Polishing circuit

- Check that all plugs are properly connected.
- Check that the stainless specimen holder contact part (Fig. 7.5) is intact and connected to the platinum strip (Fig. 7.4).
- Check that the specimen holder is connected to the polishing chamber spring contact (Fig. 4.4) when mounted in the chamber.
- Check the voltage over the spring contact and the jet holder banana plugs (or the platinum wire in the jets).

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
The drawings are not to scale. Some of the drawings may contain position numbers not used in connection with this manual.

Drawing

Fig. 1

Pos.	Description	Spare Part No.
1.1 Tenupol, complete		
2	Reservoir, uninsulated	14080125
2	Reservoir, insulated	408 MP 001
	Lid for reservoir	408 MP 019
6	Polishing cell cover	408 MP 002
	Cabinet, base	431 MP 001
	Cabinet, upper plate	431 MP 002
7	Cable with connector, 13 poles male	408 MP 051


Drawing

Fig. 2
Drawing 14080045A


Pos. Description Spare Part No.

1.2 Polishing unit, complete

	Polishing Unit, without accessories	14080001
1	Mounting plate	408 MP 006
4	Cooling coil	408 MP 008
7	Handle, 1 pc.	408 MP 005
8	Cooling coil fixture	408 MP 007
9	Y-tube	408 MP 009
	Vinyl tubes, Portex H18	408 MP 026
	ø8 x 2 x 30 mm, 1 pc.	
	ø8 x 2 x 50 mm, 2 pcs.	


(Fig. 2, Polishing unit, complete, Drawing 14080045A)

Drawing

Fig. 3
 Drawing 14080111C


Pos.	Description	Spare Part No.
1.3 Pump		
1	Pump casing w/outlet tube, glued	408 MP 017
2	Motor	392 MP 042
3	Motor screen	392 MP 049
4	Motor screen lid	392 MP 050
5	V-ring for pump motor, VS 008, 1 pc.	2IV00008
6	Motor flange	392 MP 044
7	Gasket flange	392 MP 047
8	Clamp flange and O-ring, 5.3 x 2.4	392 MP 045
9	Sling disc	392 MP 046
10	Pump axle, mounted	392 MP 048
11	Axle house	392 MP 035
12	Vane wheel + vane roller, 3 pcs.	392 MP 039
13	Glass disc	392 MP 037
14	Glass disc spring	392 MP 038

Drawing

Fig. 4
 Drawing 14080070B


Pos.	Description	Spare Part No.
1.4	Polishing cell	
	Polishing cell, complete	408 MP 010
4	Connecting spring, with cable	408 MP 014
6	Foot for polishing cell	408 MP 011
7	Inlet tubes, 2 pcs.	408 MP 012
8	Outlet tubes, 2 pcs.	408 MP 013
9	IR-transmitter CT 10 G-2	408 MP 020
9	IR-receiver CR 10 G-2	408 MP 021
10/11	Set of sealings (VITON)	408 MP 027
	O-ring for transm./receiv., 9.3x2.4 mm	
	O-ring for end cover, 49.5x3 mm	

Drawing

Fig. 5

Pos.	Description	Spare Part No.
1.5	Connection of the jets	
1	Angle plug w/2 m black cable	408 MP 015
	Angle plug spring	408 MP 016


Drawing

Pos. Description

Spare Part No.

Fig. 6

1.6 Placing the jets


Drawing

Fig. 7

Pos.	Description	Spare Part No.
1.7	Specimen holder	
	Specimen holder, complete $\varnothing 0$	04086907
	Specimen holder, complete $\varnothing 2.3$ (TETTO)	04086902
	Specimen holder, complete $\varnothing 3$ (TETRI)	04086901
	Specimen holder, complete $\varnothing 10$ (TETMA)	04086903
1	Specimen holder plate	14080146
2	Specimen holder plate, contact side	14080145
3	Displaceable Insert, $\varnothing 0$	14080200
3	Displaceable Insert, $\varnothing 2.3$	14080170
3	Displaceable Insert, $\varnothing 3$	14080148
3	Displaceable Insert, $\varnothing 10$	14080180
*4	Platinum Contact, $\varnothing 0$, $\varnothing 3$	12550460
*4	Platinum Contact, $\varnothing 2.3$	12550461
*4	Platinum Contact, $\varnothing 10$	12550073
5	Slide Contact, A4	14080142
6	Insert, contact side, $\varnothing 0$	14080149
6	Insert, contact side, $\varnothing 2.3$	14080147
6	Insert, contact side, $\varnothing 3$	14080143
6	Insert, contact side, $\varnothing 10$	14080150
7	Spring for platinum contact	12550059

*** Please Note!**


These contacts are made of pure platinum, which is a rare and valuable metal. Please recycle!

Drawing

Fig. 8

Pos.	Description	Spare Part No.
1.8	Specimen holder, seen from the side	


Drawing

Fig. 9

Pos. Description

Spare Part No.

1.9 Placing the specimen holder


Drawing

Fig. 10

Pos.	Description	Spare Part No.
1.10 Front plate, controls		
	Front plate foil	14310055
1	Switch ZW X11/2840	404 MP 024
2	Indication light 169523013.1002, 24 V, white	431 MP 016
3	Switch APR 646/2E Rubber jacket U 851	213 MP 002
4/6/8	ELMA button 020.4425,	
10/12	ELMA wing disc 041.4025, ELMA lid 040.4025)	392 MP 005
13	Push button distance bush	431 MP 008
13	Push button cap, U 224	391 MP 006
16	Thermo fuse, 2 A	427 MP 081B
17	Indication light, 24 V, yellow	431 MP 017

(Fig. 10, Front plate, controls)


Drawing

Fig. 11

Pos. Description

Spare Part No.

1.11 Placing the polishing unit


Drawing

Fig. 12

Pos. Description Spare Part No.

1.12 Placing the polishing cell cover

1	Wire cover	408 MP 003
2	Thermometer	255 MP 056
3	Thermometer pocket	408 MP 004
	Jet holder with ascending tube	408 MP 022
	Return tube	408 MP 023


Drawing

Fig. 13
Drawing 14310070

Pos.	Description	Spare Part No.
1.13 Cabinet		
2	Mounting plate for plug	431 MP 007
5	Apparatus plug STF63A3	130 MP 017
3	Mains voltage selector plug, TS 8565 w/serigraphy	
4	Mains-voltage selector socket 6168TJ	431 MP 020
	Relay print board A3	431 MP 021
6	Connector, 13 poles, female	2XM00720
6	Connector, 13 poles, male	2XM10717
A2	Power board, complete	431 MP 018
F4	Aut. fuse W28XQ1A, 5 A	431 MP 011
L2	Choke coil L2/10802	431 MP 012
T1	Transformer 11183/2	431 MP 010
Z1	Noise filter	431 MP 009

(Fig. 13, Cabinet, Drawing 14310070)


Drawing

Fig. 14
Drawing 14310008C

Pos.	Description	Spare Part No.
1.14 Front plate, control board		
1	Front plate	431 MP 004
2	Indication lamp bracket	431 MP 005
3	Voltage selector	431 MP 006
A1	Control circuit board 205P331011	431 MP 019
A1	Replacement kit for s/n up to 4310217	14310100
A1	Pot. Meter (500 Ω /2W R50)	14312901
A1	Pot. Meter (10k Ω /2W, Contelec R6, R32, R66)	R4310081
P1	Ammeter (grey scale)	14310054
P2	Voltmeter (grey scale)	14310058

(Fig. 14, Front plate, control board, Drawing 14310008C)

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












Struers A/S
Valhøjs Allé 176
DK-2610 Rødovre, Denmark
Telephone 36 70 35 00

Herewith declares that

<i>Product Name:</i>	Tenupol-3
<i>Type No:</i>	431, 408
<i>Machine Type:</i>	Electrolytic polishing apparatus

is in conformity with the provisions of the following directives:

EMC-Directive 89/336/EEC and 92/31/EEC according to the following standard(s):
EN50081-1:1991, EN50082-1:1992.

Low Voltage Directive 73/23/EEC and 93/68/EEC according to the following standard(s):
EN60204-1:1992.

Supplementary Information The equipment complies with the American standards:
FCC part 15, subpart B and UL508.

The above has been declared according to the global method, module A

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Klaus Kisbøll, Entwicklungsleiter, Struers A/S
Klaus Kisbøll, Chef du service developpement, Struers A/S