SELECTION GUIDE FOR GRINDING AND POLISHING

Plane Grinding

Disc	Application areas	Hardness range	FEPA P grit size equivalent	ANSI (US) grit size equivalent	Grain size equivalent (µm)	Abrasive (embedded)	Abrasive Bond	Indicative lifetime*	Pre-dressed/ Ready-to-use
MD-Primo 120 220	Non-ferrous metals + soft materials -	40 - 250 HV -	120 220	120 220	125 68	SiC -	Resin -	50 -	Yes -
MD-Molto 220	Aluminium alloys, and hard materials containing aluminium	50 - 2000 HV	220	220	68	Diamond	Resin	100	Yes
MD-Mezzo 220	Titanium alloys	150 - 450 HV	220	220	68	Diamond	Resin	100	Yes
MD-Piano 80 120 220	Ferrous metals + hard materials - -	150 - 2000 HV - -	80 120 220	80 120 220	200 125 68	Diamond - -	Resin - -	100 - -	Yes - -


Fine Grinding

Disc	Application areas	Hardness range	FEPA P grit size equivalent	ANSI (US) grit size equivalent	Grain size equivalent (µm)	Abrasive (embedded)	Abrasive Bond	Indicative lifetime*	Pre-dressed/ Ready-to-use
MD-Piano 500 1200 2000 4000	Ferrous metals + hard materials - - -	150-2000 HV - - -	500 1200 2000 4000**	360 600 800 1000**	30 15 10 5	Diamond - - -	Resin - - -	100 - 300 -	Yes - - -

Disc	Application areas	Hardness range	Abrasive (added)	Abrasive (µm)	Indicative lifetime*	Pre-dressed/ Ready-to-use
MD-Allegro	Fine grinding of materials harder than 150 HV	>150 HV	Suspension/ spray	15-6	150	Yes
MD-Largo	Soft materials, composites with soft matrix	40-250 HV	-	9-3	-	-

Polishing cloth - please see separate brochure for further details

Disc	Application areas	Resilience	Hardness range	Abrasive (µm)	Material
MD-Plan	Fine grinding of soft material Pre-polishing of hard materials	Very low	Hard	15-3	Coated, woven polyester
MD-Pan	Fine grinding of soft metals. Pre-polishing and polishing of hard and brittle materials	Very low	Hard	15-1	Impregn. non-woven technical textile
MD-Sat	Fine grinding and polishing of ferrous metals, non-ferrous metals, coatings and plastics	Medium	Hard	9-3	Woven acetate
MD-Dur	Fine grinding and polishing of ferrous metals, non-ferrous metals, coatings and plastics	Medium	Hard	9-1	Satin woven natural silk
MD-Dac	Polishing of all materials	Medium	Hard	6-3	Satin woven acetate
MD-Mol MD-Mol APS	Polishing of ferrous and non-ferrous metals and polymers	High	Soft	≤ 3	Taffeta woven 100 % wool APS for Automatic Preparation Systems
MD-Plus	One step polishing for sintered carbides and steels	High	Soft	≤ 3	Synthetic nap
MD-Floc	Polishing of all materials	High	Soft	≤ 3	Synthetic nap
MD-Nap	Final polishing of all materials	Very high	Very soft	≤ 1	Synthetic short nap
MD-Chem MD-Chem NonStick	Final polishing of all materials	High	Soft	<1	Porous neoprene NonStick for large specimens to avoid suction


Maintenance

The MD-Consumables are practically maintenance-free which makes them very easy to use.


MD-Primo	A brief dressing of the SiC layer is necessary from time to time
MD-Molto / MD-Mezzo	A brief dressing of the diamond layer can be required from time to time
MD-Piano	A brief dressing of the diamond layer is necessary from time to time
MD-Allegro	Are virtually maintenance-free though occasional cleaning is recommended
MD-Largo	Are virtually maintenance-free though occasional cleaning is recommended
MD-Cloths	No maintenance required

* Indicative values corresponding to preparation of number of specimen holders, 6 x ø 30 mm. Actual lifetime depends on application (prepared material, size, shape, number etc).


** Not part of FEPA P or ANSI standardization

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



CANADA
Struers Ltd.
7275 West Credit Avenue
Mississauga, Ontario L5N 5M9
Phone +1 905-814-8855
Fax +1 905-814-1440
info@struers.com

CHINA
Struers Ltd.
No. 1696 Zhang Heng Road
Zhang Jiang Hi-Tech Park
Shanghai 201203, P.R. China
Phone +86 (21) 6035 3900
Fax +86 (21) 6035 3999
struers@struers.cn

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For more information, please see Polishing Cloth brochure.

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Struers Inc.