

Technical Data

Power Supply

50/60Hz (max. load 9.8A)	1 x 100V / 120V
50/60Hz (max. load 4.9A)	1 x 220V / 240V

Output voltage/current

Polishing	0 - 100V	(1V steps) / 6A
Etching	0 - 25V	(0.5 V steps) / 6A
External Etching	0 - 15V	(0.5 V steps) / 1.5A

Software and Electronics

Display	128 x 240 dots (16 x 40 characters)
Controls	Touch pad
Database	10 Struers methods + 20 user definable methods (non volatile)

Dimensions and Weight

	Control unit	Polishing unit
Width	385 mm / 15.2"	220 mm / 8.7"
Depth	350 mm / 13.8"	350 mm / 13.8"
Height	160 mm / 6.3"	160 mm / 6.3"
Weight	18 kg / 40 lbs	6.3 kg / 14 lbs

Specifications

LectroPol-5

Fully automatic, microprocessor controlled electrolytic polishing and etching apparatus. Complete with control unit, polishing unit (04936201) and external etching unit (02250121). One set of masks (04936915) and two electrolyte containers (04936902) included.

	Cat. no.
Control unit with external etching unit	04936233
Polishing unit with one set of masks (04936915) and two electrolyte containers (04936902)	04936201
Polishing unit for low temperature, with one set of masks (04936915), one cathode (04936908) and two electrolyte containers (04936902). External cooling unit and tubing for connection to the external cooling unit are not included and have to be sourced alternatively. Internal dia. of tubing to be used, 6 mm	04936202
Extension Arm, for the preparation of specimens up to 100 mm high. Can be mounted on pol. unit 04936202 and on pol.unit 04936201 from serial no. 4935250	04936909
1 set of masks (0,5, 1, 2, 5 cm ² and 1 without aperture.)	04936915
1 set of 5 masks 0,5 cm ²	04936916
1 set of 5 masks 1 cm ²	04936917
1 set of 5 masks 2 cm ²	04936918
1 set of 5 masks 5 cm ²	04936919
1 set of 5 masks without aperture	04936920
Electrolyte container with lid	04936902
External etching unit	02250121

List of Electrolytes

A2	A3	D2
Aluminium	Martensitic Steel	Brass
Iron	Molybdenum	Copper
Nickel	Stainless Steel	Gold
Silver	Titanium	
Stainless Steel	Vanadium	
Steels	Zirconium	
Tin		
Titanium		

Struers' equipment is in conformity with the provisions of the applicable International Directives and their appurtenant Standards. (Please contact your local supplier for details)

Struers' products are subject to constant product development. Therefore, we reserve ourselves the right to introduce changes in our products without notice.

Struers A/S
 Pederstrupvej 84
 DK-2750 Ballerup, Denmark
 Phone +45 44 600 800
 Fax +45 44 600 801
 struers@struers.dk
 www.struers.com

AUSTRALIA & NEW ZEALAND

Struers Australia
 27 Mayneview Street
 Milton QLD 4064
 Australia
 Phone +61 7 3512 9600
 Fax +61 7 3369 8200
 info.au@struers.dk

BELGIQUE (Wallonie)

Struers S.A.S.
 370, rue du Marché Rollay
 F-94507 Champigny
 sur Marne Cedex
 Téléphone +33 1 5509 1430
 Télécopie +33 1 5509 1449
 struers@struers.fr

BELGIUM (Flanders)

Struers GmbH Nederland
 Elekraweg 5
 3144 CB Maassluis
 Telefoon +31 (10) 599 7209
 Fax +31 (10) 5997201
 netherlands@struers.de

CANADA

Struers Ltd.
 7275 West Credit Avenue
 Mississauga, Ontario L5N 5M9
 Phone +1 905-814-8855
 Fax +1 905-814-1440
 info@struers.com

CHINA

Struers Ltd.
 No. 1696 Zhang Heng Road
 Zhang Jiang Hi-Tech Park
 Shanghai 201203, P.R. China
 Phone +86 (21) 6035 3900
 Fax +86 (21) 6035 3999
 struers@struers.cn

CZECH REPUBLIC

Struers GmbH
 Organizační složka
 Havlíčkova 361
 CZ-252 63 Roztoky u Prahy
 Phone +420 233 312 625
 Fax +420 233 312 640
 czechrepublic@struers.de

DEUTSCHLAND

Struers GmbH
 Carl-Friedrich-Benz-Straße 5
 D-47877 Willich
 Telefon +49 (0) 2154 486-0
 Fax +49 (0) 2154 486-222
 verkauf@struers.de

FRANCE

Struers S.A.S.
 370, rue du Marché Rollay
 F-94507 Champigny
 sur Marne Cedex
 Téléphone +33 1 5509 1430
 Télécopie +33 1 5509 1449
 struers@struers.fr

HUNGARY

Struers GmbH
 Magyarországi Fióktelep
 Tatai ut 53
 2821 Gyermely
 Phone +36 (34) 880546
 Fax +36 (34) 880547
 hungary@struers.de

IRELAND

Struers Ltd.
 Unit 11 Evolution @ AMP
 Whittle Way, Catcliffe
 Rotherham S60 5BL
 Tel. +44 0845 604 6664
 Fax +44 0845 604 6651
 info@struers.co.uk

ITALY

Struers Italia
 Via Monte Grappa 80/4
 20020 Aresè (MI)
 Tel. +39-02/38236281
 Fax +39-02/38236274
 struers.it@struers.it

JAPAN

Marumoto Struers K.K.
 Takara 3rd Building
 18-6, Higashi Ueno 1-chome
 Taiko-ku, Tokyo 110-0015
 Phone +81 3 5688 2914
 Fax +81 3 5688 2927
 struers@struers.co.jp

NETHERLANDS

Struers GmbH Nederland
 Elekraweg 5
 3144 CB Maassluis
 Telefoon +31 (10) 599 7209
 Fax +31 (10) 5997201
 netherlands@struers.de

NORWAY

Struers A/S, Norge
 Sjoskogenveien 44C
 1407 Vinterbro
 Telefon +47 970 94 285
 info@struers.no

ÖSTERREICH

Struers GmbH
 Zweigniederlassung Österreich
 Betriebsgebiet Puch Nord 8
 5412 Puch
 Telefon +43 6245 70567
 Fax +43 6245 70567-78
 austria@struers.de

POLAND

Struers Sp. z o.o.
 Oddział w Polsce
 ul. Jasnogórska 44
 31-358 Kraków
 Phone +48 12 661 20 60
 Fax +48 12 626 01 46
 poland@struers.de

ROMANIA

Struers GmbH
 Sucursala Sibiu
 Str.Scoala de Inot, nr. 18
 RO-550005 Sibiu
 Phone +40 269 244 558
 Fax +40 269 244 559
 romania@struers.de

SCHWEIZ

Struers GmbH
 Zweigniederlassung Schweiz
 Weissenbrunnstraße 41
 CH-8903 Birmensdorf
 Telefon +41 44 777 63 07
 Fax +41 44 777 63 09
 switzerland@struers.de

SINGAPORE

Struers Singapore
 627A Aljunied Road,
 #07-08 BizTech Centre
 Singapore 389842
 Phone +65 6299 2268
 Fax +65 6299 2661
 struers.sg@struers.dk

SPAIN

Struers A/S, Spain
 Camino Cerro de los Gamos 1
 Building 1 - Pozuelo de Alarcón
 CP 28224 Madrid
 Teléfono +34 917 901 204
 Fax +34 917 901 112
 struers.es@struers.es

SUOMI

Struers Suomi
 Kalevankatu 43
 00180 Helsinki
 Puhelin +358 (0)207 919 430
 Faksi +358 (0)207 919 431
 finland@struers.fi

SWEDEN

Struers Sverige
 Ekbacksvägen 22
 168 69 Bromma
 Telefon +46 (0)8 447 53 90
 Telefax +46 (0)8 447 53 99
 info@struers.se

UNITED KINGDOM

Struers Ltd.
 Unit 11 Evolution @ AMP
 Whittle Way, Catcliffe
 Rotherham S60 5BL
 Tel. +44 0845 604 6664
 Fax +44 0845 604 6651
 info@struers.co.uk

USA

Struers Inc.
 24766 Detroit Road
 Westlake, OH 44145-1598
 Phone +1 440 871 0071
 Fax +1 440 871 8188
 info@struers.com