



TECHNICAL DATA
DiaDuo-2 All-in-One diamond solutions

	DiaDuo-2	Grain size	Diamonds	Water or alcohol based	Container size
	DiaDuo-2	15 µm	Optimized with Monocrystalline	Water	500 ml 5 l
	DiaDuo-2	9 µm	Optimized with Monocrystalline	Water	500 ml 5 l
	DiaDuo-2	6 µm	Optimized with Monocrystalline	Water	500 ml 5 l
	DiaDuo-2	3 µm	Optimized with Monocrystalline	Water	500 ml 5 l
	DiaDuo-2	1 µm	Optimized with Monocrystalline	Water	500 ml 5 l

Find Safety Data Sheets (SDSs) on [struers.com](https://www.struers.com)



DiaDuo-2 All-in-One diamond solutions
that combine abrasive and lubricant