



CitoPress – **Increase** your **Capacity**



The new CitoPress increases your capacity. CitoPress gives you high speed and maximum user-friendliness, which means faster turnaround and increased productive capacity.


### The fastest way to mounting

We are pleased to present the fastest hot mounting press in the market. CitoPress offers more power and ultra-short heating and cooling times, reducing mounting times by half for some resins, compared to competing presses. Speed is further improved by using the optional dosing system. The unique automatic dosing system enables dosing of a pre-set amount of mounting resin. The advantages are faster filling of the cylinders and less spilling of resins.

### ...and increased productive capacity

The large display and the turn/push-knob facilitate a fast selection of mounting parameters. Intuitive operation secures faster turnaround, a minimum of errors and shortest possible set-up time. We know that turnaround time is one of the main criteria of our customers. The end result is nothing less than the world's fastest mounting press and thereby offering increased productive capacity!

### CitoPress mounting times for different resins, compared to competing presses – competing presses = 100%



*CitoPress can reduce mounting times by half, compared to competing presses. (applies for 40 mm dia. mounts).*

# CitoPress the **Fastest Press** in the **Market**

### On-screen Hot Mounting Application Guide

Everything you need to know about the resins and their usage – right at your fingertips! Select the Struers resin to be used, and CitoPress automatically checks the cylinder diameter used, adapts the appropriate method, after which it is displayed – ready for use.

### Available in 3 versions

The new CitoPress comes in 3 versions: CitoPress-20, dual cylinder, CitoPress-10, single cylinder, CitoPress-1, single cylinder and all manual operation.

# CitoPress – Hot Mounting Presses

## Automatic dosing

The unique automatic dosing system enables fast dosing of a pre-set amount of mounting resin.

## Top closure elevator

All models come with a top closure elevator, for easy positioning and fast handling of top closure.

## Less spilling

Fast and easy dosing means saving valuable time and less spilling of resin.

## Intuitive operation

Push a single key to activate the automatic dosing system.



## Remote signalling (CitoPress-10/-20 only)

No time wasted checking on the press: A flashing LED bar indicates the end of the mounting process.



## Low-mass mounting units

The compact, low-mass mounting units have several advantages. Since there is less mass to be heated, the heating time is significantly reduced.

**Available cylinder diameters:** 25, 30, 40, 50 mm, 1 1/2", 1 1/4"





**Struers A/S**  
Pederstrupvej 84  
DK-2750 Ballerup, Denmark  
Phone +45 44 600 800  
Fax +45 44 600 801  
struers@struers.dk  
www.struers.com



**DEUTSCHLAND**  
**Struers GmbH**  
Karl-Arnold-Strasse 13 B  
D-47877 Willich  
Telefon +49(02154) 486-0  
Telefax +49(02154) 486-222  
verkauf.struers@struers.de

**ÖSTERREICH**  
**Struers GmbH**  
Zweigniederlassung Österreich  
Ginzkeyplatz 10  
A-5020 Salzburg  
Telefon +43 662 625 711  
Telefax +43 662 625 711 78  
stefan.lintschinger@struers.de

**SCHWEIZ**  
**Struers GmbH**  
Zweigniederlassung Schweiz  
Weissenbrunnenstrasse 41  
CH-8903 Birmensdorf  
Telefon +41 44 777 63 07  
Telefax +41 44 777 63 09  
rudolf.weber@struers.de

**CZECH REPUBLIC**  
**Struers GmbH**  
Organizační složka  
Havlíčkova 361  
CZ-252 63 Roztoky u Prahy  
Tel: +420 233 312 625  
Fax: +420 233 312 640  
david.cernicky@struers.de

**POLAND**  
**Struers Sp. z o.o.**  
Oddział w Polsce  
ul. Lirowa 27  
PL-02-387 Warszawa  
Tel. +48 22 824 52 80  
Fax +48 22 882 06 43  
grzegorz.uszynski@struers.de

**HUNGARY**  
**Struers GmbH**  
Magyarországi fióktelep  
Puskás Tivadar u. 4  
H-2040 Budaörs  
Phone +36 (23) 428-742  
Fax +36 (23) 428-741  
zoltan.kiss@struers.de

**SINGAPORE**  
**Struers A/S**  
10 Eunos Road 8,  
#12-06 North Lobby  
Singapore Post Centre  
Singapore 408600  
Phone +65 6299 2268  
Fax +65 6299 2661  
struers.sg@struers.dk

**USA and CANADA**  
**Struers Inc.**  
24766 Detroit Road  
Westlake, OH 44145-1598  
Phone +1 440 871 0071  
Fax +1 440 871 8188  
info@struers.com

**SWEDEN**  
**Struers A/S**  
Smältvägen 1  
P.O. Box 11085  
SE-161 11 Bromma  
Telefon +46 (0)8 447 53 90  
Telefax +46 (0)8 447 53 99  
info@struers.dk

**FRANCE**  
**Struers S.A.S.**  
370, rue du Marché Rollay  
F- 94507 Champigny  
sur Marne Cedex  
Téléphone +33 1 5509 1430  
Télécopie +33 1 5509 1449  
struers@struers.fr

**NEDERLAND/BELGIE**  
**Struers GmbH Nederland**  
Electraweg 5  
NL-3144 CB Maassluis  
Tel. +31 (0) 10 599 72 09  
Fax +31 (0) 10 599 72 01  
glen.vugt@struers.de

**BELGIQUE (Wallonie)**  
**Struers S.A.S.**  
370, rue du Marché Rollay  
F- 94507 Champigny  
sur Marne Cedex  
Téléphone +33 1 5509 1430  
Télécopie +33 1 5509 1449  
struers@struers.fr

**UNITED KINGDOM**  
**Struers Ltd.**  
Unit 25a  
Monkspath Business Park  
Sollihull  
B90 4NZ  
Phone +44 0121 745 8200  
Fax +44 0121 733 6450  
info@struers.co.uk

**JAPAN**  
**Marumoto Struers K.K.**  
Takara 3rd Building  
18-6, Higashi Ueno 1-chome  
Taito-ku, Tokyo 110-0015  
Phone +81 3 5688 2914  
Fax +81 3 5688 2927  
struers@struers.co.jp

**CHINA**  
**Struers Ltd.**  
Office 702 Hi-Shanghai  
No. 970 Dalian Road  
Shanghai 200092, P.R. China  
Phone +86 (21) 5228 8811  
Fax +86 (21) 5228 8821  
struers.cn@struers.dk