

Product Data Sheet



Struers ApS
Pederstrupvej 84
2750 Ballerup, Denmark

Phone +45 44 600 800
Fax +45 44 600 803
struers@struers.dk
www.struers.com

Product Name **Gekko PSA**

Cat. no. 49900053 Gekko PSA, 200 mm (8") dia.
 49900054 Gekko PSA, 250 mm (10") dia.
 49900055 Gekko PSA, 300 mm (12") dia.

Revision 24-01-2017

Description Self adhesive Gekko Foil for use with SiC Foil or self adhesive consumables. To be glued on an aluminium disc.

Material Round sheets of self adhesive foil with high friction surface.

Properties Equipped with a self adhesive backing for secure fixation on an aluminium disc.
Re-usable.

Application Fixation of polishing cloth and SiC foil. The user must follow the Instructions for Use for correct work procedure.

Health and Safety - Whilst the Gekko PSA is rotating, ensure hands are kept well clear of the edge and out of the splash bowl.

 - Do not exceed a rotational speed of 600 rpm.

To the best of Struers' knowledge, the products listed above do not contain any Substances of Very High Concern (SVHC) according to REACH regulations (amended on 20th June 2016) at or above 0.1%.

According to OSHA's Hazard Communication Standard (HCS 2012) 29 CFR 1910.1200 paragraph (c) this product is considered an article and defined as *"a manufactured item other than a fluid or particle: (i) which is formed to a specific shape or design during manufacture; (ii) which has end use function(s) dependent in whole or in part upon its shape or design during end use; and (iii) which under normal conditions of use does not release more than very small quantities, e.g., minute or trace amounts of a hazardous chemical (as determined under paragraph (d) of this section), and does not pose a physical hazard or health risk to employees"*.

Therefore this product is exempted from HCS 2012 obligations on classification and labelling under stated normal conditions of use (29 CFR 1910.1200 (b)(6)(v)).


See the Instructions for Use for more information.

Skin protection Wear protective gloves during use.

Respiratory Protection It is recommended that the grinding machine is adequately ventilated during use of Gekko PSA.

Environmental precautions Dispose of waste and residues as metal waste or landfill in accordance with local regulations.

Approved by:



Krista Huggins

Chemical Safety Specialist
Struers ApS.