

Product Data Sheet



Struers ApS
Pederstrupvej 84
2750 Ballerup, Denmark

Phone +45 44 600 800
Fax +45 44 600 803
struers@struers.dk
www.struers.com

Product Name **MD/DP-Sat**

Cat. no.	40500406	MD-Sat, 200 mm (8") dia.
	40500407	MD-Sat, 250 mm (10") dia.
	40500408	MD-Sat, 300 mm (12") dia.
	40500216	DP-Sat, 200 mm (8") dia.
	40500217	DP-Sat, 250 mm (10") dia.
	40500218	DP-Sat, 300 mm (12") dia.

Revision 21-03-2017

Description Polishing cloth with stiff support with a very high degree of planeness.

Material Woven acetate with an adhesive backing / on a metal plate.

Properties DP-Sat to be applied on a preparation disc. DP-Sat is self adhesive.
MD-Sat for magnetic fixation on MD-Disc. Black screen printing on the rear.

Re-usable.

Application Polishing cloth for fine grinding and polishing of ferrous metals, non-ferrous metals, coatings and plastics. To be used with Struers diamond suspension with grain size 9 – 3 µm.

Health and Safety

- Whilst the MD/DP-Sat is rotating, ensure hands are kept well clear of the edge and out of the splash bowl.
- Do not exceed a rotational speed of 600 rpm

To the best of Struers' knowledge, the products listed above do not contain any Substances of Very High Concern (SVHC) according to REACH regulations (amended on 20th June 2016) above 0.1%.


According to Article 3(3) of the REACH Regulation (EC) No 1907/2006, this product is considered an article and is defined as *"an object which during production is given a special shape, surface or design which determines its function to a greater degree than its chemical composition"*.

Therefore this product is exempted from CLP obligations and is outside the scope of the UN Globally Harmonized System of Classification and Labelling of Chemicals (GHS).

See Instructions for Use for more information.

Skin protection	Wear protective gloves during polishing.
Respiratory Protection	It is recommended that the polishing machine is adequately ventilated during use of MD/DP-Sat.
Personal Protection	No requirements.
Environmental precautions	Dispose of waste and residues as metal waste in accordance with local authority requirements.

Approved by:



Krista Huggins

Chemical Safety Specialist
Struers ApS.