

# Product Data Sheet



**Struers ApS**  
Pederstrupvej 84  
2750 Ballerup, Denmark

Phone +45 44 600 800  
Fax +45 44 600 803  
struers@struers.dk  
www.struers.com

**Product Name** MD-Fuga

**Cat. no.** 49900021 MD-Fuga, 200 mm (8") dia.  
49900022 MD-Fuga, 250 mm (10") dia.  
49900023 MD-Fuga, 300 mm (12") dia.

**Revision** 21-03-2017

---

**Description** Metal disc with an adhesive layer for holding plain grinding papers (silicon carbide papers). Can be used repeatedly.

**Material** Transfer tape with protective liner, laid on epoxy-coated metal plate.

**Properties** Equipped with special, anti-slip coating on the back.  
For magnetic fixation on MD-Disc. White screen printing on the rear.  
Re-usable.

**Application** For fixation of silicon carbide grinding papers.  
  
The user must follow the Instructions for Use for correct work procedure.

**Health and Safety**

- Whilst the MD-Fuga is rotating, ensure hands are kept well clear of the edge and out of the splash bowl.
- Do not exceed a rotational speed of 600 rpm.

To the best of Struers' knowledge, the products listed above do not contain any Substances of Very High Concern (SVHC) according to REACH regulations (amended on 20th June 2016) at or above 0.1%.

According to Article 3(3) of the REACH Regulation (EC) No 1907/2006, this product is considered an article and is defined as *"an object which during production is given a special shape, surface or design which determines its function to a greater degree than its chemical composition"*.

Therefore this product is exempted from CLP obligations and is outside the scope of the UN Globally Harmonized System of Classification and Labelling of Chemicals (GHS).


See Instructions for Use for more information.

**Respiratory Protection** It is recommended that the grinding machine is adequately ventilated during use of MD-Fuga.

**Environmental  
precautions**

Dispose of waste and residues as metal waste in accordance with local regulations.

Approved by:

A handwritten signature in purple ink that reads "Krista Huggins".

Krista Huggins

Chemical Safety Specialist  
Struers ApS.