

Product Data Sheet



Struers ApS
Pederstrupvej 84
2750 Ballerup, Denmark

Phone +45 44 600 800
Fax +45 44 600 803
struers@struers.dk
www.struers.com

Product Name MD-Rondo

Cat. no. 40503000 MD-Rondo, 200 mm (8") dia.
40503001 MD-Rondo, 250 mm (10") dia.
40503002 MD-Rondo, 300 mm (12") dia.
40503082 MD-Rondo, 350 mm (14") dia.

Revision 11-01-2018

Description Adapter disc for application of self-adhesive consumables to MD-Disc.

Material Tinned metal sheet with lacquered patterned surface and high friction back.

Properties Adapter surface to allow easy fixation and removal of self-adhesive grinding/polishing products.

Re-usable.

Application Fixation of self-adhesive backed grinding/polishing consumables

Health and Safety - Whilst the disc is rotating, keep hands well clear of the edge and out of the splash bowl.

To the best of Struers' knowledge, the products listed above do not contain any Substances of Very High Concern (SVHC) according to REACH regulations (amended on 20th June 2016) at or above 0.1%.

According to OSHA's Hazard Communication Standard (HCS 2012) 29 CFR 1910.1200 paragraph (c) this product is considered an article and defined as *"a manufactured item other than a fluid or particle: (i) which is formed to a specific shape or design during manufacture; (ii) which has end use function(s) dependent in whole or in part upon its shape or design during end use; and (iii) which under normal conditions of use does not release more than very small quantities, e.g., minute or trace amounts of a hazardous chemical (as determined under paragraph (d) of this section), and does not pose a physical hazard or health risk to employees"*.

Therefore this product is exempted from HCS 2012 obligations on classification and labelling under stated normal conditions of use (29 CFR 1910.1200 (b)(6)(v)).


See the Instructions for Use for more information

Personal Protection Wear protective gloves during use.

**Environmental
precautions**

Dispose of waste and residues as metal waste, combustible waste, or landfill in accordance with local regulations.

Approved by:

A handwritten signature in black ink that reads "Krista Huggins". The script is cursive and fluid.

Krista Huggins

Chemical Safety Specialist
Struers ApS