

Product Data Sheet



Struers ApS
Pederstrupvej 84
2750 Ballerup, Denmark

Phone +45 44 600 800
Fax +45 44 600 803
struers@struers.dk
www.struers.com

Product Name **Saw Blade SAW13**

Cat. no. 40000057 125 mm (5") dia.

Revision 02-02-2017

Description Toothed saw blade.

Material Steel.

Properties Ability to cut plastics and very soft metals.

Re-usable.

Application Cutting of plastics and very soft metals (< HV 100).
For Secotom-10/-15/-50, Accutom-10/-100 and Accutom-50/-5.

Health and Safety - Whilst the saw blade is rotating, ensure hands are kept well clear of cutting machine.

- Do not exceed a rotational speed of 1200 rpm.

To the best of Struers' knowledge, the products listed above do not contain any Substances of Very High Concern (SVHC) according to REACH regulations (amended on 20th June 2016) at or above 0.1%.

According to Article 3(3) of the REACH Regulation (EC) No 1907/2006, this product is considered an article and is defined as *"an object which during production is given a special shape, surface or design which determines its function to a greater degree than its chemical composition"*.

Therefore this product is exempted from CLP obligations and is outside the scope of the UN Globally Harmonized System of Classification and Labelling of Chemicals (GHS).

See Instructions for Use for more information.


Personal Protection

Not required.

**Environmental
precautions**

Dispose of waste and residues as metal waste, combustible waste, or landfill in accordance with local regulations.

Approved by:

A handwritten signature in purple ink that reads "Krista Huggins".

Krista Huggins

Chemical Safety Specialist
Struers ApS.