

# Product Data Sheet



**Struers ApS**  
Pederstrupvej 84  
2750 Ballerup, Denmark

Phone +45 44 600 800  
Fax +45 44 600 803  
struers@struers.dk  
www.struers.com

**Product Name**                      **Saw Blade SAW13**

**Cat. no.**                              40000057      125 mm (5") dia.

**Revision**                            02-02-2017

---

**Description**                        Toothed saw blade.

**Material**                             Steel.

**Properties**                            Ability to cut plastics and very soft metals.  
  
Re-usable.

**Application**                        Cutting of plastics and very soft metals (< HV 100).  
For Secotom-10/-15/-50, Accutom-10/-100 and Accutom-50/-5.

**Health and Safety**                - Whilst the saw blade is rotating, ensure hands are kept well clear of cutting machine.  
  
- Do not exceed a rotational speed of 1200 rpm.

To the best of Struers' knowledge, the products listed above do not contain any Substances of Very High Concern (SVHC) according to REACH regulations (amended on 20th June 2016) at or above 0.1%.

According to OSHA's Hazard Communication Standard (HCS 2012) 29 CFR 1910.1200 paragraph (c) this product is considered an article and defined as "a manufactured item other than a fluid or particle: (i) which is formed to a specific shape or design during manufacture; (ii) which has end use function(s) dependent in whole or in part upon its shape or design during end use; and (iii) which under normal conditions of use does not release more than very small quantities, e.g., minute or trace amounts of a hazardous chemical (as determined under paragraph (d) of this section), and does not pose a physical hazard or health risk to employees".

Therefore this product is exempted from HCS 2012 obligations on classification and labelling under stated normal conditions of use (29 CFR 1910.1200 (b)(6)(v)).

See Instructions for Use for more information.


**Personal Protection**

Not required.

**Environmental precautions**

Dispose of waste and residues as metal waste, combustible waste, or landfill in accordance with local regulations.

Approved by:

A handwritten signature in purple ink that reads "Krista Huggins".

Krista Huggins

Chemical Safety Specialist  
Struers ApS.