

Product Data Sheet



Struers ApS
Pederstrupvej 84
2750 Ballerup, Denmark

Phone +45 44 600 800
Fax +45 44 600 803
struers@struers.dk
www.struers.com

Product Name **Transcopy Replica Foil**

Cat. no. 40900091 Transcopy Replica Foils

Revision 21-03-2017

Description Reflective metallographic replica foils.

Material Top film of cellulose acetate metalized with aluminium.

Properties Adhesive layer glued to white plastic lining. Self-adhesive.

Application Replication of metal samples. Best used with Transcopy liquid (40900092) from Transcopy Kit.

Health and Safety To the best of Struers' knowledge, the products listed above do not contain any Substances of Very High Concern (SVHC) according to REACH regulations (amended on 20th June 2016) at or above 0.1%.

According to Article 3(3) of the REACH Regulation (EC) No 1907/2006, this product is considered an article and is defined as *"an object which during production is given a special shape, surface or design which determines its function to a greater degree than its chemical composition"*.

Therefore, this product is exempted from CLP obligations and is outside the scope of the UN Globally Harmonized System of Classification and Labelling of Chemicals (GHS).

Respiratory Protection Not required.

Environmental precautions Dispose of waste and residues as metal waste in accordance with local authority requirements.

Approved by:

A handwritten signature in purple ink that reads 'Krista Huggins'.

Krista Huggins

Chemical Safety Specialist
Struers ApS.